

# XENON<sup>2</sup> LOWER LEG LENGTH RANGES



## Xenon<sup>2</sup> FIXED FRONT & HYBRID DUAL TUBE Frame Style LOWER LEG LENGTH RANGES

FRAME HEIGHT	EXTENSION TUBE LENGTH	Platform Aluminum Angle Adjustable	Platform Composite, Carbon Fiber & Ultra Lite	Platform Composite High Mount (high mount clamps)	Auto-Folding Platform Carbon Fiber/Composite	Platform Aluminum Angle Adjustable w/ S/A Mech.
LOW	SHORT	34cm-38cm (13.3"-15")	35cm-39cm (13.8"-15.4")	N/A	33cm-40cm (13"-15.7")	N/A
	STANDARD	39cm-46cm (15.4"-18.1")	40cm-46cm (15.7"-18.1")	22cm-33cm (8.7"-13")	41cm-46cm (16.1"-18.1")	40cm-46cm (15.7"-18.1")
LOW w/ 6" Caster (Auto-Fold ONLY)	SHORT	N/A	N/A	N/A	40cm (15.7") ONLY	N/A
	STANDARD	N/A	N/A	N/A	41cm-46cm (16.1"-18.1")	N/A
HIGH	SHORT	39cm-41cm (15.4"-16.1")	37cm-43cm (14.5"-16.9")	N/A	37cm-44cm (14.6"-17.3")	N/A
	STANDARD	42cm-50cm (16.5"-19.7")	44cm-50cm (17.3"-19.7")	22cm-36cm (8.7"-14.2")	45cm-50cm (17.7"-19.7")	42cm-49cm (16.5"-19.3")

## Xenon<sup>2</sup> SWING AWAY Frame Style LOWER LEG LENGTH RANGES

EXTENSION TUBE LENGTH	Platform Aluminum Angle Adjustable	Platform Composite & Carbon Fiber	Platform Ultra Lite	Z-Finity Multi AA & Platform Alum AA High Mount (Ext. Tube is LONG)	Platform Aluminum Angle Adjustable w/ S/A Mech.	Composite & Composite Angle Adj.	Aluminum Angle Adjustable
STANDARD	32cm-40cm (12.6"-15.7")	35cm-39cm (13.8"-15.4")	30cm-39cm (11.8"-15.4")	22cm-36cm (8.7"-14.2")	35cm-40cm (13.8"-15.7")	42cm-50cm (16.5"-19.7")	39cm-50cm (15.4"-19.7")
*DEPENDANT ON Front Seat Height (FSH)	41cm-50cm (16.1"-19.7")	40cm-46cm (15.7"-18.1")	40cm-48cm (15.7"-18.9")	37cm-39cm (14.6"-15.4")	41cm-50cm (16.1"-19.7")	N/A	N/A

**\*NOTE: Dependent on required 2" (5cm) difference from the FSH**

## Xenon<sup>2</sup> SWING AWAY Frame Style LLL RANGES "CONTINUED"

EXTENSION TUBE LENGTH	ELR w/ Composite & Composite AA	ELR w/ Aluminum AA
STANDARD	45cm-50cm (17.7"-19.7")	42cm-50cm (16.5"-19.7")