



# OCTANE

## SUB4

TITANIUM  
TECHNOLOGY:  
A NEW LEVEL OF DURABILITY





\*Excluding wheels and cushion

# OCTANE

The SUB4 logo consists of the word "SUB4" in white, bold, sans-serif capital letters, set against a dark grey rectangular background. To the left of this background is a solid blue square.

THE WORLDS FIRST MADE TO MEASURE  
TITANIUM WHEELCHAIR GUARANTEED  
UNDER 4KG\*

The Octane Sub4 elevates the daily wheelchair to a whole new level of performance and engineering down to the finest detail. We took a classic design and redefined it even further by combining Sub4 technology with a hand-crafted made to measure titanium grade 9 frame. Resulting, in an ultra-lightweight and rigid wheelchair that is perfectly designed to fit you and your lifestyle.

## TITANIUM GRADE 9:

- STRENGTH TO WEIGHT RATIO
- SMOOTH RIDING
- NO CORROSION
- DURABILITY







Frame finish options: polished or brushed



Sub4 technology –  
guaranteed under 4 kg\*

## KEY FEATURES

- Made to measure
- Less than 4 kg guaranteed (excluding wheels and cushion)
- Rigid and lightweight
- Titanium grade 9 frame
- Low profile caster arm design
- Airtech Sub4 technology
- Low profile cross braces
- Optimized single tube design
- Max user weight 110 kg



Low profile caster arm design



Single post axle



# OCTANE

## TECHNICAL DATA

### OCTANE SUB4 // Titanium Wheelchair

Max. user weight	110 kg	Fully welded open frame		Max. pivot width	34.6"	Total height	35.4"	Total length	34.6"
		Custom-made		Max. safe slope	7° (for use of parking brakes)	Front wheel	3" / 4" / 5"	Rear wheel	24" / 25"
Frame angle		Custom-made							
Seat width	9.8" - 17.3"								
Seat depth	9.8" - 17.3"			Wheelchair weight	From 6.2 kg				
Seat height front	15" - 20.5"			Dismanteled weight	From 3.6 kg				
Seat height rear	14.2" - 19.7"			Mass of heaviest part	2.1 kg				
Backrest height	3.9" - 15.7"			Removable parts	Rear wheels (for transport only)				
Backrest angle	Custom-made			Colours	- Frame finish options: Brushed frame or polished frame - 5 accent colours (decal on frame), 6 backrest upholstery stitching colours, 9 Spinnery colour options				
Centre of gravity (COG)	1" - 4.9"			Intended use as seat in motor vehicle	Not released as seat in motor vehicles				
Wheel camber	1° / 2° / 3°			Transport without user	Fix the wheelchair safe in the vehicle (e.g. car or plane)				
Total width	Max. 26.8" (SW + 9.5")			Intended use and environment	For people with limited mobility (children & adults); self-propelling or being pushed by an assistant, indoor and outdoor use.				



© 2019 Sunrise Medical Canada Inc.

237 Romina Drive, Unit 3

Concord, ON

Canada L4K 4V3

Phone: 800.263.3390

Fax: 800.561.5834

E-mail: [cscanada@sunmed.com](mailto:cscanada@sunmed.com)

[www.SunriseMedical.ca](http://www.SunriseMedical.ca)



For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical Canada Inc. with any queries you may have. If you are visually impaired, this document can be viewed in PDF format at [www.SunriseMedical.ca](http://www.SunriseMedical.ca).

