

INNOVATIVE S-SERIES **SCOOTERS**


For over 30 years, Fortress has designed and built mobility scooters to the highest standards, providing the comfort and safety every customer needs.

From the shopping centre to the country, the new luxury S-Series scooters by Fortress are ready for any journey.

Packed full of exciting features with superb driving dynamics, stylish looks from every angle, there's an S-Series mobility scooter to suit everyone.

Experience luxury in motion with the Fortress S-Series scooters.









ADVANCED LED LIGHTING

Powerful all-around directional LED lighting is designed for long-term performance.

Technology with style on the Fortress S-Series.






INNOVATION IN DESIGN

With advanced LED lighting technology, the S-Series scooters are much more energy efficient than using standard bulbs, meaning more power and range on every trip.

Safety as standard, S-Series models all feature high visibility rear LED indicators, position and brake lights.







CONTROLS DESIGNED FOR YOU

Finding your perfect driving position is a breeze with the infinitely adjustable tiller. Clever design means the headlight automatically tilts as the tiller is adjusted.







COMFORT ON EVERY TRIP

The S-Series advanced seating solution takes comfort to the highest level, from the seat height, depth and recline adjustments, to the flip-up, width, angle and depth adjustable comfort armrests — every feature has been designed for comfort — for you.







→ Infinite seat depth adjustment



→ Quick backrest angle adjustment



→ Seat rotation enables comfortable transfers

COMPACT OUTDOOR SCOOTER

The luxury S425 has all of the features you want plus great power!

Taking the front/rear LED lights, comfortable all-round suspension, superior legroom, and combining it with a maximum speed of 12 km/h (7.5 mph), the S425 is suitable for indoor and outdoor use.

The increased wheel height will help to tackle those outdoor obstacles, and longer journeys can be handled with ease as the larger 55Ah batteries help to power your way.

*Baskets shown in all lifestyle pictures are smaller than what is included in the final product.









Max forward speed: 12 km/h (7.5 mph)



Colour: Choice of Urban Grey or Red



Batteries: 55Ah (Group 22NF)



Wheel size: 30 cm / 12"



Maximum user weight:

150 kg (330 lbs)

PERFORMANCE OUTDOOR SCOOTERS

The luxury S700, with superior outdoor performance is your perfect companion for the ultimate outdoor adventures, to take you where you want to go.

There's no need to worry about distance either as the large 75Ah battery size will power an extended range and a lively top speed of 15 km/h (9.3 mph).

The high performance design means gradients and higher curbs are tackled with ease every time.

The large tires, generous legroom and advanced 4-wheel suspension provides a smooth, rewarding and comfortable drive each and every time.









Max forward speed: 15 km/h (9.3 mph)



Colour: Urban Grey



Batteries: 75Ah (Group 24)



Wheel size: 33 cm / 13"



Maximum user weight: 160 kg (350 lbs)

TECHNICAL DATA 5700 *5*2125 **5-SERIES** 12 km/h (7.5 mph) 15 km/h (9.3 mph) Maximum speed 150 kg (330 lbs) Maximum userweight 160 kg (350 lbs) Battery size 55Ah x 2 75Ah x 2 8 Amp 8 Amp Charger Maximum potential range* 42 km (26 miles) 55 km (34 miles) Maximum power (peak) 2700 W 1200 W All around LED lighting Yes Yes Easy adjust tiller & wig wag control Yes Yes Suspension Yes Yes 300 mm (12") pneumatic 330 mm (13") pneumatic Tires Overall width / length 655 mm (~26") / 1345 mm (~53") 655 mm (~26") / 1500 mm (~59") Seat width x depth 490 mm (~19") x 490 mm (~19") 490 mm (~19") x 490 mm (~19") Seat height range 450 mm - 500 mm (~18" to ~20") 450 mm - 500 mm (~18" to ~20") 8 degrees at 150 kg (330 lbs) 10 degrees at 160 kg (350 lbs) Maximum Safe Slope 3030 mm (~119") 2950 mm (~116") Minimum turning radius Colour Red and Urban Grey **Urban Grey** Scooter weight inc. batteries 146 kg 117 kg

For further information on product options, please check the order forms on our website.





www.facebook.com/sunrisemedical

© 2020 Sunrise Medical Canada Inc. 237 Romina Drive, Unit #3 Concord, Ontario Canada L4K 4V3

Tel: +1.800.263.3390 Fax: +1.800.561.5834

E-mail: cscanada@sunmed.com www.SunriseMedical.ca 723487-CAN, Rev. 6

^{*} Max potential range based on industry benchmark range figures for similar products and battery sizes.