

Spex SuperShape Back Support

June 2026



MARK FOR: _____

Date: _____

Dealer Account: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ PV: _____ PC: _____

Dealer Phone: _____ Fax: _____

Confirm via: Fax E-mail

SUBMITTING FOR: Quote Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship to City: _____ PV: _____ PC: _____

Ship to Phone: _____ Fax: _____

PACKAGE INCLUDES:

Aluminium Shell
Foam System
Outer Cover
(NB hardware must be selected)



1. SELECT SIZE

✓	WIDTH	LENGTH	ITEM CODE	\$
		12"	1256-1012-OSS	\$3,507.00
	10"	14"	1256-1014-OSS	\$3,507.00
		16"	1256-1016-OSS	\$3,507.00
		12"	1256-1112-OSS	\$3,507.00
	11"	14"	1256-1114-OSS	\$3,507.00
		16"	1256-1116-OSS	\$3,507.00
		12"	1256-1212-OSS	\$3,507.00
	12"	14"	1256-1214-OSS	\$3,507.00
		16"	1256-1216-OSS	\$3,507.00
		18"	1256-1218-OSS	\$3,507.00
		12"	1256-1312-OSS	\$3,507.00
	13"	14"	1256-1314-OSS	\$3,507.00
		16"	1256-1316-OSS	\$3,507.00
		18"	1256-1318-OSS	\$3,507.00
		14"	1256-1414-OSS	\$3,737.00
		16"	1256-1416-OSS	\$3,737.00
	14"	18"	1256-1418-OSS	\$3,737.00
		20"	1256-1420-OSS	\$3,737.00
		XX	1256-14XX-OSS	\$3,737.00
		14"	1256-1514-OSS	\$3,737.00
		16"	1256-1516-OSS	\$3,737.00
	15"	18"	1256-1518-OSS	\$3,737.00
		20"	1256-1520-OSS	\$3,737.00
		XX	1256-15XX-OSS	\$3,737.00
		16"	1256-1616-OSS	\$3,737.00
		18"	1256-1618-OSS	\$3,737.00
	16"	20"	1256-1620-OSS	\$3,737.00
		23"	1256-1623-OSS	\$3,737.00
		XX	1256-16XX-OSS	\$3,737.00
		16"	1256-1716-OSS	\$3,737.00
		18"	1256-1718-OSS	\$3,737.00
	17"	20"	1256-1720-OSS	\$3,737.00
		23"	1256-1723-OSS	\$3,737.00
		XX	1256-17XX-OSS	\$3,737.00
		16"	1256-1816-OSS	\$3,737.00
		18"	1256-1818-OSS	\$3,737.00
	18"	20"	1256-1820-OSS	\$3,737.00
		23"	1256-1823-OSS	\$3,737.00
		XX	1256-18XX-OSS	\$3,737.00
		16"	1256-1916-OSS	\$4,325.00
		18"	1256-1918-OSS	\$4,325.00
	19"	20"	1256-1920-OSS	\$4,325.00
		23"	1256-1923-OSS	\$4,325.00
		XX	1256-19XX-OSS	\$4,325.00
		16"	1256-2016-OSS	\$4,325.00
		18"	1256-2018-OSS	\$4,325.00
	20"	20"	1256-2020-OSS	\$4,325.00
		23"	1256-2023-OSS	\$4,325.00
		XX	1256-20XX-OSS	\$4,325.00
	Up to 20"	XX	1256-2929-OSS	By Quote Only
	Up to 25"	XX	1256-2939-OSS	By Quote Only

2. SELECT COLOUR NC

Your colour selection applies to any extra covers ordered

- | | | |
|--|---|---|
| <input type="radio"/> Black
(default) | <input type="radio"/> Granite
1209-2491-090 | <input type="radio"/> Calypso
1209-2491-028 |
| <input type="radio"/> Chilli
1209-2491-022 | <input type="radio"/> Mango
1209-2491-030 | <input type="radio"/> Ocean
1209-2491-060 |
| <input type="radio"/> Leaf
1209-2491-054 | <input type="radio"/> Bumblebee
1209-2491-044 | <input type="radio"/> Royal
1209-2491-077 |

3. SELECT HARDWARE

✓ ITEM	IMAGE	\$
<input type="radio"/> Quad Mount (10° angle adjustable) 1259-2400-100		NC
<input type="radio"/> Dual Mount (30° angle adjustable, for 14" backs up) 1259-2200-100		\$526.00
<input type="radio"/> Quad Mount - Heavy Duty (10° angle adjustable) 1259-2405-100		\$915.00
<input type="radio"/> Width Adjustment Bracket Set (for Quad & Single mount only) Allows a wider back installed to a narrower wheelchair 1259-2490-100		\$357.00
<input type="radio"/> Quad Extension Brackets To add an extra 10° of angle (total 20°) 1259-2450-100		\$166.00

4. ADD POSITIONING WEDGES

✓ ACCESSORY	CODE	IMAGE	\$
<input type="radio"/> Additional backrest lumbar pad and lateral stability wedges	1209-2713-300		\$134.00
<input type="radio"/> Asymmetric narrow contour wedge 200mm high x 100mm wide x 40mm deep	1209-2722-300		\$134.00
<input type="radio"/> Asymmetric bi-angular contour wedge 200mm high x 100mm wide x 40mm deep	1209-2723-300		\$134.00
<input type="radio"/> Asymmetric contour kit	1209-2720-300		\$196.00

5. SELECT ACCESSORIES Input your back dimensions

✓ ACCESSORY	IMAGE	\$
<input type="radio"/> Extra Outer Cover 1256-____-0SS-C		\$361.00
<input checked="" type="radio"/> Positioning Kit (20 pink cubes + 15 grey cubes) 1209-2715-300-BK SEBND0015		\$150.00
<input type="radio"/> Extra Kit Positioning Cubes (35 Cubes) 1256-____-300 (State width and height of the backrest)		\$150.00



Lateral Trunk Support

1. SELECT HARDWARE (All hardware scripted as each, not in pairs.)

Standard Swing Away

L	R	STYLE	IMAGE	CODE	\$
		Swing-Away, straight		1266-2201-100	\$373.00
		Swing-Away, straight extend		1266-2202-100	\$373.00
		Swing-Away, offset 20mm		1266-2203-100	\$373.00
		Swing-Away, offset 40mm		1266-2204-100	\$373.00

Swing Away Axial

L	R	STYLE	IMAGE	CODE	\$
		Kids Swing Away Hardware*		1266-2623-400	\$662.00
		Axial Rotation Hardware		1266-2301-100	\$471.00

*Kids Swing Away Hardware is not Left or Right specific

Fixed

L	R	STYLE	CODE	\$
		Lateral Bracket	1266-2101-100	\$202.00
		Kidz Lateral Bracket	1266-2111-100	\$195.00

Kidz and standard hardware options are compatible with all lateral pad sizes

2. SELECT PAD WITH COVER (Lateral pads are no charge when purchased with a back.)

	STYLE	HEIGHT (mm)	DEPTH (mm)	STANDARD COVER			COMFI COVER				
				L	R	CODE	\$	L	R	CODE	\$
Narrow pad shape 	Narrow Straight 0	85	115			1266-3010-400	\$245.00			1266-3010-800	\$245.00
	Narrow Straight 1	95	145			1266-3011-400	\$245.00			1266-3011-800	\$245.00
	Narrow Straight 2	95	180			1266-3012-400	\$245.00			1266-3012-800	\$245.00
	Narrow Straight 3	115	210			1266-3013-400	\$245.00			1266-3013-800	\$245.00
	Narrow Straight 4	115	240			1266-3014-400	\$245.00			1266-3014-800	\$245.00
	Narrow Curved 1	95	145			1266-3021-400	\$245.00			1266-3021-800	\$245.00
	Narrow Curved 2	95	180			1266-3022-400	\$245.00			1266-3022-800	\$245.00
	Narrow Curved 3	115	210			1266-3023-400	\$245.00			1266-3023-800	\$245.00
	Narrow Curved 4	115	240			1266-3024-400	\$245.00			1266-3024-800	\$245.00
Broad pad shape 	Broad 1	105	130			1266-3031-400	\$245.00			1266-3031-800	\$245.00
	Broad 2	125	140			1266-3032-400	\$245.00			1266-3032-800	\$245.00
	Broad 3	140	150			1266-3033-400	\$245.00			1266-3033-800	\$245.00
	Broad 4	160	155			1266-3034-400	\$245.00			1266-3034-800	\$245.00
	Broad 5	195	165			1266-3035-400	\$245.00			1266-3035-800	\$245.00
	Broad 6	230	175			1266-3036-400	\$245.00			1266-3036-800	\$245.00


EXTRA BRACKET OPTIONS

L	R	STYLE	IMAGE	CODE	\$	L	R	STYLE	IMAGE	CODE	\$
		Left Cane Mount (attachment directly to wheelchair canes)		1266-2411-100	\$245.00			Deep Contour Extended Bracket Kit: extends depth up to 115mm Recommended for SuperShape when deep contour is created for user		1266-2430-100	\$245.00
		Right Cane Mount (attachment directly to wheelchair canes)		1266-2412-100	\$245.00			Kidz Lateral Mounting Bracket: 20mm shorter along shell for closer positioning of laterals		1266-2442-100	\$106.00
		Reinforcement plate (forceful side leaning or clients over 100kg)		1266-2450-100	\$148.00			Kidz Lateral Mounting Bracket: 20mm longer along shell for increased width growth		1266-2444-100	\$245.00
		Deep Contour Bracket: extends depth by 30mm		1266-2420-100	\$924.00						

Axial Biangular Lateral Trunk Support Order Form

1. SELECT HARDWARE

(Laterals are scripted as each, not in pairs)



MSRP \$556.00 EA

L	R
<input type="checkbox"/>	<input type="checkbox"/>
1266-2501-100	

2. SELECT BIANGULAR PAD SIZE

AXIAL PAD



SIZES AVAILABLE

4 (165mm x 175mm)

-XX44-

6 (215mm x 195mm)

-XX46-

8 (265mm x 195mm)

-XX48-

L R

<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>

5. CUSTOMISE




Left Cane Mount
(attachment directly to wheelchair canes)
1266-2411-100

\$245.00



Right Cane Mount
(attachment directly to wheelchair canes)
1266-2412-100

\$245.00



Reinforcement plate
(forceful side leaning or clients over 100kg)
1266-2450-100

\$148.00



Deep Contour Bracket
(extends depth by 30mm)
1266-2420-100

\$126.00



Deep Contour Extended Bracket Kit
(extends depth up to 115mm)
1266-2430-100

\$245.00

3. SELECT PADDING



12mm Pressure-Reduction
Overall Pad Thickness 30mm
-X00



24mm Pressure-Reduction
Overall Pad Thickness 45mm
-X10

6. BUILD CODE FOR LATERAL PAD

REFER TO THE LATERAL TRUNK SUPPORT CODES ABOVE

1	2	6	6	-	3	0	<input type="checkbox"/>	-	<input type="checkbox"/>	<input type="checkbox"/>		
							ADD PAD CODES		ADD COVER CODES		ADD PADDING CODES	

4. SELECT COVER




Standard (-4XX)



Comfi (-8XX)

Contact Customer Service for a quote
1-(800) 263-3390 cscanada@sunmed.com

LIFE-SIZE LATERAL TEMPLATE



Spex® provides the tools to make decisions for faster results, including this lateral pad size guide to scale. When deciding which lateral support pad style and size will suit your clients, place this page against their torso to make an informed decision on the spot!