



Q300M with SEDEO® PRO™ and JAY Syncra Seating

AADL - W577A

EIPW40

Account:	
Date:	
PO#:	
Buyer Name:	
Marked For:	


Ship To:	
Name:	
Address:	
City, Prov., PC:	
Phone:	

Sunrise Medical Inc.
1000 Creditstone Road, Unit 2
Concord, ON L4K 4P8
1-800-263-3390
www.sunrisemedical.ca



Client Measurements

a)	Top of head to seat	Inch
b)	Top of shoulder to seat	Inch
c)	Underarm to seat	Inch
d)	Elbow to seat	Inch
e)	Top of leg to seat	Inch
f)	Shoulder width	Inch
g)	Chest width	Inch
i)	Hip width	Inch
j)	Actual sitting depth	Inch
m)	Foot length	Inch
(n)	Knee to foot	Inch



User Weight

Notes:

Chair Specifications

		Seat Type	SEDEO PRO / JAY Syncra
Weights starting at: (with Batteries)	NPC: 287 lbs SPC: 312 lbs MPC: 322 lbs.	Seat Size Range	Width 16" to 22" Depth 14" to 22"
Power Seating	50° of CG Tilt and 150° Power Recline with 4" mechanical sheer reduction / 12" seat elevator	STF Height Range	16.5" to 19"

Mid-Wheel Drive Base Specifications

Drive Wheel Position	Mid-wheel Drive	Drive Wheels	12" drive wheels w/flat free inserts (black) 14" drive wheels w/flat free inserts (black)
User Weight Capacity	300 lbs.	Tires/Casters	7" Front / 7" Rear (Black)
Min. Turning Radius	17" (excludes footplates)	Charger	8A Dual Mode (Compact)
Base Size	Width 20.5", Length 33"	Battery Size	50Ah Group 22
Drive Control	R-net 90A / VR-2 90A / LiNX 120A	Color Accents	9 Shroud Colors and 7 Color Accents
Driving Range	Up to 17 miles	Ground Clearance	3"
Drive Motors	4-pole 6mph		

BASE				
OPTION		Description		Retail
Group 3 Standard Duty Mid Wheel Drive Base - End Item = EIPW40				
☐	400PC291		Quickie Q300 M NPC No Power Option for SEDEO Pro Seat Included in Base Pricing: Q300 M Base, Non-Expandable Electronics, No Power Seat Functions, Six-Wheel Suspension, 12" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5". (AADL Contract Price)	\$8,072.00
☐	400PC292		Quickie Q300 M SPC Single Power Option for SEDEO Pro Seat Included in Base Pricing: Q300 M Base, Non-Expandable Electronics, No Power Seat Functions, Six-Wheel Suspension, 12" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5". (AADL Contract Price)	\$8,072.00
☐	400PC293		Quickie Q300 M MPC Tilt/Recline For SEDEO Pro Seat Included in Base Pricing: Q300 M Base, Non-Expandable Electronics, No Power Seat Functions, Six-Wheel Suspension, 12" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5". (AADL Contract Price)	\$8,072.00
User Weight				
☐	400WG23		Up to 300 lbs.	NC
Motor Package				
☐	400PF99	-	4-pole 10kph / 6mph	NC
Tracking				
☐	400MI01	R-net	SureTrac Module & Motor Encoders <i>(Improves Drive Tracking)</i> for R-net	\$1,696.00
☐	400MI04	LiNX	Traction Control System for LiNX (GYRO Technology)	\$925.00
Drive Wheels				
☐	400DW13		Drive Wheel, 12", Solid, Knobby, Black	NC
☐	400DW58		Drive Wheel, 14", Solid, Knobby, Black	NC
Caster Wheels				
☐	400CW59		Caster Wheels 7", Solid, Black	NC
Base Shroud Color NEW Base Shroud and Back Shroud are painted to match				
☐	400C164		Black	NC
☐	400C165		Candy Apple Red	NC
☐	400C166		Candy Blue	NC
☐	400C167		Candy Purple	NC
☐	400C176		Green Apple	NC
☐	400C178		Pearl Pink	NC
☐	400C182		Silver	NC
☐	400C125		White	NC

BASE			
OPTION		Description	Retail
Color Accents <i>Note: Accent colors and shroud colors may not match</i>			
☐	400EZ20	Apple Green	NC
☐	400EZ22	Black	NC
☐	400EZ21	Blue	NC
☐	400EZ23	Orange	NC
☐	400EZ24	Red	NC
☐	400EZ19	White	NC
☐	400EZ25	Silver	NC
Service Decals			
☐	400D10	Army	NC
☐	400D11	Air Force	NC
☐	400D13	Marines	NC
☐	400D14	Navy	NC
☐	400D15	Purple Heart	NC
ISO Transit Tie Downs			
☐	400TO50	ISO Transit Package	NC
LED Lights (Client Cost)			
☐	400LI08	LED Fender Lights	\$646.00
☐	400LI09	LED Fender Lights <i>(Available w/ Multi-Power Only)</i>	NC
☐	400LI02	Omit LED Lights	NC
Seat Heights			
<i>Note: Available Seat-to-Floor heights are dependant on Powered Seating. Use the matrix below to determine available ranges for your desired options.</i>			
<i>Measured from front of seat pan to floor.</i>			
☐	400HT103	16.5"	NC
☐	400HT55	17"	NC
☐	400HT104	17.5"	NC
☐	400HT56	18"	NC
☐	400HT105	18.5"	NC
☐	400HT57	19"	NC
Batteries & Chargers (AADL Contract Pricing)			
Group 22 Battery			
☐	400AG35	GP22/50Ah Flat Top Pair	\$678.00
Omit Batteries.			
☐	400AG24	Battery Harnesses Only <i>(Omit Batteries)</i>	NC
Chargers			
☐	400BJ30	10A Dual Mode	\$162.00
☐	400BJ19	8A Dual Mode	NC

POWERED SEATING			
OPTION		Description	Retail
Power Seating Options			
Seat Elevator (Client Cost)			
☐	400PI07	Seat Elevator, Powered, 300 mm / 12"	\$3,477.00
☐	400PI08	C-Me Seat Elevator, Powered, 300 mm / 12" (C-Me speed up to 3.0 mph)	\$3,477.00
Seat Tilt (AADL W992 For Power Tilt)			
☐	400PA48	Posterior 50° Tilt	\$6,496.00
Seat Recline (AADL W982TR for Powr Tilt and Recline)			
☐	400PA37	150° Powered Recline w/ 4" Shear Reduction	\$8,001.00
Seat Tilt/Recline			
☐	400PA41	Posterior 50° Tilt and 150° Recline Powered, w/ 4" Shear Reduction	\$13,275.00
Power Recline Backrest Preset Home Angle			
☐	400QK1	90°	NC
☐	400QK2	95° (Default Position)	NC
Pre Set Seat Frame Angle (No Power Seating Options, Increases Seat to Floor Height)			
☐	400FA10	0° (Adds 0" to STFH)	NC
☐	400FA3	3° (Adds 1" to STFH)	NC
Backrest Frame Type			
☐	400BF16	Tool Angle Adjustable Back 90°-120°	NC
☐	400BF14	User Adjustable Manual Semi-Recline Back, 90°-130°, Right (Recline lever mounted on opposite side of joystick)	\$633.00
☐	400BF15	User Adjustable Manual Semi-Recline Back, 90°-130°, Left (Recline lever mounted on opposite side of joystick)	\$633.00
Backrest Frame Option			
☐	400BF13	Flip Down, Quick Release	\$179.00
Back Cushion Type			
☐	400BT134	SEDEO Back Cushion (Choose Seating in Section A)	NC
☐	400BT138	JAY Syncra (Choose Seating in section B)	NC
☐	400BT135	JAY J3 (Choose Seating in section C)	NC
☐	400BT137	SEDEO PRO Aftermarket Back Adapter (Height Adjustable Round Tubes for Traditional Mounting)	\$633.00
☐	400BT133	OMIT (For Custom Mounting. No Headrest Mount)	NC

POWERED SEATING

OPTION	Description	Retail
--------	-------------	--------

A - SEDEO Backrest Dimensions (Seating Benefits)

Back Width

☐	400BN28	400mm / 16"	NC
☐	400BN29	460mm / 18"	NC
☐	400BN30	500mm / 20"	NC
☐	400BN31	560mm / 22"	NC

Back Shell Height

			Tool Angle Adjust.	150° Recline	Installed Back Height	
☐	400B255	460mm / 18"	22"-23"	23"-25"		NC
☐	400B256	560mm / 22"	24"-27"	26"-29"		NC
☐	400B257	660mm / 26"	28"-30"	30"		NC

Installed Back Height - (From Pan to Top of Back)

☐	400JR32	22"	NC
☐	400JR33	23"	NC
☐	400JR34	24"	NC
☐	400JR35	25"	NC
☐	400JR36	26"	NC
☐	400JR37	27"	NC
☐	400JR38	28"	NC
☐	400JR39	29"	NC
☐	400JR40	30"	NC

POWERED SEATING

OPTION	Description	Retail
--------	-------------	--------

B - JAY Syncra Seating

B - JAY Syncra Back (Seating Benefits)

☐	400JP14	JAY Syncra Back 16-20" Wide with 3" Contour Depth	\$1,339.00
☐	400JP15	JAY Syncra Back 22" Wide with 3" Contour Depth	\$1,339.00

B - JAY Syncra Back Width

Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.

☐	400JQ07	16" Wide (Adjustment range 16-17" Wide)*	NC
☐	400JQ08	18" Wide (Adjustment range 18-19" Wide)*	NC
☐	400JQ09	20" Wide (Adjustment range 20-21" Wide)*	NC
☐	400JQ10	22" Wide (Adjustment range 22-23" Wide)	NC

*Additional +1" of shell width adjustment is achievable, requires purchase of 1" larger foam and upholstery separately

B - JAY Syncra Back Shell Options and Height Ranges

JAY Syncra Back shells are mounted with a minimim 5.5" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field with-in the ranges shown based on the back frame type.

JAY Syncra Shell Height

☐	400JR42	16" Tall Shell (22-24" Back Height)	NC
☐	400JR43	18" Tall Shell (24-26" Back Height)	NC
☐	400JR44	20" Tall Shell (26-28" Back Height)	NC
☐	400JR45	22" Tall Shell (28-30" Back Height)	NC

JAY Syncra Back Installed Back Height (From Pan to Top of Back)

☐	400JR50	22" (16" shell height)	NC
☐	400JR51	23" (16" shell height)	NC
☐	400JR52	24" (16" shell height)	NC
☐	400JR53	24" (18" shell height)	NC
☐	400JR54	25" (18" shell height)	NC
☐	400JR55	26" (18" shell height)	NC
☐	400JR56	26" (20" shell height)	NC
☐	400JR57	27" (20" shell height)	NC
☐	400JR58	28" (20" shell height)	NC
☐	400JR67	28" (22" shell height)	NC
☐	400JR68	29" (22" shell height)	NC
☐	400JR69	30" (22" shell height)	NC

B - JAY Syncra Back Upholstery

☐	400BU96	Syncra Premium Leatherette	NC
☐	400BU97	Syncra Stretch	NC

POWERED SEATING

OPTION	Description	Retail
--------	-------------	--------

C - JAY J3 Back

C - JAY J3 Back - Contour (Seating Benefits)

☐	400JP01	J3 PA Back for 14"-20" Wide	\$912.00
☐	400JP02	J3 PA Back for 22" Wide	\$1,941.00
☐	400JP11	J3 PL Back for 14"-20" Wide	\$986.00
☐	400JP13	J3 PD Back for 14"-20" Wide (18" and 21" shell height only)	\$1,068.00

C - JAY J3 Back Width

Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.

☐	400JQ02	16"	NC
☐	400JQ03	18"	NC
☐	400JQ04	20"	NC
☐	400JQ05	22"	NC

C - JAY J3 Back Shell Options and Height Ranges.

JAY J3 shells are mounted with a minimim 4" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field within the ranges shown.

JAY J3 Shell Height

			Tool Angle Adjust.	150° Recline		
☐	400JR41	15" tall shell (Upper Thoracic Med) for 22" Wide	22"-23"	22"-23"	Installed Back Height	NC
☐	400JR14	16.5" tall shell (Upper Thoracic Short)	21"-23"	22"-24"		NC
☐	400JR15	18" tall shell (Upper Thoracic Med)	22"-24"	23"-25"		NC
☐	400JR16	19.5" tall shell (Upper Thoracic Tall)	24"-26"	24"-26"		NC
☐	400JR18	18" tall shell (Shoulder Height) for 22" wide	24"	24"		NC
☐	400JR17	21" tall shell (Shoulder Height)	24"-26"	25"-28"		NC
☐	400JR27	22.5" tall shell (Shoulder Height)	26"-28"	27"-30"		NC
☐	400JR28	24" tall shell (Shoulder Height)	27"-29"	28"-30"		NC

JAY J3 Installed Back Height (From Pan to Top of Back)

☐	400JR29	21"	NC
☐	400JR19	22"	NC
☐	400JR20	23"	NC
☐	400JR21	24"	NC
☐	400JR22	25"	NC
☐	400JR23	26"	NC
☐	400JR24	27"	NC
☐	400JR25	28"	NC
☐	400JR30	29"	NC
☐	400JR31	30"	NC

POWERED SEATING

OPTION	Description	Retail
--------	-------------	--------

C - JAY Anterior Thoracic Support (Seating Benefits)

Classic-Structure Pull Down Harness (4 point strap, Only available with JAY backrest)

☐	400K10	Medium (ATS8302)	\$325.00
☐	400K11	Large (ATS8303)	\$325.00

C - JAY Lateral Thoracic Swing-Away Supports (Requires Pad, Upholstery and Mount)

JAY Lateral Thoracic Support Pads

☐	400JW35	Lateral Thoracic (LxH) 100 mm x 127 mm / 4"x 5" Pad - Right	\$73.00
☐	400JW38	Lateral Thoracic (LxH) 100 mm x 127 mm / 4"x 5" Pad - Left	\$73.00
☐	400JW36	Lateral Thoracic (LxH) 100 mm x 152 mm / 4" x 6" Pad - Right	\$73.00
☐	400JW39	Lateral Thoracic (LxH) 100 mm x 152 mm / 4" x 6" Pad - Left	\$73.00
☐	400JW37	Lateral Thoracic (LxH) 152 mm x 100 mm / 6" x 4" Pad - Right	\$73.00
☐	400JW40	Lateral Thoracic (LxH) 152 mm x 100 mm / 6" x 4" Pad - Left	\$73.00

JAY Lateral Thoracic Support Pad Upholstery

☐	400NU01	Syncra Premium Leatherette	NC
☐	400NU02	Syncra Stretch	NC

JAY Lateral Thoracic Swing-Away Mount

☐	400JW41	Lateral Thoracic Swing Away Bracket - Right	\$325.00
☐	400JW42	Lateral Thoracic Swing Away Bracket - Left	\$325.00

C - JAY Syncra Fixed Lateral Support Accessory

JAY Fixed Lateral Support Accessory Hardware

☐	400JW33	Fixed Lateral Support Accessory (LxH) 90 mm x 140 mm / 3.5" x 8" - Right	\$190.00
☐	400JW34	Fixed Lateral Support Accessory (LxH) 90 mm x 140 mm / 3.5" x 8" - Left	\$190.00

*JAY Fixed Lateral Support is adjustable

POWERED SEATING

OPTION	Description	Retail
--------	-------------	--------

Headrests (Requires Pad, Upholstery and Mount) (Seating Benefits)

Headrest Pad

☐	400HE36	SEDEO Pad Small With Cover (L x H) 240 x 150 mm / 10" x 6" (Only Available with SEDEO Multi-Position Headrest Mount and SEDEO Upholstery)	\$227.00
☐	400HE38	Whitmyer Plush Headrest Pad 203 mm / 8"	\$377.00
☐	400HE39	Whitmyer Plush Headrest Pad 254 mm / 10"	\$377.00
☐	400HE44	Whitmyer Plush Headrest Pad 356 mm / 14"	\$377.00
☐	400HE46	Whitmyer Adjust-A-Plush Headrest Pad Medium (Adjustable from 7" to 16" Wide x 4" Tall)	\$445.00
☐	400HE48	Omit Headrest Pad	NC

Headrest Pad Upholstery

☐	400NU10	SEDEO Upholstery	NC
☐	400NU11	Whitmyer Lycra / Stretch	NC

Headrest Mount

☐	400MH32	SEDEO Multi-Position Headrest Mount with Quick Release (Only available with SEDEO Pad)	\$557.00
☐	400MH33	Whitmyer Linx Headrest Mount with Quick Release	\$247.00
☐	400MH34	Whitmyer Linx-2 Headrest Mount with Quick Release	\$298.00
☐	400MH35	Whitmyer AXYS Headrest Mount with Quick Release	\$247.00
☐	400MH17	Omit Headrest Mount	NC

Push Handles

☐	400PH27	Height Adjustable Push Handles	\$61.00
--------------------------	---------	--------------------------------	---------

Seat Frame Rails

☐	400SN14	Seat Rail Short, For Effective Seat Depth 14" - 18"	NC
☐	400SN15	Seat Rail Long, For Effective Seat Depth 17" - 22"	NC

Seat Cushion Type

☐	400ST81	JAY Syncra Cushion (See Section E)	NC
☐	400ST80	JAY (See Section F)	NC
☐	400ST75	Omit	NC

POWERED SEATING

OPTION	Description	Retail
--------	-------------	--------

Seat Size

Seat Pan Width

☐	400R53	16"	NC
☐	400R54	17"	NC
☐	400R55	18"	NC
☐	400R56	19"	NC
☐	400R57	20"	NC
☐	400R58	21" (N/A for JAY Syncra Cushions)	NC
☐	400R59	22"	NC

Seat Pan Depth (Front of seat pan to front of back cushion - effective depth)

☐	400RF52	14"	NC
☐	400RF53	15"	NC
☐	400RF51	16"	NC
☐	400RF45	17"	NC
☐	400RF46	18"	NC
☐	400RF47	19"	NC
☐	400RF48	20"	NC
☐	400RF49	21" (N/A for JAY Syncra Cushions)	NC
☐	400RF50	22"	NC

E - JAY Syncra Cushion (Jay Your Way Modifications Available Through Customer Service (800) 263-3390) (Seating Benefits)

☐	400CU89	JAY Syncra Cushion featuring JAY Fluid with Cryo Technology (16"-20" wide)	\$1,017.00
☐	400CU90	JAY Syncra Cushion featuring JAY Fluid with Cryo Technology (22" wide)	\$1,125.00

JAY Syncra Cushion Cover Upholstery

☐	400YR73	Syncra Premium Leatherette	NC
☐	400YR74	Syncra Stretch	NC

F - JAY Seat Cushion (Seating Benefits)

☐	400CU59	JAY Ion (Sizes 16" - 21" Wide)	\$505.00
☐	400CU60	JAY Ion (Size 22" Wide)	\$698.00
☐	400CU61	JAY Union (Sizes 16" - 21" Wide)	\$764.00
☐	400CU62	JAY Union (Size 22" Wide)	\$806.00
☐	400CU64	JAY Fusion (Sizes 16" - 21" Wide)	\$1,045.00
☐	400CU65	JAY Fusion (Size 22" Wide)	\$1,172.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	HPCP	Description	Retail
--------	------	-------------	--------

Knee and Thigh Swing Away Support Pads (Seating Benefits)

Lateral Knee Support Pads and Hardware (Must pick pad and mount)

☐	400XY05	Lateral Knee Support Pad, Right (4" round)	\$247.00
☐	400XY06	Lateral Knee Support Pad, Left (4" round)	\$247.00
☐	400XX08	Lateral Swing Away Bracket Knee - Right	\$325.00
☐	400XX09	Lateral Swing Away Bracket Knee - Left	\$325.00

Lateral Thigh Support Pads and Hardware (Must pick pad and mount)

☐	400KF22	Lateral Thigh Support Pad Short, Right (4.8" H x 4.8" L / 120mm x 120mm)	\$146.00
☐	400KF23	Lateral Thigh Support Pad Short, Left (4.8" H x 4.8" L / 120mm x 120mm)	\$146.00
☐	400KF24	Lateral Thigh Support Pad Long, Right (4.8" H x 7" L / 120 mm x 180mm)	\$146.00
☐	400KF25	Lateral Thigh Support Pad Long, Left (4.8" H x 7" L / 120mm x 180mm)	\$146.00
☐	400KF26	Lateral Thigh Removable Mount - Right	\$325.00
☐	400KF27	Lateral Thigh Removable Mount - Left	\$325.00

JAY Lateral Thigh Supports (Must Pick Pad, Upholstery and Mount) (Seating Benefits)

JAY Lateral Thigh Support Pads

☐	400KF16	Lateral Thigh Support Pad Standard, Right (8.5" x 3.25")	\$146.00
☐	400KF17	Lateral Thigh Support Pad Standard, Left (8.5" x 3.25")	\$146.00
☐	400KF18	Lateral Thigh Support Pad Large, Right (12" x 3.25")	\$146.00
☐	400KF19	Lateral Thigh Support Pad Large, Left (12" x 3.25")	\$146.00

JAY Lateral Thigh Support Pad Upholstery

☐	400NU07	Syncra Premium Leatherette	NC
☐	400NU08	Syncra Stretch	NC

JAY Lateral Thigh Removable Mount

☐	400KF20	Lateral Thigh Removable Mount - Right	\$325.00
☐	400KF21	Lateral Thigh Removable Mount - Left	\$325.00

Armrest Structure

☐	400EG35	Reclining, Flip Back, Long (Height Adjustable 9" to 15")	NC
☐	400EG41	Reclining, Flip Back, Short (Best used with shorter seat depths) (Height Adjustable 9" to 15")	NC
☐	400EG36	Reclining, Flip Back with 2-point post locking arm (Height Adjustable 11" to 15" x .625")	\$351.00
☐	400EG51	Reclining, Flip Back, Dual Link Armrest (No Charge for Multi-Power) (Height Adjustable 8" to 15" including length and lateral adjustment)	NC
☐	400EG52	Reclining, Flip Back, Dual Link Armrest (Height Adjustable 8" to 15" including length and lateral adjustment)	\$350.00
☐	400EG32	Single Post, Flip Back (N/A w/ thigh support)(Height Adjustable 10" to 15")	NC
☐	400EG50	Dual Post, Flip Back, Removeable (Height Adjustable)	NC

POSITIONING SUPPORTS & ACCESSORIES

OPTION	Description	Retail
--------	-------------	--------

Arm Pads

SEDEO Armrest Pad

☐	400AB43	300 mm / 12", Integral Foam, Right	NC
☐	400AE43	300 mm / 12", Integral Foam, Left	NC
☐	400AB44	400 mm / 16", Integral Foam, Right	NC
☐	400AE44	400 mm / 16", Integral Foam, Left	NC

Dual Post Armrest Pad

☐	400AB74	250 mm / 10", Classic Desk Foam, Right	NC
☐	400AE74	250 mm / 10", Classic Desk Foam, Left	NC
☐	400AB76	250 mm / 10", Waterfall Desk, Right	NC
☐	400AE76	250 mm / 10", Waterfall Desk, Left	NC
☐	400AB75	250 mm / 14", Classic Full Foam, Right	NC
☐	400AE75	250 mm / 14", Classic Full Foam, Left	NC
☐	400AB77	250 mm / 10", Waterfall Full, Right	NC
☐	400AE77	250 mm / 10", Waterfall Full, Left	NC

Comfort Arm Pad (Foam)

☐	400AB46	2" Wide Right	\$113.00
☐	400AE46	2" Wide Left	\$113.00

Comfort Arm Pad (Gel)

☐	400AB48	2" Wide Right	\$213.00
☐	400AE48	2" Wide Left	\$213.00

Ottobock Arm Troughs

☐	400AE06	Ottobock Arm Trough, Med (12"), Left	\$139.00
☐	400AB06	Ottobock Arm Trough, Med (12"), Right	\$139.00
☐	400AE07	Ottobock Arm Trough, Large (14"), Left	\$139.00
☐	400AB07	Ottobock Arm Trough, Large (14"), Right	\$139.00
☐	400OL01	Ottobock Hand Pad, Med (7"), Left	\$70.00
☐	400OR01	Ottobock Hand Pad, Med (7"), Right	\$70.00
☐	400OL02	Ottobock Hand Pad, Large (9"), Left	\$70.00
☐	400OR02	Ottobock Hand Pad, Large (9"), Right	\$70.00
☐	400AE78	Ottobock 1 Piece Arm Trough, (18"), Left	\$198.00
☐	400AB78	Ottobock 1 Piece Arm Trough, (18"), Right	\$198.00

Side Guards (Client Cost)

☐	400SG44	Side Guards Integral Foam (11"L x 5.5"H) (N/A with Reclining Armrests)	\$356.00
--------------------------	---------	--	----------

Centermount Legrest

☐	400H184	Fixed Centermount Legrest, Short (Includes Calf Pads)(11.5" to 14.5" Leg Length)	NC
☐	400H185	Fixed Centermount Legrest, Long (Includes Calf Pads)(13.5" to 18.5" Leg Length)	NC
☐	400H136	Angle Adjustable Centermount (Includes Calf Pads and Single Footplate)(11.5" to 18.5" Leg Length)	\$633.00

Power Centermount Legrest (Client Cost)

☐	400H109	SEDEO Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension (Powered Articulation) (3" to 17.5" Leg Length)	\$3,435.00
--------------------------	---------	---	------------

POSITIONING SUPPORTS & ACCESSORIES

OPTION	Description	Retail
--------	-------------	--------

SEDEO Centermount Legrest Accessories (Client Cost)

SEDEO Power Centermount Footplates

☐	400F127	Dual, Height, Width and Angle Adjustable, Standard (14"W x 10"D)	\$284.00
☐	400F128	Dual, Height, Width and Angle Adjustable, Large (16"W x 12"D)	\$284.00

SEDEO Power Footplate Accessories

☐	400DI68	Extra Large Centermount Footplate	\$315.00
☐	400DI69	Footplate Connector Link	\$96.00
☐	400DI57	Gel Footpads	\$401.00

SEDEO Power Centermount Calf Pads

☐	400DI34	Calf Pads, Pair, Standard Integral Foam	NC
--------------------------	---------	---	----

SEDEO PRO Swing-Away Legrest (Pair) (400H179 is Client Cost)

☐	400H177	Legrests, Swing Away, 70°, SEDEO PRO (14" to 21" Leg Length)	NC
☐	400H178	Legrests, Manual Elevating, Swing Away, SEDEO PRO (Includes Padded Calf Pads) (14".5 to 21.5" Leg Length)	\$380.00
☐	400H179	Legrests, Power Elevating, Swing Away, SEDEO PRO (Includes Padded Calf Pads) (15.5" to 22.5" Leg Length)	\$3,435.00
☐	400H176	Omit Legrests	NC

Swing Away Legrest Footplates

☐	400F369	Angle Adjustable Flip Up Footplates	\$284.00
--------------------------	---------	-------------------------------------	----------

Swing Away Legrest Calf Supports

☐	400DI42	Calf Support, Depth Adjustable, Padded (6.5" to 8" Wide, Depending On Chair Width, x 9.5" Tall)	\$512.00
--------------------------	---------	---	----------

Swing Away Legrest Options

☐	400DI44	Heel Loops	\$0.00
☐	400KG4	Knee Pads, Integral Foam (2.5" Wide x 4.5" Tall)	\$282.00

SEDEO Lower Extremity Swing Away Support

☐	400VD01	Lower Extremity Support, Padded, Swing Away, SEDEO, Right	\$512.00
☐	400VD02	Lower Extremity Support, Padded, Swing Away, SEDEO, Left	\$512.00

Pelvic Positioning Belts

☐	400A837	Positioning Belt, Carseat Buckle 50mm x 1600mm / 2"	NC
☐	400A839	Positioning Belt, Aircraft Buckle 50mm x 1300mm / 2"	NC
☐	400A840	Positioning Belt, Aircraft Buckle, Padded 50mm x 1300mm / 2"	\$103.00
☐	400A845	Bodypoint - Evoflex Pelavic Stabilizer Belt	\$216.00
☐	400A846	Bodypoint - Evoflex Pelavic Stabilizer Belt with Extension	\$271.00

POSITIONING SUPPORTS & ACCESSORIES			
OPTION		Description	Retail
Accessories			
Ventilator Option			
☐	400VT12	Vent Tray Fixed (<i>Pulmonetic LTV Models</i>)	\$609.00
O2 Tank Holder			
☐	400A718	O2 Tank Holder Mounted to Backrest (<i>D and E Styles</i>)	\$391.00
Wheelchair Tray Table			
☐	4007A09	Tray, Swing Away	\$381.00
Medial Thigh Support			
☐	4008A80	Medial Thigh Support	\$244.00
☐	4008A81	Pommel Swing Away Hardware	\$358.00
Backpack			
☐	400A1	Black Utility	\$112.00
Misc. Accessories			
☐	4008A59	Transfer Bar (<i>Pair</i>)	\$493.00
☐	4008A52	Ram Mount Phone	\$77.00
☐	4008A53	Ram Mount Tablet	\$308.00
☐	4008A54	Ram Mount Cup Holder	\$77.00
☐	4008A58	Rear View Mirror	\$77.00
☐	400A611	Crutch Holder	\$115.00
☐	4008A73	Bag For Armrest	\$115.00
☐	4008A75	Reflectors (<i>Backrest shroud mounted</i>)	\$115.00
☐	4009A10	Utility Hooks	\$58.00
☐	4008A93	Wheelchair Cover	#N/A

JOYSTICK / ELECTRONICS				
OPTION		Type	Description	Retail
Electronics Package				
☐	400EL01		VR-2 or R-net Electronics	NC
☐	400EL02		LiNX Electronics	NC
VR-2 Joystick				
☐	400J127	VR-2	VR-2 Joystick	NC
R-net Joystick				
☐	400J121	R-net	Std LED Joystick for R-net Controls	NC
☐	400J116	R-net	Color LCD (CJSM1)	NC
☐	400J118	R-net	High Resolution Color LCD, (CJSM2+BT) 3.5" Color Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control (Client Cost)	\$1,162.00
☐	400J119	R-net	Omit Remote Control (Required with arm mounted specialty controls)	NC
LiNX Joystick				
☐	400J143	LiNX	REM110 - Drive Only Remote Joystick Adaptive Load Compensation	NC
☐	400J144	LiNX	REM200 Series - LED Joystick Adaptive Load Compensation, LED icon-based menu, Two Rocker Switches	NC
☐	400J145	LiNX	REM400 - Touch Screen Color Display Joystick 3.5" Color Display, Highly Configurable Touch Interface, Two Multi-Purpose Buttons, Two Stereo Jacks, Audible Cues, Bluetooth, & Adaptive Load Compensation	\$1,339.00
☐	400J148	LiNX	REM400 - Touch Screen Color Display Joystick - N/C for multipower 3.5" Color Display, Highly Configurable Touch Interface, Two Multi-Purpose Buttons, Two Stereo Jacks, Audible Cues, Bluetooth, & Adaptive Load Compensation	NC
☐	400J146	LiNX	REM400 - Touch Screen Color Display Joystick with Toggle Switches 3.5" Color Display, Highly Configurable Touch Interface, Two Multi-Purpose Buttons, Programmable Toggle Switches, Two Stereo Jacks, Audible Cues, Bluetooth, & Adaptive Load Compensation	\$1,695.00
☐	400J149	LiNX	REM400 - Touch Screen Color Display Joystick with Toggle Switches - Multi-Power 3.5" Color Display, Highly Configurable Touch Interface, Two Multi-Purpose Buttons, Programmable Toggle Switches, Two Stereo Jacks, Audible Cues, Bluetooth, & Adaptive Load Compensation	NC
☐	400J147	LiNX	Omit Remote Control (Required with arm mounted specialty controls)	NC
Joystick Mounting				
☐	400JM34		Fixed on Armrest	NC
☐	400JM35		Parallel Swing Away	\$447.00
☐	400JM36		Link-It Mid Line Mount on Armrest (Swing Away)	\$447.00
☐	400JM49		Attendant (Available with R-net Std LED or LCD, or LiNX REM 100, 200, or 400)	NC
Joystick Mount (Location)				
☐	400HM16		Right Hand	NC
☐	400HM17		Left Hand	NC
Display Language Selection				
☐	400L1		English	NC
☐	400L3		French	NC
☐	400L7		Spanish	NC
Joystick Handle (Optional) N/A with HD Joystick				
☐	400JH13	Both	Ball Handle	\$110.00
☐	400JH15	R-net	T-Handle	\$110.00
☐	400JH24	Both	Bodypoint - Rubber Dome Joystick Handle	\$216.00
☐	400JH25	Both	Bodypoint - U-Shaped Joystick Handle (3 in)	\$216.00
☐	400JH26	Both	Bodypoint - Grooved Mushroom Joystick	\$216.00
☐	400JH30	LiNX	Chin Cup	\$179.00
☐	400JH31	LiNX	Golf Ball	\$179.00
Electronic Programmers				
☐	400PR13		R-net PC Wireless Programmer	\$660.00
☐	400PR14		LiNX Access Key (LAK)	\$240.00

Power Modules				
☐	400ER12	VR-2	Power Module, VR-2, 90A, Non-Expandable	NC
☐	400ER22	R-net	Power Module, R-net, 90A, Non-Expandable <i>(Only Available with Single Power)</i>	NC
☐	400ER13	R-net	Power Module, R-net, 90A, Expandable (Required for specialty controls or 2+ actuators)	\$102.00
	400EX11		Expandable 90A Controller - R-net	\$0.00
	400EX12		Harness for 90A Expandable Controller - R-net	\$0.00
☐	400ER24	LiNX	Power Module, LiNX, 120A, Non-Expandable	NC
☐	400ER25	LiNX	Power Module, LiNX, 120A, Expandable (Required for specialty controls or 2+ actuators)	\$1,760.00
	400EX19		Expandable 120A Controller - LiNX	\$850.00
	400EX20		Harness for 120A Expandable Controller - LiNX	\$695.00

JOYSTICK / ELECTRONICS

OPTION	Type	Description	Retail
--------	------	-------------	--------

Seating Control

Powered Seating Control

☐	400EO40		Omit Thru-drive	NC
☐	400EO37		Thru-Drive Control for 1 Actuator.	\$1,443.00
☐	400EO38		Thru-Drive Control for 2+ Actuators.	\$2,903.00

Switch Control Only *(Available with Omit Thru-Drive Only) (1 Actuator for Tilt or Elevate only)*

☐	400EO33	R-net	Dual Toggle for 1 Actuator	NC
☐	400EO34	Both	2.5" Single Button Switch for 1 Actuator <i>(Button Switches are not Mounted)</i>	NC
☐	400EO41	R-net	Ctrl+5 Switch Box with Voice Activation for 2 Actuators or more	NC
☐	400EO42	LiNX	10 Way Switch Box with Voice Activation for 1 Actuator or more	NC

Additional Power Seating Control *(Available with Thru-drive Only)*

☐	400QG11	R-net	Ctrl+5 Switch Box with Voice Activation	\$1,078.00
☐	400QG12	R-net	Ctrl+5 Switch Box with Voice Activation <i>(Multi-Power Only)</i>	NC
☐	400QG16	LiNX	10 Way Switch Box for Actuator Control	NC

Dual Toggle Mounting Location

☐	400TN1		Toggle Mount, Right	NC
☐	400TN2		Toggle Mount, Left	NC

Power Seating Control Mounting Location *(Ctrl+5 or 10 Way) (If mounted on same side as joystick, it will swing away)*

☐	400EY19		Right Armrest	NC
☐	400EY20		Left Armrest	NC
☐	400EY21		Back Cane, Right	\$0.00
☐	400EY24		Back Cane, Left	\$0.00

Ctrl+5 Switch Type

☐	400EY22	R-net	Push Buttons	NC
☐	400EY23	R-net	Toggle Switches	\$356.00

Attendant Controls *(Prior Approval Required)*

☐	4000A05	R-net	Dual Control Module for Attendant	\$1,225.00
☐	4000A06	R-net	R-net Color LCD (CJSM1), Attendant Control Only	\$1,225.00
☐	4000A11	LiNX	LiNX Attendant Control Unit (ACU200)	\$1,195.00

Attendant Mounting Location

☐	400DN5		Right Back Cane	NC
☐	400DN6		Left Back Cane	NC

JOYSTICK / ELECTRONICS

OPTION	Type	Description	Retail
--------	------	-------------	--------

Programming

Assignable Buttons

Assignable Button feature allows for customization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

☐	400PM12	R-net	Include Assignable buttons	NC
--------------------------	---------	-------	----------------------------	----

Power Centermount Footplate Extension Option

☐	400IE03	BOTH	Include 7.5" Powered footplate extension for transfer assist (Electronics programmed to allow footplates to be extended with or without legrest being raised.) Only Available with Powered Articulation	\$592.00
--------------------------	---------	------	---	----------

Memory Seating

SitSmarter Memory Seating allows for positions to be saved and then can be accessed with a single command, moving all related actuators at once to the preset location. Positions can be locked for clinical goals and have others unlocked for easily adjustable positions for comfort and functionality that can be adjusted by the user. Angles can be shown on screen to provide live feedback.

☐	400PM19	LiNX	SitSmarter Memory Seating Positions	NC
--------------------------	---------	------	-------------------------------------	----

SitSmarter not compatible with Recline Only configurations

R-net Omni2 Display Options (For Expandable Controls, choose All That Apply) (AADL W100 Grant)

☐	400ON14	R-net	R-Net Omni2 Only Input Control Module and 3.5" LCD Display w/ Built-In Infrared and Bluetooth Choose Mount Below	\$2,384.00
☐	400ON15	R-net	R-Net Omni2 Bundled Input Control Module And 3.5" LCD Display w/ Built-In Infrared and Bluetooth. (Bundled Omni2 available w/ Switch-It Input Devices)	NC

LiNX REM500 Display Options (For Expandable Controls, choose All That Apply)

☐	400AI01	LiNX	LiNX REM500 Only Input Control Module and 3.5" Touch LCD Display with 2 Stereo Jacks, 2 Buttons, Audible Cues and Bluetooth (Choose Mount Below)	\$2,250.00
☐	400AI02	LiNX	LiNX REM500 Bundled Input Control Module and 3.5" Touch LCD Display with 2 Stereo Jacks, 2 Buttons, Audible Cues and Bluetooth (Bundled REM500 available w/ Switch-It Input Devices)	NC

Display Mount (If Used with Standard Joystick, SunDisk, or Bullet Hand Tray, Must Be Mounted Opposite Side) (AADL W100 Grant)

☐	400DS14		Link-it with Parallel Swing Away Mount, Right (Omni2 or REM500)	\$698.00
☐	400DS15		Link-It with Parallel Swing Away Mount, Left (Omni2 or REM500)	\$698.00

SPECIALTY INPUT DEVICES

OPTION	Description	Retail
--------	-------------	--------

Hand Control Packages (AADL W100 Grant)

(Requires Standard Joystick Omit / Only one hand control device can be ordered)

Proportional Input Devices

☐	400EU14	Switch-It MicroGuide	\$4,574.00
☐	400EU17	Switch-It VersaGuide	\$2,557.00
☐	400EU18	Switch-It VersaGuide EZ	\$3,023.00

Mounting Position (Required)

☐	400EM04	Link-It to the Backrest (Not Available w/ Omit Back Cushion Type)	\$1,338.00
☐	400EM05	Link-It to the Fixed Joystick Mount	\$416.00
☐	400EM06	Link-It to the Parallel Swing Away Joystick Mount	\$416.00
☐	400EM07	Link-It to the Midline Swing Away	\$1,480.00
☐	400EM08	Omit Mount	NC

Mounting Platforms (Optional) Not available with Link-It to the Backrest

☐	400MG16	SunDisk Hand Tray (Only available w/ MicroSeries joysticks)	\$416.00
☐	400MG15	Bullet Hand Tray	\$416.00

Mounting Location (Required)

☐	400IM09	Right	NC
☐	400IM10	Left	NC

Switch Options

See Switch Selections (Maximum 2)

Head Control Packages (Select headrest mount on next page) (Seating Benefits)

Only one head control device can be ordered.

Proportional Input Devices

☐	400EU25	Switch-It DualPro Proportional Head Array - Facial Pads	\$7,057.00
☐	400EU26	Switch-It DualPro Proportional Head Array - Spot Pads	\$7,057.00
☐	400EU51	Vigo Head Drive	\$6,293.00

Switched Input Devices

☐	400EU28	Switch-It 3-Switch Head Array - Short Curved Wings	\$6,154.00
☐	400EU30	Switch-It 3-Switch Head Array - Long Curved Wings	\$6,154.00

Switch Options

See Switch Selections (maximum 2)

SPECIALTY INPUT DEVICES			
OPTION		Description	Retail
Sip 'N Puff (AADL W100 Grant)			
☐	400ID73	Switch-It Sip 'N Puff Combo Head Array with Link-It Hardware <i>(Requires headrest mount)</i>	\$7,037.00
☐	400ID74	Switch-It OptiStop with Sip 'N Puff	\$2,741.00
Headrest Mounts (Applies to Head/Sip 'N Puff control packages) (AADL W100 Grant)			
☐	400IJ13	Whitmyer Linx Headrest Mount with Quick Release	\$247.00
☐	400IJ14	Whitmyer Linx-2 Headrest Mount with Quick Release	\$298.00
☐	400IJ15	Whitmyer AXYS Headrest Mount with Quick Release	\$247.00
Chin Control Packages (AADL W100 Grant)			
Proportional Input Devices			
☐	400EU34	Switch-It MicroGuide <i>(Switches mounted to the headrest only)</i>	\$4,574.00
☐	400EU36	Switch-It VersaGuide	\$2,557.00
☐	400EU37	Switch-It VersaGuide EZ	\$3,023.00
Mounting Position			
☐	400EM09	Link-It To Chin Harness, Small <i>(150mm / 6")</i>	\$1,760.00
☐	400EM10	Link-It To Chin Harness, Large <i>(220mm / 8")</i>	\$1,760.00
☐	400EM11	Link-It To Headrest <i>(Swing Away)</i>	\$1,760.00
☐	400EM16	Omit Mount	NC
Input Device Mounting Location			
☐	400IM13	Right	NC
☐	400IM14	Left	NC
Switch Options			
-		See Switch Selections <i>(Minimum 1 for the Power Chin Boom, Maximum 3)</i>	-
Safety Devices (AADL W100 Grant)			
☐	400SA06	Switch-It Remote Stop Switch - XLR Connector	\$1,317.00
Switched Driving Packages (AADL W100 Grant)			
☐	400SW44	Switch-It Cool Cube with DB9 connector for proximity / Proportional switch driving. <i>(Pro Spot Proportional Switch's are required if proportional driving is desired)</i> <i>(For Non-Proportional driving, Single Switches can be used)(For Future Development)</i>	\$3,792.00

Specialty Input Devices

OPTION	Description	Retail
--------	-------------	--------

Switches

Single switches can be programmed to activate any function of the wheelchair including: On/Off, Mode, Profile, Speeds, Seating, Memory Positions, Horn, Lights & Indicators, IR functions, Bluetooth mouse clicks, Environmental Controls when plugged into assignable port.

Switches can connect via: Joysticks (Standard, MicroGuide, Versa series) OMNI, Cool Cube, Head controls

Switch ordering instructions:

1. Select Switch type and quantity in section "A" using columns below. (Max 2 switches. see limitations below)
2. Select Switch Mounting location from Section "B" using corresponding switch column.
3. Select Switch Mounting Side from Section "C" using corresponding switch column.

Section A - Switch Selections

Single Function Switches - Mechanical Switches

2.5" Single Switch - Black

2.5" Single Switch - Blue

2.5" Single Switch - Green

2.5" Single Switch - Red

2.5" Single Switch - Yellow

Egg Switch - Black

Egg Switch - Blue

Egg Switch - Green

Egg Switch - Red

Egg Switch - Yellow

MicroLite (Includes color labels)

Mini Button 1 - Blue (Velcro)

Mini Button 1 - Green (Velcro)

Mini Button 1 - Red (Velcro)

1" Single Switch - Black

1" Single Switch - Blue

1" Single Switch - Green

1" Single Switch - Orange

1" Single Switch - Red

1" Single Switch - Yellow

Electronic Switches

Flex Switch (Headrest mount only)

Pro Spot - Proportional Switches

Pro Spot Proportional Switch - Red

Pro Spot Proportional Switch - Green

Pro Spot Proportional Switch - Blue

Pro Spot Proportional Switch - Yellow

Pro Spot Proportional Switch - Orange

Section B - Switch Mounting Position

Link-It to the Backrest (Up to 1 per side)

Headrest Mount - Gooseneck & Rigid (up to 2 per side)

Link-It to the Fixed Joystick Mount (Up to 2 per side)

Link-It to the Parallel Swing Away Joystick Mount (Up to 2 per side)

Link-It to the Seat Rail (Up to 2 per side)
(N/A w/ Knee & Thigh Supports)

Omit Mount (includes Velcro)

		Switch / Quantity			
		Switch 1	\$	Switch 2	\$
☐	400IR31	\$110.00	☐	400IR56	\$110.00
☐	400IR32	\$110.00	☐	400IR57	\$110.00
☐	400IR33	\$110.00	☐	400IR58	\$110.00
☐	400IR34	\$110.00	☐	400IR59	\$110.00
☐	400IR35	\$110.00	☐	400IR60	\$110.00
☐	400IR36	\$110.00	☐	400IR61	\$110.00
☐	400IR37	\$110.00	☐	400IR62	\$110.00
☐	400IR38	\$110.00	☐	400IR63	\$110.00
☐	400IR39	\$110.00	☐	400IR64	\$110.00
☐	400IR40	\$110.00	☐	400IR65	\$110.00
☐	400IR41	\$115.00	☐	400IR66	\$115.00
☐	400IR42	\$308.00	☐	400IR67	\$308.00
☐	400IR43	\$308.00	☐	400IR68	\$308.00
☐	400IR44	\$308.00	☐	400IR69	\$308.00
☐	400IR45	\$388.00	☐	400IR70	\$388.00
☐	400IR46	\$388.00	☐	400IR71	\$388.00
☐	400IR47	\$388.00	☐	400IR72	\$388.00
☐	400IR48	\$388.00	☐	400IR73	\$388.00
☐	400IR49	\$388.00	☐	400IR74	\$388.00
☐	400IR50		☐	400IR75	
☐	400IR54	\$542.00	☐	400IR79	\$542.00
☐	400IR181	\$542.00	☐	400IR186	\$542.00
☐	400IR182	\$542.00	☐	400IR187	\$542.00
☐	400IR183	\$542.00	☐	400IR188	\$542.00
☐	400IR184	\$542.00	☐	400IR189	\$542.00
☐	400IR185	\$542.00	☐	400IR190	\$542.00
☐	400EM17	\$1,338.00	☐	400EM23	\$1,338.00
☐	400EM18	\$1,760.00	☐	400EM24	\$1,760.00
☐	400EM19	\$416.00	☐	400EM25	\$416.00
☐	400EM20	\$416.00	☐	400EM26	\$416.00
☐	400EM21	\$556.00	☐	400EM27	\$556.00
☐	400EM22	NC	☐	400EM28	NC

Switches Continued

Section C - Switch Mounting Location

Right side of chair
Left side of chair

☐	400IM17	NC	☐	400IM19	NC
☐	400IM18	NC	☐	400IM20	NC

Joystick Integrated Switches (Pair)

☐	400JO03	-	Proximity Switches on Flex Mount, Pair (Used with VersaGuide/VersaGuide EZ)	\$1,403.00
--------------------------	---------	---	---	------------

Syncra Seating - Measuring Guide

JAY Syncra Back



Backrest Width
Max User Trunk Width

16"	18"	20"	22"
13.5"	15.5"	17.5"	19.5"


JAY Swing-Away Lateral Thoracic Adjustment Range

JAY Syncra Back Width	Width Range	Height range from bottom of shell to seat pan
16"	10" - 24"	12" - 19.5"
18"	12" - 26"	12.5" - 21.5"
20"	14" - 28"	12.5" - 22"
22"	16" - 30"	13.5" - 23"

Syncra Fixed Lateral Support Accessory Range

Syncra Back Width	Width Range (loaded)	Height range from bottom of shell to seat pan	Contour Depth Range
16"	16.4"-16.8"	6"-8"	3"-5"
18"	18.4"-18.8"	6"-8"	3"-5"
20"	20.4"-20.8"	6"-8"	3"-5"
22"	22.4"-22.8"	6"-8"	3"-5"

Syncra Back Mounting Height Diagram



JAY Syncra Cushion



JAY Swing-Away Lateral Thoracic Adjustment Range

JAY Syncra Cushion Width	Width Range
16"	8.5" - 16.0"
17"	9.5" - 17.0"
18"	10.5" - 18.0"
19"	11.5" - 19.0"
20"	12.5" - 20.0"
22"	14.5" - 22.0"



1000 Creditstone Road, Unit 2, Concord ON L4K 4P8

Customer Service: (800) 263-3390
E-mail: cscanada@sunmed.com
www.sunrisemedical.ca

Specifications are subject to change without notice.

© 10/2026 Sunrise Medical Canada Inc. MK-129808-CAN Rev 10