



ORDER FORM

October 2024 Rev 1

Q700M


with SEDEO[®] LITE[™]

AADL -W581
EIPW25

| Customer | | Ship To | |
|-----------|--|----------|-------------|
| Account # | | Name | |
| Date | | Address | |
| PO# | | City | |
| Buyer | | Province | Postal Code |
| Mark For | | Phone | |

Client Measurements

| | |
|----------------------------|----|
| a) Top of head to seat | in |
| b) Top of shoulder to seat | in |
| c) Underarm to seat | in |
| d) Elbow to seat | in |
| e) Top of leg to seat | in |
| f) Shoulder width | in |
| g) Chest width | in |
| i) Hip width | in |
| j) Actual sitting depth | in |
| m) Foot length | in |
| n) Knee to foot | in |
| o) Seat width | in |



Notes

Chair Specifications

| | |
|--|--------------------------------------|
| Weights Starting at: (with Batteries) | NPC: 294 lbs. SPC: 318 lbs. |
| Power Seating | 50° of CG Tilt 12" seat elevator |
| Seat Type | SEDEO® LITE™ |
| Seat Size Range | Width 15" to 20" Depth 16" to 22" |
| Seat-to-Floor Height Range | 16.5" to 19" |

Mid-Wheel Drive Base Specifications

| | |
|--------------------------|--|
| Drive Wheel Position | Mid-wheel drive |
| User Weight Capacity | 300 lbs. 350 lbs. HD |
| Base Size | Width 25", Length 36.5" |
| Minimum Turning Radius | 19.8" (excludes footplates) |
| Drive Control | R-net PM 120A / R-net PM 90A |
| Driving Range | Up to 18 miles (w/ Group 24 Batteries) |
| Drive Motors | 4-pole 6 mph, 5 mph HD, 8 mph |
| Drive Wheels | 14" w/flat-free inserts (black) |
| Casters | 6" front / 6" rear (black) |
| Charger | 10A Dual Mode (Compact) |
| Battery Size | Group 24 / Group 34 |
| Colour | 9 shroud colours and 7 colour accents |
| Maximum Ground Clearance | 3" |

| BASE | | | |
|---|--|--|-------------|
| OPTION | | DESCRIPTION | RETAIL |
| Standard Duty Mid Wheel Drive Base <i>End Item Number EIPW25 (LITE™ Seat)</i> | | | \$ |
| 250PC267 | | Quickie Q700 M NPC No Power Option for SEDEO® LITE™ Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5". (AADL Contract Price) | \$8,340.00 |
| 250PC268 | | Quickie Q700 M SPC Single Power Option for SEDEO® LITE™ Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5". (AADL Contract Price) | \$8,340.00 |
| Heavy Duty Mid Wheel Drive Base <i>End Item Number EIPW25 (LITE™ Seat)</i> | | | \$ |
| 250PC269 | | Quickie Q700 M NPC HD No Power Option for SEDEO® LITE™ Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 5mph, Weight Capacity = 350lbs., STFH = 16.5". (AADL Contract Price) | \$9,822.00 |
| 250PC270 | | Quickie Q700 M SPC HD Single Power Option for SEDEO® LITE™ Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 5mph, Weight Capacity = 350lbs., STFH = 16.5". (AADL Contract Price) | \$9,822.00 |
| Standard Duty Mid Wheel Drive Base <i>End Item Number EIPW25 (LITE™ Seat)</i> | | | \$ |
| 250PC271 | | Quickie Q700 M NPD No Power Option for SEDEO® LITE™ Seat Included in Base Pricing : Q700 M Base, R-net 120amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5". (Client Cost) | \$10,622.00 |
| 250PC272 | | Quickie Q700 M SPD Single Power Option for SEDEO® LITE™ Seat Included in Base Pricing : Q700 M Base, R-net 120amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5". (Client Cost) | \$10,622.00 |

| BASE | | | |
|---|--|--|------------|
| OPTION | | DESCRIPTION | RETAIL |
| User Weight | | | \$ |
| 250WG23 | | Up to 300 lbs. | NC |
| 250WG15 | | 301 to 350 lbs. | NC |
| Motor Package | | | \$ |
| 250PF98 | | 4-pole 8 kph / 5 mph | NC |
| 250PF99 | | 4-pole 10 kph / 6 mph | NC |
| 250PF101 | | 4-pole 12.5 kph / 8 mph (Client Cost) | \$1,249.00 |
| Tracking | | | \$ |
| 250MI01 | | SureTrac® Module & Motor Encoders <i>(Improves drive tracking)</i> | \$1,671.00 |
| Drive Wheels | | | \$ |
| 250DW58 | | Drive Wheel, 14", Dual Density, Knobby, Black | NC |
| Caster Wheels | | | \$ |
| 250CW55 | | Caster Wheels 6", Solid, Black | NC |
| Shroud Colour | | | \$ |
| 250C119 | | Matte Black | NC |
| 250C125 | | Pearl White | NC |
| 250C165 | | Candy Apple Red | NC |
| 250C166 | | Candy Blue | NC |
| 250C167 | | Candy Purple | NC |
| 250C176 | | Green Apple | NC |
| 250C178 | | Pearl Pink | NC |
| 250C182 | | Silver | NC |
| Colour Accents | | | \$ |
| <i>Note: Colour Accents and Shroud Colours may not match.</i> | | | |
| 250EZ20 | | Apple Green | NC |
| 250EZ22 | | Black | NC |
| 250EZ21 | | Blue | NC |
| 250EZ23 | | Orange | NC |
| 250EZ24 | | Red | NC |
| 250EZ19 | | White | NC |
| 250EZ25 | | Silver | NC |

| BASE | | | |
|--------------------------|--|---------------------|----------|
| OPTION | | DESCRIPTION | RETAIL |
| ISO Transit Tie Downs | | | \$ |
| 250T050 | | ISO Transit Package | NC |
| LED Lights (Client Cost) | | | \$ |
| 250LI08 | | LED Fender Lights | \$636.00 |
| 250LI02 | | Omit LED Lights | NC |
| Seat Heights | | | \$ |

Note: Available seat-to-floor heights are dependent on 2 options: **1- Battery Size** and **2- Powered Seating**.
Use the matrix below to determine available ranges for your desired options. Measured from front of seat pan to floor with 0° frame angle.

| | Seat Height | Power Seating | No Power Seating | |
|----------|-------------|---------------|------------------|----|
| 250HT103 | 16.5" | Group 34 | Group 34 | NC |
| 250HT55 | 17" | Group 34 | GP24 or GP34 | NC |
| 250HT104 | 17.5" | Group 34 | GP24 or GP34 | NC |
| 250HT56 | 18" | GP24 or GP34 | GP24 or GP34 | NC |
| 250HT105 | 18.5" | GP24 or GP34 | GP24 or GP34 | NC |
| 250HT57 | 19" | GP24 or GP34 | GP24 or GP34 | NC |

Batteries and Chargers

Group 24 Battery Options (Recommended for max daily range)

| | | | |
|---------|--|-----------------------------------|----------|
| 250AG25 | | GP24 Flat Top Pair w/Strap MK Gel | \$941.00 |
| 250AG26 | | GP24 Flat Top Pair w/Strap AGM | \$901.00 |

Group 34 Battery Options (Required for lower seat-to-floor heights)

| | | | |
|---------|--|------------------------------------|----------|
| 250AG27 | | GP34 Flat Top Pair w/Straps MK Gel | \$941.00 |
| 250AG29 | | GP34 Flat Top Pair w/Straps AGM | \$850.00 |

Omit Batteries

| | | | |
|---------|--|---|----|
| 250AG24 | | Battery Harnesses Only (Omit Batteries) | NC |
|---------|--|---|----|

Chargers

| | | | |
|---------|--|-------------------------|----|
| 250BJ28 | | 10A Dual Mode (Compact) | NC |
|---------|--|-------------------------|----|

POWERED SEATING

| OPTION | DESCRIPTION | RETAIL |
|--------|-------------|--------|
|--------|-------------|--------|

Power Seating Options \$

| Seat Elevator | | (Client Cost) |
|----------------------|---|-----------------------------------|
| 250PI08 | C-Me® Seat Elevator, Powered, 300 mm / 12" (Not available with Group 4) <i>C-Me speed up to 3 mph.</i> | \$3,426.00 |
| Seat Tilt | | (AADL W992 for Power Tilt) |
| 250PA35 | Posterior 50° Tilt | \$8,535.00 |
| 250PA36 | Posterior 50° HD Tilt | \$8,791.00 |

Pre Set Seat Frame Angle (No Power Seating Options, Increases Seat to Floor Height) \$

| | | |
|---------|----------------------|----|
| 250FA10 | 0° (Adds 0" to STFH) | NC |
| 250FA3 | 3° (Adds 1" to STFH) | NC |

SEDEO® LITE™ Backrest Frame (Angle Adjustable 86° - 106°) \$

| Back Frame Width | | |
|-------------------------|---|----|
| 250BF21 | Width Adjustable, 385 - 485 mm / 15 - 19" | NC |
| 250BF22 | Width Adjustable, 460 - 510 mm / 18 - 20" | NC |

| Back Cane Height | | |
|-------------------------|--|----|
| 250B263 | Standard, 18" - 20" | NC |
| 250B264 | Tall, 21"-23" | NC |
| 250PH30 | Height Adjustable, Short w/Trailing Push Handles, 20" to 22" | NC |
| 250PH31 | Height Adjustable, Tall w/Trailing Push Handles, 23" to 26" | NC |

Backrest Frame Option \$

| | | |
|---------|--------------------------|----------|
| 250BF23 | Flip Down, Quick Release | \$176.00 |
| 250BF24 | Removeable | NC |

Back Cushion Type \$

| | | |
|----------|--|----|
| 250BT148 | SEDEO® LITE™ Sling Back with Tension Adjustable Straps (Choose Seating in Section A) | NC |
| 250BT133 | OMIT (For personalized mounting. No headrest mount) | NC |

A - SEDEO® LITE™ Sling Back Dimensions \$

| Installed Back Height - (From Pan to Top of Back) | | |
|--|-----|----|
| 250JR71 | 18" | NC |
| 250JR72 | 19" | NC |
| 250JR73 | 20" | NC |
| 250JR74 | 21" | NC |
| 250JR75 | 22" | NC |
| 250JR76 | 23" | NC |

POWERED SEATING

| OPTION | DESCRIPTION | RETAIL |
|--------|-------------|--------|
|--------|-------------|--------|

Seat Size \$

Seat Pan Width

| | | |
|--------|-----|----|
| 250R78 | 15" | NC |
| 250R53 | 16" | NC |
| 250R54 | 17" | NC |
| 250R55 | 18" | NC |
| 250R56 | 19" | NC |
| 250R57 | 20" | NC |

Seat Pan Depth (Front of seat pan to front of back cushion – effective depth)

| | | |
|---------|-----|----|
| 250RF51 | 16" | NC |
| 250RF45 | 17" | NC |
| 250RF46 | 18" | NC |
| 250RF47 | 19" | NC |
| 250RF48 | 20" | NC |
| 250RF49 | 21" | NC |
| 250RF50 | 22" | NC |

POSITIONING SUPPORTS & ACCESSORIES

| OPTION | DESCRIPTION | RETAIL |
|---|---|-----------|
| Knee and Thigh Swing-Away Support Pads <i>Seating Benefits</i> | | \$ |
| Lateral Knee Support Pads and Hardware (<i>must pick pad and mount</i>) | | |
| 250XY05 | Lateral Knee Support Pad, Right (4" round) | \$243.00 |
| 250XY06 | Lateral Knee Support Pad, Left (4" round) | \$243.00 |
| 250XX08 | Lateral Swing-Away Bracket Knee - Right | \$320.00 |
| 250XX09 | Lateral Swing-Away Bracket Knee - Left | \$320.00 |
| Lateral Thigh Support Pads and Hardware (<i>must pick pad and mount</i>) | | |
| 250KF22 | Lateral Thigh Support Pad Short, Right (5"x5") | \$144.00 |
| 250KF23 | Lateral Thigh Support Pad Short, Left (5"x5") | \$144.00 |
| 250KF24 | Lateral Thigh Support Pad Long, Right (5"x7") | \$144.00 |
| 250KF25 | Lateral Thigh Support Pad Long, Left (5"x7") | \$144.00 |
| 250KF26 | Lateral Thigh Removable Mount - Right | \$320.00 |
| 250KF27 | Lateral Thigh Removable Mount - Left | \$320.00 |
| Armrest Structure | | \$ |
| 250EG35 | Reclining, Flip Back (<i>Height Adjustable 8" to 13"</i>) | NC |
| 250EG36 | Reclining, Flip Back with 2-point post locking arm (<i>Height Adjustable 11" to 13" x .625"</i>) | \$346.00 |
| 250EG34 | Single Post, SEDEO LITE, Removable (<i>N/A with thigh support</i>), (<i>Height Adjustable 10" to 15"</i>) | NC |
| Arm Pads | | \$ |
| SEDEO® Armrest Pad | | |
| 250AB62 | 300 mm / 12", Integral Foam, Right | NC |
| 250AE62 | 300 mm / 12", Integral Foam, Left | NC |
| 250AB63 | 400 mm / 16", Integral Foam, Right | NC |
| 250AE63 | 400 mm / 16", Integral Foam, Left | NC |
| 250AB61 | SEDEO LITE Pad, Right | NC |
| 250AE61 | SEDEO LITE Pad, Left | NC |
| Comfort Arm Pad (Foam) | | |
| 250AB46 | 2" Wide Right | \$111.00 |
| 250AE46 | 2" Wide Left | \$111.00 |
| Comfort Arm Pad (Gel) | | |
| 250AB48 | 2" Wide Right | \$210.00 |
| 250AE48 | 2" Wide Left | \$210.00 |
| Ottobock® Arm Troughs | | |
| 250AE06 | Ottobock Arm Trough, Med (12"), Left | \$137.00 |
| 250AB06 | Ottobock Arm Trough, Med (12"), Right | \$137.00 |
| 250AE07 | Ottobock Arm Trough, Large (14"), Left | \$137.00 |
| 250AB07 | Ottobock Arm Trough, Large (14"), Right | \$137.00 |
| 250OL01 | Ottobock Hand Pad, Med (7"), Left | \$69.00 |
| 250OR01 | Ottobock Hand Pad, Med (7"), Right | \$69.00 |
| 250OL02 | Ottobock Hand Pad, Large (9"), Left | \$69.00 |
| 250OR02 | Ottobock Hand Pad, Large (9"), Right | \$69.00 |
| 250AE78 | Ottobock One-Piece Arm Trough (18"), Left | \$195.00 |
| 250AB78 | Ottobock One-Piece Arm Trough (18"), Right | \$195.00 |

POSITIONING SUPPORTS & ACCESSORIES

| OPTION | DESCRIPTION | RETAIL |
|---|--|------------|
| Fixed Centermount Legrest \$ | | |
| 250H183 | SEDEO® Fixed Centermount Legrest (3" to 17.5" Leg Length) | \$346.00 |
| 250H184 | Fixed Centermount Legrest, Short (Includes Calf Pads & Single Footplate) (11.5" to 14.5" Leg Length) | NC |
| 250H185 | Fixed Centermount Legrest, Long (Includes Calf Pads & Single Footplate) (13.5" to 18.5" Leg Length) | NC |
| 250H136 | Angle Adjustable Centermount Legrest (11.5" to 18.5" Leg Length) | \$596.00 |
| Power Centermount Legrest (Client Cost) \$ | | |
| 250H175 | Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension (Mechanical Articulation) (5" to 17.5" Leg Length) | \$2,450.00 |
| 250H109 | Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension (Powered Articulation) (3" to 17.5" Leg Length) | \$3,384.00 |
| Power Centermount Legrest Accessories (Client Cost) \$ | | |
| SEDEO Centermount Footplates | | |
| 250F127 | Dual, Height, Width and Angle Adjustable, Standard (14"W x 10"D) | \$280.00 |
| 250F128 | Dual, Height, Width and Angle Adjustable, Large (16"W x 12"D) | \$280.00 |
| Angle Adjustable Centermount Footplates | | |
| 250F45 | Single / Dual Footplate | \$280.00 |
| 250F44 | Single Platform Footplate | \$280.00 |
| SEDEO Centermount Calf Pads | | |
| 250DI34 | Calf Pads, Pair, Standard Integral Foam | NC |
| 250DI35 | Calf Pads, Pair, Padded | \$351.00 |
| Angle Adjustable Centermount Calf Pads | | |
| 250DI49 | Calf Pads, Pair, Standard Integral Foam | NC |
| 250DI50 | Extra Large Single Calf Pad | \$351.00 |
| SEDEO® PRO™ Swing-Away Legrest (Pair) (250 H179 is Client Cost) \$ | | |
| 250H177 | Legrests, Swing-Away, 70°, SEDEO® PRO™ (14" to 21" Leg Length) | NC |
| 250H178 | Legrests, Manual Elevating, Swing-Away, SEDEO® PRO™ (Includes Padded Calf Pads) (14".5 to 21.5" Leg Length) | \$596.00 |
| 250H179 | Legrests, Power Elevating, Swing-Away, SEDEO® PRO™ (Includes Padded Calf Pads) (15.5" to 22.5" Leg Length) | \$3,384.00 |
| 250DI76 | Omit Legrests | NC |
| Swing-Away Legrest Footplates | | |
| 250F369 | Angle Adjustable Flip-Up Footplates | \$280.00 |
| Swing-Away Legrest Calf Supports | | |
| 250DI42 | Calf Support, Depth Adjustable, Padded (6.5" to 8" Wide, Depending On Chair Width, x 9.5" Tall) | \$504.00 |
| Swing-Away Legrest Options | | |
| 250DI44 | Heel Loops | NC |
| 250KG4 | Knee Pads, Integral Foam (2.5" Wide x 4.5" Tall) | \$278.00 |
| SEDEO® Lower Extremity Swing-Away Support | | |
| 250VD01 | Lower Extremity Support, Padded, Swing-Away, SEDEO®, Right | \$504.00 |
| 250VD02 | Lower Extremity Support, Padded, Swing-Away, SEDEO®, Left | \$512.00 |

POSITIONING SUPPORTS & ACCESSORIES

| OPTION | DESCRIPTION | RETAIL |
|--|--|----------|
| Pelvic Positioning Belts | | |
| 250A837 | Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2" | NC |
| 250A839 | Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2" | NC |
| 250A840 | Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2" | \$101.00 |
| 250A845 | Bodypoint – Evoflex Pelvic Stabilizer Belt | \$213.00 |
| 250A846 | Bodypoint – Evoflex Pelvic Stabilizer Belt with Extension | \$267.00 |
| Accessories | | |
| | | \$ |
| Ventilator Option | | |
| 250VT12 | Vent Tray Fixed (<i>Pulmonetic LTV Models</i>) | \$600.00 |
| O2 Tank Option | | |
| 250A718 | O2 Tank Holder Mounted to Backrest (<i>D and E Styles</i>) | \$385.00 |
| Wheelchair Tray Table | | |
| 2507A09 | Tray, Swing-Away | \$375.00 |
| Backpack (Client Cost) | | |
| 250A1 | Black Utility | \$110.00 |
| Medial Thigh Support | | |
| 2508A80 | Medial Thigh Support | \$267.00 |
| 2508A81 | Pommel Swing-Away Hardware | \$384.00 |
| Misc. Accessories (Client Cost) | | |
| 2508A59 | Transfer Bar (<i>Pair</i>) | \$486.00 |
| 2508A52 | Ram Mount Phone | \$76.00 |
| 2508A53 | Ram Mount Tablet | \$303.00 |
| 2508A54 | Ram Mount Cup Holder | \$76.00 |
| 2508A58 | Rear View Mirror | \$76.00 |
| 250A611 | Crutch Holder | \$171.00 |
| 2508A73 | Bag for Armrest | \$76.00 |

| JOYSTICK / ELECTRONICS | | | |
|----------------------------------|---------|--|------------|
| OPTION | | DESCRIPTION | RETAIL |
| R-net Joystick (<i>Remote</i>) | | | \$ |
| 250J121 | | Standard LED Joystick (<i>w/o Screen Display</i>) | NC |
| 250J116 | | Colour Joystick (<i>2.3" Colour LCD Screen Display, CJSM1</i>) | NC |
| 250J118 | | High Resolution Colour Joystick,(CJSM2+BT) <i>3.5" Colour LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control (Client Cost)</i> | \$1,245.00 |
| 250J119 | | Omit Remote Control (<i>Required with arm mounted specialty controls</i>) | NC |
| Joystick Mounting | | | \$ |
| 250JM34 | | Fixed on Armrest | NC |
| 250JM35 | | Parallel Swing-Away | \$440.00 |
| 250JM36 | | Link-It Mid Line Mount on Armrest (<i>Swing-Away</i>) | \$440.00 |
| Joystick Mount (Location) | | | \$ |
| 250HM16 | | Right Hand | NC |
| 250HM17 | | Left Hand | NC |
| Display Language Selection R-net | | | \$ |
| 250L1 | | English | NC |
| 250L3 | | French | NC |
| 250L7 | | Spanish | NC |
| Joystick Handle (Optional) | | | \$ |
| 250JH13 | | Ball Handle | \$108.00 |
| 250JH15 | | T-Handle | \$108.00 |
| 250JH24 | | Bodypoint - Rubber Dome Joystick Handle | \$213.00 |
| 250JH25 | | Bodypoint - U-Shaped Joystick Handle (<i>4 inch</i>) | \$213.00 |
| 250JH26 | | Bodypoint - Grooved Mushroom Joystick Handle | \$213.00 |
| Electronic Programmings | | | \$ |
| 250PR13 | | R-net PC Wireless Programmer | \$696.00 |
| R-net Power Modules | | | \$ |
| 250ER11 | | Power Module, R-net, 120A, Non-Expandable | NC |
| 250ER10 | | Power Module, R-net, 120A, Expandable <i>Required for specialty controls or 3+ actuators</i> | NC |
| <input type="checkbox"/> | 250EX09 | Expandable 120A Controller - R-net | NC |
| <input type="checkbox"/> | 250EX10 | Harness for 120A Expandable Controller | NC |
| | 250ER18 | Power Module, R-net, 90A, Non-Expandable | NC |
| | 250ER17 | Power Module, R-net, 90A, Expandable <i>Required for specialty controls or 3+ actuators</i> | \$102.00 |
| <input type="checkbox"/> | 250EX15 | Expandable 90A Controller - R-net | NC |
| <input type="checkbox"/> | 250EX16 | Harness for 90A Expandable Controller | NC |

JOYSTICK / ELECTRONICS

| OPTION | DESCRIPTION | RETAIL |
|--|---|------------|
| Seating Control \$ | | |
| Powered Seating Control | | |
| 250E040 | Omit Thru-Drive <i>(Requires Switch Control)</i> | NC |
| 250E037 | Thru-Drive Control for 1 Actuator | \$1,422.00 |
| 250E038 | Thru-Drive Control for 2+ Actuators | \$2,860.00 |
| Switch Control <i>(Available with Omit Thru-Drive only)</i> | | |
| 250E033 | Dual Toggle for 1 Actuator | NC |
| 250E034 | 2.5" Single Button Switch for 1 Actuator <i>(Button switches are not mounted)</i> | NC |
| 250E041 | Ctrl+5 Switch Box with Voice Activation | NC |
| Additional Power Seating Control <i>(Available with Thru-Drive only)</i> | | |
| 250QG11 | Ctrl+5 Switch Box with Voice Activation | \$1,062.00 |
| Dual Toggle Mounting Location | | |
| 250TN1 | Toggle Mount, Right | NC |
| 250TN2 | Toggle Mount, Left | NC |
| Ctrl+5 Mounting Location <i>(Will swing away if mounted in same location as joystick)</i> | | |
| 250EY19 | Right Armrest | NC |
| 250EY20 | Left Armrest | NC |
| 250EY21 | Back Cane, Right | NC |
| 250EY24 | Back Cane, Left | NC |
| Ctrl+5 Switch Type | | |
| 250EY22 | Push Buttons | NC |
| 250EY23 | Toggle Switches | \$351.00 |
| Attendant Controls (Prior Approval Required) \$ | | |
| 2500A05 | Dual Control Module for Attendant | \$1,207.00 |
| 2500A06 | R-net Colour LCD (CJSM1), Attendant Control Only | \$1,207.00 |
| Attendant Mounting Location \$ | | |
| 250DN5 | Right Back Cane | NC |
| 250DN6 | Left Back Cane | NC |
| R-net Programming \$ | | |
| Assignable Buttons | | |
| <i>Assignable Button feature allows for personalization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.</i> | | |
| 250PM12 | Include Assignable Buttons | NC |
| Power Centermount Footplate Extension Option | | |
| 250IE03 | Include 7.5" Powered Footplate Extension for Transfer Assist <i>(Electronics programmed to allow footplates to be extended with or without legrest being raised.) Only available with Powered Articulation.</i> | \$583.00 |



Q700 M with SEDEO® LITE™

\$

| | |
|-------------------------------|--|
| Base Price | |
| Total of additional equipment | |
| Packaging and transport | |
| Total Price | |

Signature

