



ORDER FORM

October 2024 REV 1

Q700M


with SEDEO® LITE™

**AADL - W581
EIPW25**

Customer	Ship To
Account #	Name
Date	Address
PO#	City
Buyer	Province Postal Code
Mark For	Phone

Client Measurements

a) Top of head to seat	in
b) Top of shoulder to seat	in
c) Underarm to seat	in
d) Elbow to seat	in
e) Top of leg to seat	in
f) Shoulder width	in
g) Chest width	in
i) Hip width	in
j) Actual sitting depth	in
m) Foot length	in
n) Knee to foot	in
o) Seat width	in



Notes

Chair Specifications

Weights Starting at: (with Batteries)	NPC: 294 lbs. SPC: 318 lbs.
Power Seating	50° of CG Tilt 12" seat elevator
Seat Type	SEDEO® LITE™
Seat Size Range	Width 15" to 20" Depth 16" to 22"
Seat-to-Floor Height Range	16.5" to 19"

Mid-Wheel Drive Base Specifications

Drive Wheel Position	Mid-wheel drive
User Weight Capacity	300 lbs. 350 lbs. HD
Base Size	Width 25", Length 36.5"
Minimum Turning Radius	19.8" (excludes footplates)
Drive Control	R-net PM 120A / R-net PM 90A
Driving Range	Up to 18 miles (w/ Group 24 Batteries)
Drive Motors	4-pole 6 mph, 5 mph HD, 8 mph
Drive Wheels	14" w/flat-free inserts (black)
Casters	6" front / 6" rear (black)
Charger	10A Dual Mode (Compact)
Battery Size	Group 24 / Group 34
Color	9 shroud colors and 7 color accents
Maximum Ground Clearance	3"

BASE			
OPTION		DESCRIPTION	RETAIL
Standard Duty Mid Wheel Drive Base <i>End Item Number EIPW25 (LITE™ Seat)</i>			\$
250PC267		Quickie Q700 M NPC No Power Option for SEDEO® LITE™ Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5". (AADL Contract Price)	\$8,340.00
250PC268		Quickie Q700 M SPC Single Power Option for SEDEO® LITE™ Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5". (AADL Contract Price)	\$8,340.00
Heavy Duty Mid Wheel Drive Base <i>End Item Number EIPW25 (LITE™ Seat)</i>			\$
250PC269		Quickie Q700 M NPC HD No Power Option for SEDEO® LITE™ Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 5mph, Weight Capacity = 350lbs., STFH = 16.5". (AADL Contract Price)	\$9,822.00
250PC270		Quickie Q700 M SPC HD Single Power Option for SEDEO® LITE™ Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 5mph, Weight Capacity = 350lbs., STFH = 16.5". (AADL Contract Price)	\$9,822.00
High Performance Mid Wheel Drive Base <i>End Item Number EIPW25 (LITE™ Seat)</i>			\$
250PC271		Quickie Q700 M NPD No Power Option for SEDEO® LITE™ Seat Included in Base Pricing : Q700 M Base, R-net 120amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5". (Client Cost)	\$10,622.00
250PC272		Quickie Q700 M SPD Single Power Option for SEDEO® LITE™ Seat Included in Base Pricing : Q700 M Base, R-net 120amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5". (Client Cost)	\$10,622.00

BASE			
OPTION		DESCRIPTION	RETAIL
User Weight			\$
250WG23		Up to 300 lbs.	NC
250WG15		301 to 350 lbs.	NC
Motor Package			\$
250PF98		4-pole 8 kph / 5 mph	NC
250PF99		4-pole 10 kph / 6 mph	NC
250PF101		4-pole 12.5 kph / 8 mph (Client Cost)	\$1,249.00
Tracking			\$
250MI01		SureTrac® Module & Motor Encoders	\$1,671.00
Drive Wheels  P. 15			\$
250DW58		Drive Wheel, 14", Dual Density, Knobby, Black	NC
Caster Wheels  P. 15			\$
250CW55		Caster Wheels 6", Solid, Black	NC
Shroud Color			\$
250C119		Matte Black	NC
250C125		Pearl White	NC
250C165		Candy Apple Red	NC
250C166		Candy Blue	NC
250C167		Candy Purple	NC
250C176		Green Apple	NC
250C178		Pearl Pink	NC
250C182		Silver	NC
Color Accents			\$
<i>Note: Color Accents and Shroud Colors may not match.</i>			
250EZ20		Apple Green	NC
250EZ22		Black	NC
250EZ21		Blue	NC
250EZ23		Orange	NC
250EZ24		Red	NC
250EZ19		White	NC
250EZ25		Silver	NC

BASE				
OPTION		DESCRIPTION		RETAIL
ISO Transit Tie Downs				\$
250T050		ISO Transit Package		NC
LED Lights  P. 15 (Client Cost)				\$
250LI04		LED Lights and Indicators		\$1,207.00
250LI08		LED Fender Lights		\$636.00
250LI02		Omit LED Lights		NC

Seat Heights \$

Note: Available seat-to-floor heights are dependent on 2 options: **1- Battery Size** and **2- Powered Seating**.

Use the matrix below to determine available ranges for your desired options. Measured from front of seat pan to floor with 0° frame angle.

		Seat Height	Power Seating	No Power Seating	
250HT103		16.5"	Group 34	Group 34	NC
250HT55		17"	Group 34	GP24 or GP34	NC
250HT104		17.5"	Group 34	GP24 or GP34	NC
250HT56		18"	GP24 or GP34	GP24 or GP34	NC
250HT105		18.5"	GP24 or GP34	GP24 or GP34	NC
250HT57		19"	GP24 or GP34	GP24 or GP34	NC

Batteries and Chargers P. 15 \$

Group 24 Battery Options (Recommended for max daily range)

250AG25		GP24 Flat Top Pair w/Strap MK Gel		\$941.00
250AG26		GP24 Flat Top Pair w/Strap Brightway AGM		\$901.00

Group 34 Battery Options (Required for lower seat-to-floor heights)

250AG27		GP34 Flat Top Pair w/Straps MK Gel		\$941.00
250AG29		GP34 Flat Top Pair w/Straps Brightway AGM		\$850.00

Omit Batteries

250AG24		Battery Harnesses Only (Omit Batteries)		NC
---------	--	---	--	----

Chargers

250BJ28		10A Dual Mode (Compact)		NC
---------	--	-------------------------	--	----

POWERED SEATING

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

Power Seating Options **P. 15** \$
Seat Elevator (Client Cost)

250PI08	C-Me® Seat Elevator, Powered, 300 mm / 12" (Not available with HD or Group 4) <i>C-Me speed up to 3 mph.</i>	\$3,426.00
---------	---	------------

Seat Tilt (AADL W992 for Power Tilt)

250PA35	Posterior 50° Tilt	\$8,535.00
250PA36	Posterior 50° HD Tilt	\$8,791.00

Pre Set Seat Frame Angle (No Power Seating Options, Increases Seat to Floor Height) \$

250FA10	0° (Adds 0" to STFH)	NC
250FA3	3° (Adds 1" to STFH)	NC

SEDEO® LITE™ Backrest Frame (Angle Adjustable 86° - 106°) \$
Back Frame Width

250BF21	Width Adjustable, 385 - 485 mm / 15 - 19"	NC
250BF22	Width Adjustable, 460 - 510 mm / 18 - 20"	NC

Back Cane Height

250B263	Standard, 18" - 20"	NC
250B264	Tall, 21"-23"	NC
250PH30	Height Adjustable, Short w/Trailing Push Handles, 20" to 22"	NC
250PH31	Height Adjustable, Tall w/Trailing Push Handles, 23" to 26"	NC

Backrest Frame Option \$

250BF23	Flip Down, Quick Release	\$176.00
250BF24	Removeable	NC

Back Cushion Type \$

250BT148	SEDEO® LITE™ Sling Back with Tension Adjustable Straps (Choose Seating in Section A)	NC
250BT133	OMIT (For personalized mounting. No headrest mount)	NC

A - SEDEO® LITE™ Sling Back Dimensions \$
Installed Back Height - (From Pan to Top of Back)

250JR71	18"	NC
250JR72	19"	NC
250JR73	20"	NC
250JR74	21"	NC
250JR75	22"	NC
250JR76	23"	NC

POWERED SEATING

OPTION		DESCRIPTION	RETAIL
--------	--	-------------	--------

Seat Size

\$

Seat Pan Width

250R78		15"	NC
250R53		16"	NC
250R54		17"	NC
250R55		18"	NC
250R56		19"	NC
250R57		20"	NC

Seat Pan Depth (*Front of seat pan to front of back cushion – effective depth*)

250RF51		16"	NC
250RF45		17"	NC
250RF46		18"	NC
250RF47		19"	NC
250RF48		20"	NC
250RF49		21"	NC
250RF50		22"	NC

POSITIONING SUPPORTS & ACCESSORIES

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

Knee and Thigh Swing-Away Support Pads **Seating Benefits** \$

Lateral Knee Support Pads and Hardware *(must pick pad and mount)*

250XY05	Lateral Knee Support Pad, Right (4" round)	\$243.00
250XY06	Lateral Knee Support Pad, Left (4" round)	\$243.00
250XX08	Lateral Swing-Away Bracket Knee - Right	\$320.00
250XX09	Lateral Swing-Away Bracket Knee - Left	\$320.00

Lateral Thigh Support Pads and Hardware *(must pick pad and mount)*

250KF22	Lateral Thigh Support Pad Short, Right (5"x5")	\$144.00
250KF23	Lateral Thigh Support Pad Short, Left (5"x5")	\$144.00
250KF24	Lateral Thigh Support Pad Long, Right (5"x7")	\$144.00
250KF25	Lateral Thigh Support Pad Long, Left (5"x7")	\$144.00
250KF26	Lateral Thigh Removable Mount - Right	\$320.00
250KF27	Lateral Thigh Removable Mount - Left	\$320.00

Armrest Structure P. 16 \$

250EG35	Reclining, Flip Back (Height Adjustable 8" to 13")	NC
250EG36	Reclining, Flip Back with 2-point post locking arm (Height Adjustable 11" to 13" x .625")	\$346.00
250EG34	Single Post, SEDEO LITE, Removable (N/A with thigh support), (Height Adjustable 10" to 15")	NC

Arm Pads P. 16 \$

SEDEO® Armrest Pad

250AB62	300 mm / 12", Integral Foam, Right	NC
250AE62	300 mm / 12", Integral Foam, Left	NC
250AB63	400 mm / 16", Integral Foam, Right	NC
250AE63	400 mm / 16", Integral Foam, Left	NC
250AB61	SEDEO LITE Pad, Right	NC
250AE61	SEDEO LITE Pad, Left	NC

Comfort Arm Pad (Foam)

250AB46	2" Wide Right	\$111.00
250AE46	2" Wide Left	\$111.00

Comfort Arm Pad (Gel)

250AB48	2" Wide Right	\$210.00
250AE48	2" Wide Left	\$210.00

Ottobock® Arm Troughs

250AE06	Ottobock Arm Trough, Med (12"), Left	\$137.00
250AB06	Ottobock Arm Trough, Med (12"), Right	\$137.00
250AE07	Ottobock Arm Trough, Large (14"), Left	\$137.00
250AB07	Ottobock Arm Trough, Large (14"), Right	\$137.00
250OL01	Ottobock Hand Pad, Med (7"), Left	\$69.00
250OR01	Ottobock Hand Pad, Med (7"), Right	\$69.00
250OL02	Ottobock Hand Pad, Large (9"), Left	\$69.00
250OR02	Ottobock Hand Pad, Large (9"), Right	\$69.00
250AE78	Ottobock One-Piece Arm Trough (18"), Left	\$195.00
250AB78	Ottobock One-Piece Arm Trough (18"), Right	\$195.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	DESCRIPTION	RETAIL
Fixed Centermount Legrest		\$
250H183	SEDEO® Fixed Centermount Legrest (3" to 17.5" Leg Length)	\$346.00
250H184	Fixed Centermount Legrest, Short (Includes Calf Pads & Single Footplate) (11.5" to 14.5" Leg Length)	NC
250H185	Fixed Centermount Legrest, Long (Includes Calf Pads & Single Footplate) (13.5" to 18.5" Leg Length)	NC
250H136	Angle Adjustable Centermount Legrest (11.5" to 18.5" Leg Length)	\$596.00
Power Centermount Legrest  P. 16 (Client Cost)		\$
250H175	Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension (Mechanical Articulation) (5" to 17.5" Leg Length)	\$2,450.00
250H109	Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension (Powered Articulation) (3" to 17.5" Leg Length)	\$3,384.00
Power Centermount Legrest Accessories (Client Cost)		\$
SEDEO Centermount Footplates		
250F127	Dual, Height, Width and Angle Adjustable, Standard (14"W x 10"D)	\$280.00
250F128	Dual, Height, Width and Angle Adjustable, Large (16"W x 12"D)	\$280.00
Angle Adjustable Centermount Footplates		
250F45	Single / Dual Footplate	\$280.00
250F44	Single Platform Footplate	\$280.00
SEDEO Centermount Calf Pads		
250DI34	Calf Pads, Pair, Standard Integral Foam	NC
250DI35	Calf Pads, Pair, Padded	\$351.00
Angle Adjustable Centermount Calf Pads		
250DI49	Calf Pads, Pair, Standard Integral Foam	NC
250DI50	Extra Large Single Calf Pad	\$351.00
SEDEO® PRO™ Swing-Away Legrest (Pair) (250H179 is Client Cost)		\$
250H177	Legrests, Swing-Away, 70°, SEDEO® PRO™ (14" to 21" Leg Length)	NC
250H178	Legrests, Manual Elevating, Swing-Away, SEDEO® PRO™ (Includes Padded Calf Pads) (14".5 to 21.5" Leg Length)	\$596.00
250H179	Legrests, Power Elevating, Swing-Away, SEDEO® PRO™ (Includes Padded Calf Pads) (15.5" to 22.5" Leg Length)	\$3,384.00
250DI76	Omit Legrests	NC
Swing-Away Legrest Footplates		
250F369	Angle Adjustable Flip-Up Footplates	\$280.00
Swing-Away Legrest Calf Supports		
250DI42	Calf Support, Depth Adjustable, Padded (6.5" to 8" Wide, Depending On Chair Width, x 9.5" Tall)	\$504.00
Swing-Away Legrest Options		
250DI44	Heel Loops	NC
250KG4	Knee Pads, Integral Foam (2.5" Wide x 4.5" Tall)	\$278.00
SEDEO® Lower Extremity Swing-Away Support		
250VD01	Lower Extremity Support, Padded, Swing-Away, SEDEO®, Right	\$504.00
250VD02	Lower Extremity Support, Padded, Swing-Away, SEDEO®, Left	\$504.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION		DESCRIPTION	RETAIL
Pelvic Positioning Belts  P. 16			
250A837		Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
250A839		Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"	NC
250A840		Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$101.00
250A845		Bodypoint – Evoflex Pelvic Stabilizer Belt	\$213.00
250A846		Bodypoint – Evoflex Pelvic Stabilizer Belt with Extension	\$267.00
Accessories			\$
Ventilator Option			
250VT12		Vent Tray Fixed (<i>Pulmonetic LTV Models</i>)	\$600.00
O2 Tank Option			
250A718		O2 Tank Holder Mounted to Backrest (<i>D and E Styles</i>)	\$385.00
Wheelchair Tray Table			
2507A09		Tray, Swing-Away	\$375.00
Backpack (Client Cost)			
250A1		Black Utility	\$110.00
Medial Thigh Support			
2508A80		Medial Thigh Support	\$267.00
2508A81		Pommel Swing-Away Hardware	\$384.00
Misc. Accessories (Client Cost)			
2508A59		Transfer Bar (<i>Pair</i>)	\$486.00
2508A52		Ram Mount Phone	\$76.00
2508A53		Ram Mount Tablet	\$303.00
2508A54		Ram Mount Cup Holder	\$76.00
2508A58		Rear View Mirror	\$76.00
250A611		Crutch Holder	\$171.00
2508A73		Bag for Armrest	\$76.00

JOYSTICK / ELECTRONICS

OPTION		DESCRIPTION	RETAIL
R-net Joystick <i>(Remote)</i>  P. 17			\$
250J121		Standard LED Joystick <i>(w/o Screen Display)</i>	NC
250J116		Color Joystick <i>(2.3" Color LCD Screen Display, CJS1M1)</i>	NC
250J118		High Resolution Color Joystick,(CJS1M2+BT) <i>3.5" Color LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control (Client Cost)</i>	\$1,245.00
250J119		Omit Remote Control <i>(Required with arm mounted specialty controls)</i>	NC
Joystick Mounting			\$
250JM34		Fixed on Armrest	NC
250JM35		Parallel Swing-Away	\$440.00
250JM36		Link-It Mid Line Mount on Armrest <i>(Swing-Away)</i>	\$440.00
Joystick Mount (Location)			\$
250HM16		Right Hand Armrest	NC
250HM17		Left Hand Armrest	NC
Display Language Selection R-net			\$
250L1		English	NC
250L3		French	NC
250L7		Spanish	NC
Joystick Handle (Optional)  P. 17			\$
250JH13		Ball Handle	\$108.00
250JH15		T-Handle	\$108.00
250JH24		Bodypoint - Rubber Dome Joystick Handle	\$213.00
250JH25		Bodypoint - U-Shaped Joystick Handle <i>(4 inch)</i>	\$213.00
250JH26		Bodypoint - Grooved Mushroom Joystick Handle	\$213.00
Electronic Programmers  P. 18			\$
250PR10		R-net PC Programmer w/ Dongle	\$696.00
R-net Power Modules  P. 18			\$
250ER11		Power Module, R-net, 120A, Non-Expandable <i>(Only available with Group 4)</i>	NC
250ER10		Power Module, R-net, 120A, Expandable <i>(Only available with Group 4)</i> <i>Required for specialty controls or 3+ actuators</i>	NC
250EX09		Expandable 120A Controller - R-net	NC
250EX10		Harness for 120A Expandable Controller	NC
250ER18		Power Module, R-net, 90A, Non-Expandable	NC
250ER17		Power Module, R-net, 90A, Expandable <i>Required for specialty controls or 3+ actuators</i>	\$100.00
250EX15		Expandable 90A Controller - R-net	NC
250EX16		Harness for 90A Expandable Controller	NC

JOYSTICK / ELECTRONICS

OPTION		DESCRIPTION	RETAIL
--------	--	-------------	--------

Seating Control P. 18

\$

Powered Seating Control

250E040		Omit Thru-Drive <i>(Requires Switch Control)</i>	NC
250E037		Thru-Drive Control for 1 Actuator	\$1,422.00
250E038		Thru-Drive Control for 2+ Actuators	\$2,860.00

Switch Control *(Available with Omit Thru-Drive only) (1 actuator for Tilt or Elevate only)*

250E033		Dual Toggle for 1 Actuator	NC
250E034		2.5" Single Button Switch for 1 Actuator <i>(Button switches are not mounted)</i>	NC
250E041		Ctrl+5 Switch Box for 2 Actuators or More	NC

Additional Power Seating Control *(Available with Thru-Drive only)*

250QG11		Ctrl+5 Switch Box with Free Assignable Functions	\$1,062.00
---------	--	--	------------

Dual Toggle Mounting Location

250TN1		Toggle Mount, Right	NC
250TN2		Toggle Mount, Left	NC

Ctrl+5 Mounting Location

250EY19		Right Armrest	NC
250EY20		Left Armrest	NC
250EY21		Back Cane, Right	NC
250EY24		Back Cane, Left	NC

Ctrl+5 Switch Type P. 18

250EY22		Push Buttons	NC
250EY23		Toggle Switches	\$351.00

Attendant Controls P. 18 (Prior Approval Required)

\$

2500A05		Dual Control Module for Attendant	\$1,207.00
2500A06		R-net Color LCD (CJSM1), Attendant Control Only	\$1,207.00

Attendant Mounting Location

\$

250DN5		Right Back Cane	NC
250DN6		Left Back Cane	NC

R-net Programming

\$

Assignable Buttons

Assignable Button feature allows for personalization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

250PM12		Include Assignable Buttons	NC
---------	--	----------------------------	----

Power Centermount Footplate Extension Option

250IE03		Include 7.5" Powered Footplate Extension for Transfer Assist <i>(Electronics programmed to allow footplates to be extended with or without legrest being raised.) Only available with Powered Articulation.</i>	\$583.00
---------	--	---	----------

built4me

Q700 M with SEDEO® LITE™

\$

Base Price

Total of additional equipment

Packaging and transport

Total Price

Date

Signature

Month

Day

Year



Customer Service: (800) 263-3390
Fax: (800) 561-5834
E-mail: cscanada@sunmed.com
www.sunrisemedical.ca
Specifications are subject to change without notice.

© 2024 Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2
Concord, ON
L4K 4P8

MK-129765-AADL CAN Rev. 11