



October 2026

**Q700M** with SEDEO LITE  
**AADL - W581**

EIPW25

|             |  |
|-------------|--|
| Account:    |  |
| Date:       |  |
| PO#:        |  |
| Buyer Name: |  |
| Marked For: |  |


|                     |  |
|---------------------|--|
| Ship To:            |  |
| Name:               |  |
| Address:            |  |
| City, Prov.,<br>PC: |  |
| Phone:              |  |

Sunrise Medical Inc.  
1000 Creditstone Road, Unit 2  
Concord, ON L4K 4P8  
1-800-263-3390



## Client Measurements

|     |                         |      |
|-----|-------------------------|------|
| a)  | Top of head to seat     | Inch |
| b)  | Top of shoulder to seat | Inch |
| c)  | Underarm to seat        | Inch |
| d)  | Elbow to seat           | Inch |
| e)  | Top of leg to seat      | Inch |
| f)  | Shoulder width          | Inch |
| g)  | Chest width             | Inch |
| i)  | Hip width               | Inch |
| j)  | Actual sitting depth    | Inch |
| m)  | Foot length             | Inch |
| (n) | Knee to foot            | Inch |



## User Weight

Notes:

## Chair Specifications

|                                              |                                    |                         |                                      |
|----------------------------------------------|------------------------------------|-------------------------|--------------------------------------|
|                                              |                                    | <b>Seat Type</b>        | SEDEO LITE                           |
| <b>Weights Starting at: (With Batteries)</b> | NPC: 294 lbs.                      | <b>Seat Size Range</b>  | Width 15" to 20"<br>Depth 16" to 22" |
|                                              | SPD 318 lbs.                       |                         |                                      |
| <b>Power Seating</b>                         | 50° of CG Tilt / 12" seat elevator | <b>STF Height Range</b> | 16.5" to 19"                         |

## Mid-Wheel Drive Base Specifications

|                               |                                        |                         |                                              |
|-------------------------------|----------------------------------------|-------------------------|----------------------------------------------|
| <b>Drive Wheel Position</b>   | Mid-wheel Drive                        | <b>Drive Wheels</b>     | 14" drive wheels w/flat free inserts (black) |
| <b>User Weight Capacity</b>   | 300 lbs. / 350 lbs. HD                 | <b>Tires/Casters</b>    | 6" Front / 6" Rear (Black)                   |
| <b>Minimum Turning Radius</b> | 19.8" (excludes footplates)            |                         |                                              |
| <b>Base Size</b>              | Width 24.5", Length 36.5"              | <b>Charger</b>          | 10A Dual Mode (Compact)                      |
| <b>Drive Control</b>          | R-net 120A / R-net 90A                 | <b>Battery Size</b>     | Group 24 / Group 34                          |
| <b>Driving Range</b>          | Up to 18 miles (w/ Group 24 Batteries) | <b>Color Accents</b>    | 8 Shroud Colors and 7 Color Accents          |
| <b>Drive Motors</b>           | 4-pole 6mph, 8mph                      | <b>Ground Clearance</b> | 3"                                           |

| BASE                                                        |  |                                                                                                                                                                                                                                                                                                                                                                                                                     |             |
|-------------------------------------------------------------|--|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------|
| OPTION                                                      |  | Description                                                                                                                                                                                                                                                                                                                                                                                                         | Retail      |
| Standard Duty Mid Wheel Drive Base = End Item Number EIPW25 |  |                                                                                                                                                                                                                                                                                                                                                                                                                     |             |
| <input type="checkbox"/> 250PC267                           |  | Quickie Q700 M NPC No Power Option for SEDEO LITE Seat<br>Included in Base Pricing: Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5". <b>(AADL Contract Price)</b>        | \$8,993.00  |
| <input type="checkbox"/> 250PC268                           |  | Quickie Q700 M SPC Single Power Option for SEDEO LITE Seat<br>Included in Base Pricing: Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5". <b>(AADL Contract Price)</b>    | \$8,993.00  |
| Heavy Duty Mid Wheel Drive Base = End Item Number EIPW25    |  |                                                                                                                                                                                                                                                                                                                                                                                                                     |             |
| <input type="checkbox"/> 250PC269                           |  | Quickie Q700 M NPC HD No Power Option for SEDEO LITE Seat<br>Included in Base Pricing: Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 5mph, Weight Capacity = 350lbs., STFH = 16.5". <b>(AADL Contract Price)</b>     | \$9,969.00  |
| <input type="checkbox"/> 250PC270                           |  | Quickie Q700 M SPC HD Single Power Option for SEDEO LITE Seat<br>Included in Base Pricing: Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 5mph, Weight Capacity = 350lbs., STFH = 16.5". <b>(AADL Contract Price)</b> | \$9,969.00  |
| Standard Duty Mid Wheel Drive Base = End Item Number EIPW25 |  |                                                                                                                                                                                                                                                                                                                                                                                                                     |             |
| <input type="checkbox"/> 250PC271                           |  | Quickie Q700 M NPD No Power Option for SEDEO LITE Seat<br>Included in Base Pricing: Q700 M Base, R-net 120amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5". <b>(Client Cost)</b>               | \$10,781.00 |
| <input type="checkbox"/> 250PC272                           |  | Quickie Q700 M SPD Single Power Option for SEDEO LITE Seat<br>Included in Base Pricing: Q700 M Base, R-net 120amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5". <b>(Client Cost)</b>           | \$10,781.00 |

| BASE   |  |             |        |
|--------|--|-------------|--------|
| OPTION |  | Description | Retail |

#### User Weight

|                          |         |                |    |
|--------------------------|---------|----------------|----|
| <input type="checkbox"/> | 250WG23 | Up to 300 lbs. | NC |
| <input type="checkbox"/> | 250WG15 | 301 - 350 lbs. | NC |

#### Motor Package

|                          |          |                                              |            |
|--------------------------|----------|----------------------------------------------|------------|
| <input type="checkbox"/> | 250PF98  | 4-pole 8 kph / 5 mph                         | NC         |
| <input type="checkbox"/> | 250PF99  | 4-pole 10 kph / 6 mph                        | NC         |
| <input type="checkbox"/> | 250PF101 | 4-pole 12.5 kph / 8 mph <b>(Client Cost)</b> | \$1,268.00 |

#### Tracking

|                          |         |                                                                   |            |
|--------------------------|---------|-------------------------------------------------------------------|------------|
| <input type="checkbox"/> | 250MI01 | SureTrac Module & Motor Encoders <i>(Improves Drive Tracking)</i> | \$1,696.00 |
|--------------------------|---------|-------------------------------------------------------------------|------------|

#### Drive Wheels

|                          |         |                                               |    |
|--------------------------|---------|-----------------------------------------------|----|
| <input type="checkbox"/> | 250DW58 | Drive Wheel, 14", Dual Density, Knobby, Black | NC |
|--------------------------|---------|-----------------------------------------------|----|

#### Caster Wheels

|                          |         |                                |    |
|--------------------------|---------|--------------------------------|----|
| <input type="checkbox"/> | 250CW55 | Caster Wheels 6", Solid, Black | NC |
|--------------------------|---------|--------------------------------|----|

#### Shroud Color

|                          |         |                 |    |
|--------------------------|---------|-----------------|----|
| <input type="checkbox"/> | 250C119 | Matte Black     | NC |
| <input type="checkbox"/> | 250C165 | Candy Apple Red | NC |
| <input type="checkbox"/> | 250C166 | Candy Blue      | NC |
| <input type="checkbox"/> | 250C167 | Candy Purple    | NC |
| <input type="checkbox"/> | 250C176 | Green Apple     | NC |
| <input type="checkbox"/> | 250C178 | Pearl Pink      | NC |
| <input type="checkbox"/> | 250C182 | Silver          | NC |
| <input type="checkbox"/> | 250C125 | Pearl White     | NC |

#### Color Accents

|                          |         |             |    |
|--------------------------|---------|-------------|----|
| <input type="checkbox"/> | 250EZ20 | Apple Green | NC |
| <input type="checkbox"/> | 250EZ22 | Black       | NC |
| <input type="checkbox"/> | 250EZ21 | Blue        | NC |
| <input type="checkbox"/> | 250EZ23 | Orange      | NC |
| <input type="checkbox"/> | 250EZ24 | Red         | NC |
| <input type="checkbox"/> | 250EZ19 | White       | NC |
| <input type="checkbox"/> | 250EZ25 | Silver      | NC |

#### BASE

| OPTION |  | Description | Retail |
|--------|--|-------------|--------|
|--------|--|-------------|--------|

#### ISO Transit Tie Downs

|                          |         |                     |    |
|--------------------------|---------|---------------------|----|
| <input type="checkbox"/> | 250TO50 | ISO Transit Package | NC |
|--------------------------|---------|---------------------|----|

#### LED Lights **(Client Cost)**

|                          |         |                   |          |
|--------------------------|---------|-------------------|----------|
| <input type="checkbox"/> | 250LI08 | LED Fender Lights | \$646.00 |
| <input type="checkbox"/> | 250LI02 | Omit LED Lights   | NC       |

#### Seat Heights

Note: Available seat-to-Floor heights are dependant on 2 options: **1 . Battery size and 2 . Powered Seating.**

Use the matrix below to determine available ranges for your desired options. Measured from front of seat pan to floor with 0° frame angle.

|                          |          | Seat Height | Powered Seating<br>Without Elevate | No Powered<br>Seating |    |
|--------------------------|----------|-------------|------------------------------------|-----------------------|----|
| <input type="checkbox"/> | 250HT103 | 16.5"       | Group 34                           | Group 34              | NC |
| <input type="checkbox"/> | 250HT55  | 17"         | Group 34                           | GP24 or GP34          | NC |
| <input type="checkbox"/> | 250HT104 | 17.5"       | Group 34                           | GP24 or GP34          | NC |
| <input type="checkbox"/> | 250HT56  | 18"         | GP24 or GP34                       | GP24 or GP34          | NC |
| <input type="checkbox"/> | 250HT105 | 18.5"       | GP24 or GP34                       | GP24 or GP34          | NC |
| <input type="checkbox"/> | 250HT57  | 19"         | GP24 or GP34                       | GP24 or GP34          | NC |

#### Batteries & Chargers

##### Group 24 Battery Options *(Recommended for max daily range)*

|                          |         |                                          |          |
|--------------------------|---------|------------------------------------------|----------|
| <input type="checkbox"/> | 250AG25 | GP24 Flat Top Pair w/Strap MK Gel        | \$955.00 |
| <input type="checkbox"/> | 250AG26 | GP24 Flat Top Pair w/Strap Brightway AGM | \$915.00 |

##### Group 34 Battery Options *(Required for lower seat-to-floor heights)*

|                          |         |                                           |          |
|--------------------------|---------|-------------------------------------------|----------|
| <input type="checkbox"/> | 250AG27 | GP34 Flat Top Pair w/Straps MK Gel        | \$955.00 |
| <input type="checkbox"/> | 250AG29 | GP34 Flat Top Pair w/Straps Brightway AGM | \$863.00 |

#### Omit Batteries

|                          |         |                                                |    |
|--------------------------|---------|------------------------------------------------|----|
| <input type="checkbox"/> | 250AG24 | Battery Harnesses Only <i>(Omit Batteries)</i> | NC |
|--------------------------|---------|------------------------------------------------|----|

#### Chargers

|                          |         |                                |    |
|--------------------------|---------|--------------------------------|----|
| <input type="checkbox"/> | 250BJ28 | 10A Dual Mode <i>(Compact)</i> | NC |
|--------------------------|---------|--------------------------------|----|

## Powered Seating

| OPTION | Description | Retail |
|--------|-------------|--------|
|--------|-------------|--------|

### Power Seating Options

#### Seat Elevator (Client Cost)

|                          |         |                                                                                                                 |            |
|--------------------------|---------|-----------------------------------------------------------------------------------------------------------------|------------|
| <input type="checkbox"/> | 250PI08 | C-Me Seat Elevator, Powered, 300 mm / 12" (Not available with HD model or Group 4)<br>C-Me speed up to 3.0 mph. | \$3,477.00 |
|--------------------------|---------|-----------------------------------------------------------------------------------------------------------------|------------|

#### Seat Tilt (AADL W992 for Power Tilt)

|                          |         |                       |            |
|--------------------------|---------|-----------------------|------------|
| <input type="checkbox"/> | 250PA35 | Posterior 50° Tilt    | \$8,663.00 |
| <input type="checkbox"/> | 250PA36 | Posterior 50° HD Tilt | \$8,923.00 |

### Pre Set Seat Frame Angle (No Power Seating Options, Increases Seat to Floor Height)

|                          |         |                      |    |
|--------------------------|---------|----------------------|----|
| <input type="checkbox"/> | 250FA10 | 0° (Adds 0" to STFH) | NC |
| <input type="checkbox"/> | 250FA3  | 3° (Adds 1" to STFH) | NC |
| <input type="checkbox"/> | 250FA4  | 6° (Adds 2" to STFH) | NC |
| <input type="checkbox"/> | 250FA5  | 9° (Adds 3" to STFH) | NC |

### SEDEO Lite Backrest Frame (Angle Adjustable 86° - 106°)

#### Back Frame Width

|                          |         |                                           |    |
|--------------------------|---------|-------------------------------------------|----|
| <input type="checkbox"/> | 250BF21 | Width Adjustable, 385 - 485 mm / 15 - 19" | NC |
| <input type="checkbox"/> | 250BF22 | Width Adjustable, 460 - 510 mm / 18 - 20" | NC |

#### Back Cane Height

|                          |         |                                                                                                     |    |
|--------------------------|---------|-----------------------------------------------------------------------------------------------------|----|
| <input type="checkbox"/> | 250B263 | Standard (Lite Sling Height 18" - 20")                                                              | NC |
| <input type="checkbox"/> | 250B264 | Tall (Lite Sling Height 21"-23")                                                                    | NC |
| <input type="checkbox"/> | 300PH30 | Height Adjustable, Short w/ Trailing Push Handles (Push Handle Height 20" - 22", N/A w/ Sling Back) | NC |
| <input type="checkbox"/> | 300PH31 | Height Adjustable, Tall w/ Trailing Push Handles (Push Handle Height 23" - 26", N/A w/ Sling Back)  | NC |

### Backrest Frame Option

|                          |         |                          |    |
|--------------------------|---------|--------------------------|----|
| <input type="checkbox"/> | 250BF23 | Flip Down, Quick Release | NC |
|--------------------------|---------|--------------------------|----|

### Back Cushion Type

|                          |          |                                                                               |    |
|--------------------------|----------|-------------------------------------------------------------------------------|----|
| <input type="checkbox"/> | 250BT148 | LITE Sling Back, with Tension Adjustable Straps (Choose Seating in Section A) | NC |
| <input type="checkbox"/> | 250BT135 | JAY J3 (Choose Seating in Section C)                                          | NC |
| <input type="checkbox"/> | 250BT133 | OMIT (For Custom Mounting. No Headrest Mount)                                 | NC |

### A - SEDEO Lite Sling Back Dimensions

#### Installed Back Height - (From Pan to Top of Back)

|                          |         |     |    |
|--------------------------|---------|-----|----|
| <input type="checkbox"/> | 250JR71 | 18" | NC |
| <input type="checkbox"/> | 250JR72 | 19" | NC |
| <input type="checkbox"/> | 250JR73 | 20" | NC |
| <input type="checkbox"/> | 250JR74 | 21" | NC |
| <input type="checkbox"/> | 250JR75 | 22" | NC |
| <input type="checkbox"/> | 250JR76 | 23" | NC |

| Powered Seating                                                                                                                                                                                                                          |         |                                                             |            |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------|-------------------------------------------------------------|------------|
| OPTION                                                                                                                                                                                                                                   |         | Description                                                 | Retail     |
| <b>B - JAY J3 Back</b>                                                                                                                                                                                                                   |         |                                                             |            |
| <b>B - JAY J3 Back - Contour</b>                                                                                                                                                                                                         |         |                                                             |            |
| <input type="checkbox"/>                                                                                                                                                                                                                 | 250JP01 | J3 PA Back for 16"-20" Wide                                 | \$912.00   |
| <input type="checkbox"/>                                                                                                                                                                                                                 | 250JP11 | J3 PL Back for 16"-20" Wide                                 | \$986.00   |
| <input type="checkbox"/>                                                                                                                                                                                                                 | 250JP13 | J3 PD Back for 16"-20" Wide (18" and 21" shell height only) | \$1,068.00 |
| <b>B - JAY J3 Back - Mounting Hardware</b>                                                                                                                                                                                               |         |                                                             |            |
| <input type="checkbox"/>                                                                                                                                                                                                                 | 250JP20 | TH                                                          | NC         |
| <b>B - JAY J3 Back Width</b>                                                                                                                                                                                                             |         |                                                             |            |
| <i>Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.</i>                                                                                                                              |         |                                                             |            |
| <input type="checkbox"/>                                                                                                                                                                                                                 | 250JQ02 | 16"                                                         | NC         |
| <input type="checkbox"/>                                                                                                                                                                                                                 | 250JQ03 | 18"                                                         | NC         |
| <input type="checkbox"/>                                                                                                                                                                                                                 | 250JQ04 | 20"                                                         | NC         |
| <b>B - JAY J3 Back Shell Options and Height Ranges.</b>                                                                                                                                                                                  |         |                                                             |            |
| <i>JAY J3 shells are mounted with a minimim 4" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field with-in the ranges shown, based on back frame type.</i> |         |                                                             |            |
| <b>JAY J3 Shell Height</b>                                                                                                                                                                                                               |         |                                                             |            |
| <input type="checkbox"/>                                                                                                                                                                                                                 | 250JR14 | 16.5" tall shell (Upper Thoracic Short)                     | NC         |
| <input type="checkbox"/>                                                                                                                                                                                                                 | 250JR15 | 18" tall shell (Upper Thoracic Med)                         | NC         |
| <input type="checkbox"/>                                                                                                                                                                                                                 | 250JR16 | 19.5" tall shell (Upper Thoracic Tall)                      | NC         |
| <input type="checkbox"/>                                                                                                                                                                                                                 | 250JR17 | 21" tall shell (Shoulder Height)                            | NC         |
| <input type="checkbox"/>                                                                                                                                                                                                                 | 250JR27 | 22.5" tall shell (Shoulder Height)                          | NC         |
| <input type="checkbox"/>                                                                                                                                                                                                                 | 250JR28 | 24" tall shell (Shoulder Height)                            | NC         |

#### JAY J3 Installed Back Height (From Pan to Top of Back)

|                          |         |     |    |
|--------------------------|---------|-----|----|
| <input type="checkbox"/> | 250JR29 | 21" | NC |
| <input type="checkbox"/> | 250JR19 | 22" | NC |
| <input type="checkbox"/> | 250JR20 | 23" | NC |
| <input type="checkbox"/> | 250JR21 | 24" | NC |
| <input type="checkbox"/> | 250JR22 | 25" | NC |
| <input type="checkbox"/> | 250JR23 | 26" | NC |
| <input type="checkbox"/> | 250JR24 | 27" | NC |
| <input type="checkbox"/> | 250JR25 | 28" | NC |
| <input type="checkbox"/> | 250JR30 | 29" | NC |
| <input type="checkbox"/> | 250JR31 | 30" | NC |

| Powered Seating                                                                              |         |                                                                |          |
|----------------------------------------------------------------------------------------------|---------|----------------------------------------------------------------|----------|
| <b>B - JAY Anterior Thoracic Support</b>                                                     |         |                                                                |          |
| <i>Classic-Structure Pull Down Harness (4 point strap, Only available with JAY backrest)</i> |         |                                                                |          |
| <input type="checkbox"/>                                                                     | 250K10  | Medium (ATS8302)                                               | \$322.00 |
| <input type="checkbox"/>                                                                     | 250K11  | Large (ATS8303)                                                | \$322.00 |
| <b>B - JAY Lateral Thoracic Swing-Away Supports (Requires Pad, Upholstery and Mount)</b>     |         |                                                                |          |
| <b>JAY Lateral Thoracic Support Pads</b>                                                     |         |                                                                |          |
| <input type="checkbox"/>                                                                     | 250JW35 | Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right | \$75.00  |
| <input type="checkbox"/>                                                                     | 250JW38 | Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left  | \$75.00  |
| <input type="checkbox"/>                                                                     | 250JW36 | Lateral Thoracic (L x H) 100 mm x 152 mm / 4" x 6" Pad - Right | \$75.00  |
| <input type="checkbox"/>                                                                     | 250JW39 | Lateral Thoracic (L x H) 100 mm x 152 mm / 4" x 6" Pad - Left  | \$75.00  |
| <input type="checkbox"/>                                                                     | 250JW37 | Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right | \$75.00  |
| <input type="checkbox"/>                                                                     | 250JW40 | Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left  | \$75.00  |
| <b>JAY Lateral Thoracic Support Pad Upholstery</b>                                           |         |                                                                |          |
| <input type="checkbox"/>                                                                     | 250NU01 | Premium Leatherette                                            | NC       |
| <input type="checkbox"/>                                                                     | 250NU02 | Stretch                                                        | NC       |
| <b>JAY Lateral Thoracic Swing-Away Mount</b>                                                 |         |                                                                |          |
| <input type="checkbox"/>                                                                     | 250JW41 | Lateral Thoracic Swing Away Bracket - Right                    | \$334.00 |
| <input type="checkbox"/>                                                                     | 250JW42 | Lateral Thoracic Swing Away Bracket - Left                     | \$334.00 |

| OPTION | Description | Page 7 | Retail |
|--------|-------------|--------|--------|
|--------|-------------|--------|--------|

### Headrests (Requires Pad, Upholstery and Mount)

#### Headrest Pad

|                          |         |                                                                                                                                           |          |
|--------------------------|---------|-------------------------------------------------------------------------------------------------------------------------------------------|----------|
| <input type="checkbox"/> | 250HE36 | SEDEO Pad Small With Cover (L x H) 240 x 150 mm / 10" x 6" (Only Available with SEDEO Multi-Position Headrest Mount and SEDEO Upholstery) | \$227.00 |
| <input type="checkbox"/> | 250HE38 | Whitmyer Plush Headrest Pad 203 mm / 8"                                                                                                   | \$377.00 |
| <input type="checkbox"/> | 250HE39 | Whitmyer Plush Headrest Pad 254 mm / 10"                                                                                                  | \$377.00 |
| <input type="checkbox"/> | 250HE44 | Whitmyer Plush Headrest Pad 356 mm / 14"                                                                                                  | \$377.00 |
| <input type="checkbox"/> | 250HE46 | Whitmyer Adjust-A-Plush Headrest Pad Medium (Adjustable from 7" to 16" Wide x 4" Tall)                                                    | \$445.00 |
| <input type="checkbox"/> | 250HE48 | OMIT Headrest Pad                                                                                                                         | NC       |

#### Headrest Pad Upholstery

|                          |         |                  |    |
|--------------------------|---------|------------------|----|
| <input type="checkbox"/> | 250NU10 | SEDEO Upholstery | NC |
| <input type="checkbox"/> | 250NU11 | Whitmyer Lycra   | NC |

#### Headrest Mount

|                          |         |                                                                                        |          |
|--------------------------|---------|----------------------------------------------------------------------------------------|----------|
| <input type="checkbox"/> | 250MH32 | SEDEO Multi-Position Headrest Mount with Quick Release (Only available with SEDEO Pad) | \$557.00 |
| <input type="checkbox"/> | 250MH33 | Whitmyer Linx Headrest Mount with Quick Release                                        | \$247.00 |
| <input type="checkbox"/> | 250MH34 | Whitmyer Linx-2 Headrest Mount with Quick Release                                      | \$298.00 |
| <input type="checkbox"/> | 250MH35 | Whitmyer AXYS Headrest Mount with Quick Release                                        | \$247.00 |
| <input type="checkbox"/> | 250MH17 | OMIT Headrest Mount                                                                    | NC       |

### Powered Seating

#### Seat Cushion Type

|                          |         |                     |    |
|--------------------------|---------|---------------------|----|
| <input type="checkbox"/> | 250ST80 | JAY (See Section C) | NC |
| <input type="checkbox"/> | 250ST75 | Omit                | NC |

#### Seat Size

##### Seat Pan Width (Odd sizes not available on SEDEO cushions)

|                          |        |     |    |
|--------------------------|--------|-----|----|
| <input type="checkbox"/> | 250R78 | 15" | NC |
| <input type="checkbox"/> | 250R53 | 16" | NC |
| <input type="checkbox"/> | 250R54 | 17" | NC |
| <input type="checkbox"/> | 250R55 | 18" | NC |
| <input type="checkbox"/> | 250R56 | 19" | NC |
| <input type="checkbox"/> | 250R57 | 20" | NC |

##### Seat Pan Depth (Front of Seat Pan to front of back cushion. Effective Depth)

|                          |         |     |    |
|--------------------------|---------|-----|----|
| <input type="checkbox"/> | 250RF51 | 16" | NC |
| <input type="checkbox"/> | 250RF45 | 17" | NC |
| <input type="checkbox"/> | 250RF46 | 18" | NC |
| <input type="checkbox"/> | 250RF47 | 19" | NC |
| <input type="checkbox"/> | 250RF48 | 20" | NC |
| <input type="checkbox"/> | 250RF49 | 21" | NC |
| <input type="checkbox"/> | 250RF50 | 22" | NC |

#### F - JAY Seat Cushion

|                          |          |                                                        |            |
|--------------------------|----------|--------------------------------------------------------|------------|
| <input type="checkbox"/> | 250CU59  | JAY Ion (Sizes 16" - 21" Wide)                         | \$657.00   |
| <input type="checkbox"/> | 250CU61  | JAY Union (Sizes 16" - 21" Wide)                       | \$764.00   |
| <input type="checkbox"/> | 250CU64  | JAY Fusion (Sizes 16" - 21" Wide)                      | \$1,045.00 |
| <input type="checkbox"/> | 250CU104 | JAY Fusion with Cryo Technology (Sizes 16" - 21" Wide) | \$1,152.00 |

## Positioning Supports & Accessories

| OPTION | Description | Retail |
|--------|-------------|--------|
|--------|-------------|--------|

### Knee and Thigh Swing Away Support Pads (Seating Benefits)

#### Lateral Knee Support Pads and Hardware (Must pick pad and mount)

|                          |         |                                            |          |
|--------------------------|---------|--------------------------------------------|----------|
| <input type="checkbox"/> | 250XY05 | Lateral Knee Support Pad, Right (4" round) | \$247.00 |
| <input type="checkbox"/> | 250XY06 | Lateral Knee Support Pad, Left (4" round)  | \$247.00 |
| <input type="checkbox"/> | 250XX08 | Lateral Swing Away Bracket Knee - Right    | \$325.00 |
| <input type="checkbox"/> | 250XX09 | Lateral Swing Away Bracket Knee - Left     | \$325.00 |

#### Lateral Thigh Support Pads and Hardware (Must pick pad and mount)

|                          |         |                                                      |          |
|--------------------------|---------|------------------------------------------------------|----------|
| <input type="checkbox"/> | 250KF22 | Lateral Thigh Support Pad Short, Right (5" H x 5" L) | \$146.00 |
| <input type="checkbox"/> | 250KF23 | Lateral Thigh Support Pad Short, Left (5" H x 5" L)  | \$146.00 |
| <input type="checkbox"/> | 250KF24 | Lateral Thigh Support Pad Long, Right (5" H x 7" L)  | \$146.00 |
| <input type="checkbox"/> | 250KF25 | Lateral Thigh Support Pad Long, Left (5" H x 7" L)   | \$146.00 |
| <input type="checkbox"/> | 250KF26 | Lateral Thigh Removable Mount - Right                | \$325.00 |
| <input type="checkbox"/> | 250KF27 | Lateral Thigh Removable Mount - Left                 | \$325.00 |

### Armrest Structure

|                          |         |                                                                                           |          |
|--------------------------|---------|-------------------------------------------------------------------------------------------|----------|
| <input type="checkbox"/> | 250EG35 | Reclining, Flip Back, Short (Height Adjustable 8" to 13")                                 | NC       |
| <input type="checkbox"/> | 250EG36 | Reclining, Flip Back with 2-point post locking arm (Height Adjustable 11" to 13" x .625") | \$351.00 |
| <input type="checkbox"/> | 250EG34 | Single Post, SEDEO LITE, Removable (N/A w/ thigh support)(Height Adjustable 10" to 15")   | NC       |

### Arm Pads

#### SEDEO Armrest Pad

|                          |         |                                    |          |
|--------------------------|---------|------------------------------------|----------|
| <input type="checkbox"/> | 250AB62 | 300 mm / 12", Integral Foam, Right | \$113.00 |
| <input type="checkbox"/> | 250AE62 | 300 mm / 12", Integral Foam, Left  | \$113.00 |
| <input type="checkbox"/> | 250AB63 | 400 mm / 16", Integral Foam, Right | \$113.00 |
| <input type="checkbox"/> | 250AE63 | 400 mm / 16", Integral Foam, Left  | \$113.00 |
| <input type="checkbox"/> | 250AB61 | SEDEO Lite Pad, Right              | NC       |
| <input type="checkbox"/> | 250AE61 | SEDEO Lite Pad, Left               | NC       |

#### Comfort Arm Pad (Foam)

|                          |         |               |          |
|--------------------------|---------|---------------|----------|
| <input type="checkbox"/> | 250AB46 | 2" Wide Right | \$113.00 |
| <input type="checkbox"/> | 250AE46 | 2" Wide Left  | \$113.00 |

#### Comfort Arm Pad (Gel)

|                          |         |               |          |
|--------------------------|---------|---------------|----------|
| <input type="checkbox"/> | 250AB48 | 2" Wide Right | \$213.00 |
| <input type="checkbox"/> | 250AE48 | 2" Wide Left  | \$213.00 |

### Ottobock Arm Troughs

|                          |         |                                           |          |
|--------------------------|---------|-------------------------------------------|----------|
| <input type="checkbox"/> | 250AE06 | Ottobock Arm Trough, Med (12"), Left      | \$139.00 |
| <input type="checkbox"/> | 250AB06 | Ottobock Arm Trough, Med (12"), Right     | \$139.00 |
| <input type="checkbox"/> | 250AE07 | Ottobock Arm Trough, Large (14"), Left    | \$139.00 |
| <input type="checkbox"/> | 250AB07 | Ottobock Arm Trough, Large (14"), Right   | \$139.00 |
| <input type="checkbox"/> | 250OL01 | Ottobock Hand Pad, Med (7"), Left         | \$70.00  |
| <input type="checkbox"/> | 250OR01 | Ottobock Hand Pad, Med (7"), Right        | \$70.00  |
| <input type="checkbox"/> | 250OL02 | Ottobock Hand Pad, Large (9"), Left       | \$70.00  |
| <input type="checkbox"/> | 250OR02 | Ottobock Hand Pad, Large (9"), Right      | \$70.00  |
| <input type="checkbox"/> | 250AE78 | Ottobock 1 Piece Arm Trough, (18"), Left  | \$198.00 |
| <input type="checkbox"/> | 250AB78 | Ottobock 1 Piece Arm Trough, (18"), Right | \$198.00 |

### Elbow Support

|                          |         |                                                        |          |
|--------------------------|---------|--------------------------------------------------------|----------|
| <input type="checkbox"/> | 2503A02 | Elbow Support Pad, RH 3" x 10" (Includes side support) | \$228.00 |
| <input type="checkbox"/> | 2503A03 | Elbow Support Pad, LH 3" x 10" (Includes side support) | \$228.00 |
| <input type="checkbox"/> | 2503A06 | Elbow Support Pad, RH 3" x 3"                          | \$167.00 |
| <input type="checkbox"/> | 2503A07 | Elbow Support Pad, LH 3" x 3"                          | \$167.00 |
| <input type="checkbox"/> | 2503A04 | Swing Away Elbow Support Hardware RH                   | \$279.00 |
| <input type="checkbox"/> | 2503A05 | Swing Away Elbow Support Hardware LH                   | \$279.00 |

| Positioning Supports & Accessories                                  |         |                                                                                                                     |            |
|---------------------------------------------------------------------|---------|---------------------------------------------------------------------------------------------------------------------|------------|
| OPTION                                                              |         | Description                                                                                                         | Retail     |
| <b>Fixed Centermount Legrest</b>                                    |         |                                                                                                                     |            |
| <input type="checkbox"/>                                            | 250H184 | Fixed Centermount, Short ( <i>Includes Calf Pads and Single Footplate</i> )(11.5" to 14.5" Leg Length)              | NC         |
| <input type="checkbox"/>                                            | 250H185 | Fixed Centermount, Long ( <i>Includes Calf Pads and Single Footplate</i> )(13.5" to 18.5" Leg Length)               | NC         |
| <input type="checkbox"/>                                            | 250H183 | Sedeo Fixed Centermount (3" to 17.5" Leg Length)                                                                    | \$351.00   |
| <input type="checkbox"/>                                            | 250H136 | Angle Adjustable Centermount ( <i>Includes Calf Pads and Single Footplate</i> )(11.5" to 18.5" Leg Length)          | \$605.00   |
| <b>Power Centermount Legrest (Client Cost)</b>                      |         |                                                                                                                     |            |
| <input type="checkbox"/>                                            | 250H175 | Sedeo Power Centermount, 90° - 0° Angle, 7.5" Extension ( <i>Mechanical Articulation</i> ) (5" to 17.5" Leg Length) | \$2,487.00 |
| <input type="checkbox"/>                                            | 250H109 | Sedeo Power Centermount, 90° - 0° Angle, 7.5" Extension ( <i>Powered Articulation</i> ) (3" to 17.5" Leg Length)    | \$3,435.00 |
| <b>Centermount Legrest Accessories (Client Cost)</b>                |         |                                                                                                                     |            |
| <b>Sedeo Centermount Footplates</b>                                 |         |                                                                                                                     |            |
| <input type="checkbox"/>                                            | 250F127 | Dual, Height, Width and Angle Adjustable, Standard (14"W x 10"D)                                                    | \$284.00   |
| <input type="checkbox"/>                                            | 250F128 | Dual, Height, Width and Angle Adjustable, Large (16"W x 12"D)                                                       | \$284.00   |
| <input type="checkbox"/>                                            | 250DI71 | Flat Footplate ( <i>Vertical Edges removed from PCM Footplates, ideal for mounting footboxes</i> )                  | \$274.00   |
| <b>SEDEO Centermount Footplate Accessories</b>                      |         |                                                                                                                     |            |
| <input type="checkbox"/>                                            | 250DI68 | Extra Large Centermount Footplate                                                                                   | \$315.00   |
| <input type="checkbox"/>                                            | 250DI69 | Footplate Connector Link                                                                                            | \$96.00    |
| <input type="checkbox"/>                                            | 250DI57 | Gel Footpads                                                                                                        | \$401.00   |
| <b>Sedeo Centermount Calf Pads</b>                                  |         |                                                                                                                     |            |
| <input type="checkbox"/>                                            | 250DI34 | Calf Pads, Pair, Standard Integral Foam                                                                             | NC         |
| <input type="checkbox"/>                                            | 250DI35 | Calf Pads, Pair, Padded ( <i>N/A w/ Std footplate</i> )                                                             | \$356.00   |
| <b>SEDEO PRO Swing-Away Legrest (Pair) (250H179 is Client Cost)</b> |         |                                                                                                                     |            |
| <input type="checkbox"/>                                            | 250H177 | Legrests, Swing Away, 70°, SEDEO PRO (14" to 21" Leg Length)                                                        | NC         |
| <input type="checkbox"/>                                            | 250H178 | Legrests, Manual Elevating, Swing Away, SEDEO PRO ( <i>Includes Padded Calf Pads</i> ) (14".5 to 21.5" Leg Length)  | \$605.00   |
| <input type="checkbox"/>                                            | 250H179 | Legrests, Power Elevating, Swing Away, SEDEO PRO ( <i>Includes Padded Calf Pads</i> ) (15.5" to 22.5" Leg Length)   | \$3,435.00 |
| <input type="checkbox"/>                                            | 250H176 | Omit Legrests                                                                                                       | NC         |
| <b>Swing Away Legrest Footplates</b>                                |         |                                                                                                                     |            |
| <input type="checkbox"/>                                            | 250F369 | Angle Adjustable Flip-Up Footplates                                                                                 | \$284.00   |
| <b>Swing Away Legrest Calf Supports</b>                             |         |                                                                                                                     |            |
| <input type="checkbox"/>                                            | 250DI42 | Calf Support, Depth Adjustable, Padded (6.5" to 8" Wide, Depending On Chair Width, x 9.5" Tall)                     | \$512.00   |
| <b>Swing Away Legrest Options</b>                                   |         |                                                                                                                     |            |
| <input type="checkbox"/>                                            | 250DI44 | Heel Loops                                                                                                          | \$0.00     |
| <input type="checkbox"/>                                            | 250KG4  | Knee Pads, Integral Foam (2.5" Wide x 4.5" Tall)                                                                    | \$282.00   |
| <b>SEDEO Lower Extremity Swing Away Support</b>                     |         |                                                                                                                     |            |
| <input type="checkbox"/>                                            | 250VD01 | Lower Extremity Support, Padded, Swing Away, SEDEO, Right                                                           | \$512.00   |
| <input type="checkbox"/>                                            | 250VD02 | Lower Extremity Support, Padded, Swing Away, SEDEO, Left                                                            | \$512.00   |

| Positioning Supports & Accessories     |         |                                                                |          |
|----------------------------------------|---------|----------------------------------------------------------------|----------|
| OPTION                                 |         | Description                                                    | Retail   |
| <b>Pelvic Positioning Belts</b>        |         |                                                                |          |
| <input type="checkbox"/>               | 250A837 | Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"        | NC       |
| <input type="checkbox"/>               | 250A839 | Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"        | NC       |
| <input type="checkbox"/>               | 250A840 | Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2" | \$103.00 |
| <input type="checkbox"/>               | 250A845 | Bodypoint - Evoflex Pelvic Stabilizer Belt                     | \$216.00 |
| <input type="checkbox"/>               | 250A846 | Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension      | \$271.00 |
| <b>Accessories</b>                     |         |                                                                |          |
| <b>Ventilator Option</b>               |         |                                                                |          |
| <input type="checkbox"/>               | 250VT12 | Vent Tray Fixed ( <i>Pulmonetic LTV Models</i> )               | \$609.00 |
| <b>O2 Tank Holder</b>                  |         |                                                                |          |
| <input type="checkbox"/>               | 250A718 | O2 Tank Holder Mounted to Backrest ( <i>D and E Styles</i> )   | \$391.00 |
| <b>Wheelchair Tray Table</b>           |         |                                                                |          |
| <input type="checkbox"/>               | 2507A09 | Tray, Swing Away                                               | \$381.00 |
| <b>Backpack (Client Cost)</b>          |         |                                                                |          |
| <input type="checkbox"/>               | 250A1   | Black Utility                                                  | \$112.00 |
| <b>Medial Thigh Support</b>            |         |                                                                |          |
| <input type="checkbox"/>               | 2508A80 | Medial Thigh Support                                           | \$271.00 |
| <input type="checkbox"/>               | 2508A81 | Pommel Swing Away Hardware                                     | \$390.00 |
| <b>Misc. Accessories (Client Cost)</b> |         |                                                                |          |
| <input type="checkbox"/>               | 2508A59 | Transfer Bar ( <i>Pair</i> )                                   | \$493.00 |
| <input type="checkbox"/>               | 2508A52 | Ram Mount Phone                                                | \$77.00  |
| <input type="checkbox"/>               | 2508A53 | Ram Mount Tablet                                               | \$308.00 |
| <input type="checkbox"/>               | 2508A54 | Ram Mount Cup Holder                                           | \$77.00  |
| <input type="checkbox"/>               | 2508A58 | Rear View Mirror                                               | \$77.00  |
| <input type="checkbox"/>               | 250A611 | Crutch Holder                                                  | \$174.00 |
| <input type="checkbox"/>               | 2508A73 | Bag For Armrest                                                | \$77.00  |
| <input type="checkbox"/>               | 2508A93 | Wheelchair Cover                                               | \$240.00 |

| Joystick / Electronics                                 |         |                                                                                                                                                                                                    |            |
|--------------------------------------------------------|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|
| OPTION                                                 |         | Description                                                                                                                                                                                        | Retail     |
| <b>R-net Joystick (Remote)</b>                         |         |                                                                                                                                                                                                    |            |
| <input type="checkbox"/>                               | 250J121 | Standard LED Joystick ( <i>w/o Screen Display</i> )                                                                                                                                                | NC         |
| <input type="checkbox"/>                               | 250J116 | Color Joystick ( <i>2.3" Color LCD Screen Display, CJS1</i> )                                                                                                                                      | NC         |
| <input type="checkbox"/>                               | 250J118 | High Resolution Color Joystick, (CJS2+BT)<br><i>3.5" Color Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control (Client Cost)</i> | \$1,162.00 |
| <input type="checkbox"/>                               | 250J119 | Omit Remote Control ( <i>Required with arm mounted specialty controls</i> )                                                                                                                        | NC         |
| <b>Joystick Mounting</b>                               |         |                                                                                                                                                                                                    |            |
| <input type="checkbox"/>                               | 250JM34 | Fixed on Armrest                                                                                                                                                                                   | NC         |
| <input type="checkbox"/>                               | 250JM35 | Parallel Swing Away ( <i>Flip-back - Swing In/Out, Single Post - Swing Out</i> )                                                                                                                   | \$447.00   |
| <input type="checkbox"/>                               | 250JM36 | Link-It Mid Line Mount on Armrest ( <i>Swing Away</i> )                                                                                                                                            | \$447.00   |
| <b>Joystick Mount (Location)</b>                       |         |                                                                                                                                                                                                    |            |
| <input type="checkbox"/>                               | 250HM16 | Right Hand armrest                                                                                                                                                                                 | NC         |
| <input type="checkbox"/>                               | 250HM17 | Left Hand armrest                                                                                                                                                                                  | NC         |
| <b>Display Language Selection R-net</b>                |         |                                                                                                                                                                                                    |            |
| <input type="checkbox"/>                               | 250L1   | English                                                                                                                                                                                            | NC         |
| <input type="checkbox"/>                               | 250L3   | French                                                                                                                                                                                             | NC         |
| <input type="checkbox"/>                               | 250L7   | Spanish                                                                                                                                                                                            | NC         |
| <b>Joystick Handle (Optional) N/A with HD Joystick</b> |         |                                                                                                                                                                                                    |            |
| <input type="checkbox"/>                               | 250JH13 | Ball Handle                                                                                                                                                                                        | \$110.00   |
| <input type="checkbox"/>                               | 250JH15 | T-Handle                                                                                                                                                                                           | \$110.00   |
| <input type="checkbox"/>                               | 250JH24 | Bodypoint - Rubber Dome Joystick Handle                                                                                                                                                            | \$216.00   |
| <input type="checkbox"/>                               | 250JH25 | Bodypoint - U-Shaped Joystick Handle ( <i>4 in</i> )                                                                                                                                               | \$216.00   |
| <input type="checkbox"/>                               | 250JH26 | Bodypoint - Grooved Mushroom Joystick                                                                                                                                                              | \$216.00   |
| <b>Electronic Programmers</b>                          |         |                                                                                                                                                                                                    |            |
| <input type="checkbox"/>                               | 250PR13 | R-net PC Wireless Programmer                                                                                                                                                                       | \$660.00   |
| <b>R-net Power Modules</b>                             |         |                                                                                                                                                                                                    |            |
| <input type="checkbox"/>                               | 250ER11 | Power Module, R-net, 120A, Non-Expandable ( <i>Only available with Group 4</i> )                                                                                                                   | NC         |
| <input type="checkbox"/>                               | 250ER10 | Power Module, R-net, 120A, Expandable ( <i>Only available with Group 4</i> )                                                                                                                       | \$102.00   |
|                                                        | 250EX09 | Expandable 120A Controller - R-net                                                                                                                                                                 | \$0.00     |
|                                                        | 250EX10 | Harness for 120A Expandable Controller                                                                                                                                                             | \$0.00     |
| <input type="checkbox"/>                               | 250ER18 | Power Module, R-net, 90A, Non-Expandable ( <i>Available with Powered Seating</i> )                                                                                                                 | NC         |
| <input type="checkbox"/>                               | 250ER17 | Power Module, R-net, 90A, Expandable ( <i>Required for specialty controls or 3+ actuators</i> )                                                                                                    | \$102.00   |
|                                                        | 250EX15 | Expandable 90A Controller - R-net                                                                                                                                                                  | \$0.00     |
|                                                        | 250EX16 | Harness for 90A Expandable Controller                                                                                                                                                              | \$0.00     |

| Joystick / Electronics                                                                                                                                                                                                                                                                                      |         |                                                                                                                                                                                                                |            |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|
| OPTION                                                                                                                                                                                                                                                                                                      |         | Description                                                                                                                                                                                                    | Retail     |
| <b>Seating Control</b>                                                                                                                                                                                                                                                                                      |         |                                                                                                                                                                                                                |            |
| <b>Powered Seating Control</b>                                                                                                                                                                                                                                                                              |         |                                                                                                                                                                                                                |            |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250EO40 | Omit Thru-drive <i>(Requires Switch Control)</i>                                                                                                                                                               | NC         |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250EO37 | Thru-drive Control for 1 Actuator                                                                                                                                                                              | \$1,443.00 |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250EO38 | Thru-drive Control for 2+ Actuators                                                                                                                                                                            | \$2,903.00 |
| <b>Switch Control <i>(Available with Omit Thru-Drive Only) (1 actuator for Tilt or Elevate only)</i></b>                                                                                                                                                                                                    |         |                                                                                                                                                                                                                |            |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250EO33 | Dual Toggle for 1 actuator                                                                                                                                                                                     | NC         |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250EO34 | 2.5" Single Button Switch 1 actuator <i>(Button switch are not mounted)</i>                                                                                                                                    | NC         |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250EO35 | Switch-It Mini 1 button for 1 actuator                                                                                                                                                                         | NC         |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250EO36 | Switch-It Mini 2 button for 1 actuator                                                                                                                                                                         | NC         |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250EO41 | Ctrl+5 Switch Box with Voice Activation for 2 actuators or more                                                                                                                                                | NC         |
| <b>Additional Power Seating Control <i>(Available with Thru-Drive Only)</i></b>                                                                                                                                                                                                                             |         |                                                                                                                                                                                                                |            |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250QG11 | Ctrl+5 Switch Box With Voice Activation                                                                                                                                                                        | \$1,078.00 |
| <b>Dual Toggle Mounting Location</b>                                                                                                                                                                                                                                                                        |         |                                                                                                                                                                                                                |            |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250TN1  | Toggle Mount, Right                                                                                                                                                                                            | NC         |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250TN2  | Toggle Mount, Left                                                                                                                                                                                             | NC         |
| <b>Ctrl+5 Location <i>(If mounted on same side as joystick, it will Swing Away)</i></b>                                                                                                                                                                                                                     |         |                                                                                                                                                                                                                |            |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250EY19 | Right Armrest                                                                                                                                                                                                  | NC         |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250EY20 | Left Armrest                                                                                                                                                                                                   | NC         |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250EY21 | Back Cane, Right                                                                                                                                                                                               | NC         |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250EY24 | Back Cane, Left                                                                                                                                                                                                | NC         |
| <b>Ctrl+5 Switch Type</b>                                                                                                                                                                                                                                                                                   |         |                                                                                                                                                                                                                |            |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250EY22 | Push Buttons                                                                                                                                                                                                   | NC         |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250EY23 | Toggle Switches                                                                                                                                                                                                | \$356.00   |
| <b>Attendant Controls <i>(Prior Approval Required)</i></b>                                                                                                                                                                                                                                                  |         |                                                                                                                                                                                                                |            |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 2500A05 | Dual Control Module for Attendant                                                                                                                                                                              | \$1,225.00 |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 2500A06 | R-net Color LCD (CJSM1), Attendant Control Only                                                                                                                                                                | \$1,225.00 |
| <b>Attendant Mounting Location</b>                                                                                                                                                                                                                                                                          |         |                                                                                                                                                                                                                |            |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250DN5  | Right Back Cane                                                                                                                                                                                                | NC         |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250DN6  | Left Back Cane                                                                                                                                                                                                 | NC         |
| <b>R-net Programming</b>                                                                                                                                                                                                                                                                                    |         |                                                                                                                                                                                                                |            |
| <b>Assignable Buttons</b>                                                                                                                                                                                                                                                                                   |         |                                                                                                                                                                                                                |            |
| Assignable Button feature allows for customization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle. |         |                                                                                                                                                                                                                |            |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250PM12 | Include Assignable buttons                                                                                                                                                                                     | NC         |
| <b>Power Center Mount Footplate Extension Option</b>                                                                                                                                                                                                                                                        |         |                                                                                                                                                                                                                |            |
| <input type="checkbox"/>                                                                                                                                                                                                                                                                                    | 250IE03 | Include 7.5" Powered footplate extension for transfer assist <i>(Electronics programmed to allow footplates to be extended with or without legrest being raised.) Only Available with Powered Articulation</i> | \$592.00   |



1000 Creditstone Road, Unit 2, Concord ONL4K 4P8

Customer Service: (800) 263-3390, E-mail: [cscanada@sunmed.com](mailto:cscanada@sunmed.com), [www.sunrisemedical.ca](http://www.sunrisemedical.ca)  
Specifications subject to change without notice.

© 10/2026 Sunrise Medical Canada Inc. MK-129765-CAN Rev 12.