

Design Template

OCTANE



EIRGK03



Rider Information

Name: \_\_\_\_\_

Phone: \_\_\_\_\_

Email: \_\_\_\_\_

Height: \_\_\_\_\_ Weight: \_\_\_\_\_

Current Chair: \_\_\_\_\_ Current Cushion: \_\_\_\_\_

Ability/Diagnosis: \_\_\_\_\_

RGK Contact Authorization:

RGK would like to send you exclusive information, news and promotional offers on RGK.

I would like to be contacted by RGK:

Signature: \_\_\_\_\_ Date: \_\_\_\_\_

\*We do not pass your personal information onto any third party.

Account Information

Account Name : \_\_\_\_\_

Account #: \_\_\_\_\_ PO#: \_\_\_\_\_

Contact Name : \_\_\_\_\_

Phone: \_\_\_\_\_

Email: \_\_\_\_\_

| Ship to: | Account | Rider |
|----------|---------|-------|
|----------|---------|-------|

|       |       |
|-------|-------|
| Name: | _____ |
|-------|-------|

|          |       |
|----------|-------|
| Address: | _____ |
|----------|-------|

|       |
|-------|
| _____ |
|-------|


|             |       |
|-------------|-------|
| RGK Genius: | _____ |
|-------------|-------|

## Frame Styles


703XZ010 **Octane** (Max User Weight 110kg/242lbs.)  
 703XA070 **Open design**  
 703XA080 **Dual-Tube design/Closed frame design**

703XW009 **Quote Only** (No CAD for Approval)  
 703XW011 **Order Only** (No CAD for Approval)  
 703XW012 **Quote CAD drawing for approval** (1st CAD is NCO)  
 703XW010 **Order CAD drawing for approval** (1st CAD is NCO)

CAD approved by: 703XW013 Customer 703XW014 RGK Genius



Open Design



Dual-Tube Design

\$8,628  
 STD  
 NCO  
 STD  
 STD  
 \$550  
 \$300

## Measurements

**Please Note:** All RGK chairs are hand built and have a manufacturing tolerance of +/- 5mm (.25"). All measurements can be done in 5mm (.25") increments.

**Please select either Millimeters or Inches.**

**Measurements:** ☐ **Millimeters (mm.)** ☐ **Inches (in.)**

703XY010 **Rear Seat Width** Measured from inside of side guard to inside of opposite side guard. Or outside seat tube to outside of opposite seat tube near backrest. (Min. 250mm/9.84" and Max. 500mm/19.69")

Current Chair:  mm.  in. New Chair:  mm.  in.

703XY020 **Seat Depth/Seat Length** Measured from front of back post to the front edge of seat sling. Frame bend will start .5" in front of seat sling unless Extra Frame Length requested. (Min. 250mm/9.84" and Max. 500mm/19.69")

Current Chair:  mm.  in. New Chair:  mm.  in.

703XY030 **Extra Frame Length** Measured from front edge of seat sling to beginning of Front frame bend. (Seat length + Extra seat length)



Current Chair:  mm.  in. New Chair:  mm.  in.

703XY050 **Front Frame Height** Measured from floor to apex of seat tube. (Min. 380mm/14.96" and Max. 550mm/21.65")

Current Chair:  mm.  in. New Chair:  mm.  in.

703XY060 **Rear Frame Height** Measured from floor to top of seat tube at back edge of seat sling. (Min. 360mm/14.17" and Max. 550mm/21.65")

Current Chair:  mm.  in. New Chair:  mm.  in.



# Design Template

(New Options in Blue) **OCTANE**

**703XY070 Occupied Frame Length/Frame Length** Measured from the front of the backrest tube to the front edge of the footrest. (Min. 390mm/15.35" and Max. 750mm/29.53")

Current Chair:   New Chair:

**703XY072 Caster Stem Position** Measured from the front edge of the footrest to the centerline of the caster barrel. **703XY073** RGK to Optimize.

Current Chair:   New Chair:

**703XA020 Ergonomic Seat** Measured from front of back post to beginning of upward bend of seat tube. **\$121**

Current Chair:   New Chair:

**703XY080 Front Seat Taper/Seat Taper Design** Measured from inside of front frame tube at the front frame bend to inside of opposite front frame tube at the front frame bend.

Current Chair:   New Chair:

**703XY091 Footrest Width/Parallel Taper** Measured from inside of front frame tube at footrest adjustment to inside of opposite front frame tube at footrest adjustment. (Tip: when ordered together 703XY091 is usually 20mm (1") less than 703XY080).

Current Chair:   New Chair:

**703XY115 Front Footrest Height** Measured from ground up to the top of the Front footrest tube or cover.

Current Chair:   New Chair:

**703XY120 Rear Footrest Height** Measured from ground up to the top of the Rear footrest tube or cover. Standard is 20mm less than Front Footrest Height.




Current Chair:   New Chair:

**703XY125 Footrest Length** Measured from front edge of footrest to back edge of footrest. (Max: 140mm for certified weight)

Current Chair:   New Chair:

**703XB200 High Mount Footrest/High Level Footrest** Measured from ground up to the top of the Front footrest tube or cover. **\$201**

Current Chair:   New Chair:



# Design Template

(New Options in Blue) **OCTANE**

**703XY330 Stem Bolt/Caster Pin Width** Measured from centerline of caster pin to the centerline of the opposite caster pin. 703XY331 RGK to Optimize:

Current Chair:   New Chair:

**703XY231 Backrest Height** Measured from top of seat tube at rear of frame to top of back-post.

Current Chair:   New Chair:

**703XY360 Backrest Angle** Measured as the angle from top front of the backrest tube to the floor. (Max. -10° to +10°) -5° -3° 0° 3° 5°

Current Chair:  New Chair:

**703XG110 Inward Taper** Measured from inside backrest tube to inside of opposite backrest tube at the top of the backrest. \$146

Current Chair:   New Chair:

**703XG100 Ergonomic Backrest** Measured from top of seat tube at rear of frame to the start of the rearward bend of back-post. \$146

Current Chair:   New Chair:

**703XG120 Outward Flare** Measured from inside backrest tube to inside of opposite backrest tube at the top of the backrest. \$146






Current Chair:   New Chair:

**703XY390 Camber Angle** Place a straight edge across the outside of rear wheel and use an angle finder to show the degree of camber.

|                |    |            |          |    |
|----------------|----|------------|----------|----|
| Current Chair: | 0° | New Chair: | 703XY396 | 0° |
|                | 1° |            | 703XY397 | 1° |
|                | 2° |            | 703XY398 | 2° |
|                | 3° |            | 703XY399 | 3° |
|                | 5° |            | 703XY400 | 5° |

**703XH165 Center of Gravity** Measured from front of back post at the seat tube to the centerline of the rear axle.

Current Chair:   New Chair:



# Design Template


(New Options in Blue) **OCTANE**

**703XY460 Rear Wheel Gap** Measured from outer most point of chair (seat back post, side guard or armrest) to inside of rear tire.

|                |          |            |          |          |
|----------------|----------|------------|----------|----------|
| Current Chair: | .6"(15)  | New Chair: | 703XY464 | .6"(15)  |
|                | .78"(20) |            | 703XY465 | .78"(20) |
|                | 1"(25)   |            | 703XY466 | 1"(25)   |

**703XR060 Welded Calf Strap Mounting Loops** Measured from the ground to the top of where the calf strap is required. **703XR061 RGK to Optimize:** \$154

|                |                          |                          |            |                          |                          |
|----------------|--------------------------|--------------------------|------------|--------------------------|--------------------------|
| Current Chair: | <input type="checkbox"/> | <input type="checkbox"/> | New Chair: | <input type="checkbox"/> | <input type="checkbox"/> |
|----------------|--------------------------|--------------------------|------------|--------------------------|--------------------------|



## Footrest Design

|          |                      |     |
|----------|----------------------|-----|
| 703XY100 | <b>Round Design</b>  | STD |
| 703XY105 | <b>Square Design</b> | NCO |

## Footrest Styles

|          |  |       |
|----------|--|-------|
| 703XB020 | <b>Tubular Titanium Footrest</b>                                 | STD   |
| 703XB070 | <b>Angle adjustable aluminum footrest</b>                        | \$207 |
| 703XB080 | <b>Tubular footrest welded into frame</b> (No height adjustment) | \$460 |
| 703XB030 | <b>Flip up plastic footrest</b> (Dual-tube frame design only)    | \$352 |
| 703XB150 | <b>Heavy duty tubular footrest</b>                               | \$180 |

## Footrest Cover

|          |   |       |
|----------|---|-------|
| 793XB050 | <b>Lightweight 3k Carbon Fiber Footrest cover</b>                     | \$248 |
| 703XB155 | <b>Lightweight Spinerger CX6 Cast Carbon Fiber Footrest cover</b>     | \$245 |
| 703XB040 | <b>Plastic Cover</b> (Not available with 703XB070 & 703XB030)         | \$98  |
| 703XB170 | <b>Aluminum cover</b> (Not available with 703XB070 & 703XB030)        | \$98  |
| 703XB165 | <b>Rubber grip tape</b> (3M Safety walk. Not available with 703XB030) | \$86  |

## Front Wheels

|           |  |       |
|-----------|--|-------|
| 703XC060  | <b>Aluminum Center 4"x1.3" (Black Hub only)</b>  | STD   |
| 703XC110  | <b>Aluminum Center 4"x1" (Silver hub only)</b>   | \$232 |
| 703XC135  | <b>Aluminum Center 4"x1.25"</b>  | \$232 |
| 703XC1131 | Sliver hub & black tire  |       |
| 703XC1132 | Black hub & black tire   |       |
| 703XC1133 | Purple hub & black tire  |       |
| 703XC1134 | Green hub & black tire   |       |
| 703XC1135 | Red hub & black tire   |       |
| 703XC1136 | Blue hub & black tire  |       |
| 703XC155  | <b>4" x 1.25" Spinerger CX6 Cast Carbon fiber front wheel</b>  | \$235 |
| 703XC156  | <b>Cerakoted Spinerger CX6 Cast Carbon fiber front wheel</b> (Hub only. See Cerakot color options on page 8) | \$121 |
| 703XC070  | <b>Soft roll 5"x1.3"</b>   | NCO   |
| 703XC100  | <b>Pneumatic Aluminum Center 5"x1"</b>   | \$248 |
| 703XC120  | <b>Aluminum Center 5"x1"</b>   | \$232 |
| 703XC140  | <b>Aluminum Center 5"x1.4"</b> (Sliver hub & black tire)   | \$232 |

# Design Template

(New Options in Blue) **OCTANE**

|             |  |                 |              |               |         |
|-------------|--|-----------------|--------------|---------------|---------|
| 703XC145    | <b>Aluminum Center 6"x1.4"</b> (Sliver hub & black tire)       |                 |              |               | \$232   |
| 703XC150    | <b>Plastic Center 6"x1.4"soft roll</b> (Black hub & grey tire) |                 |              |               | \$232   |
| Forks       |  |                 |              |               |         |
| 703XC192    | <b>High Trail fork</b>   |                 |              |               | NCO     |
| 703XC210    | <b>Froglegs suspension fork</b>                                |                 |              |               | \$633   |
|             | 703XC211 Black fork  |                 |              |               | STD     |
| Wheel Lock  |  |                 |              |               |         |
| 703XD020    | <b>Underseat Scissor</b>                                       |                 |              |               | NCO     |
| 703XD030    | <b>Compact Scissor</b>   |                 |              |               | \$266   |
| 703XD035    | <b>Compact Scissor with Direct Mount</b>                       |                 |              |               | \$352   |
| 703XD050    | <b>Push To Lock</b>  |                 |              |               | NCO     |
| 703XD085    | <b>Eagle</b> (Indicate color)                                  |                 |              |               | \$468   |
| 703XD086    | <b>Eagle with Direct Mount</b> (Indicate color)                |                 |              |               | \$509   |
|             | 703XD088 Black   | 703XD089 Purple | 703XD090 Red | 703XD091 Blue |         |
| 703XD010    | <b>Omit Wheel Lock</b>   |                 |              |               | NCO     |
| Side Guards |  |                 |              |               |         |
| 703XE030    | <b>Screw on aluminum</b>                                       |                 |              |               | \$327   |
| 703XE080    | <b>Lift out aluminum</b>                                       |                 |              |               | \$532   |
| 703XE130    | <b>Fold flat aluminum</b>                                      |                 |              |               | \$670   |
| 703XE275    | <b>Screw on Spinergy CX6 Cast Carbon fiber</b>                 |                 |              |               | \$729   |
| 703XE100    | <b>Lift Out 3k Carbon fiber</b>                                |                 |              |               | \$766   |
| 703XE106    | <b>Lift Out Spinergy CX6 Cast Carbon fiber</b>                 |                 |              |               | \$1,015 |
| 703XE150    | <b>Fold flat 3k carbon fiber</b>                               |                 |              |               | \$1,052 |
| 703XE155    | <b>Fold flat Spinergy CX6 Cast Carbon fiber</b>                |                 |              |               | \$1,052 |
| 703XE070    | <b>Screw on ABS plastic</b>                                    |                 |              |               | \$167   |
| 703XE120    | <b>Lift out ABS plastic</b>                                    |                 |              |               | \$487   |
| 703XE170    | <b>Fold flat ABS plastic</b>                                   |                 |              |               | \$830   |
|             | 703XEM05 25mm above tire (Follow tire)                         |                 |              |               | NCO     |
|             | 703XEM06 Individual above tire (Follow tire - max. 25mm)       |                 |              |               | NCO     |
|             | 703XEM02 Level with tire                                       |                 |              |               | STD     |
|             | 703XEM07 Side guard length                                     |                 |              |               | STD     |
|             | 703XEM08 Side guard height                                     |                 |              |               | STD     |
| 703XE040    | <b>Screw on aluminum with fender</b>                           |                 |              |               | \$279   |
| 703XE090    | <b>Lift out aluminum with fender</b>                           |                 |              |               | \$579   |
| 703XE140    | <b>Fold flat aluminum with fender</b>                          |                 |              |               | \$707   |
| 703XE060    | <b>Screw on 3k carbon fiber with fender</b>                    |                 |              |               | \$952   |
| 703XE110    | <b>Lift out 3k carbon fiber with fender</b>                    |                 |              |               | \$1,132 |
| 703XE160    | <b>Fold flat 3k carbon fiber with fender</b>                   |                 |              |               | \$1,248 |
|             | 703XEM03 10mm above tire                                       |                 |              |               | NCO     |
|             | 703XEM01 Individual above tire (Follow tire - max. 25mm)       |                 |              |               | NCO     |
|             | 703XEM09 Center of tire  |                 |              |               | STD     |
|             | 703XEM10 Outside of tire                                       |                 |              |               | NCO     |

# Design Template

(New Options in Blue) **OCTANE**

| Armrest           |   |   |       |
|-------------------|---|---|-------|
| 703XF010          | <b>Lift out armrest with plastic sideguard</b>  | <b>Height</b> (Top of seatrail to top of armrest) | \$592 |
| 703XF020          | <b>Lift out armrest with aluminum sideguard</b>   | <b>Height</b> (Top of seatrail to top of armrest) | \$467 |
| Backrest Height   |   |   |       |
| 703XY230          | <b>Fixed height</b>   |   | STD   |
| 703XY245          | <b>Adjustable height</b> (Range: +/- 25mm)  |   | \$248 |
| Backrest Design   |   |   |       |
| 703XG010          | <b>Fixed backrest</b> (non-folding)   |   | STD   |
| 703XG020          | <b>Folding backrest with 2 position snap lock and angle adjustable</b> (Range: +/- 5°)  |   | \$692 |
| Push Handles      |   |   |       |
| 703XG070          | <b>Screw in push handles</b> (Not available with 703XY245)  |   | \$225 |
| 703XG080          | <b>Flip down push handles</b> (Not available with 703XY245)   |   | \$441 |
| 703XG085          | <b>Flip down push handles height adjustable</b> (Not available with 703XY245)   |   | \$720 |
| 703XG090          | <b>Stroller push handles height adjustable</b> (Measure height from floor to the top of the push handle)                          |   | \$441 |
| 703XG097          | <b>Integrated push handles</b>  |   | \$676 |
| Backrest Options  |   |   |       |
| 703XG040          | <b>Heavy duty backrest upgrade</b> (Not available with 703XY245)  |   | \$207 |
| 703XG130          | <b>Extra deep brace/rigidizer bar</b> (Recommended for JAY backs)   |   | \$115 |
| Center of Gravity |   |   |       |
| 703XH110          | <b>Square back adjustable of gravity</b> (Range: 25mm - 150mm)  |   | STD   |
| 703XH120          | <b>Round back adjustable of gravity</b> (Range: 60mm - 150mm)   |   | NCO   |
| 703XH100          | <b>Single post fixed center of gravity</b> (Range: 25mm - 150mm)  |   | \$727 |
| 703XH090          | <b>V-Style fixed center of gravity</b> (Range: 25mm - 150mm)  |   | \$727 |
| Axle Pin Options  |   |   |       |
| 703XK227          | <b>Standard axle pin</b>  |   | STD   |
| 703XK220          | <b>Titanium 130mm quick release axle pin</b>  |   | \$160 |
| 703XK240          | <b>Quad release axle pin</b>  |   | \$241 |
| Anti-Tip          |   |   |       |
| 703XJ010          | <b>Swing-away anti tip LEFT side</b> (Only available with square back adjustable center of gravity. Not available with 703XH140)  |   | \$167 |
| 703XJ020          | <b>Swing-away anti tip RIGHT side</b> (Only available with square back adjustable center of gravity. Not available with 703XH140) |   | \$167 |
| 703XJ100          | <b>Single sport anti-tip</b> (Centered mounted)   |   | \$504 |



## Frame and Footrest Finishes

**Non painted finishes**

703XN030 Bead Blast frame

STD

703XN050 Polished frame

\$1,332

**Painted finishes**703XN500 **CeraKote frame**

\$534

703XN600 **CeraKote components** (Component color will match frame if not otherwise specified. Components include: Forks, Wheel Lock clamps, and Camber Tube clamps.)

\$534

703XN800 **CeraKote Spinergy wheels** (See Rear Wheel Type on page 8. Hubs and Rims are CeraKoted on Rear Wheels, Hubs on Front Wheels.)

\$655

| Frame    | Components | Spinergy | Colors                | Frame    | Components | Spinergy | Colors                  |
|----------|------------|----------|-----------------------|----------|------------|----------|-------------------------|
| 703XN501 | 703XN601   | 703XN801 | <b>Midnight Black</b> | 703XN509 | 703XN609   | 703XN809 | <b>Rose Gold</b>        |
| 703XN502 | 703XN602   | 703XN802 | <b>Java Bronze</b>    | 703XN510 | 703XN610   | 703XN810 | <b>Retro Blue</b>       |
| 703XN503 | 703XN603   | 703XN803 | <b>Pure Titanium</b>  | 703XN511 | 703XN611   | 703XN811 | <b>Hot Pink</b>         |
| 703XN504 | 703XN604   | 703XN804 | <b>Fire Red</b>       | 703XN512 | 703XN612   | 703XN812 | <b>Wicked Green</b>     |
| 703XN505 | 703XN605   | 703XN805 | <b>Prime Blue</b>     | 703XN513 | 703XN613   | 703XN813 | <b>Bright White</b>     |
| 703XN506 | 703XN606   | 703XN806 | <b>Sahara Gold</b>    | 703XN514 | 703XN614   | 703XN814 | <b>Tangerine Orange</b> |
| 703XN507 | 703XN607   | 703XN807 | <b>Black Cherry</b>   | 703XN515 | 703XN615   | 703XN815 | <b>Sunbeam Yellow</b>   |
| 703XN508 | 703XN608   | 703XN808 | <b>Fresh Green</b>    | 703XN516 | 703XN616   | 703XN816 | <b>Matte Navy</b>       |

## Accent Color

|          |                                |       |
|----------|--------------------------------|-------|
| 703XN251 | <b>Black</b>                   | STD   |
| 703XN253 | <b>Red</b>                     | NCO   |
| 703XN255 | <b>Purple</b>                  | NCO   |
| 703XN252 | <b>Silver</b>                  | NCO   |
| 703XN254 | <b>Blue</b>                    | NCO   |
| 703XN250 | <b>Individual Accent Color</b> | \$318 |

## Seat Upholstery

|          |                                    |       |
|----------|------------------------------------|-------|
| 703XP040 | <b>Seat Sling Adjustable Strap</b> | STD   |
| 703XR050 | <b>Carbon Fiber Seat Plate</b>     | \$480 |

## Back Upholstery

|          |   |       |
|----------|---|-------|
| 703XP035 | <b>Sub4 Airtech Adjustable Strap</b>  | STD   |
| 703XP050 | <b>Airtech Adjustable Strap</b>   | NCO   |
| 703XP070 | Backrest upgrade to leather corners (Only with Airtech upholstery 703XP050) | \$121 |
| 703XP110 | Backrest zip pocket (Only with Airtech upholstery 703XP050)                 | \$127 |
|          | <b>Upholstery Stitching Color</b>   |       |
| 703XP550 | Black   | STD   |
| 703XP551 | Silver  | NCO   |
| 703XP552 | Red   | NCO   |
| 703XP553 | Orange  | NCO   |
| 703XP554 | Blue  | NCO   |
| 703XP555 | Purple  | NCO   |
| 703XP556 | Individual Pantone color number   | NCO   |
|          | <b>RGK Backrest Logo Color</b>  |       |
| 703XP557 | White (Standard if nothing selected)  | STD   |
| 703XP558 | Red   | STD   |



# Design Template

(New Options in Blue) **OCTANE**

|                            |  |  |                           |         |
|----------------------------|--|--|---------------------------|---------|
| 703XP559                   | Black  |  |                           | NCO     |
| 703XP560                   | Blue   |  |                           | NCO     |
| 703XP561                   | Purple   |  |                           | NCO     |
| 703XP562                   | Individual   | Pantone color number   |                           | NCO     |
| Rear Wheel Size            |  |  |                           |         |
| 703XY420                   | 24" (540) Rear wheel size  | 703XY430   | 25" (559) Rear wheel size | STD     |
| Rear Wheel Type            |  |  |                           |         |
| 703XK015                   | Proton (Not available with Spinergy anodized color options.)   |  |                           | STD     |
| 703XK320                   | Spinergy XLX (12 spoke, cross radial spokes.)  |  |                           | \$1,067 |
| 703XK050                   | Spinergy LX (12 spoke)   |  |                           | \$1,000 |
| 703XK060                   | Spinergy Spox (18 spoke. Not available with Spinergy anodized color options or spoke colors, purple, orange, pink and green.)  |  |                           | \$665   |
| 703XK070                   | Spinergy BLXL (18 bladed spoke.)   |  |                           | \$1,067 |
| 703XK150                   | Spinergy CLX (18 bladed spokes, carbon fiber rim. Spinergy color option available on hub only. Not available with silver hub)  |  |                           | \$3,139 |
| 703XK155                   | RGK EXCLUSIVE Spinergy CX6 Cast Carbon fiber (Cast Carbon Fiber Hub, Matte Black Carbon Fiber Rim and Black Bladed Spokes with 2 White Bladed Spokes Around Valve Stem. Not available with Spinergy anodized color options.) |  |                           | \$3,592 |
|                            | Spinergy spoke colors  |  |                           |         |
| 703XK650                   | Black  |  |                           | STD     |
| 703XK651                   | Red  |  |                           | NCO     |
| 703XK652                   | Purple   |  |                           | NCO     |
| 703XK653                   | White  |  |                           | NCO     |
| 703XK654                   | Blue   |  |                           | NCO     |
| 703XK655                   | Orange   |  |                           | NCO     |
| 703XK656                   | Yellow   |  |                           | NCO     |
| 703XK657                   | Pink   |  |                           | NCO     |
| 703XK658                   | Green  |  |                           | NCO     |
| 703XK180                   | Individual spoke design  |  |                           | \$213   |
|                            | Spinergy anodized colors (Rims and Hubs are anodized not spokes, on Spinergy BLXL (703XK075), Hubs only on 703XK150)   |  |                           |         |
| 703XK190                   | Black  |  |                           | STD     |
| 703XK202                   | Purple   |  |                           | \$409   |
| 703XK203                   | Blue   |  |                           | \$409   |
| 703XK201                   | Red  |  |                           | \$409   |
| 703XK204                   | Gold   |  |                           | \$409   |
| 703XK205                   | Individual hub/rims  |  |                           | \$568   |
| Supply Without Rear Wheels |  |  |                           |         |
| 703K040                    | Omit Rear Wheels   |  |                           | NCO     |
|                            | 703XK041 Wheel Type:   | 703XK042 Tire Type:  |                           |         |
| Handrims Style             |  |  |                           |         |
| 703XY485                   | Tab Length Wide  |  |                           | STD     |
| 703XY480                   | Tab Length Narrow  |  |                           | NCO     |
|                            | 703XL090   | Remove excess mounting tabs (Recommended for small wheel gaps) |                           | \$146   |
| 703XL010                   | Aluminum Hard Anodized 6 Tab   |  |                           | STD     |
| 703XL160                   | Elipse Ergonomic 6 Tab   |  |                           | \$631   |
| 703XL020                   | QGrip 6 Tab  |  |                           | \$400   |

# Design Template

(New Options in Blue) **OCTANE**

|  |  |                         |                 |                |         |
|--|--|-------------------------|-----------------|----------------|---------|
| 703XL040                                     | Titanium 6 Tab                                     |                         |                 |                | \$446   |
| 703XL035                                     | Natural-Fit LT                                     |                         |                 |                | \$743   |
| 703XL150                                     | Carbolife Gekko                                    |                         |                 |                | \$795   |
| 703XL170                                     | Carbolife Curve L                                  |                         |                 |                | \$695   |
| 703XL180                                     | Carbolife Curve L Tetra                            |                         |                 |                | \$885   |
| 703XL190                                     | Carbolife Quadro                                   |                         |                 |                | \$885   |
| 703XL200                                     | Spinergy TC  |                         |                 |                | \$736   |
| 703XL210                                     | Spinergy TC Max                                    |                         |                 |                | \$864   |
| Tires  |  |                         |                 |                |         |
| 703XM060                                     | Schwalbe One (Lightest option)                     |                         |                 |                | \$133   |
| 703XM150                                     | IRC EXERACER                                       |                         |                 |                | \$188   |
| 703XM010                                     | Kenda  | 703XM011 Grey           | 703XM012 Blue   | 703XM013 Red   | STD     |
| 703XM020                                     | Schwalbe Marathon Plus                             |                         |                 |                | \$98    |
| 703XM040                                     | Shox Solid Tire                                    |                         |                 |                | \$160   |
| 703XM050                                     | Kenda K885 Mountain Bike Tire                      |                         |                 |                | \$40    |
| Inner Tubes                                  |  |                         |                 |                |         |
| 703XY500                                     | Schraeder Valve Inner Tube                         |                         |                 |                | STD     |
| Cushions                                     |  |                         |                 |                |         |
| 703XP080                                     | Cushion foam                                       |                         |                 |                | \$121   |
|  | Density  | 703XP001 Soft           | 703XP002 Medium |                |         |
|  | Thickness  | 703XP003 50mm           | 703XP004 75mm   | 703XP005 100mm |         |
| Cushion Covers                               |  |                         |                 |                |         |
| 703XP085                                     | Nylon Cushion Cover (Black)                        |                         |                 |                | STD     |
| 703XP090                                     | Airlite Cushion Cover (Black)                      |                         |                 |                | \$98    |
| Jay Products - Please Include Jay Order Form |  |                         |                 |                |         |
| 703XP900                                     | Jay Cushion  | Refer to Jay Order Form |                 |                |         |
| 703XP910                                     | Jay Backrest                                       | Refer to Jay Order Form |                 |                |         |
| Accessories                                  |  |                         |                 |                |         |
| 703XR080                                     | Frontwheel Assembly                                |                         |                 |                | \$1,053 |
| 703XP120                                     | Zip Black Frame Protectors Neoprene Fasten (Black) |                         |                 |                | \$84    |
| 703XP420                                     | Velcro Black Frame Protectors Nylon Fasten (Black) |                         |                 |                | \$81    |
| 703XR010                                     | Black Seat Pouch Velcro Attach (Black)             |                         |                 |                | \$59    |
| 703XR200                                     | Black Rucksack with RGK Logo (Black)               |                         |                 |                | \$151   |
| 703XR090                                     | Tool Kit   |                         |                 |                | STD     |
| 703XS060                                     | Calf Strap   |                         |                 |                | STD     |
| Ordering Information                         |  |                         |                 |                |         |

**Estimated MSRP\***

**Quote**

**Order**

(\*Price does not include any one off requested options)

Specification Notes

