



Empulse R20

October 2022



EISP9

Mark For:

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ PV: _____ PC: _____

Dealer Phone: _____ Fax: _____

Confirmation Email: _____

Submitting for:

Quote Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ PV: _____ PC: _____

Ship To Phone: _____

Empulse R20 Power Assist

Each Empulse R20 includes the following components:

- Power Assist Motor
- Adjustable attachment hardware
- Lithium Battery
- Modern Backpack
- Attendant controller
- Device Charger

R20

Selection Required.

	Item Code	Description		Retail
<input type="checkbox"/>	169PC64	Empulse R20		\$4,995.00

Wheelchair Attachment Type

Selection Required. See Images

	Item Code	Description		Retail
<input type="checkbox"/>	9V49	Folding Ultralight Wheelchair <i>Round tube hardware for diameter between 0.8" - 1.3" (22mm - 32mm)</i>		NC
<input type="checkbox"/>	9V50	Tilt-in-Space Wheelchair <i>Flat surface mount hardware</i>		NC

Wheelchair Seat Width

Selection Required.

	Item Code	Description		Retail
<input type="checkbox"/>	169S4	13 inches		NC
<input type="checkbox"/>	169S5	14 inches		NC
<input type="checkbox"/>	169S6	15 inches		NC
<input type="checkbox"/>	169S7	16 inches		NC
<input type="checkbox"/>	169S8	17 inches		NC
<input type="checkbox"/>	169S9	18 inches		NC
<input type="checkbox"/>	169S10	19 inches		NC
<input type="checkbox"/>	169S11	20 inches		NC
<input type="checkbox"/>	169S40	21 inches		NC
<input type="checkbox"/>	169S27	22 inches		NC
<input type="checkbox"/>	169S41	23 inches		NC

Wheelchair Rear Seat Height

Selection Required.

	Item Code	Description		Retail
<input type="checkbox"/>	169RH15	15 inches		NC
<input type="checkbox"/>	169RH16	15.5 inches		NC
<input type="checkbox"/>	169RH17	16 inches		NC
<input type="checkbox"/>	169RH18	16.5 inches		NC
<input type="checkbox"/>	169RH19	17 inches		NC
<input type="checkbox"/>	169RH26	17.5 inches		NC
<input type="checkbox"/>	169RH20	18 inches		NC
<input type="checkbox"/>	169RH27	18.5 inches		NC
<input type="checkbox"/>	169RH21	19 inches		NC
<input type="checkbox"/>	169RH37	19.5 inches		NC
<input type="checkbox"/>	169RH22	20 inches		NC
<input type="checkbox"/>	169RH38	20.5 inches		NC

Accessories

	Item Code	Description		Retail
<input type="checkbox"/>	169AG33	Additional Lithium Battery (5.8 Ah)		\$1,150.00

Compatibility Information

For additional information reference the Owner's Manual or Mounting Instruction Sheet

R20 Arm Size		Rear Seat Height										
		15.5	16	16.5	17	17.5	18	18.5	19	19.5	20	20.5
Seat width	13*	XS	XS	XS	XS	XS	XS	S	S	S	S	S
	14*	XS	XS	XS	XS	XS	S	S	S	S	S	S
	15	XS	XS	XS	XS	XS	S	S	S	S	S	S
	16	XS	XS	XS	XS	XS	S	S	S	S	S	S
	17	XS	S	S	S	S	S	S	S	S	S	S
	18	S	S	S	S	S	L	L	L	L	L	L
	19	S	S	S	S	S	L	L	L	L	L	L
	20	S	S	S	S	L	L	L	L	L	L	L
	21	S	L	L	L	L	L	L	L	L	L	L
	22	S	L	L	L	L	L	L	L	L	L	L
	23	S	L	L	L	L	L	L	L	L	L	L

* = Not recommended for Zippie TS 12" Wheels

 High likelihood of interference. May require adjustment of wheelchair components.



The Empulse R20 is not compatible with the following:

- Folding Push Handles
- Omit Push Handles
- Tilt-in-Space Foot-Activated Tilt Lever
- Tilt-in-Space Powered Tilt
- Tilt-in-Space Battery & Vent Tray
- Crutch Holder

The Empulse R20 may have interference with:

- O2 Holder
- Backpack
- Adj. and Extension Stroller Handles
- Reclining Backrest
- Backrest Stabilizer Bar
- Foot Tilt Mechanism
- Foot Operated Wheel Lock

Tilt-in-Space Flat surface Clamp



Folding Ultralight Round Clamp

