



ORDER FORM

Everyday Activity Seat

EILK03

Activity chair for children or young adults with moderate needs who still require pelvic stability during the day.

Customer

Account			
Account #		Date	
☐ Quote	☐ Order	PO #	
Buyer			
Mark For			

Shipping Address

Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

Everyday Activity Seat

Size 1 Weight Capacity: 66 lbs | Size 2 Weight Capacity: 110 lbs | Size 3 Weight Capacity: 220 lbs

✓ OPTION	DESCRIPTION	RETAIL
Base Frame		\$
Selection Required.		
☐ 151BF04	Size 1 Seat Base, Backrest and Hi-Low Chassis with Gas Spring	\$3,537
☐ 151BF05	Size 2 Seat Base, Backrest and Hi-Low Chassis with Gas Spring	\$3,642
☐ 151BF06	Size 3 Seat Base, Backrest and Hi-Low Hydraulic Assist Chassis	\$7,179
Chair Cover		\$
Selection Required. Must match base frame size.		
☐ 151IU16	Size 1 Vinyl Seat and Back Covers - Blue	\$409
☐ 151IU17	Size 2 Vinyl Seat and Back Covers - Blue	\$467
☐ 151IU18	Size 3 Vinyl Seat and Back Covers - Blue	\$773
Lateral Support		\$
Must match base frame size.		
☐ 151ZY19	Size 1 Rigid Chest Laterals - Black Vinyl (Pair)	\$434
☐ 151ZY20	Size 2 Rigid Chest Laterals - Black Vinyl (Pair)	\$434
☐ 151ZY21	Size 3 Rigid Chest Laterals - Black Vinyl (Pair)	\$585
☐ 151ZY22	Size 1 Swing-Away Chest Laterals - Black Vinyl (Pair)	\$537
☐ 151ZY23	Size 2 Swing-Away Chest Laterals - Black Vinyl (Pair)	\$537
☐ 151ZY24	Size 3 Swing-Away Chest Laterals - Black Vinyl (Pair)	\$719
☐ 151KE47	Size 1 and 2 Seat Protraction Pads - Black Vinyl (Pair)	\$527
☐ 151KE48	Size 3 Seat Protraction Pads - Black Vinyl (Pair)	\$563
☐ 151ZY25	Chest Support - Size 1 & 2 (Requires rigid or swing-away laterals)	\$204
☐ 151ZY26	Chest Support - Size 3 (Requires rigid or swing-away laterals)	\$204
Harness		\$
☐ 151NH105	Butterfly Harness (Small)	\$415
☐ 151NH106	Butterfly Harness (Medium)	\$415
☐ 151NH107	Butterfly Harness (Large)	\$415
Upper Extremity Positioning		\$
Must match base frame size. Tray requires armrests.		
☐ 151AA167	Size 1 Armrests (Pair)	\$365
☐ 151AA168	Size 2 Armrests (Pair)	\$365
☐ 151AA169	Size 3 Armrests (Pair)	\$473
☐ 151A217	Size 1 Tray - Plastic Gray	\$399
☐ 151A218	Size 2 Tray - Plastic Gray	\$399
☐ 151A219	Size 3 Tray - Plastic Black	\$399



Flip Away Laterals



**Flip Away Laterals and
Chest Support**



Protraction Pads



Armrest

Everyday Activity Seat

✓ OPTION		DESCRIPTION	RETAIL
Head Supports			\$
☐	151Y053	Size 1 and 2 Flat Headrest Support and Cushion - Black Vinyl Must match base frame size	\$369
☐	151Y052	Size 1 and 2 Laterals for Flat Headrest - Black Vinyl (Pair) Must match base frame size. Includes width adjustable lateral pads & covers. Requires Flat Headrest Support and Cushion.	\$314
☐	151Y055	Size 3 Contoured Headrest - Black Neoprene Must match base frame size	\$543
☐	151Y056	Size 1 and 2 Contoured Headrest Adapter Bracket	\$85
☐	151AW24	Size 1 and 2 Contoured Headrest Support Hardware	\$396
☐	151Y027	Size 1 and 2 Contoured Headrest Cushion - Black	\$188
Foot Positioning			\$
☐	151F130	Size 1 One Piece Footplate - Metal	\$531
☐	151F131	Size 2 One Piece Footplate - Metal	\$531
☐	151F132	Size 3 One Piece Footplate - Metal	\$531
☐	1518A107	Small Sandals w/ Straps (Pair)	\$273
☐	1518A108	Medium Sandals w/ Straps (Pair)	\$291
☐	1518A109	Large Sandals w/ Straps (Pair)	\$315
☐	1519A25	Size 1 and 2 Sandal Raiser 1" (Each) Requires Sandal Selection, must match Sandal size.	\$63
☐	1519A26	Size 3 Sandal Raiser 1" (Each) Requires Sandal Selection, must match Sandal size.	\$70
	QTY		
	QTY		
Pelvic Fixation			\$
☐	151ZZ36	Size 1 Pelvic Belt 4-Point	\$164
☐	151ZZ37	Size 2 Pelvic Belt 4-Point	\$177
☐	151ZZ38	Size 3 Pelvic Belt 4-Point	\$188
☐	151KE49	Size 1 Side Pads - Black Vinyl (Pair)	\$574
☐	151KE50	Size 2 Side Pads - Black Vinyl (Pair)	\$682
☐	151KE51	Size 3 Side Pads - Black Vinyl (Pair)	\$870
☐	151KE52	Size 2 Adduction Side Pad - Black Vinyl	\$324
☐	151KE53	Size 3 Adduction Side Pad - Black Vinyl	\$324
☐	1518A111	Size 1 Wedge Shape Pommel - Black Vinyl	\$245
☐	1518A112	Size 2 Wedge Shape Pommel - Black Vinyl	\$301
☐	1518A113	Size 3 Wedge Shape Pommel - Black Vinyl	\$308



Flat Headrest



Contoured Headrest



Footplate



Sandals



Pelvic Belt 4-Point



Wedge Shape Pommel



Side Pads

Technical Specifications - Imperial			
Frame Sizing Specification	Size 1	Size 2	Size 3
Max Load	66 lbs	110 lbs	220 lbs
Weight (Seat and Hi-Low Chassis)	65 lbs	75 lbs	108 lbs
Seat Height (Hi-Low Chassis)	9.8" - 23.6"	12.6" - 23.6"	17.7" - 23.6"
Seat Width	6.3" - 11"	10.2" - 15"	13.4" - 19.7"
Seat Depth	5.9" - 11"	9" - 15"	13" - 18.9"
Lateral Distance (Chest)	5.9" - 8.7"	6.7" - 9.4"	9.4" - 14.6"
Back Height	11.8" - 15.7"	15.7" - 20.5"	17.7" - 24.8"
Back Angle	30° Supine / 10° Prone	30° Supine / 10° Prone	30° Supine / 10° Prone
Armrest Height	5.9" - 9.8"	7.9" - 11.8"	7" - 12.7"
Lower Leg Length	5.9" - 9.8"	9.8" - 17.7"	9.8" - 17.7"
Lower Leg Angle	-10° to +50°	-10° to +50°	-10° to +50°
Footplate Angle	-15° to +10°	-15° to +10°	-15° to +10°
Tilt-In-Space	30° Supine / 10° Prone	30° Supine / 10° Prone	30° Supine / 10° Prone

Technical Specifications - Metric			
Frame Sizing Specification	Size 1	Size 2	Size 3
Max Load	30 kg	50 kg	100 kg
Weight (Seat and Hi-Low Chassis)	29.5 kg	34 kg	49 kg
Seat Height (Hi-Low Chassis)	250 mm - 600 mm	320 mm - 600 mm	450 mm - 600 mm
Seat Width	160 mm - 280 mm	260 mm - 380 mm	340 mm - 500 mm
Seat Depth	150 mm - 280 mm	230 mm - 380 mm	330 mm - 480 mm
Lateral Distance (Chest)	150 mm - 220 mm	170 mm - 240 mm	240 mm - 370 mm
Back Height	300 mm - 400 mm	400 mm - 520 mm	450 mm - 630 mm
Back Angle	30° Supine / 10° Prone	30° Supine / 10° Prone	30° Supine / 10° Prone
Armrest Height	150 mm - 250 mm	200 mm - 300 mm	180 mm - 325 mm
Lower Leg Length	150 mm - 250 mm	250 mm - 450 mm	250 mm - 450 mm
Lower Leg Angle	-10° to +50°	-10° to +50°	-10° to +50°
Footplate Angle	-15° to +10°	-15° to +10°	-15° to +10°
Tilt-In-Space	30° Supine / 10° Prone	30° Supine / 10° Prone	30° Supine / 10° Prone

Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2
Concord, ON L4K 4P8, Canada



MK-129846-CAN Rev. 14 © 04/2025

