



ORDER FORM

October 2024 REV 1


Q700M^{HD}

AADL- W581
EIPW25

Customer	Ship To	
Account #	Name	
Date	Address	
PO#	City	
Buyer	Province	Postal Code
Mark For	Phone	

Client Measurements

a) Top of head to seat	in
b) Top of shoulder to seat	in
c) Underarm to seat	in
d) Elbow to seat	in
e) Top of leg to seat	in
f) Shoulder width	in
g) Chest width	in
i) Hip width	in
j) Actual sitting depth	in
m) Foot length	in
n) Knee to foot	in
o) Seat width	in



Notes

Chair Specifications

Typical Weight w/out Batteries	330 lbs.
Typical Weight with Batteries	410 lbs.
Power Seating	50° of CG Tilt and Power Recline with 4" mechanical sheer reduction
Seat Type	SEDEO® PRO™ / JAY® Syncra™
Seat Size Range	Width 16" to 22" Depth 14" to 22"
Seat-to-Floor Height Range	16.5" to 19"

Mid-Wheel Drive Base Specifications

Drive Wheel Position	Mid-wheel drive
User Weight Capacity	350 lbs. HD / 450 lbs. HD
Base Size	Width 24.5", Length 36.5"
Minimum Turning Radius	19.8" (excludes footplates)
Drive Control	R-net PM 120A / R-net PM 90A
Driving Range	Up to 18 miles (w/ Group 24 Batteries)
Drive Motors	4-pole 5 mph HD
Drive Wheels	14" w/flat-free inserts (black)
Casters	6" front / 6" rear (black)
Charger	10A Dual Mode (Compact)
Battery Size	Group 24 / Group 34
Color	8 shroud colors and 7 color accents
Maximum Ground Clearance	3"

BASE			
OPTION		DESCRIPTION	RETAIL
Heavy Duty Mid Wheel Drive Base <i>End Item Number EIPW25 (PRO™ Seat)</i>			\$
250PC199		Quickie Q700 M NPC HD No Power Option for SEDEO® PRO™ Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 5 mph, Weight Capacity = 450 lbs., STFH = 16.5". (AADL Contract Price)	\$9,822.00
250PC200		Quickie Q700 M SPC HD Single Power Option for SEDEO® PRO™ Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 5 mph, Weight Capacity = 450lbs., STFH = 16.5". (AADL Contract Price)	\$9,822.00
250PC204		Quickie Q700 M MPC HD Tilt/Recline for SEDEO® PRO™ Seat Included in Base Pricing: Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 5 mph, Weight Capacity = 350lbs., STFH = 16.5". (AADL Contract Price)	\$9,822.00
User Weight			\$
250WG15		301 to 350 lbs.	NC
250WG25		351 to 450 lbs.	NC
Motor Package			\$
250PF98		4-pole 8 kph / 5 mph	NC
Tracking			\$
250MI01		SureTrac® Module & Motor Encoders <i>(Improves drive tracking)</i>	\$1,671.00
Drive Wheels  P. 23			\$
250DW58		Drive Wheel, 14", Dual Density, Knobby, Black	NC
Caster Wheels  P. 23			\$
250CW55		Caster Wheels 6", Solid, Black	NC
Shroud Color			\$
250C119		Matte Black	NC
250C165		Candy Apple Red	NC
250C166		Candy Blue	NC
250C167		Candy Purple	NC
250C176		Green Apple	NC
250C178		Pearl Pink	NC
250C182		Silver	NC
250C125		Pearl White	NC

BASE

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

Color Accents

\$

Note: Color Accents and Shroud Colors may not match.

250EZ20	Green	NC
250EZ22	Black	NC
250EZ21	Blue	NC
250EZ23	Orange	NC
250EZ24	Red	NC
250EZ25	Silver	NC
250EZ19	White	NC

ISO Transit Tie Downs

\$

250T050	ISO Transit Package	NC
---------	---------------------	----

LED Lights



P. 23

(Client Cost)

\$

250LI04	LED Lights and Indicators	\$1,207.00
250LI08	LED Fender Lights	\$636.00
250LI09	LED Fender Lights (Available with Multi-Power only))	NC
250LI02	Omit LED Lights	NC

Seat Heights

\$

Note: Available seat-to-floor heights are dependent on 2 options: **1-** Battery Size and **2-** Powered Seating.

Use the matrix below to determine available ranges for your desired options. Measured from front of seat pan to floor with 0° frame angle.

					PRO ADVANCED Recline		
		Seat Height	Power Seating	No Power Seating	w/o Elevate	w/ Elevate	
250HT103		16.5"	Group 34	Group 34	N/A	N/A	NC
250HT55		17"	Group 34	GP24 or GP34	GP34	GP34	NC
250HT104		17.5"	Group 34	GP24 or GP34	GP34	GP34	NC
250HT56		18"	GP24 or GP34	GP24 or GP34	GP24 / GP34	GP34	NC
250HT105		18.5"	GP24 or GP34	GP24 or GP34	GP24 / GP34	GP24 / GP34	NC
250HT57		19"	GP24 or GP34	GP24 or GP34	GP24 / GP34	GP24 / GP34	NC

Batteries and Chargers



P. 23

\$

Group 24 Battery Options (Recommended for max daily range)

250AG25	GP24 Flat Top Pair w/Strap MK Gel	\$941.00
250AG26	GP24 Flat Top Pair w/Strap Brightway AGM	\$901.00

BASE			
OPTION		DESCRIPTION	RETAIL
Group 34 Battery Options (Required for lower seat-to-floor heights)			
250AG27		GP34 Flat Top Pair w/Straps MK Gel	\$941.00
250AG29		GP34 Flat Top Pair w/Straps Brightway AGM	\$850.00
Omit Batteries			
250AG24		Battery Harnesses Only <i>(Omit Batteries)</i>	NC
Chargers			
250BJ28		10A Dual Mode (Compact)	NC
POWERED SEATING			
OPTION		DESCRIPTION	RETAIL
Power Seating Options  P. 23			\$
Seat Tilt (AADL W992 for Power Tilt)			
250PA36		Posterior 50° HD Tilt	\$8,791.00
Seat Recline (Available up to 350 lbs.) (AADL W992TR for Power Recline)			
250PA39		150° HD Powered Recline w/ 4" Shear Reduction	\$8,791.00
Seat Tilt/Recline (Available up to 350 lbs.) (AADL W992TR for Power Tilt/Recline)			
250PA43		Posterior 50° HD Tilt and 150° HD Recline Powered, w/ 4" Shear Reduction	\$14,179.00
250PA63		Posterior 50° HD Tilt and PRO ADVANCED 165° HD Recline Powered, w/ 7" Shear Reduction	\$14,678.00
Power Recline Backrest Preset Home Angle			
250QK1		90°	NC
250QK2		95° (Default Position)	NC
Pre Set Seat Frame Angle (No Power Seating Options, Increases Seat to Floor Height)			\$
250FA10		0° (Adds 0" to STFH)	NC
250FA3		3° (Adds 1" to STFH)	NC
Backrest Frame Type (N/A with Powered Recline)			\$
250BF10		Tool Angle Adjustable Back, 90°-120°	NC
250BF11		User Adjustable Manual Semi-Recline Back, 90°-130°, Right <i>(Adjustment lever mounted on opposite side of joystick)</i>	\$812.00
250BF12		User Adjustable Manual Semi-Recline Back, 90°-130°, Left <i>(Adjustment lever mounted on opposite side of joystick)</i>	\$812.00
Backrest Frame Option (N/A with PRO ADVANCED)			\$
250BF13		Flip Down, Quick Release	\$176.00
Back Cushion Type			\$
250BT134		SEDEO® Back Cushion <i>(Choose Seating in Section A) (N/A with PRO ADVANCED)</i>	NC
250BT138		JAY® Syncra™ <i>(Choose Seating in Section B)</i>	NC
250BT135		JAY® J3® <i>(Choose Seating in Section C)</i>	NC
250BT137		SEDEO® Aftermarket Back Adapter <i>(Height adjustable round tubes for traditional mounting)</i>	\$624.00
250BT133		OMIT <i>(For personalized mounting. No headrest mount)</i>	NC

POWERED SEATING

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

A - SEDEO® Backrest Dimensions (N/A with PRO ADVANCED) Seating Benefits \$
Back Width

250BN28	400mm / 16"	NC
250BN29	460mm / 18"	NC
250BN30	500mm / 20"	NC
250BN31	560mm / 22"	NC

Back Shell Height

Back Shell Height			Tool Angle Adjust.	150° Recline	Installed Back Height	
250B255		460mm / 18"	22"-23"	23"-25"		NO
250B256		560mm / 22"	24"-27"	26"-29"		NO
250B257		660mm / 26"	28"-30"	30"		NO

Installed Back Height - (From Pan to Top of Back)

250JR32	22"	NC
250JR33	23"	NC
250JR34	24"	NC
250JR35	25"	NC
250JR36	26"	NC
250JR37	27"	NC
250JR38	28"	NC
250JR39	29"	NC
250JR40	30"	NC

B - JAY® Syncra™ Seating Seating Benefits
B - JAY® Syncra™ Back P. 22 \$

250JP14	JAY® Syncra™ Back 16-20" Wide with 3" Contour Depth	\$1,319.00
250JP15	JAY® Syncra™ Back 22" Wide with 3" Contour Depth	\$1,319.00

B - JAY® Syncra™ Back Width \$

Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.

250JQ07	16" Wide (Adjustment range 16 to 17" Wide)*	NC
250JQ08	18" Wide (Adjustment range 18 to 19" Wide)*	NC
250JQ09	20" Wide (Adjustment range 20 to 21" Wide)*	NC
250JQ10	22" Wide (Adjustment range 22 to 23" Wide)	NC

*Additional +1" of shell width adjustment is achievable. Requires purchase of 1" larger foam and upholstery separately.

POWERED SEATING

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

B - JAY® Syncra™ Back Shell Options and Height Ranges \$

JAY® Syncra™ Back shells are mounted with a minimim 5.5" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field within the ranges shown based on the back frame type.

JAY® Syncra™ Shell Height

250JR42	16" Tall Shell (22" to 24" Back Height)	NC
250JR43	18" Tall Shell (24" to 26" Back Height)	NC
250JR44	20" Tall Shell (26" to 28" Back Height)	NC
250JR45	22" Tall Shell (28" to 30" Back Height)	NC

JAY® Syncra™ Back Installed Back Height (From Pan to Top of Back)

250JR50	22" (16" Shell Height)	NC
250JR51	23" (16" Shell Height)	NC
250JR52	24" (16" Shell Height)	NC
250JR53	24" (18" Shell Height)	NC
250JR54	25" (18" Shell Height)	NC
250JR55	26" (18" Shell Height)	NC
250JR56	26" (20" Shell Height)	NC
250JR57	27" (20" Shell Height)	NC
250JR58	28" (20" Shell Height)	NC
250JR67	28" (22" Shell Height)	NC
250JR68	29" (22" Shell Height)	NC
250JR69	30" (22" Shell Height)	NC

B - JAY® Syncra™ Back Upholstery \$

250BU96	Syncra™ Premium Leatherette	NC
250BU97	Syncra™ Stretch	NC

C - JAY® J3® Back **Seating Benefits**

C - JAY® J3® Back - Contour P. 23 \$

250JP01	J3® PA Back for 16"-20" Wide	\$899.00
250JP02	J3® PA Back for 22" Wide	\$1,912.00
250JP11	J3® PL Back for 16"-20" Wide	\$971.00
250JP13	J3® PD Back for 16"-20" Wide (18" and 21" shell height only)	\$1,052.00
250JP18	J3® Plus MT Back for 22"-26" Wide	\$1,918.00

POWERED SEATING

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

C - JAY® J3® Back Width

\$

Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.

250JQ02	16"	NC
250JQ03	18"	NC
250JQ04	20"	NC
250JQ05	22"	NC
250JQ15	24"	NC
250JQ16	26"	NC

C - JAY® J3® Back Shell Options and Height Ranges

\$

JAY® J3® shells are mounted with a minimim 4" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field within the ranges shown, based on the back frame type and the amount of usable hardware adjustability.

JAY® J3® Shell Height

Tool Angle Adjust.

150° Recline

PRO ADVANCED

250JR41	15" tall shell (Upper Thoracic Med) for 22" Wide	22"-23"	22"-23"	21"-24"	Installed Back Height	NC
250JR14	16.5" tall shell (Upper Thoracic Short)	21"-23"	22"-24"	21"-24"		NC
250JR15	18" tall shell (Upper Thoracic Med)	22"-24"	23"-25"	22"-26"		NC
250JR16	19.5" tall shell (Upper Thoracic Tall)	24"-26"	24"-26"	23"-27"		NC
250JR18	18" tall shell (Shoulder Height) for 22" wide	24"	24"	21"-25"		NC
250JR79	20" tall shell (Mid Thoracic) for 22" - 26" wide	21"-22"	21"-23"	N/A		NC
250JR17	21" tall shell (Shoulder Height)	24"-26"	25"-28"	25"-28"		NC
250JR27	22.5" tall shell (Shoulder Height)	26"-28"	27"-30"	26"-29"		NC
250JR80	23" tall shell (Mid Thoracic) for 22" - 26" wide	23"-24"	23"-25"	N/A		NC
250JR28	24" tall shell (Shoulder Height)	27"-29"	28"-30"	28"-30"		NC

JAY® J3® Installed Back Height (From Pan to Top of Back)

250JR29	21"	NC
250JR19	22"	NC
250JR20	23"	NC
250JR21	24"	NC
250JR22	25"	NC
250JR23	26"	NC
250JR24	27"	NC
250JR25	28"	NC
250JR30	29"	NC
250JR31	30"	NC

POWERED SEATING

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

C - JAY® Anterior Thoracic Support **P. 23** **Seating Benefits** **\$**

Classic-Structure Pull Down Harness (4 point strap, only available with JAY® backrest)

250K10	Medium (ATS8302)	\$317.00
250K11	Large (ATS8303)	\$317.00

JAY® Lateral Thoracic Swing-Away Supports (Requires Pad, Upholstery and Mount) **\$**
JAY® Lateral Thoracic Support Pads

250JW35	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right	\$74.00
250JW38	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left	\$74.00
250JW36	Lateral Thoracic (L x H) 100 mm x 152 mm / 4" x 6" Pad - Right	\$74.00
250JW39	Lateral Thoracic (L x H) 100 mm x 52 mm / 4" x 6" Pad - Left	\$74.00
250JW37	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right	\$74.00
250JW40	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left	\$74.00

JAY® Lateral Thoracic Support Pad Upholstery

250NU01	Syncra™ Premium Leatherette	NC
250NU02	Syncra™ Stretch	NC

JAY® Lateral Thoracic Swing-Away Mount

250JW41	Lateral Thoracic Swing Away Bracket - Right	\$329.00
250JW42	Lateral Thoracic Swing Away Bracket - Left	\$329.00

JAY® Syncra™ Fixed Lateral Support Accessory **Seating Benefits** **\$**

JAY® Fixed Lateral Support Accessory Hardware (JAY® Fixed Lateral Support is height adjustable)

250JW33	Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm / 3.5" x 8" - Right	\$329.00
250JW34	Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm / 3.5" x 8" - Left	\$329.00

Headrests (Requires Pad, Upholstery, & Mount) **P. 25** **Seating Benefits** **\$**
Headrest Pad

250HE36	SEDEO® Pad Small With Cover (LxH) 240 x 150 mm / 10" x 6" <i>Can only be used with Multi Position Headrest Mount with Quick Release and SEDEO® Upholstery</i>	\$224.00
250HE38	Whitmyer® Plush Headrest Pad 203 mm / 8"	\$371.00
250HE39	Whitmyer® Plush Headrest Pad 254 mm / 10"	\$371.00
250HE44	Whitmyer® Plush Headrest Pad 356 mm / 14"	\$371.00
250HE46	Whitmyer® Adjust-A-Plush Headrest Pad Medium (Adjustable from 7 to 16" Wide x 4" Tall)	\$438.00
250HE48	OMIT Headrest Pad	NC

POWERED SEATING			
OPTION		DESCRIPTION	RETAIL
Headrest Pad Upholstery			
250NU10		SEDEO® Upholstery	NC
250NU04		Syncra™ Premium Leatherette	NC
250NU05		Syncra™ Stretch	NC
250NU11		Whitmyer® Lycra	NC
Headrest Mount  P. 25			
250MH32		SEDEO® Multi-Position Headrest Mount with Quick Release <i>(Only available with SEDEO® Pad)</i>	\$549.00
250MH33		Whitmyer® Linx Headrest Mount with Quick Release	\$243.00
250MH34		Whitmyer® Linx-2 Headrest Mount with Quick Release	\$294.00
250MH35		Whitmyer® AXYS Headrest Mount with Quick Release	\$243.00
250MH17		OMIT Headrest Mount	NC
Push Handles \$			
250PH27		Height Adjustable Push Handles	\$60.00
Seat Frame Rails \$			
250SN16		HD Seat Rail <i>(For effective seat depth 18" to 24") (N/A with PRO ADVANCED)</i>	NC
250SN15		Seat Rail Long <i>(For effective seat depth 18" to 22")</i>	NC
Seat Cushion Type Seating Benefits \$			
250ST79		SEDEO® (See Section D) <i>(N/A with PRO ADVANCED)</i>	NC
250ST81		JAY® Syncra™ Cushion (See Section E)	NC
250ST80		JAY® (See Section F)	NC
250ST75		Omit	NC
C - Seat Size \$			
Seat Pan Width <i>(Odd sizes not available on SEDEO® cushions)</i>			
250R53		16"	NC
250R54		17"	NC
250R55		18"	NC
250R56		19"	NC
250R57		20"	NC
250R58		21" <i>(Not available with JAY® Syncra™ Cushions)</i>	NC
250R59		22"	NC
250R79		23" <i>(Not available for SEDEO, JAY® Syncra™ Cushions, or PRO ADVANCED)</i>	NC
250R80		24" <i>(Not available for SEDEO, JAY® Syncra™ Cushions, or PRO ADVANCED)</i>	NC
250R81		25" <i>(Not available for SEDEO, JAY® Syncra™ Cushions, or PRO ADVANCED)</i>	NC
250R82		26" <i>(Not available for SEDEO, JAY® Syncra™ Cushions, or PRO ADVANCED)</i>	NC

POWERED SEATING			
OPTION		DESCRIPTION	RETAIL
Seat Pan Depth (<i>Front of seat pan to front of back cushion – effective depth</i>)			
250RF46		18"	NC
250RF47		19"	NC
250RF48		20"	NC
250RF49		21" (<i>Not available with JAY® Syncra™ Cushions</i>)	NC
250RF50		22"	NC
250RF65		23" (<i>N/A with PRO ADVANCED</i>)	NC
250RF66		24" (<i>N/A with PRO ADVANCED</i>)	NC

D - SEDEO® Cushions (<i>N/A with PRO ADVANCED</i>)  P. 23 Seating Benefits \$			
250CU82		SEDEO® Cushion (<i>Sizes 16", 18", & 20" Wide</i>) (<i>Not available with 14" Depth</i>)	\$315.00
250CU83		SEDEO® Cushion (<i>Size 22" Wide</i>)	\$529.00

E - JAY® Syncra™ Cushion  P. 22 Seating Benefits \$			
250CU89		JAY® Syncra™ Cushion featuring JAY Flow® Fluid with Cryo™ Technology (<i>16" to 20" wide</i>)	\$1,002.00
250CU90		JAY® Syncra™ Cushion featuring JAY Flow® Fluid with Cryo™ Technology (<i>22" wide</i>)	\$1,108.00

JAY® Syncra™ Cushion Cover Upholstery

250YR73		Syncra™ Premium Leatherette	NC
250YR74		Syncra™ Stretch	NC

F - JAY® Seat Cushion  P. 24 Seating Benefits \$			
250CU59		JAY® Ion® (<i>Sizes 16" - 21" Wide</i>)	\$647.00
250CU60		JAY® Ion® (<i>Size 22" Wide</i>)	\$688.00
250CU61		JAY® Union® (<i>Sizes 16" - 21" Wide</i>)	\$753.00
250CU62		JAY® Union® (<i>Size 22" Wide</i>)	\$794.00
250CU64		JAY® Fusion® (<i>Sizes 16" - 21" Wide</i>)	\$1,030.00
250CU65		JAY® Fusion® (<i>Size 22" Wide</i>)	\$1,170.00
250CU95		JAY® J2® Plus, Fluid (<i>Sizes 24" - 26" Wide</i>)	\$1,297.00
250CU96		JAY® J2® Plus (<i>Sizes 24" - 26" Wide</i>)	\$983.00

F - JAY® Solid Seat Seating Benefits \$			
250CU84		JAY® Solid Seat	\$343.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION		DESCRIPTION	RETAIL
Knee and Thigh Swing-Away Support Pads Seating Benefits			\$
Lateral Knee Support Pads and Hardware (must pick pad and mount)			
250XY05		Lateral Knee Support Pad, Right (4" round)	\$243.00
250XY06		Lateral Knee Support Pad, Left (4" round)	\$243.00
250XX08		Lateral Swing-Away Bracket Knee - Right	\$320.00
250XX09		Lateral Swing-Away Bracket Knee - Left	\$320.00
Lateral Thigh Support Pads and Hardware (must pick pad and mount)			
250KF22		Lateral Thigh Support Pad Short - Right (5"x5")	\$144.00
250KF23		Lateral Thigh Support Pad Short - Left (5"x5")	\$144.00
250KF24		Lateral Thigh Support Pad Long - Right (5"x7")	\$144.00
250KF25		Lateral Thigh Support Pad Long - Left (5"x7")	\$144.00
250KF26		Lateral Thigh Removable Mount - Right	\$320.00
250KF27		Lateral Thigh Removable Mount - Left	\$320.00
JAY® Lateral Thigh Supports (Must pick pad, upholstery, and mount) Seating Benefits			\$
JAY® Lateral Thigh Support Pads			
250KF16		Lateral Thigh Support Pad Standard, Right (8.25" x 3.25")	\$144.00
250KF17		Lateral Thigh Support Pad Standard, Left (8.25" x 3.25")	\$144.00
250KF18		Lateral Thigh Support Pad Large, Right (12" x 3.25")	\$144.00
250KF19		Lateral Thigh Support Pad Large, Left (12" x 3.25")	\$144.00
JAY® Lateral Thigh Support Pad Upholstery			
250NU07		Syncra™ Premium Leatherette	NC
250NU08		Syncra™ Stretch	NC
JAY® Lateral Thigh Removable Mount			
250KF20		Lateral Thigh Removable Mount - Right	\$320.00
250KF21		Lateral Thigh Removable Mount - Left	\$320.00
Armrest Structure  P. 24			\$
250EG35		Reclining, Flip Back, Long (Height adjustable 9" to 15")	NC
250EG41		Reclining, Flip Back, Short (Best used with shorter seat depths) (Height adjustable 9" to 15")	NC
250EG36		Reclining, Flip Back with 2-point post locking arm (Height adjustable 11" to 15" x .625")	\$346.00
250EG51		Reclining, Flip Back, Dual Link Armrest (Required for PRO ADVANCED) (Height adjustable 8" to 15" including length and lateral adjustment)	NC
250EG32		Single Post, Flip Back (N/A with thigh support), (Height adjustable 10" to 15")	NC

POSITIONING SUPPORTS & ACCESSORIES		Seating Benefits	
OPTION		DESCRIPTION	RETAIL
Arm Pads  P. 24			\$
SEDEO® Armrest Pad			
250AB43		300 mm / 12", Integral Foam, Right	NC
250AE43		300 mm / 12", Integral Foam, Left	NC
250AB44		400 mm / 16", Integral Foam, Right	NC
250AE44		400 mm / 16", Integral Foam, Left	NC
Comfort Arm Pad (Foam)			
250AB46		2" Wide Right	\$111.00
250AE46		2" Wide Left	\$111.00
Comfort Arm Pad (Gel)			
250AB48		2" Wide Right	\$210.00
250AE48		2" Wide Left	\$210.00
Ottobock® Arm Troughs			
250AE06		Ottobock Arm Trough, Med (12"), Left	\$137.00
250AB06		Ottobock Arm Trough, Med (12"), Right	\$137.00
250AE07		Ottobock Arm Trough, Large (14"), Left	\$137.00
250AB07		Ottobock Arm Trough, Large (14"), Right	\$137.00
250OL01		Ottobock Hand Pad, Med (7"), Left	\$69.00
250OR01		Ottobock Hand Pad, Med (7"), Right	\$69.00
250OL02		Ottobock Hand Pad, Large (9"), Left	\$69.00
250OR02		Ottobock Hand Pad, Large (9"), Right	\$69.00
250AE78		Ottobock One Piece Arm Trough (18"), Left	\$195.00
250AB78		Ottobock One Piece Arm Trough (18"), Right	\$195.00
Side Guards (Pair) (Client Cost)			\$
250SG44		Side Guards Integral Foam (11"L x 5.5"H) (N/A with Reclining Armrests)	\$351.00
Fixed Centermount Legrest			\$
250H183		SEDEO® Fixed Centermount Legrest (3" to 17.5" Leg Length)	\$346.00
250H184		Fixed Centermount Legrest, Short (Includes Calf Pads & Single Footplate) (11.5" to 14.5" Leg Length)	NC
250H185		Fixed Centermount Legrest, Long (Includes Calf Pads & Single Footplate) (13.5" to 18.5" Leg Length)	NC

POSITIONING SUPPORTS & ACCESSORIES

OPTION		DESCRIPTION	RETAIL
Power Centermount Legrest  P. 24 (Client Cost)			\$
250H175		Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension <i>(Mechanical Articulation) (5" to 17.5" Leg Length)</i>	\$2,450.00
250H109		Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension <i>(Powered Articulation) (3" to 17.5" Leg Length) (Required for PRO ADVANCED)</i>	\$3,384.00
Power Centermount Legrest Accessories (Client Cost)			\$
Footplates			
250F127		Dual, Height, Width and Angle Adjustable, Standard (14"W x 10"D)	\$280.00
250F128		Dual, Height, Width and Angle Adjustable, Large (16"W x 12"D)	\$280.00
Calf Pad Options			
250DI34		Calf Pads, Pair, Standard	NC
250DI35		Extra Large Pair Calf Pads for Centermount Legrests. Fully adjustable <i>(N/A w/ Std footplate)</i>	\$351.00
SEDEO® PRO™ Swing-Away Legrest <i>(Pair)</i> (250H191 is Client Cost)			\$
250H190		HD Legrests, Swing-Away, 70°, SEDEO® PRO™ <i>(14" to 21" Leg Length)</i>	NC
250H199		HD Legrests, Manual Elevating, Swing-Away, SEDEO® PRO™ <i>(14".5 to 21.5" Leg Length) (Includes Padded Calf Pads)</i>	\$596.00
250H191		HD Legrests, Power Elevating, Swing-Away, SEDEO® PRO™ <i>(15.5" to 22.5" Leg Length) (Includes Padded Calf Pads)</i>	\$3,384.00
250DI76		Omit Legrests	NC
Swing-Away Legrest Footplates			
250F369		Angle Adjustable Flip-Up Footplates	\$280.00
Swing-Away Legrest Calf Supports			
250DI42		Calf Support, Depth Adjustable, Padded <i>(6.5" to 8" Wide, Depending On Chair Width, x 9.5" Tall)</i>	\$504.00
Swing-Away Legrest Options			
250DI44		Heel Loops	NC
250KG4		Knee Pads, Integral Foam <i>(2.5" Wide x 4.5" Tall)</i>	\$278.00
SEDEO® Lower Extremity Swing-Away Support			
250VD01		Lower Extremity Support, Padded, Swing-Away, SEDEO®, Right	\$504.00
250VD02		Lower Extremity Support, Padded, Swing-Away, SEDEO®, Left	\$504.00
Pelvic Positioning Belts  P. 24			
250A837		Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
250A839		Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"	NC
250A840		Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$101.00
250A845		Bodypoint - Evoflex Pelvic Stabilizer Belt	\$213.00
250A846		Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	\$267.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION		DESCRIPTION	RETAIL
Accessories			\$
Ventilator Option			
250VT12		Vent Tray Fixed <i>(Pulmonetic LTV Models)</i>	\$600.00
O2 Tank Option			
250A718		O2 Tank Holder Mounted to Backrest <i>(D and E Styles)</i>	\$385.00
Wheelchair Tray Table			
2507A09		Tray, Swing-Away	\$375.00
Backpack		(Client Cost)	
250A1		Black Utility	\$110.00
Medial Thigh Support (N/A with Centermount Legrests)			
2508A80		Medial Thigh Support	\$267.00
2508A81		Pommel Swing-Away Hardware	\$384.00
Misc. Accessories		(Client Cost)	
2508A59		Transfer Bar <i>(Pair)</i>	\$486.00
2508A52		Ram Mount Phone	\$76.00
2508A53		Ram Mount Tablet	\$303.00
2508A54		Ram Mount Cup Holder	\$76.00
2508A58		Rear View Mirror	\$76.00
250A611		Crutch Holder <i>(N/A w/Power Recline)</i>	\$171.00
2508A73		Bag for Armrest	\$76.00
2508A75		Reflectors <i>(Backrest shroud mounted)</i>	\$113.00
2509A10		Utility Hooks	\$29.00

JOYSTICK / ELECTRONICS

OPTION		DESCRIPTION	RETAIL
R-net Joystick (Remote)  P. 25			\$
250J121		Standard LED Joystick (w/o Screen Display)	NC
250J116		Color Joystick (2.3" Color LCD Screen Display, CJS1M1)	NC
250J118		High Resolution Color Joystick,(CJS2M2+BT) 3.5" Color LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereos Jack, IR and dual (IOS/Android) Bluetooth control (Client Cost)	\$1,145.00
250J119		Omit Remote Control (Required with arm mounted specialty controls)	NC
Joystick Mounting			\$
250JM34		Fixed on Armrest	NC
250JM35		Parallel Swing-Away	\$440.00
250JM36		Link-It Mid Line Mount on Armrest (Swing-Away)	\$440.00
Joystick Mount (Location)			\$
250HM16		Right Hand Armrest	NC
250HM17		Left Hand Armrest	NC
Display Language Selection R-net			\$
250L1		English	NC
250L3		French	NC
250L7		Spanish	NC
Joystick Handle (Optional)  P. 25			\$
250JH13		Ball Handle	\$108.00
250JH15		T-Handle	\$108.00
250JH24		Bodypoint – Rubber Dome Joystick Handle	\$213.00
250JH25		Bodypoint – U-Shaped Joystick Handle (4 inch)	\$213.00
250JH26		Bodypoint – Grooved Mushroom Joystick	\$213.00
Electronic Programmers  P. 26			\$
250PR10		R-net PC Programmer w/ Dongle	\$696.00
R-net Power Modules  P. 26			\$
250ER18		Power Module, R-net, 90A, Non-Expandable	NC
250ER17		Power Module, R-net, 90A, Expandable Required for specialty controls or 3+ actuators	\$100.00
250EX15		Expandable 90A Controller - R-net	NC
250EX16		Harness for 90A Expandable Controller	NC

JOYSTICK / ELECTRONICS

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

Seating Control \$

Powered Seating Control

250E040	Omit Thru-Drive <i>(Requires Switch Control)</i>	NC
250E037	Thru-Drive Control for 1 Actuator	\$1,422.00
250E038	Thru-Drive Control for 2+ Actuators	\$2,860.00
250E032	CxSM, Thru-Drive Control for 2+ Actuators <i>(Required for PRO ADVANCED)</i>	\$2,860.00

Switch Control *(Available with Omit Thru-Drive only) (1 actuator for Tilt or Elevate only)*

250E033	Dual Toggle for 1 Actuator	NC
250E034	2.5" Single Button Switch for 1 Actuator <i>(Button switches are not mounted)</i>	NC
250E041	Ctrl+5 Switch Box for 2 Actuators or More	NC

Additional Power Seating Control *(Available with Thru-Drive only)*

250QG11	Ctrl+5 Switch Box with Free Assignable Functions	\$1,062.00
250QG12	Ctrl+5 Switch Box with Free Assignable Functions <i>(Multi-power only)</i>	NC

Dual Toggle Mounting Location

250TN1	Toggle Mount, Right	NC
250TN2	Toggle Mount, Left	NC

Ctrl+5 Mounting Location

250EY19	Right Armrest	NC
250EY20	Left Armrest	NC
250EY21	Back Cane, Right	NC
250EY24	Back Cane, Left	NC

Ctrl+5 Switch Type P. 26

250EY22	Push Buttons	NC
250EY23	Toggle Switches	\$351.00

Attendant Controls P. 25 (Prior Approval Required) \$

2500A05	Dual Control Module for Attendant	\$1,207.00
2500A06	R-net Color LCD (CJSM1), Attendant Control Only	\$1,207.00

Attendant Mounting Location \$

250DN5	Right Back Cane	NC
250DN6	Left Back Cane	NC

JOYSTICK / ELECTRONICS

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

R-net Programming \$

Assignable Buttons

Assignable Button feature allows for personalization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

250PM12	Include Assignable Buttons	NC
---------	----------------------------	----

Power Centermount Footplate Extension Option (Client Cost)

250IE03	Include 7.5" Powered Footplate Extension for Transfer Assist (Electronics programmed to allow footplates to be extended with or without legrest being raised.) Only available with Powered Articulation.	\$583.00
---------	--	----------

R-net Omni2 Options P. 26 For Expandable Controls, choose all that apply. \$

2500N14	R-net Omni2 Only Input Control Module and 3.5" LCD Display w/ Built-In Infrared & Bluetooth Choose Mount Below	\$2,349.00
2500N15	R-net Omni2 Bundled Input Control Module and 3.5" LCD Display w/ Built-In Infrared & Bluetooth Bundled Omni2 Available w/ Switch-It Input Devices	NC

Omni2 Mount If used w/ Standard Joystick, it must be mounted opposite side. \$

250DS14	Link-It with Parallel Swing-Away Omni2 Mount w/ Rear Mounted Interface, Right	\$688.00
250DS15	Link-It with Parallel Swing-Away Omni2 Mount w/ Rear Mounted Interface, Left	\$688.00

PRO ADVANCED Programming \$

Biometric Repositioning

250PM03	Include Biometric Repositioning (Requires Memory Position #1, #2, #3 below)	NC
---------	---	----

Memory Seating

"Memory Seating" feature allows for up to six custom seating positions that are accessed either through the joystick, Ctrl+5, joystick shortcut keys or programmable jacks on Omni2 and joystick. Choose desired preset settings below. Default positions are dealer adjustable.

250PM06	Include #1 Memory Seating (Seated) (Required)	NC
250PM07	Include #2 Memory Seating (Supine)	NC
250PM08	Include #3 Memory Seating (Pressure Relief)	NC
250PM09	Include #4 Memory Seating (Feet Above Heart)	NC

PRAR Pressure Relief Alert and Reminders

PRAR Pressure Relief Alert and Record feature provides programmable alerts to be set into the R-net electronics that will prompt the user when it is time to perform a pressure relief. The R-net system will record how often the pressure relief is completed or ignored.

250PM14	Include PRAR (Minimum requirement Memory Position #1 and #3)	NC
---------	--	----

SPECIALTY INPUT DEVICES (AADL W100 Grant)

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

Hand Control Packages  P. 26 *Requires Standard Joystick Omit / Only one hand control device can be ordered.* (AADL W100 Grant) \$

Proportional Input Devices

250EU14	Switch-It MicroGuide	\$4,574.00
250EU15	Switch-It MicroPilot	\$4,574.00
250EU17	Switch-It VersaGuide	\$2,519.00
250EU18	Switch-It VersaGuide EZ	\$2,978.00

Mounting Position (Required)

250EM04	Link-It to the Backrest <i>(Not Available w/ Omit Back Cushion Type)</i>	\$1,318.00
250EM05	Link-It to the Fixed Joystick Mount	\$410.00
250EM06	Link-It to the Parallel Swing-Away Joystick Mount	\$410.00
250EM07	Link-It to the Midline Swing-Away	\$1,458.00
250EM08	Omit Mount	NC

Mounting Platforms (Optional) *Not available w/ Link-It to Backrest*

250MG16	SunDisk Hand Tray <i>(Only available w/ MicroSeries Joysticks)</i>	\$410.00
250MG15	Bullet Hand Tray	\$410.00

Mounting Location (Required)

250IM09	Right	NC
250IM10	Left	NC

Switch Options

—	See Switch Selections <i>(Maximum 2)</i>	—
---	--	---

Head Control Packages  P. 27 *(Select headrest mount on next page)* (AADL W100 Grant) \$

Only one head control device can be ordered.

Proportional Input Devices

250EU25	Switch-It DualPro® Proportional Head Array - Facial Pads	\$6,953.00
250EU26	Switch-It DualPro® Proportional Head Array - Spot Pads	\$6,953.00

Switched Input Devices

250EU28	Switch-It 3-Switch Head Array - Short Curved Wings	\$6,063.00
250EU30	Switch-It 3-Switch Head Array - Long Curved Wings	\$6,063.00

Switch Options

—	See Switch Selections <i>(Maximum 2)</i>	—
---	--	---

SPECIALTY INPUT DEVICES

OPTION		DESCRIPTION	RETAIL
Sip 'N Puff (AADL W100 Grant)			\$
250ID73		Switch-It Sip 'N Puff Combo Head Array with Link-It Hardware <i>(Requires headrest mount)</i>	\$6,933.00
250ID74		Switch-It OptiStop with Sip 'N Puff	\$2,700.00
Headrest Mounts <i>(Applies to Head/Sip 'N Puff control packages)</i> (AADL W100 Grant)			\$
250IJ13		Whitmyer® Linx Headrest Mount w/ Quick Release	\$243.00
250IJ14		Whitmyer® Linx-2 Headrest Mount w/ Quick Release	\$294.00
250IJ15		Whitmyer® AXYS Headrest Mount with Quick Release	\$243.00
Chin Control Packages  P. 27 (AADL W100 Grant)			\$
Proportional Input Devices			
250EU34		Switch-It MicroGuide <i>(Switches mounted to the headrest only.)</i>	\$4,574.00
250EU37		Switch-It VersaGuide EZ	\$2,978.00
250EU36		Switch-It VersaGuide	\$2,519.00
250EU37		Switch-It VersaGuide EZ	\$2,978.00
Mounting Position			
250EM09		Link-It To Chin Harness, Small (150mm / 6")	\$1,734.00
250EM10		Link-It To Chin Harness, Large (220mm / 8")	\$1734.00
250EM11		Link-It To Headrest (Swing-Away)	\$1,734.00
250EM53		Link-It To Manual Chin Boom - Right Side	\$5,267.00
250EM54		Link-It To Manual Chin Boom - Left Side	\$5,267.00
250EM16		Omit Mount	NC
Input Device Mounting Location			
250IM13		Right	NC
250IM14		Left	NC
Switch Options			
—		See Switch Selections	—
Safety Devices (AADL W100 Grant)			\$
250SA06		Switch-It Remote Stop Switch - XLR Connector	\$1,298.00
Switched Driving Packages (AADK W100 Grant)			\$
250SW44		Switch-It™ Cool Cube with DB9 connector for proximity/proportional switch driving. <i>(Pro Spot proportional switches are required if proportional driving is desired.)</i> <i>(For non-proportional driving, single switches can be used.)</i>	\$3,891.00

SPECIALTY INPUT DEVICES

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

Switches P. 27

Single switches can be programmed to activate any function of the wheelchair including: On/Off, Mode, Profile, Speeds, Seating, Memory Positions, Horn, Lights & Indicators, IR functions, Bluetooth mouse clicks, Environmental Controls when plugged into assignable port.

Switches can connect via: Joysticks (Standard, Micro series, Versa series) Omni2, SJX6 Hub, head controls

Switch ordering instructions:

1. Select switch type and quantity in Section "A" using columns below. (Max 6 switches. See limitations below.)
2. Select switch mounting position from Section "B" using corresponding switch column.
3. Select switch mounting location from Section "C" using corresponding switch column.

Section A – Switch Selections

Switch / Quantity

Single Function Switches - Mechanical Switches	Switch 1	Switch 2	Switch 3	Switch 4	Switch 5	Switch 6	\$
2.5" Single Switch - Black							\$108.00
2.5" Single Switch - Blue							\$108.00
2.5" Single Switch - Green							\$108.00
2.5" Single Switch - Red							\$108.00
2.5" Single Switch - Yellow							\$108.00
Egg Switch - Black							\$108.00
Egg Switch - Blue							\$108.00
Egg Switch - Green							\$108.00
Egg Switch - Red							\$108.00
Egg Switch - Yellow							\$108.00
MicroLite <i>(Includes color labels)</i>							\$113.00
Mini Button 1 - Blue <i>(Velcro)</i>							\$303.00
Mini Button 1 - Green <i>(Velcro)</i>							\$303.00
Mini Button 1 - Red <i>(Velcro)</i>							\$303.00
1" Single Switch - Black							\$382.00
1" Single Switch - Blue							\$382.00
1" Single Switch - Green							\$382.00
1" Single Switch - Orange							\$382.00
1" Single Switch - Red							\$382.00
1" Single Switch - Yellow							\$382.00
Electronic Switches							
Flex Switch (Headrest mount only)							\$534.00
ProSpot - Proportional Switches							
ProSpot Proportional Switch - Red							\$534.00
ProSpot Proportional Switch - Green							\$534.00
ProSpot Proportional Switch - Blue							\$534.00
ProSpot Proportional Switch - Yellow							\$534.00
ProSpot Proportional Switch - Orange							\$534.00

Section B – Switch Mounting Position

Link-It to the Backrest <i>(Up to 1 per side)</i>							\$1,388.00
Link-It to the Headrest - Gooseneck & Rigid <i>(Up to 2 per side)</i>							\$1,734.00
Link-It to the Fixed Joystick Mount <i>(Up to 2 per side)</i>							\$410.00
Link-It to the Parallel Swing-Away Joystick Mount <i>(Up to 2 per side)</i>							\$410.00
Link-It to the Seat Rail <i>(Up to 2 per side)</i> <i>(Not avail w/knee & thigh supports)</i>							\$548.00
Omit Mount <i>(includes Velcro)</i>							NC

Section C – Switch Mounting Location

Right side of chair							NC
Left side of chair							NC

Joystick Integrated Switches *(Pair)*

250J003	–	Proximity Switches on Flex Mount, Pair <i>(Used with VersaGuide/VersaGuide EZ)</i>	\$1,382.00
---------	---	--	------------

JAY® SYNCRA™ SEATING MEASURING GUIDE

JAY® SYNCRA™ BACK



Backrest Width	16"	18"	20"	22"
Max User Trunk Width	13.5"	15.5"	17.5"	19.5"


JAY® Swing-Away Lateral Thoracic Adjustment Range

Syncra Back Width	Width Range	Height range from bottom of shell to seat pan
16"	10" to 24"	12" to 19.5"
18"	12" to 26"	12.5" to 21.5"
20"	14" to 28"	12.5" to 22"
22"	16" to 30"	13.5" to 23"

Syncra™ Fixed Lateral Support Accessory Range

Syncra Back Width	Width Range	Height range from bottom of shell to seat pan	Contour Depth Range
16"	16.4" to 16.8"	6" to 8"	3" to 5"
18"	18.4" to 18.8"	6" to 8"	3" to 5"
20"	20.4" to 20.8"	6" to 8"	3" to 5"
22"	22.4" to 22.8"	6" to 8"	3" to 5"

Syncra™ Back Shell Height Diagram



JAY® SYNCRA™ CUSHION



JAY® Lateral Thigh Support Adjustment Range

Syncra Cushion Width	Width Range
16"	8.5" to 16.0"
17"	9.5" to 17.0"
18"	10.5" to 18.0"
19"	11.5" to 19.0"
20"	12.5" to 20.0"
22"	14.5" to 22.0"

built4me

Q700 M with SEDEO® PRO™

\$

Base Price

Total of additional equipment

Packaging and transport

Total Price

Date

Signature

Month

Day

Year



Customer Service: (800) 263-3390
Fax: (800) 561-5834
E-mail: cscanada@sunmed.com
www.sunrisemedical.ca
Specifications are subject to change without notice.

© 2024 Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2
Concord, ON
L4K 4P8

MK-129763-AADL Rev. 10