

JAY[®]

UNION BACK



☰ **Comfort and Support all in One**

JAY UNION BACK

≡ The JAY Union Back is the complementary addition to any manual, positioning and power wheelchairs.

Innovative Hardware Options:

The Jay Union back comes with the standard JAY hardware. It can also be ordered with the 4 Point hardware to support greater force to the back.

Multiple Foams:

The Jay Union back is made of 3 foam layers for unparalleled comfort

- Structural foam
- Dispersion foam
- Comfort foam

This combination maximizes support and comfortably contours to the user's anatomy.

Built in Lateral Support:

The Jay Union shell is molded to provide natural lateral support, forgiving contour and excellent lateral trunk stability.

Whitmyer Union Headrest Package:

Our newest concept in contemporary head support systems, the Whitmyer Union headrest package is available as an option to provide infinite range of support surface placement.

Lateral Support Option:

To increase trunk control, optional lateral supports can be added.

Extension:

The Jay Union can also be ordered as a tall back with the addition of an extension. This extension gives further pelvic support. The extension is adjustable in length and can be aligned to accommodate a kyphosis or a lordosis.



Specifications

JAY Union Back	
Weight (without brackets and hardware)	4 lbs.
Product Width	16" – 22"
Product Depth	6"
Back Height	Shell: 18"; Regular: 20.5"; Tall: 23" – 24"
Weight Capacity	300 lbs.
Cover Options	Soft 4-way Stretch Lycra

JAY®

