



EIRGK02



### Rider Information

**Name:** \_\_\_\_\_

**Phone:** \_\_\_\_\_

**Email:** \_\_\_\_\_

**Height:** \_\_\_\_\_ **Weight:** \_\_\_\_\_

**Current Chair:** \_\_\_\_\_ **Current Cushion:** \_\_\_\_\_

**Ability/Diagnosis:** \_\_\_\_\_

### RGK Contact Authorization:

RGK would like to send you exclusive information, news and promotional offers on RGK.

I would like to be contacted by RGK:

**Signature:** \_\_\_\_\_ **Date:** \_\_\_\_\_

\*We do not pass your personal information onto any third party.

### Account Information

**Account Name :** \_\_\_\_\_

**Account #:** \_\_\_\_\_ **PO#:** \_\_\_\_\_

**Contact Name :** \_\_\_\_\_

**Phone:** \_\_\_\_\_

**Email:** \_\_\_\_\_

**Ship to:** **Account** **Rider**

**Name:** \_\_\_\_\_

**Address:** \_\_\_\_\_


**RGK Genius:** \_\_\_\_\_

## Frame Styles

702XZ010 **Octane FX** (Max User Weight 110kg/242lbs.)  
702XA070 **Open design**

702XW009 **Quote Only** (No CAD for Approval)  
702XW011 **Order Only** (No CAD for Approval)  
702XW012 **Quote CAD drawing for approval** (1st CAD is NCO)  
702XW010 **Order CAD drawing for approval** (1st CAD is NCO)

CAD approved by: 702XW013 Customer 702XW014 RGK Genius



Open Design

\$13,690  
STD

STD

STD

\$550

\$300

## Measurements

**Please Note:** All RGK chairs are hand built and have a manufacturing tolerance of +/- 5mm (.25"). All measurements can be done in 5mm (.25") increments.

**Please select either Millimeters or Inches.**

**Measurements:** ☐ **Millimeters (mm.)** ☐ **Inches (in.)**

**702XY010 Rear Seat Width** Measured from inside of side guard to inside of opposite side guard. Or outside seat tube to outside of opposite seat tube near backrest. (Min. 250mm/9.84" and Max. 460mm/18.11")

mm. in. mm. in.  
Current Chair:   New Chair:

**702XY020 Seat Depth/Seat Length** Measured from front of back post to the front edge of seat sling. Frame bend will start .5" in front of seat sling unless Extra Frame Length requested. (Min. 250mm/9.84" and Max. 460mm/18.11")

Current Chair:   New Chair:

**702XY030 Extra Frame Length** Measured from front edge of seat sling to beginning of Front frame bend. (Max. 460mm/18.11" Seat length + Extra seat length)



Current Chair:   New Chair:

**702XY050 Front Frame Height** Measured from floor to apex of seat tube. (Min. 480mm/18.89" and Max. 550mm/21.65")

Current Chair:   New Chair:

**702XY060 Rear Frame Height** Measured from floor to top of seat tube at back edge of seat sling. (Min. 360mm/14.17" and Max. 550mm/21.65")

Current Chair:   New Chair:



# Design Template

**702XY070 Occupied Frame Length/Frame Length** Measured from the front of the backrest tube to the front edge of the footrest. (Min. 390mm/15.35" and Max. 625mm/24.6")

Current Chair:   New Chair:

**702XY072 Caster Stem Position** Measured from the front edge of the footrest to the centerline of the caster barrel. **702XY073** RGK to Optimize:

Current Chair:   New Chair:

**702XA020 Ergonomic Seat** Measured from front of back post to beginning of upward bend of seat tube. \$121

Current Chair:   New Chair:

**702XY080 Front Seat Taper/Seat Taper Design** Measured from inside of front frame tube at the front frame bend to inside of opposite front frame tube at the front frame bend. **702XY081** RGK to Optimize.

Current Chair:   New Chair:

**702XY115 Front Footrest Height** Measured from ground up to the top of the Front footrest tube or cover.

Current Chair:   New Chair:

**702XY120 Rear Footrest Height** Measured from ground up to the top of the Rear footrest tube or cover. Standard is 20mm less than Front Footrest Height.

Current Chair:   New Chair:

**702XY125 Footrest Length** Measured from front edge of footrest to back edge of footrest. (160mm and over requires heavy duty upgrade)


Current Chair:   New Chair:

**702XB200 High Mount Footrest/High Level Footrest** Measured from ground up to the top of the Front footrest tube or cover. \$201

Current Chair:   New Chair:

**702XY330 Stem Bolt/Caster Pin Width** Measured from centerline of caster pin to the centerline of the opposite caster pin. **702XY331** RGK to Optimize:

Current Chair:   New Chair:



# Design Template

702XY231 **Backrest Height** Measured from top of seat tube at rear of frame to top of back-post.

Current Chair:   New Chair:

702XY360 **Backrest Angle** Measured as the angle from top front of the backrest tube to the floor. (Max. -10° to +10°) -5° -3° 0° 3° 5°

Current Chair:  New Chair:

702XG110 **Inward Taper** Measured from inside backrest tube to inside of opposite backrest tube at the top of the backrest. \$146

Current Chair:   New Chair:

702XG100 **Ergonomic Backrest** Measured from top of seat tube at rear of frame to the start of the rearward bend of back-post. \$146

Current Chair:   New Chair:

702XG120 **Outward Flare** Measured from inside backrest tube to inside of opposite backrest tube at the top of the backrest. \$146

Current Chair:   New Chair:

702XY390 **Camber Angle** Place a straight edge across the outside of rear wheel and use an angle finder to show the degree of camber.






Current Chair:	0°	New Chair:	702XY396 0°
	1°		702XY397 1°
	2°		702XY398 2°
	3°		702XY399 3°
	5°		702XY400 5°

702XH165 **Center of Gravity** Measured from front of back post at the seat tube to the centerline of the rear axle.

Current Chair:   New Chair:

702XY460 **Rear Wheel Gap** Measured from outer most point of chair (seat back post, side guard or armrest) to inside of rear tire.

Current Chair:	.6"(15)	New Chair:	702XY464 .6"(15)
	.78"(20)		702XY465 .78"(20)
	1"(25)		702XY466 1"(25)



## Footrest Design

702XY100	<b>Round Design</b>	STD
702XY105	<b>Square Design</b>	NCO

## Footrest Styles

702XB080	<b>Tubular footrest welded into frame aluminum plate</b> (No height adjustment)	STD
702XB110	<b>Tubular footrest welded into frame with carbon fiber footplate</b> (No height adjustment)	\$326
702XB050	<b>Lightweight 3k Carbon Fiber Footrest cover</b>	STD
702XB155	<b>Lightweight Spinergy CX6 Cast Carbon Fiber Footrest cover</b>	NCO
702XB120	<b>Clamp-on aluminum footrest plate</b> (Angle and height adjustment)	NCO
702XB165	<b>Rubber grip tape</b> (3M Safety walk. Not available with 702XB030)	\$86

## Front Wheels

702XC060	<b>Aluminum Center 4"x1.3" (Black hub only)</b>	STD
702XC110	<b>Aluminum Center 4"x1" (Silver hub only)</b>	\$232
702XC135	<b>Aluminum Center 4"x1.25"</b>	\$232
702XC1131	Sliver hub & black tire	
702XC1132	Black hub & black tire	
702XC1133	Purple hub & black tire	
702XC1134	Green hub & black tire	
702XC1135	Red hub & black tire	
702XC1136	Blue hub & black tire	
702XC155	<b>4" x 1.25" Spinergy CX6 Cast Carbon fiber front wheel</b>	\$235
702XC156	<b>Cerakoted Spinergy CX6 Cast Carbon fiber front wheel</b> (Hub only. See Cerakot color options on page 6)	\$121
702XC070	<b>Soft roll 5"x1.3"</b>	NCO
702XC100	<b>Pneumatic Aluminum Center 5"x1"</b>	\$248
702XC120	<b>Aluminum Center 5"x1"</b>	\$232
702XC140	<b>Aluminum Center 5"x1.4"</b> (Sliver hub & black tire)	\$232
702XC145	<b>Aluminum Center 6"x1.4"</b> (Sliver hub & black tire)	\$232

## Forks

702XC190	<b>Low Ride Fork</b>	STD
702XC192	<b>High Trail fork</b>	NCO
702XC210	<b>Froglegs suspension fork</b>	\$633
702XC211	Black fork	STD

## Wheel Lock

702XD020	<b>Underseat Scissor</b>	NCO
702XD030	<b>Compact Scissor</b>	\$266
702XD035	<b>Compact Scissor with Direct Mount</b>	\$352
702XD050	<b>Push To Lock</b>	NCO
702XD085	<b>Eagle</b> (Indicate color)	\$468
702XD086	<b>Eagle with Direct Mount</b> (Indicate color)	\$508
702XD088	Black	
702XD089	Purple	
702XD090	Red	
702XD091	Blue	
702XD010	<b>Omit Wheel Lock</b>	NCO

## Side Guards

702XE080	<b>Lift out aluminum</b>	\$532
702XE130	<b>Fold flat aluminum</b>	\$670
702XE100	<b>Lift Out 3k Carbon Fiber</b>	\$1,015
702XE106	<b>Lift Out Spinergy CX6 Cast Carbon fiber</b>	\$1,015
702XE150	<b>Fold flat 3k carbon fiber</b>	\$1,052
702XE155	<b>Fold flat Spinergy CX6 Cast Carbon fiber</b>	\$1,052
702XE120	<b>Lift out ABS plastic</b>	\$487
702XE170	<b>Fold flat ABS plastic</b>	\$830
702XEM05	25mm above tire (Follow tire)	NCO
702XEM06	Individual above tire (Follow tire - max. 25mm)	NCO
702XEM02	Level with tire	STD
702XEM07	Side guard length	STD
702XEM08	Side guard height	STD
702XE090	<b>Lift out aluminum with fender</b>	\$579
702XE140	<b>Fold flat aluminum with fender</b>	\$707
702XE110	<b>Lift Out 3k Carbon Fiber with fender</b>	\$1,132
702XE211	<b>Lift Out Spinergy CX6 Cast Carbon fiber with fender</b>	\$1,132
702XEM03	10mm above tire	NCO
702XEM01	Individual above tire (Follow tire - max. 25mm)	NCO
702XEM09	Center of tire	STD
702XEM10	Outside of tire	NCO

## Armrest

702XF010	<b>Lift out armrest with plastic sideguard</b>	<b>Height</b> (Top of seatrail to top of armrest)	\$592
702XF020	<b>Lift out armrest with aluminum sideguard</b>	<b>Height</b> (Top of seatrail to top of armrest)	\$467

## Backrest Height

702XY230	<b>Fixed height</b>	STD
702XY245	<b>Adjustable height</b> (Maximum height 15.75(400) Range: +/- 25mm)	\$248

## Backrest Design

702XG060	<b>Folding backrest with Q Lock</b>	STD
----------	-------------------------------------	-----

## Push Handles

702XG070	<b>Screw in push handles</b> (Not available with 702XY245)	\$225
702XG080	<b>Flip down push handles</b> (Not available with 702XY245)	\$441
702XG085	<b>Flip down push handles height adjustable</b> (Not available with 702XY245)	\$720
702XG097	<b>Integrated push handles</b>	\$676

## Center of Gravity

702XH110	<b>Square back adjustable of gravity</b> (Range: 25mm - 150mm)	STD
702XH120	<b>Round back adjustable of gravity</b> (Range: 60mm - 150mm)	NCO
702XH090	<b>V-Style fixed center of gravity</b> (Range: 25mm - 125mm)	\$727

## Axle Pin Options

702XK227	<b>Standard axle pin</b>	STD
702XK226	<b>Titanium 100mm quick release axle pin</b>	\$453
702XK220	<b>Titanium 130mm quick release axle pin</b>	\$160
702XK240	<b>Quad release axle pin</b>	\$241

## Anti-Tip

702XJ010	<b>Swing-away anti tip LEFT side</b> (Only available with square back adjustable center of gravity. Not available with 702XH140)	\$167
702XJ020	<b>Swing-away anti tip RIGHT side</b> (Only available with square back adjustable center of gravity. Not available with 702XH140)	\$167
702XJ100	<b>Single sport anti-tip</b> (Centered mounted)	\$504

## Frame and Footrest Finishes

## Non painted finishes

702XN030	Bead blast frame	STD
702XN050	Polished frame	\$1,332

## Painted finishes

702XN500	<b>CeraKote frame</b>	\$534
702XN600	<b>CeraKote components</b> (Component color will match frame if not otherwise specified. Components include: Forks, Wheel Lock clamps, and Camber Tube clamps.)	\$534
702XN800	<b>CeraKote Spinergy wheels</b> (See Rear Wheel Type on page 8. Hubs and Rims are CeraKoted on Rear Wheels, Hubs on Front Wheels.)	\$655

Frame	Components	Spinergy	Colors	Frame	Components	Spinergy	Colors
702XN501	702XN601	702XN801	<b>Midnight Black</b>	702XN509	702XN609	702XN809	<b>Rose Gold</b>
702XN502	702XN602	702XN802	<b>Java Bronze</b>	702XN510	702XN610	702XN810	<b>Retro Blue</b>
702XN503	702XN603	702XN803	<b>Pure Titanium</b>	702XN511	702XN611	702XN811	<b>Hot Pink</b>
702XN504	702XN604	702XN804	<b>Fire Red</b>	702XN512	702XN612	702XN812	<b>Wicked Green</b>
702XN505	702XN605	702XN805	<b>Prime Blue</b>	702XN513	702XN613	702XN813	<b>Bright White</b>
702XN506	702XN606	702XN806	<b>Sahara Gold</b>	702XN514	702XN614	702XN814	<b>Tangerine Orange</b>
702XN507	702XN607	702XN807	<b>Black Cherry</b>	702XN515	702XN615	702XN815	<b>Sunbeam Yellow</b>
702XN508	702XN608	702XN808	<b>Fresh Green</b>	702XN516	702XN616	702XN816	<b>Matte Navy</b>

## Accent Color

702XN251	<b>Black</b>	STD
702XN253	<b>Red</b>	NCO
702XN255	<b>Purple</b>	NCO
702XN252	<b>Silver</b>	NCO
702XN254	<b>Blue</b>	NCO
702XN250	<b>Individual Accent Color</b>	\$318

## Seat Upholstery

702XP040	<b>Seat Sling Adjustable Strap</b>	STD
702XR050	<b>Carbon Fiber Seat Plate</b>	\$480

## Back Upholstery

702XP050	<b>Airtech Adjustable Strap</b>	NCO
702XP060	Backrest upgrade to black full leather	\$500
702XP070	Backrest upgrade to black leather corners	\$121
702XP110	Backrest zip pocket	\$127



# Design Template

<b>RGK Backrest Logo Color</b>			
702XP557	White (Standard if nothing selected)		STD
702XP558	Red		STD
702XP559	Black		NCO
702XP560	Blue		NCO
702XP561	Purple		NCO
702XP562	Individual	Pantone color number	NCO
<b>Rear Wheel Size</b>			
702XY420	24" (540) Rear wheel size	702XY430 25" (559) Rear wheel size	STD
<b>Rear Wheel Type</b>			
702XK015	<b>Proton</b> (Not available with Spinergy anodized color options.)		STD
702XK320	<b>Spinergy XLX</b> (12 spoke, cross radial spokes.)		\$1,067
702XK050	<b>Spinergy LX</b> (12 spoke)		\$1,000
702XK060	<b>Spinergy SPOX</b> (18 spoke. Not available with Spinergy anodized color options or spoke colors, purple, orange, pink and green.)		\$665
702XK070	<b>Spinergy BLXL</b> (18 bladed spoke.)		\$1,067
702XK150	<b>Spinergy CLX</b> (18 bladed spokes, carbon fiber rim. Spinergy color option available on hub only. Not available with silver hub)		\$3,139
702XK330	<b>RGK EXCLUSIVE Spinergy CX6 Cast Carbon Fiber</b> (Cast Carbon Fiber Hub, Matte Black Carbon Fiber Rim and Black Bladed Spokes with 2 White Bladed Spokes Around Valve Stem. Not available with Spinergy anodized color options.)		\$3,592
<b>Spinergy spoke colors</b>			
702XK650	Black		STD
702XK651	Red		NCO
702XK652	Purple		NCO
702XK653	White		NCO
702XK654	Blue		NCO
702XK655	Orange		NCO
702XK656	Yellow		NCO
702XK657	Pink		NCO
702XK658	Green		NCO
702XK180	Individual spoke design		\$213
<b>Spinergy anodized colors</b> (Rims and Hubs are anodized not spokes, on Spinergy BLXL (XK075), Hubs only on 702XK150)			
702XK190	Black		STD
702XK202	Purple (LX, XLX & BLXL only)		\$409
702XK203	Blue (LX, XLX & BLXL only)		\$409
702XK201	Red (LX, XLX & BLXL only)		\$409
702XK204	Gold (LX, XLX & BLXL only)		\$409
702XK205	Individual hub/rims		\$568
<b>Supply Without Rear Wheels</b>			
702XK040	<b>Omit Rear Wheels</b>		NCO
	702XK041 Wheel Type:	702XK042 Tire Type:	
<b>Handrims Style</b>			
702XY485	<b>Tab Length Wide</b>		STD
702XY480	<b>Tab Length Narrow</b>		NCO
	702XL090 <b>Remove excess mounting tabs</b> (Recommended for small wheel gaps XY460)		\$146



# Design Template

702XL010	Aluminum Hard Anodized 6 Tab				STD
702XL160	Elipse Ergonomic 6 Tab				\$631
702XL020	QGrip 6 Tab				\$400
702XL040	Titanium 6 Tab				\$446
702XL035	Natural-Fit LT				\$743
702XL150	Carbolife Gekko				\$795
702XL170	Carbolife Curve L				\$695
702XL180	Carbolife Curve L Tetra				\$885
702XL190	Carbolife Quadro				\$885
702XL200	Spinergy TC				\$736
702XL210	Spinergy TC Max				\$864
Tires					
702XM060	Schwalbe One (Lightest option)				\$133
702XM150	IRC EXERACER				\$188
702XM010	Kenda	702XM011 Grey	702XM012 Blue	702XM013 Red	STD
702XM020	Schwalbe Marathon Plus				\$98
702XM040	Shox Solid Tire				\$160
702XM050	Kenda K885 Mountain Bike Tire				\$40
Inner Tubes					
702XY500	Schraeder Valve Inner Tube				STD
Cushions					
702XP080	Cushion foam				\$121
	Density	702XP001 Soft	702XP002 Medium		
	Thickness	702XP003 50mm	702XP004 75mm	702XP005 100mm	
Cushion Covers					
702XP085	Nylon Cushion Cover (Black)				STD
702XP090	Airlite Cushion Cover (Black)				\$98
Jay Products - Please Include Jay Order Form					
702XP900	Jay Cushion	Refer to Jay Order Form			
702XP910	Jay Backrest	Refer to Jay Order Form			
Accessories					
702XR080	Frontwheel Assembly (Only with 702XB010, 702XB080, 702XB150)				\$1,053
702XP120	Black Frame Protectors Neoprene Zip Fasten (Black)				\$84
702XP420	Black Frame Protectors Nylon Velcro Fasten (Black)				\$81
702XR010	Black Seat Pouch Velcro Attach (Black)				\$59
702XR200	Black Rucksack with RGK Logo (Black)				\$151
702XR100	FX carry bag				\$334
702XR090	Tool Kit				STD
702XS060	Calf Strap				STD

Ordering Information

**Estimated MSRP\***  
(\*Price does not include any one off requested options)

Quote

Order

Specification Notes

