



Zippie GS - AADL W138P

October 2026 EIZGS

Date: _____	SUBMITTING FOR: ☐ Quote ☐ Order
Dealer Acct #: _____	
Dealer: _____	PO#: _____
Dealer Contact: _____	Ship To: _____
Dealer Address: _____	Attention: _____
Dealer City: _____ Prov: _____ PC: _____	Address: _____
Dealer Phone: () _____ Fax: () _____	Address: _____
Confirm Via: ☐ Fax ☐ Email	Ship To City: _____ Prov: _____ PC: _____
Mark For: _____	Ship To Phone: () _____ Fax: () _____

✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

1. FRAME

All Zippie GS models grow in width and depth.

☐	087WM70	Folding Growing Zippie GS Transit Included Includes True Fit Growth	\$3,202	☐	087WM72	Folding Growing Zippie GS Transit Optional Includes True Fit Growth	\$3,202
☐	087WM71	Rigid Growing Zippie GS Transit Included Includes True Fit Growth	\$3,202	☐	087WM73	Rigid Growing Zippie GS Transit Optional Includes True Fit Growth	\$3,202

Frame Options

☐	087FC4	Reverse Configuration (N/A w/ Swing-Away hangers non-rem. in reverse config.; lift off hangers recommended.)	NC	☐	087FC3	Standard Configuration	NC
--------------------------	--------	---	----	--------------------------	--------	------------------------	----

2. TRANSIT SYSTEM (For use when Transit optional is chosen above)

Tie Down Requires: Solid Seat and Solid Back or, Growing Seat Pan and JAY Fit or Conventional Back or, Growing Seat Pan and Solid Back or, JAY Fit and/or Conventional Seat & Back

☐	087TO19	Transit Optional Std. Seating (Req. Selection of Seat & Back)	NC	☐	087TO13	Transit Optional w/o Seating	NC
				☐	087TO15	Transit Omit	NC

2. SIZE/SEATING

Frame Width

Height Adjustable Arms NA with 10"-13" Frame Widths

☐	087BW4	Small Cross Brace	NC
	☐ 10" ☐ 11" ☐ 12"		
☐	087BW5	Medium Cross Brace	NC
	☐ 12" ☐ 13" ☐ 14"		
☐	087BW6	Large Cross Brace	NC
	☐ 14" ☐ 15" ☐ 16"		
☐	087BW7	X-Large Cross Brace	\$89
	☐ 16" ☐ 17" ☐ 18"		

Frame Depth and Back Post Location

☐	087FL3	Short Frame	NC
	☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14"		
☐	087FL9	Regular Frame	NC
	☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16"		
☐	087FL18	Long Frame	NC
	☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18"		
☐	087FL94	X-Long Frame	\$139
	☐ 16" ☐ 17" ☐ 18" ☐ 19" ☐ 20"		

2. SIZE/SEATING CONT.

Seats			
☐	087ST55	Solid Seat (Req. Solid Back) (N/A w/ 10-11, 19-20 inch frame depth or 17-18 inch frame width.)	\$320
☐	087PN3	Growing Seat Pan Depth (Fig. 1)	\$278
☐	087ST27	Omit Seat	STD

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list. Cushion name or part number must be filled out with desired width and depth.

☐	087JC20	JAY Bundled Cushion	Price List
		Name / Part #	
		Width: Depth:	

REFERENCE ONLY

Seat Depth Pan Size	Depth Range	
	Minimum	Maximum
Short 10-12"	10"	14"
Medium 13"-14"	11 3/4"	16 1/4"
Long 15"-17"	13 3/4"	18 1/4"
X-Long 18"-20"	15 3/4"	20 1/4"




Figure 7.1 - Frame Depth Growth Range Matrix

3. CASTERS/SEAT HEIGHT

Caster Wheel			
☐	087CW30	3" Micro Light Up	\$73
☐	087CW31	4" Micro Light Up	\$110
☐	087CW51	5" Micro Light Up	\$110
☐	087CW5	5" Low-Profile Polyurethane	NC
☐	087CW7	5" x 1.5" Semi-Pneumatic	\$110
☐	087CW33	6" Pneumatic	\$110
☐	087CW34	6" Polyurethane	\$0.00
☐	087CW35	6" x 1.5" Semi-Pneumatic	\$110
☐	087CW21	8" Pneumatic	\$146
☐	087CW32	8" Polyurethane	NC
☐	087CW39	8" x 1.5" Semi-Pneumatic	\$146
☐	087CW36	8" x 2" Pneumatic	\$247
☐	087CW37	8" x 2" Pneumatic w/Airless Insert	\$393

Caster Options			
☐	087CS5	3/4" Longer Fork Stem Bolt	NC
☐	087CS9	1 1/2" Longer Fork Stem Bolt	NC
☐	087CS3	Caster Pin Locks (N/A with 8" x 2" casters or Transit Opt.)	\$152
☐	087CW38	Omit Caster, Wheels & Tires	NC

4. REAR WHEELS/SEAT HEIGHT

Axle Plate			
☐	087AL17	One Arm Drive (Avail w/ 20", 22", 24" Mag Rear Wheels & 13"-18" wide)	\$1,382
☐	087OA2	Standard Handrim	NC
☐	087HM1	Right	NC
☐	087HM2	Left	NC
☐	087OA1	Small Handrim	NC
☐	087HM1	Right	NC
☐	087HM2	Left	NC

Axle			
☐	087AX26	One Arm Drive Axle (Only avail. w/One Arm Drive Option)	STD
☐	087AX1	Quick Release	NC
☐	087AX11	Quad Quick Release (N/A with Lite Spokes)	\$182

4. REAR WHEELS/SEAT HEIGHT CONT.

Rear Wheel Size			
☐	087RW48	12" Mag (N/A with handrims)	\$47
☐	087RW14	16" Mag (N/A with handrims)	NC
☐	087RW3	20" Mag	NC
☐	087RW8	22" Mag	NC
☐	087RW1	24" Mag	NC
☐	087RW74	20" Lite Spoke	NC
☐	087RW75	22" Lite Spoke	NC
☐	087RW54	24" Lite Spoke	\$320
☐	087RW25	Omit Rear Wheel & Axle	NC

Handrim			
☐	087HR1	Aluminum Anodized	NC
☐	087HR4	Plastic Coated	\$146
☐	087HR7	Omit (N/A with wheel/axle credit)	NC

Projections (N/A on Plastic Coated)

**Only available with 24" wheels & N/A with spoke guards.*

☐	087HR8	8 Oblique	\$146
--------------------------	--------	-----------	-------

Rear Wheel Tire			
☐	087RT1	Pneumatic	NC
☐	087RT2	Pneumatic Airless Insert	\$247
☐	087RT3	Full Polyurethane (N/A w/ 12" and 16" rear wheels)	NC
☐	087RT4	Low Polyurethane (Only available with 12" and 24" wheels)	NC
☐	087RT6	High Pressure Clincher (Avail. w/ 24" spoke wheels, N/A with mag wheels)	\$105

5. WHEEL LOCK/HUB LOCK

Wheel Lock			
<i>N/A with Hub Locks</i>			
☐	087WL1	High Mount Push	NC
☐	087WL2	High Mount Pull	NC
☐	087WL25	Foot Lock (may limit folding width, N/A in rev. config.)	NC
☐	087WL10	Omit	NC
☐	087WL12	Do Not Mount	NC

Hub Lock			
<i>*Only available with 12" or 16" wheels. N/A with One Arm Drive.</i>			
☐	087LK01	12" Hub Lock	\$162
☐	087LK02	16" Hub Lock	\$162

Extension Handle			
☐	087EH1	6" Pair	\$73

6. HANGER and FOOTRESTS

Hanger Mounted Footrest (Footplates Mount to Hanger Tubing)			
Hanger Type			
☐	087H89	70° Swing In/Out Front Mount	NC
☐	087H90	80° Swing In/Out Front Mount	NC
☐	087H6	90° Swing Away Front Mount	NC
☐	087H9	70° HD Front Mount Lift-Off	NC
☐	087H10	80° HD Front Mount Lift-Off	NC
☐	087H11	90° HD Front Mount Lift-Off	NC
☐	087H23	90° Elevating (N/A 10"-11" wide)	\$339
Extensions (may interfere with some seat heights & casters)			
☐	087E23	2" (not mounted)	\$53
☐	087E24	4" (not mounted)	\$53
Footplate			
☐	087F9	Kids Angle Adjustable (N/A w/17"-18" widths)	NC
☐	087F8	Adult Angle Adj. (Available with 10"-18" widths)	\$176
☐	087F62	Footboard (Available with 10"-16" widths)	\$230
Ext. Tube Mounted Footrest (Footplate Mount to Ext. Tube)			
Hanger Type			
☐	087H85	60° Swing In/Out	NC
☐	087H86	70° Swing In/Out	STD
☐	087H116	80° Swing In/Out (N/A with 6" or 8" Casters)	STD
☐	087H24	Articulating (N/A w/composite footplates or 10"-13")	\$526
☐	087H14	Elevating (N/A w/ composite or platform footplates or < 16" width)	\$339
Extensions			
☐	087E1	Short (11.5" - 14")	NC
☐	087E2	Med (14" - 17") (N/A w/ ELR)	NC
☐	087E5	Long (17" - 19")	NC
☐	087E16	High Mount	\$97
Footplate			
☐	087F1	Composite (N/A In 10", 11" or 12" width)	STD
☐	087F3	Aluminum Angle Adjustable	\$176
☐	087F4	Kids Angle Adjustable (N/A w/17"-18" width)	NC
☐	087F5	Aluminum (Only available with elevating leg rests 14"-18")	NC
☐	087F48	Platform Flip-Up (Avail on 60° 70° and 80° ext. mount hangers only)	\$230

Front Seat to Floor Height Guide (Reference Only)							
Rear Wheel Size	CASTER SIZE						
	3"	4"	5"	6"	8"X1"	8"X1.5"	8"x2"
12"	14.0"	14.0"	14.0"	N/A	N/A	N/A	N/A
	14.5"	14.5"	14.5"				
	15.0"	15.0"	15.0"				
	15.5"	15.5"	15.5"				
	16.0"	16.0"	16.0"	16.0"			
12" Low Poly	14.0"	14.0"	14.0"	15.0"	N/A	N/A	N/A
	14.5"	14.5"	14.5"				
	15.0"	15.0"	15.0"				
16"	14.0"	14.0"	14.0"	N/A	N/A	N/A	N/A
	14.5"	14.5"	14.5"				
	15.0"	15.0"	15.0"				
	15.5"	15.5"	15.5"	15.5"			
	16.0"	16.0"	16.0"	16.0"			
	16.5"	16.5"	16.5"	16.5"	16.5"	16.5"	16.5"
20"	17.0"	17.0"	17.0"	17.0"	17.0"	17.0"	N/A
	14.0"	14.0"	14.0"	N/A	N/A	N/A	N/A
	14.5"	14.5"	14.5"				
	15.0"	15.0"	15.0"				
	15.5"	15.5"	15.5"	16.0"			
	16.0"	16.0"	16.0"	16.5"			
	16.5"	16.5"	16.5"	17.0"			
	17.0"	17.0"	17.0"	17.5"			
	17.5"	17.5"	17.5"	18.0"			
	N/A	18.0"	18.0"	18.5"			
	N/A	N/A	18.5"	19.0"			
	N/A	N/A	N/A	15.0"			
22"	15.0"	15.0"	15.0"	15.5"	16.5"	16.5"	N/A
	15.5"	15.5"	15.5"	16.0"	17.0"	17.0"	
	16.0"	16.0"	16.0"	16.5"	17.5"	17.5"	
	16.5"	16.5"	16.5"	17.0"	18.0"	18.0"	N/A
	17.0"	17.0"	17.0"	17.5"	18.5"	18.5"	17.5"
	17.5"	17.5"	17.5"	18.0"	19.0"	19.0"	18.0"
	N/A	18.0"	18.0"	18.5"	19.5"	19.5"	18.5"
	N/A	N/A	18.5"	19.0"	20.0"	20.0"	19.0"
24"	15.5"	15.5"	15.5"	15.5"	N/A	N/A	N/A
	16.0"	16.0"	16.0"	16.0"	16.5"	16.5"	16.5"
	16.5"	16.5"	16.5"	16.5"	17.0"	17.0"	N/A
	17.0"	17.0"	17.0"	17.0"	17.5"	17.5"	17.5"
	17.5"	17.5"	17.5"	17.5"	18.0"	18.0"	18.0"
	N/A	18.0"	18.0"	18.0"	18.5"	18.5"	18.5"
	N/A	N/A	18.5"	18.5"	19.0"	19.0"	19.0"
Based on 1.5" seat angle. To achieve desired seat height see matrices on last page.							

6. FOOTRESTS CONT.

Footrest Hanger Options

☐	087IG2	Neoprene Covers (Avail. on 60°/ 70° / 80° ext. mount hangers only)	\$61
☐	087HL1	Heel Loops	NC
☐	087A110	Toe Loop (Only avail. w/ adult angle adjust. footplates)	\$43

☐	087HL5	Leg Strap	NC
☐	087A54	Extra Single Leg Strap	\$61
☐	087H33	Omit Hanger	NC
☐	087F7	Omit Footplate & Ext. Tube	NC
☐	087HL4	Omit Heel Loop	NC
☐	087HL6	Omit Single Leg Strap	NC

7. BACKREST/SEATING

Back Height

(Solid Back Depth = 2" Less Than Frame Depth)

☐	087B240	Short Solid Back 12"	\$411
☐	087B241	Med. Solid Back 14"	\$411
☐	087B242	Tall Solid Back 16"	\$411
☐	087BU14	Omit Backrest	NC

JAY Backs (Optional)

Ordering JAY backs on-chair, See the JAY Price List. Requires Part #

☐	087JB02	JAY Bundled Backrest	Price List
		Backrest Name or Part #	
		Width: Height:	

Back Posts

Fixed Height

☐	087BP9	Low (19" Ht)	NC
☐	087BP10	Medium (22" Ht)	NC
☐	087BP11	Tall (25" Ht)	NC
☐	087BP13	Removable Extension	\$176
☐	087RG1	Folding Adjustable Rigidizer	\$207


Height Adjustable

Rigid Frame w/ Angle Adjust Push Handle

☐	087B202	Std. (16"-22" Height.)	NC
☐	087PH8	Std. Angle Adj. Push Handle	NC
☐	087PH9	Long Angle Adj. Push Handle	NC
☐	087B203	Tall (19"-25" Height.)	NC
☐	087PH8	Std. Angle Adj. Push Handle	NC
☐	087PH9	Long Angle Adj. Push Handle	NC

Folding Frame with Push Handle

☐	087B204	Std. (17"-23" Height.)	NC
☐	087BP47	Swing Away Handle (N/A with dual triggers)	\$85
☐	087B205	Tall (20"-26" Height.)	NC
☐	087BP47	Swing Away Handle (N/A with dual triggers)	\$85



8. ARMRESTS

Side Guard

(N/A with height adjustable armrests or T-Post)

☐	087SG33	Plastic Std	\$188
☐	087SG29	Plastic Low	\$188

Armrest and Armpad

☐	087AA12	T-Post Regular w/Foam	NC
☐	087AA13	T-Post Tall w/Foam	NC
☐	087AA3	Flip-Up	\$0.00
☐	087AA1	Adjustable Locking Flip-Up	\$339
☐	087AA7	Length Adjust. Locking Flip-Up	\$339
☐	087AM23	Desk Classic (10") (N/A w/T Post)	NC
☐	087AM24	Full Classic (14") (N/A w/T Post)	NC
☐	087AM9	Foam Short (N/A w/T Post)	NC
☐	087AM10	Foam Long (N/A w/T Post)	NC

Armrest and Armpad - Height Adjustable

☐	087AA11	Single-Post Height Adj. Std.	\$339
☐	087AM23	Desk Classic	NC
☐	087AM21	Desk Contour	NC
☐	087AM24	Full Classic	NC
☐	087AM22	Full Contour	NC
☐	087AA29	Single-Post Height Adj. Low	\$339
☐	087AM23	Desk Classic	NC
☐	087AM21	Desk Contour	NC
☐	087AM24	Full Classic	NC
☐	087AM22	Full Contour	NC

Armrest Receiver for Height Adjustable Armrest

☐	087AR17	Low	NC
☐	087AR18	Std	NC

9. KEY OPTIONS

Anti-Tip

☐	087AO2	Front Anti-Tip (N/A with caster pin locks or transit option)	\$157
☐	087AT17	Rear Anti-Tip	\$121

Positioning Belt (w/clamps)

☐	087A112	1" Dual Pull Side Squeeze Padded	\$132
☐	087A117	Belt 2" Aircraft Buckle	\$139
☐	087A118	Belt 2" Aircraft Buckle Padded	\$204
☐	087A141	Belt 1.5" Auto Buckle	\$97
☐	087A814	1" Comfort Fit APB5116	\$146

Upper Extremity Support Tray

☐	087A29	X-Small (Suggested 10"-12" chair width)	\$381
☐	087A30	Small (Suggested 13"-14" chair width)	\$381
☐	087A31	Medium (Suggested 15"-17" chair width)	\$381
☐	087A32	Large (Suggested 18"-20" chair width)	\$381

Spoke Guard

(N/A with 12" and 16" rear wheels & Projections Handrims)

☐	087A23	Black	\$157
☐	087A303	Clear Zippie	\$157

Canopy (SunShade)

(Need 2" back post clearance to mount canopy hardware, may require next height based on seating; please note some interference may occur with some headrest configurations)

☐	0878A06	Canopy (Sunshade) 10"-14" Width	\$428
☐	0878A07	Canopy (Sunshade) 15"-20" Width	\$428

Back Pack (Kids Backpack included)

☐	087A1	Utility Backpack - Black	\$112
☐	087A2	Additional Kids Black Backpack	\$105

Miscellaneous

☐	087A83	Tool Kit	\$43
☐	087A8	Black Seat Pouch	\$83

10. COLOR

Standard Colors

☐	087C211	Celestial Sapphire	NC
☐	087C212	Desert Sage	NC
☐	087C213	Desert Tan	NC
☐	087C214	Gloss Black	NC
☐	087C215	Matte Black	NC
☐	087C216	Metallic Shadow	NC
☐	087C217	Midnight Blue	NC
☐	087C218	Mystic Purple	NC
☐	087C219	Petal Pink	NC
☐	087C220	Pink Power	NC
☐	087C221	Red Pepper	NC
☐	087C222	Royal Blue	NC
☐	087C223	Salute Green	NC
☐	087C224	Scarlet Noir	NC
☐	087C225	Silver Ore	NC
☐	087C226	Solar Orange	NC
☐	087C227	Starlight White	NC
☐	087C228	Storm Gray	NC
☐	087C229	Turbo Green	NC
☐	087C230	Velocity Yellow	NC

Personalized Colors

☐	087C43	Mossy Oak Camouflage	\$368
☐	087C42	Stars and Stripes	\$368
☐	087C108	Zebra (Choose base colour). NA w/ Gloss Black or Matte Black.	\$368

Additional Options Through Built-4-Me

To inquire about any unique personalization for Zippie products, contact: Phone: (800) 263-3390, Email cscanada@sunmed.com




Commonly Requested Options Include:

- ✓ Dynamic Back
- ✓ Reclining Backrest
- ✓ Personalized Anodizing
- ✓ Personalized Paint
- ✓ O2 Holder with Transit
- ✓ IV Pole with Transit


Zippie GS Seat Height Setup

Caster Fork Size	Caster Stem	Caster Fork Hole Position	Caster Size						
			3"	4"	5"	6"	8"x1"	8"x1.5"	8"x2"
5-1/4" (N/A with 8"x2" caster)	Std	1	N/A	N/A	14.0"	N/A	N/A	N/A	N/A
		2		14.0"	14.5"	15.0"			
		3	14.0"	14.5"	15.0"	15.5"	17.0"	16.5"	
		4	14.5"	15.0"	15.5"	16.0"			
	+3/4"	1	N/A	N/A	14.5"	N/A	N/A	N/A	N/A
		2		14.5"	15.0"	15.5"			
		3	14.5"	15.0"	15.5"	16.0"	17.5"	17.5"	
		4	15.0"	15.5"	16.0"	16.5"			
	+1-1/2"	1	N/A	N/A	15.5"	N/A	N/A	N/A	N/A
		2		15.5"	16.0"	16.5"			
		3	15.5"	16.0"	16.5"	17.0"	18.5"	18.0"	
		4	16.0"	16.5"	17.0"	17.5"			
6" (N/A with 8"x2" caster)	Std	1	N/A	N/A	14.5"	15.0"	N/A	N/A	N/A
		2		14.5"	15.0"	15.5"	16.5"	16.5"	
		3	14.5"	15.0"	15.5"	16.0"	17.0"	17.0"	
		4	15.0"	15.5"	16.0"	16.5"	17.5"	17.5"	
	+3/4"	1	N/A	N/A	15.5"	16.0"	N/A	N/A	N/A
		2		15.5"	16.0"	16.5"	17.5"	17.0"	
		3	15.5"	16.0"	16.5"	17.0"	18.0"	17.5"	
		4	16.0"	16.5"	17.0"	17.5"	18.5"	18.0"	
	+1-1/2"	1	N/A	N/A	16.0"	16.5"	N/A	N/A	N/A
		2		16.0"	16.5"	17.0"	18.0"	18.0"	
		3	16.0"	16.5"	17.0"	17.5"	18.5"	18.5"	
		4	16.5"	17.0"	17.5"	18.0"	19.0"	19.0"	
7" (N/A with 8"x2" caster)	Std	1	N/A	N/A	15.5"	16.0"	17.0"	17.0"	N/A
		2		15.5"	16.0"	16.5"	17.5"	17.5"	
		3	15.5"	16.0"	16.5"	17.0"	18.0"	18.0"	
		4	16.0"	16.5"	17.0"	17.5"	18.5"	18.5"	
	+3/4"	1	N/A	N/A	16.0"	16.5"	18.0"	17.5"	N/A
		2		16.0"	16.5"	17.0"	18.5"	18.0"	
		3	16.0"	16.5"	17.0"	17.5"	19.0"	18.5"	
		4	16.5"	17.0"	17.5"	18.0"	19.5"	19.0"	
	+1-1/2"	1	N/A	N/A	17.0"	17.5"	18.5"	18.5"	N/A
		2		17.0"	17.5"	18.0"	19.0"	19.0"	
		3	17.0"	17.5"	18.0"	18.5"	19.5"	19.5"	
		4	17.5"	18.0"	18.5"	19.0"	20.0"	20.0"	
8" Wide (Must use with 8"x2" caster)	Std	1	N/A						16.5"
		2	N/A						17.5"
	+3/4"	1	N/A						17.5"
		2	N/A						18.5"
	+1-1/2"	1	N/A						18.0"
		2	N/A						19.0"

5-1/4", 6" & 7" Fork



8" Wide Fork



1000 Creditstone Road, Unit 2, Concord ON L4K 4P8