



ORDER FORM

# Freestander

# EILK02

Upright accessible standing frame  
from Early Intervention to adult.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

## Freestander

Size 1 Weight Capacity: 132 lbs | Size 2 Weight Capacity: 176 lbs | Size 3 Weight Capacity: 198 lbs

✓ OPTION	DESCRIPTION	RETAIL
----------	-------------	--------

### Freestander

\$

Selection Required. Includes Chest Pad, Flexible Lateral, PU Knee Supports, Footplate, Tray Cover and Activity Bowl.

<input type="checkbox"/>	141QE31	Size 1 Freestander - Meadow	\$2,643
<input type="checkbox"/>	141QE32	Size 2 Freestander - Mondeo Blue	\$2,913
<input type="checkbox"/>	141QE33	Size 3 Freestander - Mondeo Blue	\$3,277

### Additional Positioning

\$

If Solid Pelvic Plate selected, pelvic band selection required.

<input type="checkbox"/>	141A800	Size 1 Hip Belt - Meadow	\$281
<input type="checkbox"/>	141A854	Size 2 Hip Belt - Mondeo Blue	\$291
<input type="checkbox"/>	141A855	Size 3 Hip Belt - Mondeo Blue	\$316
<input type="checkbox"/>	141ZZ33	Size 1 Solid Pelvic Plate - Meadow	\$454
<input type="checkbox"/>	141ZZ34	Size 2 Solid Pelvic Plate - Mondeo Blue	\$513
<input type="checkbox"/>	141ZZ35	Size 3 Solid Pelvic Plate - Mondeo Blue	\$573
<input type="checkbox"/>	141ZZ30	Size 1 Pelvic Band - Meadow	\$265
<input type="checkbox"/>	141ZZ31	Size 2 Pelvic Band - Mondeo Blue	\$291
<input type="checkbox"/>	141ZZ32	Size 3 Pelvic Band - Mondeo Blue	\$316

### Additional Knee Supports

\$

<input type="checkbox"/>	141XX12	Size 1 Upholstered Knee Cups - Meadow (Pair)	\$182
<input type="checkbox"/>	141XX13	Size 2 Upholstered Knee Cups - Mondeo Blue (Pair)	\$182
<input type="checkbox"/>	141XX14	Size 3 Upholstered Knee Cups - Mondeo Blue (Pair)	\$182
<input type="checkbox"/>	141XX15	Size 4 Upholstered Knee Cups - Mondeo Blue (Pair)	\$182

### Additional Foot Positioning

\$

Sandal selection required with Sandal Raiser.

<input type="checkbox"/>	1418A107	Small Sandals with Straps (Pair)	\$273
<input type="checkbox"/>	1418A108	Medium Sandals with Straps (Pair)	\$291
<input type="checkbox"/>	1418A109	Large Sandals with Straps (Pair)	\$315
<input type="checkbox"/>	1418A110	QTY Sandal Raiser 1" (Each) Requires Sandal Selection.	\$107

### Technical Specifications - Imperial

Frame Sizing Specification	Size 1	Size 2	Size 3
Age (approximate)	1 1/2 - 5 years	4 - 10 years	9 - 18 years
User Body Height (Min / Max)	31.5" / 45.3"	43.3" / 57"	55.1" / 72.8"
Maximum User Weight	132 lbs	176 lbs	198 lbs
Chest Pad Height (Min / Max)	20.7" / 33.5"	26.8" / 43.3"	38.4" / 56.1"
Width Between Chest Supports (Min / Max)	6" / 10"	8" / 12"	10" / 14"
Activity Tray	16" x 21"	18" x 23"	20" x 25"

### Technical Specifications - Metric

Frame Sizing Specification	Size 1	Size 2	Size 3
Age (approximate)	1 1/2 - 5 years	4 - 10 years	9 - 18 years
User Body Height (Min / Max)	800 mm / 1151 mm	1100 mm / 1448 mm	1400 mm / 1849 mm
Maximum User Weight	60 kgs	80 kgs	90 kgs
Chest Pad Height (Min / Max)	470 mm / 680 mm	680 mm / 975 mm	975 mm / 1350 mm
Width Between Chest Supports (Min / Max)	152 mm / 254 mm	203 mm / 305 mm	254 mm / 356 mm
Activity Tray	400 mm x 505 mm	450 mm x 575 mm	500 mm x 625 mm

Sunrise Medical Canada Inc.  
 1000 Creditstone Road, Unit 2  
 Concord, ON L4K 4P8, Canada



MK-129847-CAN Rev. 8 © 04/2025

