



Squiggles Plus Stander

February 2021



Mark For: _____ Date: _____ Dealer Acct #: _____ Dealer: _____ Dealer Contact: _____ Dealer Address: _____ Dealer City: _____ PV: _____ PC: _____ Dealer Phone: _____ Fax: _____ Confirmation Email: _____ Confirm Via: <input type="checkbox"/> Fax <input type="checkbox"/> Email	Submitting for: <input type="checkbox"/> Quote <input type="checkbox"/> Order PO#: _____ <div style="text-align: center;">ADDITIONAL SHIPPING INFORMATION</div> Ship To: _____ Attention: _____ Address: _____ Address: _____ Ship To City: _____ PV: _____ PC: _____ Ship To Phone: _____ Fax: _____
---	--

Squiggles + Stander

Weight Capacity: 48.4 lbs

SECTION A - PACKAGES

Package includes Squiggles Standard Model with Support Shell, Pivot Base, Positioning Pads, Flat Headrest Hardware, Cushion and Laterals, Tray, Small Sandals with Straps and Product Assembly.

	Item Code	Description	Retail
<input type="checkbox"/>	PSPSTB	Squiggles Plus Stander Package - Blue	\$6,533.00
<input type="checkbox"/>	PSPSTP	Squiggles Plus Stander Package - Purple	\$6,533.00
<input type="checkbox"/>	PSPSTR	Squiggles Plus Stander Package - Red	\$6,533.00
<input type="checkbox"/>	PSPSTRV	Squiggles Plus Stander Package - Red Vinyl	\$6,533.00

SECTION B - MADE-TO-ORDER *(Required if you did not choose a package in Section A)*

Standard Model

Selection Required. Includes hip and chest support with sternum prompt; chest and hip lateral supports; size 2 knee cups, footplates and tray attachment.

	Item Code	Description	Retail
<input type="checkbox"/>	184-600	Standard Model with Support Shell	\$2,938.00
<input type="checkbox"/>	184-764	Pivot Base - Pneumatic Angle Adjustment for 3-in-1 Positioning	\$1,581.00

Positioning Pads

Selection Required. Includes chest and hip pad covers and knee support covers.

	Item Code	Description	Retail
<input type="checkbox"/>	184-1800-02	Squiggles Plus Stander Positioning Pads - Blue	\$655.00
<input type="checkbox"/>	184-1800-04	Squiggles Plus Stander Positioning Pads - Purple	\$655.00
<input type="checkbox"/>	184-1800-09	Squiggles Plus Stander Positioning Pads - Red	\$655.00
<input type="checkbox"/>	184-1800-09VIN	Squiggles Plus Stander Positioning Pads - Red Vinyl	\$655.00

Accessories

	Item Code	Description	Retail
<input type="checkbox"/>	118-619	Tray - Gray	\$371.00
<input type="checkbox"/>	120-270	Padded Tray Insert - Gray	\$305.00
<input type="checkbox"/>	117-769	Grab Rail / Toy Bar for Tray	\$301.00
<input type="checkbox"/>	152-1600	Small Sandals w/ Straps (Pair)	\$252.00
<input type="checkbox"/>	152-2600	Medium Sandals w/ Straps (Pair)	\$269.00

Headrest

Cushion selection required with hardware

Item Code	Description	Retail
<input checked="" type="checkbox"/> 117-697	Flat Headrest Support Hardware	\$228.00
<input type="checkbox"/> 184-1812-02	Flat Headrest Cushions - Blue	\$116.00
<input type="checkbox"/> 184-1812-04	Flat Headrest Cushions - Purple	\$116.00
<input type="checkbox"/> 184-1812-09	Flat Headrest Cushions - Red	\$116.00
<input type="checkbox"/> 184-1812-09VIN	Flat Headrest Cushions - Red Vinyl (Pair)	\$116.00
<input checked="" type="checkbox"/> 118-667	Lateral Bracket Assembly Hardware (Pair) <i>Requires Flat Headrest Support Hardware 117-697 and Cushion</i>	\$179.00
<input type="checkbox"/> 184-1814-02	Flat Headrest and Lateral Cushions - Blue	\$197.00
<input type="checkbox"/> 184-1814-04	Flat Headrest and Lateral Cushions - Purple	\$197.00
<input type="checkbox"/> 184-1814-09	Flat Headrest and Lateral Cushions - Red	\$197.00
<input type="checkbox"/> 184-1814-09VIN	Flat Headrest and Lateral Cushions - Red Vinyl (Pair)	\$197.00
<input type="checkbox"/> 184-1838-02	Spine Cap Cover - Blue	\$42.00
<input type="checkbox"/> 184-1838-04	Spine Cap Cover - Purple	\$42.00
<input type="checkbox"/> 184-1838-09	Spine Cap Cover - Red	\$42.00
<input type="checkbox"/> 184-1838-09VIN	Spine Cap Cover - Red Vinyl	\$42.00

Assembly

Item Code	Description	Retail
<input type="checkbox"/> ASSY-LKY	Product Assembly Fee	\$150.00

Technical Specifications

Specification	Imperial	Metric
Age (approximate)	1 - 5 years	1 - 5 years
Maximum User Weight	48.4 lbs	22 kgs
User Standing Height	28.5" - 43.3"	730 mm - 1100 mm
Chest Width Between Laterals	6.3" - 9"	160 mm - 230 mm
Chest Support Height	21.7" - 32"	550 mm - 820 mm
Chest Support Angle Adjustment	+/- 30°	+/- 30°
Chest Support Depth Adjustment	2"	50 mm
Hip Width Between Hip Guides	6.3" - 9"	160 mm - 230 mm
Hip Support Height from Footplate to Mid Buttocks	11.8" - 20.8"	300 mm - 530 mm
Distance Between Midline of Knees	5.5" - 8.4"	140 mm - 215 mm
Footplate Angle	10° Plantarflexion / 10° Dorsiflexion	10° Plantarflexion / 10° Dorsiflexion
Small Sandal Size	6.3" L x 3" W	175 mm L x 75 mm W
Medium Sandal Size	8.3" L x 3.2" W	210 mm L x 80 mm W
Tray Height: Prone	21.9" - 32.5"	555 mm - 825 mm
Tray Height: Supine	18.5" - 31.5"	480 mm - 800 mm
Tray Angle Adjustment	Prone 30° / Supine 40°	Prone 30° / Supine 40°
Stander Support Frame Weight	15.4 lbs	7 kgs
Stander Support Frame Dimensions	26.8" L x 15.4" W x 15.4" H	680 mm L x 390 mm W x 390 mm H
Pivot Base Weight	14.3 lbs	6.5 kgs
Pivot Base Dimensions (Unfolded)	37.4" L x 23.2" W x 20.3" H	950 mm L x 590 mm W x 515 mm H
Pivot Base Dimensions (Folded)	26" L x 23.2" W x 20.3" H	660 mm L x 590 mm W x 515 mm H
Pivot Base Angle Range	Prone 90-70° / Supine 90-160°	Prone 90-70° / Supine 90-160°

