

CUSTOMER				SHIP TO			
Company				Name			
Account #		Date		Address			
<input type="checkbox"/> Quote	<input type="checkbox"/> PO #			City		Prov	
Buyer				Postal Code		Phone	
Mark For				Email			
<input type="checkbox"/> Quote has Additional Sunrise Products with different Quote / Order #				Quote / Order #			

Technical Data Overview			
Maximum User Weight	120 Kg (265 lbs)	Motor	150 W
Wheelchair Rear Wheel Sizes	24" (540)	Operating Voltage	24V
Total Device Weight	12 kg (26.4 lbs)	Maximum Motor Noise	60Db
Per Battery Weight (4.0Ah)	0.6 kg (1.3 lbs)	Battery Type	Lithium Ion
Optimal Range Per Charge	15 km (9.3 mi)	Device Charger Type	28V DC 5.7A 160W
Maximum Speed	6 kph (3.7 mph)	Docking Charger	28V DC 1.5A IP44
Maximum Slope	6°	Control Box	3.7V 1Ah 3.7Wh
Motor Wheel Size	24" (540) Tires	Bluetooth Standard	BLE (125Kbps Coded PHY)

BASE			
✓ OPTION	Description		RETAIL
Empulse M90 Option			\$
<input type="checkbox"/> 139MW94	<b>Empulse M90 Wheels</b> M90 Wheels, Dual Controller, Mounting Hardware, Quick Release Axles, and Device Charger		\$ 11,995.00
	<input checked="" type="checkbox"/> <b>Battery</b> 15km (9.3 mi) range		\$ 500.00
Options and Accessories			\$
<input type="checkbox"/> 139AG44	Additional Batteries	Includes two interchangeable M90 batteries, usable with either wheel.	\$ 650.00
<input type="checkbox"/> 139BJ39	Docking Charger	Includes docking charger, capable of charging two M90 batteries.	\$ 600.00
<input type="checkbox"/> 139JM50	Swing Away Mount Left Hand	Left hand swing away mount for M90 Dual Controller.	\$ 550.00
<input type="checkbox"/> 139JM51	Swing Away Mount Right Hand	Right hand swing away mount for M90 Dual Controller.	\$ 550.00
<input type="checkbox"/> 139NP02	Wrist support	Additional wrist and hand support for the M90 Dual Controller.	\$ 200.00
<input type="checkbox"/> 139DI56	Thigh Strap	Includes thigh strap and hardware, to mount the Dual Controller to the riders leg.	\$ 75.00
<input type="checkbox"/> 1399A33	Docking Tray	Docking tray for the M90 Dual Controller	\$ 50.00

## WHEELCHAIR

✓ OPTION	DESCRIPTION	DESCRIPTION	
----------	-------------	-------------	--

### Sunrise Medical Wheelchair Models - Specific mounting hardware will be provide for the chair selected \$

Select the specific model of the wheelchair the M90 Wheels will be installed on.

#### Quickie Rigid Chairs NC

<input type="checkbox"/>	1399V58	Quickie Nitrum	<input type="checkbox"/>	1399V56	Quickie 5R	<input type="checkbox"/>	1399V59	Quickie Qri
<input type="checkbox"/>	1399V80	Quickie GP	<i>RGK -The Empulse M90 works with select RGK models. Contact your local RGK representative for details.</i>					
<input type="checkbox"/>	1399V66	Quickie 7R						

#### Quickie Folding Chairs NC

<input type="checkbox"/>	1399V69	Quickie QS5X	<input type="checkbox"/>	1399V68	Quickie 2	<input type="checkbox"/>	1399V60	Quickie Qxi
<input type="checkbox"/>	1399V77	Quickie Qx	<input type="checkbox"/>	1399V78	Quickie 2 Lite			

#### Quickie Tilt In Space Chairs NC

<input type="checkbox"/>	1399V05	Quickie Iris						
--------------------------	---------	--------------	--	--	--	--	--	--

**Zippie Pediatric** - The Empulse M90 is only available in a 24" (540) wheel.

<input type="checkbox"/>	1399V04	Zippie GS	<input type="checkbox"/>	1399V64	Zippie Zone	Zippie X'Cape	Coming Soon
<input type="checkbox"/>	1399V33	Zippie Iris					

### Other Wheelchair Models - Specific mounting hardware will be provide for the chair selected \$

Specific model of the wheelchair the M90 Wheels by Empluse will be used on

Sunrise Medical cannot guarantee fitment on all configurations of the below chairs, please contact your dealer for guidance.

#### TiLite Chairs - Coming Soon NC

<i>TiLite ZR</i>	<i>TiLite TRA</i>	<i>TiLite X</i>					
<i>TiLite ZRA</i>	<i>TiLite Aero Z</i>	<i>TiLite Aero T</i>					
<i>TiLite TR</i>	<i>TiLite 2GX</i>	<i>TiLite AeroX</i>					

#### Ki Mobility Chairs - Coming soon NC

<i>Ki Mobility Catalyst ST</i>	<i>Ki Mobility Liberty FT</i>						
<i>Ki Mobility Catalyst 5</i>	<i>Ki Mobility Focus CR</i>						

#### Motion Composites - Coming soon NC

<i>Apex A</i>	<i>Helio C2</i>	<i>Helio A6</i>					
<i>Veloce</i>	<i>Helio A7</i>						

### Wheelchair Wheel Size \$

M90 Wheels are **only** offered in 24" (540) wheels

M90 Wheels are **only** compatible with 0° of camber

<input type="checkbox"/>	139WS6	24" Wheel Size ( 540 Wheels)	NC
--------------------------	--------	------------------------------	----

### Handrim \$

Select only one Handrim option (24" Wheels Only)

<input type="checkbox"/>	139HR45	Powder Coated Aluminum Black	NC	<input type="checkbox"/>	139HR39	Ellipse 3R	\$	550.00	
<input type="checkbox"/>	139HR4	Plastic Coated	\$	146.00	<input type="checkbox"/>	139HR22	Natural Fit - Standard Grip	\$	374.00
<input type="checkbox"/>	139HR7	Omit Handrim	NC						

### Tire \$

Select only one Tire option (24" Wheels Only)

<input type="checkbox"/>	139RT21	Pneumatic Schwalbe Marathon Plus	\$	157.00	<input type="checkbox"/>	139RT25	Primo Orion Poly	\$	271.00
--------------------------	---------	----------------------------------	----	--------	--------------------------	---------	------------------	----	--------