



FORM. FUNCTION. FIT

JAY® Fit Seating simplifies the world of custom seating by offering a “semi-custom” seat and back solution. A great option for individuals who require more complex seating than off-the-shelf options can provide, while still allowing for adjustments and growth.



FIT

Adjustability

The JAY Fit's modular and interchangeable components allow the system to be grown and modified over time.

- Well Inserts
- Lateral Well Reducers
- Lateral Pelvic Supports
- Lateral Thigh Supports
- Flip-Down Medial Thigh Supports
- Posterior Contour Strips
- Lateral Thoracic Supports
- Lateral/Posterior Pelvic Supports
- Lateral Pyramids
- Thoracic Pads

Available with Whitmyer Headrest and JAY Positioning Supports

The JAY Fit system is available with the most widely used positioning support options. The most commonly-used options are now on the order form and can be easily mounted at the factory.

- JAY Pelvic Positioning Belts
- JAY Anterior Trunk Supports
- JAY Shoe Holders
- Whitmyer® Headrests
- JAY Spine Align Components

Easy to Order

The JAY Fit's online ordering system allows you to spec, quote and order full seating systems. Quickly determine pricing and option compatibility 24/7 from the comfort of the home or office.

- New easy-to-follow order form includes more options.
- Similar process to ordering the JAY J3 back system; It's as simple as selecting your back and options.



New Swing-Away Lateral Design

The JAY Fit system now features a new swing-away lateral design with increased adjustability and improved stability.

- Adjustable for a wide range of mounting options.
- Secure locked position for improved feel when engaged.

JAY FIT SEAT SPECS

HCPCS Code
E2609

Weight Capacity
10" to 18" Widths: 250 lbs.

Seat Width Range
10" to 18"

Seat Depth Range
8" to 20"

JAY FIT BACK SPECS

HCPCS Code
E2617

Weight Capacity
10" to 18" Widths: 250 lbs.

Back Width Range
10" to 22"

Back Height Range
10" to 22"