





THE ULTIMATE
MID-WHEEL DRIVE
EXPERIENCE!





Designed to meet your demanding needs. The Ω stands for quality.



RELIABLE

Every Q700 M is built to keep you moving. The chassis is formed from a single sheet of metal, and its steel components receive automotive undercoat plating. The industry's best motors undergo additional end-of-line testing before leaving our factory.

MANEUVERABLE

True mid-wheel drive allows the Q700 to turn on its own axis, so you enjoy intuitive driving and excellent maneuverability in tight, compact spaces.

PERFORMANCE

With its independent movement of all six wheels, the Q700 M ProGrip suspension gives you the ability to safely crawl up and down steep obstacles, keeping all six wheels on the ground.

ELEVATE AND DRIVE WITH C-ME®

C-Me enables you to elevate up to 12" while traveling up to 3 mph! You can lift or descend at the touch of a button with no need to stop while transitioning. Our suspension doesn't lock out when you drive with C-Me, giving you a smooth and comfortable ride from start to finish.





SEDEO® PRO ADVANCED™

Designed for comfort and aesthetics, the SEDEO PRO ADVANCED covers the most demanding clinical needs. With Biometric Repositioning and a full range of programmable seating positions, it is the best ergonomic seating from QUICKIE.



SEDEO® LITE™

The SEDEO LITE is an adjustable seat frame with a sleek look. Its minimal design alleviates weight capacity concerns when using aftermarket equipment. Traditional back canes are width-adjustable from 15" to 20" and easily fit JAY® backs. An array of options cover mild to moderate clinical needs, and intuitive adjustability meets individual requirements.



SEDEO® PRO™

Truly modular in design, the SEDEO PRO provides a new level of functionality, comfort and design. For moderate to complex clinical needs, the PRO combines a highly adjustable seat frame with an excellent choice of postural support options.



Captain's Seat

The Captain's Seat provides a relaxed sitting position and features a semi-recline, fold-down back, adjustable seat depth, and convertible seat pan that allows you to use your preferred cushion. Available with a 50° tilt and power centermount leg support.

Q700 M SEAT FRAME COMPARISON

	sedeo PRO ADVANCED	SEDEO PRO	SEDEO LITE	CAPTAIN SEAT
Cryo™ Fluid	Available on JAY Syncra Seating	Available on JAY Syncra Seating	Available on JAY Fusion	Available on JAY Fusion
Seat Width	16" to 22"	16" to 26"	15" to 20"	16", 18", 20"
Seat Depth	17" to 22"	14" to 22"	16" to 22"	16" to 22"
Back Support Height	21" to 30"	21" to 30"	18" to 30"	25"
Biometric Repositioning	•			
Memory Seating	•			
Powered Recline	165°	150°		
Anti-shear Back Support	7"	4" (optional)		
Pressure Relief Alert Reminders (PRAR)	•			
Tilt Angle	50°	50°	50°	50°
Anterior Tilt	30°			
Elevate with C-Me	12"	12"	12"	
JAY & WHITMYER	•	•	•	•

Q700 M	See order form as specification	order form as specifications may vary based upon configuration				
HCPC codes	No Power Seating Options	Single Power Seating Options	Multi Power Seating Options			
Group 3 Standard Duty	K0848/K0849	K0856/K0857	K0861			
Group 3 Heavy Duty	K0850	K0858	K0862			
Group 4 Standard Duty	K0868/K0869	K0877/K0878	K0884			
Max. user weight		300 lbs* (450 lbs for HD)				
Seat-to-floor height		16.5" to 19"				
Total width		25"				
Total length		36.5"				
Turning radius		19.8" (Excludes footplates)				
Max. speed	6 mph (Group 3 Standard)	5 mph (Group 3 HD)	8 mph (optional)			
Battery size		Group 34 or Group 24				
Overall weight	Starting at 294 lbs	Starting at 318 lbs	Starting at 366 lbs			
Drive wheel position		Mid				
Ground clearance		3"				
Wheel options		14" drive wheels; 6" casters				
Max. range	(w/G	Up to 18 miles (w/Group 24 batteries and 6 mph motors)				
Transit Brackets	I:	ISO 7176-19:2008 Transit Approved				
Warranty		Lifetime Limited				

COLOR OPTIONS

SHROUD COLORS



COLOR ACCENTS









Order Form



