

<b>Mark For:</b> _____ Date: _____ Dealer Acct #: _____ Dealer: _____ Dealer Contact: _____ Dealer Address: _____ Dealer City: _____ PV: _____ PC: _____ Dealer Phone: ( ) _____ Fax: ( ) _____ Confirmation Email: _____ Confirm Via: <input type="checkbox"/> Fax <input type="checkbox"/> Email	<b>Submitting for:</b> <input type="checkbox"/> Quote <input type="checkbox"/> Order PO#: _____  <p style="text-align: center;"><b>ADDITIONAL SHIPPING INFORMATION</b></p> Ship To: _____ Attention: _____ Address: _____ Address: _____ Ship To City: _____ PV: _____ PC: _____ Ship To Phone: ( ) _____ Fax: ( ) _____
---	--



**S410**

**S410 HP**

300 lbs Std/350 HP

**MAXIMUM USER WEIGHT:**  **kg**

Mandatory information to ensure correct configuration

**ADP DEVICE CODE WP0002088**

<b>Speed</b>	<b>Colour Options</b>	<b>Wheels</b>
● 10 km/hour Std 10.5km/hour HP	● Red, Blue & Black	● 10" / 25 cm black pneumatic
<b>LED Lightning</b>	<b>Tiller</b>	<b>Scooter Control</b>
● Headlight, front/rear indicators, brake, position	● Gas strut easy angle adjust, delta handle bar	● Easy touch control pad with LED battery gauge
<b>User Control</b>	<b>Braking</b>	<b>Suspension</b>
● Wig wag, finger thumb control	● Motor brake lever with detent button safety feature	● Front/rear adjustable
<b>Batteries</b>	<b>Charger</b>	<b>Controller</b>
● Groups 22 (not included)	● 8 A, tiller mounted charging socket	● Standard Model: 90 Amp ● High Performance Model: 120 Amp
<b>Seating, size 49 cm x 49 cm</b>		
● Comfort seat with sports upholstery	● Armrests width, depth, angle adjustable	● Headrest
● Easy adjust recline	● Comfort arm pads, flip up	● Seat slider with rotation
● Height adjustable 45-50 cm 3 positions	● Buckle type lap strap	● Accessory mount with handle
<b>Safety</b>		
● Front and rear bumpers	● Side and rear reflectors	● Anti-tip wheels
● Carbon effect mirrors (Left/Right)	● Key switch	

**Fortress S410**

Part Number	Model	MSRP	Qty	Total
<input type="checkbox"/> 22200410-001	S410 Std Red	\$ 4,715.00		
<input type="checkbox"/> 22200410-002	S410 Std Blue	\$ 4,715.00		
<input type="checkbox"/> 22200410-006	S410 Std Black	\$ 4,715.00		
<input type="checkbox"/> 22200411-001	S410 HP Red	\$ 5,222.00		
<input type="checkbox"/> 22200411-002	S410 HP Blue	\$ 5,222.00		
<input type="checkbox"/> 22200411-006	S410 HP Black	\$ 5,222.00		
<b>Options &amp; Accessories</b>				
<input type="checkbox"/> N14751001	Large Basket	\$ 102.00		
<input type="checkbox"/> 22219060	Seat Belt	\$ 81.00		
<input type="checkbox"/> 22219025	Mirrors (Right & Left)	\$ 65.00		
<input type="checkbox"/> BW-12550-C	Group 22 NF Battery 12V 55AH (Each)	\$ 460.00		
<input type="checkbox"/> 22219107	Flat Free foam Filled tires- set of 4	\$ 687.00		
<input type="checkbox"/> SP-000 228 229	Safety Flag	\$ 61.00		
<input type="checkbox"/> 22219094	Crutch Holder	\$ 192.00		