



AADL ORDER FORM

K285-1
K285-2
K285-3

October 1, 2026

MyWay+

EILK01

Unique, open frame gait trainer designed to support children in an upright, hands-free position.

Customer			
Account			
Account #		Date	
☐ Quote	☐ Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

MyWay+		
✓ OPTION	DESCRIPTION	RETAIL
Frame		\$
Selection Required. Includes multi-functional casters.		
☐ 122FT17 K285-1	Size 1 MyWay+ Frame	\$4,629
☐ 122FT18 K285-2	Size 2 MyWay+ Frame	\$5,339
☐ 122FT19 K285-3	Size 3 MyWay+ Frame	\$5,561
Harness		\$
Selection Required.		
☐ 122NH34	Size 1 Harness - Purple	\$398
☐ 122NH35	Size 2 Harness - Green	\$402
☐ 122NH36	Size 3 Harness - Blue	\$437
☐ 122NH37	Size 4 Harness - Red	\$477
☐ 122NH38	Size 5 Harness - Orange	\$523
☐ 122NH39	Size 6 Harness - Grey	\$542








Size 1 - Purple	Size 2 - Green	Size 3 - Blue	Size 4 - Red	Size 5 - Orange	Size 6 - Grey
-----------------	----------------	---------------	--------------	-----------------	---------------

For ease of ordering configured units,
use the NEW [MyWay+ configurator](#).

MyWay+			
✓ OPTION	DESCRIPTION	RETAIL	
User Handles		\$	
User Handles MUST match Frame Size. Complex User Handles Ball Grip standard; Cylinder Grip optional.			
☐ 122HY1	Size 1 User Handles - Pair	\$364	Standard User Handles 
☐ 122HY2	Size 2 User Handles - Pair	\$364	
☐ 122HY3	Size 3 User Handles - Pair	\$364	Complex User Handles 
☐ 122HY17	Size 1 Complex User Handles	\$1,873	
☐ 122HY18	Size 2 Complex User Handles	\$1,951	
☐ 122HY19	Size 3 Complex User Handles	\$2,003	
☐ 122HY22	Cylinder Grip - Pair	\$78	
Guide Handles		\$	
Use Size 1 Adult Handles for Size 1 Frame. Use Size 2 Adult Handles for Size 2 & 3 Frame.			
☐ 122HY20	Size 1 Adult Handles - Fits Size 1 Frame	\$448	Adult Handles 
☐ 122HY21	Size 2 Adult Handles - Fits Size 2 & 3 Frame	\$448	
Headrests		\$	
☐ 122HE49	Headrest - One Size	\$247	Standard  Complex 
☐ 122HE30	Complex Headrest - One Size	\$729	
Saddles		\$	
Saddle Size 1 recommended for Size 1 & 2 Frames. Saddle size 2 recommended for Size 2 & 3 Frames.			
☐ 1229A19	Size 1 Saddle	\$387	Standard Saddle 
☐ 1229A20	Size 2 Saddle	\$390	
☐ 1229A21	Size 1 Narrow Saddle	\$448	Narrow Saddle 
☐ 1229A22	Size 2 Narrow Saddle	\$448	
Foot and Leg Positioning Accessories		\$	
Ankle Guides MUST match Frame Size.			
☐ 122DI60	Size 1 Drop - Fits Size 1 Narrow Saddle ONLY	\$240	Saddle Drop 
☐ 122DI61	Size 2 Drop - Fits Size 2 Narrow Saddle ONLY	\$260	
☐ 122DI62	Size 3 Drop - Fits Size 2 Narrow Saddle ONLY	\$291	Ankle Guides 
☐ 122DI63	Size 1 Ankle Guides - Pair	\$277	
☐ 122DI64	Size 2 Ankle Guides - Pair	\$277	Ankle Guides 
☐ 122DI65	Size 3 Ankle Guides - Pair	\$277	
☐ 122DI66	Hoist Strap - Pair	\$67	
MyWay Pedal		\$	
MyWay Pedal Compatible with Size 2 & Size 3 MyWay+ Only			
MyWay Pedal Selection Requires Sandal Selection			
☐ 1223A08	Smart Pedal Walker with Smart Tech	\$1,827	MyWay Pedal 
☐ 1228A107	Small Sandals with Straps (Pair)	\$273	
☐ 1228A108	Medium Sandals with Straps (Pair)	\$291	Ankle Wraps 
☐ 1228A109	Large Sandals with Straps (Pair)	\$315	
☐ 122DI72	Small Ankle Wraps	\$273	
☐ 122DI73	Large Ankle Wraps	\$315	

To order individual MyWay+ accessories, see parts listed below.

MyWay+ Accessories

✓ PART	DESCRIPTION	RETAIL
Harness		\$
☐ 173-4100-04	Size 1 Harness - Purple	\$398
☐ 173-4200-07	Size 2 Harness - Green	\$402
☐ 173-4300-02	Size 3 Harness - Blue	\$437
☐ 173-4400-09	Size 4 Harness - Red	\$477
☐ 173-4500-06	Size 5 Harness - Orange	\$523
☐ 173-4600-08	Size 6 Harness - Grey	\$542
Handles		\$
User Handles MUST match Frame Size. Complex User Handles Ball Grip standard; Cylinder Grip optional.		
☐ 209-1667	Size 1 User Handles - Pair	\$364
☐ 209-2667	Size 2 User Handles - Pair	\$364
☐ 209-3667	Size 3 User Handles - Pair	\$364
☐ 209-1740	Size 1 Complex User Handles	\$1,873
☐ 209-2740	Size 2 Complex User Handles	\$1,951
☐ 209-3740	Size 3 Complex User Handles	\$2,003
☐ 209-1741	Cylinder Grip - Pair	\$78
Guide Handles		\$
Use Size 1 Adult Handles for Size 1 Frame. Use Size 2 Adult Handles for Size 2 & 3 Frame.		
☐ 209-1691	Size 1 Adult Handles - Fits Size 1 Frame	\$448
☐ 209-2691	Size 2 Adult Handles - Fits Size 2 & 3 Frame	\$448
Headrests		\$
☐ 173-1680	Headrest - One Size	\$247
☐ 209-1680	Complex Headrest - One Size	\$729
Saddles		\$
Saddle Size 1 recommended for Size 1 & 2 Frames. Saddle size 2 recommended for Size 2 & 3 Frames.		
☐ 173-1671	Size 1 Saddle	\$387
☐ 173-2671	Size 2 Saddle	\$390
☐ 209-1671	Size 1 Narrow Saddle	\$572
☐ 209-2671	Size 2 Narrow Saddle	\$598
Foot and Leg Positioning		\$
Ankle Guides MUST match Frame Size.		
☐ 209-1679	Size 1 Drop - Fits Size 1 Narrow Saddle ONLY	\$240
☐ 209-2679	Size 2 Drop - Fits Size 2 Narrow Saddle ONLY	\$260
☐ 209-3679	Size 3 Drop - Fits Size 2 Narrow Saddle ONLY	\$291
☐ 209-1700	Size 1 Ankle Guides - Pair	\$443
☐ 209-2700	Size 2 Ankle Guides - Pair	\$443
☐ 209-3700	Size 3 Ankle Guides - Pair	\$443
☐ 173-1800	Hoist Strap - Pair	\$67

MyWay+ Accessories

✓ PART	DESCRIPTION	RETAIL
MyWay Pedal		\$
☐ 173-1750	Smart Pedal Walker with Smart Tech	\$1,827
☐ 152-1600	Small Sandals with Straps (Pair)	\$273
☐ 152-2600	Medium Sandals with Straps (Pair)	\$291
☐ 152-3600	Large Sandals with Straps (Pair)	\$315
☐ 173-2305	Small Ankle Wraps	\$273
☐ 173-3305	Large Ankle Wraps	\$315

Technical Specifications - Imperial

Frame Sizing Specification	Size 1		Size 2		Size 3	
Age (approximate)	1 - 5 years		4 - 11 years		9 - 16 years	
User Shoulder Height (to floor)	23 - 34"		34 - 47"		42 - 59"	
User Saddle Height (inside leg)	9.5 - 16"		16 - 24"		21.5 - 32"	
Maximum User Weight	55 lbs.		110 lbs.		176 lbs.	
Frame Weight	26.5 lbs.		14.5 lbs.		21 lbs.	
Upper Support Unit Weight	Size 1 Frame and Upper Support Unit do not separate		23 lbs.		24 lbs.	
Frame Width	22.5"		26.5"		28.5"	
Frame Length	28.5"		33.5"		43.5"	
Saddle	9" L x 3.25" W		10" L x 3.25" W			
Harness Sizing Specification	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6
Age (approximate)	1 - 4 years	2 - 6 years	5 - 9 years	8 - 12 years	12 - 15 years	16 - 18 years
Waist at Navel Level	18 - 22"	19.5 - 23.5"	21.5 - 25.5"	23.5 - 27.5"	27.5 - 31.5"	31.5 - 35.5"
Color	Purple	Green	Blue	Red	Orange	Gray
Maximum User Weight	110 lbs.	110 lbs.	110 lbs.	176 lbs.	176 lbs.	176 lbs.

Technical Specifications - Metric

Frame Sizing Specification	Size 1	Size 2	Size 3			
Age (approximate)	1 - 5 years	4 - 11 years	9 - 16 years			
User Shoulder Height (to floor)	635 - 920 mm	870 - 1190 mm	1070 - 1490 mm			
User Saddle Height (inside leg)	220 - 440 mm	400 - 610 mm	550 - 810 mm			
Maximum User Weight	25 kgs	50 kgs	80 kgs			
Frame Weight	12 kgs	6.5 kgs	9.5 kgs			
Upper Support Unit Weight	Size 1 Frame and Upper Support Unit do not separate	10.5 kgs	11 kgs			
Frame Width	570 mm	670 mm	720 mm			
Frame Length	720 mm	850 mm	1100 mm			
Saddle	225 mm L x 80 mm W	255 mm L x 80 mm W				
Harness Sizing Specification	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6
Age (approximate)	1 - 4 years	2 - 6 years	5 - 9 years	8 - 12 years	12 - 15 years	16 - 18 years
Waist at Navel Level	460 - 560 mm	500 - 600 mm	550 - 650 mm	600 - 700 mm	700 - 800 mm	800 - 900 mm
Color	Purple	Green	Blue	Red	Orange	Gray
Maximum User Weight	50 kgs	50 kgs	50 kgs	80 kgs	80 kgs	80 kgs

Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2,
Concord, ON L4K 4P8



MK-130103-CAN Rev. 4 © 10/2026

