



Q300M Mini

AADL -W577P

EIPW40

October 2026

Account:	
Date:	
PO#:	
Buyer Name:	
Marked For:	


Ship To:	
Name:	
Address:	
City, Prov, PC:	
Phone:	

Sunrise Medical Inc.
1000 Creditston Road, Unit 2
Concord, ON L4K 4P8
1-800-263-3390
www.sunrisemedical.ca



Client Measurements

a)	Top of head to seat	Inch
b)	Top of shoulder to seat	Inch
c)	Underarm to seat	Inch
d)	Elbow to seat	Inch
e)	Top of leg to seat	Inch
f)	Shoulder width	Inch
g)	Chest width	Inch
i)	Hip width	Inch
j)	Actual sitting depth	Inch
m)	Foot length	Inch
(n)	Knee to foot	Inch



User Weight

Notes:

Chair Specifications

		Seat Type	SEDEO PRO / JAY Syncra
Weights starting at: (with Batteries)	NPC: 250 lbs SPC: 274 lbs	Seat Size Range	Width 10" to 16" Depth 10" to 19"
Power Seating	50° of CG Tilt / 12" seat elevator	STF Height Range	17" to 19"

Mid-Wheel Drive Base Specifications

Drive Wheel Position	Mid-wheel Drive	Drive Wheels	12" drive wheels w/flat free inserts (black) 14" drive wheels w/flat free inserts (black)
User Weight Capacity	165 lbs.	Tires/Casters	7" Front / 7" Rear (Black)
Min. Turning Radius	17" (excludes footplates)	Charger	8A Dual Mode (Compact)
Base Size	Width 20.5", Length 33"	Battery Size	50Ah Group 22
Drive Control	R-net 90A / VR-2 90A	Color Accents	9 Shroud Colors and 7 Color Accents
Driving Range	Up to 17 miles	Ground Clearance	3"
Drive Motors	4-pole 6mph		

BASE			
OPTION		Description	Retail
Standard Duty Mid Wheel Drive Base - End Item = EIPW40 (AADL Contract Price)			
☐	400PC308	Quickie Q300 M NPC No Power Option for Zippie Seat Included in Base Pricing: Q300 M Base, VR-2 90amp Electronics for No Power Seat Functions, Six-Wheel Suspension, 12" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 165 lbs., STFH = 17".	\$7,877.00
☐	400PC309	Quickie Q300 M SPC Single Power Option for Zippie Seat Included in Base Pricing: Q300 M Base, VR-2 90amp Electronics for No Power Seat Functions, Six-Wheel Suspension, 12" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 165 lbs., STFH = 17".	\$7,877.00
User Weight			
☐	400WG32	Up to 165 lbs.	NC
Motor Package			
☐	400PF99	4-pole 10kph / 6mph	NC
Tracking			
☐	400MI01	SureTrac Module & Motor Encoders (<i>Improves Drive Tracking</i>)	\$1,696.00
Drive Wheels			
☐	400DW13	Drive Wheel, 12", Solid, Knobby, Black	NC
☐	400DW58	Drive Wheel, 14", Solid, Knobby, Black	NC
Caster Wheels			
☐	400CW59	Caster Wheels 7", Solid, Black	NC
Shroud Color			
☐	400C164	Black	NC
☐	400C165	Candy Apple Red	NC
☐	400C166	Candy Blue	NC
☐	400C167	Candy Purple	NC
☐	400C176	Green Apple	NC
☐	400C178	Pearl Pink	NC
☐	400C182	Silver	NC
☐	400C125	White	NC
Color Accents			
Note: Accent colors and shroud colors may not match			
☐	400EZ20	Apple Green	NC
☐	400EZ22	Black	NC
☐	400EZ21	Blue	NC
☐	400EZ23	Orange	NC
☐	400EZ24	Red	NC
☐	400EZ19	White	NC
☐	400EZ25	Silver	NC

BASE

OPTION	Description	Retail
--------	-------------	--------

ISO Transit Tie Downs

☐	400TO50	ISO Transit Package	NC
--------------------------	---------	---------------------	----

LED Lights (Client Cost)

☐	400LI08	LED Fender Lights (<i>Joystick Activation</i>)	\$646.00
☐	400LI02	Omit LED Lights	NC

Seat Heights

Note: Available Seat-to-Floor heights are dependant on Powered Seating. Use the matrix below to determine available ranges for your desired options.
Measured from front of seat pan to floor.

Seat Height				12" Drive Wheel			14" Drive Wheel				
				No Power	Tilt	Elevate	No Power	Tilt	Elevate		
☐	400HT55		17"		N/A	N/A	N/A	N/A	N/A		NC
☐	400HT104		17.5"	N/A		N/A		N/A	N/A		NC
☐	400HT56		18"		N/A						NC
☐	400HT105		18.5"			N/A	N/A		N/A		NC
☐	400HT57		19"								NC

Batteries & Chargers (AADL Contract Pricing)

Group 22 Battery

☐	400AG35	GP22/50Ah Flat Top Pair	\$678.00
--------------------------	---------	-------------------------	----------

Omit Batteries.

☐	400AG24	Battery Harnesses Only (<i>Omit Batteries</i>)	NC
--------------------------	---------	--	----

Chargers

☐	400BJ30	10A Dual Mode	\$162.00
☐	400BJ19	8A Dual Mode	NC

POWERED SEATING

OPTION	Description	Retail
--------	-------------	--------

Power Seating Options (Client Cost)

Seat Elevator

☐	400PI07	Seat Elevator, Powered, 300 mm / 12"	\$3,477.00
☐	400PI08	C-Me Seat Elevator, Powered, 300 mm / 12" (C-Me speed up to 3.0 mph) (N/A with VR-2)	\$3,477.00

Seat Tilt

☐	400PA48	Posterior 50° Tilt (W993 for Power Tilt)	\$6,496.00
--------------------------	---------	---	------------

Pre Set Seat Frame Angle (No Power Seating Options, Increases Seat to Floor Height)

☐	400FA10	0° (Adds 0" to STFH)	NC
☐	400FA3	3° (Adds 1" to STFH)	NC

Backrest Frame

Backrest Type

☐	400BF28	Angle Adjustable	NC
☐	400BF29	40° Manual Recline	\$617.00

Back Canes

☐	400B177	Medium, 17.5"	NC
☐	400B178	Tall, 19.5"	NC
☐	400B179	Height Adjustable, 20" - 26"	\$0.00

Backrest Frame Option

☐	400BF13	Flip Down, Quick Release (N/A w/ Manual Recline)	\$179.00
--------------------------	---------	--	----------

Push Handles

☐	400PH28	Trailing Push Handles	NC
☐	400PH29	Stroller Handle	\$67.00

Towel Bar

☐	400QT03	Shallow	NC
☐	400QT04	Deep (Extends 4" from Back Canes)	NC

Back Cushion Type

☐	400BT95	SPOT Back Cushion (Choose Seating in Section A)	NC
☐	400BT135	JAY J3 (Choose Seating in section C)	NC
☐	400BT133	OMIT (For Custom Mounting. No Headrest Mount)	NC

A - SPOT Backrest Dimensions

Installed Back Height - (From Pan to Top of Back)

☐	400JR46	12	NC
☐	400JR47	14	NC
☐	400JR48	16	NC
☐	400JR71	18	NC
☐	400JR73	20	NC

POWERED SEATING

OPTION	Description	Retail
--------	-------------	--------

B - JAY J3 Back

B - JAY J3 Back - Contour (Seating Benefits)

☐	400JP01	J3 PA Back for 12"-16" Wide	\$912.00
☐	400JP11	J3 PL Back for 14"-16" Wide	\$986.00
☐	400JP13	J3 PD Back for 14"-16" Wide (18" Shell Height Only)	\$1,068.00
☐	400JP39	JAY Zip Back for 10"-14" Wide	\$1,015.00

B - JAY J3 Back Shell Options and Height Ranges.

JAY J3 shells are mounted with a minimim 4" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field within the ranges shown.

JAY J3 Shell Height

☐	400JR81	8" tall shell (JAY Zip Only)	NC
☐	400JR82	10" tall shell (JAY Zip Only)	NC
☐	400JR83	12" tall shell (J3 Mid Thoracic Short or JAY Zip)	NC
☐	400JR84	13.5" tall shell (Mid Thoracic Med)	NC
☐	400JR85	15" tall shell (Mid Thoracic Tall)	NC
☐	250JR14	16.5" tall shell (Upper Thoracic Short)	NC
☐	250JR15	18" tall shell (Upper Thoracic Med)	NC

JAY J3 Installed Back Height (From Pan to Top of Back)

☐	400JR86	12	NC
☐	400JR87	14	NC
☐	400JR88	16	NC
☐	400JR89	18	NC
☐	400JR90	20	NC
☐	400JR19	22	NC

B - JAY Anterior Thoracic Support (Seating Benefits)

Classic-Structure Pull Down Harness (4 point strap, Only available with JAY backrests)

☐	400K15	Extra-Small (JZ8300)	\$322.00
☐	400K16	Small (JZ8301)	\$322.00
☐	400K17	Medium (JZ8302)	\$322.00

POWERED SEATING

OPTION	Description	Retail
--------	-------------	--------

JAY Lateral Thoracic Swing Away Supports *(Requires Pad, Upholstery and Mount)* (Seating Benefits)

JAY Lateral Thoracic Support Pads

☐	400JW35	Lateral Thoracic (LxH) 100 mm x 127 mm / 4"x 5" Pad - Right	\$73.00
☐	400JW38	Lateral Thoracic (LxH) 100 mm x 127 mm / 4"x 5" Pad - Left	\$73.00
☐	400JW36	Lateral Thoracic (LxH) 100 mm x 152 mm / 4" x 6" Pad - Right	\$73.00
☐	400JW39	Lateral Thoracic (LxH) 100 mm x 152 mm / 4" x 6" Pad - Left	\$73.00
☐	400JW37	Lateral Thoracic (LxH) 152 mm x 100 mm / 6" x 4" Pad - Right	\$73.00
☐	400JW40	Lateral Thoracic (LxH) 152 mm x 100 mm / 6" x 4" Pad - Left	\$73.00

JAY Lateral Thoracic Support Pad Upholstery

☐	400NU01	Premium Leatherette	NC
☐	400NU02	Syncra Stretch	NC

JAY Lateral Thoracic Swing-Away Mount

☐	400JW41	Lateral Thoracic Swing Away Bracket - Right	\$325.00
☐	400JW42	Lateral Thoracic Swing Away Bracket - Left	\$325.00

Headrests *(Requires Pad, Upholstery and Mount)* (Seating Benefits)

Headrest Pad

☐	400HE38	Whitmyer Plush Headrest Pad 203 mm / 8"	\$377.00
☐	400HE39	Whitmyer Plush Headrest Pad 254 mm / 10"	\$377.00
☐	400HE45	Whitmyer Adjust-A-Plush Headrest Pad Small <i>(Narrow) (Adjustable from 6" to 14" Wide x 3" Tall)</i>	\$445.00
☐	400HE46	Whitmyer Adjust-A-Plush Headrest Pad Medium <i>(Adjustable from 7" to 16" Wide x 4" Tall)</i>	\$445.00
☐	400HE48	Omit Headrest Pad	NC

Headrest Mount

☐	400MH33	Whitmyer Linx Headrest Mount with Quick Release	\$247.00
☐	400MH34	Whitmyer Linx-2 Headrest Mount with Quick Release	\$298.00
☐	400MH35	Whitmyer AXYS Headrest Mount with Quick Release	\$247.00
☐	400MH17	Omit Headrest Mount	NC

Seat Pan

☐	400OP08	Small <i>(10"-13"W and 10"-16"D)</i>	NC
☐	400OP09	Standard <i>(12"-16"W and 12"-19"D)</i>	NC

Seat Cushion Type

☐	400ST80	JAY <i>(See Section F)</i>	NC
☐	400ST75	Omit	NC

POWERED SEATING

OPTION	Description	Retail
--------	-------------	--------

Seat Size

Seat Pan Width

☐	400R83	10"	NC
☐	400R84	11"	NC
☐	400R85	12"	NC
☐	400R86	13"	NC
☐	400R87	14"	NC
☐	400R78	15"	NC
☐	400R53	16"	NC

Seat Pan Depth (Front of Seat Pan to Back Structure)

☐	400RF75	10"	NC
☐	400RF76	11"	NC
☐	400RF77	12"	NC
☐	400RF78	13"	NC
☐	400RF52	14"	NC
☐	400RF53	15"	NC
☐	400RF51	16"	NC
☐	400RF45	17"	NC
☐	400RF46	18"	NC
☐	400RF47	19"	NC

C - JAY Seat Cushion (Seating Benefits)

☐	400CU110	Jay ZIP Cushion (Sizes 10" - 16" Wide)	\$711.00
☐	400CU59	JAY Ion (Sizes 14" - 16" Wide)	\$505.00
☐	400CU61	JAY Union (Sizes 14" - 16" Wide)	\$764.00
☐	400CU64	JAY Fusion (Sizes 14" - 16" Wide)	\$1,045.00
☐	400CU104	JAY Fusion with Cryo Technology (Sizes 14"-16" Wide)	\$1,152.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	Description	Retail
--------	-------------	--------

Knee and Thigh Swing Away Support Pads (Seating Benefits)

Lateral Knee Support Pads and Hardware (Must pick pad and mount)

☐	400XY07	Lateral Knee Support Pad, Right (3" round)	\$247.00
☐	400XY08	Lateral Knee Support Pad, Left (3" round)	\$247.00
☐	400XY09	Lateral Knee Support Pad, Right (4" round)	\$247.00
☐	400XY10	Lateral Knee Support Pad, Left (4" round)	\$247.00
☐	400XX10	Lateral Swing Away Bracket Knee - Right	\$325.00
☐	400XX11	Lateral Swing Away Bracket Knee - Left	\$325.00

Lateral Thigh Support Pads and Hardware (Must pick pad and mount)

☐	400KF22	Lateral Thigh Support Pad Short, Right (5" H x 5" L)	\$146.00
☐	400KF23	Lateral Thigh Support Pad Short, Left (5" H x 5" L)	\$146.00
☐	400KF24	Lateral Thigh Support Pad Long, Right (5" H x 7" L)	\$146.00
☐	400KF25	Lateral Thigh Support Pad Long, Left (5" H x 7" L)	\$146.00
☐	400KF26	Lateral Thigh Removable Mount - Right	\$325.00
☐	400KF27	Lateral Thigh Removable Mount - Left	\$325.00

Armrest Structure

☐	400EG44	Standard Cantilever Armrest (Height Adjustable 9" to 12")	NC
☐	400EG45	Low Cantilever Armrest (Height Adjustable 7" to 10")(Available w/ 12" Armpads only)	NC
☐	400EG46	Dual Post Short, Removeable, Flip-Back (Height Adjustable (10" to 15"))	\$434.00
☐	400EG47	Single Post, Removable (N/A w/ thigh support) (Height Adjustable 11" to 15")	\$333.00

Arm Pads

SEDEO Armrest Pad (Available w/ Cantilever Armrest)

☐	400AB43	12", Integral Foam, Right	NC
☐	400AE43	12", Integral Foam, Left	NC
☐	400AB44	16", Integral Foam, Right	NC
☐	400AE44	16", Integral Foam, Left	NC

Waterfall Armpad (Available w/ Dual Post and Single Post Armrest)

☐	400AB01	10" Waterfall Armpad, Right	NC
☐	400AE01	10" Waterfall Armpad, Left	NC
☐	400AB02	14" Waterfall Armpad, Right	NC
☐	400AE02	14" Waterfall Armpad, Left	NC

Side Guards

☐	400SG46	Side Guards Mounted On Seat Rail Bracket (Only with Reclining Armrests) (N/A w/ Thigh Supports)	\$356.00
--------------------------	---------	---	----------

POSITIONING SUPPORTS & ACCESSORIES

OPTION		Description	Retail
--------	--	-------------	--------

Centermount Legrest (400H36 is Client Cost)

☐	400H34	Fixed Centermount Legrest, Short (Includes Calf Pads (12.5" to 15.5" Leg Length)	NC
☐	400H35	Fixed Centermount Legrest, Long (Includes Calf Pads (14.5" to 19.5" Leg Length)	NC
☐	400H123	Angle Adjustable Centermount Legrest, Small (Includes Calf Pads)(3" to 7" Leg Length)	\$438.00
☐	400H124	Angle Adjustable Centermount Legrest, Medium (Includes Calf Pads)(3" to 13" Leg Length)	\$438.00
☐	400H36	Power Centermount Legrest, 90° - 0°, 7.5" Extension (Powered articulation) (3" to 17.5" Leg Length)	\$3,435.00

Angle Adjustable Centermount Footplates

☐	400F32	Dual Angle Adjustable Footplates, Small (5.75" x 6")	\$235.00
☐	400F33	Dual Angle Adjustable Footplates, Large (5.75" x 10")	\$235.00

Power Centermount Footplates

☐	400F127	Dual, Height, Width and Angle Adjustable, Standard (14"W x 10"D)	\$284.00
☐	400F128	Dual, Height, Width and Angle Adjustable, Large (16"W x 12"D)	\$284.00

Power Centermount Footplate Accessories

☐	400DI68	Extra Large Centermount Footplate	\$315.00
☐	400DI69	Footplate Connector Link (Available with Sedeo Fixed Centermount)	\$96.00
☐	400DI57	Gel Footpads	\$401.00

Power Centermount Calf Pads

☐	400DI34	Calf Pads, Pair, Standard Intergral Foam	NC
--------------------------	---------	--	----

Swing-Away Legrest (Pair)

☐	400H196	70° Swing Away (13.5" to 16" Leg Length)	NC
☐	400H197	80° Swing In-Out (Front Mount) (5.5" to 9" Leg Length)	NC
☐	400H198	90° Swing In-Out (Front Mount)(5.5" to 9" Leg Length)	NC
☐	400H51	Manual Elevating Legrests, Swing Away (12.5" to 17" Leg Length)	\$406.00
☐	400H176	Omit Legrests	NC

Swing Away Legrest Footplates

☐	400F115	Small Angle Adjustable Flip Up Footplates (Requires Leg Strap)	\$235.00
☐	400F116	Standard Angle Adjustable Flip Up Footplates (Requires Heel Loop)	\$235.00

Swing Away Legrest Options

☐	400DI09	Leg Strap	NC
☐	400DI10	Heel Loops	NC
☐	400DI53	2" Hanger Extension (7.5" to 11" Leg Length)(Available w/ Swing In-Out)	\$97.00
☐	400DI54	High Mount Extension (7" to 14" Leg Length)(Available w/ 70° Swing Away)	\$97.00

Pelvic Positioning Belts

☐	400A837	Positioning Belt, Carseat Buckle 50mm x 1600mm / 2"	NC
☐	400A839	Positioning Belt, Aircraft Buckle 50mm x 1300mm / 2"	NC
☐	400A840	Positioning Belt, Aircraft Buckle, Padded 50mm x 1300mm / 2"	\$103.00
☐	400A845	Bodypoint - Evoflex Pelvic Stabilizer Belt	\$216.00
☐	400A846	Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	\$271.00

POSITIONING SUPPORTS & ACCESSORIES			
OPTION		Description	Retail
Accessories			
Ventilator Option			
☐	400VT12	Vent Tray Fixed (<i>Pulmonetic LTV Models</i>)	\$609.00
O2 Tank Holder			
☐	400A718	O2 Tank Holder (<i>D and E Styles</i>)	\$391.00
Wheelchair Tray Table			
☐	4007A09	Tray, Swing Away	\$381.00
Backpack (Client Cost)			
☐	400A1	Black Utility	\$112.00
Misc. Accessories (Client Cost)			
☐	4008A59	Transfer Bar (<i>Pair</i>)	\$493.00
☐	4008A52	Ram Mount Phone	\$77.00
☐	4008A53	Ram Mount Tablet	\$308.00
☐	4008A54	Ram Mount Cup Holder	\$77.00
☐	4008A58	Rear View Mirror	\$77.00
☐	400A610	Walker Holder Holder	\$191.00
☐	400A608	Cane/Crutch Holder	\$308.00
☐	4009A10	Utility Hooks	\$58.00
☐	4008A93	Wheelchair Cover	\$240.00

JOYSTICK / ELECTRONICS			
OPTION		Description	Retail
VR-2 Joystick			
☐	400J127	VR-2 Joystick	NC
R-net Joystick			
☐	400J121	Std LED Joystick for R-net Controls	NC
☐	400J116	Color LCD (CJSM1)	NC
☐	400J118	High Resolution Color LCD, (CJSM2+BT) (Client Cost) <i>3.5" Color Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control</i>	\$1,162.00
☐	400J119	Omit Remote Control <i>(Required with arm mounted specialty controls)</i>	NC
Joystick Mounting			
☐	400JM34	Fixed Mount	NC
☐	400JM35	Parallel Swing-Away Mount	\$447.00
☐	400JM36	Link-it Mid Line Mount <i>(Swing Away)</i>	\$447.00
☐	400JM49	Attendant <i>(Available with LED or CJSM1 only, Mounted to Push Handles)</i>	NC
Joystick Mount (Location)			
☐	400HM16	Right Hand armrest	NC
☐	400HM17	Left Hand armrest	NC
Display Language Selection R-net			
☐	400L1	English	NC
☐	400L3	French	NC
☐	400L7	Spanish	NC

JOYSTICK / ELECTRONICS

OPTION	Description	Retail
--------	-------------	--------

Joystick Handle (Optional) N/A with HD Joystick

☐	400JH13	Ball Handle	\$110.00
☐	400JH15	T-Handle	\$110.00
☐	400JH24	Bodypoint - Rubber Dome Joystick Handle	\$216.00
☐	400JH25	Bodypoint - U-Shaped Joystick Handle (3 in)	\$216.00
☐	400JH26	Bodypoint - Grooved Mushroom Joystick	\$216.00

Electronic Programmers

☐	400PR13	R-net PC Wireless Programmer	\$660.00
--------------------------	---------	------------------------------	----------

Power Modules

☐	400ER12	Power Module, VR-2, 90A, Non-Expandable	NC
☐	400ER22	Power Module, R-net, 90A, Non-Expandable (Only Available with Single Power)	NC
☐	400ER13	Power Module, R-net, 90A, Expandable (Required for specialty controls or 2+ actuators)	\$102.00
	400EX11	Expandable 90A Controller - R-net	NC
	400EX12	Harness for 90A Expandable Controller	NC

Seating Control (Client Cost)

Powered Seating Control

☐	400EO40	Omit Thru-drive	NC
☐	400EO37	Thru-Drive Control for 1 Actuator.	\$1,443.00
☐	400EO38	Thru-Drive Control for 2+ Actuators.	\$2,903.00

Switch Control Only (Available with Omit Thru-Drive Only) (1 Actuator for Tilt or Elevate only)

☐	400EO33	Dual Toggle for 1 actuator	NC
☐	400EO34	2.5" Single Button Switch for 1 Actuator (Button Switches are not Mounted)	NC
☐	400EO41	Ctrl+5 Switch Box with Voice Activation for 2 actuators or more	NC

Additional Power Seating Control (Available with Thru-drive Only)

☐	400QG11	Ctrl+5 Switch Box With Voice Activation	\$1,078.00
--------------------------	---------	---	------------

Dual Toggle Mounting Location

☐	400TN1	Toggle Mount, Right	NC
☐	400TN2	Toggle Mount, Left	NC

Ctrl+5 Location (If mounted on same side as joystick, it will Swing Away)

☐	400EY19	Right Armrest	NC
☐	400EY20	Left Armrest	NC
☐	400EY21	Back Cane, Right	NC
☐	400EY24	Back Cane, Left	NC

Ctrl+5 Switch Type

☐	400EY22	-	Push Buttons	NC
☐	400EY23	K0108	Toggle Switches	\$356.00

JOYSTICK / ELECTRONICS			
OPTION		Description	Retail
Attendant Controls (Prior Approval Required)			
☐	4000A05	Dual Control Module for Attendant	\$1,225.00
☐	4000A06	R-net Color LCD (CJSM1), Attendant Control Only	\$1,225.00
Attendant Mounting Location			
☐	400DN5	Right Back Cane	NC
☐	400DN6	Left Back Cane	NC
R-net Programming			
Assignable Buttons			
Assignable Button feature allows for customization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.			
☐	400PM12	Include Assignable buttons	NC
Power Centermount Footplate Extension Option			
☐	400IE03	Include 7.5" Powered footplate extension for transfer assist (Electronics programmed to allow footplates to be extended with or without legrest being raised.) Only Available with Powered Articulation	\$592.00
R-net Omni2 Options (For Expandable Controls, Choose All That Apply)			
☐	400ON14	R-Net Omni2 Only Input Control Module and 3.5" LCD Display w/ Built-In Infrared and Bluetooth Choose Mount Below	\$2,384.00
☐	400ON15	R-Net Omni2 Bundled Input Control Module And 3.5" LCD Display w/ Built-In Infrared and Bluetooth. (Bundled Omni2 available w/ Switch-It Input Devices)	NC
Omni2 Mount (If Used with Standard Joystick, SunDisk, or Bullet Hand Tray, Must Be Mounted Opposite Side)			
☐	400DS21	Swing Away Omni2 Mount w/ Rear Mounted Interface, Right	\$698.00
☐	400DS22	Swing Away Omni2 Mount w/ Rear Mounted Interface, Left	\$698.00

SPECIALTY INPUT DEVICES (AADL W100 Grant)

OPTION	HCPC	Description	Retail
--------	------	-------------	--------

Hand Control Packages (Requires Standard Joystick Omit/Only one hand control device can be ordered) (AADL W100 Grant)**Proportional Input Devices**

☐	400EU14	Switch-It MicroGuide	\$4,574.00
☐	400EU17	Switch-It VersaGuide	\$2,557.00
☐	400EU18	Switch-It VersaGuide EZ	\$3,023.00

Mounting Position (Required)

☐	400EM05	Link-It to the Fixed Joystick Mount	\$416.00
☐	400EM06	Link-It to the Parallel Swing Away Joystick Mount	\$416.00
☐	400EM07	Link-It to the Midline Swing Away	\$1,480.00
☐	400EM08	Omit Mount	NC

Mounting Platforms (Optional) (Only available w/ MicroSeries joysticks)

☐	400MG16	SunDisk Hand Tray	\$416.00
☐	400MG15	Bullet Hand Tray	\$416.00

Mounting Location (Required)

☐	400IM09	Right	NC
☐	400IM10	Left	NC

Switch Options

☐	See Switch Selections (Maximum 2)
--------------------------	-----------------------------------

Head Control Packages (Select headrest mount on next page) (AADL W100 Grant)

Only one head control device can be ordered.

Proportional Input Devices

☐	400EU25	Switch-It DualPro Proportional Head Array - Facial Pads	\$7,057.00
☐	400EU26	Switch-It DualPro Proportional Head Array - Spot Pads	\$7,057.00
☐	400EU51	Vigo Head Drive	\$6,293.00

Switched Input Devices

☐	400EU28	Switch-It 3-Switch Head Array - Short Curved Wings	\$6,154.00
☐	400EU30	Switch-It 3-Switch Head Array - Long Curved Wings	\$6,154.00

Switch Options

☐	See Switch Selections (maximum 2)
--------------------------	-----------------------------------

Sip 'N Puff (AADL W100 Grant)

☐	400ID73	Switch-It Sip 'N Puff Combo Head Array with Link-It Hardware (Requires headrest mount)	\$7,037.00
☐	400ID74	Switch-It OptiStop (Requires Sip 'N Puff)	\$2,741.00

SPECIALTY INPUT DEVICES

OPTION	Description	Retail
--------	-------------	--------

Headrest Mounts *(Applies to Head/Sip 'N Puff control packages)* (Seating Benefits)

☐	400IJ13	Whitmyer Linx Headrest Mount with Quick Release	\$247.00
☐	400IJ14	Whitmyer Linx-2 Headrest Mount with Quick Release	\$298.00
☐	400IJ15	Whitmyer AXYS Headrest Mount with Quick Release	\$247.00

Chin Control Packages (AADL W100 Grant)

Proportional Input Devices

☐	400EU34	Switch-It MicroGuide <i>(Switches mounted to the headrest only)</i>	\$4,574.00
☐	400EU36	Switch-It VersaGuide	\$2,557.00
☐	400EU37	Switch-It VersaGuide EZ	\$3,023.00

Mounting Position

☐	400EM09	Link-It To Chin Harness, Small <i>(150mm / 6")</i>	\$1,760.00
☐	400EM10	Link-It To Chin Harness, Large <i>(220mm / 8")</i>	\$1,760.00
☐	400EM11	Link-It To Headrest <i>(Swing Away)</i>	\$1,760.00
☐	400EM16	Omit Mount	NC

Input Device Mounting Location

☐	400IM13	Right	NC
☐	400IM14	Left	NC

Switch Options

-	See Switch Selections <i>(Minimum 1 for the Power Chin Boom, Maximum 3)</i>	-
---	---	---

Safety Devices (AADL W100 Grants)

☐	400SA06	Switch-It Remote Stop Switch - XLR Connector	\$1,317.00
--------------------------	---------	--	------------

Switched Driving Packages (AADL W100 Grants)

☐	400SW44	Switch-It Cool Cube with DB9 connector for proximity / Proportional switch driving. <i>(Pro Spot Proportional Switch's are required if proportional driving is desired)</i> <i>(For Non-Proportional driving, Single Switches can be used)(For Future Development)</i>	\$3,792.00
--------------------------	---------	--	------------

Specialty Input Devices (AADL W100 Grant)

Switches

Single switches can be programmed to activate any function of the wheelchair including: On/Off, Mode, Profile, Speeds, Seating, Memory Positions, Horn, Lights & Indicators, IR functions, Bluetooth mouse clicks, Environmental Controls when plugged into assignable port.

Switches can connect via: Joysticks (Standard, Micro series, Versa series) OMNI, SJX6 Hub, Cool Cube, Head controls

Switch ordering instructions:

1. Select Switch type and quantity in section "A" using columns below. (Max 2 switches. see limitations below)
2. Select Switch Mounting location from Section "B" *using corresponding switch column.*
3. Select Switch Mounting Side from Section "C" *using corresponding switch column.*

Section A - Switch Selections

Single Function Switches - Mechanical Switches

2.5" Single Switch - Black
 2.5" Single Switch - Blue
 2.5" Single Switch - Green
 2.5" Single Switch - Red
 2.5" Single Switch - Yellow
 Egg Switch - Black
 Egg Switch - Blue
 Egg Switch - Green
 Egg Switch - Red
 Egg Switch - Yellow
 MicroLite (Includes color labels)
 Mini Button 1 - Blue (Velcro)
 Mini Button 1 - Green (Velcro)
 Mini Button 1 - Red (Velcro)
 1" Single Switch - Black
 1" Single Switch - Blue
 1" Single Switch - Green
 1" Single Switch - Orange
 1" Single Switch - Red
 1" Single Switch - Yellow

Switch / Quantity

Switch 1	\$	Switch 2	\$
☐ 400IR31	\$110.00	☐ 400IR56	\$110.00
☐ 400IR32	\$110.00	☐ 400IR57	\$110.00
☐ 400IR33	\$110.00	☐ 400IR58	\$110.00
☐ 400IR34	\$110.00	☐ 400IR59	\$110.00
☐ 400IR35	\$110.00	☐ 400IR60	\$110.00
☐ 400IR36	\$110.00	☐ 400IR61	\$110.00
☐ 400IR37	\$110.00	☐ 400IR62	\$110.00
☐ 400IR38	\$110.00	☐ 400IR63	\$110.00
☐ 400IR39	\$110.00	☐ 400IR64	\$110.00
☐ 400IR40	\$110.00	☐ 400IR65	\$110.00
☐ 400IR41	\$115.00	☐ 400IR66	\$115.00
☐ 400IR42	\$308.00	☐ 400IR67	\$308.00
☐ 400IR43	\$308.00	☐ 400IR68	\$308.00
☐ 400IR44	\$308.00	☐ 400IR69	\$308.00
☐ 400IR45	\$388.00	☐ 400IR70	\$388.00
☐ 400IR46	\$388.00	☐ 400IR71	\$388.00
☐ 400IR47	\$388.00	☐ 400IR72	\$388.00
☐ 400IR48	\$388.00	☐ 400IR73	\$388.00
☐ 400IR49	\$388.00	☐ 400IR74	\$388.00
☐ 400IR50	\$388.00	☐ 400IR75	\$388.00

Electronic Switches

Flex Switch (Headrest mount only)

☐ 400IR54	\$542.00	☐ 400IR79	\$542.00
----------------------------------	----------	----------------------------------	----------

Pro Spot - Proportional Switches

Pro Spot Proportional Switch - Red
 Pro Spot Proportional Switch - Green
 Pro Spot Proportional Switch - Blue
 Pro Spot Proportional Switch - Yellow
 Pro Spot Proportional Switch - Orange

☐ 400IR181	\$542.00	☐ 400IR186	\$542.00
☐ 400IR182	\$542.00	☐ 400IR187	\$542.00
☐ 400IR183	\$542.00	☐ 400IR188	\$542.00
☐ 400IR184	\$542.00	☐ 400IR189	\$542.00
☐ 400IR185	\$542.00	☐ 400IR190	\$542.00

Section B - Switch Mounting Position

Link-It to the Backrest (Up to 1 per side)
 Headrest Mount - Gooseneck & Rigid (up to 2 per side)
 Link-It to the Fixed Joystick Mount (Up to 2 per side)
 Link-It to the Parallel Swing Away Joystick Mount (Up to 2 per side)
 Link-It to the Seat Rail (Up to 2 per side)
 (N/A w/ Knee & Thigh Supports)
 Omit Mount (includes Velcro)

☐ 400EM17	\$1,338.00	☐ 400EM23	\$1,338.00
☐ 400EM18	\$1,760.00	☐ 400EM24	\$1,760.00
☐ 400EM19	\$416.00	☐ 400EM25	\$416.00
☐ 400EM20	\$416.00	☐ 400EM26	\$416.00
☐ 400EM21	\$556.00	☐ 400EM27	\$556.00
☐ 400EM22	NC	☐ 400EM28	NC

Section C - Switch Mounting Location

Right side of chair	☐	400IM17	NC	☐	400IM19	\$698.00
Left side of chair	☐	400IM18	NC	☐	400IM20	\$698.00

Joystick Integrated Switches (Pair)

☐	400JO03	Proximity Switches on Flex Mount, Pair (Used with VersaGuide/VersaGuide EZ)	\$1,403.00
--------------------------	---------	---	------------



1000 Creditstone Road, Unit 2, Concord ON L4K 4P8