



Q300M with SEDEO LITE
AADL - W577A

EIPW40

Account:	
Date:	
PO#:	
Buyer Name:	
Marked For:	


Ship To:	
Name:	
Address:	
City, Prov., PC:	
Phone:	

Sunrise Medical Inc.
1000 Creditstone Road, Unit 2
Concord, ON L4K 4P8
1-800-263-3390



Client Measurements

a)	Top of head to seat	Inch
b)	Top of shoulder to seat	Inch
c)	Underarm to seat	Inch
d)	Elbow to seat	Inch
e)	Top of leg to seat	Inch
f)	Shoulder width	Inch
g)	Chest width	Inch
i)	Hip width	Inch
j)	Actual sitting depth	Inch
m)	Foot length	Inch
(n)	Knee to foot	Inch



User Weight

Notes:

Chair Specifications

		Seat Type	SEDEO LITE
Weights starting at: (with Batteries)	NPC: 250 lbs SPC: 274 lbs	Seat Size Range	Width 15" to 20" Depth 16" to 22"
Power Seating	50° of CG Tilt / 12" seat elevator	STF Height Range	16.5" to 19"

Mid-Wheel Drive Base Specifications

Drive Wheel Position	Mid-wheel Drive	Drive Wheels	12" drive wheels w/flat free inserts (black) 14" drive wheels w/flat free inserts (black)
User Weight Capacity	300 lbs.	Tires/Casters	7" Front / 7" Rear (Black)
Min. Turning Radius	17" (excludes footplates)	Charger	8A Dual Mode (Compact)
Base Size	Width 20.5", Length 33"	Battery Size	50Ah Group 22
Drive Control	R-net 90A / VR-2 90A / LiNX 120A	Color Accents	8 Shroud Colors and 7 Color Accents
Driving Range	Up to 17 miles	Ground Clearance	3"
Drive Motors	4-pole 6mph		

BASE				
OPTION		Description		Retail
Standard Duty Mid Wheel Drive Base - End Item = EIPW40				
☐	400PC294		Quickie Q300 M NPC No Power Option for SEDEO Lite Seat Included in Base Pricing: Q300 M Base, Non-Expandable Electronics, No Power Seat Functions, Six-Wheel Suspension, 12" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	\$8,072.00
☐	400PC295		Quickie Q300 M SPC Single Power Option for SEDEO Lite Seat Included in Base Pricing: Q300 M Base, Non-Expandable Electronics, No Power Seat Functions, Six-Wheel Suspension, 12" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	\$8,072.00
User Weight				
☐	400WG23		Up to 300 lbs.	NC
Motor Package				
☐	400PF99		4-pole 10kph / 6mph	NC
Tracking				
☐	400MI01		SureTrac Module & Motor Encoders <i>(Improves Drive Tracking)</i>	\$1,696.00
☐	400MI04		Traction Control System for LINUX (GYRO Technology)	\$925.00
Drive Wheels				
☐	400DW13		Drive Wheel, 12", Solid, Knobby, Black	NC
☐	400DW58		Drive Wheel, 14", Solid, Knobby, Black	NC
Caster Wheels				
☐	400CW59		Caster Wheels 7", Solid, Black	NC
Shroud Color				
☐	400C164		Black	NC
☐	400C165		Candy Apple Red	NC
☐	400C166		Candy Blue	NC
☐	400C167		Candy Purple	NC
☐	400C176		Green Apple	NC
☐	400C178		Pearl Pink	NC
☐	400C182		Silver	NC
☐	400C125		White	NC
Color Accents				
Note: Accent colors and shroud colors may not match				
☐	400EZ20		Apple Green	NC
☐	400EZ22		Black	NC
☐	400EZ21		Blue	NC
☐	400EZ23		Orange	NC
☐	400EZ24		Red	NC
☐	400EZ19		White	NC
☐	400EZ25		Silver	NC

BASE			
OPTION		Description	Retail
ISO Transit Tie Downs			
☐	400TO50	ISO Transit Package	NC
LED Lights			
☐	400LI08	LED Fender Lights	\$646.00
☐	400LI02	Omit LED Lights	NC
Seat Heights			
Note: Available Seat-to-Floor heights are dependant on Powered Seating. Use the matrix below to determine available ranges for your desired options. Measured from front of seat pan to floor.			
Seat Height			
☐	400HT103	16.5"	NC
☐	400HT55	17"	NC
☐	400HT104	17.5"	NC
☐	400HT56	18"	NC
☐	400HT105	18.5"	NC
☐	400HT57	19"	NC
Batteries & Chargers			
Batteries			
☐	400AG35	GP22/50Ah Flat Top Pair	\$678.00
Omit Batteries.			
☐	400AG24	Battery Harnesses Only (Omit Batteries)	NC
Chargers			
☐	400BJ30	10A Dual Mode	\$162.00
☐	400BJ19	8A Dual Mode	NC

POWERED SEATING			
OPTION		Description	Retail
Power Seating Options			
Seat Elevator			
☐	400PI07	Seat Elevator, Powered, 300 mm / 12"	\$3,477.00
☐	400PI08	C-Me Seat Elevator, Powered, 300 mm / 12"	\$3,477.00
Seat Tilt			
☐	400PA48	Posterior 50° Tilt (W993 for Power Tilt)	\$6,496.00
Pre Set Seat Frame Angle (No Power Seating Options, Increases Seat to Floor Height)			
☐	400FA10	0° (Adds 0" to STFH)	NC
☐	400FA3	3° (Adds 1" to STFH)	NC
SEDEO Lite Backrest Frame Type (Angle Adjustable 86° to 106°)			
Back Frame Width			
☐	400BF21	Width Adjustable, 385 - 485 mm / 15 - 19"	NC
☐	400BF22	Width Adjustable, 460 - 510 mm / 18 - 20"	NC
Back Cane Height			
☐	400B263	Standard, 18" - 20"	NC
☐	400B264	Tall, 21"-23"	NC
☐	400PH30	Height Adjustable, Short w/ Trailing Push Handles, 20" - 22"	NC
☐	400PH31	Height Adjustable, Tall w/ Trailing Push Handles, 23" -26"	NC
Backrest Frame Option			
☐	400BF23	Flip Down, Quick Release	NC
Back Cushion Type			
☐	400BT148	LITE Sling Back, with Tension Adjustable Straps	NC
☐	400BT135	JAY J3	NC
☐	400BT133	OMIT (For Custom and Jay Your Way Mounting. No Headrest Mount)	NC

POWERED SEATING

OPTION	Description	Retail
--------	-------------	--------

A - SEDEO LITE

A - SEDEO Lite Sling Backrest Dimensions

Installed Back Height - (From Pan to Top of Back)

☐	400JR71	18"	NC
☐	400JR72	19"	NC
☐	400JR73	20"	NC
☐	400JR74	21"	NC
☐	400JR75	22"	NC
☐	400JR76	23"	NC

B - JAY J3 Back

B - JAY J3 Back - Contour

☐	400JP01	J3 PA Back for 14"-20" Wide	\$912.00
☐	400JP11	J3 PL Back for 14"-20" Wide	\$986.00
☐	400JP13	J3 PD Back for 14"-20" Wide (18" and 21" shell height only)	\$1,068.00

B - JAY J3 Back - Mounting Hardware

☐	400JP20	TH	NC
--------------------------	---------	----	----

B - JAY J3 Back Shell Options and Height Ranges

JAY J3 shells are mounted with a minimim 4" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field within the ranges shown.

JAY J3 Shell Height

☐	400JR14	16.5" tall shell (Upper Thoracic Short)	NC
☐	400JR15	18" tall shell (Upper Thoracic Med)	NC
☐	400JR16	19.5" tall shell (Upper Thoracic Tall)	NC
☐	400JR17	21" tall shell (Shoulder Height)	NC
☐	400JR27	22.5" tall shell (Shoulder Height)	NC
☐	400JR28	24" tall shell (Shoulder Height)	NC

POWERED SEATING

OPTION	Description	Retail
--------	-------------	--------

JAY J3 Installed Back Height (From Pan to Top of Back)

☐	400JR29	21"	NC
☐	400JR19	22"	NC
☐	400JR20	23"	NC
☐	400JR21	24"	NC
☐	400JR22	25"	NC
☐	400JR23	26"	NC
☐	400JR24	27"	NC
☐	400JR25	28"	NC
☐	400JR30	29"	NC
☐	400JR31	30"	NC

C - JAY Anterior Thoracic Support

Classic-Structure Pull Down Harness (4 point strap, Only available with Jay J3 backrest)

☐	400K10	Medium (ATS8302)	\$325.00
☐	400K11	Large (ATS8303)	\$325.00

JAY Lateral Thoracic Swing-Away Supports (Requires Pad, Upholstery and Mount)

JAY Lateral Thoracic Support Pads

☐	400JW35	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right	\$73.00
☐	400JW38	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left	\$73.00
☐	400JW36	Lateral Thoracic (L x H) 100 mm x 152 mm / 4" x 6" Pad - Right	\$73.00
☐	400JW39	Lateral Thoracic (L x H) 100 mm x 152 mm / 4" x 6" Pad - Left	\$73.00
☐	400JW37	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right	\$73.00
☐	400JW40	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left	\$73.00

JAY Lateral Thoracic Support Pad Upholstery

☐	400NU01	Premium Leatherette	NC
☐	400NU02	Stretch	NC

JAY Lateral Thoracic Swing-Away Mount

☐	400JW41	Lateral Thoracic Swing Away Bracket - Right	\$325.00
☐	400JW42	Lateral Thoracic Swing Away Bracket - Left	\$325.00

POWERED SEATING

OPTION	Description	Retail
--------	-------------	--------

Headrests *(Requires Pad, Upholstery and Mount)*

Headrest Pad

☐	400HE36	SEDEO Pad Small With Cover (L x H) 240 x 150 mm / 10" x 6" <i>(Only Available with SEDEO Multi-Position Headrest Mount and SEDEO Upholstery)</i>	\$227.00
☐	400HE38	Whitmyer Plush Headrest Pad 203 mm / 8"	\$377.00
☐	400HE39	Whitmyer Plush Headrest Pad 254 mm / 10"	\$377.00
☐	400HE44	Whitmyer Plush Headrest Pad 356 mm / 14"	\$377.00
☐	400HE46	Whitmyer Adjust-A-Plush Headrest Pad Medium <i>(Adjustable from 7" to 16" Wide x 4" Tall)</i>	\$445.00
☐	400HE48	Omit Headrest Pad	NC

Headrest Pad Upholstery

☐	400NU10	SEDEO Upholstery	NC
☐	400NU11	Whitmyer Lycra / Stretch	NC

Headrest Mount

☐	400MH32	SEDEO Multi-Position Headrest Mount with Quick Release <i>(Only available with SEDEO Pad)</i>	\$557.00
☐	400MH33	Whitmyer Linx Headrest Mount with Quick Release	\$247.00
☐	400MH34	Whitmyer Linx-2 Headrest Mount with Quick Release	\$298.00
☐	400MH35	Whitmyer AXYS Headrest Mount with Quick Release	\$247.00
☐	400MH17	Omit Headrest Mount	NC

Seat Cushion Type

☐	400ST80	JAY <i>(See Section F)</i>	NC
☐	400ST75	Omit	NC

Seat Size

Seat Pan Width

☐	400R78	15"	NC
☐	400R53	16"	NC
☐	400R54	17"	NC
☐	400R55	18"	NC
☐	400R56	19"	NC
☐	400R57	20"	NC

Seat Pan Depth *(Front of Seat Pan to Back Structure - Effective Depth)*

☐	400RF51	16"	NC
☐	400RF45	17"	NC
☐	400RF46	18"	NC
☐	400RF47	19"	NC
☐	400RF48	20"	NC
☐	400RF49	21"	NC
☐	400RF50	22"	NC

POWERED SEATING

OPTION	Description	Retail
--------	-------------	--------

C - JAY Seat Cushion

☐	400CU59	JAY Ion (Sizes 16" - 20" Wide)	\$505.00
☐	400CU61	JAY Union (Sizes 16" - 20" Wide)	\$764.00
☐	400CU64	JAY Fusion (Sizes 16" - 20" Wide)	\$1,045.00

Knee and Thigh Swing Away Support Pads

Lateral Knee Support Pads and Hardware (Must pick pad and mount)

☐	400XY05	Lateral Knee Support Pad, Right (4" round)	\$247.00
☐	400XY06	Lateral Knee Support Pad, Left (4" round)	\$247.00
☐	400XX08	Lateral Swing Away Bracket Knee - Right	\$325.00
☐	400XX09	Lateral Swing Away Bracket Knee - Left	\$325.00

Lateral Thigh Support Pads and Hardware (Must pick pad and mount)

☐	400KF22	Lateral Thigh Support Pad Short, Right (5" H x 5" L)	\$146.00
☐	400KF23	Lateral Thigh Support Pad Short, Left (5" H x 5" L)	\$146.00
☐	400KF24	Lateral Thigh Support Pad Long, Right (5" H x 7" L)	\$146.00
☐	400KF25	Lateral Thigh Support Pad Long, Left (5" H x 7" L)	\$146.00
☐	400KF26	Lateral Thigh Removable Mount - Right	\$325.00
☐	400KF27	Lateral Thigh Removable Mount - Left	\$325.00

Armrest Structure

☐	400EG41	Reclining, Flip Back, Short (Best used with shorter seat depths) (Height Adjustable 9" to 15")	NC
☐	400EG36	Reclining, Flip Back with 2-point post locking arm (Height Adjustable 11" to 15" x .625")	\$351.00
☐	400EG34	Single Post SEDEO Lite, Removable (N/A w/ thigh support)(Height Adjustable 10" to 15")	NC
☐	400EG50	Dual Post, Flip Back, Removeable (Height Adjustable)	NC

Arm Pads

SEDEO Armrest Pad

☐	400AB62	300 mm / 12", Integral Foam, Right	\$54.00
☐	400AE62	300 mm / 12", Integral Foam, Left	\$54.00
☐	400AB63	400 mm / 16", Integral Foam, Right	\$54.00
☐	400AE63	400 mm / 16", Integral Foam, Left	\$54.00
☐	400AB61	Lite Pad, Right	NC
☐	400AE61	Lite Pad, Left	NC

Dual Post Armrest Pad

☐	400AB74	250 mm / 10", Classic Desk Foam, Right	NC
☐	400AE74	250 mm / 10", Classic Desk Foam, Left	NC
☐	400AB76	250 mm / 10", Waterfall Desk, Right	NC
☐	400AE76	250 mm / 10", Waterfall Desk, Left	NC
☐	400AB75	250 mm / 14", Classic Full Foam, Right	NC
☐	400AE75	250 mm / 14", Classic Full Foam, Left	NC
☐	400AB77	250 mm / 10", Waterfall Full, Right	NC
☐	400AE77	250 mm / 10", Waterfall Full, Left	NC

Comfort Arm Pad (Foam)

☐	400AB46	2" Wide Right	\$113.00
☐	400AE46	2" Wide Left	\$113.00

Comfort Arm Pad (Gel)

☐	400AB48	2" Wide Right	\$213.00
☐	400AE48	2" Wide Left	\$213.00

Ottobock Arm Troughs

☐	400AE06	Ottobock Arm Trough, Med (12"), Left	\$139.00
☐	400AB06	Ottobock Arm Trough, Med (12"), Right	\$139.00
☐	400AE07	Ottobock Arm Trough, Large (14"), Left	\$139.00
☐	400AB07	Ottobock Arm Trough, Large (14"), Right	\$139.00
☐	400OL01	Ottobock Hand Pad, Med (7"), Left	\$70.00
☐	400OR01	Ottobock Hand Pad, Med (7"), Right	\$70.00
☐	400OL02	Ottobock Hand Pad, Large (9"), Left	\$70.00
☐	400OR02	Ottobock Hand Pad, Large (9"), Right	\$70.00
☐	400AE78	Ottobock 1 Piece Arm Trough, (18"), Left	\$198.00
☐	400AB78	Ottobock 1 Piece Arm Trough, (18"), Right	\$198.00

Elbow Support

☐	4003A02	Elbow Support Pad, RH 3" x 10" <i>(Includes side support)</i>	\$203.00
☐	4003A03	Elbow Support Pad, LH 3" x 10" <i>(Includes side support)</i>	\$203.00
☐	4003A06	Elbow Support Pad, RH 3" x 3"	\$167.00
☐	4003A07	Elbow Support Pad, LH 3" x 3"	\$167.00
☐	4003A04	Swing Away Elbow Support Hardware RH	\$279.00
☐	4003A05	Swing Away Elbow Support Hardware LH	\$279.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	Description	Retail
--------	-------------	--------

Centermount Legrest

☐	400H184	Fixed Centermount, Short <i>(Includes Calf Pads)(11.5" to 14.5" Leg Length)</i>	NC
☐	400H185	Fixed Centermount, Long <i>(Includes Calf Pads)(13.5" to 18.5" Leg Length)</i>	NC
☐	400H136	Angle Adjustable Centermount <i>(Includes Calf Pads and Single Footplate)(11.5" to 18.5" Leg Length)</i>	\$633.00

Power Centermount Legrest

☐	400H109	SEDEO Power Centermount, 90° - 0° Angle, 7.5" Extension <i>(Powered Articulation) (3" to 17.5" Leg Length)</i>	\$3,435.00
--------------------------	---------	--	------------

Centermount Legrest Accessories

SEDEO Power Centermount Footplates

☐	400F127	Dual, Height, Width and Angle Adjustable, Standard <i>(14"W x 10"D)</i>	\$284.00
☐	400F128	Dual, Height, Width and Angle Adjustable, Large <i>(16"W x 12"D)</i>	\$284.00

SEDEO Power Footplate Accessories

☐	400DI68	Extra Large Centermount Footplate	\$315.00
☐	400DI69	Footplate Connector Link <i>(Available with Sedeo Fixed Centermount)</i>	\$96.00
☐	400DI57	Gel Footpads	\$401.00

Power Centermount Calf Pads

☐	400DI34	Calf Pads, Pair, Standard Integral Foam	NC
☐	400DI35	Calf Pads, Pair, Padded <i>(N/A w/ Std footplate)</i>	\$356.00

SEDEO Swing-Away Legrest (Pair)

☐	400H177	Legrests, Swing Away, 70°, SEDEO PRO <i>(14" to 21" Leg Length)</i>	NC
☐	400H178	Legrests, Manual Elevating, Swing Away, SEDEO PRO <i>(Includes Padded Calf Pads) (14".5 to 21.5" Leg Length)</i>	\$380.00
☐	400H179	Legrests, Power Elevating, Swing Away, SEDEO PRO <i>(Includes Padded Calf Pads) (15.5" to 22.5" Leg Length)</i>	\$3,435.00
☐	400H176	Omit Legrests	NC

Swing Away Legrest Footplates

☐	400F369	Angle Adjustable Flip Up Footplates	\$284.00
--------------------------	---------	-------------------------------------	----------

Swing Away Legrest Calf Supports *(Sedeo Pro Only)*

☐	400DI42	Calf Support, Depth Adjustable, Padded <i>(6.5" to 8" Wide, Depending On Chair Width, x 9.5" Tall)</i>	\$512.00
--------------------------	---------	--	----------

Swing Away Legrest Options

☐	400DI44	Heel Loops	NC
☐	400KG4	Knee Pads, Integral Foam <i>(2.5" Wide x 4.5" Tall)</i>	\$282.00

SEDEO Lower Extremity Swing Away Support

☐	400VD01	Lower Extremity Support, Padded, Swing Away, SEDEO, Right	\$512.00
☐	400VD02	Lower Extremity Support, Padded, Swing Away, SEDEO, Left	\$512.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	Description	Retail
--------	-------------	--------

Pelvic Positioning Belts

☐	400A837	Positioning Belt, Autostyle Buckle 50mm x 1600mm / 2"	NC
☐	400A839	Positioning Belt, Aircraft Buckle 50mm x 1300mm / 2"	NC
☐	400A840	Positioning Belt, Aircraft Buckle, Padded 50mm x 1300mm / 2"	\$103.00
☐	400A845	Bodypoint - Evoflex Pelavic Stabilizer Belt	\$216.00
☐	400A846	Bodypoint - Evoflex Pelavic Stabilizer Belt with Extension	\$271.00

Accessories

Ventilator Option

☐	400VT12	Vent Tray Fixed (<i>Pulmonetic LTV Models</i>)	\$609.00
--------------------------	---------	--	----------

O2 Tank Holder

☐	400A718	O2 Tank Holder Mounted to Backrest (<i>D and E Styles</i>)	\$391.00
--------------------------	---------	--	----------

Wheelchair Tray Table

☐	4007A09	Tray, Swing Away	\$381.00
--------------------------	---------	------------------	----------

Backpack

☐	400A1	Black Utility	\$112.00
--------------------------	-------	---------------	----------

Medial Thigh Support

☐	4008A80	Medial Thigh Support	\$244.00
☐	4008A81	Pommel Swing Away Hardware	\$358.00

Misc. Accessories

☐	4008A59	Transfer Bar (<i>Pair</i>)	\$493.00
☐	4008A52	Ram Mount Phone	\$77.00
☐	4008A53	Ram Mount Tablet	\$308.00
☐	4008A54	Ram Mount Cup Holder	\$77.00
☐	4008A58	Rear View Mirror	\$77.00
☐	400A611	Crutch Holder	\$115.00
☐	4008A73	Bag For Armrest	\$115.00
☐	4009A10	Utility Hooks	\$58.00
☐	4008A93	Wheelchair Cover	\$240.00

JOYSTICK / ELECTRONICS				
OPTION		Type	Description	Retail
Electronics Package				
☐	400EL01		VR-2 or R-net Electronics	NC
☐	400EL02		LiNX Electronics	NC
VR-2 Joystick (Remote)				
☐	400J127	VR-2	VR-2 Joystick	NC
R-net Joystick (Remote)				
☐	400J121	R-net	Std LED Joystick for R-net Controls	NC
☐	400J116	R-net	Color LCD (<i>CJSM1</i>)	NC
☐	400J118	R-net	High Resolution Color LCD, (<i>CJSM2+BT</i>) <i>3.5" Color Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control</i>	\$1,162.00
☐	400J119	R-net	Omit Remote Control (<i>Required with arm mounted specialty controls</i>)	NC
LiNX Joystick				
☐	400J143	LiNX	REM110 - Drive Only Remote Joystick <i>Adaptive Load Compensation</i>	NC
☐	400J144	LiNX	REM200 Series - LED Joystick <i>Adaptive Load Compensation, LED icon-based menu, Two Rocker Switches</i>	NC
☐	400J145	LiNX	REM400 - Touch Screen Color Display Joystick <i>3.5" Color Display, Highly Configurable Touch Interface, Two Multi-Purpose Buttons, Two Stereo Jacks, Audible Cues, Bluetooth, & Adaptive Load Compensation</i>	\$1,339.00
☐	400J146	LiNX	REM400 - Touch Screen Color Display Joystick with Toggle Switches <i>3.5" Color Display, Highly Configurable Touch Interface, Two Multi-Purpose Buttons, Programmable Toggle Switches, Two Stereo Jacks, Audible Cues, Bluetooth, & Adaptive Load Compensation</i>	\$1,695.00
☐	400J147	LiNX	Omit Remote Control (<i>Required with arm mounted specialty controls</i>)	NC
Joystick Mounting				
☐	400JM34		Fixed on Armrest	NC
☐	400JM35		Parallel Swing Away Capsulated on armrest	\$447.00
☐	400JM36		Link-It Mid Line Mount on Armrest (<i>Swing Away</i>)	\$447.00
Joystick Mount (Location)				
☐	400HM16		Right Hand armrest	NC
☐	400HM17		Left Hand armrest	NC
Display Language Selection R-net				
☐	400L1		English	NC
☐	400L3		French	NC
☐	400L7		Spanish	NC
Joystick Handle (Optional) N/A with HD Joystick				
☐	400JH13	Both	Ball Handle	\$110.00
☐	400JH15	R-net	T-Handle	\$110.00
☐	400JH24	Both	Bodypoint - Rubber Dome Joystick Handle	\$216.00
☐	400JH25	Both	Bodypoint - U-Shaped Joystick Handle (<i>3 in</i>)	\$216.00
☐	400JH26	Both	Bodypoint - Grooved Mushroom Joystick	\$216.00
☐	400JH30	LiNX	Chin Cup	\$179.00
☐	400JH31	LiNX	Golf Ball	\$179.00
Electronic Programmers				
☐	400PR13		R-net PC Wireless Programmer	\$660.00
☐	400PR14		LiNX Access Key (LAK)	\$240.00
Power Modules				
☐	400ER12	VR-2	Power Module, VR-2, 90A, Non-Expandable	NC
☐	400ER13	R-net	Power Module, R-net, 90A, Expandable (<i>Required for specialty controls or 2+ actuators</i>)	\$102.00
	400EX13		Expandable 90A Controller - R-net	\$0.00
	400EX14		Harness for 90A Expandable Controller - R-net	\$0.00
☐	400ER24	LiNX	Power Module, LiNX, 120A, Non-Expandable	NC
☐	400ER25	LiNX	Power Module, LiNX, 120A, Expandable (<i>Required for specialty controls or 2+ actuators</i>)	\$1,760.00
	400EX19		Expandable 120A Controller - LiNX	\$850.00
	400EX20		Harness for 120A Expandable Controller - LiNX	\$695.00

JOYSTICK / ELECTRONICS

OPTION	Type	Description	Retail
--------	------	-------------	--------

Seating Control

Powered Seating Control

☐	400EO40	Omit Thru-drive (Requires Switch Control)	NC
☐	400EO37	Thru-Drive Control for 1 Actuator.	\$1,443.00
☐	400EO38	Thru-Drive Control for 2+ Actuators.	\$2,903.00

Switch Control Only (Available with Omit Thru-Drive Only) (Available with Tilt or Elevate Only)

☐	400EO33	R-net	Dual Toggle for 1 Actuator	NC
☐	400EO34	Both	2.5" Single Button Switch for 1 Actuator (Button Switches are not Mounted)	NC
☐	400EO41	R-net	Ctrl+5 Switch Box with Voice Activation for 2 Actuators or more	NC
☐	400EO42	LiNX	10 Way Switch Box for Actuator Control	NC

Additional Power Seating Control (Available with Thru-drive Only)

☐	400QG11	R-net	Ctrl+5 Switch Box with Voice Activation	\$1,078.00
☐	400QG12	R-net	Ctrl+5 Switch Box with Voice Activation (Multi-Power Only)	NC
☐	400QG16	LiNX	10 Way Switch Box for Actuator Control	NC

Dual Toggle Mounting Location

☐	400TN1	Toggle Mount, Right	NC
☐	400TN2	Toggle Mount, Left	NC

Power Seating Control Mounting Location (Ctrl+5 or 10 Way) (If mounted on same side as joystick, it will Swing Away)

☐	400EY19	Right Armrest	NC
☐	400EY20	Left Armrest	NC
☐	400EY21	Back Cane, Right	\$0.00
☐	400EY24	Back Cane, Left	\$0.00

Ctrl+5 Switch Type

☐	400EY22	R-net	Push Buttons	NC
☐	400EY23	R-net	Toggle Switches	\$356.00

Attendant Controls

☐	4000A05	R-net	Dual Control Module for Attendant	\$1,225.00
☐	4000A06	R-net	R-net Color LCD (CJSM1), Attendant Control Only	\$1,225.00
☐	4000A11	LiNX	LiNX Attendant Control Unit (ACU200)	\$1,195.00

Attendant Mounting Location

☐	400DN5	Right Back Cane	NC
☐	400DN6	Left Back Cane	NC

R-net Programming Selection

Assignable Buttons

Assignable Button feature allows for customization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

☐	400PM12	Include Assignable buttons	NC
--------------------------	---------	----------------------------	----

Power Centermount Footplate Extension Option

☐	400IE03	Include 7.5" Powered footplate extension for transfer assist (Electronics programmed to allow footplates to be extended with or without legrest being raised.) Only Available with Powered Articulation	\$592.00
--------------------------	---------	---	----------

Memory Seating

SitSmarter Memory Seating allows for positions to be saved and then can be accessed with a single command, moving all related actuators at once to the preset location. Positions can be locked for clinical goals and have others unlocked for easily adjustable positions for comfort and functionality that can be adjusted by the user. Angles can be shown on screen to provide live feedback.

☐	400PM19	LiNX	SitSmarter Memory Seating Positions	NC
--------------------------	---------	------	-------------------------------------	----

SitSmarter not compatible with Recline Only



1000 Creditstone Road, Unit 2, Concord ON L4K 4P8

Customer Service: (800) 2363-3390, Email: cscanada@sunmed.com, www.sunrisemedical.ca . Specifications are subject to change without notice. © 10/2026 Sunrise Medical Canada Inc. MK-129794-CAN Rev.10