

JAY[®]

SUREFIT™

≡ Simply Custom Seating



JAY SUREFIT

The JAY SureFit custom seating solution offers a simple and smart option mix to meet a wide range of client requirements and dealer reimbursement needs. It's easy to order, offering one of the first online configurators available for the custom seating market.

Simple Option Mix

- Options you need and know in an easy to navigate, easy to order format.

≡ BACKREST SHAPES

Standard Planar



Contour Modification



I-Back Modification



T-Back Modification



Biangular Modification



Grid Back



Lumbar



≡ SEAT SHAPES

Standard Planar



Asymmetrical Seat



Antithrust Modification



Antithrust with Leg Troughs Modification



≡ ATTACHING HARDWARE



Snap Tite

Hook and latch surround the tubing completely providing a solid, secure mounting. To use, simply drop the seat or back onto the rail and latch into place.



Quickie TS Snap 2 Pin

Durable and economical attaching system that clips around the back or seat rails. No height adjustment capabilities.



Universal Hardware

Adjusts 3" above rail and saddle cutout to mount on double seat rail tubes.



J3 Back Style

Lightweight 2-point attachment hardware that's great for reduced mounting space and need for anterior adjustment capabilities. Only available with back.



Adjustable J&L

Combination can accommodate both angle adjustability and fixed 90° orientation.

WHITMYER HEAD SUPPORTS



Plush Pad
Dual-density foam with contoured, curved base. Available in 6", 8", 10", 14" or 19" widths.



Contoured Cradle Pad
Single layer of foam with contoured, curved base. Available in pediatric, small, medium and large sizes.



Adjust-A-Plush Pad
Lateral opening adjusts at two tapered locking joints under zippered cover. Available in Standard or Narrow widths.



Onyx Mount
Available with fixed, detachable or detachable w/offset vertical bar bracket. Adjustable at vertical/horizontal rod interface.



Axys Mount
Available with detachable bracket only. Adjustable at vertical/horizontal rod interface.

UPHOLSTERY

- Rubatex (black only), Vinyl (Naugahyde), Dartex or Reverse Dartex

Vinyl (Naugahyde) Upholstery Color Options



Note: Print processes do not always accurately depict color.

POSITIONING SUPPORTS

Available with variety of bracket styles, contours and sizes dependent on support.



Lateral Thoracic Supports
Fixed bracket,
curved pad shown



Lateral Pelvis and Thigh Supports
Fixed bracket shown.
Removable also available.



Medial Thigh Supports
Multi-ply solid core and
1/2" pressure relieving
foam. Available in oval or
wedge.



Shoe Holders
ABS Shown. Also
available with Tendon
Relief and/or padded
straps.



Anterior Trunk Supports
Classic Comforfit shown.
Also available with Contour
Comforfit. Both available
in structured or dynamic
style

Online Configurator

- Place an order 24/7
- Determine pricing immediately
- Place the seating on the same PO as a mobility base

Single Location Assembly of Seating with Mobility Base

- Ensures compatibility before delivery
- Efficient delivery time

Jay SureFit
Products configured online are for standard options only. For any customization please call customer service.

Specs Back Backrest Seat Seatrest Head Supports Lateral Lateralrest Lateralrest Cost Accessories

Chair Model: None

State Option: Ca, Oh, Va

Wheelchair Specs: WHEELCHAIR SPECS

Custom Embroidery: Custom Embroidery

Embroidery Print: Curase Lowercase, Block/Print Uppercase, Block/Print Lowercase

Embroidery Color: None

T-Nuts: Do not mark t-nut locations

Insert: Insert with chair, Insert with transit chair, Insert w/o chair

Chair Width: 10 inch, 11 inch, 12 inch, 13 inch, 14 inch, 15 inch, 16 inch, 17 inch, 18 inch, 19 inch, 20 inch, 21 inch, 22 inch, Other

Seat Rail Dimensions: 7/8 inch, 1 inch

N/A - Price Dependent On Other Selections. N/C - No Charge

VIEW SELECTIONS INTERNATIONAL OPTIONS

SPECIFICATIONS

JAY[®]

HCPCS Code	E2609 & E2617
Product Weight	Dependent on system size and componentry
Seat Width	Minimum 10", all other dimensions variable
Foam Options	High Resiliency, Sunmate, Pudgee
Cover Materials	Rubatex, Nagahyde, Dartex, Reverse Dartex
Product Configuration	Seat Only, Back Only, Seat and Back – all with available accessories
User Weight Capacity	Non-transit: 250 lbs. Non-transit HD: 300 lbs.
Transit	Minimum 48 lbs. – Max. 200 lbs.
Warranty	Sunrise Limited Warranty