



Quickie ACCESS
July 2025
EIR14

Mark For:	Submitting for: ☐ Quote ☐ Order
Date: _____	PO#: _____
Dealer Acct #: _____	ADDITIONAL SHIPPING INFORMATION
Dealer: _____	
Dealer Contact: _____	
Dealer Address: _____	
Dealer City: _____ Prov: _____ PC: _____	
Dealer Phone: () _____ Fax: () _____	
Confirmation Email: _____	
Confirm Via: ☐ Fax ☐ Email	Ship To: _____
	Attention: _____
	Address: _____
	Address: _____
	Ship To City: _____ Prov _____ PC: _____
	Ship To Phone: () _____ Fax: () _____

✓	Option #	Definition	Price	✓	Option #	Definition	Price
1. FRAME							

Base Model			
☐	208WM105	Quickie Access	\$3,953.00

Transit			
☐	208TO45	Transit Option	\$374.00

Pre-Set Tilt Angle			
Certain wheel combinations will limit the seat angle			
☐	208TA11	-5° to 20° Degree Tilt	NC
☐	208TA12	0° to 25° Degree Tilt	NC

Frame Options			
☐	208FS23	Standard Base and Seat Frame (Maximum Weight: 300 lbs)	STD

Tilt Release Options			
☐	208TT8	Foot Release NA with 19"-22" Depth and Flip Back HA Armrest.	NC
		☐ Left ☐ Right (Facing from rear)	
☐	208TT12	Backrest trigger cable release	NC

User's Weight			
Gas spring selection depends User's weight selection.			
☐	208WG26	Up to 150 lbs	NC
☐	208WG27	151 - 180 lbs	NC
☐	208WG28	181 - 210 lbs	NC
☐	208WG29	211 - 240 lbs	NC
☐	208WG30	241 - 270 lbs	NC
☐	208WG31	271 - 300 lbs	NC




Option #	Definition	Price	Option #	Definition	Price		
2. SIZE / SEATING							
Frame Width			Seat Type				
21" & 22" NA w/ removable stroller handle extension or back upholsteries.							
Small Frame 14" - 17"							
☐	208W5	14"	NC	☐	208ST35	Standard Seat Pan For growth ranges, see Figure 1.5	\$234.00
☐	208W6	15"	NC	☐	208ST54	Omit seat Pan	NC
☐	208W7	16"	NC				
☐	208W8	17"	NC				
Large Frame 18" - 22"							
☐	208W9	18"	NC				
☐	208W10	19"	NC				
☐	208W11	20"	\$308.00				
☐	208W12	21"	\$308.00				
☐	208W13	22"	\$308.00				
Seat Depth			JAY Cushions (Optional)				
See Figure 1 Dimension A for measurement			To order any JAY cushions on-chair, please see the JAY price list.				
Foot release tilt option NA w/19"-22"			Cushion name or part number must be filled out with desired width and depth.				
☐	208SD7	14"	NC	☐	208JC20	JAY Bundled Cushion Cushion Name/Part #:	Price List
☐	208SD8	15"	NC	Width: Depth:			
☐	208SD9	16"	NC				
Medium 16" - 18"							
☐	208SD9	16"	NC				
☐	208SD10	17"	NC				
☐	208SD11	18"	NC				
Long 18" - 20"							
☐	208SD11	18"	NC				
☐	208SD12	19"	NC				
☐	208SD13	20"	\$308.00				
Extra Long 20" - 22"							
☐	208SD13	20"	\$308.00				
☐	208SD14	21"	\$308.00				
☐	208SD33	22"	\$308.00				

Custom Seating			
☐	208QC02	JAY Surefit	See OF
☐	208QC03	JAY Fit	See OF

Figure 1

Frame Width	Seat Pan Growth
14"	14" - 16"
15"	15" - 16"
16"	16" - 18"
17"	17" - 18"
18"	18" - 20"
19"	19" - 20"
20"	20" - 22"
21"	21" - 22"
22"	NA



Frame Width	Seat Pan Growth
14"	14" - 16"
15"	15" - 16"
16"	16" - 18"
17"	17" - 18"
18"	18" - 20"
19"	19" - 20"
20"	20" - 22"
21"	21" - 22"
22"	NA

3. CASTERS / SEAT HEIGHTS / REAR WHEELS

Front Seat Height							
12"-12.5" with 4" casters are limited to a 1.5" floor clearance.							
Rear Wheel Size	CASTER SIZE						
	4"	5"	6"	7"	8" x 1"	8"x 1.5"	8" x 2"
12" - 24"	☐ 12.0"	☐ 12.0"	☐ 12.0"				
	☐ 12.5"	☐ 12.5"	☐ 12.5"				
	☐ 13.0"	☐ 13.0"	☐ 13.0"	☐ 13.0"			
	☐ 13.5"	☐ 13.5"	☐ 13.5"	☐ 13.5"			
	☐ 14.0"	☐ 14.0"	☐ 14.0"	☐ 14.0"			
	☐ 14.5"	☐ 14.5"	☐ 14.5"	☐ 14.5"		☐ 14.5"	☐ 14.5"
	☐ 15.0"	☐ 15.0"	☐ 15.0"	☐ 15.0"	☐ 15.0"	☐ 15.0"	☐ 15.0"
	☐ 15.5"	☐ 15.5"	☐ 15.5"	☐ 15.5"	☐ 15.5"	☐ 15.5"	☐ 15.5"
	☐ 16.0"	☐ 16.0"	☐ 16.0"	☐ 16.0"	☐ 16.0"		
	☐ 16.5"	☐ 16.5"	☐ 16.5"	☐ 16.5"		☐ 16.5"	☐ 16.5"
	☐ 17.0"	☐ 17.0"	☐ 17.0"	☐ 17.0"	☐ 17.0"	☐ 17.0"	☐ 17.0"
	☐ 17.5"	☐ 17.5"	☐ 17.5"	☐ 17.5"	☐ 17.5"	☐ 17.5"	☐ 17.5"
	☐ 18.0"	☐ 18.0"	☐ 18.0"	☐ 18.0"	☐ 18.0"	☐ 18.0"	
	☐ 18.5"	☐ 18.5"	☐ 18.5"	☐ 18.5"	☐ 18.5"	☐ 18.5"	
	☐ 19.0"	☐ 19.0"	☐ 19.0"	☐ 19.0"	☐ 19.0"		
	☐ 19.5"	☐ 19.5"	☐ 19.5"	☐ 19.5"			
	☐ 20.0"	☐ 20.0"	☐ 20.0"	☐ 20.0"		☐ 20.0"	

Rear Tire & Wheel Matrix				
Tire	Rear Wheel			
	5 Spoke Mag	Spoke	Mag Wheel	Drum Brake
Pneumatic	20,22,24	24	12,16	24
Pneumatic w/ Airless Insert	20,22,24	24	12,16	24
Full Poly	20,22,24	24	N/A	24
Low Poly	N/A	N/A	12	N/A
12" = 203mm 16" = 305mm 20" = 451mm 22" = 501mm 24" = 540mm				

Axle Plate			
☐	208AL33	Standard Axle Plate	NC
Casters			
☐	208CW3	4" Low Profile Polyurethane N/A w/ HD package	NC
☐	208CW72	4" x 1" Aluminum Polyurethane	\$129.00
☐	208CW31	4" Micro Lighted	\$82.00
☐	208CW53	5" x 1" Polyurethane	NC
☐	208CW7	5" x 1.5" Semi Pneumatic	\$100.00
☐	208CW45	6" x 1" Polyurethane	NC
☐	208CW35	6" x 1.5" Semi Pneumatic	\$100.00
☐	208CW16	7" Polyurethane	NC
☐	208CW46	8" x 1" Polyurethane	NC
☐	208CW39	8"x 1.5" Semi Pneumatic	\$100.00
☐	208CW37	8"x 2" Pneumatic w/Insert	\$263.00
Rear Wheel			
☐	208RW13	12" Mag Wheel	NC
☐	208RW14	16" Mag Wheel	NC
☐	208RW92	24" Lite Spoke Wheel	NC
☐	208RW103	20" Mag Wheel 5 Spoke	NC
☐	208RW104	22" Mag Wheel 5 Spoke	NC
☐	208RW102	24" Mag Wheel 5 Spoke	NC
☐	208RW150	24" Drum brake wheel - 36 spokes	NC
Rear Wheel Tire			
☐	208RT3	Full Polyurethane	NC
☐	208RT4	Low Polyurethane	NC
☐	208RT1	Pneumatic	NC
☐	208RT2	Pneumatic Airless Insert	\$100.00

✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

3. CASTERS / SEAT HEIGHTS / REAR WHEELS (CONT)

Handrim			
<i>12" and 16" wheels not available with hand rims. See matrix for availability</i>			
☐	208HR1	Aluminum Anodized	NC
☐	208HR4	Plastic Coated	\$146.00
☐	208HR22	Natural Fit - Standard Grip	\$374.00
☐	208HR23	Natural Fit - Super Grip	\$434.00
☐	208HR39	Ellipse 3R	\$550.00
☐	208HR40	Carbolife - Gekko	\$575.00
☐	208HR41	Carbolife - Curve L Tetra Grip	\$625.00
☐	208HR7	Omit Handrim	NC

Axle			
☐	208AX1	Quick Release Axles	NC
☐	208AX17	Threaded Axles <i>N/A w/ 24" spoke wheels & 5 Spoke Mags</i>	NC

Rear Wheel & Handrim Matrix

Handrim	Rear Wheel		
	5 Spoke Mag	Spoke	Drum Brake
Aluminum Anodized	20,22,24	24	24
Plastic Coat	20,22,24	24	24
Natural Fit - Standard Grip	24	24	
Natural Fit - No Thumb Grip	24	24	
Natural Fit - Super Grip	24	24	
Natural Fit LT - Standard Grip	24	24	
Natural Fit LT - Super Grip	24	24	
Natural Fit LT - No Thumb Grip	24	24	

4. WHEEL LOCKS

Wheel Lock			
☐	208WL8	Push to Lock	NC
☐	208WL41	Attendant Wheel Lock (Foot Lock) <i>NA w/ ext handles. NA w/19"-22" Depth and Flip Back HA armrest</i>	\$89.00
☐	208WL79	Drum Brake	\$698.00
☐	208WL10	Omit Wheel Lock	NC

Extension Handle			
☐	208EH7	6" Pair	\$88.00

Hub Lock			
<i>Only available with smaller wheels</i>			
☐	208LK01	12" Hub Lock	\$116.00
☐	208LK02	16" Hub Lock	\$116.00

5. HANGER AND FOOTREST

Hangers			
<i>Footrests mount to extension tubes</i>			
☐	208H145	60° Swing In/Swing Out	NC
☐	208H139	70° Swing In/Swing Out	NC
☐	208H141	80° Swing In/Swing Out <i>NA w/ 8" casters</i>	NC
☐	208H160	Elevating Hanger - Ext Mount.	\$339.00
☐	208H159	Articulating Hanger	\$526.00
☐	208H32	Omit Hanger	NC

Extension Tubes			
☐	208E1	Short	NC
☐	208E2	Medium <i>NA with 12.5" seat height</i>	NC
☐	208E5	Long <i>NA with 12.5" seat height</i>	NC
☐	208E16	High Mount <i>N/A with ELR/ALR hangers. Only available with adult footplates</i>	\$97.00


Footrest Footplates			
<i>Heel loop or leg strap not included, See Footplate Options Section</i>			
☐	208F1	Composite	STD
☐	208F271	Composite Angle Adj-(16"-22"w)	\$176.00
☐	208F5	Aluminum <i>Only available w/ELR & ALR</i>	NC
☐	208F3	Adult Angle Adjustable - Ext Mount	\$176.00
☐	208F7	Omit Footplate and Ext. Tube	NC

Footrest Options			
☐	208HL1	Heel Loop	\$47.00
☐	208HL5	Leg Strap	\$47.00
☐	208IG2	Impact Guard Neoprene <i>Only available w/ 60° or 70° ext. mount hanger</i>	\$105.00

JAY Positioning			
☐	208K3	Left JAY Lower Extremity S/A support	\$460.00
	☐	208K23 Pad 6" x 8" Left	NC
	☐	208K25 Pad 8" x 10" Left	NC
☐	208K4	Right JAY Lower Extremity S/A support	\$460.00
	☐	208K22 Pad 6" x 8" Right	NC
	☐	208K24 Pad 8" x 10" Right	NC
☐	208K2	Omit Hanger ☐ Left ☐ Right	NC


Articulating / Elevating Hanger Cover			
☐	208DI47	Neoprene Cover	\$397.00
☐	208DI48	Gel Pad Cover	\$467.00

✓	Option #	Definition	Price	✓	Option #	Definition	Price
6. BACKREST / SEATING							
Backpost Style - (Choose ONE Option - A, B, C)							
A. Stroller Handle Backpost (Fixed)				OR			
☐	208BP9	Low (18")	NC	C. Reclining Backrest (0-50°)			
☐	208BP10	Medium (21")	NC	NA w/ 22" Deep Frame and Flip Back HA Armrests.			
☐	208BP11	Tall (24")	NC	☐	208BT21	Reclining Backrest	\$1,032.00
A				☐	208B5	Medium	NC
				☐	208B9	Tall	NC
OR				☐	208BP64	Angle Adj stroller handles Height adj range 19"- 23"	\$175.00
B. Height Adjustable Back w/ Adj. Stroller Handles				☐	208PH8	Standard push handle	NC
Select desired Push Handle				☐	208PH9	Long push handle	NC
☐	208PH8	Standard Push Handle	NC	Backrest Options and Accessories			
☐	208PH9	Long Push Handle	NC	☐	208BP13	Removable Extension NA w/ 21"-22" frame width, height adj back with adj stroller handles	\$100.00
AND select desired Back Cane				☐	208BT101	Dynamic Backrest Adds 1.25" to seat depth. NA w/ Reclining Backrest or 19"-22" seat depths and Flip Back HA Armrests.	\$1,744.00
☐	208B200	Standard 15"-21"	NC	JAY J3 Back			
☐	208B201	Tall 18"-24"	NC	☐	208JB01	JAY J3 Back NA w/ O2 Holder	See OF
B				Backrest Upholstery			
				ONLY Available with Low Fixed, Medium Fixed, Standard Height Adjustable (15"-21" back height) stroller handles and in 14"-20" widths. N/A with Dynamic, Reclining Backrest, O2 holder, IV Pole, JAY Anterior Thoracic Support and JAY Padded Trunk Strap Accessories. Towel bar included.			
				EXO Tension Adj upholstery binding is <u>Standard black</u> .			
				☐	208BU70	Backrest Tension Adjust Ballistic *Reference 186BU70 for billing	\$367.00
				☐	208BU89	Backrest Tension Adj EXO *Reference 186BU70 for billing	\$387.00
				Colors: ☐ Silver ☐ Orange ☐ Red ☐ Blue			\$26.00



7. ARMREST

Armrest				Armpad			
☐	208AA65	Flip Back Height Adj (LITE 8.5"-13.5") NA w/Foot release tilt, Reclining Back, Dynamic Back and Attendant Foot WL.	\$381.00	Ergo armpad only available with Single post Height adj armrest			
☐	208AA11	Single Post Height Adjustable - Std	\$339.00	☐	208AM23	Desk Length Classic	NC
☐	208AA29	Single Post Height Adjustable - Low	\$339.00	☐	208AM21	Desk Length Contour	NC
☐	208AA7	Length Adjustable Locking Flip-Up NA w/ Contour Armpad	\$339.00	☐	208AM24	Full Length Classic	NC
☐	208AA1	Adjustable Locking Flip-Up	\$339.00	☐	208AM22	Full Length Contour	NC
☐	208AA14	Omit Armrest	NC	☐	208AM9	Foam Short Option for Length Adj Locking Flip Up	NC
				☐	208AM10	Foam Long Option for Length Adj Locking Flip Up	NC
				☐	208AM27	Ergo armpad	NC

Armrest Matrix for Single Post						
Armrest Size	Transfer Tube Length D	Receiver	Classic Pad			
			H Range	Pad length L		
				Desk	Full	
Std	8"	Std	10-3/4" to 15-3/4"	10-1/4"	14"	
		Low	8-3/4" to 13-3/4"			
Low	5-1/4"	Std	9" to 13"	10-1/4"	14"	
		Low	7-1/4" to 11-1/4"			
			Contour Pad			
			Receiver	H Range	Pad length L	
					Desk	Full
			Std	10" to 15"	10-1/2"	14-1/4"
			Low	8-1/4" to 13-1/4"		
			Std	8-1/2" to 12-1/2"	10-1/2"	14-1/4"
			Low	6-3/4" to 10-3/4"		
			Ergo Pad			
			Receiver	H Range	Pad length L	
					11-3/4"	
			Std	10-1/2" to 15-1/2"	11-3/4"	
			Low	8-3/4" to 13-3/4"		
Std	9" to 13"					
Low	7-1/4" to 11-1/4"					


✓	Option #	Definition	Price	✓	Option #	Definition	Price
8. KEY OPTIONS AND ACCESSORIES							
Anti Tips				JAY Anterior Thoracic Support			
Not ordered as separate line item.				Style			
☐	208AT1	Anti Tips (Pair)	\$121.00	☐	208JY01	Classic Dynamic	\$177.00
Positioning Belt				☐	208JY02	Classic Structure	\$177.00
☐	208ZZ25	JAY 2PT Belt 2" Push Btn No Pad	\$64.00	☐	208JY03	Contour Dynamic	\$177.00
☐	208ZZ26	JAY 2PT Belt 2" Push Btn 9.5" Pad	\$126.00	☐	208JY04	Contour Structure	\$177.00
☐	208ZZ14	JAY 2PT Belt 1.5" Push Btn No Pad	\$51.00	☐	208JY05	Center Opening with Buckle	\$177.00
☐	208ZZ15	JAY 2PT Belt 1.5" Push Btn 7.5" Pad	\$80.00	Size			
☐	208A117	2" Aircraft Buckle w/Clamps	\$139.00	☐	208JZ01	Small	NC
☐	208A118	2" Aircraft Buckle Pad w/Clamp	\$204.00	☐	208JZ02	Medium	NC
☐	208A120	Positioning Belt Clamps Only	\$61.00	☐	208JZ03	Large	NC
Wheelchair Tray Table				☐	208JZ04	X-Large	NC
Only available if armrests are included. NA with JAY Arm Troughs.				JAY Padded Trunk Strap			
☐	208A30	Small (13"-14")	\$381.00	☐	208JA01	Padded Trunk Strap - Small	\$165.00
☐	208A31	Medium (15"-17")	\$381.00	☐	208JA02	Padded Trunk Strap - Medium	\$165.00
☐	208A32	Large (18"-20")	\$381.00	☐	208JA03	Padded Trunk Strap - Large	\$165.00
☐	208A33	X-Large (21"-22")	\$381.00	☐	208JA04	Padded Trunk Strap - XL	\$165.00
				Miscellaneous Accessories			
				☐	208A1	Black Utility Backpack	\$88.00
				☐	208A8	Black Seat Pouch	\$58.00
				☐	208LS14	Extra Leg Strap Black	\$36.00
				☐	208A601	O2 Holder NA w/ Height Adj Stroller Handles. NA w/ Reclining Backrest	\$187.00
				☐	208A701	IV Pole	\$321.00
				☐	208A999	Quickie Water Bottle & Holder	\$63.00

9. COLOR							
Standard Frame Color				Anodized Color Package			
Color options only to available to hangers, seat tubes on seat. Front and back extensions on the base				Anodize package include Forks, axle plate, back plate.			
☐	208C211	Celestial Sapphire	NC	☐	208EB07	Black	NC
☐	208C212	Desert Sage	NC	☐	208EB08	Blue	\$187.00
☐	208C213	Desert Tan	NC	☐	208EB13	Gold	\$187.00
☐	208C214	Gloss Black	NC	☐	208EB09	Lime Green	\$187.00
☐	208C215	Matte Black	NC	☐	208EB10	Purple	\$187.00
☐	208C216	Metallic Shadow	NC	☐	208EB11	Red	\$187.00
☐	208C217	Midnight Blue	NC	☐	208EB12	Ti Grey	\$187.00
☐	208C218	Mystic Purple	NC	Additional Options Through Built for Me			
☐	208C219	Petal Pink	NC	To inquire about any unique personalization for Quickie products, contact:			
☐	208C220	Pink Power	NC	Phone: (800) 263-3390, Email: cscanada@sunmed.com			
☐	208C221	Red Pepper	NC				
☐	208C222	Royal Blue	NC				
☐	208C223	Salute Green	NC				
☐	208C224	Scarlet Noir	NC				
☐	208C225	Silver Ore	NC				
☐	208C226	Solar Orange	NC				
☐	208C227	Starlight White	NC				
☐	208C228	Storm Green	NC				
☐	208C229	Turbo Green	NC				
☐	208C230	Velocity Yellow	NC				
Personalized Frame Color							
☐	208C116	Desert Camouflage	NC				
☐	208C151	Mossy Oak Camouflage	NC				
☐	208C150	Stars & Stripes	NC				
☐	208C149	Zebra (choose above base color) Std w/ Starlight White. N/A w/ Matte Black and Gloss Black	NC				

Lower Leg Length Matrix

Front seat height must be at least 2" taller than lower leg length on all swing in/out hangers. Front seat height must be at least 2" taller than maximum lower leg length on all front mount hangers.

Hanger	Extension	Seat Pan to Top of Footrest	
		Min	Max
60° Swing In/Out	Short	8 1/2	14
	Med	8 1/2	15 1/2
	Long	9 1/2	17 1/2
70° Swing In/Out	Short	9 1/2	14 1/2
	Med	10	16 1/2
	Long	11 1/2	18 1/2
80° Swing In/Out	Short	10	14
	Med	13	16
	Long	13	18
ALR	Short	9	13
	Med	9	14 1/2
	Long	10	16 1/2
ELR	Short	9	13
	Medium	9	14 1/2
	Long	10	16 1/2





JAY® Signature Series Backrest and Accessories

Note: Not all sizes and configurations shown below, for full offering visit www.sunrisemedical.ca

When ordering mark the part number desired and it will be added to the order as a separate line item .

JAY BACKREST


JAY J3 Posterior (PA) BACK - 2.2" Contour Depth												
	Shell Height		Shell Width									
Mid Thoracic (MT)	12.0" (S)	J3PAMTSxxFS	☐ 12	\$888	☐ 14	\$888	☐ 16	\$888	☐ 18	\$888	☐ 20	\$888
		J3PAMTSxxTH	☐ 12	\$888	☐ 14	\$888	☐ 16	\$888	☐ 18	\$888	☐ 20	\$888
Upper Thoracic (UT)	16.5" (S)	J3PAUTSxxFS	☐ 12	\$991	☐ 14	\$991	☐ 16	\$991	☐ 18	\$991	☐ 20	\$991
		J3PAUTSxxTH	☐ 12	\$991	☐ 14	\$991	☐ 16	\$991	☐ 18	\$991	☐ 20	\$991
	18.0" (M)	J3PAUTMxxFS			☐ 14	\$991	☐ 16	\$991	☐ 18	\$991	☐ 20	\$991
		J3PAUTMxxTH			☐ 14	\$991	☐ 16	\$991	☐ 18	\$991	☐ 20	\$991
	19.5" (T)	J3PAUTTxxFS			☐ 14	\$991	☐ 16	\$991	☐ 18	\$991	☐ 20	\$991
		J3PAUTTxxTH			☐ 14	\$991	☐ 16	\$991	☐ 18	\$991	☐ 20	\$991
Shoulder Height (SH)	21.0" (S)	J3PASHSxxFS			☐ 14	\$1,021	☐ 16	\$1,021	☐ 18	\$1,021	☐ 20	\$1,021
		J3PASHSxxTH			☐ 14	\$1,021	☐ 16	\$1,021	☐ 18	\$1,021	☐ 20	\$1,021
	22.5" (M)	J3PASHMxxFS					☐ 16	\$1,021	☐ 18	\$1,021	☐ 20	\$1,021
		J3PASHMxxTH					☐ 16	\$1,021	☐ 18	\$1,021	☐ 20	\$1,021
	24.0" (T)	J3PASHTxxFS							☐ 18	\$1,021	☐ 20	\$1,021
		J3PASHTxxTH							☐ 18	\$1,021	☐ 20	\$1,021

FS Mounting Hardware



JAY J3 Posterior Lateral (PL) BACK - 3" Contour Depth												
Mid Thoracic (MT)	15.0" (T)	J3PLMTTxxFS	☐ 12	\$979	☐ 14	\$979	☐ 16	\$979	☐ 18	\$979	☐ 20	\$979
		J3PLMTTxxTH	☐ 12	\$979	☐ 14	\$979	☐ 16	\$979	☐ 18	\$979	☐ 20	\$979
Upper Thoracic (UT)	16.5" (S)	J3PLUTSxxFS	☐ 12	\$1,063	☐ 14	\$1,063	☐ 16	\$1,063	☐ 18	\$1,063	☐ 20	\$1,063
		J3PLUTSxxTH	☐ 12	\$1,063	☐ 14	\$1,063	☐ 16	\$1,063	☐ 18	\$1,063	☐ 20	\$1,063
	18.0" (M)	J3PLUTMxxFS		☐ 14	\$1,063	☐ 16	\$1,063	☐ 18	\$1,063	☐ 20	\$1,063	
		J3PLUTMxxTH		☐ 14	\$1,063	☐ 16	\$1,063	☐ 18	\$1,063	☐ 20	\$1,063	
	19.5" (T)	J3PLUTTxxFS		☐ 14	\$1,063	☐ 16	\$1,063	☐ 18	\$1,063	☐ 20	\$1,063	
		J3PLUTTxxTH		☐ 14	\$1,063	☐ 16	\$1,063	☐ 18	\$1,063	☐ 20	\$1,063	
Shoulder Height (SH)	21.0" (S)	J3PLSHSxxFS		☐ 14	\$1,092	☐ 16	\$1,092	☐ 18	\$1,092	☐ 20	\$1,092	
		J3PLSHSxxTH		☐ 14	\$1,092	☐ 16	\$1,092	☐ 18	\$1,092	☐ 20	\$1,092	
	22.5" (M)	J3PLSHMxxFS			☐ 16	\$1,092	☐ 18	\$1,092	☐ 20	\$1,092		
		J3PLSHMxxTH			☐ 16	\$1,092	☐ 18	\$1,092	☐ 20	\$1,092		
	24.0" (T)	J3PLSHTxxFS				☐ 18	\$1,092	☐ 20	\$1,092			
		J3PLSHTxxTH				☐ 18	\$1,092	☐ 20	\$1,092			

TH Mounting Hardware



JAY J3 Posterior Deep (PD) BACK - 6" Contour Depth												
Mid Thoracic (MT)	15.0" (T)	J3PDMTTxxFS	☐ 12	\$1,073	☐ 14	\$1,073	☐ 16	\$1,073	☐ 18	\$1,073	☐ 20	\$1,073
		J3PDMTTxxTH	☐ 12	\$1,073	☐ 14	\$1,073	☐ 16	\$1,073	☐ 18	\$1,073	☐ 20	\$1,073
Upper Thoracic (UT)	18.0" (M)	J3PDUTMxxFS			☐ 14	\$1,161	☐ 16	\$1,161	☐ 18	\$1,161	☐ 20	\$1,161
		J3PDUTMxxTH			☐ 14	\$1,161	☐ 16	\$1,161	☐ 18	\$1,161	☐ 20	\$1,161
Shoulder Height (SH)	21.0" (S)	J3PDSHSxxFS			☐ 14	\$1,237	☐ 16	\$1,237	☐ 18	\$1,237	☐ 20	\$1,237
		J3PDSHSxxTH			☐ 14	\$1,237	☐ 16	\$1,237	☐ 18	\$1,237	☐ 20	\$1,237

JAY J3 Back Whitmyer Headrest

AXYS Mount and Pad

☐	J3HEDX8	AXYS Headrest with 8" Plush Pad	\$403.00
☐	J3HEDX10	AXYS Headrest with 10" Plush Pad	\$403.00
☐	J3HEDX14	AXYS Headrest with 14" Plush Pad	\$403.00
☐	J3HEDXAAPN	AXYS with Adj-a-Plush Narrow Pad	\$551.00
☐	J3HEDXAAPS	AXYS with Adj-a-Plush Std Pad	\$551.00

LINX Mount and Pad

☐	J3HEDL8	LINX Headrest with 8" Plush Pad	\$528.00
☐	J3HEDL10	LINX Headrest with 10" Plush Pad	\$528.00
☐	J3HEDL14	LINX Headrest with 14" Plush Pad	\$528.00
☐	J3HEDLAAPN	LINX with Adj-a-Plush Narrow Pad	\$508.00
☐	J3HEDLAAPS	LINX with Adj-a-Plush Std Pad	\$508.00

JAY J3 Back Laterals Swing Away Laterals With Pad - Pair

☐	JLATSA44PR	SwingAway Lat with 4"x4" Pad (PR)	\$668.00	☐	JLATSA64PR	SwingAway Lat with 6"x4" Pad (PR)	\$668.00
☐	JLATSA45PR	SwingAway Lat with 4"x5" Pad (PR)	\$668.00	☐	JLATSA65PR	SwingAway Lat with 6"x5" Pad (PR)	\$668.00
☐	JLATSA46PR	SwingAway Lat with 4"x6" Pad (PR)	\$668.00	☐	JLATSA66PR	SwingAway Lat with 6"x6" Pad (PR)	\$668.00



JAY® Signature Series Cushions

Note: Not all sizes and configurations shown below, for full offering visit www.sunrisemedical.ca
When ordering mark the part number desired and it will be added to the order as a separate line item .

JAY CUSHIONS

JAY Fusion Cushion

☐	JFUSION1416	JAY Fusion 14" x 16"	\$1,045.00
☐	JFUSION1616	JAY Fusion 16" x 16"	\$1,045.00
☐	JFUSION1618	JAY Fusion 16" x 18"	\$1,045.00
☐	JFUSION1620	JAY Fusion 16" x 20"	\$1,045.00
☐	JFUSION1717	JAY Fusion 17" x 17"	\$1,045.00
☐	JFUSION1816	JAY Fusion 18" x 16"	\$1,045.00
☐	JFUSION1818	JAY Fusion 18" x 18"	\$1,045.00
☐	JFUSION1820	JAY Fusion 18" x 20"	\$1,045.00
☐	JFUSION2018	JAY Fusion 20" x 18"	\$1,045.00
☐	JFUSION2020	JAY Fusion 20" x 20"	\$1,045.00
☐	JFUSION2220	JAY Fusion 22" x 20"	\$1,172.00

JAY Fusion Cushion w/ Single Valve Air Insert

☐	JFUSION1416AS	JAY Fusion 14" x 16"	\$1,045.00
☐	JFUSION1616AS	JAY Fusion 16" x 16"	\$1,045.00
☐	JFUSION1618AS	JAY Fusion 16" x 18"	\$1,045.00
☐	JFUSION1620AS	JAY Fusion 16" x 20"	\$1,045.00
☐	JFUSION1717AS	JAY Fusion 17" x 17"	\$1,045.00
☐	JFUSION1816AS	JAY Fusion 18" x 16"	\$1,045.00
☐	JFUSION1818AS	JAY Fusion 18" x 18"	\$1,045.00
☐	JFUSION1820AS	JAY Fusion 18" x 20"	\$1,045.00
☐	JFUSION2018AS	JAY Fusion 20" x 18"	\$1,045.00
☐	JFUSION2020AS	JAY Fusion 20" x 20"	\$1,045.00
☐	JFUSION2220AS	JAY Fusion 22" x 20"	\$1,172.00

JAY Fusion Cushion w/ Dual Valve Air Insert

☐	JFUSION1416AD	JAY Fusion 14" x 16"	\$1,063.00
☐	JFUSION1616AD	JAY Fusion 16" x 16"	\$1,063.00
☐	JFUSION1618AD	JAY Fusion 16" x 18"	\$1,063.00
☐	JFUSION1620AD	JAY Fusion 16" x 20"	\$1,063.00
☐	JFUSION1717AD	JAY Fusion 17" x 17"	\$1,063.00
☐	JFUSION1816AD	JAY Fusion 18" x 16"	\$1,063.00
☐	JFUSION1818AD	JAY Fusion 18" x 18"	\$1,063.00
☐	JFUSION1820AD	JAY Fusion 18" x 20"	\$1,063.00
☐	JFUSION2018AD	JAY Fusion 20" x 18"	\$1,063.00
☐	JFUSION2020AD	JAY Fusion 20" x 20"	\$1,063.00
☐	JFUSION2220AD	JAY Fusion 22" x 20"	\$1,191.00

JAY Union Cushion

☐	JU1416	JAY UNION 14" x 16"	\$764.00
☐	JU1616	JAY UNION 16" x 16"	\$764.00
☐	JU1618	JAY UNION 16" x 18"	\$764.00
☐	JU1620	JAY UNION 16" x 20"	\$764.00
☐	JU1717	JAY UNION 17" x 17"	\$764.00
☐	JU1816	JAY UNION 18" x 16"	\$764.00
☐	JU1818	JAY UNION 18" x 18"	\$764.00
☐	JU1820	JAY UNION 18" x 20"	\$806.00
☐	JU2018	JAY UNION 20" x 18"	\$806.00
☐	JU2020	JAY UNION 20" x 20"	\$806.00
☐	JU2220	JAY UNION 22" x 20"	\$806.00

JAY Union Reduced Profile Cushion

☐	JU1416RP	JAY Union RP 14" x 16"	\$764.00
☐	JU1616RP	JAY Union RP 16" x 16"	\$764.00
☐	JU1618RP	JAY Union RP 16" x 18"	\$764.00
☐	JU1620RP	JAY Union RP 16" x 20"	\$764.00
☐	JU1717RP	JAY Union RP 17" x 17"	\$764.00

JAY J2 Cushion

☐	2108-MJ	JAY J2 14" x 16"	\$832.00
☐	2105-MJ	JAY J2 15.5" x 16"	\$832.00
☐	2106-MJ	JAY J2 15.5" x 18"	\$832.00
☐	2107-MJ	JAY J2 15.5" x 20"	\$958.00
☐	2117-MJ	JAY J2 17" x 17"	\$832.00
☐	2100-MJ	JAY J2 18" x 16"	\$832.00
☐	2101	JAY J2 18" x 18"	\$832.00
☐	2103-MJ	JAY J2 18" x 20"	\$958.00
☐	2118-MJ	JAY J2 20" x 18"	\$941.00
☐	2102-MJ	JAY J2 20" x 20"	\$1,112.00

JAY J2 Deep Contour Cushion

☐	2446	JAY J2 Deep Contour 14" x 16"	\$870.00
☐	2466	JAY J2 Deep Contour 16" x 16"	\$870.00
☐	2468	JAY J2 Deep Contour 16" x 18"	\$870.00
☐	2460	JAY J2 Deep Contour 16" x 20"	\$945.00
☐	2477	JAY J2 Deep Contour 17" x 17"	\$870.00
☐	2486	JAY J2 Deep Contour 18" x 16"	\$870.00
☐	2488	JAY J2 Deep Contour 18" x 18"	\$870.00
☐	2480	JAY J2 Deep Contour 18" x 20"	\$945.00
☐	2408	JAY J2 Deep Contour 20" x 18"	\$945.00
☐	2400-MJ	JAY J2 Deep Contour 20" x 20"	\$1,093.00
☐	2428	JAY J2 Deep Contour 22" x 18"	\$1,093.00

JAY X2 Cushion

☐	JX21416	JAY X2 14" x 16"	\$944.00
☐	JX21616	JAY X2 16" x 16"	\$944.00
☐	JX21618	JAY X2 16" x 18"	\$944.00
☐	JX21620	JAY X2 16" x 20"	\$944.00
☐	JX21717	JAY X2 17" x 17"	\$944.00
☐	JX21816	JAY X2 18" x 16"	\$944.00
☐	JX21818	JAY X2 18" x 18"	\$944.00
☐	JX21820	JAY X2 18" x 20"	\$944.00
☐	JX22020	JAY X2 20" x 20"	\$944.00

JAY Ion Cushion

☐	JION1416	JAY ION 14" x 16"	\$657.00
☐	JION1616	JAY ION 16" x 16"	\$657.00
☐	JION1618	JAY ION 16" x 18"	\$657.00
☐	JION1620	JAY ION 16" x 20"	\$657.00
☐	JION1717	JAY ION 17" x 17"	\$657.00
☐	JION1816	JAY ION 18" x 16"	\$657.00
☐	JION1818	JAY ION 18" x 18"	\$657.00
☐	JION1820	JAY ION 18" x 20"	\$657.00
☐	JION2018	JAY ION 20" x 18"	\$698.00
☐	JION2020	JAY ION 20" x 20"	\$698.00
☐	JION2220	JAY ION 22" x 20"	\$698.00

☐	JU1816RP	JAY Union RP 18" x 16"	\$764.00
☐	JU1818RP	JAY Union RP 18" x 18"	\$764.00
☐	JU1820RP	JAY Union RP 18" x 20"	\$764.00
☐	JU2018RP	JAY Union RP 20" x 18"	\$764.00
☐	JU2020RP	JAY Union RP 20" x 20"	\$764.00



Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2, Concord, ON L4K 4P8