



Breezy

BREEZY 600

The Breezy 600 is the lightest weight feature-packed chair in its class, making it the ideal solution for foot propellers and those needing environmental access. Ten different seat height adjustments allow this chair to fit every client.



Features:

- Seat to Floor Height Range
Standard 15"-20"
Ultra Hemi 13"-18"
- Seat Width 14"-20"
- Lightweight aluminum sideframes
- Angle adjustable back
- Adjustable caster housing for better stability and maximum rolling efficiency
- Padded seat sling for extra comfort
- Locking crossbrace hinge limits flex in the frame for more efficient propelling
- Weight Capacity: 250 lbs
- Lifetime warranty on the frame & 2 years on parts (See owner's manual for details)



One-Arm Drive positive lock strut system increases propelling efficiency and is easily removed for folding. Option Fits Breezy 600 models.



Standard, hemi and ultra-hemi configurations



www.SunriseMedical.ca

© 2018 Sunrise Medical Canada Inc.

Concord Head Office
Phone (Toll-Free): +1.800.263.3390
Fax (Toll-Free): +1.800.561.5834
E-Mail: cscanada@sunmed.com

Maple Ridge Office
Phone (Toll-Free): +1.800.888.285.2488
Fax (Toll-Free): +1.888.830.6684
E-Mail: cscanada@sunmed.com

Lasalle Office
Phone (Toll-Free): +1.888.548.8688
Fax (Toll-Free): +1.888.909.5090
E-Mail: cquebec@sunmed.com