



ORDER FORM

Advance Bath Chair

EILK04

Versatile bath chair for all ages
ensuring safety and comfort for your child

Customer

Account			
Account #		Date	
☐ Quote	☐ Order	PO #	
Buyer			
Mark For			

Shipping Address

Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

Advance Bath Chair

Weight Capacity: 159 lbs

✓ OPTION	DESCRIPTION	RETAIL
Packages		\$
Includes Chair Frame, Chair Cover, Hip Belt, and Chest Belt. Substitutions are not available for packages.		
☐ 152QE23	Size 1 Advance Bath Chair Package - Pink	\$1,080
☐ 152QE24	Size 1 Advance Bath Chair Package - Royal Blue	\$1,080
☐ 152QE25	Size 2 Advance Bath Chair Package - Pink	\$1,080
☐ 152QE26	Size 2 Advance Bath Chair Package - Royal Blue	\$1,080
☐ 152QE27	Size 3 Advance Bath Chair Package - Pink	\$1,156
☐ 152QE28	Size 3 Advance Bath Chair Package - Royal Blue	\$1,156
☐ 152QE29	Size 4 Advance Bath Chair Package - Pink	\$1,156
☐ 152QE30	Size 4 Advance Bath Chair Package - Royal Blue	\$1,156
Chair Frame		\$
Selection Required if an Advance Bath Chair Package was not selected.		
☐ 152BS29	Size 1 Chair Frame	\$628
☐ 152BS30	Size 2 Chair Frame	\$628
☐ 152BS31	Size 3 Chair Frame	\$688
☐ 152BS32	Size 4 Chair Frame	\$719
Chair Cover		\$
Selection Required. Must match chair frame size.		
☐ 152IU23	Size 1 Cover - Pink	\$211
☐ 152IU24	Size 1 Cover - Royal Blue	\$211
☐ 152IU25	Size 2 Cover - Pink	\$211
☐ 152IU26	Size 2 Cover - Royal Blue	\$211
☐ 152IU27	Size 3 Cover - Pink	\$227
☐ 152IU28	Size 3 Cover - Royal Blue	\$227
☐ 152IU29	Size 4 Cover - Pink	\$227
☐ 152IU30	Size 4 Cover - Royal Blue	\$227
Accessories		\$
Abduction Belts must match chair frame size.		
☐ 152A876	Hip Belt - Pink	\$118
☐ 152A877	Hip Belt - Royal Blue	\$118
☐ 152EJ07	Chest Belt - Pink	\$123
☐ 152EJ08	Chest Belt - Royal Blue	\$123
☐ 152A878	Size 1 and 2 Abduction Belt - Pink	\$123
☐ 152A879	Size 1 and 2 Abduction Belt - Royal Blue	\$123
☐ 152A880	Size 3 and 4 Abduction Belt - Pink	\$123
☐ 152A881	Size 3 and 4 Abduction Belt - Royal Blue	\$123
☐ 152HQ08	Head Laterals - Pink	\$123
☐ 152HQ09	Head Laterals - Royal Blue	\$123
☐ 1528A119	Shower Trolley	\$1,327

Technical Specifications - Imperial

Specification	Size 1	Size 2	Size 3	Size 4
Age	1 - 5 years	4 - 9 years	8 - 14 years	12 - 18 years
User Height (Max / Min)	30" / 41"	37" / 49"	45" / 57"	53" / 65"
Head Support	6"	6"	8.5"	8.5"
Back Support Height	14"	16"	18"	22"
Seat Depth	9"	11"	13"	16"
Calf Support	10"	10"	13"	13"
Seat Height	2" - 17"	2" - 17"	2" - 17"	2" - 17"
Maximum Weight	159 lbs	159 lbs	159 lbs	159 lbs
Total Width	17"	17"	17"	17"

Technical Specifications - Metric

Specification	Size 1	Size 2	Size 3	Size 4
Age	1 - 5 years	4 - 9 years	8 - 14 years	12 - 18 years
User Height (Max / Min)	750 mm / 1050 mm	950 mm / 1250 mm	1150 mm / 1450 mm	1350 mm / 1650 mm
Head Support	150 mm	150 mm	215 mm	215 mm
Back Support Height	350 mm	400 mm	460 mm	550 mm
Seat Depth	225 mm	275 mm	335 mm	400 mm
Calf Support	260 mm	260 mm	340 mm	340 mm
Seat Height	50 mm - 425 mm	50 mm - 425 mm	50 mm - 425 mm	50 mm - 425 mm
Maximum Weight	72 kgs	72 kgs	72 kgs	72 kgs
Total Width	430 mm	430 mm	430 mm	430 mm

Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2
Concord, ON L4K 4P8, Canada



MK-129844-CAN Rev. 6 © 04/2025

