



October 2026

Q700M with SEDEO[®] PRO[™], PRO[™] ADVANCED and JAY[®]
Syncra[™] Seating
AADL - W581

EIPW25

Account:	
Date:	
PO#:	
Buyer Name:	
Marked For:	


Ship To:	
Name:	
Address:	
City, Prov., PC:	
Phone:	

Sunrise Medical Inc.
1000 Creditstone Road, Unit 2
Concord, ON L4K 4P8
1-800-263-3390



Client Measurements

a)	Top of head to seat	Inch
b)	Top of shoulder to seat	Inch
c)	Underarm to seat	Inch
d)	Elbow to seat	Inch
e)	Top of leg to seat	Inch
f)	Shoulder width	Inch
g)	Chest width	Inch
i)	Hip width	Inch
j)	Actual sitting depth	Inch
m)	Foot length	Inch
(n)	Knee to foot	Inch



User Weight

Notes:

Chair Specifications

		Seat Type	SEDEO PRO / JAY Syncra
Typical Weight w/out Batteries	330 lbs.	Seat Size Range	Width 16" to 22" Depth 14" to 22"
Typical Weight with Batteries	410 lbs.		
Power Seating	50° of CG Tilt and 150° Power Recline with 4" mechanical sheer reduction / 12" seat elevator	STF Height Range	16.5" to 19"

Mid-Wheel Drive Base Specifications

Drive Wheel Position	Mid-wheel Drive	Drive Wheels	14" drive wheels w/flat free inserts (black)
User Weight Capacity	300 lbs.	Tires/Casters	6" Front / 6" Rear (Black)
Minimum Turning Radius	19.8" (excludes footplates)		
Base Size	Width 24.5", Length 36.5"	Charger	10A Dual Mode (Compact)
Drive Control	R-net 120A / R-net 90A	Battery Size	Group 24 / Group 34
Driving Range	Up to 18 miles (w/ Group 24 Batteries)	Color Accents	9 Shroud Colors and 7 Color Accents
Drive Motors	4-pole 6mph, 8mph	Ground Clearance	3"

BASE			
OPTION		Description	Retail
Standard Duty Mid Wheel Drive Base = End Item Number EIPW25			
☐	250PC197	Quickie Q700 M NPC No Power Option for SEDEO PRO Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5". (AADL Contract Price)	\$8,993.00
☐	250PC198	Quickie Q700 M SPC Single Power Option for SEDEO PRO Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5". (AADL Contract Price)	\$8,993.00
☐	250PC203	Quickie Q700 M MPC Tilt/Recline For SEDEO PRO Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5". (AADL Contract Price)	\$8,993.00
Standard Duty Mid Wheel Drive Base = End Item Number EIPW25			
☐	250PC201	Quickie Q700 M NPD No Power Option for SEDEO PRO Seat Included in Base Pricing : Q700 M Base, R-net 120amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5". (Client Cost)	\$10,781.00
☐	250PC202	Quickie Q700 M SPD Single Power Option for SEDEO PRO Seat Included in Base Pricing : Q700 M Base, R-net 120amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5". (Client Cost)	\$10,781.00
☐	250PC205	Quickie Q700 M MPD Tilt/Recline For SEDEO PRO Seat Included in Base Pricing : Q700 M Base, R-net 120amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5". (Client Cost)	\$10,781.00

BASE			
OPTION		Description	Retail
User Weight			
☐	250WG23	Up to 300 lbs.	NC
☐	250WG15	301 - 350 lbs.	NC
☐	250WG25	351 - 450 lbs.	NC
Motor Package			
☐	250PF99	4-pole 10 kph / 6 mph	NC
☐	250PF100	4-pole 12.5 kph / 8 mph	NC
☐	250PF101	4-pole 12.5 kph / 8 mph (Client Cost)	\$1,268.00
Tracking			
☐	250MI01	SureTrac Module & Motor Encoders <i>(Improves Drive Tracking)</i>	\$1,696.00
Drive Wheels			
☐	250DW58	Drive Wheel, 14", Dual Density, Knobby, Black	NC
Caster Wheels			
☐	250CW55	Caster Wheels 6", Solid, Black	NC
Base Shroud Color			
NEW: Base Shroud and Backrest Shroud Painted to Match			
☐	250C119	Matte Black	NC
☐	250C165	Candy Apple Red	NC
☐	250C166	Candy Blue	NC
☐	250C167	Candy Purple	NC
☐	250C176	Green Apple	NC
☐	250C178	Pearl Pink	NC
☐	250C182	Silver	NC
☐	250C125	Pearl White	NC
Color Accents			
☐	250EZ20	Apple Green	NC
☐	250EZ22	Black	NC
☐	250EZ21	Blue	NC
☐	250EZ23	Orange	NC
☐	250EZ24	Red	NC
☐	250EZ25	Silver	NC
☐	250EZ19	White	NC

BASE									
OPTION		Description						Retail	
Service Decals									
☐	250D10		Army						NC
☐	250D11		Air Force						NC
☐	250D13		Marines						NC
☐	250D14		Navy						NC
☐	250D15		Purple Heart						NC
ISO Transit Tie Downs									
☐	250TO50		ISO Transit Package						NC
LED Lights <i>(Requires R-net Electronics and operated through joystick)</i>									
☐	250LI08		LED Fender Lights						\$646.00
☐	250LI09		LED Fender Lights <i>(Available w/ Multi-Power Only)</i>						NC
☐	250LI02		Omit LED Lights						NC
Seat Heights									
<i>Note: Available seat-to-Floor heights are dependant on 2 options: 1. Battery size and 2. Powered Seating.</i>									
<i>Use the matrix below to determine available ranges for your desired options. Measured from front of seat pan to floor with 0° frame angle.</i>									
			Seat Height	Powered Seating Without Elevate	No Powered Seating	Pro Advanced Recline			
						w/o Elevate	w/ Elevate		
☐	250HT103		16.5"	Group 34	Group 34	N/A	N/A		NC
☐	250HT55		17"	Group 34	GP24 or GP34	GP34	GP34		NC
☐	250HT104		17.5"	Group 34	GP24 or GP34	GP34	GP34		NC
☐	250HT56		18"	GP24 or GP34	GP24 or GP34	GP24 / GP34	GP34		NC
☐	250HT105		18.5"	GP24 or GP34	GP24 or GP34	GP24 / GP34	GP24 / GP34		NC
☐	250HT57		19"	GP24 or GP34	GP24 or GP34	GP24 / GP34	GP24 / GP34		NC
Batteries & Chargers									
Group 24 Battery Options <i>(Recommended for max daily range)</i>									
☐	250AG25		GP24 Flat Top Pair w/Strap MK Gel						\$955.00
☐	250AG26		GP24 Flat Top Pair w/Strap Brightway AGM						\$915.00
Group 34 Battery Options <i>(Required for lower seat-to-floor heights)</i>									
☐	250AG27		GP34 Flat Top Pair w/Straps MK Gel						\$955.00
☐	250AG29		GP34 Flat Top Pair w/Straps Brightway AGM						\$863.00
Omit Batteries									
☐	250AG24		Battery Harnesses Only <i>(Omit Batteries)</i>						NC
Chargers									
☐	250BJ28		10A Dual Mode						NC

Powered Seating			
OPTION		Description	Retail
Power Seating Options			
Seat Elevator (Client Cost)			
☐	250PI08	C-Me Seat Elevator, Powered, 12" (C-Me speed up to 3.0 mph)	\$3,477.00
☐	250PI07	Seat Elevator, Powered, 12" (Required for Anterior Tilt)	\$3,477.00
Seat Tilt (AADL W992 For Power Tilt)			
☐	250PA35	Posterior 50° Tilt	\$8,663.00
Seat Recline (AADL W992TR for Power Recline)			
☐	250PA37	150° Powered Recline w/ 4" Shear Reduction	\$8,001.00
Seat Tilt/Recline (AADL W992TR for Tilt and Recline)			
☐	250PA41	Posterior 50° Tilt and Pro 150° Recline Powered, w/ 4" Shear Reduction	\$13,275.00
☐	250PA62	Posterior 50° Tilt and Pro Advanced 165° Recline Powered, w/ 7" Shear Reduction	\$13,785.00
Seat Anterior Tilt (Requires Pro Advanced)			
☐	250TG02	Anterior tilt 0° - 30° Tilt	\$3,329.00
Power Recline Backrest Preset Home Angle			
☐	250QK1	90°	NC
☐	250QK2	95° (Default Position)	NC
Pre Set Seat Frame Angle (No Power Seating Options, Increases Seat to Floor Height)			
☐	250FA10	0° (Adds 0" to STFH)	NC
☐	250FA3	3° (Adds 1" to STFH)	NC
Backrest Frame Type (N/A with Powered Recline)			
☐	250BF10	Tool Angle Adjustable Back, 90°-120°	NC
☐	250BF11	User Adjustable Manual Semi-Recline Back, 90°-130°, Right (Adjustment Lever Mounted on Opposite Side of Joystick)	\$824.00
☐	250BF12	User Adjustable Manual Semi-Recline Back, 90°-130°, Left (Adjustment Lever Mounted on Opposite Side of Joystick)	\$824.00
Backrest Frame Option (N/A with Pro Advanced)			
☐	250BF13	Flip Down, Quick Release	\$179.00
Back Cushion Type			
☐	250BT134	SEDEO Back Cushion (Choose Seating in Section A) (N/A with Pro Advanced)	NC
☐	250BT138	JAY Syncra (Choose Seating in Section B)	NC
☐	250BT135	JAY J3 (Choose Seating in Section C)	NC
☐	250BT136	JAY J3 Interface Only (Allows JAY Your Way modifications)	NC
☐	250BT137	SEDEO Aftermarket Back Adapter (Height Adjustable Round Tubes for Traditional Mounting)	\$633.00
☐	250BT133	OMIT (For Custom Mounting. No Headrest Mount)	NC

Powered Seating

OPTION	Description	Retail
--------	-------------	--------

A- SEDEO Backrest Dimensions (N/A with Pro Advanced)

Back Width

☐	250BN28	400mm / 16"	NC
☐	250BN29	460mm / 18"	NC
☐	250BN30	500mm / 20"	NC
☐	250BN31	560mm / 22"	NC

Back Shell Height

☐	250B255		460mm / 18"	Height Back	22"-23"	23"-25"	NC
☐	250B256		560mm / 22"		24"-27"	26"-29"	NC
☐	250B257		660mm / 26"		28"-30"	30"	NC

Installed Back Height (From Pan to Top of Back)

☐	250JR32	22"	NC
☐	250JR33	23"	NC
☐	250JR34	24"	NC
☐	250JR35	25"	NC
☐	250JR36	26"	NC
☐	250JR37	27"	NC
☐	250JR38	28"	NC
☐	250JR39	29"	NC
☐	250JR40	30"	NC

Powered Seating

OPTION	Description	Retail
--------	-------------	--------

B - JAY Syncra Seating (Seating Benefits)

B - JAY Syncra Back

☐	250JP14	JAY Syncra Back 16-20" Wide with 3" Contour Depth	\$1,339.00
☐	250JP15	JAY Syncra Back 22" Wide with 3" Contour Depth	\$1,339.00

B - JAY Syncra Back Width

Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.

☐	250JQ07	16" Wide (Adjustment range 16 to 17" Wide)*	NC
☐	250JQ08	18" Wide (Adjustment range 18 to 19" Wide)*	NC
☐	250JQ09	20" Wide (Adjustment range 20 to 21" Wide)*	NC
☐	250JQ10	22" Wide (Adjustment range 22 to 23" Wide)	NC

*Additional +1" of shell width adjustment is achievable, requires purchase of 1" larger foam and upholstery separately

B - JAY Syncra Back Shell Options and Height Ranges

JAY Syncra Back shells are mounted with a minimim 5.5" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field with-in the ranges shown based on the back frame type.

JAY Syncra Shell Height

☐	250JR42	16" Tall Shell (22 to 24" Back Height)	NC
☐	250JR43	18" Tall Shell (24 to 26" Back Height)	NC
☐	250JR44	20" Tall Shell (26 to 28" Back Height)	NC
☐	250JR45	22" Tall Shell (28 to 30" Back Height)	NC

OPTION	Description	Retail
B - JAY Syncra Back Shell Options and Height Ranges Continued		
JAY Syncra Back Installed Back Height (From Pan to Top of Back)		
☐ 250JR50	22" (16" shell height)	NC
☐ 250JR51	23" (16" shell height)	NC
☐ 250JR52	24" (16" shell height)	NC
☐ 250JR53	24" (18" shell height)	NC
☐ 250JR54	25" (18" shell height)	NC
☐ 250JR55	26" (18" shell height)	NC
☐ 250JR56	26" (20" shell height)	NC
☐ 250JR57	27" (20" shell height)	NC
☐ 250JR58	28" (20" shell height)	NC
☐ 250JR67	28" (22" shell height)	NC
☐ 250JR68	29" (22" shell height)	NC
☐ 250JR69	30" (22" shell height)	NC

B - JAY Syncra Back Upholstery		
☐ 250BU96	Syncra Premium Leatherette	NC
☐ 250BU97	Syncra Stretch	NC

Powered Seating

C - JAY J3 Back (Seating Benefits)

C - JAY J3 Back - Contour		
☐ 250JP01	J3 PA Back for 14"-20" Wide	\$912.00
☐ 250JP02	J3 PA Back for 22" Wide	\$1,941.00
☐ 250JP11	J3 PL Back for 14"-20" Wide	\$986.00
☐ 250JP13	J3 PD Back for 14"-20" Wide (18" and 21" shell height only)	\$1,068.00

C - JAY J3 Back Width		
<i>Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.</i>		
☐ 250JQ02	16"	NC
☐ 250JQ03	18"	NC
☐ 250JQ04	20"	NC
☐ 250JQ05	22"	NC

C - JAY J3 Back Shell Options and Height Ranges.

JAY J3 shells are mounted with a minimim 4" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field with-in the ranges shown, based on back frame type.

JAY J3 Shell Height			
	Tool Angle Adjust.	150° Recline	Pro Advanced
☐ 250JR41	22"-23"	22"-23"	21"-24"
☐ 250JR14	21"-23"	22"-24"	21"-24"
☐ 250JR15	22"-24"	23"-25"	22"-26"
☐ 250JR16	24"-26"	24"-26"	23"-27"
☐ 250JR18	24"	24"	21"-25"
☐ 250JR17	24"-26"	25"-28"	25"-28"
☐ 250JR27	26"-28"	27"-30"	26"-29"
☐ 250JR28	27"-29"	28"-30"	28"-30"

Installed Back Height

JAY J3 Installed Back Height (From Pan to Top of Back)		
☐ 250JR29	21"	NC
☐ 250JR19	22"	NC
☐ 250JR20	23"	NC
☐ 250JR21	24"	NC
☐ 250JR22	25"	NC
☐ 250JR23	26"	NC
☐ 250JR24	27"	NC
☐ 250JR25	28"	NC
☐ 250JR30	29"	NC
☐ 250JR31	30"	NC

Powered Seating

OPTION	Description	Retail
--------	-------------	--------

C - JAY Anterior Thoracic Support (Seating Benefits)

Classic-Structure Pull Down Harness (4 point strap, Only available with JAY backrest)

☐	250K10	Medium (ATS8302)	\$322.00
☐	250K11	Large (ATS8303)	\$322.00

JAY Lateral Thoracic Swing-Away Supports (Requires Pad, Upholstery and Mount)

JAY Lateral Thoracic Support Pads

☐	250JW35	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right	\$75.00
☐	250JW38	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left	\$75.00
☐	250JW36	Lateral Thoracic (L x H) 100 mm x 152 mm / 4" x 6" Pad - Right	\$75.00
☐	250JW39	Lateral Thoracic (L x H) 100 mm x 152 mm / 4" x 6" Pad - Left	\$75.00
☐	250JW37	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right	\$75.00
☐	250JW40	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left	\$75.00

JAY Lateral Thoracic Support Pad Upholstery

☐	250NU01	Premium Leatherette	NC
☐	250NU02	Stretch	NC

JAY Lateral Thoracic Swing-Away Mount

☐	250JW41	Lateral Thoracic Swing Away Bracket - Right	\$334.00
☐	250JW42	Lateral Thoracic Swing Away Bracket - Left	\$334.00

JAY Syncra Fixed Lateral Support Accessory

JAY Fixed Lateral Support Accessory Hardware (JAY® Fixed Lateral Support is height adjustable)

☐	250JW33	Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm / 3.5" x 8" - Right	\$334.00
☐	250JW34	Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm / 3.5" x 8" - Left	\$334.00

*JAY Fixed Lateral Support is height adjustable

Headrests (Requires Pad, Upholstery and Mount)

Headrest Pad

☐	250HE36	SEDEO Pad Small With Cover (L x H) 240 x 150 mm / 10" x 6" <i>Can only be used with Multi-Position Headrest Mount and SEDEO Upholstery</i>	\$227.00
☐	250HE38	Whitmyer Plush Headrest Pad 203 mm / 8"	\$377.00
☐	250HE39	Whitmyer Plush Headrest Pad 254 mm / 10"	\$377.00
☐	250HE44	Whitmyer Plush Headrest Pad 356 mm / 14"	\$377.00
☐	250HE46	Whitmyer Adjust-A-Plush Headrest Pad Medium (Adjustable from 7" to 16" Wide x 4" Tall)	\$445.00
☐	250HE48	OMIT Headrest Pad	NC

Headrest Pad Upholstery

☐	250NU10	SEDEO Upholstery	NC
☐	250NU11	Whitmyer Lycra	NC

Powered Seating

OPTION	Description	Retail
--------	-------------	--------

Headrest Mount

☐	250MH32	SEDEO Multi-Position Headrest Mount with Quick Release (<i>Only available with SEDEO Pad</i>)	\$557.00
☐	250MH33	Whitmyer Linx Headrest Mount with Quick Release	\$247.00
☐	250MH34	Whitmyer Linx-2 Headrest Mount with Quick Release	\$298.00
☐	250MH35	Whitmyer AXYS Headrest Mount with Quick Release	\$247.00
☐	250MH17	OMIT Headrest Mount	NC

Push Handles

☐	250PH27	Height Adjustable Push Handles	\$61.00
--------------------------	---------	--------------------------------	---------

Seat Frame Rails

☐	250SN14	Seat Rail Short (<i>For Effective Seat Depth 14" - 18"</i>)	NC
☐	250SN15	Seat Rail Long (<i>For Effective Seat Depth 17" - 22"</i>)	NC

Seat Cushion Type

☐	250ST79	SEDEO (<i>See Section D</i>)	NC
☐	250ST81	JAY Syncra Cushion (<i>See Section E</i>)	NC
☐	250ST80	JAY (<i>See Section F</i>)	NC
☐	250ST75	Omit	NC

C - Seat Size

Seat Pan Width (*Odd sizes not available on SEDEO cushions*)

☐	250R53	16"	NC
☐	250R54	17"	NC
☐	250R55	18"	NC
☐	250R56	19"	NC
☐	250R57	20"	NC
☐	250R58	21" (<i>N/A for JAY Syncra Cushions</i>)	NC
☐	250R59	22"	NC

Seat Pan Depth (*Front of Seat Pan to front of back cushion. Effective Depth*)

☐	250RF52	14" (<i>N/A with Pro Advanced</i>)	NC
☐	250RF53	15" (<i>N/A with Pro Advanced</i>)	NC
☐	250RF51	16" (<i>N/A with Pro Advanced</i>)	NC
☐	250RF45	17"	NC
☐	250RF46	18"	NC
☐	250RF47	19"	NC
☐	250RF48	20"	NC
☐	250RF49	21" (<i>N/A for JAY Syncra Cushions</i>)	NC
☐	250RF50	22"	NC

Powered Seating

OPTION	Description	Retail
--------	-------------	--------

D - JAY Syncra Cushion (*Jay Your Way Modifications Available Through Customer Service (800) 263-3390*)

☐	250CU89	JAY Syncra Cushion featuring JAY Cryo/Flow Fluid (<i>16"-20" wide</i>)	\$1,017.00
☐	250CU90	JAY Syncra Cushion featuring JAY Cryo/Flow Fluid (<i>22" wide</i>)	\$1,125.00

JAY Syncra Cushion Cover Upholstery

☐	250YR73	Syncra Premium Leatherette	NC
☐	250YR74	Syncra Stretch	NC

F - JAY Seat Cushion (*Seating Benefits*)

☐	250CU59	JAY Ion (<i>Sizes 16" - 21" Wide</i>)	\$657.00
☐	250CU60	JAY Ion (<i>Size 22" Wide</i>)	\$698.00
☐	250CU61	JAY Union (<i>Sizes 16" - 21" Wide</i>)	\$764.00
☐	250CU62	JAY Union (<i>Size 22" Wide</i>)	\$806.00
☐	250CU64	JAY Fusion (<i>Sizes 16" - 21" Wide</i>)	\$1,045.00
☐	250CU65	JAY Fusion (<i>Size 22" Wide</i>)	\$1,188.00

Positioning Supports & Accessories

OPTION	Description	Retail
--------	-------------	--------

Knee and Thigh Swing Away Support Pads (Seating Benefits)

Lateral Knee Support Pads and Hardware (Must pick pad and mount)

☐	250XY05	Lateral Knee Support Pad, Right (4" round)	\$247.00
☐	250XY06	Lateral Knee Support Pad, Left (4" round)	\$247.00
☐	250XX08	Lateral Swing Away Bracket Knee - Right	\$325.00
☐	250XX09	Lateral Swing Away Bracket Knee - Left	\$325.00

Lateral Thigh Support Pads and Hardware (Must pick pad and mount)

☐	250KF22	Lateral Thigh Support Pad Short, Right (5" H x 5" L)	\$146.00
☐	250KF23	Lateral Thigh Support Pad Short, Left (5" H x 5" L)	\$146.00
☐	250KF24	Lateral Thigh Support Pad Long, Right (5" H x 7" L)	\$146.00
☐	250KF25	Lateral Thigh Support Pad Long, Left (5" H x 7" L)	\$146.00
☐	250KF26	Lateral Thigh Removable Mount - Right	\$325.00
☐	250KF27	Lateral Thigh Removable Mount - Left	\$325.00

JAY Lateral Thigh Supports (Must pick pad, upholstery and mount) (Seating Benefits)

JAY Lateral Thigh Support Pads

☐	250KF16	Lateral Thigh Support Pad Standard, Right (8.5" x 3.25")	\$146.00
☐	250KF17	Lateral Thigh Support Pad Standard, Left (8.5" x 3.25")	\$146.00
☐	250KF18	Lateral Thigh Support Pad Large, Right (12" x 3.25")	\$146.00
☐	250KF19	Lateral Thigh Support Pad Large, Left (12" x 3.25")	\$146.00

JAY Lateral Thigh Support Pad Upholstery

☐	250NU07	Syncra Premium Leatherette	NC
☐	250NU08	Syncra Stretch	NC

JAY Lateral Thigh Removable Mount

☐	250KF20	Lateral Thigh Removable Mount - Right	\$325.00
☐	250KF21	Lateral Thigh Removable Mount - Left	\$325.00

Armrest Structure

☐	250EG35	Reclining, Flip Back, Long (Height Adjustable 9" to 15")	NC
☐	250EG41	Reclining, Flip Back, Short (Best used with shorter seat depths) (Height Adjustable 9" to 15")	NC
☐	250EG36	Reclining, Flip Back with 2-point post locking arm (Height Adjustable 11" to 15" x .625")	\$351.00
☐	250EG51	Reclining, Flip Back, Dual Link Armrest (Required for PRO ADVANCED) (No Charge for Multi-Power) (Height Adjustable 8" to 15" including length and lateral adjustment)	NC
☐	250EG52	Reclining, Flip Back, Dual Link Armrest (Height Adjustable 8" to 15" including length and lateral adjustment)	\$450.00
☐	250EG32	Single Post, Flip Back (N/A w/ thigh support)(Height Adjustable 10" to 15")	NC

Arm Pads

SEDEO Armrest Pad

☐	250AB43	300 mm / 12", Integral Foam, Right	NC
☐	250AE43	300 mm / 12", Integral Foam, Left	NC
☐	250AB44	400 mm / 16", Integral Foam, Right	NC
☐	250AE44	400 mm / 16", Integral Foam, Left	NC

Comfort Arm Pad (Foam)

☐	250AB46	2" Wide Right	\$113.00
☐	250AE46	2" Wide Left	\$113.00

Comfort Arm Pad (Gel)

☐	250AB48	2" Wide Right	\$213.00
☐	250AE48	2" Wide Left	\$213.00

Positioning Supports & Accessories

OPTION		Description	Retail
Ottobock Arm Troughs			
☐	250AE06	Ottobock Arm Trough, Med (12"), Left	\$139.00
☐	250AB06	Ottobock Arm Trough, Med (12"), Right	\$139.00
☐	250AE07	Ottobock Arm Trough, Large (14"), Left	\$139.00
☐	250AB07	Ottobock Arm Trough, Large (14"), Right	\$139.00
☐	250OL01	Ottobock Hand Pad, Med (7"), Left	\$70.00
☐	250OR01	Ottobock Hand Pad, Med (7"), Right	\$70.00
☐	250OL02	Ottobock Hand Pad, Large (9"), Left	\$70.00
☐	250OR02	Ottobock Hand Pad, Large (9"), Right	\$70.00
☐	250AE78	Ottobock 1 Piece Arm Trough, (18"), Left	\$198.00
☐	250AB78	Ottobock 1 Piece Arm Trough, (18"), Right	\$198.00
Elbow Support			
☐	2503A02	Elbow Support Pad, RH 3" x 10" <i>(Includes side support)</i>	\$228.00
☐	2503A03	Elbow Support Pad, LH 3" x 10" <i>(Includes side support)</i>	\$228.00
☐	2503A06	Elbow Support Pad, RH 3" x 3"	\$149.00
☐	2503A07	Elbow Support Pad, LH 3" x 3"	\$149.00
☐	2503A04	Swing Away Elbow Support Hardware RH	\$245.00
☐	2503A05	Swing Away Elbow Support Hardware LH	\$245.00
Side Guards (Pair) (Client Cost)			
☐	250SG44	Side Guards Integral Foam (11"L x 5.5"H) <i>(N/A with Reclining Armrests)</i>	\$356.00
Fixed Centermount Legrest			
☐	250H184	Fixed Centermount Legrest, Short <i>(Includes Calf Pads and Single Footplate)(11.5" to 14.5" Leg Length)</i>	NC
☐	250H185	Fixed Centermount Legrest, Long <i>(Includes Calf Pads and Single Footplate)(13.5" to 18.5" Leg Length)</i>	NC
☐	250H183	Sedeo Fixed Centermount Legrest (3" to 17.5" Leg Length)	\$351.00
☐	250H136	Angle Adjustable Centermount Legrest <i>(Includes Calf Pads and Single Footplate)(11.5" to 18.5" Leg Length)</i>	\$605.00
Power Centermount Legrest			
☐	250H175	SEDEO Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension <i>(Mechanical Articulation) (5" to 17.5" Leg Length) (Required for Pro Advanced)</i>	\$2,487.00
☐	250H109	SEDEO Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension <i>(Powered Articulation) (3" to 17.5" Leg Length)</i>	\$3,435.00
SEDEO Centermount Legrest Accessories			
SEDEO Centermount Footplates			
☐	250F127	Dual, Height, Width and Angle Adjustable, Standard (14"W x 10"D) - With Side Supports	\$284.00
☐	250F128	Dual, Height, Width and Angle Adjustable, Large (16"W x 12"D)	\$284.00
☐	250DI71	Flat Footplate <i>(Vertical Edges removed from SEDEO PCM Footplates, ideal for mounting footboxes)</i>	\$274.00
SEDEO Footplate Accessories			
☐	250DI68	Extra Large Centermount Footplate	\$315.00
☐	250DI69	Footplate Connector Link	\$96.00
☐	250DI57	Gel Footpads	\$401.00
SEDEO Centermount Calf Pad Options			
☐	250DI34	Calf Pads, Pair, Standard Integral Foam	NC
☐	250DI35	Calf Pads, Pair, Padded	\$356.00
Knee Block (Required for Anterior Tilt)			
☐	250DI37	Knee Blocks: Dual, Height, Width and Angle Adjustable	\$704.00
☐	250DI55	Knee Block Quick Release Removable Hardware	\$179.00

OPTION	Description	Retail
SEDEO PRO Swing-Away Legrest (Pair) (N/A with Pro Advanced) (250H179 Client Cost)		
☐ 250H177	Legrests, Swing Away, 70°, SEDEO PRO (14" to 21" Leg Length)	NC
☐ 250H178	Legrests, Manual Elevating, Swing Away, SEDEO PRO (Includes Padded Calf Pads) (14".5 to 21.5" Leg Length)	\$605.00
☐ 250H179	Legrests, Power Elevating, Swing Away, SEDEO PRO (Includes Padded Calf Pads) (15.5" to 22.5" Leg Length)	\$3,435.00
☐ 250H176	Omit Legrests	NC

Swing Away Legrest Footplates

☐ 250F369	Angle Adjustable Flip-Up Footplates	\$284.00
----------------------------------	-------------------------------------	----------

Positioning Supports & Accessories

Swing Away Legrest Calf Supports

☐ 250DI42	Calf Support, Depth Adjustable, Padded (6.5" to 8" Wide, Depending On Chair Width, x 9.5" Tall)	\$512.00
----------------------------------	---	----------

Swing Away Legrest Options

☐ 250DI44	Heel Loops	\$0.00
☐ 250KG4	Knee Pads, Integral Foam (2.5" Wide x 4.5" Tall)	\$282.00

SEDEO Lower Extremity Swing Away Support

☐ 250VD01	Lower Extremity Support, Padded, Swing Away, SEDEO, Right	\$512.00
☐ 250VD02	Lower Extremity Support, Padded, Swing Away, SEDEO, Left	\$512.00

Pelvic Positioning Belts

☐ 250A837	Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
☐ 250A839	Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"	NC
☐ 250A840	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$103.00
☐ 250A845	Bodypoint - Evoflex Pelvic Stabilizer Belt	\$216.00
☐ 250A846	Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	\$271.00

Chest Straps

☐ 250ZR06	Chest Strap, 2-Point, Padded (Required)	NC
☐ 250ZR07	Additional Chest Strap, 2-Point, Padded	\$351.00

Chest Bars

☐ 250ZR08	Chest Bar, Detachable, Round Pad	\$704.00
☐ 250ZR09	Chest Bar, Detachable, Flat Pad, 200 mm x 100mm	\$748.00

Accessories

Ventilator Option

☐ 250VT12	Vent Tray Fixed (Pulmonetic LTV Models)	\$609.00
----------------------------------	---	----------

O2 Tank Holder

☐ 250A718	O2 Tank Holder Mounted to Backrest (D and E Styles)	\$391.00
----------------------------------	---	----------

Wheelchair Tray Table

☐ 2507A09	Tray, Swing Away	\$381.00
----------------------------------	------------------	----------

Backpack (Client Cost)

☐ 250A1	Black Utility	\$112.00
--------------------------------	---------------	----------

Medial Thigh Support

☐ 2508A80	Medial Thigh Support	\$271.00
☐ 2508A81	Pommel Swing Away Hardware	\$390.00

Misc. Accessories (Client Cost)

☐ 2508A59	Transfer Bar (Pair)	\$493.00
☐ 2508A52	Ram Mount Phone	\$77.00
☐ 2508A53	Ram Mount Tablet	\$308.00
☐ 2508A54	Ram Mount Cup Holder	\$77.00
☐ 2508A58	Rear View Mirror	\$77.00
☐ 250A611	Crutch Holder (N/A w/ Power Recline)	\$174.00
☐ 2508A73	Bag For Armrest	\$77.00
☐ 2508A75	Reflectors (Backrest shroud mounted)	\$115.00
☐ 2509A10	Utility Hooks	\$30.00
☐ 2508A93	Wheelchair Cover	#N/A

Joystick / Electronics			
OPTION		Description	Retail
R-net Joystick (Remote)			
☐	250J121	Standard LED Joystick (w/o Screen Display)	NC
☐	250J116	Color Joystick (2.3" Color LCD Screen Display, CJSM1)	NC
☐	250J118	High Resolution Color Joystick,(CJSM2+BT) 3.5" Color LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control (Client Cost)	\$1,162.00
☐	250J119	Omit Remote Control (Required with arm mounted specialty controls)	NC
Joystick Mounting			
☐	250JM34	Fixed on Armrest	NC
☐	250JM35	Parallel Swing Away	\$447.00
☐	250JM36	Link-It Mid Line Mount on Armrest (Swing Away)	\$447.00
☐	250JM49	Attendant (Available with LED or CJSM1 only, Mounted to Push Handles)	NC
Joystick Mount (Location)			
☐	250HM16	Right Hand	NC
☐	250HM17	Left Hand	NC
Display Language Selection R-net			
☐	250L1	English	NC
☐	250L3	French	NC
☐	250L7	Spanish	NC
Joystick Handle (Optional) N/A with HD Joystick			
☐	250JH13	Ball Handle	\$110.00
☐	250JH15	T-Handle	\$110.00
☐	250JH24	Bodypoint - Rubber Dome Joystick Handle	\$216.00
☐	250JH25	Bodypoint - U-Shaped Joystick Handle (4 in)	\$216.00
☐	250JH26	Bodypoint - Grooved Mushroom Joystick	\$216.00
Electronic Programmers			
☐	250PR13	R-net PC Wireless Programmer	\$660.00
R-net Power Modules			
☐	250ER11	Power Module, R-net, 120A, Non-Expandable (Only available with Group 4)	NC
☐	250ER10	Power Module, R-net, 120A, Expandable (Only available with Group 4)	\$102.00
	250EX09	Expandable 120A Controller - R-net	\$0.00
	250EX10	Harness for 120A Expandable Controller	\$0.00
☐	250ER18	Power Module, R-net, 90A, Non-Expandable (Available with No Power and Single Power Options)	NC
☐	250ER17	Power Module, R-net, 90A, Expandable (Required for specialty controls or 3+ actuators)	\$102.00
	250EX15	Expandable 90A Controller - R-net	\$0.00
	250EX16	Harness for 90A Expandable Controller	\$0.00

Joystick / Electronics

OPTION	Description	Retail
--------	-------------	--------

Seating Control

Powered Seating Control

☐	250EO40	Omit Thru-drive <i>(Requires Switch Control)</i>	NC
☐	250EO37	Thru-drive Control for 1 Actuator	\$1,443.00
☐	250EO38	Thru-drive Control for 2+ Actuators	\$2,903.00
☐	250EO32	CxSM, Actuator control module, Thru Drive 2+ Actuators <i>(Required for Pro Advanced)</i>	\$2,903.00

Switch Control *(Available with Omit Thru-Drive Only) (1 actuator for Tilt or Elevate Only)*

☐	250EO33	Dual Toggle for 1 actuator	NC
☐	250EO34	2.5" Single Button Switch for 1 actuator <i>(Button Switches are not mounted)</i>	NC
☐	250EO41	Ctrl+5 Switch Box with Voice Activation for 2 actuators or more	NC

Additional Power Seating Control *(Available with Thru-Drive Only)*

☐	250QG11	Ctrl+5 Switch Box With Voice Activation	\$1,078.00
☐	250QG12	Ctrl+5 Switch Box With Voice Activation <i>(Multi-Power only)</i>	NC

Dual Toggle Mounting Location

☐	250TN1	Toggle Mount, Right	NC
☐	250TN2	Toggle Mount, Left	NC

Ctrl+5 Location *(If mounted on same side as joystick, it will Swing Away)*

☐	250EY19	Right Armrest	NC
☐	250EY20	Left Armrest	NC
☐	250EY21	Back Cane, Right	\$0.00
☐	250EY24	Back Cane, Left	\$0.00

Ctrl+5 Switch Type

☐	250EY22	Push Buttons	NC
☐	250EY23	Toggle Switches	\$356.00

Attendant Controls *(Prior Approval Required)*

☐	2500A05	Dual Control Module for Attendant	\$1,225.00
☐	2500A06	R-net Color LCD (CJSM1), Attendant Control Only	\$1,225.00

Attendant Mounting Location

☐	250DN5	Right Back Cane	NC
☐	250DN6	Left Back Cane	NC

R-net Programming

Assignable Buttons

Assignable Button feature allows for customization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

☐	250PM12	Include Assignable buttons	NC
--------------------------	---------	----------------------------	----

Power Center Mount Footplate Extension Option *(Client Cost)*

☐	250IE03	Include 7.5" Powered footplate extension for transfer assist <i>(Electronics programmed to allow footplates to be extended with or without legrest being raised.) Only Available with Powered Articulation</i>	\$592.00
--------------------------	---------	--	----------

Joystick / Electronics

R-net Omni2 Options *(For Expandable Controls, Choose All That Apply)*

☐	250ON14	R-net Omni2 Only Input Control Module and 3.5" LCD Display w/ Built-In Infrared & Bluetooth Choose Mount Below	\$2,384.00
☐	250ON15	R-net Omni2 Bundled Input Control Module and 3.5" LCD Display w/ Built-In Infrared & Bluetooth Bundled Omni2 Available w/ Switch-It Input Devices	NC

Omni2 Mount *If used w/ Standard Joystick, it must be mounted opposite side.*

☐	250DS14	Link-It with Parallel Swing-Away Omni2 Mount w/ Rear Mounted Interface, Right	\$698.00
☐	250DS15	Link-It with Parallel Swing-Away Omni2 Mount w/ Rear Mounted Interface, Left	\$698.00

OPTION	Description	Retail
Pro Advanced Programming		
Biometric Repositioning		
☐ 250PM03	Include Biometric Repositioning (Requires Memory Position #1, #2, #3 Below)	NC
Memory Seating		
<i>"Memory Seating" feature allows for up to 6 custom seating positions that are accessed either thru the joystick, Ctrl+5, Joystick Shortcut keys or programmable jacks on Omni2 and joystick. Choose desired preset settings below. Default positions are dealer adjustable.</i>		
☐ 250PM06	Include #1 Memory Seating (Seated) (Required)	NC
☐ 250PM07	Include #2 Memory Seating (Supine)	NC
☐ 250PM08	Include #3 Memory Seating (Pressure Relief)	NC
☐ 250PM09	Include #4 Memory Seating (Feet Above Heart)	NC
☐ 250PM10	Include #5 Memory Seating (Transfer)	NC
☐ 250PM11	Include #6 Memory Seating (C-Me)	NC
PRAR Pressure Relief Alert and Recording		
<i>PRAR Pressure Relief Alert and Record feature provides programmable alerts to be set into the Rnet electronics that will prompt the user when it is time to perform a pressure relief. The R-net system will record how often the pressure relief is completed or ignored. Requires Memory Seating function.</i>		
☐ 250PM14	Include PRAR (Minimum Requirement Memory Position #1 and #3)	NC

Specialty Input Devices

Hand Control Packages (Requires Standard Joystick Omit / Only one hand control device can be ordered) (AADL W100 Grants)

Proportional Input Devices

☐ 250EU14	Switch-It MicroGuide	\$4,574.00
☐ 250EU17	Switch-It VersaGuide	\$2,557.00
☐ 250EU18	Switch-It VersaGuide EZ	\$3,023.00

Mounting Position (Required)

☐ 250EM04	Link-It to the Backrest (Not Available w/ Omit Back Cushion Type)	\$1,338.00
☐ 250EM05	Link-It to the Fixed Joystick Mount	\$416.00
☐ 250EM06	Link-It to the Parallel Swing Away Joystick Mount	\$416.00
☐ 250EM07	Link-It to the Midline Swing Away	\$1,480.00
☐ 250EM08	Omit Mount	NC

Mounting Platforms (Optional) Not available with Link-It to the Backrest

☐ 250MG16	SunDisk Hand Tray (Only available w/ MicroSeries Joysticks)	\$416.00
☐ 250MG15	Bullet Hand Tray	\$416.00

Mounting Location (Required)

☐ 250IM09	Right	NC
☐ 250IM10	Left	NC

Switch Options

See Switch Selections (Maximum 2)

OPTION	Description	Retail
Head Control Packages (Select headrest mount on next page) (AADL W100 Grants)		
Only one head control device can be ordered.		
Proportional Input Devices		
☐ 250EU25	Switch-It DualPro Proportional Head Array - Facial Pads	\$7,057.00
☐ 250EU26	Switch-It DualPro Proportional Head Array - Spot Pads	\$7,057.00
☐ 250EU51	Vigo Head Drive	\$6,293.00
Switched Input Devices		
☐ 250EU28	Switch-It 3-Switch Head Array - Short Curved Wings	\$6,154.00
☐ 250EU30	Switch-It 3-Switch Head Array - Long Curved Wings	\$6,154.00
Switch Options		
☐ -	See Switch Selections (Maximum 2)	-
Sip 'N Puff (AADL W100 Grants)		
☐ 250ID73	Switch-It Sip 'N Puff Combo Head Array with Link-It Hardware (Requires headrest mount)	\$7,037.00
☐ 250ID74	Switch-It OptiStop with Sip 'N Puff	\$2,741.00
Headrest Mounts (Applies to Head/Sip 'N Puff control packages) (AADL W100 Grants)		
☐ 250IJ13	Whitmyer Linx Headrest Mount w/ Quick Release	\$247.00
☐ 250IJ14	Whitmyer Linx-2 Headrest Mount w/ Quick Release	\$298.00
☐ 250IJ15	Whitmyer AXYS Headrest Mount w/ Quick Release	\$247.00
Specialty Input Devices		
Chin Control Packages (AADL W100 Grant)		
Proportional Input Devices		
☐ 250EU34	Switch-It MicroGuide (Switches mounted to the headrest only.)	\$4,574.00
☐ 250EU36	Switch-It VersaGuide	\$2,557.00
☐ 250EU37	Switch-It VersaGuide EZ	\$3,023.00
Mounting Position		
☐ 250EM09	Link-It To Chin Harness, Small (150mm / 6")	\$1,760.00
☐ 250EM10	Link-It To Chin Harness, Large (220mm / 8")	\$1,760.00
☐ 250EM11	Link-It To Headrest (Swing Away)	\$1,760.00
☐ 250EM16	Omit Mount	NC
Input Device Mounting Location		
☐ 250IM13	Right	NC
☐ 250IM14	Left	NC
Switch Options		
☐ -	See Switch Selections (Maximum 2)	-
Safety Devices		
☐ 250SA06	Switch-It Remote Stop Switch - XLR Connector	\$1,317.00
Switched Driving Packages		
☐ 250SW44	Switch-It Cool Cube with DB9 connector for proximity / Proportional switch driving. (Pro Spot Proportional Switch's are required if proportional driving is desired) (For Non-Proportional driving, Single Switches can be used)(For Future Development)	\$3,949.00

Specialty Input Devices

OPTION	Description	Retail
--------	-------------	--------

Switches

Single switches can be programmed to activate any function of the wheelchair including: On/Off, Mode, Profile, Speeds, Seating, Memory Positions, Horn, Lights & Indicators, IR functions, Bluetooth mouse clicks, Environmental Controls when plugged into assignable port.

Switches can connect via: Joysticks (Standard, Micro series, Versa series) OMNI, SJX6 Hub, Cool Cube, Head controls

Switch ordering instructions:

1. Select Switch type and quantity in section "A" using columns below. (Max 6 switches. see limitations below)
2. Select Switch Mounting location from Section "B" using corresponding switch column.
3. Select Switch Mounting Side from Section "C" using corresponding switch column.

Section A - Switch Selections

Single Function Switches - Mechanical Switches

2.5" Single Switch - Black

2.5" Single Switch - Blue

2.5" Single Switch - Green

2.5" Single Switch - Red

2.5" Single Switch - Yellow

Egg Switch - Black

Egg Switch - Blue

Egg Switch - Green

Egg Switch - Red

Egg Switch - Yellow

MicroLite (Includes color labels)

Mini Button 1 - Blue (Velcro)

Mini Button 1 - Green (Velcro)

Mini Button 1 - Red (Velcro)

1" Single Switch - Black

1" Single Switch - Blue

1" Single Switch - Green

1" Single Switch - Orange

1" Single Switch - Red

1" Single Switch - Yellow

Switch / Quantity

Switch 1

\$

Switch 2

\$

☐

250IR31

\$110.00

☐

250IR56

\$110.00

☐

250IR32

\$110.00

☐

250IR57

\$110.00

☐

250IR33

\$110.00

☐

250IR58

\$110.00

☐

250IR34

\$110.00

☐

250IR59

\$110.00

☐

250IR35

\$110.00

☐

250IR60

\$110.00

☐

250IR36

\$110.00

☐

250IR61

\$110.00

☐

250IR37

\$110.00

☐

250IR62

\$110.00

☐

250IR38

\$110.00

☐

250IR63

\$110.00

☐

250IR39

\$110.00

☐

250IR64

\$110.00

☐

250IR40

\$110.00

☐

250IR65

\$110.00

☐

250IR41

\$115.00

☐

250IR66

\$115.00

☐

250IR42

\$308.00

☐

250IR67

\$308.00

☐

250IR43

\$308.00

☐

250IR68

\$308.00

☐

250IR44

\$308.00

☐

250IR69

\$308.00

☐

250IR45

\$388.00

☐

250IR70

\$388.00

☐

250IR46

\$388.00

☐

250IR71

\$388.00

☐

250IR47

\$388.00

☐

250IR72

\$388.00

☐

250IR48

\$388.00

☐

250IR73

\$388.00

☐

250IR49

\$388.00

☐

250IR74

\$388.00

☐

250IR50

\$388.00

☐

250IR75

\$388.00

Electronic Switches

Flex Switch (Headrest mount only)

☐

250IR54

\$303.00

☐

250IR79

\$303.00

Pro Spot - Proportional Switches

Pro Spot Proportional Switch - Red

☐

250IR181

\$542.00

☐

250IR186

\$542.00

Pro Spot Proportional Switch - Green

☐

250IR182

\$542.00

☐

250IR187

\$542.00

Pro Spot Proportional Switch - Blue

☐

250IR183

\$542.00

☐

250IR188

\$542.00

Pro Spot Proportional Switch - Yellow

☐

250IR184

\$542.00

☐

250IR189

\$542.00

Pro Spot Proportional Switch - Orange

☐

250IR185

\$542.00

☐

250IR190

\$542.00

Section B - Switch Mounting Position

Link-It to the Backrest (Up to 1 per side)

☐

250EM17

\$1,338.00

☐

250EM23

\$1,338.00

Headrest Mount - Gooseneck & Rigid (up to 2 per side)

☐

250EM18

\$1,760.00

☐

250EM24

\$1,760.00

Link-It to the Fixed Joystick Mount (Up to 2 per side)

☐

250EM19

\$416.00

☐

250EM25

\$416.00

Link-It to the Parallel Swing Away Joystick Mount (Up to 2 per side)

☐

250EM20

\$416.00

☐

250EM26

\$416.00

Link-It to the Seat Rail (Up to 2 per side)
(N/A w/ Knee & Thigh Supports)

☐

250EM21

\$556.00

☐

250EM27

\$556.00

Omit Mount (includes Velcro)

☐

250EM22

NC

☐

250EM28

NC

Switches Continued			
OPTION	Description		Retail

Section C - Switch Mounting Location

Right side of chair



250IM17 NC



250IM19 NC

Left side of chair



250IM18 NC



250IM20 NC

Joystick Integrated Switches (Pair)

☐ 250JO03 Proximity Switches on Flex Mount, Pair (Used with VersaGuide/VersaGuide EZ)

\$1,403.00

Syncra Seating - Measuring Guide

JAY Syncra Back



Backrest Width

Max User Trunk Width

16"	18"	20"	22"
13.5"	15.5"	17.5"	19.5"

JAY Syncra Cushion



JAY Swing-Away Lateral Thoracic Adjustment Range

JAY Syncra Back Width	Width Range	Height range from bottom of shell to seat pan
16"	10" - 24"	12" - 19.5"
18"	12" - 26"	12.5" - 21.5"
20"	14" - 28"	12.5" - 22"
22"	16" - 30"	13.5" - 23"


JAY Swing-Away Lateral Thoracic Adjustment Range

JAY Syncra Cushion Width	Width Range
16"	8.5" - 16.0"
17"	9.5" - 17.0"
18"	10.5" - 18.0"
19"	11.5" - 19.0"
20"	12.5" - 20.0"
22"	14.5" - 22.0"

Syncra Fixed Lateral Support Accessory Range

Syncra Back Width	Width Range (loaded)	Height range from bottom of shell to seat pan	Contour Depth Range
16"	16.4"-16.8"	6"-8"	3"-5"
18"	18.4"-18.8"	6"-8"	3"-5"
20"	20.4"-20.8"	6"-8"	3"-5"
22"	22.4"-22.8"	6"-8"	3"-5"

Syncra Back Mounting Height Diagram



Customer Service: (800) 263-3390
E-mail: cscanada@sunmed.com
www.sunrisemedical.ca
Specifications are subject to change without notice.

© 10/2026 Sunrise Medical Canada Inc. MK129630-CAN Rev 27