



Quickie GP/GPV - AADL W133B

October 2026

EIGP70

Mark For:

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ Prov: _____ PC: _____

Dealer Phone: () _____ Fax: () _____

Confirmation Email: _____

Confirm Via: ☐ Fax ☐ Email

Submitting for:

☐ Quote ☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ Prov: _____ PC: _____

Ship To Phone: () _____ Fax: () _____

✓ Option # Definition Price ✓ Option # Definition Price

1. FRAME

Base Model

☐	024WM23	Quickie GP-70°	\$2,828.00
☐	024WM23	Quickie GPV-70° (NA 12"-13")	\$2,828.00

Frame Angle

☐	024FA6	60° (N/A 18 1/4 and 19 1/4 frame length)	NC
☐	024FA7	70°	NC
☐	024FA1	75° (N/A with rigid footrest)	NC
☐	024FA8	80° (N/A with rigid footrest)	NC

2. SIZE / SEATING

Seat Width

Seat Width 1/2" Narrower

☐	024W3	12"	NC
☐	024W4	13"	NC
☐	024W5	14"	NC
☐	024W6	15"	NC
☐	024W7	16"	NC
☐	024W8	17"	NC
☐	024W9	18"	NC
☐	024W41	19"	\$308.00
☐	024W11	20"	\$308.00
☐	024W12	21"	\$308.00
☐	024W13	22"	\$308.00

Seat Depth

21 & 22" seat depths only available with 16" or wider chairs

☐	024SD5	12"	NC
☐	024SD6	13"	NC
☐	024SD7	14"	NC
☐	024SD8	15"	NC
☐	024SD9	16"	NC
☐	024SD10	17"	NC
☐	024SD11	18"	NC
☐	024SD72	19"	\$308.00
☐	024SD13	20"	\$308.00
☐	024SD14	21"	\$308.00
☐	024SD33	22"	\$308.00

Frame Length

☐	024FL30	18-1/4" (Avail. on 12" Sling Depth)	NC
☐	024FL31	19-1/4" (Avail. on 13" Sling Depth)	NC
☐	024FL32	20-1/4" (Avail. on 14", 15" Sling Depth)	NC
☐	024FL33	21-1/4" (Avail. on 14", 15", 16" Sling Depth)	NC
☐	024FL34	22-1/4" (Avail. on 16", 17" Sling Depth)	NC
☐	024FL35	23-1/4" (Avail. on 17", 18" Sling Depth)	NC
☐	024FL36	24-1/4" (Avail. on 18" Sling Depth)	NC
☐	024FL37	25-1/4" (Avail. on 19" Sling Depth)	NC
☐	024FL38	26-1/4" (Avail. on 20" Sling Depth)	NC
☐	024FL39	27-1/4" (Avail. on 21" Sling Depth)	NC
☐	024FL40	28-1/4" (Avail. on 22" Sling Depth)	NC

Frame Height

See Appendix 11 for height chart. Short frame N/A w/Xtender

☐	024FH2	Short (front 9.25"-rear 8.25") N/A with Xtender Power Assist	\$230.00
☐	024FH5	Standard (front 10" - rear 9")	NC
☐	024FH7	Tall (front 11" - rear 10")	\$230.00

Cushion

☐	024CU9	Reinforced Seat Sling	\$53.00
☐	024CU4	Omit Seat Sling	NC

JAY Bundled Cushions - Optional

To order any JAY cushions on-chair, please see JAY Price List. Cushion name or part number must be filled out with desired width & depth.

☐	024JC20	JAY Bundled Cushion Cushion Name or Part #: Width: _____ Depth: _____	Price List
--------------------------	---------	---	------------

✓	Option #	Definition	Price	✓	Option #	Definition	Price
3. CASTERS							
Caster Wheels				Caster Wheels - Continued			
☐	024CW1	3" Micro	NC	☐	024CW32	8" Polyurethane <i>Only available on 60° frame option</i>	NC
☐	024CW3	4" Low Profile Polyurethane	NC	☐	024CW36	8" x 2" Pneumatic <i>Only available on 60° frame option</i>	\$247.00
☐	024CW31	4" Micro Lighted	\$110.00	☐	024CW37	8" x 2" Pneumatic w/Airless Insert <i>Only available on 60° frame option</i>	\$393.00
☐	024CW6	5" Soft Roll Semi-Pneumatic	\$126.00	Caster Fork			
☐	024CW5	5" Low Profile Polyurethane	NC	☐	024CS5	3/4" Longer Fork Stem Bolt	NC
☐	024CW7	5" x 1.5" Semi-Pneumatic <i>only avail. For use in bottom fork holes</i>	\$110.00	☐	024CS9	1-1/2" Longer Fork Stem Bolt	NC
☐	024CW33	6" Pneumatic	\$110.00	☐	024CS3	Caster Pin Locks	\$152.00
☐	024CW34	6" Polyurethane	NC				
☐	024CW35	6" x 1.5" Semi-Pneumatic <i>only avail. For use in bottom fork holes</i>	\$110.00				
☐	024CW21	8" Pneumatic <i>Only available on 60° frame option</i>	\$146.00				

4. REAR WHEELS							
Xtender Power Assist				Rear Wheel Tire			
<i>RSH 15"-19" with 22" wheels, RSH 16"-20" with 24" wheels. Spare batteries must match set up. Anti-Tips included. Dual Drive mode available w/B4ME</i>				☐	024RT1	Pneumatic	NC
☐	024XT6	22" Xtender Power Assist JWX-2 NICKEL	\$9,996.00	☐	024RT2	Pneumatic Airless Insert <i>N/A with 26" wheels & Xtender</i>	\$247.00
☐	024XT7	22" Xtender Power Assist JWX-2 LITHIUM	\$11,604.00	☐	024RT3	Full Poly <i>N/A with 26" wheels</i>	NC
☐	024XT4	24" Xtender Power Assist JWX-2 NICKEL	\$9,996.00	☐	024RT6	High-Pressure Clincher	\$105.00
☐	024XT5	24" Xtender Power Assist JWX-2 LITHIUM	\$11,604.00	Hand Rim			
<i>Please select one of the following handrim options</i>				☐	024HR1	Aluminum Anodized	NC
☐	024XH5	Plastic Coated	NC	☐	024HR4	Plastic Coated - <i>N/A on all rims</i>	\$146.00
☐	024XH6	Stainless Steel	NC	☐	024HR7	Omit Hand Rim	NC
Xtender - Other Options				Projection Hand Rim			
☐	024XA10	Spare NiMH Power Pack - NICKEL	\$1,848.00	<i>Projections N/A w/ plastic coated handrims. N/A w/ Spoke guards</i>			
☐	024XA11	Spare Li-ion Power Pack - LITHIUM	\$2,379.00	☐	024HR8	8 Oblique (24"/26")	\$146.00
Rear Wheel Size				Axle			
<i>Spinergy N/A w/ Xtender, spoke guards. Contact Built 4 Me for Spinergy Color Options</i>				☐	024AX1	Quick Release Steel	NC
☐	024RW8	22" Mag	NC	☐	024AX11	Quad Quick Release	\$182.00
☐	024RW9	22" Spoke	\$320.00	Axle Plate			
☐	024RW1	24" Mag	NC	☐	024AL2	Standard Adjustable	NC
☐	024RW2	24" Spoke	NC	☐	024AL13	Double Slotted Adjustable	NC
☐	024RW11	26" Spoke (<i>N/A with Xtender Power Assist</i>)	\$89.00	☐	024AL6	Amputee Adjustable <i>N/A with padded swing-away armrests or stroller handles</i>	\$146.00
☐	024RW42	24" Spinergy (540) SPOX Black	\$864.00				
☐	024RR25	Omit Rear Wheel and Axle Intended Size for Set Up :	NC				

5. WHEEL LOCKS							
Wheel Lock				Extension Handle			
☐	024WL1	High Mount Push	NC	☐	024EH1	6" Pair	\$73.00
☐	024WL2	High Mount Pull	NC	Grade Aid			
☐	024WL5	High Mount Scissor <i>N/A w/ extension handles, grade aides</i>	NC	☐	024GA1	Grade Aid <i>N/A w/ poly, high-pressure clincher tires, or Xtender</i>	\$157.00
☐	024WL12	Do Not Mount Wheel Locks	NC				
☐	024WL10	Omit	NC				

✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

6. FOOTREST

Footplates			
<i>Heel loop or leg strap not included. See Footplate Options Section</i>			
1. GP FOOTREST OPTIONS			
☐	024F56	Rigid N/A w/ Ext. Tubes; N/A with 75° or 80° frame angle	NC
☐	024F50	Platform N/A w/ Ext. Tubes	\$152.00
☐	024F1	Composite Footplate	NC
☐	024F3	Adult Angle Adj.	\$176.00
☐	024F48	Platform Flip Up	\$230.00
2. GPV FOOTREST OPTIONS			
☐	024F56	Rigid N/A w/ Ext. Tubes; N/A with 75° or 80° frame angle	NC
☐	024F50	Platform N/A w/ Ext. Tubes	\$152.00
☐	024F48	Platform Flip Up	\$230.00

Footrest Extension Tubes			
☐	024E1	Short (14" to 16-1/2")	NC
☐	024E2	Medium (16-1/2" to 19")	NC
☐	024E5	Long (19" to 21-1/2")	NC
☐	024E16	High Mount Only avail. with adult angle adjustable footplates	\$97.00
Footplate Options (Ext Tube & Direct Mount)			
☐	024HL1	Heel Loop	\$0.00
☐	024HL5	Leg Strap	\$0.00
☐	024HL4	Omit Heel Loop	NC
☐	024HL6	Omit Leg Strap	NC
☐	024F7	Omit Footplate & Ext Tube	NC
☐	024IG2	Impact Guard Neoprene	\$61.00

7. BACKREST

Back Type			
☐	024BT3	Non-Folding Backrest (Fixed) N/A with push handles or adj. upholstery	
☐	024BB1	8" Back Height	NC
☐	024BB2	9" Back Height	NC
☐	024BB3	10" Back Height	NC
☐	024BB4	11" Back Height	NC
☐	024BB5	12" Back Height	NC
☐	024BB6	13" Back Height	NC
☐	024BB7	14" Back Height	NC
☐	024BB8	15" Back Height	NC
☐	024BB9	16" Back Height	NC
☐	024BT1	Folding Backrest (Adjustable)	
☐	024BA1	11" - 14" Back Height	\$121.00
☐	024BA2	12" - 16" Back Height	\$121.00
☐	024BA3	14" - 18" Back Height	\$121.00
☐	024BA4	16" - 20" Back Height	\$121.00

Back Upholstery			
☐	024BU1	Standard	NC
☐	024BU8	Adjustable	\$204.00
☐	024BU16	Reinforced (N/A w/ adj. upholstery)	\$53.00
☐	024BU14	Omit Back Upholstery	NC
Backrest Options			
☐	024BT33	Folding w/ Integral Push Handles	NC
☐	024A85	Bolt-On Push Handle (Folding only)	\$139.00
☐	024SR2	Regular Stroller Handle N/A with amputee axle plate	\$193.00
☐	024SR3	Tall Stroller Handle N/A with amputee axle plate	\$193.00
JAY Back (Optional)			
☐	024JB01	J3 Back (See J3 Order Form)	Price List

8. ARMREST / SIDEGUARDS

Single Post			
☐	024AA9	Single Post Ht. Adjust. Low Must be removed to fold back down	\$339.00
☐	024AM1	Desk Length - Classic (10")	NC
☐	024AM2	Full Length - Classic (14")	NC
☐	024AA8	Single Post Ht. Adjust. Std Must be removed to fold back down	\$339.00
☐	024AM1	Desk Length - Classic (10")	NC
☐	024AM2	Full Length - Classic (14")	NC

Padded Swing Away Armrest			
<i>All Padded Swing-Away N/A with amputee axle plate or Xtender Power Assist</i>			
☐	024AA19	Padded S-A Medium	\$314.00
☐	024AA16	Padded S-A Adult	\$314.00
Side Guard			
☐	024SG3	Plastic Regular	\$188.00
☐	024SG4	Plastic Kids	\$188.00

9. KEY OPTIONS AND ACCESSORIES

Rear Anti-Tip			
☐	024AT1	Rear Anti-Tip	\$121.00
Backpack			
☐	024A1	Quickie Backpack Black	\$112.00
☐	024A2	Quickie Kids Backpack Black	\$105.00
Seat Pouch			
☐	024A8	Quickie Seat Pouch Black	\$83.00
Positioning Belts			
☐	024A142	Belt 1.5" Auto Buckle w/Clamps	\$97.00
☐	024A12	Belt 2" Aircraft Buckle w/Clamps	\$139.00
☐	024A13	Belt 2" Aircraft Buckle Pad w/Clamp	\$204.00

Spoke Guards			
<i>N/A w/ Xtender Power Assist; Only available w/ Mag and Spoke wheels</i>			
☐	024A23	Black Spoke Guards	\$157.00
☐	024A24	Clear Spoke Guards	\$157.00
Crutch Holder			
☐	024A53	Left Side	\$115.00
☐	024A52	Right Side	\$115.00
Extra Leg Strap			
☐	024A55	Extra Leg Strap Black	\$61.00
Miscellaneous			
☐	024A83	Tool Kit	\$43.00
☐	024A999	Quickie Water Bottle & Holder	\$81.00
Tray Table			
☐	024A31	Tray Table	\$381.00

10. COLOR

Standard Frame Colors			
☐	024C211	Celestial Sapphire	NC
☐	024C212	Desert Sage	NC
☐	024C213	Desert Tan	NC
☐	024C214	Gloss Black	NC
☐	024C215	Matte Black	NC
☐	024C216	Metallic Shadow	NC
☐	024C217	Midnight Blue	NC
☐	024C218	Mystic Purple	NC
☐	024C219	Petal Pink	NC
☐	024C220	Pink Power	NC
☐	024C221	Red Pepper	NC
☐	024C222	Royal Blue	NC
☐	024C223	Salute Green	NC
☐	024C224	Scarlet Noir	NC
☐	024C225	Silver Ore	NC
☐	024C226	Solar Orange	NC
☐	024C227	Starlight White	NC
☐	024C228	Storm Gray	NC
☐	024C229	Turbo Green	NC
☐	024C230	Velocity Yellow	NC

Personalized Frame Color			
☐	024C112	Desert Camouflage	\$368.00
☐	024C43	Mossy Oak Camouflage	\$368.00
☐	024C42	Stars & Stripes	\$368.00
☐	024C108	Zebra (Choose base color) N/A w/ Gloss Black and Matte Black.	\$368.00

11. APPENDIX

GP Series Rear STF Chart

*All measurements given are an estimate and have been rounded to the nearest .5"

SHORT Frame (9 1/4" Front / 8 1/4" Rear)			Rear Axle Plate Hole Position					
			Standard			Inverted		
Wheel	Rim	Tire	Top	Middle	Bottom	Top	Middle	Bottom
22"	Mag / Spoke	Full Poly	14.5	15	15.5	16	16.5	17
22"	Mag / Spoke	Pneum / Insert	14.5	15	16	16	17	17.5
24"	Mag / Spoke	Full Poly	15	15.5	16.5	16.5	17	17.5
24"	Spoke / Spiner	Clincher	15	16	16.5	16.5	17	18
24"	Mag / Spoke / Spiner	Pneum / Insert	15.5	16	17	17	17.5	18
26"	Spoke / Spiner	Clincher	16	17	17.5	17.5	18	19
26"	Spoke / Spiner	Pneumatic	16.5	17	17.5	18	18.5	19

TALL Frame (11" Front / 10" Rear)			Rear Axle Plate Hole Position					
			Standard			Inverted		
Wheel	Rim	Tire	Top	Middle	Bottom	Top	Middle	Bottom
22"	Mag / Spoke	Full Poly	16	16.5	17.5	17.5	18	19
22"	Mag / Spoke	Pneum / Insert	16.5	17	17.5	18	18.5	19
24"	Mag / Spoke	Full Poly	17	17.5	18	18.5	18.5	19.5
24"	Spoke / Spiner	Clincher	17	17.5	18	18.5	19	19.5
24"	Mag / Spoke / Spiner	Pneum / Insert	17	17.5	18.5	19	19	20
26"	Spoke	Clincher	18	18.5	19	19.5	20	20.5
26"	Spoke	Pneumatic	18	19	19.5	19.5	20	21

STANDARD Frame (10" Front / 9" Rear)			Rear Axle Plate Hole Position					
			Standard			Inverted		
Wheel	Rim	Tire	Top	Middle	Bottom	Top	Middle	Bottom
22"	Mag / Spoke	Full Poly	15	15.5	16.5	16.5	17.5	18
22"	Mag / Spoke	Pneum / Insert	15.5	16	17	17	17.5	18
24"	Mag / Spoke	Full Poly	16	16.5	17	17	18	18.5
24"	Spoke / Spiner	Clincher	16	16.5	17	17.5	18	18.5
24"	Mag / Spoke / Spiner	Pneum / Insert	16	16.5	17.5	17.5	18	18.5
26"	Spoke / Spiner	Clincher	17	17.5	18	18.5	19	19.5
26"	Spoke / Spiner	Pneumatic	17	18	18.5	18.5	19.5	20

Maximum Dump 2 3/4"

*** Custom Frame Height Through Configurator - REAR 6 1/2 - 13" (Max frame squeeze = 3 1/2") Calculate Accordingly

GP Series Front STF Chart


*All measurements given are an estimate and have been rounded to the nearest .5"

SHORT Frame (9 1/4" Front / 8 1/4" Rear)			3"	4" Poly 4" Lighted	5" Poly 5" x 1.5 5" Soft Roll	6" Poly 6" x 1.5	6" Pneu	8" Poly	8" Pneu	8x2
Fork	Spacer	Fork Hole	15.50	16.00	16.50	17.00	17.00	x	x	x
4 1/4"	STD	Bottom	16.50	17.00	17.00	17.50	18.00	x	x	x
4 1/4"	3/4"	Bottom	17.00	17.50	18.00	18.50	18.50	x	x	x
5 1/4"	STD	Top	x	x	16.50	17.00	17.00	x	x	x
5 1/4"	STD	Bottom	16.50	17.00	17.5	18.00	18.00	19.00	20.5	x
5 1/4"	3/4"	Top	x	x	17.50	18.00	18.00	x	x	x
5 1/4"	3/4"	Bottom	17.50	18.00	18.50	18.50	19.00	19.50	20.00	x
5 1/4"	1 1/2"	Top	x	x	18.00	18.50	18.50	x	x	x
5 1/4"	1 1/2"	Bottom	18.00	18.50	19.00	19.50	19.50	20.50	20.50	x
6"	STD	Top	x	x	17.50	18.00	18.00	18.50	19.00	x
6"	STD	Bottom	17.50	18.00	18.50	19.00	19.00	19.50	20.00	x
6"	3/4"	Top	x	x	18.00	18.50	18.50	19.50	19.50	x
6"	3/4"	Bottom	18.00	18.50	19.00	19.50	19.50	20.50	20.50	x
6"	1 1/2"	Top	x	x	19.00	19.50	19.50	20.00	20.50	x
6"	1 1/2"	Bottom	19.00	19.50	20.00	20.50	20.50	21.00	21.50	x
7"	STD	Top	x	x	18.50	19.00	19.00	19.50	20.00	x
7"	STD	Bottom	18.50	19.00	19.50	20.00	20.00	20.50	21.00	x
7"	3/4"	Top	x	x	19.00	19.50	19.50	20.50	20.50	x
7"	3/4"	Bottom	19.00	19.50	20.00	20.50	20.50	21.50	21.50	x
7"	1 1/2"	Top	x	x	20.00	20.50	20.50	21.00	21.50	x
7"	1 1/2"	Bottom	20.00	20.50	21.00	21.00	21.50	22.00	22.50	x
8x2	STD	Top	x	x	x	x	x	x	x	19.00
8x2	STD	Bottom	x	x	x	x	x	x	x	20.00
8x2	3/4"	Top	x	x	x	x	x	x	x	19.50
8x2	3/4"	Bottom	x	x	x	x	x	x	x	20.50
8x2	1 1/2"	Top	x	x	x	x	x	x	x	20.50
8x2	1 1/2"	Bottom	x	x	x	x	x	x	x	21.00

STANDARD Frame (10" Front / 9" Rear) ADD 3/4" to each measurement above
TALL Frame (11" Front / 10" Rear) ADD 1 3/4" to each measurement above

*** Custom Frame Height through Configurator - Front 9 - 13" - Calculate Accordingly (Minimum Front STF = 15 1/2" / Max = 26")

Frame Dimensions



Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2, Concord, ON L4K 4P8