















A COMPREHENSIVE
GUIDE TO SELECTING
THE RIGHT QUICKIE
POWER WHEELCHAIR

POWER WHEELCHAIR

REFERENCE GUIDE



MODEL CODES	Q700-UP M					Q700 M					Q500 M				Q300 M Mini				Q200 R																																				
1st Digit	MID WHEEL DRIVE					MID WHEEL DRIVE					MID WHEEL DRIVE				MID WHEEL DRIVE				REAR WHEEL DRIVE																																				
N- No Power Seating																																																							
2nd Digit	K0861 K0884					K0861 K0862 K0884					K0848 K0850 K0868 K0856 K0858 K0877 K0861 K0862 K0884					K0848 K0850 K0856 K0858 K0868 K0877				K0848 K0849 K0868 K0869 K0856 K0857 K0877 K0878				K0848 K0856 K0861				K0848 K0856				K0848 K0849 K0856 K0857				K0848 K0856 K0861				K0848 K0856				K0848 K0849 K0856 K0857				K0848 K0856		K0848 K0849 K0856 K0857					
3rd Digit	GRP 3					GRP 4					GRP 3					GRP 3				GRP 3				GRP 3				GRP 3				GRP 3		GRP 3																					
D- Group 4	MPC					MPD					MPC					MPC HD				MPD				NPC SPC				NPC HD				NPD				MPC				SPC HD				SPD				NPC SCC		NCD SCD					
SEAT TYPE	PRO ADVANCED					PRO ADVANCED					PRO					LITE				CAPTAIN				PRO				LITE				CAPTAIN				PRO				LITE				CAPTAIN				REHAB		CAPTAIN					
STD MAX SPEED	6 mph, 8 mph					5 mph, 6 mph, 8 mph					5 mph, 6 mph, 8 mph					5 mph, 6 mph, 8 mph				6 mph, 8 mph				6 mph				6 mph				6 mph				6 mph				6 mph				6 mph		6 mph									
USER WEIGHT CAPACITY	265 lbs					300 lbs 350 lbs HD					300 lbs 350 lbs HD, 450 lbs HD					300 lbs 350 lbs HD				300 lbs				300 lbs				300 lbs				300 lbs				300 lbs				300 lbs				300 lbs		300 lbs									
ELECTRONICS	R-NET					R-NET					R-NET					R-NET				R-NET				VR-2 AND R-NET				VR-2 AND R-NET				VR-2 AND R-NET				VR-2 AND R-NET				VR-2 AND R-NET				VR-2 AND R-NET				VR-2 AND R-NET		VR-2 AND R-NET					
TURNING RADIUS ¹	19.8"					19.8"					19.8"					19.8"				19.8"				18.5"				18.5"				18.5"				17"				17"				17"				18.5"		18.5"					
SEAT WIDTH	16" to 20"					16" to 22"					16" to 26"					15" to 20"				16" to 20"				16" to 22"				15" to 20"				16" to 20"				16" to 22"				15" to 20"				16" to 20"				15" to 21"		16" to 20"					
SEAT DEPTH	17" to 22"					17" to 22"					14" to 24"					16" to 22"				16" to 22"				14" to 22"				16" to 22"				16" to 22"				14" to 22"				16" to 22"				16" to 22"				16" to 20"		16" to 22"					
SEAT-TO-FLOOR HEIGHT	18" to 19"					17" to 19"					16.5" to 19"					16.5" to 19"				17" to 19"				16.5" to 19"				16.5" to 19"				17" to 19"				16.5" to 19"				16.5" to 19"				17.5" to 19"				17" to 19"		17" to 18"					
LIGHT OPTIONS	• LED LIGHTS & INDICATORS • LED FENDER LIGHTS					• LED LIGHTS & INDICATORS • LED FENDER LIGHTS					• LED LIGHTS & INDICATORS • LED FENDER LIGHTS					• LED LIGHTS & INDICATORS • LED FENDER LIGHTS				• LED LIGHTS & INDICATORS • LED FENDER LIGHTS				• LED LIGHTS & INDICATORS • LED FENDER LIGHTS				• LED LIGHTS & INDICATORS • LED FENDER LIGHTS				• LED LIGHTS & INDICATORS • LED FENDER LIGHTS				• LED LIGHTS & INDICATORS • LED FENDER LIGHTS				• LED LIGHTS & INDICATORS • LED FENDER LIGHTS				• LED LIGHTS & INDICATORS • LED FENDER LIGHTS		• LED LIGHTS & INDICATORS • LED FENDER LIGHTS									
TILT	POSTERIOR 25° ANTERIOR 30°					POSTERIOR 50° ANTERIOR 30°					POSTERIOR 50°					POSTERIOR 50°				POSTERIOR 50°				POSTERIOR 50°				POSTERIOR 50°				POSTERIOR 50°				POSTERIOR 50°				POSTERIOR 50°				POSTERIOR 50°				POSTERIOR 22°		POSTERIOR 22°					
RECLINE	180° POWERED RECLINE w/ 7" SHEAR REDUCTION					165° POWERED RECLINE w/ 7" SHEAR REDUCTION					150° POWERED RECLINE w/ 4" SHEAR REDUCTION					X				25° MANUAL ADJUSTABLE RECLINE				150° POWERED RECLINE w/ 4" SHEAR REDUCTION				X				25° MANUAL ADJUSTABLE RECLINE				150° POWERED RECLINE w/ 4" SHEAR REDUCTION				X				25° MANUAL ADJUSTABLE RECLINE				12° ANGLE ADJUSTABLE		25° MANUAL ADJUSTABLE RECLINE					
SEAT ELEVATOR	8"					12", C-Me					12", C-Me					12", C-Me				12"				12", C-Me				12", C-Me				12"				12", C-Me				12", C-Me				12", C-Me				X		X					
LEGREST	• POWER CENTER MNT					• FIXED CENTER MNT • POWER CENTER MNT					• FIXED CENTER MNT • POWER CENTER MNT • SWING-AWAY					• FIXED CENTER MNT • POWER CENTER MNT • ANGLE-ADJUSTABLE CENTER MNT • SWING-AWAY				• FIXED CENTER MNT • POWER CENTER MNT				• FIXED CENTER MNT • POWER CENTER MNT • ANGLE-ADJUSTABLE CENTER MNT • SWING-AWAY				• FIXED CENTER MNT • POWER CENTER MNT				• FIXED CENTER MNT • POWER CENTER MNT • ANGLE-ADJUSTABLE CENTER MNT • SWING-AWAY				• FIXED CENTER MNT • POWER CENTER MNT • ANGLE-ADJUSTABLE CENTER MNT • SWING-AWAY				• FIXED CENTER MNT • POWER CENTER MNT • ANGLE-ADJUSTABLE CENTER MNT • SWING-AWAY				• SWING-AWAY		• SWING-AWAY									
VENT TRAY	✓					✓					✓					✓				X				✓				✓				X				✓				✓				X				X		X					
BATTERIES	GROUP 24					GROUP 24 GROUP 34					GROUP 24 GROUP 34					GROUP 24 GROUP 34				GROUP 24 GROUP 34				GROUP 22 GROUP 34				GROUP 22 GROUP 34				GROUP 22 GROUP 34				50 Ah (GROUP 22)				50 Ah (GROUP 22)				50 Ah (GROUP 22)				GROUP 22		GROUP 22					
DRIVE WHEELS	14"					14"					14"					14"				14"				13"				13"				13"				12", 14"				12", 14"				12", 14"				12", 14"		12", 14"					
OVERALL BASE WIDTH	24.5"					24.5"					24.5"					24.5"				25"				24"				24"				24"				20.5"				20.5"				20.5"				23"		23"					
OVERALL LENGTH	36.5"					36.5"					36.5"					36.5"				36.5"				33"				33"				33"				33"				33"				33"				33"				42"		42"	
OCCUPIED TRANSIT TESTED	✓					✓					✓					✓				✓				✓				✓				✓				✓				✓				✓				✓		✓					

1. This measurement excludes footplates.

SEDEO® PRO ADVANCED™	SEDEO® PRO™	SEDEO® LITE™	CAPTAIN SEAT
			
A highly specialized and technical seating system, the SEDEO PRO ADVANCED allows for precision fitting for complex clinical needs. Biometric Repositioning and a full range of programmable seating positions make SEDEO PRO ADVANCED the best ergonomic seating from QUICKIE.	Truly modular in design, the SEDEO PRO seating system provides a new level of function, comfort and design. For moderate to complex clinical needs, the SEDEO PRO combines a highly adjustable seat frame with an excellent selection of postural support options.	The SEDEO LITE is an adjustable seat frame with a sleek look. Its minimal design alleviates weight capacity concerns when using aftermarket equipment. Traditional back canes are width-adjustable from 15" to 20" and easily fit JAY® backs. An array of options cover mild to moderate clinical needs, and intuitive adjustability meets individual requirements.	Designed for comfort and convenience, the Captain Seat adjusts in seat depth, legrest placement, and armrest width. The convertible seat pan allows alternate cushion selection, and the backrest has both semi-recline and fold-down capability.

CLINICAL NEED	HIGH	MED TO HIGH	MED TO HIGH	LOW
Q700 M COMPATIBILITY	✓	✓	✓	✓
Q500 M COMPATIBILITY	✗	✓	✓	✓
Q300 M COMPATIBILITY	✗	✓	✓	✓
JAY® SYNCRA™ COMPATIBILITY	✓	✓	✗	✗
BIOMETRIC REPOSITIONING	✓	✗	✗	✗
REMOTE SEATING APP	✓	✗	✗	✗
MEMORY SEATING	✓	✗	✗	✗
PRESSURE RELIEF ALERT REMINDER (PRAR)	✓	✗	✗	✗
POWER RECLINE w/7" of SHEAR (165°)	✓	✗	✗	✗
POWER RECLINE w/4" of SHEAR (150°)	✗	✓	✗	✗
MANUAL RECLINE OPTION	✗	✓	✗	✓
POWER 50° CENTER GRAVITY TILT	✓	✓	✓	✓
ANTERIOR 30° TILT	✓	✗	✗	✗
POWER STANDING	✓	✗	✗	✗
POWER LEGS	✓	✓	✓	✓
12" LIFT	✓	✓	✓	✓
SEAT WIDTH	16" to 22"	16" to 26"	15" to 20"	16", 18", 20"
SEAT DEPTH	17" to 22"	14" to 24"	16" to 22"	16" to 22"
BACKREST HEIGHT	21" to 30"	21" to 30"	18" to 30"	25"

ELEVATE AND DRIVE WITH C-ME® AVAILABLE ON SEDEO PRO ADVANCED, SEDEO PRO, AND SEDEO LITE

No more waiting around or driving at reduced speed when elevated. C-Me enables you to elevate up to 12" while traveling up to 3 mph! You can lift or descend at the touch of a button with no need to stop while transitioning. Just rise-up and go at the same time!