








A COMPREHENSIVE
GUIDE TO SELECTING
THE RIGHT QUICKIE
POWER WHEELCHAIR

POWER WHEELCHAIR

REFERENCE GUIDE



MODEL CODES	Q700 M									Q500 M				Q500 H			Q300 M Mini				Q200 R											
1st Digit	MID WHEEL DRIVE									MID WHEEL DRIVE				HYBRID (MID/REAR) DRIVE			MID WHEEL DRIVE				REAR WHEEL DRIVE											
N- No Power Seating																																
S- Single Power																																
M- Multiple Power																																
2nd Digit																																
P- Rehab Seat																																
C-Captain Seat																																
3rd Digit																																
C- Group 3																																
D- Group 4																																
HCPCS CODES	K0861 K0862 K0884			K0848 K0850 K0868 K0856 K0858 K0877 K0861 K0862 K0884			K0848 K0850 K0856 K0858 K0868 K0877			K0848 K0849 K0868 K0869 K0856 K0857 K0877 K0878		K0848 K0856 K0861		K0848 K0856		K0848 K0849 K0856 K0857		K0848 K0850 K0868 K0856 K0858 K0877 K0861 K0862 K0884		K0848 K0856 K0861		K0848 K0856		K0848 K0849 K0856 K0857		K0848 K0856		K0848 K0849 K0856 K0857				
MODELS	GRP 3 GRP 3HD GRP 4			GRP 3 GRP 3HD GRP 4			GRP 3 GRP 3HD GRP 4			GRP 3 GRP 4		GRP 3		GRP 3		GRP 3		GRP 3 GRP 3HD GRP 4			GRP 3		GRP 3		GRP 3		GRP 3					
	MPC	MPC HD	MPD	NPC SPC MPC	NPC HD SPC HD MPC HD	NPD SPD MPD	NPC SPC	NPC HD SPC HD	NPD SPD	NCC SCC	NCD SCD	NPC SPC MPC	NPC SPC	NCC SCC	NPC SPC MPC	NPC HD SPC HD MPC HD	NPD SPD MPD	NPC SPC MPC	NPC SPC	NCC SCC	NPC SPC	NCC SCC	NPC SPC	NCC SCC	NPC SPC	NCC SCC	NPC SPC	NCC SCC				
SEAT TYPE	PRO ADVANCED			PRO			LITE			CAPTAIN		PRO		LITE		CAPTAIN		PRO			PRO		LITE		CAPTAIN		REHAB		CAPTAIN			
STD MAX SPEED	5 mph, 6 mph, 8 mph			5 mph, 6 mph, 8 mph			5 mph, 6 mph, 8 mph			6 mph, 8 mph		6 mph		6 mph		6 mph		5 mph, 6 mph, 8 mph			6 mph		6 mph		6 mph		6 mph		6 mph			
USER WEIGHT CAPACITY	300 lbs 350 lbs HD			300 lbs 350 lbs HD, 450 lbs HD			300 lbs 350 lbs HD			300 lbs		300 lbs		300 lbs		300 lbs 350 lbs HD			300 lbs		300 lbs		300 lbs		300 lbs		300 lbs		300 lbs			
ELECTRONICS	R-NET			R-NET			R-NET			R-NET		VR-2 AND R-NET		VR-2 AND R-NET		VR-2 AND R-NET		VR-2 AND R-NET			VR-2 AND R-NET		VR-2 AND R-NET		VR-2 AND R-NET		VR-2 AND R-NET		VR-2 AND R-NET			
TURNING RADIUS ¹	19.8"			19.8"			19.8"			19.8"		18.5"		18.5"		18.5"		25.5"			17"		17"		17"		18.5"		18.5"			
SEAT WIDTH	16" to 22"			16" to 26"			15" to 20"			16" to 20"		16" to 22"		15" to 20"		16" to 20"		16" to 22"			16" to 22"		15" to 20"		16" to 20"		15" to 21"		16" to 20"			
SEAT DEPTH	17" to 22"			14" to 24"			16" to 22"			16" to 22"		14" to 22"		16" to 22"		16" to 22"		14" to 22"			14" to 22"		16" to 22"		16" to 22"		16" to 20"		16" to 22"			
SEAT-TO-FLOOR HEIGHT	17" to 19"			16.5" to 19"			16.5" to 19"			17" to 19"		16.5" to 19"		16.5" to 19"		17" to 19"		16.5" to 19"			16.5" to 19"		16.5" to 19"		17.5" to 19"		17" to 19"		17" to 18"			
LIGHT OPTIONS	• LED LIGHTS & INDICATORS • LED FENDER LIGHTS			• LED LIGHTS & INDICATORS • LED FENDER LIGHTS			• LED LIGHTS & INDICATORS • LED FENDER LIGHTS			• LED LIGHTS & INDICATORS • LED FENDER LIGHTS		• LED LIGHTS & INDICATORS • LED FENDER LIGHTS		• LED LIGHTS & INDICATORS • LED FENDER LIGHTS		• LED LIGHTS & INDICATORS • LED FENDER LIGHTS			• LED LIGHTS & INDICATORS • LED FENDER LIGHTS		• LED LIGHTS & INDICATORS • LED FENDER LIGHTS		• LED LIGHTS & INDICATORS • LED FENDER LIGHTS		• LED LIGHTS & INDICATORS • LED FENDER LIGHTS		• LED LIGHTS & INDICATORS • LED FENDER LIGHTS					
TILT	POSTERIOR 50°			POSTERIOR 50°			POSTERIOR 50°			POSTERIOR 50°		POSTERIOR 50°		POSTERIOR 50°		POSTERIOR 50°			POSTERIOR 30°, 50°		POSTERIOR 30°, 50°		POSTERIOR 30°, 50°		POSTERIOR 30°, 50°		POSTERIOR 22°		POSTERIOR 22°			
RECLINE	165° POWERED RECLINE w/ 7" SHEAR REDUCTION			150° POWERED RECLINE w/4" SHEAR REDUCTION			X			25° MANUAL ADJUSTABLE RECLINE		150° POWERED RECLINE w/4" SHEAR REDUCTION		X		25° MANUAL ADJUSTABLE RECLINE			150° POWERED RECLINE w/4" SHEAR REDUCTION		150° POWERED RECLINE w/4" SHEAR REDUCTION		X		25° MANUAL ADJUSTABLE RECLINE		12° ANGLE ADJUSTABLE		25° MANUAL ADJUSTABLE RECLINE			
SEAT ELEVATOR	12", C-Me			12", C-Me			12", C-Me			12"		12", C-Me		12", C-Me		12"			12", C-Me		12", C-Me		12", C-Me		X		X		X			
LEGREST	• FIXED CENTER MNT • POWER CENTER MNT			• FIXED CENTER MNT • POWER CENTER MNT • SWING-AWAY			• FIXED CENTER MNT • POWER CENTER MNT • ANGLE-ADJUSTABLE CENTER MNT • SWING-AWAY			• FIXED CENTER MNT • POWER CENTER MNT		• FIXED CENTER MNT • POWER CENTER MNT • SWING-AWAY		• FIXED CENTER MNT • POWER CENTER MNT • ANGLE-ADJUSTABLE CENTER MNT • SWING-AWAY		• FIXED CENTER MNT • POWER CENTER MNT • SWING-AWAY			• SWING-AWAY		• FIXED CENTER MNT • POWER CENTER MNT • ANGLE-ADJUSTABLE CENTER MNT • SWING-AWAY		• FIXED CENTER MNT • POWER CENTER MNT • ANGLE-ADJUSTABLE CENTER MNT • SWING-AWAY		• SWING-AWAY		• SWING-AWAY		• SWING-AWAY			
VENT TRAY	✓			✓			✓			X		✓		✓		X			✓		✓		X		X		X		X			
BATTERIES	GROUP 24 GROUP 34			GROUP 24 GROUP 34			GROUP 24 GROUP 34			GROUP 24 GROUP 34		GROUP 22 GROUP 34		GROUP 22 GROUP 34		GROUP 22 GROUP 34			GROUP 24 GROUP 34			50 Ah (GROUP 22)		50 Ah (GROUP 22)		50 Ah (GROUP 22)		GROUP 22		GROUP 22		
DRIVE WHEELS	14"			14"			14"			14"		13"		13"		13"			13", 14"			12", 14"		12", 14"		12", 14"		12", 14"		12", 14"		
OVERALL BASE WIDTH	24.5"			24.5"			24.5"			25"		24"		24"		24"			24.5"			20.5"		20.5"		20.5"		23"		23"		
OVERALL LENGTH	36.5"			36.5"			36.5"			36.5"		33"		33"		33"			36"			33"		33"		33"		42"		42"		
OCCUPIED TRANSIT TESTED	✓			✓			✓			✓		✓		✓		✓			✓		✓		✓		✓		✓		✓		✓	

1. This measurement excludes footplates.

SEDEO® PRO ADVANCED™	SEDEO® PRO™	SEDEO® LITE™	CAPTAIN SEAT
			
A highly specialized and technical seating system, the SEDEO PRO ADVANCED allows for precision fitting for complex clinical needs. Biometric Repositioning and a full range of programmable seating positions make SEDEO PRO ADVANCED the best ergonomic seating from QUICKIE.	Truly modular in design, the SEDEO PRO seating system provides a new level of function, comfort and design. For moderate to complex clinical needs, the SEDEO PRO combines a highly adjustable seat frame with an excellent selection of postural support options.	The SEDEO LITE is an adjustable seat frame with a sleek look. Its minimal design alleviates weight capacity concerns when using aftermarket equipment. Traditional back canes are width-adjustable from 15" to 20" and easily fit JAY® backs. An array of options cover mild to moderate clinical needs, and intuitive adjustability meets individual requirements.	Designed for comfort and convenience, the Captain Seat adjusts in seat depth, legrest placement, and armrest width. The convertible seat pan allows alternate cushion selection, and the backrest has both semi-recline and fold-down capability.

CLINICAL NEED	HIGH	MED TO HIGH	MED TO HIGH	LOW
Q700 M COMPATIBILITY	✓	✓	✓	✓
Q500 M COMPATIBILITY	✗	✓	✓	✓
Q500 H COMPATIBILITY	✗	✓	✗	✗
Q300 M COMPATIBILITY	✗	✓	✓	✓
JAY® SYNCRA™ COMPATIBILITY	✓	✓	✗	✗
BIOMETRIC REPOSITIONING	✓	✗	✗	✗
REMOTE SEATING APP	✓	✗	✗	✗
MEMORY SEATING	✓	✗	✗	✗
PRESSURE RELIEF ALERT REMINDER (PRAR)	✓	✗	✗	✗
POWER RECLINE w/7" of SHEAR (165°)	✓	✗	✗	✗
POWER RECLINE w/4" of SHEAR (150°)	✗	✓	✗	✗
MANUAL RECLINE OPTION	✗	✓	✗	✓
POWER 50° CENTER GRAVITY TILT	✓	✓	✓	✓
POWER LEGS	✓	✓	✓	✓
12" LIFT	✓	✓	✓	✓
SEAT WIDTH	16" to 22"	16" to 26"	15" to 20"	16", 18", 20"
SEAT DEPTH	17" to 22"	14" to 24"	16" to 22"	16" to 22"
BACKREST HEIGHT	21" to 30"	21" to 30"	18" to 30"	25"

ELEVATE AND DRIVE WITH C-ME® AVAILABLE ON SEDEO PRO ADVANCED, SEDEO PRO, AND SEDEO LITE

No more waiting around or driving at reduced speed when elevated. C-Me enables you to elevate up to 12" while traveling up to 3 mph! You can lift or descend at the touch of a button with no need to stop while transitioning. Just rise-up and go at the same time!