



## JAY® FUSION® with Cryo™ Technology

Introducing CryoFluid, the new revolutionary fluid option available for the JAY Fusion cushion that excels at reducing the risk of pressure injuries. Patent-pending CryoFluid is now being offered alongside industry favorite, JAY Flow™ Fluid, to provide your patients with more options for skin protection.

CryoFluid actively cools a patient's seated skin surface for up to 8 hours\* while evenly distributing pressure, reducing shear, and lowering the risk of moisture.



Heat leaves skin surface,  
cooling the skin



A long lasting effect

Learn more at: [www.SunriseMedical.com/Cryo](http://www.SunriseMedical.com/Cryo)

\* Internal testing data at 25°C. Results may vary.



**SUPERIOR CLINICAL SEATING**

# JAY® FUSION

## CUSHION WITH CRYO TECHNOLOGY

### INNOVATIVE DESIGN TECHNOLOGY

- X-STATIC® stretch outer cover reduces tension and resists bacterial growth
- JAY CryoFluid insert is designed for skin protection, actively lowers the skin temperature and evenly distributes pressure
- Inner cover with AquaGuard® moisture resistant zipper and anti-wicking thread for moisture protection

### SIMPLE TO FIT

- Pelvic Loading Area (PLA) well and insert sized around anthropometric data, specific for each cushion width
- Pre-contoured foam base with high-resiliency posterior pelvic wall design encourages orthopedic alignment
- The CryoFluid PLA insert is divided into 5 chambers, the two chambers under the ischial tuberosities (ITs) include CryoFluid and the other three chambers include industry favorite, JAY Flow™ Fluid
- Also available in a Positioning Option and Reduced Profile Option

### SIMPLE TO ORDER

- Smart part numbers and a wide range of available modifications
- JAY-Your-Way customizations available

#### JAY TECHNOLOGY IN THE FUSION WITH CRYOFLUID



LAYERED  
SUPPORT



CONTOURING



POSITIONING



PRESSURE  
DISTRIBUTION



FOUR WAY  
STRETCH



ACTIVE  
COOLING



### Specifications:

#### HCPCS Codes

E2622, E2623,  
E2624, E2625

#### Width

16" to 22"

#### Depth

16" to 22"

#### Height

3 3/8" to 4"

#### Weight Capacity

300 lbs for 16" to 21"  
500 lbs for 22"

#### Base

Pre-contoured, closed cell  
high resiliency foam

#### Insert

JAY CryoFluid

#### Inner Cover

Moisture-resistant with  
AquaGuard zipper, anti-  
wicking seam thread

#### Outer Cover

X-STATIC silver  
thread stretch



Sunrise Medical (US) LLC, Fresno, CA 93727  
MK-129563 Rev.A ©11.2019

800.333.4000 [www.SunriseMedical.com](http://www.SunriseMedical.com)

