



ORDER FORM


Q200R with Captain's Seat

EIPW35

Customer		Ship To	
Account #		Name	
Date		Address	
PO#		City	
Buyer		Province	Postal Code
Mark For		Phone	

Client Measurements

a) Top of head to seat	in
b) Top of shoulder to seat	in
c) Underarm to seat	in
d) Elbow to seat	in
e) Top of leg to seat	in
f) Shoulder width	in
g) Chest width	in
i) Hip width	in
j) Actual sitting depth	in
m) Foot length	in
n) Knee to foot	in
o) Seat width	in



Notes

Chair Specifications

Weights starting at: (with Batteries)	260 lbs
Power Seating	22° Tilt
Seat Type	Captain Seat w/ Fold Down Back and 25° Angle Adjustment
Seat Size Range	Width 16" to 20" Depth 16" to 22"
Seat-to-Floor Height Range	17" to 18"

Mid-Wheel Drive Base Specifications

Drive Wheel Position	Rear-wheel drive
User Weight Capacity	300 lbs.
Base Size	Width 23", Length 42"
Minimum Turning Radius	18.5" (excludes footplates)
Drive Control	R-net PM 90A / VR-2 90A
Driving Range	Up to 15.5 miles
Drive Motors	2-pole, 6 mph / 10 kph
Drive Wheels	12" w/ Flat Free Inserts 14" w/ Pneumatic Inserts 14" w/ Flat Free Inserts
Front Casters	9" w/ Pneumatic Inserts 9" w/ Flat Free Inserts
Charger	8A Dual Mode (Compact), 10A
Battery Size	Group 22
Color	7 color accents
Ground Clearance	3"

BASE			
OPTION		DESCRIPTION	RETAIL
Standard Duty Rear Wheel Drive Base <i>End Item Number EIPW35 (Captain's Seat)</i>			\$
350PC304		Quickie Q200 R NCC No Power Option, Captain's Seat and Back Included in Base Pricing : Q200 R Base, VR-2 90amp Electronics for No Power Seat Functions, 4-wheel Suspension System, 12" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 10kph / 6mph, Weight Capacity = 300 lbs., STFH = 17".	\$7,579.00
350PC305		Quickie Q200 R SCC Single Power Option, Captain's Seat and Back Included in Base Pricing : Q200 R Base, VR-2 90amp Electronics for No Power Seat Functions, 4-wheel Suspension System, 12" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 10kph / 6mph, Weight Capacity = 300 lbs., STFH = 17".	\$7,579.00
350PC306		Quickie Q200 R NCC No Power Option, Captain's Seat Back and Solid Seat Included in Base Pricing : Q200 R Base, VR-2 90amp Electronics for No Power Seat Functions, 4-wheel Suspension System, 12" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 10kph / 6mph, Weight Capacity = 300 lbs., STFH = 17".	\$7,579.00
350PC307		Quickie Q200 R SCC Single Power Option, Captain's Seat Back and Solid Seat Included in Base Pricing : Q200 R Base, VR-2 90amp Electronics for No Power Seat Functions, 4-wheel Suspension System, 12" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 10kph / 6mph, Weight Capacity = 300 lbs., STFH = 17".	\$7,579.00
User Weight			\$
350WG23		Up to 300 lbs.	NC
Motor Package			\$
350PF107		2-pole 6 mph / 10 kph	NC
Drive Wheels			\$
350DW77		Drive Wheel, 12", Solid, Knobby, Black	NC
350DW79		Drive Wheel, 14", Solid, Knobby, Black	NC
Caster Wheels			\$
350CW121		Caster Wheels 9", Pneumatic, Black	NC
350CW122		Caster Wheels 9", Solid, Black	NC
Color Accents			\$
350EZ20		Apple Green	NC
350EZ22		Black	NC
350EZ21		Blue	NC
350EZ23		Orange	NC
350EZ24		Red	NC
350EZ19		White	NC
350EZ25		Silver	NC

BASE			
OPTION		DESCRIPTION	RETAIL
ISO Transit Tie Downs			\$
350T050		ISO Transit Package	NC
LED Lights			\$
350LI02		Omit LED Lights	NC
350LI08		LED Fender Lights	\$646.00
Seat Heights <i>(Seat height is measured from front of seat pan to floor with 0° frame angle.)</i>			\$
350HT55		Seat Height 17" (N/A with Tilt)	NC
350HT56		Seat Height 18"	NC
Batteries and Chargers			\$
Group 22 Battery			
350AG31		GP22 Flat Top Pair w/Straps	\$678.00
350AG35		GP22/50Ah Flat Top Pair w/Straps <i>(Required for Tilt)</i>	\$662.00
Chargers			
350BJ30		10A Dual Mode (Compact)	\$162.00
350BJ19		8A Dual Mode (Compact)	NC

SEATING

OPTION	DESCRIPTION	RETAIL
Seating Options		\$
Power Seating Options		\$
Seat Tilt		
350PA56	Posterior CG 22° Tilt	\$3,636.00
Pre Set Seat Frame Angle (No Power Seating Options, Increases Seat to Floor Height)		\$
350FA10	0° (Adds 0" to STFH)	NC
350FA3	3° (Adds 1" to STFH)	NC
350FA4	6° (Adds 2" to STFH)	NC
Backrest Frame Type		\$
350BF08	Captain's Seat Frame (Includes back cushion)	NC
Seat Type		\$
350ST80	JAY®	NC
350ST77	Captain's Seat Cushion	NC
350ST75	Omit	NC
Seat Size		\$
Seat Pan Width		
350R53	16" Available Depths	NC
350R55	18" 16" to 18"	NC
350R57	20" 16" to 20"	NC
	18" to 22"	
Seat Pan Depth (Front of Seat Pan to Back Structure)		
350RF51	16"	NC
350RF45	17"	NC
350RF46	18"	NC
350RF47	19"	NC
350RF48	20"	NC
350RF49	21"	NC
350RF50	22"	NC
JAY® Seat Cushion		\$
350CU59	JAY® Ion®	\$505.00
350CU61	JAY® Union®	\$764.00
350CU64	JAY® Fusion®	\$1,045.00
350CU104	JAY® Fusion® with Cryo® Technology	\$1,152.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	DESCRIPTION	RETAIL
Armrest Structure		\$
350EG37	Flip Back <i>(Includes Waterfall Arm Pads) (Height Adjustable 11" to 14")</i>	271.00
Swing-Away Legrest (Pair)		\$
350H177	SEDEO PRO Legrests, Swing Away, 70° <i>(14" to 21" Leg Length)</i>	NC
350H178	SEDEO PRO Legrests, Manual Elevating, Swing Away <i>(Includes Padded Calf Pads) (14".5 to 21.5" Leg Length)</i>	\$380.00
350H179	SEDEO PRO Legrests, Power Elevating, Swing Away <i>(Includes Padded Calf Pads) (15.5" to 22.5" Leg Length)</i>	\$3,435.00
350H180	SEDEO PRO Hanger Receivers Only	NC
Swing Away Legrest Calf Supports		
350DI42	Calf Support, Depth Adjustable, Padded <i>(6.5" to 8" Wide, Depending On Chair Width, x 9.5" Tall)</i>	\$512.00
Footplates		
350F370	Angle Adjustable Flip Up Footplates	\$277.00
Accessories		
350DI44	Heel Loops	NC
350KG4	Knee Pads, Integral Foam <i>(2.5" Wide x 4.5" Tall)</i>	\$282.00
SEDEO Lower Extremity Swing Away Support		
350VD01	Lower Extremity Support, Padded, Swing Away, SEDEO, Right	\$512.00
350VD02	Lower Extremity Support, Padded, Swing Away, SEDEO, Left	\$512.00
Pelvic Positioning Belts		
350A837	Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
350A839	Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"	NC
350A840	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$103.00
Accessories		\$
O2 Tank Holder		
350A718	O2 Tank Holder Mounted to Backrest (D & E Styles)	\$391.00
Backpack		
350A1	Black Utility	\$112.00
Captain's Seat Accessories		
350A608	Cane / Crutch Holder	\$115.00
350A610	Walker Holder	\$191.00
Medial Thigh Support		
3508A80	Medial Thigh Support	\$271.00
3508A81	Pommel Swing-Away Hardware	\$390.00
Misc. Accessories		
3508A52	Ram Mount Phone	\$77.00
3508A53	Ram Mount Tablet	\$308.00
3508A54	Ram Mount Cup Holder	\$77.00
3508A58	Rear View Mirror	\$77.00
3509A10	Utility Hooks	\$58.00

JOYSTICK / ELECTRONICS

OPTION	DESCRIPTION	RETAIL
VR-2 Joystick \$		
350J127	VR-2 Joystick	NC
R-net Joystick (Remote) \$		
350J121	Standard LED Joystick (w/o Screen Display)	NC
350J138	Color Joystick (2.3" Color LCD Screen Display, CJSM1)	NC
350J118	High Resolution Color Joystick, (CJSM2+BT) 3.5" Color LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control	\$1,264.00
350J119	Omit Remote Control	NC
Joystick Mounting \$		
350JM34	Fixed on Armrest	NC
350JM35	Parallel Swing-Away	\$447.00
Joystick Mount (Location) \$		
350HM16	Right Hand Armrest	NC
350HM17	Left Hand Armrest	NC
Display Language Selection R-net \$		
350L1	English	NC
350L3	French	NC
350L7	Spanish	NC
Joystick Handle (Optional) N/A with HD Joystick \$		
350JH13	Ball Handle	\$110.00
350JH15	T-Handle	\$110.00
Electronic Programmings \$		
350PR10	R-net PC Programmer w/ Dongle	\$706.00
R-net Power Modules \$		
350ER12	Power Module, VR-2 90A, Non-Expandable	NC
350ER13	Power Module, R-net, 90A, Expandable	\$100.00
<input type="checkbox"/> 350EX15	Expandable 90A Controller - R-net	NC
<input type="checkbox"/> 350EX16	Harness for 90A Expandable Controller	NC

JOYSTICK / ELECTRONICS

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

Seating Control

\$

Powered Seating Control

350E040	Omit Thru-Drive <i>(Requires Switch Control and Standard LED Joystick Requires Ctrl+5) (N/A w/ VR-2 electronics)</i>	NC
350E037	Thru-drive Control for 1 Actuator	\$1,443.00
350E038	Thru-drive Control for 2+ Actuators	\$2,903.00

Switch Control *(Available with Omit Thru-Drive only) (Requires use of Omni2 or Joystick Jacks)*

350E033	Dual Toggle for 1 Actuator	NC
350E034	2.5" Single Button Switch for 1 Actuator <i>(Button switches are not mounted)</i>	NC
350E035	Switch-It Mini 1 button for 1 Actuator	NC
350E036	Switch-It Mini 2 button for 1 Actuator	NC
350E041	Ctrl+5 Switch Box for 2 Actuators or more	NC

Additional Power Seating Control

350QG11	Ctrl+5 Switch Box with Voice Activation	\$1,078.00
---------	---	------------

Dual Toggle Mounting Location

350TN1	Toggle Mount, Right	NC
350TN2	Toggle Mount, Left	NC

Ctrl+5 Mounting Location

350EY19	Right Armrest	NC
350EY20	Left Armrest	NC

Ctrl+5 Switch Type

350EY22	Push Buttons	NC
350EY23	Toggle Switches	\$356.00

R-net Programming

\$

Assignable Buttons

Assignable Button feature allows for customization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

350PM12	Include Assignable Buttons	NC
---------	----------------------------	----



Q200 R with Captain's Seat

\$

Base Price

Total of additional equipment

Packaging and transport

Total Price

Date

Signature

Month

Day

Year



www.SunriseMedical.ca

800.263.3390

E-mail: cscanada@sunmed.com

Specifications are subject to change without notice.

© 2025 Sunrise Medical Canada Inc.

1000 Creditstone Road, Unit 2

Concord, Ontario L4K 4P8

07.2025 MK-129835 Rev.5