



Quickie 2 October 2026

EI/Q2X

**Mark For:**

Date: _____
Dealer Acct #: _____
Dealer: _____
Dealer Contact: _____
Dealer Address: _____
Dealer City: _____ PV: _____ PC: _____
Dealer Phone: () _____ Fax: () _____
Confirmation Email: _____
Confirm Via: ☐ Fax ☐ Email

Submitting for:☐ Quote ☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
Attention: _____
Address: _____
Address: _____
Ship To City: _____ PV: _____ PC: _____
Ship To Phone: () _____ Fax: () _____

✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

1. FRAME**Base model**

☐	002WM32	Quickie 2	\$3,245.00
☐	002WM109	Quickie 2 - Transit included	\$3,245.00

Transit Option

N/A on Stroller Handles (Removable Push Handles) or Xtender

☐	002TO46	Transit System WC19	NC
--------------------------	---------	---------------------	----

Frame Style

☐	002FS31	Standard Swing Away	NC
☐	002FS30	Hemi Swing Away	NC

Rear Frame Type (300 lbs. Weight Capacity)

☐	002RS1	Rear Side Frame Straight	NC
☐	002RS2	Rear Side Frame Rounded	NC
☐	002RS3	Depth Adjustable Back (DAB) 3.0" of rearward adjustment. N/A w/ swing-away or flip-back armrest.	\$368.00

Heavy Duty Package

☐	002HV1	Heavy Duty (18"-24" Wide only) (350 lbs. weight capacity). Includes Reinforced Back Upholstery & Stabilizer Bar.	\$1,262.00
--------------------------	--------	--	------------

2. SIZE / SEATING**Frame Width**

☐	002W3	12"	NC
☐	002W4	13"	NC
☐	002W5	14"	NC
☐	002W6	15"	NC
☐	002W7	16"	NC
☐	002W8	17"	NC
☐	002W9	18"	NC
☐	002W10	19"	NC
☐	002W11	20"	\$308.00
☐	002W12	21" (Only Avail. w/ Heavy Duty)	\$308.00
☐	002W13	22" (Only Avail. w/ Heavy Duty)	\$308.00
☐	002W57	B4ME - 23" Width (Heavy Duty)	\$308.00
☐	002W58	B4ME - 24" Width (Heavy Duty)	\$308.00

Frame Depth

Seat Depth can be ordered up to 3" smaller for growth. Std frame available with 13"-22" frame Depth & HD frame available with 16"-22" frame Depth. Hemi available with 16"-18" & 20"-22" frame depth. Please refer to matrix in page 8 of Appendix for frame Depth Vs seat Depth matrix. Growth is achieved to the front of the depth.

☐	002FD4	13"	NC
☐	002FD5	14"	NC
☐	002FD6	15"	NC
☐	002FD7	16"	NC
☐	002FD8	17"	NC
☐	002FD9	18"	NC
☐	002FD10	19" Steel reinforced	NC
☐	002FD36	20" Steel reinforced	NC
☐	002FD51	B4ME - 21" Frame Depth	\$308.00
☐	002FD52	B4ME - 22" Frame Depth	\$308.00

Seat Depth

12" seat depth not available with Hemi frame

☐	002SD5	12"	NC
☐	002SD6	13"	NC
☐	002SD7	14"	NC
☐	002SD8	15"	NC
☐	002SD9	16"	NC
☐	002SD10	17"	NC
☐	002SD11	18"	NC
☐	002SD72	19"	\$308.00
☐	002SD13	20"	\$308.00
☐	002SD105	B4ME - 21" Depth	\$308.00
☐	002SD106	B4ME - 22" Depth	\$308.00
☐	002SD107	B4ME - 21" Depth (Heavy Duty)	\$308.00
☐	002SD108	B4ME - 22" Depth (Heavy Duty)	\$308.00

✓	Option #	Definition	Price	✓	Option #	Definition	Price
2. SIZE / SEATING (CONT)							

Seat Sling			
☐	002SL1	Standard Seat Sling	NC
☐	002SL6	Omit Seat Sling	NC

JAY Cushions (Optional)			
To order any JAY cushions on-chair, please see the JAY price list.			
Cushion name or part number must be filled out with desired width and depth			
JAY Bundled Cushion			
☐	002JC20	Cushion Name or Part #	Price List
		Width:	Depth:

3. CASTERS / SEAT HEIGHT


Front Seat Height							
Please see matrix on Page 9 for additional information. 13.5"-14" Hemi only 21.5"-22" not available in Hemi Frame							
☐	002SH40	13.5"	NC	☐	002SH18	18.0"	NC
☐	002SH41	14.0"	NC	☐	002SH20	18.5"	NC
☐	002SH42	14.5"	NC	☐	002SH21	19.0"	NC
☐	002SH12	15.0"	NC	☐	002SH23	19.5"	NC
☐	002SH13	15.5"	NC	☐	002SH24	20.0"	NC
☐	002SH14	16.0"	NC	☐	002SH26	20.5"	NC
☐	002SH15	16.5"	NC	☐	002SH27	21.0"	NC
☐	002SH16	17.0"	NC	☐	002SH29	21.5"	NC
☐	002SH17	17.5"	NC	☐	002SH36	22.0"	NC

Caster Wheels			
See matrix on page 9.			
☐	002CW1	3" Micro	NC
☐	002CW3	4" Low Profile Polyurethane N/A w/ HD package	NC
☐	002CW31	4" Micro Lighted	\$110.00
☐	002CW72	4" x 1" Aluminum Polyurethane	\$73.00
☐	002CW77	4" x 1.25" Semi Pneumatic	NC
☐	002CW64	4" x 1.5" Aluminum Soft Roll	\$110.00
☐	002CW5	5" Low Profile Polyurethane	NC
☐	002CW51	5" Micro Lighted	\$110.00
☐	002CW6	5" Soft Roll (Grey) N/A w/ Frog Legs	\$98.00
☐	002CW65	5" x 1" Aluminum Polyurethane	\$126.00
☐	002CW78	5" x 1.5" Semi-Pneumatic	\$110.00
☐	002CW66	5" x 1.5" Aluminum Soft Roll	NC
☐	002CW34	6" Polyurethane	NC
☐	002CW33	6" Pneumatic	\$110.00
☐	002CW79	6" x 1.5" Semi Pneumatic	\$110.00
☐	002CW68	6" x 1.5" Aluminum Soft Roll	NC
☐	002CW16	7" Polyurethane	NC
☐	002CW32	8" Polyurethane	NC
☐	002CW21	8" Pneumatic	\$146.00
☐	002CW39	8" x 1.5" Semi Pneumatic	\$146.00
☐	002CW36	8" x 2" Pneumatic	\$247.00
☐	002CW37	8" x 2" Pneumatic w/Airless Insert	\$393.00

Caster Housing			
Integral bearing forks only available with Integral caster housing			
☐	002CP2	Standard	NC
☐	002CP14	Integral Caster Housing N/A with caster Pin locks	\$110.00

Caster Fork			
Please see matrix on Page 9 for additional information. Frog Legs not available with Integral caster housing			
☐	002CF46	Multi Position Fork-Short	NC
☐	002CF47	Multi Position Fork- Tall	NC
☐	002CF62	3" Fork	NC
☐	002CF12	5 1/4" Fork	NC
☐	002CF17	6 " Fork	NC
☐	002CF23	7" Fork	NC
☐	002CF29	8" x 2" Fork	NC
☐	002CF65	Multiple Position Integral caster housing only	NC
☐	002CF69	Multiple Position Single Sided Fork Integral caster housing only	\$176.00
☐	002CF53	Frog Legs Suspension (0-100 lbs)	\$574.00
☐	002CF54	Frog Legs Suspension (101-190 lbs)	\$574.00
☐	002CF55	Frog Legs Suspension (191-265 lbs)	\$574.00

Caster Options			
☐	002OC1	Caster Pin Locks (N/A w/8"x2") NA w/ Xtender, Frog Legs or Transit	\$209.00



✓	Option #	Definition	Price	✓	Option #	Definition	Price
4. REAR WHEELS / SEAT HEIGHT							

Rear Seat Height							
☐	002RH11	13.0"	NC	☐	002RH26	17.5"	NC
☐	002RH12	13.5"	NC	☐	002RH20	18.0"	NC
☐	002RH13	14.0"	NC	☐	002RH27	18.5"	NC
☐	002RH14	14.5"	NC	☐	002RH21	19.0"	NC
☐	002RH15	15.0"	NC	☐	002RH37	19.5"	NC
☐	002RH16	15.5"	NC	☐	002RH22	20.0"	NC
☐	002RH17	16.0"	NC	☐	002RH38	20.5"	NC
☐	002RH18	16.5"	NC	☐	002RH23	21.0"	NC
☐	002RH19	17.0"	NC	☐	002RH39	21.5"	NC
				☐	002RH40	22.0"	NC

Axle Plate			
☐	002AL102	Standard Adjustable (All in-one)	NC
☐	002AL6	Amputee Adjustable N/A w/ Padded Armrest	\$146.00

Axle			
☐	002AX1	Quick Release Steel	NC
☐	002AX11	Quad Quick Release 5 Spoke Mag & Mountain Wheel only	\$182.00

One Arm Drive			
N/A w/ Xtender, 12", 23" & 24" wide, amputee axle plates, and small tray table. N/A w/ combination of Padded SA arms and Anti tips. Only available w/ 20", 22", 24" mag wheels. Max weight capacity of 300 lbs. Only available in 6 spoke Mags			

☐	002AL16	One Arm Drive Adj. Axle Plate	\$1,382.00
	☐	002OA1 Small Handrim (N/A w/ 20" wheels)	NC
		☐ RH position ☐ LH position	
	☐	002OA2 Standard Handrim	NC
		☐ RH position ☐ LH position	
Select Wheel Size:	☐ 20"	☐ 22"	☐ 24"

Xtender Power Assist			
RSH of 14.5"-20" for 22" rear wheels and RSH 15.5"-21" for 24" wheels. Front seat height must be greater than rear. N/A w/ Transit Option Max weight capacity of 265 lbs. Available 12"-22" widths.			

Xtender Power Assist			
☐	002XT6	22" Xtender Power Assist JWX-2 NICKEL	\$9,996.00
☐	002XT7	22" Xtender Power Assist JWX-2 LITHIUM	\$11,604.00
☐	002XT4	24" Xtender Power Assist JWX-2 NICKEL	\$9,996.00
☐	002XT5	24" Xtender Power Assist JWX-2 LITHIUM	\$11,604.00

Please select one of the following handrim options

☐	002XH5	Plastic Coated	NC
☐	002XH6	Stainless Steel	NC

Xtender - Other Options			
☐	002XA10	Spare NiMH Power Pack - NICKEL	\$1,925.00
☐	002XA11	Spare Li-ion Power Pack - LITHIUM	\$2,479.00

Xtender Wheel Options			
☐	002XA8	SPARE 22" Lite Spoke Pneu Alum	NC
☐	002XA9	SPARE 22" Lite Spoke Pneu Plastic	\$144.00
☐	002XA2	SPARE 24" Lite Spoke Pneu Alum	NC
☐	002XA4	SPARE 24" Lite Spoke Pneu Plastic	\$144.00

Rear Wheel Type			
Not applicable w/ Xtender Power Assist Option. Ultra Lightweight Wheels only available w/ Alum Anodized Handrim. Handrims cannot be omitted on Lite Spoke Wheels.* Omit Rear wheel will default chair configuration set up to a 24" wheel			
☐	002RW103	20" Mag Wheel 5 Spoke	NC
☐	002RW104	22" Mag Wheel 5 Spoke	NC
☐	002RW102	24" Mag Wheel 5 Spoke	NC
☐	002RW158	20" Lite Spoke Wheel	\$320.00
☐	002RW159	22" Lite Spoke Wheel	\$320.00
☐	002RW92	24" Lite Spoke Wheel	NC
☐	002RW93	25" Lite Spoke Wheel	NC
☐	002RW160	26" Lite Spoke Wheel	\$126.00
☐	002RW77	22" Ultra Lightweight Spoke Wheel	\$357.00
☐	002RW78	24" Ultra Lightweight Spoke Wheel	\$357.00
☐	002RW28	22" Spinergy SPOX Black	\$864.00
☐	002RW95	24" Spinergy SPOX Black	\$864.00
☐	002RW96	25" Spinergy SPOX Black	\$864.00
☐	002RW97	26" Spinergy SPOX Black	\$864.00
☐	002RW32	22" Spinergy LX Black	\$1,111.00
☐	002RW98	24" Spinergy LX Black	\$1,111.00
☐	002RW99	25" Spinergy LX Black	\$1,111.00
☐	002RW40	26" Spinergy LX Black	\$1,111.00
☐	002RW50	24" Mountain	\$393.00
☐	002RW25	Omit Rear Wheel and Axle *	NC
☐	002RW49	Omit Rear Wheel (Use w/24" Xtender)	NC
☐	002RW111	Omit Rear Wheel (Use w/22" Xtender)	NC

Rear Wheel Tire			
☐	002RT1	Pneumatic	\$0.00
☐	002RT2	Pneumatic Airless Insert N/A with Xtender.	\$247.00
☐	002RT24	Pneumatic Airless Insert - Infinergy® Lightweight - E-TPU material	\$247.00
☐	002RT3	Full Poly	\$0.00
☐	002RT38	SHOX Performance Poly	\$146.00
☐	002RT25	Primo Orion Poly	\$271.00
☐	002RT39	High-Pressure Clincher Black	\$105.00
☐	002RT27	Knobby V-TRAK	\$121.00
☐	002RT21	Schwalbe Marathon Plus Evolution	\$157.00
☐	002RT12	Mountain Tire <i>Only available on 24" Mountain Rear Wheel</i>	NC

✓	Option #	Definition	Price	✓	Option #	Definition	Price
4. REAR WHEELS / SEAT HEIGHT (CONT)							


Rear Seat Height Availability Matrices

Rear tire, side frame, and axle plate selections dictate seat height availability. *Not all combinations are possible within acceptable seat height ranges

Rear Tire	Min RSH	Max RSH
20" Pneumatic, Poly	13"	19.0"
22" Pneumatic, Poly	14.0"	20.0"
22" HP Clincher Blk	13.5"	19.5"
22" Schwalbe	13.5"	19.5"
24" Pneumatic, Poly	15"	21.0"
24" HP Clincher Blk	14.5"	20.5"
24" Knobby	15"	21.0"
24" Schwalbe	14.5"	20.5"
24" SHOX Poly	14.5"	20.5"
24" Mountain	15.0"	21.0"
25" HP Clincher Blk	15.0"	21.0"
25" Schwalbe	15.0"	21.0"
25" SHOX Poly	14.5"	20.5"
26" Pneumatic	16.0"	22.0"
26" HP Clincher Blk	15.5"	21.5"
26" Schwalbe	15.5"	21.5"

Axle Plate	Min RSH	Max RSH
Std Adj (All in one)	13"	22.0"
Amputee	14"	22.0"

Rear Side Frame	Min RSH	Max RSH
Curved	13"	21.0"
Straight	13	22.0"
DAB	13	22.0"



Wheel	5 Spoke	Mag	Lite	Ultra	Spinergy
Tire	Mag		Spoke	Lt wt	
Pneumatic	20,22,24	20,22,24	20,22,24,26		24,26
Pneumatic w/Airless Insert	20,22,24	20,22,24	20,22,24,26		24,26
Full Poly Tire	20,22,24	20,22,24	20,22,24		
High Pressure Clinch Blk			22,24,26	22,24	22,24,26
Knobby V-TRAK			24		24
Schwalbe Marathon			24,25,26	22,24	24,25,26
SHOX Perform. Poly	20,22,24		24,25		24,25
20"= 451 mm 22" = 501 mm 24" = 540 mm 25" = 559 26" = 590					

Handrim

Natural Fits available in 22", 24", and 25". Handrim selection N/A w/ Xtender.

Ultra lightweight wheels only available w/ alum anodized handrim. Plastic coated limited to certain rims and sizes. Natural Fit w/ 5 spoke MAG only available in 24"

☐	002HR1	Aluminum Anodized	NC
☐	002HR4	Plastic Coated - N/A on every rim	\$146.00
☐	002HR22	Natural Fit - Standard Grip	\$374.00
☐	002HR23	Natural Fit - Super Grip	\$434.00
☐	002HR28	Q-Grip Neoprene Coated	\$471.00
☐	002HR39	Ellipse 3R	\$550.00
☐	002HR40	Carbolife - Gekko	\$575.00
☐	002HR41	Carbolife - Curve L Tetra Grip	\$625.00
☐	002HR7	Omit Handrim Option may not be available on all wheels	NC
☐	002HR8	8 Oblique N/A w/ Xtender. Only available with 24" or 26" 6 spokeMag or Lite Spoke rims.	NC

Camber

☐	002CI1	0°	NC
☐	00217CI	2° (default)	NC
☐	002CI3	4°	NC

Center of Gravity Preset

Amputee Axle Plate will come preset at 1.0" Center of Gravity (CG), range of adjustment is -3.0" to -1.75" and 1.0" to 3.5" CG.

☐	002CC2	1.0" (default)	NC
☐	002CC3	1.5"	NC
☐	002CC4	2.0"	NC
☐	002CC5	2.5"	NC
☐	002CC6	3.0"	NC
☐	002CC7	3.5"	NC

One Arm Drive CG

☐	002CC30	1AD Forward position	NC
☐	002CC31	1AD Reward position	NC
☐	002CC32	1AD Center position	NC

5. WHEEL LOCKS

Wheel Lock

Ergo & compact N/A w/ 13"-14" frame depth or amputee axle plate. Additional No Gos may exist depending on wheel size and position

☐	002WL46	High Mount Push	NC
☐	002WL50	High Mount Pull	NC
☐	002WL75	Ergo Scissor	\$157.00
☐	002WL64	Compact	\$162.00
☐	002WL24	Unilateral Wheel Lock Right	\$320.00
☐	002WL23	Unilateral Wheel Lock Left	\$320.00
☐	002WL10	Omit	NC

Wheel Lock Option

☐	002WL12	Do Not Mount	NC
--------------------------	---------	--------------	----

Extension Handle

N/A w/ ergo scissor or compact wheel locks

☐	002EH7	Extension Handle (Pair)	\$73.00
--------------------------	--------	-------------------------	---------

Grade Aid

Only available w/ high mount push or pull to lock wheel locks on pneumatic or pneumatic w/ airless tire options

☐	002GA1	Grade Aid	\$157.00
--------------------------	--------	-----------	----------

Anti Roll-Back

☐	002WL100	Automatic Braking system	\$2,865.00
--------------------------	----------	--------------------------	------------

PRATIKO Automatic Braking System. Installed by Sunrise Medical. Wheel Lock selection required - If no selection, Push to Lock will be provided. Weight capacity of 250 lbs (113 kg). Not available with transit, drum brake, seat pan, OneArmDrive, Reclining Back, Heavy Duty, or seat widths wider than 21". Available with the following rear wheels: 24" Lite Spoke and XCore Mags. Folding width maybe be affected.

✓	Option #	Definition	Price	✓	Option #	Definition	Price
6. HANGER & FOOTREST							

Hangers			
<i>Front mount hangers only available w/ aluminum angle adj footplates. ALR/ELR Hangers N/A w/ Push To Lock wheel locks.</i>			
☐	002H145	60° Swing In/Out Hanger	NC
☐	002H139	70° Swing In/Out Hanger	NC
☐	002H141	80° Swing In/Out Hanger	NC
☐	002H154	60° Swing In/Out Hanger Hemi	NC
☐	002H155	70° Swing In/Out Hanger Hemi	NC
☐	002H140	70° Swing In/Out Front Mount	NC
☐	002H142	80° Swing In/Out Front Mount	NC
☐	002H158	90° Swing In/Out Front Mount	NC
☐	002H159	Articulating Hanger - Ext Mount N/A w/ Platform Flip Up and Alum Angle Adj Front Mount footplates, 14"-22" width	\$526.00
☐	002H160	Elevating Hanger - Ext Mount N/A w/ Platform Flip Up and Alum Angle Adj Front Mount footplates, 14"-22" width	\$339.00
☐	002H161	Elevating Hanger - Front Mount Only available w/ front mount alum angle adj footplates and 14"-22" width	\$339.00
☐	002H32	Omit Hanger & Footplate	NC

Articulating / Elevating Hanger Cover			
☐	002DI47	Neoprene Cover	\$325.00
☐	002DI48	Gel Pad Cover	\$380.00

Footplates			
<i>Heel loop or leg strap not included, See Footplate Options Section</i>			
☐	002F1	Composite (13"-22"w)	NC
☐	002F271	Composite Angle Adj-(16"-22"w)	\$176.00
☐	002F3	Alum Angle Adj - Ext Mt (12"-22"w)	\$176.00
☐	002F8	Alum Angle Adj - Front Mt(12"-22"w)	\$176.00
☐	002F5	Alum Footplate - Fixed (14"-22"w) Available w/ Ext Mount ELR or ALR	NC
☐	002F48	Platform Flip Up Angle Adj (12"-22"w) N/A w/ Hemi Hangers,ELR/ALR,Front Mnt	\$230.00
☐	002F7	Omit Footplate & Ext Tube	NC

Footplate Options (Ext Tube & Front Mount)			
☐	002HL1	Heel Loop	NC
☐	002HL5	Leg Strap	NC

Footrest Extension Tubes			
<i>Required w/ extension tube mounted hangers. See leg length matrix</i>			
☐	002E61	Short Hemi	NC
☐	002E1	Short	NC
☐	002E2	Medium	NC
☐	002E5	Long	NC
☐	002E16	High Mount N/A w/ one-arm drive option or ELR/ALR or 90° front mount. Only available w/angle adj footplates	\$97.00
☐	002E23	2" - Front Mount Only N/A w/ 6", 7" or 8" casters with 90° hangers, or ELR's	\$53.00
☐	002E24	4" - Front Mount Only N/A w/6", 7" or 8" casters w/ 90° hangers or ELR/ALR	\$53.00

Lower Leg Length Setting (Extension Mount Only)							
Recommend 2" less than FSTFH. (*)denotes frame, 9" & 9.5" Hemi Only							
☐	002LO01	9.0"(*)	NC	☐	002LO13	15.0"	NC
☐	002LO02	9.5"(*)	NC	☐	002LO14	15.5"	NC
☐	002LO03	10.0"	NC	☐	002LO15	16.0"	NC
☐	002LO04	10.5"	NC	☐	002LO16	16.5"	NC
☐	002LO05	11.0"	NC	☐	002LO17	17.0"	NC
☐	002LO06	11.5"	NC	☐	002LO18	17.5"	NC
☐	002LO07	12.0"	NC	☐	002LO19	18.0"	NC
☐	002LO08	12.5"	NC	☐	002LO20	18.5"	NC
☐	002LO09	13.0"	NC	☐	002LO21	19.0"	NC
☐	002LO10	13.5"	NC	☐	002LO22	19.5"	NC
☐	002LO11	14.0"	NC	☐	002LO23	20.0"	NC
☐	002LO12	14.5"	NC	☐	002LO24	20.5"Std. Only	NC

Z-finity Hangers/Footrest System			
<i>(Maximum weight capacity for Z-Finity footrests is 165 lbs)</i>			
<i>Some interference may exist in certain configurations. N/A with Heavy Duty</i>			
☐	002H162	Z-finity Hanger	NC

Z-finity Footrest Mount			
☐	002DI22	Direct Hanger Mt/Fixed Frame - High	NC
☐	002DI23	Direct Hanger Mt/Fixed Frame - Low	NC
☐	002DI24	Extension Tube Mount	NC

Z-finity Footrest Pivot Arm			
☐	002DI20	Short 2"	NC
☐	002DI21	Long 3 1/4"	NC

Z-finity Footplate			
<i>Some interference may exist in certain configurations.</i>			
☐	002F320	Multi Angle Adjustable (Pair) (not avail. w/8"-9" widths)	\$247.00
☐	002F321	Angle Adjustable Platform Footboard	\$247.00
☐	002F6	Omit Footplate	NC

Z-finity ExtensionTube Length (When Choosing Ext. Tube Mount)			
☐	002E62	Hemi Short	NC
☐	002E1	Short	NC
☐	002E2	Medium	NC
☐	002E5	Long	NC

✓	Option #	Definition	Price	✓	Option #	Definition	Price
7. BACKREST / SEATING							

Back Type							
☐	002BT10	Standard (see section A)	NC	☐	002BT22	Angle Adjustable Back (5"-7" Pivot) (see section B)	\$357.00
☐	002BT20	Angle Adjustable Back (3-5" Pivot) (see section B)	\$357.00	☐	002BT6	Depth Adjustable (Required with DAB) (see section C)	NC

SECTION A

Back Height Setting							
☐	002UX4	9"	NC	☐	002UX16	15"	NC
☐	002UX6	10"	NC	☐	002UX18	16"	NC
☐	002UX8	11"	NC	☐	002UX20	17"	NC
☐	002UX10	12"	NC	☐	002UX22	18"	NC
☐	002UX12	13"	NC	☐	002UX24	19"	NC
☐	002UX14	14"	NC				

Back Height			
☐	002B1	Low (9"-12" Std Back height range)	NC
☐	002B5	Medium (12"-15" Std Back height range)	NC
☐	002B9	Tall (16"-19" Std Back height range)	NC

Back Post Options			
<i>Omit push handles N/A w/ stabilizer bar. Folding HA push handle N/A w/8 deg or Low back height.</i>			
Standard Back type options:			
☐	002BP1	Push Handle Straight	NC
☐	002BP24	Push Handle 8° Med or Tall only. N/A w/ Angle Adjustable Backrest	NC
☐	002BP2	Omit Push Handle Straight (N/A w/HD)	NC
☐	002B76	Omit Push Handle 8° (N/A w/HD) Med or Tall only	NC
☐	002PH17	Folding Push Handle with Ht. Adj. Back N/A w/ Angle Adj Back & Heavy D.	\$310.00

SECTION B


Angle Adjustable Backrest Pivot Height			
<i>Standard push handles only. Some interference w/ Padded SA arms. N/A w/ flip-back or flip-back height adjustable armrests (classic).</i>			
☐	002AH1	Angle Adjustable Pivot 3"	NC
☐	002AH2	Angle Adjustable Pivot 4"	NC
☐	002AH3	Angle Adjustable Pivot 5"	NC
☐	002AH4	Angle Adjustable Pivot 6"	NC
☐	002AH5	Angle Adjustable Pivot 7"	NC

Preset Angle Adjustable Backrest (-15° to +3°) (Optional)			
☐	002BD5	Back Angle - 15°	NC
☐	002BD4	Back Angle - 12°	NC
☐	002BD3	Back Angle - 6°	NC
☐	002BD2	Back Angle - 0°	NC
☐	002BD1	Back Angle +3°	NC

Angle Adj Back type options:			
☐	002BP5	Angle Adj. Straight Push Handles	NC
☐	002BP6	Angle Adj. Omit Push Handles, N/A w/HD	NC

Circle the height setting in the below chart:

Angle Adjustable Backrest Chart			
Back Tube Height	Low	Medium	Tall
Pivot Regular set at 3"	12" 13" 14"	15" 16" 17"	17" 18" 19"
Pivot Regular set at 4"	13" 14" 15"	16" 17" 18"	18" 19" 20"
Pivot Regular set at 5"	14" 15" 16"	17" 18" 19"	19" 20" 21"
Pivot Tall set at 5"	14" 15" 16"	17" 18" 19"	19" 20" 21"
Pivot Tall set at 6"	15" 16" 17"	18" 19" 20"	N/A
Pivot Tall set at 7"	16" 17" 18"	19" 20" 21"	N/A
Adjustable Pivot Point Measurement from Seat Sling			
If back height option is not chosen, middle option will be present			



SECTION C

DAB Back Post Options			
<i>Diagram of DAB in Appendix on page 9. SA Stroller Handle height range is 21.25"-28.5"</i>			
☐	002BP33	DAB - Swing Away Stroller Handles	\$193.00
☐	002BP27	DAB - Post Standard (19")	NC
☐	002BP28	DAB - Post Medium (22")	NC
☐	002BP29	DAB - Post Tall (25")	NC

DAB Backrest Options			
<i>DAB Solid Back NA w/23" & 24" frame widths</i>			
☐	002BU28	Omit DAB Solid Back and Hardware	NC
☐	002BU29	Omit DAB Solid Back Include Hardware	\$162.00
☐	002BU23	Depth Adjustable Solid Back	\$188.00

DAB Solid Backrest Height			
☐	002B66	Depth Adjustable 12"	NC
☐	002B67	Depth Adjustable 14"	NC
☐	002B68	Depth Adjustable 16"	NC
☐	002B69	Depth Adjustable 18"	NC
☐	002B70	Depth Adjustable 20"	NC

Back Type - Build 4 Me Options			
<i>Must choose DAB Back type with B4ME options</i>			
☐	002BT127	B4ME - Dynamic Back	\$1,558.00
		☐ Fixed Back post	NC
		☐ Std (19") ☐ Med (22") ☐ Tall (25")	
☐	002BT128	Height Adjustable Stroller handles	NC
		☐ Std (18"-23") ☐ Tall (19"-26")	
		B4ME - Reclining Backrest	
		☐ Med (19") ☐ Tall (22")	\$1,458.00

✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

7. BACKREST / SEATING (cont)

BACKREST UPHOLSTERY AND ACCESSORIES

Backrest Upholstery (Standard or Angle Adjustable)

Std EXO upholstery and Tension adj binding is always black. EXO Tension Adj Uph N/A with Heavy Duty. Max back height 20" for EXO.

Not available with DAB back type.

☐	002BU1	Standard Upholstery	\$0.00
☐	002BU88	Standard Upholstery EXO	\$75.00
Color Binding: ☐ Silver ☐ Orange ☐ Red ☐ Blue			
☐	002BU70	Backrest Tension Adjust Ballistic *Reference 186BU70 for billing	\$190.00
☐	002BU89	Backrest Tension Adj EXO *Reference XXXBU89 for billing	\$449.00
Color Binding: ☐ Silver ☐ Orange ☐ Red ☐ Blue			
☐	002BU71	Backrest Tension Adjust 3DX Vented *Reference 186BU71 for billing	\$471.00
☐	002BU82	Jay GO Backrest N/A on Heavy Duty	\$598.00
☐	002BU14	Omit Back Upholstery	NC

Standard Backrest Option

☐	002OD4	Reinforced Upholstery Only avail w/ Std Back Upholstery. Std w/HD package	\$53.00
--------------------------	--------	---	---------

Stabilizer Bar

N/A with Folding Push Handle, or with Ht. Adj. Back or DAB Stroller handles

☐	002SZ1	Stabilizer Bar Standard with HeavyDuty Package and Recline	\$168.00
--------------------------	--------	--	----------

Military Service Patches (Optional)

N/A with standard Upholstery

☐	002VU2	Army Patch	NC
☐	002VU4	Air Force Patch	NC
☐	002VU12	Coast Guard Patch	NC
☐	002VU1	Marines Patch	NC
☐	002VU3	Navy Patch	NC
☐	002VU5	Purple Heart Patch	NC
☐	002VU9	Purple Heart with Army Patch	NC
☐	002VU11	Purple Heart with Air Force Patch	NC
☐	002VU13	Purple Heart with Coast Guard Patch	NC
☐	002VU8	Purple Heart with Marines Patch	NC
☐	002VU10	Purple Heart with Navy Patch	NC


JAY Back (Optional)

☐	002JB01	J3 Back (See J3 Order Form)	Price List
--------------------------	---------	-----------------------------	------------

Stroller Handles (Removable Push Handles)

Only available w/ omit push handles. N/A w/ DAB, amputee axle plate, push handles, or when angle adjustable back is selected, flip-back or flip-back height adjustable armrests, and JAY backs.

☐	002SR2	Regular Stroller Handles	\$193.00
☐	002SR3	Tall Stroller Handles	\$193.00




NEW!

Stroller Handle Extension - Removable

Add the part number to your order corresponding to the chair width. Only available for widths 12-20". 7/8" to 1" tubing push cane conversion kit will be required except with DAB.

☐	#842052	Stroller Extension Assembly 12" wide	\$185.00
☐	#842053	Stroller Extension Assembly 13" wide	\$185.00
☐	#842054	Stroller Extension Assembly 14" wide	\$185.00
☐	#842055	Stroller Extension Assembly 15" wide	\$185.00
☐	#842056	Stroller Extension Assembly 16" wide	\$185.00
☐	#842057	Stroller Extension Assembly 17" wide	\$185.00
☐	#842058	Stroller Extension Assembly 18" wide	\$185.00
☐	#842059	Stroller Extension Assembly 19" wide	\$185.00
☐	#842063	Stroller Extension Assembly 20" wide	\$185.00
☐	#125324	Push Handle Conversion Kit 7/8" to 1"	\$146.00



✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

8. ARMREST / SIDE GUARDS

Padded Swing Away

N/A w/ DAB, Amputee axle. Some limitations may exist with integral push handles, some RSTFHs and wheel sizes. **Now available w/Transit.

☐	002AA176	Padded Swing-Away STD (9" to 13")	\$278.00
Setting ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13"			

Flip Back

N/A w/ DAB, Xtender, 12" - 14" frame depth, Classic Style armpad in highest/lowest axle RSH on rounded rear frame, or angle adjustable back. 4 Point (FS) J3 mounting bracket will cause interference with Flip Back Height Adj armrest. Flip back N/A with amputee axle plate for seat heights requiring the highest axle plate position

☐	002AA6	Flip Back Height Adj (Classic 8.5"-13.5")	\$381.00
☐	002AA5	Flip Back (incl. Desk Length Pad)	\$0.00
☐	002AA65	Flip Back Height Adj (LITE 8.5"-13.5")	\$381.00

Flip Up (For use with Depth Adjustable Back)

Requires an Armpad or Foam from the Arm Pad section below

☐	002AA7	Length Adjustable Locking Flip-up	\$339.00
☐	002AA1	Adj. Locking Flip-up (Std. w/ foam)	\$339.00

Height Adjustable

☐	002AA11	Single Post Hgt. Adj. Std (8.5"-13.5")	\$339.00
☐	002AA29	Single Post Hgt. Adj. Low (7"-11")	\$339.00

Omit Armrest

☐	002AA14	Omit Armrest	NC
--------------------------	---------	--------------	----

Armpad for Height Adjustable

Ergo armpad only available with Single post Height adj & Flip back armrest

Ergo not available with Classic Flip Back armrest

☐	002AM23	Desk Length - Classic (10")	NC
☐	002AM24	Full Length - Classic (14")	NC
☐	002AM9	Armpad Foam - Short (Flip Up armrests)	NC
☐	002AM10	Armpad Foam - Long (Flip Up armrests)	NC
☐	002AM27	Ergo Pad (9.5"-13.5")	NC

Side Guard

Plastic side guards N/A w/ height adj. arms. Side guard interference w/ Padded SA arms at some heights may exist.

☐	002SG33	Plastic Standard (8.5" High)	\$188.00
☐	002SG29	Plastic Low (6.5" High)	\$188.00

9. KEY OPTIONS & ACCESSORIES

Anti-Tip

☐	002AT1	Rear Anti-Tip (NC w/ Xtender)	\$121.00
--------------------------	--------	-------------------------------	----------

Backpack

☐	002A1	Quickie Backpack Black	\$112.00
☐	002A2	Quickie Kids Backpack Black	\$105.00

Seat Pouch

☐	0021A01	Seat Pouch Black /Charcoal	\$83.00
--------------------------	---------	----------------------------	---------

JAY Positioning

☐	002K1	JAY Lower Extremity S/A support	NC
☐ Left ☐ Right Position			
☐	002K2	Omit Hanger	NC
☐ Left ☐ Right			
☐	002K7	Pad Size	NC
☐ 6" x 8" ☐ 8" x 10"			

Additional Cover & Pack Accessories

☐	0028A02	Gel Cover Hanger/Frame	\$162.00
☐	0028A01	Gel Cover Swing Away Arm	\$162.00

Positioning Belts

☐	002A141	Belt 1.5" Auto Buckle w/Clamps	\$97.00
☐	002A117	Belt 2" Aircraft Buckle w/Clamps	\$139.00
☐	002A118	Belt 2" Aircraft Buckle Pad w/Clamp	\$204.00
☐	002A120	Positioning Belt Slide Clamps	\$61.00

Spoke Guards

Only available w/Mountain, and Lite Spoke wheels,

Not available with 5 spoke MAG wheels or Projection handrims

☐	0024A02	Quickie Clear Spoke Guards	\$157.00
☐	0024A01	Quickie Black Spoke Guards	\$157.00

✓ Option #	Definition	Price	✓ Option #	Definition	Price
9. KEY OPTIONS & ACCESSORIES (cont)					
Tray Tables			Miscellaneous		
☐ 002A30	Small (13"-14"W)	\$381.00	<i>O2 holder N/A with Xtender or 12" - 13" frame width. Crutch holder N/A with EXO Tension Adj upholstery and Omit push handle</i>		
☐ 002A31	Medium (15"-17"W)	\$381.00	☐ 0028A03	Caddy (N/A w/ front mount hangers)	\$89.00
☐ 002A32	Large (18"-20"W)	\$381.00	☐ 002A2100	Crutch Holder	\$115.00
☐ 002A33	X-Large (21"-22"W)	\$381.00	☐ 002LS14	Extra Leg Strap Black	\$37.00
			☐ 002IG2	Impact Guard Neoprene	\$61.00
			☐ 002A73	IV Holder	\$218.00
			☐ 002A145	O2 Holder	\$428.00
			☐ 0028A21	Tool Kit	\$43.00
			☐ 002A999	Quickie Water Bottle & Holder	\$81.00

10. COLOUR

Standard Frame Colors			
☐ 002C211	Celestial Sapphire	NC	
☐ 002C212	Desert Sage	NC	
☐ 002C213	Desert Tan	NC	
☐ 002C214	Gloss Black	NC	
☐ 002C215	Matte Black	NC	
☐ 002C216	Metallic Shadow	NC	
☐ 002C217	Midnight Blue	NC	
☐ 002C218	Mystic Purple	NC	
☐ 002C219	Petal Pink	NC	
☐ 002C220	Pink Power	NC	
☐ 002C221	Red Pepper	NC	
☐ 002C222	Royal Blue	NC	
☐ 002C223	Salute Green	NC	
☐ 002C224	Scarlet noir	NC	
☐ 002C225	Silver Ore	NC	
☐ 002C226	Solar Orange	NC	
☐ 002C227	Starlight White	NC	
☐ 002C228	Storm Gray	NC	
☐ 002C229	Turbo Green	NC	
☐ 002C230	Velocity Yellow	NC	

Legacy - Frame Colors			
☐ 002C141	Black Cherry	NC	
☐ 002C10	Black Opal	NC	
☐ 002C14	Candy Purple	NC	
☐ 002C142	Electric Blue	NC	

Personalized Frame Color			
☐ 002C112	Desert Camouflage	\$368.00	
☐ 002C43	Mossy Oak Camouflage	\$368.00	
☐ 002C42	Stars & Stripes	\$368.00	
☐ 002C108	Zebra (Choose base color) Std w/ Starlight White. NA w/ Gloss Black and Matte Black.	\$368.00	

Anodized Color Package			
<i>Standard Axle plate, Multi Forks and Aluminum Soft Roll casters only. N/A with Frog Legs.</i>			
☐ 002EB07	Black - Caster only	NC	
☐ 002EB08	Blue	\$212.00	
☐ 002EB13	Gold	\$212.00	
☐ 002EB09	Lime Green	\$212.00	
☐ 002EB10	Purple	\$212.00	
☐ 002EB11	Red	\$212.00	
☐ 002EB12	Ti Grey	\$212.00	

Additional Options Through Built for Me
 To inquire about any unique personalized Quickie products, contact: Phone: (800) 263 - 3390, Email cscanda@sunmed.com



Commonly Requested Options Include:


- ✓ Personalized Anodizing
- ✓ Xtender Battery on right side
- ✓ Dual Mount Wheel Locks - Attendant Foot Lock
- ✓ Hub Wheel Locks on Mag or Spoke Wheels
- ✓ Colored handrims



11. APPENDIX

Frame Depth matrix	
Frame Depth	Seat Depth
	12" 13" 14" 15" 16" 17" 18" 19" 20" 21" 22"
	13" S S
	14" S S S
	15" S S S S
	16" S S S S S
	Hemi* S S S S
	17" S S S S S
	Hemi* S S S S
	18" S S S S S
	Hemi* S S S S
	19" S S S S S
	20" S S S S S
	Hemi* S S S S S
	21" S S S S S
	Hemi* S S S S S
	22" S S S S S
	Hemi*

"S" = Standard Frame Style
 "H" = Hemi Frame Style
 *Only (6) Frame Depths avail. on Hemi Style Frame



✓	Option #	Definition	Price	✓	Option #	Definition	Price
11. APPENDIX							

Lower Leg Length Matrix

2" & 4" extensions for front mount hangers shipped separate, not mounted. Front seat height must be at least 2" taller than lower leg length on all swing in/out hangers. Front seat height must be at least 2" taller than maximum lower leg length on all front mount hangers. All values shown apply to Standard frame style, subtract 1.5" to get Hemi equivalent. More than 2" will require on field adjustment.

Hanger	Extension	Seat Sling to Top of Footrest in inches		Min Allowable FSH
		Min	Max	
60° Hemi Swing In/Out	Short Hemi	10	13	14
	Short	-	-	-
	Med	12.5	15.5	15.5
	Long	14.5	17.5	17.5
70° Hemi Swing In/Out	Short Hemi	11	14	14
	Short	-	-	-
	Med	13.5	16.5	16
	Long	15.5	18.5	18
60° Swing In/Out	Short	11.5	15	15
	Med	11.5	17	15
	Long	14	18.5	17
70° Swing In/Out	Short	12	16.5	15
	Med	12.5	18.5	15.5
	Long	14.5	20.5	17.5
80° Swing In/Out	Short	13.5	17	15.5
	Med	13.5	19	15.5
	Long	15.5	20.5	17.5
60° Swing In/Out w/Platform Flip Up	Short	11.5	14.5	14
	Long	14	18.5	16
70° Swing In/Out w/Platform Flip Up	Short	13	15.5	15
	Long	15.5	20	17.5
80° Swing In/Out w/Platform Flip Up	Short	13.5	14.5	15.5
	Long	15	19	17
70° Front Mount	No Ext	6	11.5	14
	2"	6	14	16
	4"	6	16	18
80° Front Mount	No Ext	6	12	14
	2"	6	14.5	16.5
	4"	6	16.5	18.5
90° Front Mount	No Ext	5.5	11.5	14
	Short	15	16.5	17
	Med	15	18.5	17
ALR	Long	16.5	20.5	18.5
	Short	12	16.5	15.5
	Med	12	18.5	15.5
ELR Ext. Mount	Long	13.5	20.5	17
	No Ext	5.5	9.5	13

Bckrest Height Matrix

Back Type	Ang-Adj Pivot Ht Setting	Push Handles	Backrest Height	Avail. Ht. Range
Standard - (Fixed)	N/A	Integral	Low	9"-12"
			Medium	12"-15"
			Tall	16"-19"
		Omit	Low	9"-12"
			Medium	12"-15"
Angle Adjustable (3"-5" Pivot)	3"	Integral	Tall	16"-19"
			Low	12"-14"
			Medium	15"-17"
	4"	Integral	Tall	17"-19"
			Low	13"-15"
			Medium	16"-18"
	5"	Integral	Tall	18"-20"
			Low	14"-16"
			Medium	17"-19"
	3"	Omit	Tall	19"-20"
			Low	12"-14"
			Medium	15"-17"
	4"	Omit	Tall	17"-19"
			Low	13"-15"
			Medium	16"-18"
	5"	Omit	Tall	18"-20"
			Low	14"-16"
			Medium	17"-19"
Angle Adjustable (5"-7" Pivot)	5"	Integral	Tall	19"-20"
			Low	14"-16"
			Medium	17"-19"
	6"	Integral	Tall	19"-20"
			Low	15"-17"
			Medium	18"-20"
	7"	Integral	Tall	N/A
			Low	16"-18"
			Medium	18"-20"
	5"	Omit	Tall	N/A
			Low	14"-16"
			Medium	17"-19"
	6"	Omit	Tall	19"-20"
			Low	15"-17"
			Medium	18"-20"
Folding Push Handle - (Stroller Handle w/ Knuckle)	N/A	N/A	Tall	N/A
			Low	16"-18"
			Medium	18"-20"
			Low	14"-16"
			Medium	17"-19"

✓	Option #	Definition	Price	✓	Option #	Definition	Price
11. APPENDIX							

Quickie 2 Seat To Floor Height Matrix - Standard Frame *Hemi options will result in a reduction of approx 2 1/2"

All 3" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	16.5	17.0
	1/2"	-	-	-	17.0	17.5
	1"	-	-	-	17.5	18.0
Multi Position Tall	1 1/2"	-	-	-	18.0	18.5
	Std	-	-	-	18.5	19.0
	1/2"	-	-	-	19.0	19.5
3"	1"	-	-	-	19.5	20.0
	1 1/2"	-	-	-	20.0	20.5
	Std	-	-	-	15.0	-
5 1/4"	1/2"	-	-	-	15.5	-
	1"	-	-	-	16.0	-
	1 1/2"	-	-	-	16.5	-
6"	Std	-	-	-	17.5	-
	1/2"	-	-	-	18.0	-
	1"	-	-	-	18.5	-
7"	1 1/2"	-	-	-	19.0	-
	Std	-	-	-	19.5	-
	1/2"	-	-	-	20.0	-
8" x 2"	1"	-	-	-	20.5	-
	1 1/2"	-	-	-	-	-
	Std	-	-	-	-	-
Frog Legs	1/2"	-	-	-	16.5	17.0
	1"	-	-	-	18.0	18.5
	1 1/2"	-	-	-	18.5	-
Integral Multi Pos Forks (Low)	Std	-	-	-	15.5	16.0
	1/2"	-	-	-	16.0	16.5
	1"	-	-	-	16.5	17.0
Integral Multi Pos Forks (Std)	Std	-	-	-	16.5	17.0
	1/2"	-	-	-	17.0	17.5
	1"	-	-	-	17.5	18.0

All 4" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	16.5	17.0
	1/2"	-	-	-	17.5	18.0
	1"	-	-	-	18.0	18.5
Multi Position Tall	1 1/2"	-	-	-	18.5	19.0
	Std	-	-	-	18.5	19.0
	1/2"	-	-	-	19.0	19.5
3"	1"	-	-	-	19.5	20.0
	1 1/2"	-	-	-	20.0	20.5
	Std	-	-	-	20.5	21.0
5 1/4"	1/2"	-	-	-	17.0	17.5
	1"	-	-	-	18.0	18.5
	1 1/2"	-	-	-	19.0	19.5
6"	Std	-	-	-	17.5	18.0
	1/2"	-	-	-	18.5	19.0
	1"	-	-	-	19.0	19.5
7"	1 1/2"	-	-	-	19.5	20.0
	Std	-	-	-	18.5	19.0
	1/2"	-	-	-	19.5	20.0
8" x 2"	1"	-	-	-	20.0	20.5
	1 1/2"	-	-	-	20.5	21.0
	Std	-	-	-	-	-
Frog Legs	1/2"	-	-	-	16.5	17.0
	1"	-	-	-	18.0	18.5
	1 1/2"	-	-	-	18.5	19.0
Integral Multi Pos Forks (Low)	Std	-	-	-	16.0	16.5
	1/2"	-	-	-	16.5	17.0
	1"	-	-	-	17.0	17.5
Integral Multi Pos Forks (Std)	Std	-	-	-	16.5	17.0
	1/2"	-	-	-	17.5	18.0
	1"	-	-	-	18.0	18.5

All 5" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	16.5	17.0
	1/2"	-	-	-	17.5	18.0
	1"	-	-	-	18.5	19.0
Multi Position Tall	1 1/2"	-	-	-	19.0	19.5
	Std	-	-	-	18.5	19.0
	1/2"	-	-	-	20.0	20.5
3"	1"	-	-	-	20.5	21.0
	1 1/2"	-	-	-	21.0	21.5
	Std	-	-	-	-	-
5 1/4"	1/2"	-	-	-	17.5	18.0
	1"	-	-	-	18.5	19.0
	1 1/2"	-	-	-	19.5	20.0
6"	Std	-	-	-	18.0	18.5
	1/2"	-	-	-	19.0	19.5
	1"	-	-	-	19.5	20.0
7"	1 1/2"	-	-	-	20.0	20.5
	Std	-	-	-	19.0	19.5
	1/2"	-	-	-	20.0	20.5
8" x 2"	1"	-	-	-	20.5	21.0
	1 1/2"	-	-	-	21.0	21.5
	Std	-	-	-	-	-
Frog Legs	1/2"	-	-	-	17.5	18.0
	1"	-	-	-	18.0	18.5
	1 1/2"	-	-	-	19.0	19.5
Integral Multi Pos Forks (Low)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Integral Multi Pos Forks (Std)	Std	-	-	-	17.0	17.5
	1/2"	-	-	-	18.0	18.5
	1"	-	-	-	18.5	19.0

All 6" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	17.5	18.0
	1/2"	-	-	-	18.0	18.5
	1"	-	-	-	18.5	19.0
Multi Position Tall	1 1/2"	-	-	-	19.0	19.5
	Std	-	-	-	19.5	20.0
	1/2"	-	-	-	20.0	20.5
3"	1"	-	-	-	20.5	21.0
	1 1/2"	-	-	-	21.0	21.5
	Std	-	-	-	21.5	22.0
5 1/4"	1/2"	-	-	-	18.5	19.0
	1"	-	-	-	19.0	19.5
	1 1/2"	-	-	-	19.5	20.0
6"	Std	-	-	-	20.0	20.5
	1/2"	-	-	-	20.5	21.0
	1"	-	-	-	21.0	21.5
7"	1 1/2"	-	-	-	21.5	22.0
	Std	-	-	-	19.0	19.5
	1/2"	-	-	-	20.0	20.5
8" x 2"	1"	-	-	-	20.5	21.0
	1 1/2"	-	-	-	21.0	21.5
	Std	-	-	-	21.5	22.0
Frog Legs	1/2"	-	-	-	18.5	19.0
	1"	-	-	-	19.0	19.5
	1 1/2"	-	-	-	20.0	20.5
Integral Multi Pos Forks (Low)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Integral Multi Pos Forks (Std)	Std	-	-	-	17.5	18.0
	1/2"	-	-	-	18.0	18.5
	1"	-	-	-	18.5	19.0

All 7" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	-	19.0
	1/2"	-	-	-	-	19.5
	1"	-	-	-	-	20.0
Multi Position Tall	1 1/2"	-	-	-	-	20.5
	Std	-	-	-	-	21.0
	1/2"	-	-	-	-	21.5
3"	1"	-	-	-	-	22.0
	1 1/2"	-	-	-	-	-
	Std	-	-	-	-	-
5 1/4"	1/2"	-	-	-	-	19.0
	1"	-	-	-	-	19.5
	1 1/2"	-	-	-	-	20.0
6"	Std	-	-	-	-	20.5
	1/2"	-	-	-	-	21.0
	1"	-	-	-	-	21.5
7"	1 1/2"	-	-	-	-	22.0
	Std	-	-	-	-	19.5
	1/2"	-	-	-	-	20.0
8" x 2"	1"	-	-	-	-	21.0
	1 1/2"	-	-	-	-	21.5
	Std	-	-	-	-	22.0
Frog Legs	1/2"	-	-	-	-	16.5
	1"	-	-	-	-	17.0
	1 1/2"	-	-	-	-	17.5
Integral Multi Pos Forks (Low)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Integral Multi Pos Forks (Std)	Std	-	-	-	-	19.0
	1/2"	-	-	-	-	19.5
	1"	-	-	-	-	20.0

All 8" x 1.5" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Multi Position Tall	1 1/2"	-	-	-	-	-
	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
3"	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
	Std	-	-	-	-	-
5 1/4"	1/2"	-	-	-	-	19.0
	1"	-	-	-	-	19.5
	1 1/2"	-	-	-	-	20.0
6"	Std	-	-	-	19.5	20.0
	1/2"	-	-	-	20.0	20.5
	1"	-	-	-	20.5	21.0
7"	1 1/2"	-	-	-	21.0	21.5
	Std	-	-	-	20.0	20.5
	1/2"	-	-	-	20.5	21.0
8" x 2"	1"	-	-	-	21.5	22.0
	1 1/2"	-	-	-	22.0	-
	Std	-	-	-	-	-
Frog Legs	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Integral Multi Pos Forks (Low)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Integral Multi Pos Forks (Std)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-

11. APPENDIX

Quickie 2 Seat To Floor Height Matrix - Standard Frame *Hemi options will result in a reduction of approx 2 1/2"

All 8" x 2" Casters						All Other 8" Casters							
Fork Size	Stem Size	Bolt Position					Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1			5	4	3	2	1
Multi Position Short	Std	-	-	-	-	-	Multi Position Short	Std	-	-	-	-	-
	1/2"	-	-	-	-	-		1/2"	-	-	-	-	-
	1"	-	-	-	-	-		1"	-	-	-	-	-
Multi Position Tall	Std	-	-	-	-	-	Multi Position Tall	Std	19.5	20.0	20.5	21.0	21.5
	1/2"	-	-	-	-	-		1/2"	-	-	-	21.5	22.0
	1"	-	-	-	-	-		1"	-	-	-	22.0	-
3"	Std	-	-	-	-	-	3"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-		1/2"	-	-	-	-	-
	1"	-	-	-	-	-		1"	-	-	-	-	-
5 1/4"	Std	-	-	-	-	-	5 1/4"	Std	-	-	-	-	20.0
	1/2"	-	-	-	-	-		1/2"	-	-	-	-	20.5
	1"	-	-	-	-	-		1"	-	-	-	-	21.0
6"	Std	-	-	-	-	-	6"	Std	-	-	-	20.0	20.5
	1/2"	-	-	-	-	-		1/2"	-	-	-	20.5	21.0
	1"	-	-	-	-	-		1"	-	-	-	21.0	21.5
7"	Std	-	-	-	-	-	7"	Std	-	20.0	20.5	21.0	21.5
	1/2"	-	-	-	-	-		1/2"	-	-	-	21.5	22.0
	1"	-	-	-	-	-		1"	-	-	-	22.0	-
8" x 2"	Std	-	-	-	20.0	20.5	8" x 2"	Std	-	-	-	-	-
	1/2"	-	-	-	20.5	21.0		1/2"	-	-	-	-	-
	1"	-	-	-	21.0	21.5		1"	-	-	-	-	-
Frog Legs	Std	-	-	-	21.5	22.0	Frog Legs	Std	-	-	-	-	-
	1/2"	-	-	-	-	-		1/2"	-	-	-	-	-
	1"	-	-	-	-	-		1"	-	-	-	-	-
Integral Multi Pos Forks (Low)	Std	-	-	-	-	-	Integral Multi Pos Forks (Low)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-		1/2"	-	-	-	-	-
	1"	-	-	-	-	-		1"	-	-	-	-	-
Integral Multi Pos Forks (Std)	Std	-	-	-	-	-	Integral Multi Pos Forks (Std)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-		1/2"	-	-	-	-	-
	1"	-	-	-	-	-		1"	-	-	-	-	-



Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2 • Concord, ON L4K 4P8 • Canada