

CTRL+5 WITH VOICE ACTIVATION

POWER SEATING
MADE SIMPLE —
JUST SAY THE WORD!



VOICE ACTIVATION: A FIRST IN THE INDUSTRY

Now standard, voice activation lets you control up to 10 seating movements. Just say "Hello Sunrise" to activate, then follow with "Sunrise 1" through "Sunrise 10" to tilt, recline, elevate legs, and more—no scrolling or reaching required. Each voice command moves a seat function for a set amount of time that you can personalize.

SEAMLESS CONTROL OF POWER SEATING

CTRL+5 with Voice Activation provides intuitive control of your power seating. Compatible with QUICKIE, ZIPPIE and Magic Mobility power chairs with the R-net system.

EFFORTLESSLY CHANGE POSITIONS

Switching between drive and seating modes can be time consuming, especially for riders using sip-and-puff, a head array or a mini joystick. Use Voice Activation to change positions without leaving drive mode.

SUPPORT FOR COMPLEX DIAGNOSES

For riders with spinal cord injuries, multiple sclerosis, muscular dystrophy, or limited mobility, voice commands offer a hands-free way to access seating functions, promote independence, and improve pressure relief.

FLEXIBLE SETUP, FAMILIAR FEEL

With push button or toggle switch options and USB-C connectivity, CTRL+5 adapts to your setup. If you've used Siri or Alexa, you'll feel right at home with how easy it is to speak and go. For soft speakers, a lapel microphone can be added.



PRESSURE RELIEF MADE EASY

With Memory Seating on the SEDEO PRO ADVANCED (Q700 M and Q700-UP M), CTRL+5 lets you move into your memory seat position using a single voice command. Recall positions like "pressure relief" or "transfer to bed" completely hands-free.

STAY IN CONTROL WITH A SIMPLE "STOP"

Say "stop" at any time to immediately halt movement and turn off voice control. To reactivate, just say "Hello Sunrise." You can also use the CTRL+5 buttons at any time as an alternative.



LEARN MORE:



CTRL+5 Webpage



This document can be viewed in PDF format at www.SunriseMedical.com.

SUNRISE MEDICAL CARES ABOUT YOUR SAFETY.
Please read your owner's manual for safe operating instructions and maintenance guidelines.



Sunrise Medical (US), LLC
MK-130346 Rev.A ©07.2025



800.333.4000

www.SunriseMedical.com

SWITCH-IT®
SMART TECHNOLOGIES