



Quickie QS5 X October 2026

EIF7



Mark For:

Date: _____
Dealer Acct #: _____
Dealer: _____
Dealer Contact: _____
Dealer Address: _____
Dealer City: _____ Prov: _____ PC: _____
Dealer Phone: () _____ Fax: () _____
Confirmation Email: _____
Confirm Via: ☐ Fax ☐ Email

Submitting for:

☐ Quote ☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
Attention: _____
Address: _____
Address: _____
Ship To City: _____ Prov _____ PC: _____
Ship To Phone: () _____ Fax: () _____

✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

1. FRAME

Base model

☐ 170WM112 Quickie QS5 X \$3,344.00

Frame Style

See figure 1. Measure from front frame to floor

☐	170FS31	Swing Away Frame	NC
☐	170FS43	Fixed Front 75°	NC
☐	170FS44	Fixed Front 85°	NC

Frame Inset (Fixed Front Only)

See FIGURE 2 Frame inset is the distance from the outside of the seat tube to the outside of the footrest. Footrest width = frame width minus 2. N/A for 14" and 15" Wide using composite angle adj footrest. 75° only available with 0" inset.

☐	170FF1	0"	NC
☐	170FF2	1"	NC

Transit Option

ISO 7176-19 Transit Approved

☐ 170TO49 Transit System NC

Frame Options

Ergo Touch Points is a rubber cover on the knee bend area of the frame tubes

☐	1709A11	Ergo Touch Points Swing Away hangers or Fixed front	\$198.00
☐	1709A12	Frame Fold-N-Lock mechanism	NC
☐	1709A06	Twin LED lights	\$310.00



QS5 X - Swing Away



QS5 X - Fixed Front

2. SIZE / SEATING

Frame Width

See FIGURE 2 Measured from outside of seat tube to the outside of opposite seat tube.

☐	170W4	13"	NC
☐	170W5	14"	NC
☐	170W6	15"	NC
☐	170W7	16"	NC
☐	170W8	17"	NC
☐	170W9	18"	NC
☐	170W10	19"	NC
☐	170W11	20"	\$308.00

Seat Depth

See Seat Squeeze Matrix in Appendix section. Measured from back canes to front edge of upholstery. Frame Fold-N-Lock mechanism NA with 13"-14".

☐	170SD6	13"	NC
☐	170SD7	14"	NC
☐	170SD8	15"	NC
☐	170SD9	16"	NC
☐	170SD10	17"	NC
☐	170SD11	18"	NC
☐	170SD12	19"	NC
☐	170SD13	20"	\$308.00

2. SIZE / SEATING (CONT)

Seat Sling

Standard seat sling includes utility bag sewn at the bottom of upholstery.

☐	170SL1	Standard Seat Sling	NC
☐	170SL6	Omit Seat Sling	NC




Figure 1

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.

Cushion name or part number must be filled out with desired width and depth

☐	170JC20	JAY Bundled Cushion	Price List
		Cushion Name or Part #	
		Width:	Depth:



Figure 2

3. CASTERS / SEAT HEIGHT

Front Seat Height

☐	170SH41	14.0"	NC	☐	170SH18	18.0"	NC
☐	170SH42	14.5"	NC	☐	170SH20	18.5"	NC
☐	170SH12	15.0"	NC	☐	170SH21	19.0"	NC
☐	170SH13	15.5"	NC	☐	170SH23	19.5"	NC
☐	170SH14	16.0"	NC	☐	170SH24	20.0"	NC
☐	170SH15	16.5"	NC	☐	170SH26	20.5"	NC
☐	170SH16	17.0"	NC	☐	170SH27	21.0"	NC
☐	170SH17	17.5"	NC	☐	170SH29	21.5"	NC
				☐	170SH36	22.0"	NC

Caster Wheels (option A)

☐	170CW3	4" Low Profile Polyurethane	NC
☐	170CW77	4" x 1.25" Semi Pneumatic	\$110.00
☐	170CW5	5" Low Profile Polyurethane	NC
☐	170CW78	5" x 1.5" Semi-Pneumatic	\$110.00
☐	170CW34	6" Polyurethane	NC
☐	170CW33	6" Pneumatic	\$110.00
☐	170CW79	6" x 1.5" Semi Pneumatic	\$110.00
☐	170CW16	7" Polyurethane	NC
☐	170CW32	8" Polyurethane	NC

Caster Fork


☐	170CF5	Multiple Position	NC
☐	170CF69	Multiple Position Single Sided Fork	\$176.00
☐	170CF53	Frog Legs 0 weight to 100 lbs.	\$574.00
☐	170CF54	Frog Legs 101 lbs. weight to 190 lbs.	\$574.00
☐	170CF55	Frog Legs 191 lbs. weight to 265 lbs.	\$574.00

Caster Wheels colors (option B)

Casters are silver color unless you choose color option. Refer to Anodize package section on Page 7 for Anodize color selection. Match Ano color pkg pick is N/A with Black option.

☐	170CW64	4" x 1.5" Aluminum Soft Roll	\$110.00
☐	170CW66	5" x 1.5" Aluminum Soft Roll	\$110.00
☐	170CW68	6" x 1.5" Aluminum Soft Roll	\$110.00
	Colors:	☐ Match anodize color package	Package price
		☐ Match frame color	NC
		☐ Black anodize caster only	\$68.00

Multi Position Single side Frog Leg



Seat Squeeze Matrix

Seat Depth	Max Seat Squeeze
13"-15"	2.0"
16"-19"	2.5"
20"	3.0"

✓	Option #	Definition	Price	✓	Option #	Definition	Price
4. REAR WHEELS / SEAT HEIGHT							

Rear Seat Height						
☐	170RH13	14.0"	NC	☐	170RH19	17.0"
☐	170RH14	14.5"	NC	☐	170RH26	17.5"
☐	170RH15	15.0"	NC	☐	170RH20	18.0"
☐	170RH16	15.5"	NC	☐	170RH27	18.5"
☐	170RH17	16.0"	NC	☐	170RH21	19.0"
☐	170RH18	16.5"	NC	☐	170RH37	19.5"
				☐	170RH22	20.0"

Axle			
☐	170AX1	Quick Release Steel	NC
☐	170AX11	Quad Quick Release	\$182.00

Rear Wheel Type			
-----------------	--	--	--

Max weight capacity on Proton wheels is 220 lbs.

☐	170RW103	20" Mag Wheel 5 Spoke	NC
☐	170RW104	22" Mag Wheel 5 Spoke	NC
☐	170RW102	24" Mag Wheel 5 Spoke	NC
☐	170RW158	20" Lite Spoke Wheel	\$320.00
☐	170RW159	22" Lite Spoke Wheel	\$320.00
☐	170RW92	24" Lite Spoke Wheel	NC
☐	170RW93	25" Lite Spoke Wheel	NC
☐	170RW160	26" Lite Spoke Wheel	\$126.00
☐	170RW95	24" Spinergy SPOX Black (silver hub)	\$864.00
☐	170RW96	25" Spinergy SPOX Black (silver hub)	\$864.00
☐	170RW98	24" Spinergy LX Black (silver hub)	\$1,111.00
☐	170RW99	25" Spinergy LX Black (silver hub)	\$1,111.00
☐	170RW148	24" Proton Wheel	\$218.00
☐	170RW149	25" Proton Wheel	\$218.00
☐	170RW25	Omit Rear Wheel and Axle *	NC

Spinergy Spoke Colors (Optional)			
----------------------------------	--	--	--

Please select one of the following Spinergy wheel options first. Includes black hub.

Spinergy Spox spokes color:

☐ Black ☐ Blue ☐ Red ☐ White ☐ Yellow

Spinergy LX spokes color:

☐ Black ☐ Blue ☐ Green ☐ Orange ☐ Pink

☐ Purple ☐ Red ☐ White ☐ Yellow

Rear Wheel Tire			
☐	170RT1	Pneumatic	NC
☐	170RT24	Pneumatic Airless Insert - Infinergy® Lightweight - E-TPU material	\$247.00
☐	170RT3	Full Poly	NC
☐	170RT38	SHOX Performance Poly	\$146.00
☐	170RT25	Primo Orion Poly	\$271.00
☐	170RT39	High-Pressure Clincher Black	\$105.00
☐	170RT21	Schwalbe Marathon Plus Evolution	\$157.00

Handrim			
---------	--	--	--

Natural Fits available in 22", 24", and 25". Natural Fit w/ 5 spoke MAG only available in 24". Not all handrims available on all wheels. See page 6 for matrices.

☐	170HR1	Aluminum Anodized	NC
☐	170HR29	Aluminum Powder coated colored	\$95.00
☐	170HR4	Plastic Coated	\$146.00
☐	170HR22	Natural Fit - Standard Grip	\$374.00
☐	170HR23	Natural Fit - Super Grip	\$434.00
☐	170HR39	Ellipse 3R	\$550.00
☐	170HR40	Carbolife - Gekko	\$575.00
☐	170HR41	Carbolife - Curve L Tetra Grip	\$625.00
☐	170HR28	Q-Grip Neoprene Coated	\$471.00
☐	170HR7	Omit Handrim Option may not be available on all wheels	NC

Camber			
☐	170CI1	0°	NC
☐	170CI13	2° (default)	NC
☐	170CI3	4°	NC

Center of Gravity Preset			
☐	170CC20	1.0" (default)	NC
☐	170CC15	1.5"	NC
☐	170CC17	2.0"	NC
☐	170CC18	2.5"	NC
☐	170CC19	3.0"	NC

5. WHEEL LOCKS

Wheel Lock			
------------	--	--	--

No Gos may exist depending on wheel size and position. Ergo & compact N/A w/ 13"-14" frame depth

☐	170WL1	High Mount Push	NC
☐	170WL2	High Mount Pull	NC
☐	170WL45	Ergo Scissor	NC
☐	170WL88	Fold Down Push (Knee lever)	\$110.00
☐	170WL30	Compact	\$162.00
☐	170WL10	Omit	NC

Wheel Lock Option			
☐	170WL12	Do Not Mount	NC

Extension Handle			
------------------	--	--	--

N/A w/ ergo scissor or compact wheel locks

☐	170EH7	Extension Handle (Pair)	\$73.00
--------------------------	--------	-------------------------	---------

✓	Option #	Definition	Price	✓	Option #	Definition	Price
6. HANGER & FOOTREST							
Swing Away Hangers				Footrest Extension Tubes			
<i>Swing Away frame only</i>				<i>See Extension matrix for leg length range</i>			
☐	170H145	60° Swing In/Out Hanger	NC	☐	170E1	Short	NC
☐	170H139	70° Swing In/Out Hanger	NC	☐	170E2	Medium	NC
☐	170H141	80° Swing In/Out Hanger	NC	☐	170E5	Long	NC
☐	170H160	Elevating Hanger - Ext Mount.	\$339.00	Footrest Options			
☐	170H159	Articulating Hanger - Ext Mount	\$526.00	☐	170HL1	Heel Loop <i>Swing away only</i>	NC
☐	170H32	Omit Hanger and Footplate	NC	☐	170HL5	Leg Strap	NC
Footplates				☐	107DI04	Extra Leg Strap	\$37.00
<i>Heel loop or leg strap not included, See Footplate Options Section.</i>				☐	170IG2	Impact Guard Neoprene	\$61.00
<i>Platform footplates flip up (sideways).</i>				☐	170DI06	Platform Position Plates	\$68.00
☐	170F1	Composite	NC	☐	170F35	Auto-Folding Platform Carbon fiber only	\$89.00
☐	170F271	Composite Angle Adj-(16"-22"w)	\$176.00	Articulating / Elevating Hanger Cover			
☐	170F5	Aluminum (ELR & ALR)	NC	☐	170DI47	Neoprene Cover	\$325.00
☐	170F3	Alum Adult Angle Adjustable - Ext Mount	\$176.00	☐	170DI48	Gel Pad Cover	\$380.00
☐	170F259	Multi Angle Adjustable ☐ Black <i>NC</i> ☐ Std Frame color\$51	\$176.00	Lower Leg Length (LLL)			
☐	170F50	Platform Footrest ☐ Black <i>NC</i> ☐ Std Frame color\$51	\$152.00	<i>Lower Leg Length is measured from seat sling to the top rear of the footrest. See Appendix on Page 7. A minimum 2" ground clearance is required below the footrest. If the selected configuration results in less than 2" of ground clearance the footrest will be raised to 2" ground clearance, which may alter the lower leg length measurement. Available lower leg length is dependent on frame selection including adapter spacers. Reference LLL matrix in page 7</i>			
☐	107F335	Platform Ultra Lite Footrest	\$176.00	☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19"			
☐	107F336	Platform Carbon Fiber footrest	\$151.00	7. BACKREST / SEATING			
☐	170F7	Omit Footplate and Ext. Tube	NC	Back Type			
Backrest Upholstery				☐	107BT123	Standard Fix (90°)	NC
<i>EXO standard and Tension adjustable upholstery binding is always black.</i>				☐	170BT47	Angle Adjustable Back (80° - 105°)	NC
☐ 170BU95 Standard Upholstery EXO NC				Select desired setting ☐ 80° ☐ 85° ☐ 90° ☐ 100°			
Color Binding: ☐ Silver ☐ Orange ☐ Red ☐ Blue				☐	107BT124	Fixed Half Fold down (to the front)	\$246.00
☐ 170BU70 Backrest Tension Adjust Ballistic \$190.00				☐	170BT49	Lumbar curve bend	NC
☐ 170BU89 Backrest Tension Adj EXO \$449.00				Select desired setting ☐ 5° forward ☐ 5° backward			
Color Binding: ☐ Silver ☐ Orange ☐ Red ☐ Blue				Back Height Setting			
☐ 170BU90 Backrest Tension Adj EXO Pro 3DX \$449.00				<i>Half Fold Down N/A with the 10" - 13" back heights.</i>			
Color Binding: ☐ Silver ☐ Orange ☐ Red ☐ Blue				<i>Lumbar curve backrest with integral push handles available with 12" thru 20" back heights.</i>			
☐ 170BU14 Omit Back Upholstery NC				<i>Lumbar curve backrest with Fold Down or Omit push handles available with 14" thru 20" back heights.</i>			
JAY Back (Optional)				☐	170B161	Setting 10 "- 12" Back Height NC ☐ 10" ☐ 11" ☐ 12"	
<i>To order a J3 Back please select Omit Back Upholstery. See J3 Back Order Form for Part #'s</i>				☐	170B162	Setting 12 "- 14" Back Height NC ☐ 12" ☐ 13" ☐ 14"	
☐ 170JB01 J3 Back (See J3 Order Form) Price List				☐	170B163	Setting 14 "- 16" Back Height NC ☐ 14" ☐ 15" ☐ 16"	
				☐	170B164	Setting 16 "- 18" Back Height NC ☐ 16" ☐ 17" ☐ 18"	
				☐	170B165	Setting 18 "- 20" Back Height NC ☐ 18" ☐ 19" ☐ 20"	
Push Handles				Push Handles			
<i>Std EXO Uph N/A with Fold down. Height Adjustable push handles available with EXO / EXO Pro Tension Adjus</i>				☐	170PH6	Integral Push Handles	NC
☐	170PH6	Integral Push Handles	NC	☐	170PH5	Fold Down Push Handles <i>Available w/12"-20" back height</i>	\$168.00
☐	170PH5	Fold Down Push Handles <i>Available w/12"-20" back height</i>	\$168.00	☐	170PH21	Height adjustable push Handles <i>Available w/ 12"-20" back height. N/A w/Half fold down back</i>	\$201.00
☐	170PH21	Height adjustable push Handles <i>Available w/ 12"-20" back height. N/A w/Half fold down back</i>	\$201.00	☐	170PH25	Height adjustable push Handles for JAY Backs <i>N/A w/Half fold down back</i>	\$207.00
☐	170PH25	Height adjustable push Handles for JAY Backs <i>N/A w/Half fold down back</i>	\$207.00	☐	170PH22	Omit Push Handles - Straight Post	NC
☐	170PH22	Omit Push Handles - Straight Post	NC				

✓	Option #	Definition	Price	✓	Option #	Definition	Price
8. ARMREST / SIDE GUARDS							
Padded Swing Away <i>Some limitations when used w/ angle adj back in the rear most setting</i> <i>Please select a setting, will default to 11" if no selection is made.</i>				Side Guard <i>Side guards N/A w/ height adj. arms. Side guard interference w/ Padded SA arms at some heights may exist</i>			
☐	170AA162	Padded Swing-Away Fixed height (14")	NC	☐	170SG33	Plastic Standard	\$188.00
☐	170AA177	Padded Swing-Away STD (9" to 14")	\$278.00	☐	170SG29	Plastic Low	\$188.00
Setting ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14"				☐	170SG35	Aluminum	\$315.00
				☐	170SG36	Aluminum w/ Fender	\$315.00
				☐	170SG12	Carbon Fiber	\$399.00
				☐	170SG41	Carbon Fiber Long w/LG Fender	\$602.00
Armrest - Height Adjustable				Side Guard - Color Option <i>Available with Aluminum or Aluminum with Fender. Not applicable to carbon fiber</i>			
☐	170AA11	Single Post Hgt. Adj. Std (8.5"-13.5")	\$339.00	☐	170GE01	Black	NC
☐	170AA29	Single Post Hgt. Adj. Low (7"-11")	\$339.00	☐	170GE02	Matching Frame- Standard Colors	NC
☐	170AA65	Flip Back Height Adj <i>LITE</i> (8.5"-13.5")	\$381.00	☐	170GE03	Matching Frame- Personalized Color	\$68.00
Armpad - For Height Adjustable <i>Ergo armpad only available with Single post Height adj & Flip back armrest</i>							
☐	170AM23	Desk Length - Classic (10")	NC				
☐	170AM24	Full Length - Classic (14")	NC				
☐	170AM27	Ergo Pad (9.5"-13.5")	NC				
Omit Armrest							
☐	170AA14	Omit Armrest	NC				

9. KEY OPTIONS & ACCESSORIES

Anti-Tip					
☐	170AT1		Anti Tips (Pair)	\$121.00	
☐	170AT		Anti Tip swing away		
		☐	170AT88	Right Mount	\$152.00
		☐	170AT89	Left Mount	\$152.00
		☐	170AT90	Pair	\$305.00

Additional Accessories				
☐	170A1		Backpack - Adult	\$112.00
☐	170A2		Backpack - Kids	\$105.00
☐	1701A01		Seat Pouch Black	\$83.00

Additional Cover Accessories				
☐	1708A02		Gel Cover Hanger/Frame	\$162.00
☐	1708A01		Gel Cover Swing Away Arm	\$162.00

Positioning Belts				
☐	170ZZ15		JAY 2PT Belt 1.5" Push Btn 7.5" Pad	\$80.00
☐	170ZZ26		JAY 2PT Belt 2" Push Btn 9.5" Pad	\$126.00
☐	170ZZ14		JAY 2PT Belt 1.5" Push Btn No Pad	\$51.00
☐	170A118		Belt 2" Aircraft Buckle Pad w/Clamp	\$204.00
☐	170A120		Positioning Belt Slide Clamps	\$61.00

Spoke Guards				
Only available w/ Lite Spoke wheels. Not available with 5 spoke MAG wheels				
☐	1704A01		Quickie Black Spoke Guards	\$157.00
☐	1704A02		Quickie Clear Spoke Guards	\$157.00

Miscellaneous					
☐	1708A03		Caddy	\$89.00	
☐	170A2100		Crutch Holder	\$115.00	
☐	1708A21		Tool Kit	\$43.00	
☐	170A999		Quickie Water Bottle & Holder	\$81.00	
☐	1708A47		Mobile phone pocket	\$34.00	
☐	170FO4		Transit wheels (Pair)	NC	
☐	170AT		Tip Assist		
		☐	170AT83	Right Mount	\$20.00
		☐	170AT84	Left Mount	\$20.00
		☐	170AT91	Pair	\$41.00

✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

10. COLOR

Personalized X-Brace Color

Option has an extended leadtime. Frame colors only, NA with Personalized frame colors.
If no color pick crossbrace will be Gloss Black.

☐	170GS15	Crossbrace Color package	NC
--------------------------	---------	--------------------------	----

Standard Frame Colors

Handrim color option only available if Aluminum Coated colored option was selected.

			Frame	X-brace	Handrim
170C211		Celestial Sapphire	☐	☐	☐
170C212		Desert Sage	☐	☐	☐
170C213		Desert Tan	☐	☐	☐
170C214		Gloss Black	☐	STD	☐
170C215		Matte Black	☐	☐	☐
170C216		Metallic Shadow	☐	☐	☐
170C217		Midnight Blue	☐	☐	☐
170C218		Mystic Purple	☐	☐	☐
170C219		Petal Pink	☐	☐	☐
170C220		Pink Power	☐	☐	☐
170C221		Red Pepper	☐	☐	☐
170C222		Royal Blue	☐	☐	☐
170C223		Salute Green	☐	☐	☐
170C224		Scarlet Noir	☐	☐	☐
170C225		Silver Ore	☐	☐	☐
170C226		Solar Orange	☐	☐	☐
170C227		Starlight White	☐	☐	☐
170C228		Storm Gray	☐	☐	☐
170C229		Turbo Green	☐	☐	☐
170C230		Velocity Yellow	☐	☐	☐

NEW!

Legacy - Frame Colors

☐	170C141	Black Cherry	NC
☐	170C10	Black Opal	NC
☐	170C14	Candy Purple	NC
☐	170C142	Electric Blue	NC

Personalized Frame Color

Only frame will have the personalized color; Crossbrace, handrim, footplate and Aluminum side guard not included (unless selected).

☐	170C112	Desert Camouflage	\$368.00
☐	170C43	Mossy Oak Camouflage	\$368.00
☐	170C42	Stars & Stripes	\$368.00
☐	170C108	Zebra (Choose base color). N/A w/ Gloss Black and Matte Black.	\$368.00

Anodized Color Package

Axle plate, Multi Forks and Aluminum Soft Roll casters only. N/A with Frog suspension Legs.

☐	170EB08	Blue	\$212.00
☐	170EB13	Gold	\$212.00
☐	170EB09	Lime Green	\$212.00
☐	170EB10	Purple	\$212.00
☐	170EB11	Red	\$212.00
☐	170EB12	Ti Grey	\$212.00

Frame Decal Colors

Applies to front caster link decals and aluminum & carbon fiber sideguards.

☐	170D16	Black	NC
☐	170D17	Blue	NC
☐	170D18	Orange	NC
☐	170D19	Red	NC
☐	170D20	Silver	NC

Built For Me

To inquire about any unique personalized for Quickie products, contact:
Phone: (800) 263-3390, Email: cscanada@sunmed.com



11. APPENDIX

Wheel - Handrims Matrix

Handrim	Rear Wheel			
	5 Spoke Mag	Lite Spoke	Spinergy	Proton
Aluminum Ano / Painted	20,22,24	20,22,24,25,26	24,25	24,25
Plastic Coat	20,22,24	22,24	24,25	-
Natural Fit - Standard Grip	24	22,24,25,26	24,25	-
Natural Fit - Super Grip	24	22,24,25,26	24,25	-
Natural Fit - No Thumb Grip	24	22,24,25,26	24,25	-
Surge - Oval w/Gription Strip	-	22,24,25,26	24,25	-
Surge LT - Oval w/Gription Strip	-	22,24,25,26	24,25	-
Q Grip Neoprene Coated	-	24,25	24,25	-
Omit Handrim	20,22,24	-	24,25	-

Wheel - Tire Matrix

Tire	Rear Wheel			
	5 Spoke Mag	Lite Spoke	Spinergy	Proton
Pneumatic	20,22,24	20,22,24,26	24	24
Pneumatic w/ Airless Insert	20,22,24	20,22,24,26	24	-
Full Poly	20,22,24	20,22,24	24	-
SHOX Performance Poly	20,22,24	20,22,24,25	24,25	-
Primo Orion Poly	20,22,24	20,22,24	24	-
High-Pressure Clincher Blk	24	22,24,26	24	24
Schwalbe Marathon	22, 24	24,25,26	24,25	24,25
20" = 451mm 22" = 501mm 24" = 540mm 25" = 559mm 26" = 590mm				

Note 1: Fork Hole 1+ indicates that a 10mm spacer is used to achieve additional Front Seat Height (FSH).

Note 2: FSH values are rounded estimates and assume flat seat heights (0 degree seat angle). Increasing seat angle can affect final FSH.

FSH for 4" Casters						
Caster Link	Fork Size	Hole Position				
		4	3	2	1	1+
Hemi	Multi 90	-	-	14.0	14.5	15.0
	Multi 115	-	-	14.5	15.0	15.5
	Multi 131	-	-	15.0	15.5	16.0
	Multi 158	-	-	16.0	16.5	17.0
	SS Low	-	-	14.0	14.5	15.0
	SS Std	-	-	14.5	15.0	15.5
Short	Frog Legs	-	-	15.5	16.0	16.5
	Multi 90	-	-	15.5	16.0	16.0
	Multi 115	-	-	16.0	16.5	17.0
	Multi 131	-	-	16.5	17.0	17.5
	Multi 158	-	-	17.5	18.0	18.5
	SS Low	-	-	15.5	16.0	16.0
Tail	SS Std	-	-	16.0	16.5	17.0
	Frog Legs	-	-	17.0	17.5	18.0
	Multi 90	-	-	16.5	17.0	17.5
	Multi 115	-	-	17.0	17.5	18.0
	Multi 131	-	-	17.5	18.0	18.5
	Multi 158	-	-	18.5	19.0	19.5
Hemi	SS Low	-	-	16.5	17.0	17.5
	SS Std	-	-	17.0	17.5	18.0
	Frog Legs	-	-	18.0	18.5	19.0
	Multi 90	-	-	18.0	18.5	19.0
	Multi 115	-	-	18.5	19.0	19.5
	Multi 131	-	-	19.0	19.5	20.0
Short	Multi 158	-	-	19.5	20.0	20.5
	SS Low	-	-	17.5	18.0	18.5
	SS Std	-	-	18.0	18.5	19.0
	Frog Legs	-	-	19.0	19.5	20.0
	Multi 90	-	-	19.0	19.5	20.0
	Multi 115	-	-	19.5	20.0	20.5
Tail	Multi 131	-	-	20.0	20.5	21.0
	Multi 158	-	-	20.5	21.0	21.5
	SS Low	-	-	20.0	20.5	21.0
	SS Std	-	-	20.5	21.0	21.5
	Frog Legs	-	-	21.0	21.5	22.0
	Multi 90	-	-	21.0	21.5	22.0

FSH for 7" Casters						
Caster Link	Fork Size	Hole Position				
		4	3	2	1	1+
Hemi	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	-	-	-
	Multi 158	17.0	17.5	18.0	18.5	19.0
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
Short	Frog Legs	-	-	-	-	-
	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	-	-	-
	Multi 158	18.5	19.0	19.5	20.0	20.5
	SS Low	-	-	-	-	-
Tail	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-
	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	-	-	-
	Multi 158	19.5	20.0	20.5	21.0	21.5
Hemi	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-
	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	-	-	-
Short	Multi 158	20.5	21.0	21.5	22.0	22.5
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-
	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
Tail	Multi 131	-	-	-	-	-
	Multi 158	21.5	22.0	22.5	23.0	23.5
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-
	Multi 90	-	-	-	-	-

Lower Leg Length Matrix				
Note for Swing Away Hangers only: subtract 1" from the max achievable LLL for long extension tubes if the Multi Angle Adjustable footplate is selected.				
Hanger	Extension	LLL (inches) (measured from top of seat sling to top of footrest)		Min Allowable FSH (due to 2" ground clearance requirement)
		Min	Max	
60° Swing In/Out	Short	10	14	No lower FSH limit
	Med	11	16	
	Long	12	18	
70° Swing In/Out	Short	10	16	No lower FSH limit
	Med	12	17	
	Long	13	19	
80° Swing In/Out	Short	11	15	No lower FSH limit
	Med	12	17	
	Long	14	19	
ELR	Short	11	16	No lower FSH limit
	Med	11	18	
	Long	13	19	
ALR	Short	11	16	No lower FSH limit
	Med	13	18	
	Long	15	19	
75° Fixed Front	Short	13	18	15.5
	Med	13	20	
	Long	14	20	
85° Fixed Front 0" inset	Short	13	18	15.5
	Med	13	20	
	Long	15	20	
85° Fixed Front 1" inset	Short	13	18	15.5
	Med	15	20	
	Long	17	20	