

ORDER FORM




Q700-UPM **EIPW25**
SEDEO[®] PRO ADVANCED[™]

Customer		Ship To	
Account #		Name	
Date		Address	
PO#		City	
Buyer		Province	Postal Code
Mark For		Phone	

Client Measurements

a) Top of head to seat	in
b) Top of shoulder to seat	in
c) Underarm to seat	in
d) Elbow to seat	in
e) Top of leg to seat	in
f) Shoulder width	in
g) Chest width	in
i) Hip width	in
j) Actual sitting depth	in
m) Foot length	in
n) Knee to foot	in
o) Seat width	in



Notes

Seating Specifications

Weights Starting at: (with Batteries)	435 lbs
Power Seating	25° of CG Tilt and Power Recline with 7" mechanical shear reduction / 8" seat elevation / 30° Anterior Tilt
Seat Type	SEDEO PRO ADVANCED
Seat Size Range	Width 16" to 20" Depth 17" to 22"
Seat-to-Floor Height Range	18" to 19"

Mid-Wheel Drive Base Specifications

Drive Wheel Position	Mid-wheel drive
User Weight Capacity	265 lbs. Standing
Base Size	Width 24.5", Length 36.5"
Turning Radius	19.8" (Excludes Footplates)
Drive Control	R-net PM 120A, R-net PM 90A
Driving Range	Up to 18 miles (w/ Group 24 Batteries)
Drive Motors	4-pole 6 mph, 8 mph
Drive Wheels	14" w/flat-free inserts (black)
Casters	6" front / 6" rear (black)
Charger	10A Dual Mode (Compact)
Battery Size	Group 24
Colour	9 shroud colours and 7 colour accents
Ground Clearance	3"

BASE			
OPTION		DESCRIPTION	RETAIL
Q700™ M Base <i>End Item #EIPW25</i>			\$
250PC335		QUICKIE Q700-UP M MPC Tilt/Recline for SEDEO PRO ADVANCED Included in Base Pricing: Q700 M Base, R-net 90amp electronics, ProGrip™ with 6-wheel independent suspension , 14" drive wheels with dual density inserts, 6" flat free casters, ISO Occupied Tie-Down brackets and 10amp charger. max speed = 6mph, weight capacity = 265lbs, STFH = 17", 3" ground clearance.	\$11,052.00
250PC336		QUICKIE Q700-UP M MPD Tilt/Recline for SEDEO PRO ADVANCED Included in Base Pricing: Q700 M Base, R-net 120amp electronics, ProGrip™ 6-wheel independent suspension , 14" drive wheels with dual density inserts, 6" flat free casters, ISO Occupied Tie-Down brackets and 10amp charger. max speed = 8mph, weight capacity = 265lbs, STFH = 17", 3" ground clearance.	\$17,607.00
User Weight			\$
250WG24		Up to 265 lbs.	NC
Motor Package			\$
250PF99		4-pole 10 kph / 6 mph	NC
250PF100		4-pole 12.5 kph / 8 mph	NC
250PF101		4-pole 12.5 kph / 8 mph	\$1,268.00
Tracking			\$
250MI01		SureTrac® Module & Motor Encoders (Improves drive tracking.)	\$1,696.00
Drive Wheels			\$
250DW58		Drive Wheel, 14", Dual Density, Knobby, Black	NC
Caster Wheels			\$
250CW55		Caster Wheels 6", Solid, Black	NC
Shroud Colour			\$
250C119		Matte Black	NC
250C165		Candy Apple Red	NC
250C166		Candy Blue	NC
250C167		Candy Purple	NC
250C176		Green Apple	NC
250C178		Pearl Pink	NC
250C182		Silver	NC
250C125		Pearl White	NC

BASE

OPTION	DESCRIPTION	RETAIL
Colour Accents		\$
250EZ20	Apple Green	NC
250EZ22	Black	NC
250EZ21	Blue	NC
250EZ23	Orange	NC
250EZ24	Red	NC
250EZ19	White	NC
250EZ25	Silver	NC
Service Decals		\$
250D10	Army	NC
250D11	Air Force	NC
250D13	Marines	NC
250D14	Navy	NC
250D15	Purple Heart	NC
ISO Transit Tie Downs		\$
250T050	ISO Transit Package	NC
LED Lights		\$
250LI11	LED Lights and Indicators <i>(Available with 4 Power Seat Functions)</i>	NC
250LI09	LED Fender Lights	NC
250LI02	Omit LED Lights	NC
Seat Heights		\$
<i>Measured from front of seat pan to floor.</i>		
250HT56	18"	NC
250HT105	18.5"	NC
250HT57	19"	NC
Batteries and Chargers		\$
Group 24 Battery Options (Recommended for max daily range)		
250AG25	GP24 Flat Top Pair w/Strap MK Gel	\$955.00
250AG26	GP24 Flat Top Pair w/Strap AGM	\$915.00
Chargers		
250BJ28	10A Dual Mode (Compact)	NC

POWERED SEATING

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

Power Seating Options

\$

Seat Elevator (Required)

250PI06	Seat Elevator, Powered, 8"	\$3,477.00
---------	----------------------------	------------

Seat Tilt/Recline (Required)

250PA51	Posterior 0° - 25° Tilt and 180° Recline Powered, w/ 7" Shear Reduction	\$13,275.00
---------	---	-------------

Seat Anterior Tilt / Standing (Must Choose One)

250TG03	Standing	\$7,910.00
---------	----------	------------

Back Frame Type

\$

250BF09	Standing Seat Frame	NC
---------	---------------------	----

Back Cushion Type

\$

250BT138	JAY Syncra® <i>(Choose Seating in section B)</i>	NC
250BT135	JAY® J3® <i>(Choose Seating in section C)</i>	NC
250BT137	SEDEO Aftermarket Adapter <i>(Height adjustable round tubes for traditional mounting)</i>	\$633.00
250BT133	OMIT <i>(Shear back plate will have 4 mounting posts for personalized mounting. No headrest mount.)</i>	NC

POWERED SEATING

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

B - JAY® Syncra™ Seating**B - JAY® Syncra™ Back**

\$

250JP14	JAY® Syncra™ Back 16-20" Wide with 3" Contour Depth	\$1,339.00
250JP15	JAY® Syncra™ Back 22" Wide with 3" Contour Depth	\$1,339.00

B - JAY® Syncra™ Back Width

\$

Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.

250JQ07	16" Wide (Adjustment range 16 to 17" Wide)*	NC
250JQ08	18" Wide (Adjustment range 18 to 19" Wide)*	NC
250JQ09	20" Wide (Adjustment range 20 to 21" Wide)*	NC
250JQ10	22" Wide (Adjustment range 22 to 23" Wide)	NC

*Additional +1" of shell width adjustment is achievable. Requires purchase of 1" larger foam and upholstery separately.

B - JAY® Syncra™ Back Shell Options and Height Ranges

\$

JAY® Syncra™ Back shells are mounted with a minimim 5.5" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field within the ranges shown based on the back frame type.

JAY® Syncra™ Shell Height

250JR42	16" Tall Shell (22" to 24" Back Height)	NC
250JR43	18" Tall Shell (24" to 26" Back Height)	NC
250JR44	20" Tall Shell (26" to 28" Back Height)	NC
250JR45	22" Tall Shell (28" to 30" Back Height)	NC

JAY® Syncra™ Back Installed Back Height (From Pan to Top of Back)

250JR50	22" (16" Shell Height)	NC
250JR51	23" (16" Shell Height)	NC
250JR52	24" (16" Shell Height)	NC
250JR53	24" (18" Shell Height)	NC
250JR54	25" (18" Shell Height)	NC
250JR55	26" (18" Shell Height)	NC
250JR56	26" (20" Shell Height)	NC
250JR57	27" (20" Shell Height)	NC
250JR58	28" (20" Shell Height)	NC
250JR67	28" (22" Shell Height)	NC
250JR68	29" (22" Shell Height)	NC
250JR69	30" (22" Shell Height)	NC

B - JAY® Syncra™ Back Upholstery

\$

250BU96	Syncra™ Premium Leatherette	NC
250BU97	Syncra™ Stretch	NC

POWERED SEATING

OPTION	DESCRIPTION	RETAIL
B- JAY J3 Back \$		
B- JAY J3 Back – Contour		
250JP01	J3 PA Back for 14"-20" Wide (J3PAXXXYYZZ)	\$912.00
250JP11	J3 PL Back for 14"-20" Wide (J3PLXXYYZZ)	\$986.00
250JP13	J3 PD Back for 14"-20" Wide (J3PDXXYYZZ)	\$1,068.00
B- JAY J3 Back – Width <i>Back can be ordered a different width than seat pan.</i> \$		
250JQ01	14"	NC
250JQ02	16"	NC
250JQ03	18"	NC
250JQ04	20"	NC
B- JAY J3 Shell Options and Height Ranges \$		
<i>JAY J3 shells are mounted with a minimim 4" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field within the ranges shown.</i>		
JAY J3 Shell Height		PRO ADVANCED
250JR41	15" Tall Shell (Mid Thoracic Tall)	21"-22" NC
250JR14	16.5" Tall Shell (Upper Thoracic Short)	21"-22" NC
250JR15	18" Tall Shell (Upper Thoracic Med)	22"-24" NC
250JR16	19.5" Tall Shell (Upper Thoracic Tall)	23"-25" NC
250JR17	21" Tall Shell (Shoulder Height)	25"-26" NC
250JR27	22.5" tall shell (<i>Shoulder Height</i>)	26"-27" NC
250JR28	24" tall shell (<i>Shoulder Height</i>)	28" NC
JAY J3 Installed Back Height (From Pan to Top of Back)		
250JR29	21"	NC
250JR19	22"	NC
250JR20	23"	NC
250JR21	24"	NC
250JR22	25"	NC
250JR23	26"	NC
250JR24	27"	NC
250JR25	28"	NC
250JR30	29"	NC
250JR31	30"	NC
B- JAY Anterior Thoracic Support <i>Classic-Structure Pull Down Harness (4 point strap, only avail. w/JAY backrest)</i> \$		
250K10	Medium (ATS8302)	\$322.00
250K11	Large (ATS8303)	\$322.00

POWERED SEATING

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

JAY® Lateral Thoracic Swing-Away Supports (Requires Pad, Upholstery and Mount)

\$

JAY® Lateral Thoracic Support Pads

250JW35	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right	\$75.00
250JW38	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left	\$75.00
250JW36	Lateral Thoracic (L x H) 100 mm x 152 mm / 4" x 6" Pad - Right	\$75.00
250JW39	Lateral Thoracic (L x H) 100 mm x 52 mm / 4" x 6" Pad - Left	\$75.00
250JW37	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right	\$75.00
250JW40	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left	\$75.00

JAY® Lateral Thoracic Support Pad Upholstery

250NU01	Premium Leatherette	NC
250NU02	Stretch	NC

JAY® Lateral Thoracic Swing-Away Mount

250JW41	Lateral Thoracic Swing Away Bracket - Right	\$334.00
250JW42	Lateral Thoracic Swing Away Bracket - Left	\$334.00

JAY Syncra Fixed Lateral Support Accessory**JAY Fixed Lateral Support Accessory Hardware** (JAY® Fixed Lateral Support is height adjustable)

250JW33	Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm / 3.5" x 8" - Right	\$334.00
250JW34	Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm / 3.5" x 8" - Left	\$334.00

*JAY Fixed Lateral Support is height adjustable

Headrests (Requires Pad, Upholstery, & Mount)

\$

Headrest Pad

250HE36	SEDEO® Pad Small With Cover (LxH) 240 x 150 mm / 10" x 6" <i>Can only be used with Multi Position Headrest Mount with Quick Release and SEDEO® Upholstery</i>	\$227.00
250HE38	Whitmyer® Plush Headrest Pad 203 mm / 8"	\$377.00
250HE39	Whitmyer® Plush Headrest Pad 254 mm / 10"	\$377.00
250HE44	Whitmyer® Plush Headrest Pad 356 mm / 14"	\$377.00
250HE46	Whitmyer® Adjust-A-Plush Headrest Pad Medium (Adjustable from 7 to 16" Wide x 4" Tall)	\$445.00

Headrest Pad Upholstery

250NU10	SEDEO® Upholstery	NC
250NU04	Syncra™ Premium Leatherette	NC
250NU05	Syncra™ Stretch	NC
250NU11	Whitmyer® Lycra	NC

Headrest Mount

250MH32	SEDEO® Multi-Position Headrest Mount with Quick Release (Only available with SEDEO® Pad)	\$557.00
250MH33	Whitmyer® Linx Headrest Mount with Quick Release	\$247.00
250MH34	Whitmyer® Linx-2 Headrest Mount with Quick Release	\$298.00
250MH35	Whitmyer® AXYS Headrest Mount with Quick Release	\$247.00

POWERED SEATING

OPTION	DESCRIPTION	RETAIL
Push Handles		\$
250PH23	Push Handles <i>(Not available with O2 holder)</i>	NC
Seat Frame Rails		\$
250SN14	Seat Rail Short <i>(For effective seat depth 14" to 18")</i>	\$806.00
250SN15	Seat Rail Long <i>(For effective seat depth 17" to 22")</i>	\$806.00
Seat Cushion Type		\$
250ST81	JAY® Syncra® Cushion (See Section E)	NC
250ST74	JAY® (See Section D)	NC
250ST75	Omit	NC
C - SEDEO PRO ADVANCED Seat Size and D - JAY® Seat Size		\$
Seat Pan Width		
250R53	16"	NC
250R54	17"	NC
250R55	18"	NC
250R56	19"	NC
250R57	20"	NC
Seat Pan Depth (Front of Seat Pan to Front of Back Cushion Effective Depth)		
250RF45	17"	NC
250RF46	18"	NC
250RF47	19"	NC
250RF48	20"	NC
250RF49	21" <i>(Not available with JAY® Syncra® Cushion)</i>	NC
250RF50	22"	NC
E - JAY® Syncra™ Cushion (JAY Your Way modifications available through Customer Service (1-800-263-3390))		\$
250CU89	JAY® Syncra™ Cushion featuring JAY Flow® Fluid with Cryo™ Technology <i>(16" to 20" wide)</i>	\$1,017.00
JAY® Syncra™ Cushion Cover Upholstery		
250YR73	Syncra™ Premium Leatherette	NC
250YR74	Syncra™ Stretch	NC
D- JAY Seat Cushion		\$
250CU59	JAY Ion® (Sizes 16" to 20" Wide)	\$657.00
250CU61	JAY Union® (Sizes 16" to 20" Wide)	\$764.00
250CU64	JAY Fusion® (Sizes 16" to 20" Wide)	\$1,045.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

JAY® Lateral Thigh Supports (Must pick pad, upholstery, and mount)

\$

JAY® Lateral Thigh Support Pads

250KF16	Lateral Thigh Support Pad Standard, Right (8.25" x 3.25")	\$146.00
250KF17	Lateral Thigh Support Pad Standard, Left (8.25" x 3.25")	\$146.00
250KF18	Lateral Thigh Support Pad Large, Right (12" x 3.25")	\$146.00
250KF19	Lateral Thigh Support Pad Large, Left (12" x 3.25")	\$146.00

JAY® Lateral Thigh Support Pad Upholstery

250NU07	Syncra™ Premium Leatherette	NC
250NU08	Syncra™ Stretch	NC

JAY® Lateral Thigh Removable Mount

250KF20	Lateral Thigh Removable Mount - Right	\$325.00
250KF21	Lateral Thigh Removable Mount - Left	\$325.00

Traditional Knee and Thigh Swing-Away Support Pads

\$

Lateral Knee Support Pads and Hardware (Must pick pad and mount)

250XY05	Lateral Knee Support Pad, Right (4" round)	\$247.00
250XY06	Lateral Knee Support Pad, Left (4" round)	\$247.00
250XX08	Lateral Swing-Away Bracket Knee - Right	\$325.00
250XX09	Lateral Swing-Away Bracket Knee - Left	\$325.00

Lateral Thigh Support Pads and Hardware (Must pick pad and mount)

250KF22	Lateral Thigh Support Pad Short, Right (5" H x 5" L)	\$146.00
250KF23	Lateral Thigh Support Pad Short, Left (5" H x 5" L)	\$146.00
250KF24	Lateral Thigh Support Pad Long, Right (5" H x 7" L)	\$146.00
250KF25	Lateral Thigh Support Pad Long, Left (5" H x 7" L)	\$146.00
250KF26	Lateral Thigh Removable Mount - Right	\$325.00
250KF27	Lateral Thigh Removable Mount - Left	\$325.00

Armrest Structure

\$

250EG51	Reclining, Flip-Back, Dual Link Armrest (Height adjustable 8" to 15" including length and lateral adjustment)	NC
---------	---	----

Arm Pads

\$

SEDEO® Armrest Pad

250AB43	300 mm / 12", Integral Foam, Right	NC
250AE43	300 mm / 12", Integral Foam, Left	NC
250AB44	400 mm / 16", Integral Foam, Right	NC
250AE44	400 mm / 16", Integral Foam, Left	NC

Comfort Arm Pad (Foam)

250AB46	2" Wide Right	\$113.00
250AE46	2" Wide Left	\$113.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	DESCRIPTION	RETAIL
Comfort Arm Pad (Gel)		
250AB48	2" Wide Right	\$213.00
250AE48	2" Wide Left	\$213.00
Ottobock® Arm Troughs		
250AE06	Ottobock Arm Trough, Med (12"), Left	\$139.00
250AB06	Ottobock Arm Trough, Med (12"), Right	\$139.00
250AE07	Ottobock Arm Trough, Large (14"), Left	\$139.00
250AB07	Ottobock Arm Trough, Large (14"), Right	\$139.00
2500L01	Ottobock Hand Pad, Med (7"), Left	\$70.00
2500R01	Ottobock Hand Pad, Med (7"), Right	\$70.00
2500L02	Ottobock Hand Pad, Large (9"), Left	\$70.00
2500R02	Ottobock Hand Pad, Large (9"), Right	\$70.00
250AE78	Ottobock 1 Piece Arm Trough (18"), Left	\$198.00
250AB78	Ottobock 1 Piece Arm Trough (18"), Right	\$198.00
Elbow Support \$		
2503A02	Elbow Support Pad, RH 3" x 10" <i>(Includes side support)</i>	\$228.00
2503A03	Elbow Support Pad, LH 3" x 10" <i>(Includes side support)</i>	\$228.00
2503A06	Elbow Support Pad, RH 3" x 3"	\$167.00
2503A07	Elbow Support Pad, LH 3" x 3"	\$167.00
2503A04	Swing-Away Elbow Support Hardware, RH	\$279.00
2503A05	Swing-Away Elbow Support Hardware, LH	\$279.00
Centermount Legrest \$		
250H173	Power Centermount Standing Legrest, 90° - 0° Angle, 7.5" Extension (12.5" to 17.5" Leg Length)	\$3,435.00
Footplates		
250F127	Dual, Height and Angle Adjustable, Standard (13.5" to 14.5" wide X 10" long)	\$284.00
250F128	Dual, Height and Angle Adjustable, Large (15.5" to 16.5" wide X 12" long)	\$284.00
Footplate Accessories		
250DI68	Extra Large Centermount Footplate <i>(Single Footplate Overlay 16"W x 12"D)</i>	\$315.00
250DI69	Footplate Connector Link <i>(Attaches footplates together)</i>	\$96.00
250DI57	Gel Footpads <i>(Pair, LH and RH)</i>	\$401.00
Calf Pad Options		
250DI34	Calf Pads, Pair, Standard	NC
Knee Block (Required for UP)		
250DI37	Knee Blocks: Dual, Height, Width and Angle Adjustable	\$704.00
250DI55	Knee Block Quick-Release Removable Hardware	\$179.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	DESCRIPTION	RETAIL
Pelvic Positioning Belts \$		
250A835	Positioning Belt, Aircraft Buckle 50mm x 1300mm / 2" x 50" (Included with UP)	NC
250A836	Positioning Belt, Aircraft Buckle, Padded 50mm x 1300mm / 2" x 50"	\$103.00
250A845	Bodypoint - Evoflex Pelvic Stabilizer Belt	\$216.00
250A846	Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	\$271.00
Chest Straps \$		
250ZR06	Chest Strap, 2-Point, Padded (Included with UP)	NC
250ZR07	Additional Chest Strap, 2-Point, Padded	\$351.00
Chest Bars \$		
250ZR08	Chest Bar, Detachable, Round Pad	\$704.00
250ZR09	Chest Bar, Detachable, Flat Pad, 200 mm x 100mm	\$748.00
Accessories \$		
Ventilator Option		
250VT11	Vent Tray Fixed (Pulmonetic LTV Models)	\$609.00
O2 Tank Option		
250A703	O2 Tank Holder to Backrest Styles (Not available with fixed push handles)	NC
Wheelchair Tray Table		
2507A08	Tray, Swing-Away	\$381.00
Backpack		
250A1	Black Utility	\$112.00
Misc. Accessories		
2508A52	Ram Mount Phone	\$77.00
2508A53	Ram Mount Tablet	\$308.00
2508A54	Ram Mount Cup Holder	\$77.00
2508A55	Rear View Mirror	\$77.00
2508A59	Transfer Bar (Pair)	\$493.00
2508A73	Armrest Bag	\$77.00
2509A10	Utility Hooks	\$30.00

JOYSTICK / ELECTRONICS			
OPTION		DESCRIPTION	RETAIL
R-net Joystick (Remote)			\$
250J116		Colour Joystick (2.3" Colour LCD Screen Display, CJM1)	NC
250J118		High Resolution Colour Joystick, (CJSM2+BT) 3.5" Colour LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and Dual (iOS/Android) Bluetooth Control	\$1,162.00
250J119		Omit Remote Control (Required with arm mounted specialty controls and table top control)	NC
Joystick Mounting			\$
250JM31		Fixed on Armrest	NC
250JM32		Parallel Swing-Away	\$447.00
250JM33		Link-It Mid-Line Mount on Armrest (Swing-Away)	\$447.00
250JM49		Attendant (Available with LED or CJSM1 only. Mounted to Push Handle)	NC
Joystick Mount (Location)			\$
250HM16		Right Hand	NC
250HM17		Left Hand	NC
Display Language Selection R-net			\$
250L1		English	NC
250L3		French	NC
250L7		Spanish	NC
Joystick Handle (Optional)			\$
250JH13		Ball Handle	\$110.00
250JH15		T-Handle	\$110.00
250JH24		Bodypoint – Rubber Dome Joystick Handle	\$216.00
250JH25		Bodypoint – U-Shaped Joystick Handle (4 inch)	\$216.00
250JH26		Bodypoint – Grooved Mushroom Joystick	\$216.00
Electronic Programmers			\$
250PR13		R-net PC Wireless Programmer	\$660.00
R-net Power Modules			\$
250ER17		Power Module, R-net, 90A, Expandable	\$102.00
250EX15		Expandable 90A Controller - R-net	NC
250EX16		Harness for 90A Expandable Controller	NC
250ER10		Power Module, R-net, 120A, Expandable (Only available with Group 4)	\$102.00
250EX09		Expandable 120A Controller - R-net	NC
250EX10		Harness for 120A Expandable Controller	NC
R-net Seating Module Control			\$
250E039		UP Actuator control module, R-Net, CxSM, Thru Drive 2+ Actuators	\$3,088.00
Attendant Controls			\$
2500A05		Dual Control Module for Attendant	\$1,225.00
2500A06		R-net Colour LCD (CJSM1) Control Attendant Only	\$1,225.00
Attendant Mounting Location			\$
250DN1		Right Back Cane	NC
250DN2		Left Back Cane	NC

JOYSTICK / ELECTRONICS

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

Ctrl+5® Switchbox Options

\$

250QG09	Ctrl+5 Switch Box w/ Free Assignable Functions and Switch-It™ Remote Seating™	NC
---------	---	----

Ctrl+5 Mounting Location *(Will swing away if mounted in same location as joystick)*

250EY09	Right Armrest	NC
250EY10	Left Armrest	NC
250EY11	Back Cane, Right	NC
250EY12	Back Cane, Left	NC

Ctrl+5 Switch Type

\$

250EY03	Push Buttons	NC
250EY17	Toggle Switches	\$356.00

R-net Programming Selection

\$

Biometric Repositioning

250PM03	Include Biometric Repositioning (Requires Memory Position #1, #2, #3 Below)	NC
---------	---	----

Memory Seating

Memory Seating feature allows for up to 6 personalized seating positions that are accessed either through the joystick, Ctrl+5, joystick shortcut keys or programmable jacks on Omni2 and joystick. Choose desired presets below. Default positions are dealer adjustable.

250PM06	Include #1 Memory Seating (Seated: Required as return position from all other positions)	NC
250PM07	Include #2 Memory Seating (Supine)	NC
250PM08	Include #3 Memory Seating (Pressure Relief)	NC
250PM09	Include #4 Memory Seating (Forward Tilt)	NC
250PM10	Include #5 Memory Seating (Elevate) (Required)	NC
250PM11	Include #6 Memory Seating (Stand)	NC

Assignable Buttons

250PM12	Include Assignable Buttons	NC
---------	----------------------------	----

PRAR Pressure Relief Alert and Recording

250PM14	Include PRAR (Minimum Requirement Memory Position #1 and #3)	NC
---------	--	----

Power Centermount Footplate Extension Option (Required for UP)

250IE03	Include 7.5" Powered Footplate Extension for Transfer Assist (Electronics programmed to allow footplates to be extended with or without legrest being raised.)	\$592.00
---------	--	----------

SPECIALTY INPUT DEVICES

OPTION	DESCRIPTION	RETAIL
R-net Omni2 Options <i>For Expandable Controls, choose all that apply.</i>		\$
2500N14	R-net Omni2 Only Input Control Module and 3.5" LCD Display w/ Built-In Infrared & Bluetooth, Choose Mount Below	\$2,384.00
2500N15	R-net Omni2 Bundled Input Control Module and 3.5" LCD Display w/ Built-In Infrared & Bluetooth, Bundled Omni2 Available w/ Switch-It Input Devices	NC
Omni2 Mount <i>If used w/ Standard Joystick, it must be mounted opposite side.</i>		\$
250DS12	Link-It with Parallel Swing-Away Omni2 Mount w/ Rear Mounted Interface, Right	\$698.00
250DS13	Link-It with Parallel Swing-Away Omni2 Mount w/ Rear Mounted Interface, Left	\$698.00
Hand Control Packages <i>Requires Standard Joystick Omit / Only one hand control device can be ordered.</i>		\$
Proportional Input Devices		
250EU14	Switch-It MicroGuide	\$4,529.00
250EU17	Switch-It VersaGuide	\$2,557.00
250EU18	Switch-It VersaGuide EZ	\$3,023.00
Mounting Position (Required)		
250EM04	Link-It to the Backrest (Not Available w/ Omit Back Cushion Type)	\$1,338.00
250EM05	Link-It to the Fixed Joystick Mount	\$416.00
250EM06	Link-It to the Parallel Swing-Away Joystick Mount	\$416.00
250EM07	Link-It to the Midline Swing-Away	\$1,480.00
250EM08	Omit Mount	NC
Mounting Platforms (Optional) <i>Not available w/ Link-It to Backrest</i>		
250MG16	SunDisk Hand Tray (Only available w/ MicroSeries Joysticks)	\$416.00
250MG15	Bullet Hand Tray	\$416.00
Mounting Location (Required)		
250IM09	Right	NC
250IM10	Left	NC
Switch Options		
–	See Switch Selections (Maximum 2)	–
Head Control Packages <i>Only one head control device can be ordered. Must select headrest mount on next page.</i>		\$
Proportional Input Devices		
250EU25	Switch-It DualPro® Proportional Head Array - Facial Pads	\$7,057.00
250EU26	Switch-It DualPro Proportional Head Array - Spot Pads	\$7,057.00
250EU51	Switch-It Vigo Head Drive	\$6,293.00
Switched Input Devices		
250EU28	Switch-It 3-Switch Head Array - Short Curved Wings	\$6,154.00
250EU30	Switch-It 3-Switch Head Array - Long Curved Wings	\$6,154.00
Switch Options		
–	See Switch Selections (Maximum 2)	–

SPECIALTY INPUT DEVICES

OPTION	DESCRIPTION	RETAIL
Sip 'N Puff		\$
250ID73	Switch-It Sip 'N Puff Combo Head Array with Link-It Hardware (Requires headrest mount)	\$7,037.00
250ID74	Switch-It OptiStop with Sip 'N Puff	\$2,741.00
Headrest Mounts <i>(Applies to Head/Sip 'N Puff control packages.)</i>		\$
250J13	Whitmyer® Linx Headrest Mount w/ Quick Release	\$247.00
250J14	Whitmyer Linx² Headrest Mount w/ Quick Release	\$298.00
250J15	Whitmyer Axys Headrest Mount w/ Quick Release	\$247.00
Chin Control Packages		\$
Proportional Input Devices		
250EU34	Switch-It MicroGuide (Switches mounted to the headrest only.)	\$4,529.00
250EU36	Switch-It VersaGuide	\$2,557.00
250EU37	Switch-It VersaGuide EZ	\$3,023.00
Mounting Position		
250EM09	Link-It to Chin Harness, Small (150mm / 6")	\$1,760.00
250EM10	Link-It to Chin Harness, Large (220mm / 8.5")	\$1,760.00
250EM11	Link-It to Headrest S/A	\$1,760.00
250EM16	Omit Mount	NC
Input Device Mounting Location		
250IM13	Right	NC
250IM14	Left	NC
Switch Options		
–	See Switch Selections	–
Safety Devices		\$
250SA06	Switch-It Remote Stop Switch - XLR Connector	\$1,317.00
Switched Driving Packages		\$
250SW44	Switch-It® Cool Cube® with DB9 connector for proximity/proportional switch driving <i>(Pro Spot® proportional switches are required if proportional driving is desired.)</i> <i>(For non-proportional driving, single switches can be used.)</i>	\$3,949.00

SPECIALTY INPUT DEVICES

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Switches

Single switches can be programmed to activate any function of the wheelchair including: On/Off, Mode, Profile, Speeds, Seating, Memory Positions, Horn, Lights & Indicators, IR functions, Bluetooth mouse clicks, Environmental Controls when plugged into assignable port.

Switches can connect via: Joysticks (Standard, Micro series, Versa series) Omni2, Cool Cube®, head controls

Switch ordering instructions:

1. Select switch type and quantity in Section "A" using columns below. (Max 6 switches. See limitations below.)
2. Select switch mounting position from Section "B" using corresponding switch column.
3. Select switch mounting location from Section "C" using corresponding switch column.

Section A – Switch Selections

Switch / Quantity

Single Function Switches - Mechanical Switches	Switch 1	Switch 2	Switch 3	Switch 4	Switch 5	Switch 6	\$
2.5" Single Switch - Black							\$230.00
2.5" Single Switch - Blue							\$230.00
2.5" Single Switch - Green							\$230.00
2.5" Single Switch - Red							\$230.00
2.5" Single Switch - Yellow							\$230.00
Egg Switch - Black							\$131.00
Egg Switch - Blue							\$131.00
Egg Switch - Green							\$131.00
Egg Switch - Red							\$131.00
Egg Switch - Yellow							\$131.00
MicroLite <i>(Includes colour labels)</i>							\$195.00
Mini Button 1 - Blue <i>(Velcro)</i>							\$125.00
Mini Button 1 - Green <i>(Velcro)</i>							\$125.00
Mini Button 1 - Red <i>(Velcro)</i>							\$125.00
1" Single Switch - Black							\$326.00
1" Single Switch - Blue							\$326.00
1" Single Switch - Green							\$326.00
1" Single Switch - Orange							\$326.00
1" Single Switch - Red							\$326.00
1" Single Switch - Yellow							\$326.00
Electronic Switches							
Flex Switch (Headrest mount only)							\$254.00
ProSpot - Proportional Switches							
ProSpot Proportional Switch - Red							\$562.00
ProSpot Proportional Switch - Green							\$562.00
ProSpot Proportional Switch - Blue							\$562.00
ProSpot Proportional Switch - Yellow							\$562.00
ProSpot Proportional Switch - Orange							\$562.00
Section B – Switch Mounting Position							
Link-It to the Backrest (Up to 1 per side)							\$770.00
Headrest Mount - Gooseneck & Rigid (Up to 2 per side)							\$770.00
Link-It to the Fixed Joystick Mount (Up to 2 per side)							\$587.00
Link-It to the Parallel Swing-Away Joystick Mount (Up to 2 per side)							\$587.00
Link-It to the Seat Rail (Up to 2 per side) <i>(Not avail w/knee & thigh supports)</i>							\$467.00
Omit Mount (includes Velcro)							NC
Section C – Switch Mounting Location							
Right side of chair							NC
Left side of chair							NC
Joystick Integrated Switches (Pair)							
250J003	–	Proximity Switches on Flex Mount, Pair (Used with VersaGuide/VersaGuide EZ)					\$1,178.00



Q700-UP M with SEDEO PRO ADVANCED

\$

Base Price	
Total of additional equipment	
Packaging and transport	
Total Price	

Signature

