




Q300M Mini

EIPW40

| Customer | | Ship To | |
|-----------|--|----------|-------------|
| Account # | | Name | |
| Date | | Address | |
| PO# | | City | |
| Buyer | | Province | Postal Code |
| Mark For | | Phone | |

Client Measurements

| | |
|----------------------------|----|
| a) Top of head to seat | in |
| b) Top of shoulder to seat | in |
| c) Underarm to seat | in |
| d) Elbow to seat | in |
| e) Top of leg to seat | in |
| f) Shoulder width | in |
| g) Chest width | in |
| i) Hip width | in |
| j) Actual sitting depth | in |
| m) Foot length | in |
| n) Knee to foot | in |
| o) Seat width | in |



Notes

Chair Specifications

| | |
|--|--------------------------------------|
| Weights starting at: (with Batteries) | NPC: 250 lbs SPC: 274 lbs |
| Power Seating | 50° of CG Tilt / 12" seat elevator |
| Seat Type | ASAP |
| Seat Size Range | Width 10" to 16" Depth 10" to 19" |
| Seat-to-Floor Height Range | 17" to 19" |

Mid-Wheel Drive Base Specifications

| | |
|------------------------|--|
| Drive Wheel Position | Mid-wheel drive |
| User Weight Capacity | 165 lbs. |
| Base Size | Width 20.5", Length 33" |
| Minimum Turning Radius | 17" (excludes footplates) |
| Drive Control | R-net PM 90A / VR-2 90A |
| Driving Range | Up to 17 miles |
| Drive Motors | 4-pole 6 mph / 10 kph |
| Drive Wheels | 12" w/flat-free inserts (black) 14" w/flat-free inserts (black) |
| Casters | 7" front / 7" rear (black) |
| Charger | 8A Dual Mode (Compact) |
| Battery Size | 50Ah Group 22 |
| Color | 9 shroud colors and 7 color accents |
| Ground Clearance | 3" |

| BASE | | | |
|--|--|---|------------|
| OPTION | | DESCRIPTION | RETAIL |
| Standard Duty Mid Wheel Drive Base <i>End Item Number EIPW40 (ZIPPIE Seat)</i> | | | \$ |
| 400PC308 | | ZIPPIE Q300 M NPC No Power Option for ZIPPIE Seat Included in Base Pricing : Q300 M Base, VR-2 90amp Electronics for No Power Seat Functions, Six-Wheel Suspension, 12" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 165 lbs., STFH = 17". | \$7,877.00 |
| 400PC309 | | ZIPPIE Q300 M SPC Single Power Option for ZIPPIE Seat Included in Base Pricing : Q300 M Base, VR-2 90amp Electronics for No Power Seat Functions, Six-Wheel Suspension, 12" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 165 lbs., STFH = 17". | \$7,877.00 |
| User Weight | | | \$ |
| 400WG32 | | Up to 165 lbs. | NC |
| Motor Package | | | \$ |
| 400PF99 | | 4-pole 10 kph / 6 mph | NC |
| Tracking | | | \$ |
| 400MI01 | | SureTrac™ Module & Motor Encoders (Improves drive tracking) | \$1,696.00 |
| Drive Wheels | | | \$ |
| 400DW13 | | Drive Wheel, 12", Solid, Knobby, Black | NC |
| 400DW58 | | Drive Wheel, 14", Solid, Knobby, Black | NC |
| Caster Wheels | | | \$ |
| 400CW59 | | Caster Wheels 7", Solid, Black | NC |
| Shroud Color | | | \$ |
| 400C164 | | Black | NC |
| 400C165 | | Candy Apple Red | NC |
| 400C166 | | Candy Blue | NC |
| 400C167 | | Candy Purple | NC |
| 400C176 | | Green Apple | NC |
| 400C178 | | Pearl Pink | NC |
| 400C182 | | Silver | NC |
| 400C125 | | White | NC |
| Color Accents | | | \$ |
| Note: Accent colors and shroud colors may not match | | | |
| 400EZ20 | | Apple Green | NC |
| 400EZ22 | | Black | NC |
| 400EZ21 | | Blue | NC |
| 400EZ23 | | Orange | NC |
| 400EZ24 | | Red | NC |
| 400EZ19 | | White | NC |
| 400EZ25 | | Silver | NC |

| BASE | | | | | | | | | |
|--|--|--|-----------------|------|---------|-----------------|------|---------|----------|
| OPTION | | DESCRIPTION | | | | | | | RETAIL |
| ISO Transit Tie Downs | | | | | | | | | \$ |
| 400T050 | | ISO Transit Package | | | | | | | NC |
| LED Lights | | | | | | | | | \$ |
| 400LI08 | | LED Fender Lights <i>(Joystick Activation)</i> | | | | | | | \$464.00 |
| 400LI02 | | Omit LED Lights | | | | | | | NC |
| Seat Heights | | | | | | | | | \$ |
| Seat Height <i>Measured from front of seat pan to floor with 0° frame angle. Not all seat heights are available with all seat depths and power seating. Please see online configurator for full availability.</i> | | | | | | | | | |
| | | | 12" Drive Wheel | | | 14" Drive Wheel | | | |
| | | | No Power | Tilt | Elevate | No Power | Tilt | Elevate | |
| 400HT55 | | 17" | | N/A | N/A | N/A | N/A | N/A | NC |
| 400HT104 | | 17.5" | N/A | | N/A | | N/A | N/A | NC |
| 400HT56 | | 18" | | N/A | | | | | NC |
| 400HT105 | | 18.5" | | | N/A | N/A | | N/A | NC |
| 400HT57 | | 19" | | | | | | | NC |
| Batteries and Chargers | | | | | | | | | \$ |
| Group 22 Battery | | | | | | | | | |
| 400AG35 | | GP22 / 50Ah Flat Top Pair | | | | | | | \$678.00 |
| Omit Batteries | | | | | | | | | |
| 400AG24 | | Battery Harnesses Only <i>(Omit Batteries)</i> | | | | | | | NC |
| Chargers | | | | | | | | | |
| 400BJ30 | | 10A Dual Mode (Compact) | | | | | | | \$162.00 |
| 400BJ19 | | 8A Dual Mode (Compact) | | | | | | | NC |

POWERED SEATING

| OPTION | DESCRIPTION | RETAIL |
|--------|-------------|--------|
|--------|-------------|--------|

Power Seating Options \$**Seat Elevator**

| | | |
|---------|--|------------|
| 400PI07 | Seat Elevator, Powered, 300 mm / 12" | \$3,477.00 |
| 400PI08 | C-Me® Seat Elevator, Powered, 300 mm / 12" <i>(C-Me speed up to 3 mph) (N/A with VR-2)</i> | \$3,477.00 |

Seat Tilt

| | | |
|---------|--------------------|------------|
| 400PA48 | Posterior 50° Tilt | \$6,496.00 |
|---------|--------------------|------------|

Pre Set Seat Frame Angle (No Power Seating Options, Increases Seat to Floor Height) \$

| | | |
|---------|-----------------------------|----|
| 400FA10 | 0° <i>(Adds 0" to STFH)</i> | NC |
| 400FA3 | 3° <i>(Adds 1" to STFH)</i> | NC |

Backrest Frame Type \$**Backrest Type**

| | | |
|---------|--------------------|----------|
| 400BF28 | Angle Adjustable | NC |
| 400BF29 | 40° Manual Recline | \$617.00 |

Back Canes

| | | |
|---------|-------------------------------|----|
| 400B177 | Medium, 17.5" | NC |
| 400B178 | Tall, 19.5" | NC |
| 400B179 | Height Adjustable, 20" to 26" | NC |

Backrest Frame Option \$

| | | |
|---------|---|----------|
| 400BF13 | Flip Down, Quick Release <i>(Not available with manual recline)</i> | \$179.00 |
|---------|---|----------|

Push Handles \$

| | | |
|---------|-----------------------|---------|
| 400PH28 | Trailing Push Handles | NC |
| 400PH29 | Stroller Handle | \$67.00 |

Towel Bar \$

| | | |
|---------|--|----|
| 400QT03 | Shallow | NC |
| 400QT04 | Deep <i>(Extends 4" from back canes)</i> | NC |

Back Cushion Type \$

| | | |
|----------|--|----|
| 400BT95 | SPOT Back Cushion <i>(Choose Seating in Section A)</i> | NC |
| 400BT135 | JAY® J3® <i>(Choose Seating in Section B)</i> | NC |
| 400BT133 | OMIT <i>(For personalized Mounting. No Headrest Mount)</i> | NC |

A - SPOT Backrest Dimensions \$**Installed Back Height - (From Pan to Top of Back)**

| | | |
|---------|-----|----|
| 400JR46 | 12" | NC |
| 400JR47 | 14" | NC |
| 400JR48 | 16" | NC |
| 400JR71 | 18" | NC |
| 400JR73 | 20" | NC |

POWERED SEATING

| OPTION | DESCRIPTION | RETAIL |
|--------|-------------|--------|
|--------|-------------|--------|

B - JAY® J3® Back**B - JAY® J3® Back - Contour** \$

| | | |
|---------|---|------------|
| 400JP01 | J3® PA Back for 12"-16" Wide | \$912.00 |
| 400JP11 | J3® PL Back for 14"-16" Wide | \$986.00 |
| 400JP13 | J3® PD Back for 14"-16" Wide <i>(18" Shell Height Only)</i> | \$1,068.00 |
| 400JP39 | JAY Zip Back for 10"-14" Wide | \$1,015.00 |

B - JAY® J3® Back Shell Options and Height Ranges \$

JAY® J3® shells are mounted with a minimim 4" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field within the ranges shown.

JAY® J3® Shell Height

| | | |
|---------|--|----|
| 400JR81 | 8" Tall Shell <i>(JAY Zip only)</i> | NC |
| 400JR82 | 10" Tall Shell <i>(JAY Zip only)</i> | NC |
| 400JR83 | 12" Tall Shell <i>(J3 Mid Thoracic Short or JAY Zip)</i> | NC |
| 400JR84 | 13.5" Tall Shell <i>(Mid Thoracic Med)</i> | NC |
| 400JR85 | 15" Tall Shell <i>(Mid Thoracic Tall)</i> | NC |
| 400JR14 | 16.5" Tall Shell <i>(Upper Thoracic Short)</i> | NC |
| 400JR15 | 18" Tall Shell <i>(Upper Thoracic Med)</i> | NC |

JAY® J3® Installed Back Height (From Pan to Top of Back)

| | | |
|---------|-----|----|
| 400JR86 | 12" | NC |
| 400JR87 | 14" | NC |
| 400JR88 | 16" | NC |
| 400JR89 | 18" | NC |
| 400JR90 | 20" | NC |
| 400JR19 | 22" | NC |

B - JAY® Anterior Thoracic Support \$

Classic-Structure Pull Down Harness (4 point strap, Only available with JAY® backrest)

| | | |
|--------|-----------------------------|----------|
| 400K15 | Extra-Small <i>(JZ8300)</i> | \$322.00 |
| 400K16 | Small <i>(JZ8301)</i> | \$322.00 |
| 400K17 | Medium <i>(JZ8302)</i> | \$322.00 |

POWERED SEATING

| OPTION | DESCRIPTION | RETAIL |
|--------|-------------|--------|
|--------|-------------|--------|

JAY® Lateral Thoracic Swing-Away Supports (Requires Pad, Upholstery and Mount) \$**JAY® Lateral Thoracic Support Pads**

| | | |
|---------|--|---------|
| 400JW35 | Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right | \$73.00 |
| 400JW38 | Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left | \$73.00 |
| 400JW36 | Lateral Thoracic (L x H) 100 mm x 152 mm / 4" x 6" Pad - Right | \$73.00 |
| 400JW39 | Lateral Thoracic (L x H) 100 mm x 52 mm / 4" x 6" Pad - Left | \$73.00 |
| 400JW37 | Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right | \$73.00 |
| 400JW40 | Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left | \$73.00 |

JAY® Lateral Thoracic Support Pad Upholstery

| | | |
|---------|-----------------------------|----|
| 400NU01 | Syncra™ Premium Leatherette | NC |
| 400NU02 | Syncra™ Stretch | NC |

JAY® Lateral Thoracic Swing-Away Mount

| | | |
|---------|---|----------|
| 400JW41 | Lateral Thoracic Swing Away Bracket - Right | \$325.00 |
| 400JW42 | Lateral Thoracic Swing Away Bracket - Left | \$325.00 |

JAY® Syncra™ Fixed Lateral Support Accessory \$**JAY® Fixed Lateral Support Accessory Hardware** (*JAY Fixed Lateral Support is height adjustable.*)

| | | |
|---------|--|----------|
| 400JW33 | Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm / 3.5" x 8" - Right | \$190.00 |
| 400JW34 | Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm / 3.5" x 8" - Left | \$190.00 |

Headrests (Requires Pad, Upholstery, & Mount) \$**Headrest Pad**

| | | |
|---------|--|----------|
| 400HE38 | Whitmyer® Plush Headrest Pad 203 mm / 8" | \$377.00 |
| 400HE39 | Whitmyer® Plush Headrest Pad 254 mm / 10" | \$377.00 |
| 400HE45 | Whitmyer® Adjust-A-Plush Headrest Pad Small (<i>Narrow</i>) (<i>Adjustable 6" to 14" Wide x 3" Tall</i>) | \$445.00 |
| 400HE46 | Whitmyer® Adjust-A-Plush Headrest Pad Medium (<i>Adjustable 7" to 16" Wide x 4" Tall</i>) | \$445.00 |
| 400HE48 | OMIT Headrest Pad | NC |

Headrest Mount

| | | |
|---------|--|----------|
| 400MH33 | Whitmyer® Linx Headrest Mount with Quick Release | \$247.00 |
| 400MH34 | Whitmyer® Linx-2 Headrest Mount with Quick Release | \$298.00 |
| 400MH35 | Whitmyer® AXYS Headrest Mount with Quick Release | \$247.00 |
| 400MH17 | OMIT Headrest Mount | NC |

Seat Pan \$

| | | |
|---------|---|----|
| 400OP08 | Small (<i>10" - 13" wide and 10" - 16" deep</i>) | NC |
| 400OP09 | Standard (<i>12" - 16" wide and 12" - 19" deep</i>) | NC |

Seat Cushion Type \$

| | | |
|---------|----------------------|----|
| 400ST80 | JAY® (See Section C) | NC |
| 400ST75 | Omit | NC |

POWERED SEATING

| OPTION | DESCRIPTION | RETAIL |
|--------|-------------|--------|
|--------|-------------|--------|

Seat Size

\$

Seat Pan Width

| | | |
|--------|-----|----|
| 400R83 | 10" | NC |
| 400R84 | 11" | NC |
| 400R85 | 12" | NC |
| 400R86 | 13" | NC |
| 400R87 | 14" | NC |
| 400R78 | 15" | NC |
| 400R53 | 16" | NC |

Seat Pan Depth *(Front of seat pan to back structure)*

| | | |
|---------|-----|----|
| 400RF75 | 10" | NC |
| 400RF76 | 11" | NC |
| 400RF77 | 12" | NC |
| 400RF78 | 13" | NC |
| 400RF52 | 14" | NC |
| 400RF53 | 15" | NC |
| 400RF51 | 16" | NC |
| 400RF45 | 17" | NC |
| 400RF46 | 18" | NC |
| 400RF47 | 19" | NC |

C - JAY® Seat Cushion

\$

| | | |
|----------|--|------------|
| 400CU110 | JAY Zip® <i>(Sizes 10" to 16" Wide)</i> | \$711.00 |
| 400CU59 | JAY Ion® <i>(Sizes 14" to 16" Wide)</i> | \$505.00 |
| 400CU61 | JAY Union® <i>(Sizes 14" to 16" Wide)</i> | \$764.00 |
| 400CU64 | JAY Fusion® <i>(Sizes 14" to 16" Wide)</i> | \$1,045.00 |

POSITIONING SUPPORTS & ACCESSORIES

| OPTION | DESCRIPTION | RETAIL |
|--|--|-----------|
| Knee and Thigh Swing-Away Support Pads | | \$ |
| Lateral Knee Support Pads and Hardware (must pick pad and mount) | | |
| 400XY07 | Lateral Knee Support Pad, Right (3" round) | \$247.00 |
| 400XY08 | Lateral Knee Support Pad, Left (3" round) | \$247.00 |
| 400XY09 | Lateral Knee Support Pad, Right (4" round) | \$247.00 |
| 400XY10 | Lateral Knee Support Pad, Left (4" round) | \$247.00 |
| 400XX10 | Lateral Swing-Away Bracket Knee - Right | \$325.00 |
| 400XX11 | Lateral Swing-Away Bracket Knee - Left | \$325.00 |
| Lateral Thigh Support Pads and Hardware (must pick pad and mount) | | |
| 400KF22 | Lateral Thigh Support Pad Short, Right (5" H x 5" L) | \$146.00 |
| 400KF23 | Lateral Thigh Support Pad Short, Left (5" H x 5" L) | \$146.00 |
| 400KF24 | Lateral Thigh Support Pad Long, Right (5" H x 7" L) | \$146.00 |
| 400KF25 | Lateral Thigh Support Pad Long, Left (5" H x 7" L) | \$146.00 |
| 400KF26 | Lateral Thigh Removable Mount - Right | \$325.00 |
| 400KF27 | Lateral Thigh Removable Mount - Left | \$325.00 |
| Armrest Structure | | \$ |
| 400EG44 | Standard Cantilever Armrest (Height adjustable 9" to 12") | NC |
| 400EG45 | Low Cantilever Armrest (Height adjustable 7" to 10") (Available w/12" armpads only) | NC |
| 400EG46 | Dual Post Short, Removeable, Flip-Back (Height adjustable 10" to 15") | \$434.00 |
| 400EG47 | Single Post, Removable (N/A with thigh support) (Height adjustable 11" to 15") | \$333.00 |
| Arm Pads | | \$ |
| SEDEO® Armrest Pad (Available with cantilever armrest) | | |
| 400AB43 | 300 mm / 12", Integral Foam, Right | NC |
| 400AE43 | 300 mm / 12", Integral Foam, Left | NC |
| 400AB44 | 400 mm / 16", Integral Foam, Right | NC |
| 400AE44 | 400 mm / 16", Integral Foam, Left | NC |
| Waterfall Armpad (Available with Dual Post and Single Post Armrest) | | |
| 400AB01 | 10" Waterfall Armpad, Right | NC |
| 400AE01 | 10" Waterfall Armpad, Left | NC |
| 400AB02 | 14" Waterfall Armpad, Right | NC |
| 400AE02 | 14" Waterfall Armpad, Left | NC |
| Side Guards (Pair) | | \$ |
| 400SG46 | Side Guards Mounted On Seat Rail Bracket (Only with Cantilever Armrests) (N/A with Thigh Supports) | \$356.00 |

POSITIONING SUPPORTS & ACCESSORIES

| OPTION | DESCRIPTION | RETAIL |
|--|---|------------|
| Centermount Legrest | | \$ |
| 400H34 | Fixed Centermount Legrest, Short <i>(Includes Calf Pads and Single Footplate) (12.5" to 15.5" Leg Length)</i> | NC |
| 400H35 | Fixed Centermount Legrest, Long <i>(Includes Calf Pads and Single Footplate) (14.5" to 19.5" Leg Length)</i> | NC |
| 400H123 | Angle Adjustable Centermount Legrest, Small <i>(Includes Calf Pads) (3" to 7" Leg Length)</i> | \$438.00 |
| 400H124 | Angle Adjustable Centermount Legrest, Medium <i>(Includes Calf Pads) (3" to 13" Leg Length)</i> | \$438.00 |
| 400H36 | Power Centermount Legrest, 90° to 0°, 7.5" Extension <i>(Power Articulation) (3" to 17.5" Leg Length)</i> | \$3,435.00 |
| Angle Adjustable Centermount Footplates | | |
| 400F32 | Dual Angle Adjustable Footplates, Small <i>(5.75" x 6")</i> | \$235.00 |
| 400F33 | Dual Angle Adjustable Footplates, Large <i>(5.75" x 10")</i> | \$235.00 |
| Fixed Centermount Footplates | | |
| 400F45 | Single/Dual Footplate | NC |
| Power Centermount Footplates | | |
| 400F127 | Dual, Height, Width and Angle Adjustable, Standard <i>(14"W x 10"D)</i> | \$284.00 |
| 400F128 | Dual, Height, Width and Angle Adjustable, Large <i>(16"W x 12"D)</i> | \$284.00 |
| Footplate Accessories (Available with Power Centermounts) | | |
| 400DI68 | Extra Large Centermount Footplate <i>(Single Footplate Overlay 16"W x 12"D)</i> | \$315.00 |
| 400DI69 | Footplate Connector Link <i>(Attaches footplates together)</i> | \$96.00 |
| 400DI57 | Gel Footpads <i>(Pair, LH and RH)</i> | \$401.00 |
| Power Centermount Calf Pads | | |
| 400DI34 | Calf Pads, Pair, Standard Integral Foam | NC |
| Swing-Away Legrest (Pair) | | \$ |
| 400H196 | 70° Swing-Away <i>(13.5" to 16" Leg Length)</i> | NC |
| 400H197 | 80° Swing-In/Out <i>(Front Mount) (5.5" to 9" Leg Length)</i> | NC |
| 400H198 | 90° Swing-In/Out <i>(Front Mount) (5.5" to 9" Leg Length)</i> | NC |
| 400H51 | Manual Elevating Legrests, Swing-Away <i>(12.5" to 17" Leg Length)</i> | \$406.00 |
| 400HI76 | Omit Legrests | NC |
| Swing-Away Legrest Footplates | | |
| 400F115 | Small Angle Adjustable Flip-Up Footplates <i>(Requires Leg Strap)</i> | \$235.00 |
| 400F116 | Standard Angle Adjustable Flip-Up Footplates <i>(Requires Heel Loop)</i> | \$235.00 |
| Swing-Away Legrest Options | | |
| 400DI09 | Leg Strap | NC |
| 400DI10 | Heel Loops | NC |
| 400DI53 | 2" Hanger Extension <i>(7.5" to 11" Leg Length) (Available w/Swing-In/Out)</i> | \$97.00 |
| 400DI54 | High Mount Extension <i>(7" to 14" Leg Length) (Available w/Swing-Away)</i> | \$97.00 |
| Pelvic Positioning Belts | | |
| 400A837 | Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2" | NC |
| 400A839 | Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2" | NC |
| 400A840 | Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2" | \$103.00 |
| 400A845 | Bodypoint – Evoflex Pelvic Stabilizer Belt | \$216.00 |
| 400A846 | Bodypoint – Evoflex Pelvic Stabilizer Belt with Extension | \$271.00 |

POSITIONING SUPPORTS & ACCESSORIES

| OPTION | DESCRIPTION | RETAIL |
|------------------------------|--|-----------|
| Accessories | | \$ |
| Ventilator Option | | |
| 400VT12 | Vent Tray Fixed <i>(Pulmonetic LTV Models)</i> | \$609.00 |
| O2 Tank Option | | |
| 400A718 | O2 Tank Holder Mounted to Backrest <i>(D and E Styles)</i> | \$391.00 |
| Wheelchair Tray Table | | |
| 4007A09 | Tray, Swing-Away | \$381.00 |
| Backpack | | |
| 400A1 | Black Utility | \$112.00 |
| Misc. Accessories | | |
| 4008A59 | Transfer Bar <i>(Pair)</i> | \$493.00 |
| 4008A52 | Ram Mount Phone | \$77.00 |
| 4008A53 | Ram Mount Tablet | \$308.00 |
| 4008A54 | Ram Mount Cup Holder | \$77.00 |
| 4008A58 | Rear View Mirror | \$77.00 |
| 400A610 | Walker Holder | \$191.00 |
| 400A608 | Cane/Crutch Holder | \$191.00 |
| 4009A10 | Utility Hooks | \$58.00 |

JOYSTICK / ELECTRONICS

| | | |
|----------------------------------|--|------------|
| Electronics Package | | \$ |
| 400EL01 | VR-2 or R-net Electronics | NC |
| VR-2 Joystick | | \$ |
| 400J127 | VR-2 Joystick | NC |
| R-net Joystick (Remote) | | \$ |
| 400J121 | Standard LED Joystick <i>(w/o Screen Display)</i> | NC |
| 400J116 | Color Joystick <i>(2.3" Color LCD Screen Display, CJSM1)</i> | NC |
| 400J118 | High Resolution Color Joystick, (CJSM2+BT) <i>3.5" Color LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control</i> | \$1,162.00 |
| 400J119 | Omit Remote Control <i>(Required with arm mounted specialty controls)</i> | NC |
| Joystick Mounting | | \$ |
| 400JM34 | Fixed on Armrest | NC |
| 400JM35 | Parallel Swing-Away | \$447.00 |
| 400JM36 | Link-It Mid Line Mount on Armrest <i>(Swing-Away)</i> | \$447.00 |
| 400JM49 | Attendant <i>(Available with LED or CJSM1 only. Mounted to Push Handle)</i> | NC |
| Joystick Mount (Location) | | \$ |
| 400HM16 | Right Hand | NC |
| 400HM17 | Left Hand | NC |

JOYSTICK / ELECTRONICS

| OPTION | DESCRIPTION | RETAIL |
|---|---|------------|
| Display Language Selection R-net \$ | | |
| 400L1 | English | NC |
| 400L3 | French | NC |
| 400L7 | Spanish | NC |
| Joystick Handle (Optional) <i>N/A with HD Joystick</i> \$ | | |
| 400JH13 | Ball Handle | \$110.00 |
| 400JH15 | T-Handle | \$110.00 |
| 400JH24 | Bodypoint – Rubber Dome Joystick Handle | \$216.00 |
| 400JH25 | Bodypoint – U-Shaped Joystick Handle (4") | \$216.00 |
| 400JH26 | Bodypoint – Grooved Mushroom Joystick Handle | \$216.00 |
| Electronic Programmers \$ | | |
| 400PR13 | R-net PC Wireless Programmer | \$660.00 |
| R-net Power Modules \$ | | |
| 400ER12 | Power Module, VR-2 90A, Non-Expandable | NC |
| 400ER22 | Power Module, R-net, 90A, Non-Expandable (<i>Only available with Single Power</i>) | NC |
| 400ER13 | Power Module, R-net, 90A, Expandable (<i>Required for Specialty Controls or 3+ Actuators</i>) | \$102.00 |
| 400EX11 | Expandable 90A Controller - R-net | NC |
| 400EX12 | Harness for 90A Expandable Controller | NC |
| Seating Control \$ | | |
| Powered Seating Control | | |
| 400EO40 | Omit Thru-Drive (<i>Requires Switch Control</i>) | NC |
| 400EO37 | Thru-drive Control for 1 Actuator | \$1,443.00 |
| 400EO38 | Thru-drive Control for 2+ Actuators | \$2,903.00 |
| Switch Control (<i>Available with Omit Thru-Drive only</i>) (<i>1 Actuator for Tilt or Elevate only</i>) | | |
| 400EO33 | Dual Toggle for 1 Actuator | NC |
| 400EO34 | 2.5" Single Button Switch for 1 Actuator (<i>Button switches are not mounted</i>) | NC |
| 400EO41 | Ctrl+5 Switch Box with Voice Activation for 2 Actuators or more | NC |
| Additional Power Seating Control (<i>Available with Thru-Drive only</i>) | | |
| 400QG11 | Ctrl+5 Switch Box with Voice Activation | \$1078.00 |
| 400QG12 | Ctrl+5 Switch Box with Voice Activation (<i>Multi-power only</i>) | NC |
| Dual Toggle Mounting Location | | |
| 400TN1 | Toggle Mount, Right | NC |
| 400TN2 | Toggle Mount, Left | NC |
| Ctrl+5 Mounting location (<i>If mounted on the same side as the joystick, it will swing away.</i>) | | |
| 400EY19 | Right Armrest | NC |
| 400EY20 | Left Armrest | NC |
| 400EY21 | Back Cane, Right | NC |
| 400EY24 | Back Cane, Left | NC |
| Ctrl+5 Switch Type | | |
| 400EY22 | Push Buttons | NC |
| 400EY23 | Toggle Switches | \$356.00 |

JOYSTICK / ELECTRONICS

| OPTION | DESCRIPTION | RETAIL |
|--------|-------------|--------|
|--------|-------------|--------|

Attendant Controls \$

| | | |
|---------|---|------------|
| 4000A05 | Dual Control Module for Attendant | \$1,225.00 |
| 4000A06 | R-net Color LCD (CJSM1), Control Attendant Only | \$1,225.00 |

Attendant Mounting Location \$

| | | |
|--------|-----------------|----|
| 400DN5 | Right Back Cane | NC |
| 400DN6 | Left Back Cane | NC |

R-net Programming \$**Assignable Buttons**

Assignable Button feature allows for personalization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

| | | |
|---------|----------------------------|----|
| 400PM12 | Include Assignable Buttons | NC |
|---------|----------------------------|----|

Power Centermount Footplate Extension Option

| | | |
|---------|--|----------|
| 400IE03 | Include 7.5" Powered Footplate Extension for Transfer Assist (<i>Electronics programmed to allow footplates to be extended with or without legrest being raised.</i>) Only Available with Powered Articulation | \$592.00 |
|---------|--|----------|

R-net Omni2 Options *For Expandable Controls, choose all that apply.* \$

| | | |
|---------|---|------------|
| 4000N14 | R-net Omni2 Only Input Control Module and 3.5" LCD Display w/ Built-In Infrared and Bluetooth® Choose Mount Below | \$2,384.00 |
| 4000N15 | R-net Omni2 Bundled Input Control Module and 3.5" LCD Display w/ Built-In Infrared and Bluetooth® Bundled Omni2 Available w/ Switch-It Input Devices | NC |

Omni2 Mount *If used w/ Standard Joystick, it must be mounted opposite side.* \$

| | | |
|---------|--|----------|
| 400DS14 | Link-It with Parallel Swing-Away Omni2 Mount w/Rear Mounted Interface, Right | \$698.00 |
| 400DS15 | Link-It with Parallel Swing-Away Omni2 Mount w/Rear Mounted Interface, Left | \$698.00 |

SPECIALTY INPUT DEVICES

| OPTION | DESCRIPTION | RETAIL |
|--------|-------------|--------|
|--------|-------------|--------|

Hand Control Packages *(Requires Standard Joystick Omit / Only one hand control device can be ordered.)*

\$

Proportional Input Devices

| | | |
|---------|-------------------------|------------|
| 400EU14 | Switch-It MicroGuide | \$4,574.00 |
| 400EU17 | Switch-It VersaGuide | \$2,557.00 |
| 400EU18 | Switch-It VersaGuide EZ | \$3,023.00 |

Mounting Position (Required)

| | | |
|---------|---|------------|
| 400EM05 | Link-It to the Fixed Joystick Mount | \$416.00 |
| 400EM06 | Link-It to the Parallel Swing-Away Joystick Mount | \$416.00 |
| 400EM07 | Link-It to the Midline Swing-Away | \$1,480.00 |
| 400EM08 | Omit Mount | NC |

Mounting Platforms (Optional) *Not available w/ Link-It to Backrest*

| | | |
|---------|--|----------|
| 400MG16 | SunDisk Hand Tray <i>(Only available w/ MicroSeries Joysticks)</i> | \$416.00 |
| 400MG15 | Bullet Hand Tray | \$416.00 |

Mounting Location (Required)

| | | |
|---------|-------|----|
| 400IM09 | Right | NC |
| 400IM10 | Left | NC |

Switch Options

| | | |
|---|--|---|
| – | See Switch Selections <i>(Maximum 2)</i> | – |
|---|--|---|

Head Control Packages *(Select headrest mount on next page)*

\$

*Only one head control device can be ordered.***Proportional Input Devices**

| | | |
|---------|--|------------|
| 400EU25 | Switch-It DualPro® Proportional Head Array - Facial Pads | \$7,057.00 |
| 400EU26 | Switch-It DualPro® Proportional Head Array - Spot Pads | \$7,057.00 |
| 400EU51 | Switch-It Vigo Head Drive | \$6,293.00 |

Switched Input Devices

| | | |
|---------|--|------------|
| 400EU28 | Switch-It 3-Switch Head Array - Short Curved Wings | \$6,154.00 |
| 400EU30 | Switch-It 3-Switch Head Array - Long Curved Wings | \$6,154.00 |

Switch Options

| | | |
|---|--|---|
| – | See Switch Selections <i>(Maximum 2)</i> | – |
|---|--|---|

Sip 'N Puff

\$

| | | |
|---------|---|------------|
| 400ID73 | Switch-It Sip 'N Puff Combo Head Array with Link-It Hardware <i>(Requires headrest mount)</i> | \$7,037.00 |
| 400ID74 | Switch-It OptiStop with Sip 'N Puff | \$2,741.00 |

SPECIALTY INPUT DEVICES

| OPTION | DESCRIPTION | RETAIL |
|--|---|------------|
| Headrest Mounts <i>(Applies to Head/Sip 'N Puff control packages)</i> | | \$ |
| 400IJ13 | Whitmyer® Linx Headrest Mount with Quick Release | \$247.00 |
| 400IJ14 | Whitmyer® Linx-2 Headrest Mount with Quick Release | \$298.00 |
| 400IJ15 | Whitmyer® AXYS Headrest Mount with Quick Release | \$247.00 |
| Chin Control Packages | | \$ |
| Proportional Input Devices | | |
| 400EU34 | Switch-It MicroGuide <i>(Switches mounted to the headrest only.)</i> | \$4,574.00 |
| 400EU36 | Switch-It VersaGuide | \$2,557.00 |
| 400EU37 | Switch-It VersaGuide EZ | \$3,023.00 |
| Mounting Position | | |
| 400EM09 | Link-It To Chin Harness, Small (150mm / 6") | \$1,760.00 |
| 400EM10 | Link-It To Chin Harness, Large (220mm / 8") | \$1,760.00 |
| 400EM11 | Link-It To Headrest, Swing-Away | \$1,760.00 |
| 400EM16 | Omit Mount | NC |
| Input Device Mounting Location | | |
| 400IM13 | Right | NC |
| 400IM14 | Left | NC |
| Switch Options | | |
| – | See Switch Selections <i>(Minimum 1 for the Power Chin Boom, Maximum 3)</i> | – |
| Safety Devices | | \$ |
| 400SA06 | Switch-It Remote Stop Switch - XLR Connector | \$1,317.00 |
| Switched Driving Packages | | \$ |
| 400SW44 | Switch-It™ Cool Cube with DB9 connector for proximity/proportional switch driving. <i>(Pro Spot proportional switches are required if proportional driving is desired.) (For non-proportional driving, single switches can be used.)</i> | \$3,792.00 |

SPECIALTY INPUT DEVICES

| OPTION | HCPC | DESCRIPTION | RETAIL |
|--------|------|-------------|--------|
|--------|------|-------------|--------|

Switches

Single switches can be programmed to activate any function of the wheelchair including: On/Off, Mode, Profile, Speeds, Seating, Memory Positions, Horn, Lights & Indicators, IR functions, Bluetooth mouse clicks, Environmental Controls when plugged into assignable port.

Switches can connect via: Joysticks (Standard, Micro series, Versa series) Omni2, SJX6 Hub, head controls.

Switch ordering instructions:

1. Select switch type and quantity in Section "A" using columns below. (Max 2 switches. See limitations below.)
2. Select switch mounting position from Section "B" using corresponding switch column.
3. Select switch mounting location from Section "C" using corresponding switch column.

Section A – Switch Selections

Switch / Quantity

| Single Function Switches - Mechanical Switches | Switch 1 | \$ | Switch 2 | \$ |
|--|----------|--|----------|------------|
| 2.5" Single Switch - Black | 400IR31 | \$110.00 | 400IR56 | \$110.00 |
| 2.5" Single Switch - Blue | 400IR32 | \$110.00 | 400IR57 | \$110.00 |
| 2.5" Single Switch - Green | 400IR33 | \$110.00 | 400IR58 | \$110.00 |
| 2.5" Single Switch - Red | 400IR34 | \$110.00 | 400IR59 | \$110.00 |
| 2.5" Single Switch - Yellow | 400IR35 | \$110.00 | 400IR60 | \$110.00 |
| Egg Switch - Black | 400IR36 | \$110.00 | 400IR61 | \$110.00 |
| Egg Switch - Blue | 400IR37 | \$110.00 | 400IR62 | \$110.00 |
| Egg Switch - Green | 400IR38 | \$110.00 | 400IR63 | \$110.00 |
| Egg Switch - Red | 400IR39 | \$110.00 | 400IR64 | \$110.00 |
| Egg Switch - Yellow | 400IR40 | \$110.00 | 400IR65 | \$110.00 |
| MicroLite <i>(Includes color labels)</i> | 400IR41 | \$115.00 | 400IR66 | \$115.00 |
| Mini Button 1 - Blue <i>(Velcro)</i> | 400IR42 | \$308.00 | 400IR67 | \$308.00 |
| Mini Button 1 - Green <i>(Velcro)</i> | 400IR43 | \$308.00 | 400IR68 | \$308.00 |
| Mini Button 1 - Red <i>(Velcro)</i> | 400IR44 | \$308.00 | 400IR69 | \$308.00 |
| 1" Single Switch - Black | 400IR45 | \$388.00 | 400IR70 | \$388.00 |
| 1" Single Switch - Blue | 400IR46 | \$388.00 | 400IR71 | \$388.00 |
| 1" Single Switch - Green | 400IR47 | \$388.00 | 400IR72 | \$388.00 |
| 1" Single Switch - Orange | 400IR48 | \$388.00 | 400IR73 | \$388.00 |
| 1" Single Switch - Red | 400IR49 | \$388.00 | 400IR74 | \$388.00 |
| 1" Single Switch - Yellow | 400IR50 | \$388.00 | 400IR75 | \$388.00 |
| Electronic Switches | | | | |
| Flex Switch (Headrest mount only) | 400IR54 | \$542.00 | 400IR79 | \$542.00 |
| ProSpot - Proportional Switches | | | | |
| ProSpot Proportional Switch - Red | 400IR181 | \$542.00 | 400IR186 | \$542.00 |
| ProSpot Proportional Switch - Green | 400IR182 | \$542.00 | 400IR187 | \$542.00 |
| ProSpot Proportional Switch - Blue | 400IR183 | \$542.00 | 400IR188 | \$542.00 |
| ProSpot Proportional Switch - Yellow | 400IR184 | \$542.00 | 400IR189 | \$542.00 |
| ProSpot Proportional Switch - Orange | 400IR185 | \$542.00 | 400IR190 | \$542.00 |
| Section B – Switch Mounting Position | | | | |
| Link-It to the Backrest <i>(Up to 1 per side)</i> | 400EM17 | \$1,338.00 | 400EM23 | \$1,338.00 |
| Link-It to the Headrest - Gooseneck & Rigid <i>(Up to 2 per side)</i> | 400EM18 | \$1,760.00 | 400EM24 | \$1,760.00 |
| Link-It to the Fixed Joystick Mount <i>(Up to 2 per side)</i> | 400EM19 | \$416.00 | 400EM25 | \$416.00 |
| Link-It to the Parallel Swing-Away Joystick Mount <i>(Up to 2 per side)</i> | 400EM20 | \$416.00 | 400EM26 | \$416.00 |
| Link-It to the Seat Rail <i>(Up to 2 per side)</i> <i>(Not avail w/knee & thigh supports)</i> | 400EM21 | \$556.00 | 400EM27 | \$556.00 |
| Omit Mount <i>(includes Velcro)</i> | 400EM22 | NC | 400EM28 | NC |
| Section C – Switch Mounting Location | | | | |
| Right side of chair | 400IM17 | \$698.00 | 400IM19 | \$698.00 |
| Left side of chair | 400IM18 | \$698.00 | 400IM20 | \$698.00 |
| Joystick Integrated Switches (Pair) | | | | |
| 400J003 | – | Proximity Switches on Flex Mount, Pair <i>(Used with VersaGuide/VersaGuide EZ)</i> | | \$1,403.00 |



ZIPPIE Q300 M Mini

\$

Base Price

Total of additional equipment

Packaging and transport

Total Price

Signature

