

April 2025



Customer Contact Details

Customer Name: _____

Advisor: _____

Address: _____

Phone: _____

Email: _____

Ship to Customer Address: ☐

Shipping Name: _____

Shipping Address: _____

Phone: _____

Email: _____

Additional Information

Notes: _____

Current Wheelchair: _____

Customer Weight: _____

Customer Height: _____

Customer Disability: _____

Guide:		<input type="checkbox"/> Base product		<input checked="" type="checkbox"/> Required option		<input type="checkbox"/> Optional Extra		Critical Information		Information Needed	
<input checked="" type="checkbox"/> SQCXZ010		CLUB SPORT		Max User Weight 125kg (275lbs)		<input type="checkbox"/> Yes, I am within the max user weight					
\$4,009.00											
Standard Specification Includes: <ul style="list-style-type: none"> Basketball or Tennis (option 3 only) Aluminium screw on sideguards Adjustable position footrest Adjustable centre of gravity (18 degree) Removable single sport anti-tip Adjustable height backrest 250-300mm Cross laced wire spoke wheels Aluminium hard anodized push rim/kenda tire 72mm caster wheels with extruded forks Adjustable seat and backrest upholstery 5cm medium density foam cushion (*option 4) Calf strap & foot strap 											
Key Product Measurements											
Frame Option (Choose one)											
<input type="checkbox"/> XY520 Option 1 (Junior size)			<input type="checkbox"/> XY530 Option 2 (Low Point)			<input type="checkbox"/> XY540 Option 3 (Low/Mid Point)			<input type="checkbox"/> XY550 Option 4 (Max height)		
Sport			Sport			Sport			Sport		
<input checked="" type="checkbox"/> Basketball (with wing)			<input checked="" type="checkbox"/> Basketball (with wing)			<input type="checkbox"/> Basketball (with wing) <input type="checkbox"/> Tennis (no wing) (5mm A/T)			<input checked="" type="checkbox"/> Basketball (with wing)		
Frame height			Frame height			Frame height			Frame height		
<input checked="" type="checkbox"/> Front frame height: 480mm			<input checked="" type="checkbox"/> Front frame height: 500mm			<input checked="" type="checkbox"/> Front frame height: 500mm			<input checked="" type="checkbox"/> Front frame height: 530mm		
<input checked="" type="checkbox"/> Rear frame height: 420mm			<input checked="" type="checkbox"/> Rear frame height: 400mm			<input checked="" type="checkbox"/> Rear frame height: 440mm			<input checked="" type="checkbox"/> Rear frame height: 530mm		
Wheel size			Wheel size			Wheel size			Wheel size		
<input type="checkbox"/> 22" <input type="checkbox"/> 24" <input type="checkbox"/> 25"			<input type="checkbox"/> 24" <input type="checkbox"/> 25"			<input type="checkbox"/> 24" <input type="checkbox"/> 25"			<input type="checkbox"/> 26" <input type="checkbox"/> 700c		
XY415 XY420 XY430			XY420 XY430			XY420 XY430			XY440 XY450		
Seat width (inside of sides)			Seat width (inside of sides)			Seat width (inside of sides)			Seat width (inside of sides)		
<input type="checkbox"/> 340mm <input type="checkbox"/> 370mm			<input type="checkbox"/> 370mm <input type="checkbox"/> 400mm <input type="checkbox"/> 440mm			<input type="checkbox"/> 370mm <input type="checkbox"/> 400mm <input type="checkbox"/> 440mm			<input type="checkbox"/> 370mm <input type="checkbox"/> 400mm <input type="checkbox"/> 440mm		
Seat length			Seat length			Seat length			Seat length		
<input type="checkbox"/> 300mm <input type="checkbox"/> 330mm <input type="checkbox"/> 360mm			<input type="checkbox"/> 360mm <input type="checkbox"/> 390mm <input type="checkbox"/> 430mm			<input type="checkbox"/> 360mm <input type="checkbox"/> 390mm <input type="checkbox"/> 430mm			<input type="checkbox"/> 360mm <input type="checkbox"/> 390mm <input type="checkbox"/> 430mm		
Frame length			Frame length			Frame length			Frame length		
<input checked="" type="checkbox"/> 520mm			<input checked="" type="checkbox"/> 580mm			<input checked="" type="checkbox"/> 580mm			<input checked="" type="checkbox"/> 580mm (SW = 370/400) <input checked="" type="checkbox"/> 630mm SW = 440		
Options											
Finish Options						Accessories					
<input checked="" type="checkbox"/> XN010 Powder coated frame						<input type="checkbox"/> XS020 Sport ratchet strap (click strap) (not fitted) \$140.00					
Specify Code <input type="text"/>						<input type="checkbox"/> XS010 100mm velcro adjustable strap \$86.00					
<input checked="" type="checkbox"/> XN150 Powder coated components (Black)						<input type="checkbox"/> XK080 Spinergy SLX 24 spoke (Wheel upgrade) \$1006.00					
						<input type="checkbox"/> XP170 Individual backrest embroidery (Please provide .jpeg file and positioning) \$240.00					