



# Quickie Xplore2 STD

March 2021

EIPW24



### ACCOUNT INFORMATION

Dealer Contact Name:	_____	P.O. # or Quote:	_____
E-mail Address:	_____		
Phone Number:	_____ Fax Number: _____	User Name or Number:	_____
Bill to (Name):	_____	User Weight:	_____
Bill to (Location):	_____	User Height:	_____
Ship to (Name):	_____	Seat Cushion Thickness:	_____
Address:	_____	Knee to Heel Length: Left: _____ Right: _____	
City:	_____		
Province:	_____ Postal Code: _____		

The order cannot be processed without the following:

- A user name or number, the client's measurements and weight
- Sunrise Medical is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Sunrise Medical summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

### STANDARD OPTIONS

Quickie Xplore 2 STD Power Base, 6-wheel independent wheel suspension, an additional seat suspension system, tie-down hardware for strap systems, 8A onboard charger, 8" front casters, and standard auto style buckle seat positioning strap, 6.5 mph motors.

### QUICKIE Xplore 2 STD POWER WHEELCHAIR

<input type="checkbox"/> XPLORESTD2	QUICKIE Xplore 2 STD Power Wheelchair	\$9,275.00
-------------------------------------	---------------------------------------	------------

### USER WEIGHT PACKAGES

<input type="checkbox"/> WCA3	300 lbs. Maximum user weight capacity	STD
-------------------------------	---------------------------------------	-----

### STANDARD FRAME COLORS

<input type="checkbox"/> FCBLK	Black	STD
<input type="checkbox"/> FCSBM	Matte Black	STD
<input type="checkbox"/> FCGMC	Gun Metal (Charcoal)	STD
<input type="checkbox"/> FCSVM	Silver Metallic	STD
<input type="checkbox"/> FCCBU	Candy Blue	STD
<input type="checkbox"/> FCCRE	Candy Red	STD
<input type="checkbox"/> FCMPL	Midnight Purple	STD
<input type="checkbox"/> FCRZB	Razzeberry	STD

### OPTIONAL FRAME COLORS (May Increase lead times)

<input type="checkbox"/> FCMYG	Military Green (New)	N/C
<input type="checkbox"/> FCBOR	Burnt Orange Metallic	N/C
<input type="checkbox"/> FCBAO	Black Opal Metallic	N/C
<input type="checkbox"/> FCOPM	Sparkle Orange	N/C
<input type="checkbox"/> FCMYT	Military Tan (New)	N/C
<input type="checkbox"/> FCHDD	Desert Digital Camo Hydro Graphics	\$340.00
<input type="checkbox"/> FCHWL	Woodland Camo Hydro Graphics	\$340.00
<input type="checkbox"/> FCHCF	Carbon Fiber Hydro Graphics (Matte Black base)	\$340.00

Includes matching accessories with the Hydro Graphics Option. The matching accessories option is not required with these colours.

### MATCHING ACCESSORIES TO THE FRAME COLOR

<input type="checkbox"/> FC***-M	Matching Accessories to the Frame Color (Wheels, front access door, headrest mount and back shell when selected)	\$256.00
----------------------------------	--	----------

This options is included at no charge with the Hydro Graphics Frame Color option.

### TRANSPORT OPTION

<input type="checkbox"/> TTD4	1 4-Point wheelchair transport brackets	STD
-------------------------------	---	-----

**Note 1** 4-Point wheelchair transport brackets have not yet received WC-19 crash-testing approval. Use these transport brackets only to secure an unoccupied wheelchair during transport.

### MOTOR PACKAGE

<input type="checkbox"/> MP65	6.5 MPH 4-Pole Motor Package	STD
<input type="checkbox"/> MP85	8.5 MPH 4-Pole High Speed Motor Package (Max. user weight 250 lbs - Only available with ECR120)	N/C

### ASSISTIVE DRIVING PACKAGE

<input type="checkbox"/> INTELLIDR	2,3 Intelli-Drive	\$1,694.00
------------------------------------	-------------------	------------

**Note 2** Intelli-Drive technology provides increased tracking capabilities with enhanced traction and steering precision. This is ideal for those with marginal hand function or switch control.

**Note 3** Requires the R-net PM 120 Amp Expandable Controller (ECR120)

**▼ BATTERY CHARGER**

<input type="checkbox"/>	BCOB	8A Offboard Battery Charger	N/C
--------------------------	------	-----------------------------	-----

**▼ BATTERY TRAY AND FRAME SIZE**

<input type="checkbox"/>	FS22B	4	Battery Tray for Group 22NF Batteries (Total frame width = 23 3/8")	STD
<input type="checkbox"/>	FS24B		Battery Tray for Group 24 Batteries (Total frame width = 25 3/8")	N/C
Please Select One Option: Front Cas: <input type="checkbox"/> Flared out (STD) <input type="checkbox"/> Straight (Optional)				
Note 4 The 23 3/8" compact frame is not compatible with most of the Center Mount Front Rigging Options due to Mechanical conflicts				
Note 5 The Straight Front Arm option is not available with any of the Center Mount Front Rigging Options due to Mechanical conflicts				

**▼ BATTERIES**

<input type="checkbox"/>	BAT22NFPR	SigmaTek Gel Batteries-Group 22NF (Pair - On chair)	\$658.00
<input type="checkbox"/>	BAT24PR	SigmaTek Gel Batteries-Group 24 (Pair - On Chair)	\$881.00

**▼ CHOOSE TIRE COLOR (APPLIES TO ALL TIRES OR CASTERS)**

- BLACK (New)  GREY

**▼ DRIVE WHEEL**

<input type="checkbox"/>	DW14PT		14" Cast Wheel Pneumatic Tire	STD
<input type="checkbox"/>	DW14PI		14" Cast Wheel Pneumatic with Flat-Free Inserts	\$256.00
<input type="checkbox"/>	DW144PT	6	14" x 4" Extra Wide Cast Wheel Pneumatic Tire	\$279.00
<input type="checkbox"/>	DW144PI	6	14" x 4" Extra Wide Cast Wheel Pneumatic with Flat-Free Inserts	\$279.00
<input type="checkbox"/>	DW12PT	7	12 1/2" Cast Wheel Pneumatic Tire	N/C
<input type="checkbox"/>	DW12PI	7	12 1/2" Cast Wheel Pneumatic with Flat-Free Inserts	\$256.00
Note 6 The 14" x 4" Extra Wide Wheel increases the overall width of the Base by 2 1/2"				
The option to match the wheels to the shroud colour (FC...M) is not available with the 14" x 4" Extra Wide Wheels.				
The 14" x 4" option is not available with the VR2 electronics (ECV902A)				
Note 7 12 1/2" Drive wheels reduce overall seat to floor height by 3/4", lowers ground clearance and speed, but increases motor torque				

**▼ FRONT CASTERS**

<input type="checkbox"/>	FC8PTR		8" Pneumatic Casters	STD
<input type="checkbox"/>	FC8PIR		8" Pneumatic Casters with Flat-Free Inserts	\$256.00
<input type="checkbox"/>	FC9PT	8	9" Pneumatic Casters	\$173.00
<input type="checkbox"/>	FC9PI	8	9" Pneumatic Casters with Flat-Free Inserts	\$313.00
Note 8 The 9" casters will increase the minimum seat to floor height by 1/2" and are NOT compatible with the 12 1/2" Drive wheel configuration (DW12PT and DW12PI)				
The 9" casters are NOT compatible with many front rigging options. Please see details in the front rigging section.				

**▼ SEAT TO FLOOR HEIGHT (Approximate measurement to the top of seat pan ± 1/4")**

► Options for a Seat Only

<input type="checkbox"/>	STFH1525	9, 10	15 1/4" (Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH18	18"	(Possible Seat angle of 0° and 3°)
<input type="checkbox"/>	STFH16	10	16" (Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH19	19"	(Possible Seat angle of 0°, 3° and 6°)
<input type="checkbox"/>	STFH17	10	17" (Possible Seat angle of 0° and 3°)	<input type="checkbox"/>	STFH20	20"	(Possible Seat angle of 0°, 3°, 6°, 9°)

► Options for a Power Tilt or a STD Power Elevating Seat

<input type="checkbox"/>	STFH1575	9, 10	15 1/4" (Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH18	18"	(Possible Seat angle of 0° only)
<input type="checkbox"/>	STFH165	10	16 1/2" (Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH185	18 1/2"	(Possible Seat angle of 0° only)
<input type="checkbox"/>	STFH175		17 1/2" (Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH19	19"	(Possible Seat angle of 0° only)

► Options for a Combo - Power Tilt with a STD Power Elevating Seat

<input type="checkbox"/>	STFH1775	9	17 3/4" (Possible Seat angle of 0° only)
<input type="checkbox"/>	STFH185		18 1/2" (Possible Seat angle of 0° only)
Note 9 Requires the 12 1/2" Drive wheels (DW12PT or DW12PI). This lowers ground clearance and speed, but increases motor torque.			
Note 10 The Heavy Duty Swing Away footrests (SOSA AF) are not available with this seat to floor height.			

**▼ SEAT TO FLOOR HEIGHT SPECIFICATION (If not standard)**

<input type="checkbox"/>	SOMHSS	Specific seat to floor height requirement other than standard. Upcharges may apply. Required at: _____" (± 1/4" Range)	RFQ
--------------------------	--------	---	-----

**▼ PRE-SET SEAT FRAME ANGLE**

<input type="checkbox"/>	SFE-0	0° Pre-set seat frame seat angle (Must be selected when using power options)	STD
<input type="checkbox"/>	SFE-3	3° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C
<input type="checkbox"/>	SFE-6	6° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C
<input type="checkbox"/>	SFE-9	9° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C

**▼ LIGHTING PACKAGE**

<input type="checkbox"/>	EOLP	11	LED Light Package (Includes front beams, side flashers, rear tail light with flashers & brake light)	\$1,131.00
Note 11 Requires the upgraded Joystick with Light Operation or the OMNI (ECVJ6KL, ECELJL, ECRJCL, ECRJCTPL, ECRJC35B or EDIB). The Light Package is not compatible with the EPSEL option.				

**▼ NON-EXPANDABLE ELECTRONICS - PG DRIVES VR2 CONTROLLER AND JOYSTICKS**

<input type="checkbox"/>	ECV902A	VR2 90 Amp Drive Controller	STD
<input type="checkbox"/>	ECVJ4K	12 VR2 Standard 4-Key Joystick	STD
<input type="checkbox"/>	ECVJ6K	12 VR2 6-Key Joystick (Allows for a single power option to operate through this joystick)	N/C
<input type="checkbox"/>	ECVJ9KL	12 VR2 9-Key Joystick with Lighting Operation (Allows for a single power function to operate through this joystick)	\$168.00
<b>Note 12</b> Requires the VR2 90A Drive Control (ECV902A)			
VR2 Electronics can only accommodate a single power seat function. The only choices are Power Tilt or Elevating seat			

**▼ EXPANDABLE ELECTRONICS - PG DRIVES RNET CONTROLLERS**

<input type="checkbox"/>	ECEL90	R-net EL 90 Amp Expandable Controller	N/C
<input type="checkbox"/>	ECR120	R-net PM 120 Amp Expandable Controller	N/C
<input type="checkbox"/>	ECRNH	13 Harness for Expandable Controller	N/C
<b>Note 13</b> Must be ordered with the expandable controller ECEL90 or ECR120 to interface with Joysticks or OMNI screen			

**▼ EXPANDABLE ELECTRONICS - PG DRIVES RNET JOYSTICKS**

<input type="checkbox"/>	ECO	Omit Joystick (For use with an Omni Display)	N/C
<input type="checkbox"/>	ECELJ	14 R-net (LED) Joystick	STD
<input type="checkbox"/>	ECRJJC	14 R-net Joystick with Color screen	\$168.00
<b>▶ JOYSTICKS WITH LIGHTING OPERATION</b>			
<input type="checkbox"/>	ECELJL	14 R-net (LED) Joystick with Lighting Operation	\$168.00
<input type="checkbox"/>	ECRJCL	14 R-net Joystick with Color screen and Lighting Operation	\$222.00
<b>▶ ADVANCED JOYSTICK</b>			
<input type="checkbox"/>	ECRJ35B	14 R-net Advanced Bluetooth Joystick (3.5" Screen). Lighting Operation, paddle switches, shortcut keys and stereo jacks	\$1,202.00
<b>Note 14</b> This option requires the expandable controller ECEL90 or ECR120			

**▼ JOYSTICK MOUNT**

<input type="checkbox"/>	JSMF	Fixed Joystick Mount	N/C
<input type="checkbox"/>	JMSA	Swing Away Joystick Mount - Height adjustable, Angle adjustable and Multi-Axis	\$412.00
<input type="checkbox"/> Set-up with standard installation to swing-away sideways <input type="checkbox"/> Set-up at 35° to swing-away sideways and downwards			
<input type="checkbox"/>	JSMCM	Midline Mounted Swing Away Joystick Mount	\$412.00

**▼ JOYSTICK MOUNT SUPPORT**

<input type="checkbox"/>	JSMFH	Joystick mount support Bracket	STD
--------------------------	-------	--------------------------------	-----

**▼ JOYSTICK HAND MOUNT**

<input type="checkbox"/>	HML	Install Joystick on Left hand side	<input type="checkbox"/>	HMR	Install Joystick on Right hand side
--------------------------	-----	------------------------------------	--------------------------	-----	-------------------------------------

**▼ JOYSTICK HANDLE**

<input type="checkbox"/>	JSH	Standard Joystick Gimble	STD
<input type="checkbox"/>	JBAT	Stick Handle 3" in Length	\$101.00
<input type="checkbox"/>	JBH	Ball Handle	\$101.00
<input type="checkbox"/>	JFBH	Foam Ball Handle	\$101.00
<input type="checkbox"/>	JTH	T-Handle	\$101.00
<input type="checkbox"/>	JMUSH	<b>New</b> 2.5" Mushroom Handle	\$122.00
<input type="checkbox"/>	JGP	Goalpost Handle	\$101.00

**▼ ATTENDANT DRIVE CONTROL OPTION**

<input type="checkbox"/>	EADC	Attendant Hand Control for R-net EL & PM Controls	\$1,131.00
--------------------------	------	---	------------

**▼ ATTENDANT JOYSTICK HAND MOUNT**

<input type="checkbox"/>	EADCL	Install Attendant Joystick on Rear Left hand s	<input type="checkbox"/>	EADCR	Install Attendant Joystick on Rear Right hand side
--------------------------	-------	--	--------------------------	-------	--

**▼ ENHANCED DISPLAY & DIGITAL INTERFACE BOX (Included with specialty controls)**

<input type="checkbox"/>	EDIB-F-CA	15 OMNI2 (Enhanced Display & Digital Interface Box with Infrared)	\$2,201.00
<b>Note 15</b> Allows the switch operation of alternative devices through driver controls. Requires one of the R-net fully Expandable Controllers (ECEL90 or ECR120)			

**▼ ENHANCED DISPLAY MOUNT**

<input type="checkbox"/>	EDSAM	User Swing Away Mount	<input type="checkbox"/>	EDML (Install Display on LH side)	<input type="checkbox"/>	EDMR (Install Display on RH side)	N/C
--------------------------	-------	-----------------------	--------------------------	-----------------------------------	--------------------------	-----------------------------------	-----

▼ ENVIRONMENTAL CONTROL MODULE

Table with 3 columns: Item Name, Price, Description. Includes ECU (Environment control module - (Output Module)) for \$1,131.00. Note 16: Requires one of the fully expandable R-net controllers (ECEL90 or ECR120).

▼ BLUETOOTH MOUSE EMULATOR (Includes USB Dongle)

Table with 3 columns: Item Name, Price, Description. Includes EDBT (Bluetooth Mouse Emulator (USB Dongle included)) for \$1,131.00. Note 17: Requires one of the fully expandable R-net controllers (ECEL90 or ECR120). Not required if the ECRJC35B joystick has been selected. Not limited to only one selection. Any two can be installed on the same chair, but only one can be used at a time.

- Windows and MAC OS  Android Devices  iPhone, iPad and iPod

▼ POWER OPTIONS & SEATING

Table with 3 columns: Item Name, Price, Description. Includes SAABA (\$881.00), SAPT-X (\$7,385.00), SAPRX (\$9,090.00), SAPTRX (\$13,070.00), SAPESX (\$3,663.00), ADDPESX (\$3,404.00), ADDMR (\$1,131.00). Notes 18-20 describe adjustable seats and weight capacity.

▼ ELECTRONICS - SINGLE OR DUAL (1 OR 2 ACTUATORS) POWER SEAT FUNCTION CONTROLLERS

Table with 3 columns: Item Name, Price, Description. Includes EPSV, EPSEL, SAC, DAC. Notes 21-25 describe activation methods and controller requirements.

▼ ELECTRONICS - MULTIPLE (3 to 6 ACTUATORS) POWER SEAT FUNCTION CONTROLLER

Table with 3 columns: Item Name, Price, Description. Includes MAC (\$2,720.00), SWTJSM, SWDUALACT, SWMULTI, SWMULTIT, SWPUSH4, SWPUSH5. Notes 26-27 describe multiple actuator controllers.

▼ SWITCH ACTIVATION HAND MOUNT

Switch Installation:  TSML (Install Switch on LH side)  TSMR (Install Switch on RH side) N/C

**▼ STANDARD ADJUSTABLE SEAT WIDTH RANGE BETWEEN 15" AND 20"**

<input type="checkbox"/> SOWN15-20-19	(Width adjustment from 15" to 20" with a maximum Depth of 19" )	N/C
<input type="checkbox"/> SOWN15-20-22	(Width adjustment from 15" to 20" with a maximum Depth of 22" )	N/C
<input type="checkbox"/> 15" (SOWN15)	<input type="checkbox"/> 17" (SOWN17)	<input type="checkbox"/> 19" (SOWN19S)
<input type="checkbox"/> 16" (SOWN16)	<input type="checkbox"/> 18" (SOWN18)	<input type="checkbox"/> 20" (SOWN20S)

**▼ LARGE ADJUSTABLE SEAT WIDTH RANGE BETWEEN 19" AND 24"**

<input type="checkbox"/> SOWN19-24-22	(Width adjustment from 19" to 24" and Depth range of 18" to 22" )	N/C
<input type="checkbox"/> 19" (SOWN19L)	<input type="checkbox"/> 21" (SOWN21)	<input type="checkbox"/> 23" (SOWN23)
<input type="checkbox"/> 20" (SOWN20L)	<input type="checkbox"/> 22" (SOWN22)	<input type="checkbox"/> 24" (SOWN24)

**▼ ADJUSTABLE SEAT DEPTH RANGE BETWEEN 15" AND 22"**

<input type="checkbox"/> 15" (SODN15)	<input type="checkbox"/> 16" (SODN16)	<input type="checkbox"/> 17" (SODN17)	<input type="checkbox"/> 18" (SODN18)	STD
<input type="checkbox"/> 19" (SODN19)	<input type="checkbox"/> 20" (SODN20)	<input type="checkbox"/> 21" (SODN21)	<input type="checkbox"/> 22" (SODN22)	\$334.00

**▼ BACK CANE HEIGHT (The back cane height measurement is taken from the seat pan to the top of the push handles)**

<input type="checkbox"/> 16" (SOBCH16)	<input type="checkbox"/> 22" (SOBCH22)
<input type="checkbox"/> 18" (SOBCH18)	<input type="checkbox"/> 24" (SOBCH24)
<input type="checkbox"/> 20" (SOBCH20)	<input type="checkbox"/> 26" (SOBCH26)

Note 28 22" is the minimum back cane height required with a power or manual recline.

**▼ PRE-SET SEAT ANGLE (For Non-Reclining Backs)**

<input type="checkbox"/> SOBAA	Please Select one:	<input type="checkbox"/> 86°	<input type="checkbox"/> 90°	<input type="checkbox"/> 94°	<input type="checkbox"/> 98°	<input type="checkbox"/> 102°	<input type="checkbox"/> 106°	N/C
--------------------------------	--------------------	------------------------------	------------------------------	------------------------------	------------------------------	-------------------------------	-------------------------------	-----

**▼ BACK SUPPORT**

<input type="checkbox"/> STDSB	29	Contour back support 15" to 20" Wide with a maximum cushion height of 20" (15" to 20" Wide)	\$624.00
<input type="checkbox"/> STDSB/CUSTOM		Contour back support (Non Standard Sizes) Specify: _____ Width: _____ Depth: _____	\$1,333.00

Note 29 It is recommended to use the Contour Back Support option with the power or manual recline.

**▼ BACK SUPPORT TOTAL HEIGHT (MEASURED FROM THE SEAT PAN TO THE TOP OF THE BACK)**

<input type="checkbox"/> 16" (SOBH16)	<input type="checkbox"/> 24" (SOBH24)
<input type="checkbox"/> 18" (SOBH18)	<input type="checkbox"/> 26" (SOBH26)
<input type="checkbox"/> 20" (SOBH20)	<input type="checkbox"/> 28" (SOBH28)
<input type="checkbox"/> 22" (SOBH22)	

Note 30 22" is the minimum back height required with a power or manual recline.

Note 31 The space between the bottom of the back and the seat pan is 4" with a non-reclining back and 6" with a reclining back.

**▼ SEAT CUSHIONS**

**► Contour Seat Cushion**

<input type="checkbox"/> SCPSC	Contour Seat Cushion (15" to 20" Wide, with maximum 22" seat depth)	\$624.00
<input type="checkbox"/> SCPSC/CUSTOM	Contour Seat Cushion (Non Standard Sizes) Width: _____ Depth: _____	\$1,359.00

**▼ CANTILEVER - HEIGHT ADJUSTABLE ARMRESTS - OPENS WITH A 5° ANGLE (For Non-Reclining Backs - Choose 2)**

<input type="checkbox"/> SOAFBFLH	Flip-Back Cantilever - Full Length Left Side (Height Adjustable from 7" to 15" with 1" increments)	STD
<input type="checkbox"/> SOAFBFRH	Flip-Back Cantilever - Full Length Right Side (Height Adjustable from 7" to 15" with 1" increments)	STD
<input type="checkbox"/> SOAFBDLH	Flip-Back Cantilever - Desk Length Left Side (Height Adjustable from 7" to 15" with 1" increments)	N/C
<input type="checkbox"/> SOAFBDRH	Flip-Back Cantilever - Desk Length Right Side (Height Adjustable from 7" to 15" with 1" increments)	N/C

**▼ DUAL POST - HEIGHT ADJUSTABLE ARMRESTS - INCLUDES TRANSFER GRIPS (For Non-Reclining Backs - Choose 2)**

<input type="checkbox"/> SOADPFLH	Dual Post Height Adjustable - Full Length Left Side (Height Adjustable from 10 ¼" - 15 ¼")	\$202.00
<input type="checkbox"/> SOADPFRH	Dual Post Height Adjustable - Full Length Right Side (Height Adjustable from 10 ¼" - 15 ¼")	\$202.00
<input type="checkbox"/> SOADPDLH	Dual Post Height Adjustable - Desk Length Left Side (Height Adjustable from 10 ¼" - 15 ¼")	\$202.00
<input type="checkbox"/> SOADPDRH	Dual Post Height Adjustable - Desk Length Right Side (Height Adjustable from 10 ¼" - 15 ¼")	\$202.00

**▼ RECLINING FLIP-UP, HEIGHT ADJUSTABLE ARMRESTS (For Reclining Backs - Choose 2)**

<input type="checkbox"/> SOANIFLH	Reclining armrest - Full Length Left Side (Height Adjustable from 10 ½" to 14" with ½" increments)	N/C
<input type="checkbox"/> SOANIFRH	Reclining armrest - Full Length Right Side (Height Adjustable from 10 ½" to 14" with ½" increments)	N/C
<input type="checkbox"/> SOANIDLH	Reclining armrest - Desk Length Left Side (Height Adjustable from 10 ½" to 14" with ½" increments)	N/C
<input type="checkbox"/> SOANIDRH	Reclining armrest - Desk Length Right Side (Height Adjustable from 10 ½" to 14" with ½" increments)	N/C

**▼ ARMREST HEIGHT ADJUSTMENT AT FACTORY**

Please set the Armrest Height at: \_\_\_\_\_ " (The Armrest Height is taken from the seat pan to the top of the armpad. Not required for Dual Post)

**▼ ARMPADS**

**► Foam Type Selection**

<input type="checkbox"/>	SOFAPD	Desk Length Standard (10")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	STD
<input type="checkbox"/>	SOFAPP	Full Length Standard (14")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	STD
<input type="checkbox"/>	SOWAPD	32 Desk Length Waterfall (10")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$34.00
<input type="checkbox"/>	SOWAPP	32 Full Length Waterfall (14")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$34.00

**► Hybrid Gel and Foam Type Selection (The Gel always stays soft)**

<input type="checkbox"/>	SOGELD	Desk Length Hybrid Gel and Foam Flat armpad (10" x 2.5")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$34.00
<input type="checkbox"/>	SOGELF	Full Length Hybrid Gel and Foam Flat armpad (14" x 2.5")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$34.00
<input type="checkbox"/>	SOWGELD	32 Desk Length Hybrid Gel and Foam Waterfall armpad (10")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$54.00
<input type="checkbox"/>	SOWGELF	32 Full Length Hybrid Gel and Foam Waterfall armpad (14")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$54.00

**► Specialty Type Selection**

<input type="checkbox"/>	SOCAR	33 Arm Trough with Hand Support (18" x 4 1/2")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$222.00
--------------------------	-------	--	--	--	----------

Note 32 The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports  
 Note 33 Arm Trough Armpads may interfere with Armrests Flipping back.

**▼ FOOTRESTS**

<input type="checkbox"/>	SOSAFA	37 Adult Swing-Away Footrests	Set at: <input type="checkbox"/> 60°	<input type="checkbox"/> 38 70°	\$279.00
			Select Footrest Range: <input type="checkbox"/> SOFLA10 (10 1/2" - 15 1/2")	<input type="checkbox"/> SOFLA14 (14 1/2" - 19 1/2")	
<input type="checkbox"/>	SOAAF	37 Adult Heavy Duty Swing-Away footrests	Set at: <input type="checkbox"/> 60°	<input type="checkbox"/> 38 70°	\$390.00
			Select Footrest Range: <input type="checkbox"/> SOFLHD12 (12 1/2" - 15 1/2")	<input type="checkbox"/> SOFLHD14 (14 1/2" - 19 1/2")	
<input type="checkbox"/>	ADDPL	34, 37 Adult Power Elevating & Articulating Legrests	Please Select Function: <input type="checkbox"/> Individual 35	<input type="checkbox"/> Combined 36	\$3,177.00
			Select Footrest Range: <input type="checkbox"/> SOFLPW12 (12 1/2" - 15 1/2")	<input type="checkbox"/> SOFLPW14 (14 1/2" - 19 1/2")	

37 Please set the Footrest Length at: \_\_\_\_\_ "

Note 34 This options includes calf pads  
 Note 35 Individual legs count as two power functins.  
 Note 36 Combined legs count as one power function.  
 Note 37 The Footrest Length is taken from the top of the seat pan to the top of the footplate.  
 Note 38 The 70° footrest angle is not compatible with the 9" front casters (FC9PT and FC9PI) and may interfere with the casters when selecting the straight front caster arm option

**▼ FOOTPLATES & ACCESSORIES**

<input type="checkbox"/>	SOSFAF	Flip-up Footplates		STD
<input type="checkbox"/>	SOAFP	Angle Adjustable & Flip-up Footplates - Each		\$145.00
<input type="checkbox"/>	SOAHL	Heel Loops - Each		\$44.00

▼ CENTER MOUNTS

<input type="checkbox"/>	SOAAF90	41,42,43,44 .45,46	90° Fixed Center Mount Interface Bracket (Offers Manual adjustability from 95° to 70°)	\$256.00
Select Foot Platform Range: <input type="checkbox"/> SOFBA12 (12 ½" - 16 ½") <input type="checkbox"/> SOFBA16 (16 ½" - 19 ½")				
<input type="checkbox"/>	ADDPF	39,40,41,43 .44,45	90° Power Center Mount Elevating & Articulating Foot Platform	\$3,177.00
Select Foot Platform Range: <input type="checkbox"/> SOFBPW12 (12 ½" - 16 ½") <input type="checkbox"/> SOFBPW16 (16 ½" - 20 ½")				
	44	Please set the Foot Platform Length at: _____"		
	Note 39	The power centermount counts as one power function		
	Note 40	This option includes calf pads		
	Note 41	Links and attaches the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB) to the seat		
	Note 42	With the split footplates option (SOAFPFB), this bracket offers independent height adjustability for each footplate.		
	Note 43	Please select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB)		
	Note 44	The Foot Platform length is taken from the top of the seat pan to the top of the foot platform or split footplates.		
	Note 45	Centermount Front Riggings are not available with the Straight Front Arm option due to Mechanical conflicts.		
	Note 46	This option is not compatible with the 9" Caster option.		

▼ FOOT PLATFORMS & ACCESSORIES

<input type="checkbox"/>	SOFUFP	Flip-up one piece Foot Platform	N/C
<input type="checkbox"/>	SOAFPFB	Angle Adjustable & Flip-up Split Footplates for center mount - Each	\$145.00
<input type="checkbox"/>	SOACPCM	Calf Pads for 90° Fixed Centermount (SOAAF90) with adjustable hardware - Each	\$145.00

▼ SEAT POSITIONING (NOT COMPATIBLE WITH THE WC19-ISO OPTION)

<input type="checkbox"/>	SOAUTO	55" Seat Positioning Strap	STD
--------------------------	--------	----------------------------	-----

▼ HEADREST

Please select one of the headrest options with the mounting hardware and specify where the headrest will be installed

<input type="checkbox"/>	SONUAHR	Using headrest mounting hardware other than SOMAH located below. Specify: _____	N/C
<input type="checkbox"/>	SOFHRS	Small Curved Headrest Pad 4" x 9"	\$228.00
<input type="checkbox"/>	SOFHRM	Medium Curved Headrest Pad 5" x 10"	\$228.00
<input type="checkbox"/>	SOFHRL	Large Curved Headrest Pad 6" x 11"	\$228.00
<input type="checkbox"/>	SOMAH	Multi-axis Removable Mounting Hardware	\$183.00
<input type="checkbox"/>	SOHRBC	Bracket to Mount on the Contour Back Support (Without Recline)	\$111.00
<input type="checkbox"/>	SOHRPR	Bracket to Mount on the Seat Frame (For Reclining Backs)	\$111.00

Please select one of the headrest options with the mounting hardware and specify where the headrest will be installed (The frame or back cushion)

▼ LATERAL TRUNK SUPPORTS

<input type="checkbox"/>	SOSLS3-C	Adjustable Laterals with Small Curved Pads (Each) - 3" x 5"	\$675.00
<input type="checkbox"/>	SOSLM4-C	Adjustable Laterals with Medium Curved Pads (Each) - 4" x 6"	\$675.00
<input type="checkbox"/>	SOSLL5-C	Adjustable Laterals with Large Curved Pads (Each) - 5" x 7"	\$675.00

▼ HIP & THIGH SUPPORTS

<input type="checkbox"/>	SOHTSS10	Adjustable Hip & Thigh Support Pads (Each) - Small 10" x 4"	\$256.00
<input type="checkbox"/>	SOHTSM12	Adjustable Hip & Thigh Support Pads (Each) - Medium 12" x 5"	\$256.00
<input type="checkbox"/>	SOHTSL14	Adjustable Hip & Thigh Support Pads (Each) - Large 14" x 6"	\$256.00
<input type="checkbox"/>	SOHTSARH	Adjustable and Removable Hip & Thigh Support Hardware - Each	\$485.00

Please select one of the hip & thigh support pads and the adjustable & removable hardware.

▼ ACCESSORIES

<input type="checkbox"/>	SOCH	Cup Holder (Self-Leveling)	\$101.00
<input type="checkbox"/>	SOCPH	Ram Mounts X-Grip Cell Phone Holder	\$173.00
<input type="checkbox"/>	SOTHS	Ram Mounts X-Grip Small Tablet Holder	\$196.00
<input type="checkbox"/>	SOTHL	Ram Mounts X-Grip Large Tablet Holder with Swivel	\$295.00

▼ REFLECTOR KIT

<input type="checkbox"/>	REFLECT	Red & White Reflector Kit Stickers (They will not be installed - Included in a bag with the owner's manual)	N/C
--------------------------	---------	---	-----

▼ SWITCHES AND OPTIONS

<input type="checkbox"/>	ASCPCP	PGDT R-NET Dongle and PC programmer kit	\$652
<input type="checkbox"/>	ASCPVR2	PGDT VR2 Dongle and PC programmer kit	\$652.00
<input type="checkbox"/>	SSBS-STD	Buddy Button Switch	\$122.00
<input type="checkbox"/>	SSES-STD	Egg Switch-Black	\$122.00
<input type="checkbox"/>	SSML	Ultra Light Switch	\$168.00

▼ **SPECIALTY DRIVER CONTROL OPTIONS**

▼ **Therafin Sip n' Puff Kit**

<input type="checkbox"/>	SCSNPK	47, 48	Therafin Sip n' Puff Kit			\$1,701.00
<input type="checkbox"/>	SCSNPBTk	48	Therafin Sip n' Puff Breath Tube Kit			\$552.00
<input type="checkbox"/>	SCSRMSNP	48	Stealth Rock-n-Lock Swingaway Arm for Sip n' Puff	<input type="checkbox"/> SWL (Swings on Left)	<input type="checkbox"/> SWR (Swings on Right)	\$307.00
<input type="checkbox"/>	SCSRM	48	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings on Left)	<input type="checkbox"/> SWR (Swings on Right)	\$552.00
<input type="checkbox"/>	SSES-STD	48, 50	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$122.00
<input type="checkbox"/>	SSBS-STD	48	Buddy Button Switch (Stop Switch) - Red			\$122.00
<input type="checkbox"/>	SOMAH	48	Multi-axis Headrest Removable Mounting Hardware			\$183.00
<input type="checkbox"/>	SOFHRS	49	Small Curved Headrest Pad 4" x 9"			\$228.00
<input type="checkbox"/>	SOFHRM	49	Medium Curved Headrest Pad 5" x 10"			\$228.00
<input type="checkbox"/>	SOFHRL	49	Large Curved Headrest Pad 6" x 11"			\$228.00

▼ **Switch-It Proximity Head Array**

<input type="checkbox"/>	SCSIHAA	47, 48	Adult - 3 Switch Proximity Head Array with 1/8" Mono Jack for mode/reverse (E2330)			\$5,679.00
<input type="checkbox"/>	SCSRM	48	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings on Left)	<input type="checkbox"/> SWR (Swings on Right)	\$552.00
<input type="checkbox"/>	SSES-STD	48, 50	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$122.00
<input type="checkbox"/>	SSBS-STD	48	Buddy Button Switch (Stop Switch) - Red			\$122.00
<input type="checkbox"/>	SOMAH	48	Multi-axis Headrest Removable Mounting Hardware			\$183.00

▼ **Proportional RIM Head Control Package**

<input type="checkbox"/>	SCPHC	47, 48	Proportional RIM Head Control Package			\$2,781.00
<input type="checkbox"/>	SCSRM	48	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings on Left)	<input type="checkbox"/> SWR (Swings on Right)	\$552.00
<input type="checkbox"/>	SSES-STD	48, 50	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$122.00
<input type="checkbox"/>	SOMAH	48	Multi-axis Headrest Removable Mounting Hardware			\$183.00

▼ **Switch-It MicroPilot Light Touch Proportional Mini joystick for chin / lip control**

<input type="checkbox"/>	SCSIMP	47, 48	MicroPilot Light touch mini Proportional Joystick			\$4,453.00
<input type="checkbox"/>	SCSRMMP	48	Stealth Rock-n-Lock Swingaway Arm for Chin MP	<input type="checkbox"/> SWL (Swings on Left)	<input type="checkbox"/> SWR (Swings on Right)	\$552.00
<input type="checkbox"/>	SCSRM	48	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings on Left)	<input type="checkbox"/> SWR (Swings on Right)	\$552.00
<input type="checkbox"/>	SSES-STD	48, 50	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$122.00
<input type="checkbox"/>	JCC		Chin Cup			\$57.00
<input type="checkbox"/>	SOMAH	48	Multi-axis Headrest Removable Mounting Hardware			\$183.00
<input type="checkbox"/>	SOFHRS	49	Small Curved Headrest Pad 4" x 9"			\$228.00
<input type="checkbox"/>	SOFHRM	49	Medium Curved Headrest Pad 5" x 10"			\$228.00
<input type="checkbox"/>	SOFHRL	49	Large Curved Headrest Pad 6" x 11"			\$228.00

▼ **Proportional Compact or Mini joystick for Midline Mount**

<input type="checkbox"/>	SCPGMJS	47, 48	PG Drives JC200 Proportional Compact Remote Joystick for Midline Mount			\$1,059.00
<input type="checkbox"/>	SCGMJSM	51	Gatlin Midline Swing-away Mount for JC200	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$435.00
<input type="checkbox"/>	SCSIMP	47, 48	Switch-It MicroPilot Light touch mini Proportional Joystick			\$4,453.00
<input type="checkbox"/>	SCGMPM	52	Gatlin Midline Swing-away Mount for MicroPilot	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$474.00
<input type="checkbox"/>	SCBTM	53	Bullet Hand Tray for MicroPilot	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$502.00

▼ **Proportional Compact or Mini joystick for Armrest Mount**

<input type="checkbox"/>	SCPGMJS	49, 50	PG Drives JC200 Proportional Compact Remote Joystick for Armrest Mount			\$1,059.00
<input type="checkbox"/>	SCARMJSM	51	Swing-away and Multi-Axis armrest mount	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$435.00
<input type="checkbox"/>	SCSIMP	47, 48	Switch-It MicroPilot Light touch mini Proportional Joystick			\$4,453.00
<input type="checkbox"/>	SCARMPM	52	Swing-away and Multi-Axis armrest mount	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$474.00
<input type="checkbox"/>	SCBTM	53	Bullet Hand Tray for MicroPilot	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$502.00

- Note 47 Requires the OMNI Enhanced Display and Digital Interface Box (EDIB)
- Note 48 This option must be selected for the kit
- Note 49 Choose one (1) Headrest Pad
- Note 50 The mode switch type can be changed from this selection on Page #6
- Note 51 Must be selected for use with the PG Drives JC200 Mini Joystick
- Note 52 Must be selected for use with the Switch-It MicroPilot
- Note 53 This option is only compatible with the Switch-It MicroPilot



**▼ MADE-TO-ORDER SECTION**

Note 54 Lead times will increase depending on the made-to-order pipeline at the time the purchase order is received. Please call for Quote and specifications for these options.

**▼ O<sup>2</sup> Holder**

<input type="checkbox"/>	SUPOXY	55	Non-Articulating Oxygen Tank Holder	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$435.00
<input type="checkbox"/>	SUPOXYA	56	Articulating Oxygen Tank Holder	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$2,040.00
			<input type="checkbox"/> D-Size Oxygen Tank (* 22" Tall x 4.3" Diameter) * The height includes the valve					
			<input type="checkbox"/> E-Size Oxygen Tank (* 31" Tall x 4.3" Diameter) * The height includes the valve					
	Note 55		Not available with a power tilt and/or power recline					
	Note 56		Required with a power tilt and/or power recline					

**▼ VENT TRAY (Please complete in Full)**

<input type="checkbox"/>	SOVENT	Vent Tray						
	<input type="checkbox"/>	SOVENT-A	57	Articulating Vent Box (Required with Power Tilt and/or Recline)	Model:	_____		\$1,957.00
	<input type="checkbox"/>	SOVENT-F	57	Fixed Vent Box	Model:	_____		\$491.00
<input type="checkbox"/>	SOVTBE	57, 58	Extended Battery Compartment					\$600.00
	<input type="checkbox"/>	The Vent is being powered by 1 U1 Battery and the chair by 2 Group 22NF batteries						
	<input type="checkbox"/>	The Vent is being powered by 1 U1 Battery and the chair by 2 Group 24 batteries						
	<input type="checkbox"/>	The Vent is being powered by 1 Group 22NF Battery and the chair Powered by 2 Group 22NF Batteries						
	<input type="checkbox"/>	The Vent is being powered by 1 Group 24 Battery and the chair Powered by 2 Group 22NF Batteries						
	Note 57	Sunrise Medical does not supply the battery, charger or the wiring for the Vent. These items are normally provided by the Vent Manufacturer.						
	Note 58	This option is not available with the Straight Front Caster Arms						

**▼ OTHER MADE-TO-ORDER OPTIONS (Request For Quote)**

<input type="checkbox"/>	PAT	59, 60	Power Anterior Tilt					RFQ
<input type="checkbox"/>	CUSTOM	59	Power Lateral Tilt 15° <input type="checkbox"/> Raise the Left Side only <input type="checkbox"/> Raise the Right Side only <input type="checkbox"/> Raise the Left & Right Side					RFQ
<input type="checkbox"/>	CUSTOM	59, 60	Power Elevating Stretcher (Please call CSR for availability with configuration and user weight capacity)					RFQ
<input type="checkbox"/>	CUSTOM	59, 60	Power Flip-up Foot Platform (Please call CSR for availability with configuration and user weight capacity)					RFQ
<input type="checkbox"/>	CUSTOM	59, 62	Power Vertical Lifting Foot Platform (Please call CSR for availability with configuration and user weight capacity)					RFQ
<input type="checkbox"/>	CUSTOM	59	Independent Depth Adjustability built into the Footrest Receiver that is attached to the Seat Pan					RFQ
<input type="checkbox"/>	CUSTOM	59	Endomorphic Back Canes (Pair) <input type="checkbox"/> 1" per side towards the inside = 2" in <input type="checkbox"/> 2" per side towards the inside = 4" in					RFQ
<input type="checkbox"/>	CUSTOM		Transfer Handles (Pair)					RFQ
<input type="checkbox"/>	CUSTOM		Manual Wheel Locks					RFQ
<input type="checkbox"/>	CUSTOM	59	Non Standard Frame Colors (Not all colors are available, please contact CSR for availability)					RFQ
<input type="checkbox"/>	CUSTOM # 1	59	Other made-to-order options. Specify in notes section below.					RFQ
<input type="checkbox"/>	CUSTOM # 2	59	Other made-to-order options. Specify in notes section below.					RFQ
	Note 59	Please specify details and specifications for these options.						
	Note 60	The maximum user weight capacity with this option is 250 lbs						
	Note 61	The 11" Power Elevating Seat Module and the Power Elevating Legs or Centermount options are compatible.						
	Note 62	This option is not available with the Straight Front Caster Arms.						

**▼ NOTES / SPECIAL REQUESTS / MADE-TO-ORDER SPECIFICATIONS**

---



---



---



---



---