

April 2025

Ti

Grade 2



Customer Contact Details

Customer Name: _____

Advisor: _____

Address: _____

Phone: _____

Email: _____

Ship to Customer Address: ☐

Shipping Name: _____

Shipping Address: _____

Phone: _____

Email: _____

Additional Information

Notes: _____

Current Wheelchair: _____

Customer Weight: _____

Customer Height: _____

Customer Disability: _____

Other

SERIAL NO: _____

Guide:		☒ Base product		☒ Required option	☐ Optional Extra	Critical Information	Information Needed
PRODUCT STYLE							
☒	HXTXZ010	HILITE XTR					
	\$ 7,066.00						
☐	XA070	Open design					
☐	XA080	Closed design					
☐	HXTXA085	Crash tested upgrade					
	\$348.00						
Max User Weight 115kg (253lbs)							
			XA070 - Open design		XA080 - Closed design		
PLEASE NOTE: ALL RGK CHAIRS ARE HAND BUILT AND HAVE A MANUFACTURING MEASUREMENT TOLERANCE OF +/- 5MM							
CAD APPROVAL							
☐	XW009	QUOTE ONLY (NO CAD FOR APPROVAL)	☐	XW011	ORDER (NO CAD FOR APPROVAL)		
☐	XW012	QUOTE CAD drawing for approval	☐	XW010	ORDER CAD drawing for approval		
		\$550 FLAT FREE (2 customer amendments ONLY)			\$300 FLAT FREE (2 customer amendments ONLY)		
KEY MEASUREMENTS							
☒	XY010	Seat width <i>(Inside to inside of sideguards)</i>					
		mm	(Max - 500mm)				
☒	XY020	Seat length					
		mm	(Max - 500mm)				
☐	XY030	Extra Seat Length <i>(Provides additional leg support & hand grip)</i>					
		mm					
☒	XY050	Front frame height					
		mm	(Max - 550mm) (Min - 380mm)				
☒	XY060	Rear frame height					
		mm	(Max - 550mm) (Min - 360mm)				
☒	XY070	Frame length					
		mm	(Min = XY020 + XY030 + 140mm)				
☒	XY072	Castor position					
		mm					
		☐ Optimise measurements					
☐	XA020	Ergonomic seat					
	\$121.00						
		Specify length					
		mm	(Standard = 120mm)				


FRONT FRAME AND FOOTREST

Front Frame Options

- ☐ XY085 Seat taper design
\$224.00

i Specify width mm


(When combined with XY091, dimension is inbetween the centre of the main bends)



- ☐ HXTXY091 Parallel taper design
\$159.00

i Specify width mm


(When combined with XY080, dimensions must be smaller by 20mm)




Footrest Design (Choose One)

- ☐ XY100 Round design

Castor Size	Max Height	With XB200
3"	140mm	190mm
4"	150mm	200mm
5"	160mm	210mm




- ☐ XB200 High level footrest (Only with XB010, HXTXB080 & XB150)
\$201.00




- ☐ XY105 Square design

Castor Size	Max Height	With XB200
3"	135mm	185mm
4"	145mm	195mm
5"	155mm	205mm



- ☐ XB200 High level footrest (Only with XB010, HXTXB080 & XB150)
\$ 201.00



Footrest Individual Measurements

- ☐ XY115 Front footrest height

! mm

- ☐ XY120 Rear footrest height


i mm

(Std = 20 lower than front height)

- ☐ XY125 Footrest length

i mm

(Std = 100mm)
(160mm and over requires heavy duty upgrade XB150)



Footrest Style (Choose One)

- ☐ XB010 Tubular aluminium footrest
(Height adjustable)

- ☐ XB020 Tubular titanium footrest
(Height adjustable)

- ☐ XB070 Angle adjustable aluminium footrest
(Angle/height adjustable)
\$207.00

- ☐ HXTXB080 Tubular footrest welded into frame
(No height adjustment)
\$ 248.00

- ☐ XB030 Flip up plastic footrest
(Angle/height adjustable, with closed frame only)
\$ 352.00

- ☐ XB150 Heavy duty tubular footrest
(Height adjustable, for Freewheel use)
\$ 180.00

Footrest Cover

- ☐ XB050 Carbon fibre cover
(Not with XB070 & XB030)
\$248.00

- ☐ XB170 Aluminium cover
(Not with XB070 & XB030)
\$ 98.00

- ☐ XB165 Rubber grip tape
(3M Safety walk)(not with XB030)
\$ 86.00

SERIAL NO: _____

CASTOR WHEELS, FORKS & WHEEL LOCKS

Castor Wheels 100mm (4")

- ☐ XC060 Soft roll 4"x1.3"
(Black)
- ☐ XC110 Aluminium centre 4"x1"
\$ 232.00 (Silver hub and black urethane tyre)
- ☐ XC130 Aluminium centre 4"x1.4"
\$ 232.00 (choose colour below)
- ☐ Silver hub & black tyre (Std) ☐ Red hub & black tyre ☐ Purple hub & black tyre ☐ Green hub & black tyre
- ☐ Black hub & black tyre ☐ Blue hub & black tyre ☐ Gold hub & black tyre

125mm (5")

- ☐ XC070 Soft roll 5"x1.3"
(Black)
- ☐ XC100 Pneumatic Aluminium centre 5"x1"
\$ 248.00 (Grey)
- ☐ XC120 Froglegs Aluminium centre 5"x1"
\$ 232.00 (Silver hub and black urethane tyre)
- ☐ XC140 Froglegs Aluminium centre 5"x1.4"
\$ 232.00 (Silver hub and black urethane tyre)

150mm (6")


- ☐ XC080 Soft roll 6"x1.3"
- ☐ XC145 Frog Legs aluminium centre 6"x1.4"
\$ 310.00 (Silver hub and black urethane tyre)

Castor Forks

- ☐ XC192 Anti-Flutter castor fork
(Std when no option selected)
- ☐ XC210 Froglegs suspension castor fork
\$ 633.00
- Note: Colour follows accent ☐ Black ☐ Silver







Castor Pin Width

- ☒ XY330 Castor pin width
- ☐ Optimise measurements



Wheel Lock Options (choose one)

- ☐ XD020 Underseat Scissor
- ☐ XD030 Compact scissor
\$ 266.00
- ☐ XD050 Push to lock
- ☐ XD046 Out Front push to lock
\$ 409.00
- ☐ XD085 Eagle Brakes (Indicate colour)
\$ 468.00
- ☐ XD086 Eagle Brakes with L/W fixing
\$ 508.00 (Only with fixed centre of gravity)
(Indicate colour)
- ☐ XD010 Omit wheel locks
- ☐ BLACK ☐ PURPLE ☐ RED ☐ BLUE
Black Std if no option selected

SIDEGUARDS AND ARMRESTS	
Sideguard Aluminium Options ☐ XE030 Screw on aluminium \$ 327.00 <i>(N/A with HXTXG020)</i> ☐ XE080 Lift out aluminium \$ 532.00 <i>(N/A with HXTXG020)</i> ☐ HXTXE130 Fold flat aluminium \$ 690.00 Carbon Fibre Options ☐ XE050 Screw on carbon fibre \$ 726.00 <i>(N/A with HXTXG020)</i> ☐ XE100 Lift out carbon fibre \$ 766.00 <i>(N/A with HXTXG020)</i> ☐ HXTXE150 Fold flat carbon fibre \$ 1,121.00 Welded Options ☐ HXTXE220 Welded Titanium \$ 738.00 <i>(N/A with HXTXG020)</i> Armrest Options ☐ XF010 Lift out armrest with plastic sideguard \$ 592.00 <i>(Indicate fixed height)</i> ☐ XF020 Lift out armrest with Aluminium sideguard \$ 467.00 <i>(Indicate fixed height)</i> Height (Top of seatrail to top of armrest) mm Sideguard Measurements Optimised measurements ☐ 25mm above tyre <i>(Follow tyre)</i> ☐ Individual above tyre <i>(Follow tyre)</i> mm ☐ Level with tyre <i>(Follow tyre)</i>   Individual measurements ☐ Sideguard length mm ☐ Sideguard height mm 	Sideguard with Fender Aluminium Options ☐ XE040 Screw on aluminium with fender \$ 279.00 <i>(N/A with HXTXG020)</i> ☐ XE090 Lift out aluminium with fender \$ 579.00 <i>(N/A with HXTXG020)</i> ☐ HXTXE140 Fold flat aluminium with fender \$ 749.00 Carbon Fibre Options ☐ XE060 Screw on carbon fibre with fender \$ 952.00 <i>(N/A with HXTXG020)</i> ☐ XE110 Lift out carbon fibre with fender \$1,132.00 <i>(N/A with HXTXG020)</i> ☐ HXTXE160 Fold flat carbon fibre with fender \$1,362.00 Sideguard Measurements Fender Height ☐ 10mm above tyre <i>(Standard)</i> ☐ Individual above tyre <i>(Min = 10mm)</i> mm i mm  Fender position ☐ Centre of tyre <i>(Standard if no option selected)</i> ☐ Outside of tyre  

SERIAL NO: _____

BACKREST			
Backrest Height - Fixed or Adjustable (Choose One)			
☐ XY230 Fixed height backrest <i>Specify height</i> 	☐ XY245 Adjustable height backrest \$ 248.00 (+/25mm) <i>Specify height</i> 		
Backrest Design (Choose one)			
☐ XG010 Fixed backrest (Non-folding)	☐ HXTXG021 Folding backrest \$ 725.00 (Oracing style) ☐ HXTXG020 Folding backrest \$ 725.00 (with 2 position snap lock and angle adjust +5 Deg)		
Ergonomic Options			
☒ XY360 Backrest angle ☐ -5 ☐ -3 ☐ 0 ☐ +3 ☐ +5 <i>Individual</i> Deg (Max -10/+20 Deg)	☐ XG100 Ergonomic backrest \$ 146.00 <i>Vertical height</i> mm <i>Rearward position</i> mm		
☐ XG110 \$ 146.00 Inward taper (inside of backrest posts) <i>Specify width or per side</i> 	☐ XG120 \$ 146.00 Outward flare (inside of backrest posts) <i>Specify width or per side</i> 		
Push handle Options			
☐ XG070 \$ 222.00 Screw in push handles (N/A with XY245) ☐ XG080 \$ 434.00 Flip down push handles (N/A with XY245) ☐ XG085 \$ 709.00 Flip down push handles height adjustable (N/A with XY245)	☐ HXTXG090 \$ 543.00 Stroller push handles height adjustable (Oracing style) <i>Height from floor</i> mm ☐ Optimise (1000mm from floor) ☐ XG097 \$ 676.00 Integrated push handle (quickie design)		
Backrest Options			
☐ XG040 \$ 225.00 Heavy duty backrest upgrade (Recommended for hand cycling) (Not with XY245)	☐ XG130 \$ 115.00 Extra deep brace bar (Recommended for JAY backs) (Std - 90mm) <i>Specify depth</i> 		
ILSA Backrest Options			
☐ XG140 \$ 479.00 Integrated lateral support A - Left mm B - Left mm C - Left mm D - Left mm A - Right mm B - Right mm C - Right mm D - Right mm			


CAMBER BAR

Adjustable Centre of Gravity

☐ XH110

Square back adjustable centre of gravity


(25-150mm)


 mm


☐ XH120

Round back adjustable centre of gravity

(60-150mm)


 mm


☐ HXTXH090

V-style fixed centre of gravity
(Open frame only)

(25-125mm)


 mm


☐ HXTXH105

Fixed centre of gravity
(Closed frame only)

(25-125mm)


 mm


Specify Camber Angle

☒ XY390

Camber angle


☐ 0 ☐ 1 ☐ 2 ☐ 3 ☐ 5

☐ XY390

Camber angle


☐ 0 ☐ 1 ☐ 2 ☐ 3 ☐ 5

Camber Bar Options

☐ XH040

\$ 400.00

Handbike camber bar fixed



Camber angle (Handbike camber bar)

☐ 0 ☐ 5


Castor height from floor for handbike fitting
(Standard 30mm if not indicated)

☐ 30 ☐ 50

☐ XH040

\$ 400.00

Handbike camber bar fixed



Camber angle (Handbike camber bar)

☐ 0 ☐ 5


Castor height from floor for handbike fitting
(Standard 30mm if not indicated)

☐ 30 ☐ 50

Anti-Tip Options

☐ XJ010

\$ 167.00

Swingaway anti tip LEFT side (Only with square back frame. Not with XH040)

☐ XJ020

\$ 167.00

Swingaway anti tip RIGHT side (Only with square back frame. Not with XH040)

☐ XJ100

\$ 504.00

Single sport anti tip
(centre mounted)

REAR WHEELS			
Rear Wheel Size			
☐ XY420 24" Wheels	☐ XY430 25" Wheels		
Rear Wheel Gap			
☒ XY460 Rear wheel gap 15 ☐ 20 ☐ 25 ☐ Individual mm			
Wheel Options (Choose one)			
☐ XK010 Fusion 16 spoke <i>(Silver Hub, black spokes)</i>		☐ XK320 Spinergy XLX 12 spoke \$1,067.00	
☐ XK050 Spinergy LX 12 spoke \$ 1000.00		☐ XK060 Spinergy SPOX 18 spoke \$ 665.00	
☐ XK150 Spinergy CLX 18 carbon blade spoke \$ 3,139.00 <i>N/A IN SILVER HUB</i>		☐ XK070 Spinergy BLXL 18 blade spoke \$1,067.00 <i>N/A IN SILVER HUB</i>	
Spinergy Colour Options: (i)			
☐ Black spokes	☐ White spokes	☐ Yellow spokes	☐ XK180 Individual spoke design \$ 213.00 <i>(Specify colour)</i> Inside: _____ Outside: _____
☐ Red spokes	☐ Blue spokes	☐ Pink spokes	
☐ Purple spokes	☐ Orange spokes	☐ Green spokes	
Spinergy Anodise Options			
☐ XK190 Black hub <i>(std option)</i>	☐ XK200 Silver hub \$ 127.00 <i>(n/a with BLXL/CLX)</i>	☐ XK201 Red hub/rims \$ 409.00 <i>(LX, XLX & BLXL only)</i>	☐ XK205 Individual: \$568.00 Select colour: _____
☐ XK202 Purple hub/rims \$409.00 <i>(LX, XLX & BLXL only)</i>	☐ XK203 Blue hub/rims \$ 409.00 <i>(LX, XLX & BLXL only)</i>	☐ XK204 Gold hub/rims \$ 539.00 <i>(LX, XLX & BLXL only)</i>	
Supply Without Wheels			
☐ XK040 Supply without Wheels		Wheel Type: _____ Tyre Type: _____	
Handrim Position / Options			
☐ XY480 Narrow set handrim		☐ XY485 Wide set handrim	
☐ XL090 Remove excess mounting tabs \$146.00 <i>(recommended for small wheel gaps)</i>			
Handrim Styles			
☐ XL010 Aluminium hard anodised 6 tab		☐ XL040 Titanium 6 tab \$ 446.00	
☐ XL160 Ellipse 6 tab \$631.00 <i>(replacement for Surge)</i>		☐ XL035 Natural Fit LT \$ 743.00	
☐ XL020 Tetra grip 6 tab \$400.00		☐ XL150 Carbolife Gekko ☐ Black ☐ Silver \$1,406.00	
☐ XL100 Para grip 6 tab \$400.00		☐ XL170 Carbolife Curve L (Black) \$1,393.00	
Tyre Options			
☐ XM010 Kenda ☐ Grey ☐ Blue ☐ Red <i>(Std)</i>	☐ XM060 Schwalbe One \$ 133.00 <i>(24", 25" only)</i>		
☐ XM130 Kenda black \$ 28.00	☐ XM050 Kenda K885 mountain bike tyre \$ 40.00 <i>(24", 25" only)</i>		
☐ XM020 Schwalbe Marathon Plus \$ 98.00			
☐ XM040 Shox solid airless (Black) \$160.00			
Inner Tube Valve Options			
☐ XY500 Schraeder inner tube valve <i>(Standard if no option selected)</i>		☐ XY490 Presta inner tube valve	
Wheel Accessories			
☐ XK220 Titanium quick release axle pins (pair) <i>(N/A with XK010 & XK070)</i> \$160.00		☐ XK300 Offroad wheel set 22" rim with handrim, MTB tyre and axle pins <i>(Replaces 24" wheel with XM010)</i> \$ 933.00	
☐ XK240 Quad release axle nuts (pair) \$241.00		☐ XK310 Offroad wheel set 24" rim with handrim, MTB tyre and axle pins <i>(Replaces 25" wheel with XM010)</i> \$ 933.00	

FINISH OPTIONS									
Frame & Component Finish									
☐ XN010		Powder coated frame			☐ XN150		Powder coated components (Frame colour if not specified)		
<i>Frame Comp.</i>		<u>Gloss</u>			<i>Frame Comp.</i>		<u>Matt</u>		
☐	☐	Black gloss 103			☐	☐	Black matt 105		
☐	☐	White gloss 100			☐	☐	White matt 101		
☐	☐	Grey gloss 132			☐	☐	Grey matt 131		
☐	☐	Cream gloss 135			☐	☐	Light grey matt 133		
☐	☐	Red gloss 108			☐	☐	Navy blue matt 152		
☐	☐	Yellow gloss 120							
☐	☐	Forest green gloss 143							
☐	☐	Green gloss 107							
☐	☐	Dark blue gloss 151							
☐	☐	Royal blue gloss 154							
☐	☐	Blue gloss 174							
☐	☐	Duckegg blue gloss 141							
☐	☐	Pink gloss 188							
☐	☐	Orange gloss 180							
☐	☐	Purple gloss 157							
☐ HXTXN030		Bead blast frame (Raw)			☐ XN165		Bead blast components		
\$472.00		(All parts black unless XN165 chosen)			\$ 250.00				
Accent Colour									
<i>Note: Accent colour includes forks & castor caps.</i>									
☐		Black accent			☐		Silver accent		
☐		Red accent			☐		Blue accent		
☐		Purple accent							
UPHOLSTERY									
Seat Upholstery									
☐ XP040		Seat sling adjustable tension with straps			☐ XR050		Carbon fibre seat plate		
					\$ 480.00				
Backrest Upholstery									
☐ XP050		Airtech backrest strap adjustable			☐ XP070		Backrest upgrade to leather corners (Black)		
					\$ 121.00				
Backrest Stitching Options									
☐	Silver	☐	Black	☐	Blue	☐ Individual: _____			
☐	Red	☐	Orange	☐	Purple	_____			
RGK Logo Colour									
☐	White (Std)	☐	Green	☐	Blue	☐ Individual: _____			
☐	Red	☐	Orange	☐	Purple	_____			
Model Logo Colour									
☐	White (Std)	☐	Green	☐	Blue	☐ Individual: _____			
☐	Red	☐	Orange	☐	Purple	_____			
Backrest Accessories									
☐ XP110		Backrest with zip pocket							
\$ 127.00									

SERIAL NO: _____

CUSHION

Cushion foam options

- ☐ XP080 Cushion foam
\$ 121.00 (select density and size below)
- ☐ Soft (CMS040) ☐ Medium (CMS050)
- ☐ 50mm ☐ 75mm ☐ 100mm

Cushion cover options


- ☐ XP085 Nylon cushion cover
(Black)
- ☐ XP090 Air Lite cushion cover
\$ 98.00 (Black)
- ☐ XP100 Leather cushion cover
\$ 254.00 (Black)

Jay Products - PLEASE ATTACH JAY ORDER FORM

- ☐ JAY Cushion
(Price dependant on size / option) ☐ Jay Order form attached
- ☐ JAY Backrest
(Price dependant on size / option) ☐ Jay Order form attached
- ☐ JAY Whitmyer Headrest
(Price dependant on size / option) ☐ Jay Order form attached

FRAME OPTIONS

- ☐ XR060 Welded calf strap mounting loops
\$ 154.00 (To locate calf strap on front frame)
- Specify height from floor
- ☐ Optimise measurements



ACCESSORIES

- ☐ XR080 Frontwheel assembly
\$1,053.00 (Only with Tubular style footrest)
- ☐ XP120 Frame protectors Neoprene zip fasten
\$ 84.00 (Back)
- ☐ XP420 Frame protectors Nylon velcro fasten
\$ 81.00 (Black)
- ☐ XR010 Seat pouch velcro attach
\$ 59.00 (Back)
- ☐ XR200 Rucksack with RGK logo
\$ 151.00 (Back)
- ☐ XR020 Underseat large storage bag
\$ 127.00 (Back)
- ☒ XR090 Tool kit
- ☒ XS060 Calf strap

ADDITIONAL NOTES

ADDITIONAL NOTES