

MWD Crossover Compact Electric Wheelchair

EIPW60



Customer				Ship To			
Account				Name			
Account #		Date		Address			
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #		City		Province	
Buyer				PC	Phone		
Mark For				Email / Fax			

Please submit orders to: cscanada@sunmed.com

Seating Specifications


Typical WC Weight w/ Batteries	Starting at 320 lbs. Weights vary depending on options		
Power Seating	50° of CG Tilt / Power Recline / 12" Seat Elevator	Seat Type	MPS and Rehab Style
Seat Size Range	Width: 13-23" Depth: 13-24"	STF Height Range	17", 18", 19"

Drive Base Specifications

Drive Wheel Position	Mid Wheel Drive	Drive Control	R-Net PM 120A controller
User Weight Capacity	350 lbs.	Drive Motors	4-pole 6mph
Base Size: Off-Road	Width 26" Length 39" Turning Radius 22.5"	Drive Wheels	14" Off-Road Drive Wheels 14" Urban Drive Wheels
Base Size: Urban	Width 24" Length 39" Turning Radius 22.5"	Tires / Casters	Pneumatic or Solid inserts
		Charger	Dual Mode
		Battery Size	Group 24
		Ground Clearance	3"
		Color Accents	10 colors

Client Measurements

a) Top of head to seat	_____	in
b) Top of shoulder to seat	_____	in
c) Underarm to seat	_____	in
d) Elbow to seat	_____	in
e) Top of leg to seat	_____	in
f) Shoulder width	_____	in
g) Chest width	_____	in
i) Hip width	_____	in
j) Actual sitting depth	_____	in
m) Foot length	_____	in
n) Knee to foot	_____	in
o) Seat width	_____	in



Notes

✓ OPTION	DESCRIPTION	RETAIL
Seating Options		\$
<input type="checkbox"/> 600SF4	REHAB Seat and Back - Back Canes and adjustable flat seat pan	NC
<input type="checkbox"/> 600SF5	MPS Seat and Back - Captain style seating	NC
<input type="checkbox"/> 600SF6	MPS Back and REHAB Seat - Captain style backrest and adjustable flat seat pan	NC
Base		
Magic 360 Base - Includes Transit and lights		End Item Number: EIPW60
		\$
<input type="checkbox"/> 600MW61	Magic Mobility - 360 No Power Seating Options	\$ 26,745.00
<input type="checkbox"/> 600MW62	Magic Mobility - 360 Single Power Seating Options	\$ 26,745.00
<input type="checkbox"/> 600MW63	Magic Mobility - 360 Multi Power Seating Options	\$ 26,745.00
<input type="checkbox"/> 600MW64	Magic Mobility - 360 No Power Seating Options - MPS (Captain)	\$ 26,745.00
<input type="checkbox"/> 600MW65	Magic Mobility - 360 Single Power Seating Options - MPS (Captain)	\$ 26,745.00
<input type="checkbox"/> 600MW66	Magic Mobility - 360 Multi Power Seating Options - MPS (Captain)	\$ 26,745.00
Drive Wheels		\$
<input type="checkbox"/> 600DW80	Off-Road Drive Wheels and Tires	NC
<input type="checkbox"/> 600DW81	Urban Drive Wheels and Tires - Pneumatic	NC
<input type="checkbox"/> 600DW82	Urban Drive Wheels and Tires - Solid	\$ 63.00
Optional Drive Tire Options		\$
<input type="checkbox"/> 6008A90	Tire Sealant for 2 Pneumatic Drive Tires	\$ 51.00
Base Accessories		\$
<input type="checkbox"/> 6008A91	Jack for Drive Tire	\$ 272.00
Frame Color		\$
<input type="checkbox"/> 600C185	Midnight Black	NC
<input type="checkbox"/> 600C186	Lucent Lime	NC
<input type="checkbox"/> 600C187	Sunset Red	NC
<input type="checkbox"/> 600C188	Horizon Blue	NC
<input type="checkbox"/> 600C189	Dusk Purple	NC
<input type="checkbox"/> 600C190	Vivid Pink	NC
<input type="checkbox"/> 600C191	Desert Orange	NC
<input type="checkbox"/> 600C192	Camouflage - Woods	NC
<input type="checkbox"/> 600C193	Camouflage - Desert	NC
<input type="checkbox"/> 600C171	Carbon Fiber	NC
<input type="checkbox"/> 600C194	American Flag	NC
Seat To Floor Height		\$
Seat to Floor Heights start at 17" +/- 0.25"		
<input type="checkbox"/> 600HT44	Lowest Seat To Floor Height	NC
<input type="checkbox"/> 600HT45	Seat To Floor Height +1"	NC
<input type="checkbox"/> 600HT46	Seat To Floor Height +2"	NC
Batteries		\$
<input type="checkbox"/> 600AG25	Group 24 Gel Battery (Pair) - Installed	\$ 1,202.00
Battery Charger		\$
<input type="checkbox"/> 600BJ28	10A Battery Charger	NC

Seating: Seat Functions			
✓ OPTION		DESCRIPTION	RETAIL
No Power Seat Options			\$
<input type="checkbox"/>	600ST87	No Power Seating System with Power Functionality	\$ 9,328.00
Power Tilt			\$
Utilized without Power Recline.			
<input type="checkbox"/>	600PA64	50° Power Tilt in Space	\$ 10,236.00
Power Seat Elevate			\$
Only available when Power Tilt or Power Tilt and Recline is selected.			
<input type="checkbox"/>	600PI16	12" Adjustable Seat Elevator	\$ 8,917.00
Power Recline			\$
Utilized without Power Tilt. Requires No Power Seating with Power Functionality.			
<input type="checkbox"/>	600PA65	150° Power Recline Only	\$ 7,286.00
Power Tilt and Recline			\$
<input type="checkbox"/>	600PA66	50° Tilt & 150° Power Recline	\$ 16,165.00

Seating: Seat Sizes			
MPS (Captain) Seat Cushion (width x depth)			\$
MPS Backrest can be paired with the REHAB flat pan. If flat pan is desired, omit the MPS cushion and select the desired REHAB flat pan.			
<input type="checkbox"/>	600SI24	16" x 18" Vinyl Seat	NC
<input type="checkbox"/>	600SI25	18" x 18" Vinyl Seat	NC
<input type="checkbox"/>	600SI26	18" x 20" Vinyl Seat	NC
<input type="checkbox"/>	600SI27	20" x 18" Vinyl Seat	NC
<input type="checkbox"/>	600SI28	20" x 20" Vinyl Seat	NC
<input type="checkbox"/>	600SI29	22" x 20" Vinyl Seat	\$ 617.00
REHAB Flat Pan Seat Sizes - Adjustable in width and depth for most sizes			\$

REHAB Seat Pan

Standard seat pan available 13-23" wide. Deep seat pan available 13"-23" wide.

<input type="checkbox"/>	600PN13	Standard Adjustable Flat Seat Pan - 13" to 19" deep	NC
<input type="checkbox"/>	600PN14	Deep Adjustable Flat Seat Pan - 20" to 22" deep	\$ 621.00
<input type="checkbox"/>	600PN15	Deep Adjustable Flat Seat Pan - 23" to 24" deep	\$ 621.00

REHAB Seat Width

- 13"
 14"
 15"
 16"
 17"
 18"
 19"
 20"
 21"
 22"
 23"

REHAB Seat Depth

- 13"
 14"
 15"
 16"
 17"
 18"
 19"
 20"
 21"
 22"
 23"
 24"

REHAB Seating Cushion Options			
Note cushion size will be determined by Rehab Seat size selected above			
<input type="checkbox"/>	600CU59	JAY® Ion® (Sizes 16" - 21" Wide)	\$ 657.00
<input type="checkbox"/>	600CU60	JAY® Ion® (Size 22" Wide)	\$ 698.00
<input type="checkbox"/>	600CU61	JAY® Union® (Sizes 16" - 21" Wide)	\$ 764.00
<input type="checkbox"/>	600CU62	JAY® Union® (Size 22" Wide)	\$ 806.00
<input type="checkbox"/>	600CU64	JAY® Fusion® (Sizes 16" - 21" Wide)	\$ 1,045.00
<input type="checkbox"/>	600CU65	JAY® Fusion® (Size 22" Wide)	\$ 1,188.00

Seating: Seat Sizes

✓ OPTION	DESCRIPTION	RETAIL
Support Pads and Hardware (must pick pad and mount)		\$

Lateral Thigh Support Pads

<input type="checkbox"/> 600KF22	Lateral Thigh Support Pad Short - Right (5x5)	\$ 75.00
<input type="checkbox"/> 600KF23	Lateral Thigh Support Pad Short - Left (5x5)	\$ 75.00

Lateral Thigh Support Mounts

<input type="checkbox"/> 600KF26	Lateral Thigh Removable Mount - Right	\$ 325.00
<input type="checkbox"/> 600KF27	Lateral Thigh Removable Mount - Left	\$ 325.00

Seating: Seatback Options

MPS Backrest Options - Comes standard with Manual Recline	\$
---	----

Backrest Width

<input type="checkbox"/> 600BN46	Standard Backrest Width - 16-20" Wide	NC
<input type="checkbox"/> 600BN47	Wide Backrest Width - 22-23" Wide	\$ 452.00

Backrest Height

<input type="checkbox"/> 600B80	Standard Backrest Height - 22" Tall	NC
<input type="checkbox"/> 600B81	Tall Backrest Height - 26" Tall	\$ 744.00

Backrest Material

<input type="checkbox"/> 600BU106	Cloth Backrest	NC
<input type="checkbox"/> 600BU107	Vinyl Backrest	\$ 410.00

MPS Back Accessories	\$
----------------------	----

<input type="checkbox"/> 600PH33	Push Handle - Needed for O2 holder	\$ 211.00
<input type="checkbox"/> 6009A13	Accessory Mount - Needed for aftermarket headrest	\$ 159.00

REHAB Seating Backrest Options	\$
--------------------------------	----

Backrest will have the width of the MPS seat or flat pan. Please choose your back height. Backrest will come with 2 point hardware with a 500 lb. weight capacity.

JAY J3 Back Contour

<input type="checkbox"/> 600JP01	J3® PA Back for 16"-20" Wide - 2.2" of lateral thoracic support	\$ 912.00
<input type="checkbox"/> 600JP02	J3® PA Back for 22" Wide - 2.2" of lateral thoracic support	\$ 1,941.00
<input type="checkbox"/> 600JP11	J3® PL Back for 16"-20" Wide - 3" of lateral thoracic support	\$ 986.00
<input type="checkbox"/> 600JP13	J3® PD Back for 16"-20" Wide (18" and 21" shell height only) - 6" of lateral thoracic support	\$ 1,068.00

REHAB Seating Backrest Options	\$
--------------------------------	----

JAY® J3® Back Shell Options and Height Ranges

<input type="checkbox"/> 600JR14	16.5" tall shell (Upper Thoracic Short)	NC
<input type="checkbox"/> 600JR15	18" tall shell (Upper Thoracic Med)	NC
<input type="checkbox"/> 600JR16	19.5" tall shell (Upper Thoracic Tall)	NC
<input type="checkbox"/> 600JR17	21" tall shell (Shoulder Height)	NC
<input type="checkbox"/> 600JR27	22.5" tall shell (Shoulder Height)	NC
<input type="checkbox"/> 600JR28	24" tall shell (Shoulder Height)	NC

REHAB Flip Forward Backrest Option	\$
------------------------------------	----

<input type="checkbox"/> 600OD23	REHAB Flip Forward Backrest Option	\$ 214.00
----------------------------------	------------------------------------	-----------

Seating: Seatback Options

✓ OPTION	DESCRIPTION	RETAIL
JAY® Lateral Thoracic Swing-Away Supports (Requires Pad and Mount) - J3 back only		\$
JAY® Lateral Thoracic Support Pads		
<input type="checkbox"/> 600JW35	Lateral Thoracic (L x H) 4" x 5" Pad - Right	\$ 75.00
<input type="checkbox"/> 600JW38	Lateral Thoracic (L x H) 4" x 5" Pad - Left	\$ 75.00
<input type="checkbox"/> 600JW37	Lateral Thoracic (L x H) 6" x 4" Pad - Right	\$ 75.00
<input type="checkbox"/> 600JW40	Lateral Thoracic (L x H) 6" x 4" Pad - Left	\$ 75.00
JAY® Lateral Thoracic Swing-Away Mount		
<input type="checkbox"/> 600JW41	Lateral Thoracic Swing Away Bracket - Right	\$ 334.00
<input type="checkbox"/> 600JW42	Lateral Thoracic Swing Away Bracket - Left	\$ 334.00
JAY® Anterior Thoracic Support		\$
<input type="checkbox"/> 600ZR10	Bodypoint Monoflex Medium Chest Strap - Pad 3.5" x 18.5" (total length 49.5")	\$ 227.00
<input type="checkbox"/> 600ZR11	Bodypoint Monoflex Large Chest Strap - Pad 3.5" x 21.5" (total length 52.5")	\$ 227.00
<input type="checkbox"/> 600K10	JAY® Anterior Thoracic Support - Medium ATS8302	\$ 322.00
<input type="checkbox"/> 600K11	JAY® Anterior Thoracic Support - Large ATS8303	\$ 322.00
Pelvic Positioning Options - Select One		\$
<input type="checkbox"/> 600A848	Positioning Belt, Auto Style Buckle, 2" x 60"	NC
<input type="checkbox"/> 600A849	Positioning Belt, Aircraft Buckle, Padded 2" x 50"	\$ 103.00
<input type="checkbox"/> 600A850	Bodypoint - Evoflex Pelvic Stabilizer Belt	\$ 159.00
<input type="checkbox"/> 600A851	Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	\$ 272.00
Headrest Options - MPS Style		\$
<input type="checkbox"/> 600NU12	Cloth Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 552.00
<input type="checkbox"/> 600NU13	Vinyl Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 552.00
Headrest Options - REHAB Style		\$
Headrest Pad		
<input type="checkbox"/> 600HE39	Whitmyer® Plush Headrest Pad, 10"	\$ 377.00
<input type="checkbox"/> 600HE46	Whitmyer® Adjust-A-Plush Headrest Pad Medium (Adjustable from 7 to 16" Wide x 4" Tall)	\$ 445.00
Headrest Mount		
<input type="checkbox"/> 600MH33	Whitmyer® Linx Headrest Mount with Quick Release	\$ 247.00
<input type="checkbox"/> 600MH34	Whitmyer® Linx-2 Headrest Mount with Quick Release	\$ 298.00
Armrest Style Options		\$
<input type="checkbox"/> 600EG53	Flip-Up Removable Armrests (Post Style) - Attaches to base. Height adj 3" & angle adj 90°-155°	\$ 1,788.00
<input type="checkbox"/> 600EG54	Cantilever Armrests (Reclining Flip Back) - Attaches to back. Height adj 4" & angle adj 90°	\$ 2,539.00
<input type="checkbox"/> 600EG55	Cantilever 2 Point Quick Release Post - Attaches front of armrest to seat frame	\$ 3,167.00
Armrest Pad Options		\$
<input type="checkbox"/> 600AB81	10" Waterfall Armrest Pad - Right	NC
<input type="checkbox"/> 600AE81	10" Waterfall Armrest Pad - Left	NC
<input type="checkbox"/> 600AB82	14" Waterfall Armrest Pad - Right	NC
<input type="checkbox"/> 600AE82	14" Waterfall Armrest Pad - Left	NC

Legrest Options

✓ OPTION	DESCRIPTION	RETAIL
Center Mount Fixed Legrest		\$
Center Mount Fixed Legrest		
<input type="checkbox"/> 600H210	Center Mount (14.5" - 17" Leg Length) with Flip-Up Footplate (11" x 10")	\$ 610.00
Center Mount Fixed Legrest - Accessories		
<input type="checkbox"/> 600DI59	Fixed Center Mount Calf Pads	\$ 697.00
<input type="checkbox"/> 600LS16	Foot Strap - 1.5" wide strap	\$ 120.00
Center Mount Power Elevating Legrest		\$
Center Mount Power Elevating Legrest		
<input type="checkbox"/> 600H213	Center Mount Power Elevating - 2 Piece Footplate - Includes calf pads	\$ 5,583.00
Center Mount Power Elevating Legrest Extension		
<input type="checkbox"/> 600E42	Center Mount Power Elevating Extension - Standard (15" - 17" Leg Length)	NC
<input type="checkbox"/> 600E43	Center Mount Power Elevating Extension - Short (10" - 15" Leg Length)	NC
Center Mount Power Elevating Legrest Footplates		
<input type="checkbox"/> 600F117	Center Mount Power Elevating Footplate - Standard (5" x 12" each)	NC
<input type="checkbox"/> 600F118	Center Mount Power Elevating Footplate - Short (4.5" x 7" each)	NC
Swing-Away Legrest		\$
Swing-Away Legrest		
<input type="checkbox"/> 600H211	Swing-Away Legrest (12" - 17" Leg Length) with flip-up footplates (6" x 8.5")	\$ 877.00
Swing-Away Power Elevating Legrest		\$
Swing-Away Power Elevating Legrest		
<input type="checkbox"/> 600H214	Swing-Away Power Elevating Legrest (12" - 17" Leg Length) with flip-up footplates (6" x 8.5")	\$ 6,094.00
S/A Power Legrest Positioning Options		
<input type="checkbox"/> 600HL15	Heel Loops - Attaches to the footplate. Length 11.5-13" and Width 1"	\$ 171.00
Accessories		\$
<input type="checkbox"/> 6008A93	Waterproof Wheelchair Cover	\$ 222.00
<input type="checkbox"/> 6008A94	Waterproof Joystick Cover	\$ 69.00
<input type="checkbox"/> 6008A95	Medical Necessity Bag - Velcro's to armrest (10.5" x 6" x 1.25")	\$ 51.00
<input type="checkbox"/> 600A221	Deluxe Backpack	\$ 144.00
<input type="checkbox"/> 600A615	Cup Holder - Self-leveling. Attaches to armrest opposite joystick	\$ 287.00
<input type="checkbox"/> 600A719	Articulating Oxygen Holder for MPS Seating System	\$ 1,492.00
<input type="checkbox"/> 600A612	Cane Holder for MPS Seating System - Includes push handles	\$ 957.00
<input type="checkbox"/> 6008A96	Luggage Rack for 360 - Attaches to base behind the chair (10.5" x 13")	\$ 499.00
<input type="checkbox"/> 6008A97	Charge Mate - 12 Volt/Car outlet and USB/cell phone charger with max draw of 6 amps	\$ 803.00
<input type="checkbox"/> 600A613	Cell Phone Holder - Adjustable to fit most cell phones	\$ 445.00
<input type="checkbox"/> 600A614	Fishing Rod Holder	\$ 445.00
Be Adaptive Products		
<input type="checkbox"/> 6008A98	Be Adaptive Limited Mobility Shooting Mount - Gun or crossbow	\$ 2,242.00
<input type="checkbox"/> 6008A99	Be Adaptive Trigger Mechanism	\$ 1,242.00
<input type="checkbox"/> 6008A100	Be Adaptive Limited Mobility Fishing Mount	\$ 2,242.00

Electronics

✓ OPTION	DESCRIPTION			RETAIL				
Joystick				\$				
<input type="checkbox"/> 600J121	Standard LED Joystick (w/o Screen Display)			NC				
<input type="checkbox"/> 600J116	Color Joystick (2.3" Color LCD Screen Display, CJSM1)			NC				
<input type="checkbox"/> 600J118	High Resolution Color Joystick (CJSM2+BT) 3.5" Color LCD Display			\$ 757.00				
Joystick Mount				\$				
<input type="checkbox"/> 600JM34	Fixed on Armrest - Depth adjustable 8.25"			NC				
<input type="checkbox"/> 600JM45	Swing-Away Joystick Mount - Includes latch to lock in the forward position			\$ 849.00				
<input type="checkbox"/> 600JM46	Midline Joystick Mount - 2 position mount (Rotation 0° - 180°)			\$ 2,050.00				
Joystick Mount - Location				\$				
<input type="checkbox"/> 600HM16	Right Hand Armrest			NC				
<input type="checkbox"/> 600HM17	Left Hand Armrest			NC				
Display Language Selection R-net				\$				
<input type="checkbox"/> 600L1	English	NC	<input type="checkbox"/> 600L3	French	NC	<input type="checkbox"/> 600L7	Spanish	NC
Joystick Handle Type				\$				
Comes standard with PG Cone-Shaped Joystick Handle								
<input type="checkbox"/> 600JH13	Ball Handle			\$ 110.00				
<input type="checkbox"/> 600JH26	Bodypoint Grooved Mushroom Joystick Handle			\$ 256.00				
<input type="checkbox"/> 600JH24	Bodypoint – Rubber Dome Joystick Handle			\$ 256.00				
<input type="checkbox"/> 600JH25	Bodypoint – U-Shaped Joystick Handle (4 inch)			\$ 240.00				
Expandable Power Module				\$				
Must select both items								
<input type="checkbox"/> 600EX05	Expandable 120A Controller - R-net			\$ 1,099.00				
<input type="checkbox"/> 600EX06	Harness for 120A Expandable Controller			\$ 768.00				
Power Seat Function				\$				
<input type="checkbox"/> 600EO37	Thru Drive Control 1 Actuator - Must be selected if ordering 1 actuator function			\$ 1,443.00				
<input type="checkbox"/> 600EO38	Thru Drive Control 2+ Actuators - Must be selected if ordering 2 or more actuator functions			\$ 2,903.00				
Power Seat Control Option				\$				
<input type="checkbox"/> 600QG11	Ctrl+5 Switch Box with Free Assignable Functions			\$ 1,017.00				
<input type="checkbox"/> 600QG12	Ctrl+5 Switch Box with Free Assignable Functions (Multi-power only)			NC				
Ctrl+5 Switch Type				\$				
<input type="checkbox"/> 600EY22	Push Button			NC				
<input type="checkbox"/> 600EY23	Toggle Switches			\$ 356.00				
Power Seat Control Option - Ctrl+5 Mounting location				\$				
<input type="checkbox"/> 600EY19	Right Armrest	NC	<input type="checkbox"/> 600EY20	Left Armrest	NC			
Attendant Controls				\$				
<input type="checkbox"/> 6000A05	Dual Control Module for Attendant - Includes accessory mount and push handle			\$ 2,460.00				
Attendant Control Mounting Location				\$				
<input type="checkbox"/> 600DN5	Right Back Cane	NC	<input type="checkbox"/> 600DN6	Left Back Cane	NC			
R-net Omni2 - 3.5" LCD Display w/ Built-In Infrared & Bluetooth				\$				
<input type="checkbox"/> 600ON14	R-net Omni2 Only Input Control Module			\$ 3,611.00				
<input type="checkbox"/> 600ON15	R-net Omni2 Bundled Input Control Module Available w/ Switch-It Input Devices			NC				
R-net Omni2 - Mounting Location				\$				
<input type="checkbox"/> 600DS23	Right Hand Armrest	NC	<input type="checkbox"/> 600DS24	Left Hand Armrest	NC			
Head Control Packages				\$				
<input type="checkbox"/> 600EU25	Switch-It DualPro® Proportional Head Array - Facial Pads			\$ 7,057.00				
<input type="checkbox"/> 600EU26	Switch-It DualPro® Proportional Head Array - Spot Pads			\$ 7,057.00				