



ORDER FORM

# Mygo Stander

# EILK02

Stand out from the crowd in this modular 3-in-1 adjustable standing system with support for hip, knee & ankle contractures

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

## Mygo Stander

Size 1 Weight Capacity: 110 lbs | Size 2 Weight Capacity: 132 lbs

✓ OPTION	DESCRIPTION	RETAIL
<b>Packages</b>		<b>\$</b>
Package includes Mygo Stander, Stander Covers, Shoulder Pad and Cover, Contoured Headrest Hardware and Cushion - Black, Fixed Chest Laterals, Chest Wraparound Harness, Supine Derotation Hip Belt, Knee Cup Height Extension Brackets, Basic Knee Straps, Tray, Sandals with Straps. Substitutions are not available for Packages.		
<input type="checkbox"/> 141QE11	Size 1 Mygo Stander Package - Blue (Includes Medium Sandals)	\$8,200
<input type="checkbox"/> 141QE12	Size 1 Mygo Stander Package - Pink (Includes Medium Sandals)	\$8,200
<input type="checkbox"/> 141QE13	Size 2 Mygo Stander Package - Blue (Includes Large Sandals)	\$8,327
<input type="checkbox"/> 141QE14	Size 2 Mygo Stander Package - Pink (Includes Large Sandals)	\$8,327
<b>Mygo Stander</b>		<b>\$</b>

Selection Required.

Includes manual angle adjustment, polyurethane chest pad, hip pad and knee supports, two piece footplate and chassis.

<input type="checkbox"/> <b>141NB112</b>	<b>Size 1 Mygo Stander</b>	<b>\$4,563</b>
143-1600-1	Mygo Stander Frame	
143-1600-2	Pivot Base	
143-1600-3	Chest Pad Assembly	
143-1600-4	Hip Assembly	
143-1600-5	Knee Supports	
143-1600-6	Positioning Pad Covers	
143-1600-7	Swivel Casters	
143-1600-8	Two Piece Footplate	
<input type="checkbox"/> <b>141NB113</b>	<b>Size 2 Mygo Stander</b>	<b>\$4,563</b>
143-2600-1	Mygo Stander Frame	
143-2600-2	Pivot Base	
143-2600-3	Chest Pad	
143-2600-4	Hip Pad	
143-2600-5	Knee Supports	
143-2600-6	Positioning Pad Covers	
143-2600-7	Swivel Casters	
143-2600-8	Two Piece Footplate	

<b>Cover</b>	<b>\$</b>
--------------	-----------

Selection Required. Must match Stander size.

Includes chest and hip pad covers and knee support covers.

<input type="checkbox"/> 141IU03	Size 1 Mygo Stander Covers - Blue	\$434
<input type="checkbox"/> 141IU04	Size 1 Mygo Stander Covers - Orange	\$434
<input type="checkbox"/> 141IU05	Size 1 Mygo Stander Covers - Pink	\$434
<input type="checkbox"/> 141IU06	Size 2 Mygo Stander Covers - Blue	\$434
<input type="checkbox"/> 141IU07	Size 2 Mygo Stander Covers - Orange	\$434
<input type="checkbox"/> 141IU08	Size 2 Mygo Stander Covers - Pink	\$434

<b>Shoulder Support</b>	<b>\$</b>
-------------------------	-----------

Cover selection required with Shoulder Pad. Shoulder Support requires Head Support.

Shoulder support pad must be selected for supine use and head support use.

<input type="checkbox"/> 141SS5	Mygo Stander Shoulder Pad	\$283
<input type="checkbox"/> 141SS6	Mygo Stander Shoulder Pad Cover - Blue	\$120
<input type="checkbox"/> 141SS7	Mygo Stander Shoulder Pad Cover - Orange	\$120
<input type="checkbox"/> 141SS8	Mygo Stander Shoulder Pad Cover - Pink	\$120
<input type="checkbox"/>	No Cover	N/C

## Mygo Stander

✓ OPTION	DESCRIPTION	RETAIL
----------	-------------	--------

### Headrest Support Hardware \$

Selection of Headrest Support Hardware Requires Headrest Cushion

Head Support selection requires Shoulder Support.

<input type="checkbox"/>	141AW24	Contoured Headrest Support Hardware	\$396
<input type="checkbox"/>	141AW22	Flat Headrest Support Hardware	\$247

### Headrest Cushion \$

Headrest Laterals only available with Flat Headrest

<input type="checkbox"/>	141Y030	Contoured Headrest Cushion Infection Control - Black Vinyl Includes Hardware	\$406
<input type="checkbox"/>	141Y027	Contoured Headrest Cushion - Black	\$188
<input type="checkbox"/>	141Y028	Flat Headrest Cushion - Black	\$124
<input type="checkbox"/>	141Y031	Flat Headrest Laterals with Covers - Black (Pair)	\$369

### Wraparound Harness \$

Harnesses may be used as chest or hip harnesses. If used as chest harness, chest laterals are required. If used as hip harness, hip laterals are required.

<input type="checkbox"/>	141NH42	QTY	Size 1 Chest/Hip Wraparound Harness - Blue	\$287
<input type="checkbox"/>	141NH43	QTY	Size 1 Chest/Hip Wraparound Harness - Orange	\$287
<input type="checkbox"/>	141NH86	QTY	Size 1 Chest/Hip Wraparound Harness - Pink	\$287
<input type="checkbox"/>	141NH87	QTY	Size 1 Chest/Hip Wraparound Harness Infection Control - Black	\$287
<input type="checkbox"/>	141NH88	QTY	Size 2 Chest/Hip Wraparound Harness - Blue	\$306
<input type="checkbox"/>	141NH89	QTY	Size 2 Chest/Hip Wraparound Harness - Orange	\$306
<input type="checkbox"/>	141NH99	QTY	Size 2 Chest/Hip Wraparound Harness - Pink	\$306
<input type="checkbox"/>	141NH100	QTY	Size 2 Chest/Hip Wraparound Harness Infection Control - Black	\$306
<input type="checkbox"/>	141NH101	QTY	Size 3 Chest/Hip Wraparound Harness - Blue	\$296
<input type="checkbox"/>	141NH102	QTY	Size 3 Chest/Hip Wraparound Harness - Orange	\$296
<input type="checkbox"/>	141NH103	QTY	Size 3 Chest/Hip Wraparound Harness - Pink	\$296
<input type="checkbox"/>	141NH104	QTY	Size 3 Chest/Hip Wraparound Harness Infection Control - Black	\$296

### Lateral Support \$

Chest Laterals required if ordering Wraparound Harness to be used as chest harness.

Hip Laterals required if ordering Wraparound Harness to be used as hip harness.

<input type="checkbox"/>	141ZY16	Size 1 Fixed Chest Laterals (Pair)	\$176
<input type="checkbox"/>	141ZY17	Size 2 Fixed Chest Laterals (Pair)	\$176
<input type="checkbox"/>	141ZY18	Size 1 and 2 Flip-Away Chest Laterals (Pair)	\$557
<input type="checkbox"/>	141KC10	Hip Laterals (Pair) Not available with Derotation Belt	\$176

### Hip Infill \$

<input type="checkbox"/>	141IU09	Hip Infill Pad	\$348
<input type="checkbox"/>	141IU10	Size 1 Hip Infill Pad Cover - Blue	\$214
<input type="checkbox"/>	141IU11	Size 1 Hip Infill Pad Cover - Orange	\$214
<input type="checkbox"/>	141IU12	Size 1 Hip Infill Pad Cover - Pink	\$214
<input type="checkbox"/>	141IU13	Size 2 Hip Infill Pad Cover - Blue	\$214
<input type="checkbox"/>	141IU14	Size 2 Hip Infill Pad Cover - Orange	\$214
<input type="checkbox"/>	141IU15	Size 2 Hip Infill Pad Cover - Pink	\$214

## Mygo Stander

✓ OPTION	DESCRIPTION	RETAIL
<b>Hip Belt</b>		<b>\$</b>
Not available with Hip Laterals or Hip Wraparound Harness		
<input type="checkbox"/> 141A858	Size 1 Prone Derotation Hip Belt with Hardware - Blue	\$326
<input type="checkbox"/> 141A859	Size 1 Prone Derotation Hip Belt with Hardware - Orange	\$326
<input type="checkbox"/> 141A860	Size 1 Prone Derotation Hip Belt with Hardware - Pink	\$326
<input type="checkbox"/> 141A861	Size 2 Prone Derotation Hip Belt with Hardware - Blue	\$369
<input type="checkbox"/> 141A862	Size 2 Prone Derotation Hip Belt with Hardware - Orange	\$369
<input type="checkbox"/> 141A863	Size 2 Prone Derotation Hip Belt with Hardware - Pink	\$369
<input type="checkbox"/> 141A864	Size 3 Prone Derotation Hip Belt with Hardware - Blue	\$414
<input type="checkbox"/> 141A865	Size 3 Prone Derotation Hip Belt with Hardware - Orange	\$414
<input type="checkbox"/> 141A866	Size 3 Prone Derotation Hip Belt with Hardware - Pink	\$414
<input type="checkbox"/> 141A867	Size 1 Supine Derotation Hip Belt with Hardware - Blue	\$326
<input type="checkbox"/> 141A868	Size 1 Supine Derotation Hip Belt with Hardware - Orange	\$326
<input type="checkbox"/> 141A869	Size 1 Supine Derotation Hip Belt with Hardware - Pink	\$326
<input type="checkbox"/> 141A870	Size 2 Supine Derotation Hip Belt with Hardware - Blue	\$369
<input type="checkbox"/> 141A871	Size 2 Supine Derotation Hip Belt with Hardware - Orange	\$369
<input type="checkbox"/> 141A872	Size 2 Supine Derotation Hip Belt with Hardware - Pink	\$369
<input type="checkbox"/> 141A873	Size 3 Supine Derotation Hip Belt with Hardware - Blue	\$414
<input type="checkbox"/> 141A874	Size 3 Supine Derotation Hip Belt with Hardware - Orange	\$414
<input type="checkbox"/> 141A875	Size 3 Supine Derotation Hip Belt with Hardware - Pink	\$414
<b>Knee</b>		<b>\$</b>
Knee Cup Height Extension Brackets increase distance from footplate to top of knee by 4"		
<input type="checkbox"/> 141JG13	Knee Cup Height Extension Brackets (Pair)	\$265
<input type="checkbox"/> 141XX20	Size 1 Basic Knee Straps - Blue (Pair)	\$197
<input type="checkbox"/> 141XX21	Size 1 Basic Knee Straps - Orange (Pair)	\$197
<input type="checkbox"/> 141XX22	Size 1 Basic Knee Straps - Pink (Pair)	\$197
<input type="checkbox"/> 141XX23	Size 1 Basic Knee Straps Infection Control - Black (Pair)	\$197
<input type="checkbox"/> 141XX24	Size 2 Basic Knee Straps - Blue (Pair)	\$240
<input type="checkbox"/> 141XX25	Size 2 Basic Knee Straps - Orange (Pair)	\$240
<input type="checkbox"/> 141XX26	Size 2 Basic Knee Straps - Pink (Pair)	\$240
<input type="checkbox"/> 141XX27	Size 2 Basic Knee Straps Infection Control - Black (Pair)	\$240
<input type="checkbox"/> 141XX28	Split Knee Straps - Blue (Pair)	\$369
<input type="checkbox"/> 141XX29	Split Knee Straps - Orange (Pair)	\$369
<input type="checkbox"/> 141XX30	Split Knee Straps - Pink (Pair)	\$369
<b>Accessories</b>		<b>\$</b>
<input type="checkbox"/> 141A211	Size 1 Mygo Stander Tray	\$674
<input type="checkbox"/> 141A212	Size 2 Mygo Stander Tray	\$682
<input type="checkbox"/> 141A213	Size 1 Mygo Stander Tray with Bowl Insert	\$695
<input type="checkbox"/> 141A214	Size 2 Mygo Stander Tray with Bowl Insert	\$704
<input type="checkbox"/> 141A215	Grab Rail / Toy Bar for Tray	\$306
<input type="checkbox"/> 1418A107	Small Sandals with Straps (Pair)	\$273
<input type="checkbox"/> 1418A108	Medium Sandals with Straps (Pair)	\$291
<input type="checkbox"/> 1418A109	Large Sandals with Straps (Pair)	\$315
<input type="checkbox"/> 1419A23	Sandal Raiser 1/2" (Each) Requires Sandal Selection.	\$140

### Technical Specifications - Imperial

Frame Sizing Specification	Size 1	Size 2
Age (approximate)	4 - 10 years	8 - 14 years
User Height (Min / Max)	42.5" / 57.1"	53.1" / 66.9"
Maximum User Weight	110 lbs	132 lbs
Chest Support Height	29.9" - 42.5"	38.6" - 51.2"
Chest Width	6.3" - 11"	7.9" - 12.6"
Hip Support Height	18.1" - 28.7"	24.4" - 34.4"
Knee Cup Height	11" - 15"	12.25" - 16.25"
Knee Cup Height with Extension Bracket	15" - 19"	16.25" - 20.25"
Hip Width	8.7" - 13"	11" - 15.4"
Footplate Angle: Plantarflexion / Dorsiflexion	10° / 10°	10° / 10°
Maximum Knee Contracture	25°	25°
Tray Height Range	0" - 6"	0" - 6"
Tray Angle Adjustment: Prone / Supine	45° / 45°	45° / 45°
Product Weight	48 lbs	55 lbs
Product Angle Range: Prone / Supine	0° - 80°	0° - 80°
Stander Footprint (length x width)	38.5" x 24.8"	44.5" x 25.6"

### Technical Specifications - Metric

Frame Sizing Specification	Size 1	Size 2
Age (approximate)	4 - 10 years	8 - 14 years
User Height (Min / Max)	1080 mm / 1450 mm	1350 mm / 1700 mm
Maximum User Weight	50 kgs	60 kgs
Chest Support Height	760 mm - 1080 mm	980 mm - 1300 mm
Chest Width	160 mm - 280 mm	200 mm - 320 mm
Hip Support Height	460 mm - 730 mm	620 mm - 875 mm
Knee Cup Height	280 mm - 380 mm	310mm - 410 mm
Knee Cup Height with Extension Bracket	381 mm - 481 mm	411 mm - 511 mm
Hip Width	220 mm - 330 mm	280 mm - 390 mm
Footplate Angle: Plantarflexion / Dorsiflexion	10° / 10°	10° / 10°
Maximum Knee Contracture	25°	25°
Tray Height Range	0 mm - 150 mm	0 mm - 150 mm
Tray Angle Adjustment: Prone / Supine	45° / 45°	45° / 45°
Product Weight	22 kgs	25 kgs
Product Angle Range: Prone / Supine	0° - 80°	0° - 80°
Stander Footprint (length x width)	978 mm - 630 mm	1130 mm x 650 mm

Sunrise Medical Canada Inc.  
1000 Creditstone Road, Unit 2  
Concord, ON L4K 4P8, Canada



MK-129850-CAN Rev. 14 © 05/2026

