



JAY® Fit
Backrest and Seat Assembly
May 2025

Mark For:

Date: _____
Dealer Acct #: _____
Dealer: _____
Dealer Contact: _____
Dealer Address: _____
Dealer City: _____ Prov: _____ PC: _____
Dealer Phone: () _____ Fax: () _____
Confirmation Email: _____
Confirm Via: ☐ Fax ☐ Email

Submitting for:

☐ Quote ☐ Order


PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
Attention: _____
Address: _____
Address: _____
Ship To City: _____ Prov: _____ PC: _____
Ship To Phone: () _____ Fax: () _____

CLIENT MEASUREMENTS

First Name	
Age	
Diagnosis	
Height	
Weight	
1. Top of Head to Seat	
2. Top of Shoulder to Seat	
3. Under Arm to Seat	
4. Trunk Depth	
5. Actual Sitting Depth	
6. Foot Length	
7. Width of Head	
8. Shoulder Width	
9. Chest Width	
10. Hip Width	
11. Left Lower Leg Length	
12. Right Lower Leg Length	



Client Wheelchair Description

Model	
Back Height	
Combine with Wheelchair	☐ See Quickie/Zippie OF ☐ Existing Chair Serial #
Chair Width	
Chair Depth	

NOTES

JAY Fit Back

☐ JFITBXXYY | JAY Fit Backrest All Sizes | \$512.00

NOTE: Select appropriate assembly number below.

250 lbs Weight Capacity

EACH BACK SHIPS WITH THE FOLLOWING:

- Choice of mounting hardware - TH Hardware Included with Back
- Choice of backrest cover - 3DX with Black Stretch Standard
- Colored backrest shell of your choice
- 1 Small and 1 Large Lumbar Support



BACKREST WIDTH	10"	12"	14"	16"	18"
FITS BACK CANE WIDTHS	10" - 12"	12" - 14"	14" - 16"	16" - 18"	18" - 20"
ACTUAL SHELL WIDTH	9"	11"	13"	15"	17"
DEPTH OF CONTOUR	1-1/4"	1-1/4"	1-1/4"	1-1/4"	1-1/4"

		Backrest Width				
		10	12	14	16	18
Backrest Height	10"	☐ JFITB1010				
	12"	☐ JFITB1012	☐ JFITB1212	☐ JFITB1412		
	14"	☐ JFITB1014	☐ JFITB1214	☐ JFITB1414	☐ JFITB1614	☐ JFITB1814
	16"		☐ JFITB1216	☐ JFITB1416	☐ JFITB1616	☐ JFITB1816
	18"		☐ JFITB1218	☐ JFITB1418	☐ JFITB1618	☐ JFITB1818
	20"			☐ JFITB1420	☐ JFITB1620	☐ JFITB1820
	22"				☐ JFITB1622	☐ JFITB1822

The JAY Fit Backrest has been dynamically tested and meets all criteria for use in motor vehicles according to Section 5.1 of ISO 16840-4 and Section 20, ANSI/RESNA WC/Volume 4.

Please visit www.sunparts.us for spare parts information.

JAY Fit Back - Mounting Hardware

Note: J3 TH Hardware is include with the back but if Snaptite or Adj Quick Mount Hardware is selected J3 TH hardware will automatically be omitted.

☐ Included in Back | J3 TH Hardware | NC

☐ JYW-JFTBHDW-OPT02 | Snaptite Hardware | \$77.00

Note: Includes collars to mount to 3/4", 7/8", and 1" tubes

☐ JYW-JFTBHDW-OPT03 | Adjustable Quick Mount | \$77.00

Note: Includes collars to mount to 3/4", 7/8", and 1" tubes

JAY Fit Back - Back Cover

Note: 3DX with Black Stretch cover is include with the back but if Reverse Dartex is selected the 3DX with Black Stretch Cover will automatically be omitted.

☐ Included in Back | 3DX with Black Stretch | NC

☐ JYW-JFTBCV-OPT02 | Reverse Dartex | \$74.00

JAY Fit Back - Back Cover Embroidery (Optional)

Embroidery - 12 Letter Limit

☐ JYW-JFTBCE-OPT01 | Script/Cursive | NC

☐ JYW-JFTBCE-OPT02 | Block | NC

Name:

Thread Color:

JAY Fit Back - Backrest Spinal Fluid Pad (Optional)

☐ JYW-JFTBSFP-OPT01 | Spinal Fluid Pad* |

*Sewn into cover

JAY Fit Back - Lateral/Posterior Pelvic Support Kit (Optional)

Kit includes 1 Narrow & 1 wide posterior contour strip, and 1 pair of lateral pyramids

☐ JYW-JFTBPSK-OPT01 | Sml Lateral/Posterior Plvc Support Kit | \$158.00

☐ JYW-JFTBPSK-OPT02 | Lrg Lateral/Posterior Plvc Support Kit | \$158.00

Additional Narrow Posterior Contour Strips

☐ JYW-JFTBAPCS-OPT01 | 1 Narrow (1-3/4" Wide) | \$35.00

☐ JYW-JFTBAPCS-OPT02 | 2 Narrow (1-3/4" Wide) | \$69.00

☐ JYW-JFTBAPCS-OPT03 | 3 Narrow (1-3/4" Wide) | \$103.00

Additional Wide Posterior Contour Strips

☐ JYW-JFTBAPCS-OPT04 | 1 Wide (2-1/4" Wide) | \$35.00

☐ JYW-JFTBAPCS-OPT05 | 2 Wide (2-1/4" Wide) | \$69.00

☐ JYW-JFTBAPCS-OPT06 | 3 Wide (2-1/4" Wide) | \$103.00

JAY Fit - Painted Back Shell (Shell to come black if not otherwise specified) - Required

☐	JYW-JFTBPBS-OPT27	Celestial Sapphire	NC	☐	JYW-JFTBPBS-OPT37	Salute Green	NC
☐	JYW-JFTBPBS-OPT28	Desert Sage	NC	☐	JYW-JFTBPBS-OPT38	Scarlet Noir	NC
☐	JYW-JFTBPBS-OPT29	Desert Tan	NC	☐	JYW-JFTBPBS-OPT39	Silver Ore	NC
☐	JYW-JFTBPBS-OPT25	Gloss Black	NC	☐	JYW-JFTBPBS-OPT40	Solar Orange	NC
☐	JYW-JFTBPBS-OPT26	Matte Black	NC	☐	JYW-JFTBPBS-OPT41	Starlight White	NC
☐	JYW-JFTBPBS-OPT30	Metallic Shadow	NC	☐	JYW-JFTBPBS-OPT42	Storm Gray	NC
☐	JYW-JFTBPBS-OPT31	Midnight Blue	NC	☐	JYW-JFTBPBS-OPT43	Turbo Green	NC
☐	JYW-JFTBPBS-OPT32	Mystic Purple	NC	☐	JYW-JFTBPBS-OPT44	Velocity Yellow	NC
☐	JYW-JFTBPBS-OPT33	Petal Pink	NC				
☐	JYW-JFTBPBS-OPT34	Pink Power	NC				
☐	JYW-JFTBPBS-OPT35	Red Pepper	NC				
☐	JYW-JFTBPBS-OPT36	Royal Blue	NC				

Whitmyer Headrests

STEP 1 - Select your WHITMYER HEADREST POST

Please select one box from the below **WHITMYER Headrest Mount & Bracket** choices. Note: Headrest post includes post, removable mount, and mounting hardware.


☐	JAXYS	WHITMYER - Alys Mount & Bracket	\$358.00	☐	JONYX	WHITMYER - Onyx Mount & Bracket	\$418.00
☐	JLINX	WHITMYER - Linx Mount & Bracket	\$477.00	☐	JCOBRA	WHITMYER - Cobra Mount & Bracket	\$418.00
☐	JLINX2	WHITMYER - Linx ² Mount & Bracket	\$477.00	☐	JCOBRAFB	WHITMYER - Cobra FB Mount & Bracket	\$418.00




STEP 2 - Select your WHITMYER HEADREST PAD

Please select one box from the below **WHITMYER Headrest Pad** choices.

☐	14	JPLSH14 - 14" Plush Pad	\$210.00	☐	AAPS	JAAPS - Adjust-A-Plush Standard (Lg)	\$358.00
☐	10	JPLSH10 - 10" Plush Pad	\$210.00	☐	AAPM	JAAPM - Adjust-A-Plush Medium	\$382.00
☐	8	JPLSH8 - 8" Plush Pad	\$330.00	☐	AAPN	JAAPN - Adjust-A-Plush Narrow (Sm)	\$358.00
☐	HUPL	JHUPL - HeadsUp Lrg: 5.3" Lng Arms	\$358.00	☐	CL	JCRDLLG - Lrg Contour Cradle Pad	\$298.00
☐	HUPM	JHUPM - HeadsUp Med: 4.1" Sht Arms	\$358.00	☐	CM	JCRDLMD - Med Contour Cradle Pad	\$298.00
☐	HUPS	JHUPS - HeadsUp Sml: 3.5" Sht Arms	\$358.00				




PLUSH



Heads Up



Adjust-A-Plush



Contour Cradle

STEP 1 - Headrest Post Example = JAXYS	+	STEP 2 - Headrest Pad # Example = 14	=	Completed JAY Fit Headrest # Example = JAXYS14
--	---	--	---	--

HEADREST OPTIONS

☐	JUHM	Universal Headrest Mounting Plate Only	\$68.00
☐	JSAHP	Swing Away Headrest Plate	\$74.00
☐	JDFSFA	Dynamic Forehead Strap Adult	\$256.00
☐	JDFSFP	Dynamic Forehead Strap Peds	\$256.00

Note: Universal Headrest Mounting Plate utilized for mounting non Whitmyer Headrest.

JAY Fit Laterals

FIXED LATERAL THORACIC SUPPORT

Includes fixed mounting hardware. Left and Right lateral are independent and different types and sizes can be selected for each side.

Left Side				Right Side			
☐	JFTLAT4X4FL	Left-LATERAL FIXED 4Lx4H	\$172.00	☐	JFTLAT4X4FR	Right-LATERAL FIXED 4Lx4H	\$172.00
☐	JFTLAT4X5FL	Left-LATERAL FIXED 4Lx5H	\$172.00	☐	JFTLAT4X5FR	Right-LATERAL FIXED 4Lx5H	\$172.00
☐	JFTLAT4X6FL	Left-LATERAL FIXED 4Lx6H	\$172.00	☐	JFTLAT4X6FR	Right-LATERAL FIXED 4Lx6H	\$172.00
☐	JFTLAT6X4FL	Left-LATERAL FIXED 6Lx4H	\$172.00	☐	JFTLAT6X4FR	Right-LATERAL FIXED 6Lx4H	\$172.00
☐	JFTLAT6X5FL	Left-LATERAL FIXED 6Lx5H	\$172.00	☐	JFTLAT6X5FR	Right-LATERAL FIXED 6Lx5H	\$172.00
☐	JFTLAT6X6FL	Left-LATERAL FIXED 6Lx6H	\$172.00	☐	JFTLAT6X6FR	Right-LATERAL FIXED 6Lx6H	\$172.00

FIXED LATERAL THORACIC SUPPORT with 1" INSET

Includes fixed mounting hardware. Left and Right lateral are independent and different types and sizes can be selected for each side.


Left Side				Right Side			
☐	JFTLAT4X4FL1I	Left-LATERAL FIXED 1"Inset 4Lx4H	\$172.00	☐	JFTLAT4X4FR1I	Right-LATERAL FIXED 1"Inset 4Lx4H	\$172.00
☐	JFTLAT4X5FL1I	Left-LATERAL FIXED 1"Inset 4Lx5H	\$172.00	☐	JFTLAT4X5FR1I	Right-LATERAL FIXED 1"Inset 4Lx5H	\$172.00
☐	JFTLAT4X6FL1I	Left-LATERAL FIXED 1"Inset 4Lx6H	\$172.00	☐	JFTLAT4X6FR1I	Right-LATERAL FIXED 1"Inset 4Lx6H	\$172.00
☐	JFTLAT6X4FL1I	Left-LATERAL FIXED 1"Inset 6Lx4H	\$172.00	☐	JFTLAT6X4FR1I	Right-LATERAL FIXED 1"Inset 6Lx4H	\$172.00
☐	JFTLAT6X5FL1I	Left-LATERAL FIXED 1"Inset 6Lx5H	\$172.00	☐	JFTLAT6X5FR1I	Right-LATERAL FIXED 1"Inset 6Lx5H	\$172.00
☐	JFTLAT6X6FL1I	Left-LATERAL FIXED 1"Inset 6Lx6H	\$172.00	☐	JFTLAT6X6FR1I	Right-LATERAL FIXED 1"Inset 6Lx6H	\$172.00

SWING-AWAY LATERAL THORACIC SUPPORT

Includes swing away mounting hardware. Left and Right lateral are independent and different types and sizes can be selected for each side.

Left Side				Right Side			
☐	JFTLAT4X4SAL	LAT 4L X 4H W/ S/A MOUNT	\$334.00	☐	JFTLAT4X4SAR	LAT 4L X 4H W/ S/A MOUNT	\$334.00
	JFTLAT4X4	LATERAL 4L X 4H	\$172.00		JFTLAT4X4	LATERAL 4L X 4H	\$172.00
	JFTSAMNT	LATERAL SWING AWAY MOUNT	\$162.00		JFTSAMNT	LATERAL SWING AWAY MOUNT	\$162.00
☐	JFTLAT4X5SAL	LAT 4L X 5H W/ S/A MOUNT	\$334.00	☐	JFTLAT4X5SAR	LAT 4L X 5H W/ S/A MOUNT	\$334.00
	JFTLAT4X5	LATERAL 4L X 5H	\$172.00		JFTLAT4X5	LATERAL 4L X 5H	\$172.00
	JFTSAMNT	LATERAL SWING AWAY MOUNT	\$162.00		JFTSAMNT	LATERAL SWING AWAY MOUNT	\$162.00
☐	JFTLAT4X6SAL	LAT 4L X 6H W/ S/A MOUNT	\$334.00	☐	JFTLAT4X6SAR	LAT 4L X 6H W/ S/A MOUNT	\$334.00
	JFTLAT4X6	LATERAL 4L X 6H	\$172.00		JFTLAT4X6	LATERAL 4L X 6H	\$172.00
	JFTSAMNT	LATERAL SWING AWAY MOUNT	\$162.00		JFTSAMNT	LATERAL SWING AWAY MOUNT	\$162.00
☐	JFTLAT6X4SAL	LAT 6L X 4H W/ S/A MOUNT	\$334.00	☐	JFTLAT6X4SAR	LAT 6L X 4H W/ S/A MOUNT	\$334.00
	JFTLAT6X4	LATERAL 6L X 4H	\$172.00		JFTLAT6X4	LATERAL 6L X 4H	\$172.00
	JFTSAMNT	LATERAL SWING AWAY MOUNT	\$162.00		JFTSAMNT	LATERAL SWING AWAY MOUNT	\$162.00
☐	JFTLAT6X5SAL	LAT 6L X 5H W/ S/A MOUNT	\$334.00	☐	JFTLAT6X5SAR	LAT 6L X 5H W/ S/A MOUNT	\$334.00
	JFTLAT6X5	LATERAL 6L X 5H	\$172.00		JFTLAT6X5	LATERAL 6L X 5H	\$172.00
	JFTSAMNT	LATERAL SWING AWAY MOUNT	\$162.00		JFTSAMNT	LATERAL SWING AWAY MOUNT	\$162.00
☐	JFTLAT6X6SAL	LAT 6L X 6H W/ S/A MOUNT	\$334.00	☐	JFTLAT6X6SAR	LAT 6L X 6H W/ S/A MOUNT	\$334.00
	JFTLAT6X6	LATERAL 6L X 6H	\$172.00		JFTLAT6X6	LATERAL 6L X 6H	\$172.00
	JFTSAMNT	LATERAL SWING AWAY MOUNT	\$162.00		JFTSAMNT	LATERAL SWING AWAY MOUNT	\$162.00

LATERAL MEASURING GUIDE (Fixed and Swing Away)



			Back Width				
DIM	Description		10" Wide	12" Wide	14" Wide	16" Wide	18" Wide
A	Lateral Range (Min Spacing - Max Spacing)	Fixed Lateral	7.2" - 9"	9.2" - 11"	11.2" - 13"	13.2" - 15"	15.2" - 17"
		Fixed Lat w/ 1" Inset	5.2" - 7"	7.2" - 9"	9.2" - 11"	11.2" - 13"	13.2" - 15"
		Swing Away Lateral	5.2" - 11"	7.2" - 13"	9.2" - 15"	11.2" - 17"	13.2" - 19"
B	Lateral Height Range (Measured from bottom edge of shell to top of lateral)	10" Tall Back	7.2"				
		12" Tall Back	9.2"				
		14" Tall Back	11.2"				
		16" Tall Back		13.2"			
		18" Tall Back		15.2"			
		20" Tall Back				17.2"	
		22" Tall Back					19.2"

Spine Align Positioning Components

(JAY Fit Backs come standard with a 1 Small and 1 Large lumbar.)

SPINE-ALIGN MINI FIT KIT 1

Includes: 3 wedges, 2 lateral shims, and 2 posterior/lateral pelvic positioning pads.

☐	J3SPNK112	SPINE-ALIGN MINI FIT KIT 1 (12")	\$146.00
☐	J3SPNK114	SPINE-ALIGN MINI FIT KIT 1 (14")	\$146.00
☐	J3SPNK116	SPINE-ALIGN MINI FIT KIT 1 (16")	\$146.00
☐	J3SPNK118	SPINE-ALIGN MINI FIT KIT 1 (18")	\$146.00

WEDGES (Back Width) (W" x H" x D")

☐	J3SPNWE12	WEDGES (12") (3.5" x 6.5" x 1.6")	\$47.00
☐	J3SPNWE14	WEDGES (14") (3.5" x 6.5" x 1.6")	\$47.00
☐	J3SPNWE16	WEDGES (16") (5.5" x 6.5" x 1.6")	\$48.00
☐	J3SPNWE18	WEDGES (18") (7.5" x 6.5" x 1.6")	\$48.00

	Sm Lumbar	Lg Lumbar	Wedge	Curved Slim	Flat Shim	Lateral Shim	Sm Bean Bag	Lg Bean Bag	Post/Lat Pelvic
W	1.8" - 9.5"	3.5" - 9.5"	3.5" - 9.5"	10"	10"	2.4"	6"	10"	5.3"
H	3.2" - 4.3"	7"	6.5"	4"	4"	4.3"	4.5"	3.2"	3.6"
D	0.8" - 1"	1.6"	1.6"	0.5"	0.3"	0.5"	0.75"	0.75"	0.75"

SPINE-ALIGN MINI FIT KIT 2

Includes: 5 wedges, 2 lat shims, 2 post/lat pelvic, 1 lg lumbar, 1 flat shim, and 1 curved

☐	J3SPNK212	SPINE-ALIGN MINI FIT KIT 2 (12")	\$216.00
☐	J3SPNK214	SPINE-ALIGN MINI FIT KIT 2 (14")	\$216.00
☐	J3SPNK216	SPINE-ALIGN MINI FIT KIT 2 (16")	\$216.00
☐	J3SPNK218	SPINE-ALIGN MINI FIT KIT 2 (18")	\$216.00

OTHER INDIVIDUAL SPINE-ALIGN COMPONENTS

☐	J3SPNCR	Curved Shim	\$37.00
☐	J3SPNBS	Small Bean Bag	\$42.00
☐	J3SPNFL	Flat Shim	\$37.00
☐	J3SPNBL	Large Bean Bag	\$42.00
☐	J3SPNLS	Lateral Shim	\$67.00
☐	J3SPNLP	Posterior/Lateral Pelvic	\$67.00

Chest Straps, Harnesses, & Brackets

(Straps & Harnesses includes strap guides & mounting hdw.)

Classic Style

☐	JFT8107	Classic-Dynamic Pull Down (EI)	\$239.00
☐	JFT8100	Classic-Dynamic Pull Down (XS)	\$239.00
☐	JFT8101	Classic-Dynamic Pull Down (Sml)	\$239.00
☐	JFT8102	Classic-Dynamic Pull Down (Med)	\$239.00
☐	JFT8103	Classic-Dynamic Pull Down (Lrg)	\$239.00
☐	JFT8104	Classic-Dynamic Pull Down (X-Lrg)	\$239.00
☐	JFT8300	Classic-Structured Pull Down (XS)	\$239.00
☐	JFT8301	Classic-Structured Pull Down (Sml)	\$239.00
☐	JFT8302	Classic-Structured Pull Down (Med)	\$239.00
☐	JFT8303	Classic-Structured Pull Down (Lrg)	\$239.00
☐	JFT8304	Classic-Structured Pull Down (X-Lrg)	\$239.00

Center Opening w/ Buckle Style

☐	JFT9407	Center Opening w/ Buckle (EI)	\$239.00
☐	JFT9400	Center Opening w/ Buckle (XS)	\$239.00
☐	JFT9401	Center Opening w/ Buckle (Sml)	\$239.00
☐	JFT9402	Center Opening w/ Buckle (Med)	\$239.00
☐	JFT9403	Center Opening w/ Buckle (Lrg)	\$239.00
☐	JFT9404	Center Opening w/ Buckle (X-Lrg)	\$239.00

Contoured Style


☐	JFT9107	Contoured-Dynamic Pull Down (EI)	\$239.00
☐	JFT9100	Contoured-Dynamic Pull Down (XS)	\$239.00
☐	JFT9101	Contoured-Dynamic Pull Down (Sml)	\$239.00
☐	JFT9102	Contoured-Dynamic Pull Down (Med)	\$239.00
☐	JFT9103	Contoured-Dynamic Pull Down (Lrg)	\$239.00
☐	JFT9104	Contoured-Dynamic Pull Down (X-Lrg)	\$239.00
☐	JFT9300	Contoured-Structured Pull Down (XS)	\$239.00
☐	JFT9301	Contoured-Structured Pull Down (Sml)	\$239.00
☐	JFT9302	Contoured-Structured Pull Down (Med)	\$239.00
☐	JFT9303	Contoured-Structured Pull Down (Lrg)	\$239.00
☐	JFT9304	Contoured-Structured Pull Down (X-Lrg)	\$239.00

Padded Trunk Strap


☐	JFT4050	Padded Anterior Trunk Straps (XS)	\$239.00
☐	JFT4000	Padded Anterior Trunk Straps (Sml)	\$239.00
☐	JFT4100	Padded Anterior Trunk Straps (Med)	\$239.00
☐	JFT4200	Padded Anterior Trunk Straps (Lrg)	\$239.00
☐	JFT4300	Padded Anterior Trunk Straps (X-Lrg)	\$239.00

Back Pack Style


☐	JFT9701	Classic-Back Pack Style 12"	\$217.00
☐	JFT9703	Classic-Back Pack Style 18"	\$217.00
☐	JFT9601	Dynamic-Back Pack Style 12"	\$217.00
☐	JFT9603	Dynamic-Back Pack Style 18"	\$217.00



	A	B	C	D	E	F	G	H
EI	9.8"	3"	2"	5.3"	3"	3"	1"	6.5"
XS	12.5"	3.5"	2"	7.5"	3.5"	3"	2"	7.5"
S	15"	4.5"	2"	9"	4"	4"	2"	8.5"
M	19"	5.5"	2"	11"	5"	5"	3"	9.5"
L	22.5"	6.5"	2"	13"	6"	6"	3.5"	10.5"
XL	25.5"	6.5"	2.5"	14"	7"	7"	4.5"	11.5"



	A	B	C	D	E	F	G	H
EI	9.8"	3"	1.8"	5.3"	1"	3.5"	1.5"	6.5"
XS	12.5"	3.5"	2"	7.5"	1"	3.5"	1.5"	7.5"
S	15"	4.5"	2"	9"	1"	4"	2"	8.5"
M	18.5"	5.5"	2"	11"	1"	5"	2.5"	9.5"
L	22.5"	6.5"	2"	13"	1.5"	6"	3.5"	10.5"
XL	25.5"	6.5"	2.5"	14"	1.5"	7"	4.5"	11.5"



	A	B	C
EI	9.8"	7"	5"
XS	11"	8"	5.8"
S	12.3"	9"	6.5"
M	14.3"	11"	9.8"
L	15.5"	11.5"	9.8"
XL	17.5"	12.8"	10"


GUIDE BRACKETS & HARDWARE ONLY - Only used when mounting third party harness to the JAY Fit backrest.

☐	JGPSG	JAY Fit Strap Positioner Kit	\$72.00
--------------------------	-------	------------------------------	---------

JAY FIT - SEAT

☐	JFITSXXYY	JAY Fit Seat All Sizes	\$540.00
--------------------------	-----------	------------------------	----------

NOTE: Select appropriate assembly number below.
250 lbs Weight Capacity

**EACH SEAT SHIPS WITH THE FOLLOWING:**

Quick Mount Hardware
3DX with Black Stretch Cover and Inner Cover
Choice of Variety of Foam Options or Fluid bladders
Lateral Well Reducers (1 Left & 1 Right)
2 Anterior Well Reducers
Plus other positioning components are available.

SEAT WIDTH	10"	12"	14"	16"	18"
FITS CHAIR WIDTHS	10" - 11"	12" - 13"	14" - 15"	16" - 17"	18" - 19"
ACTUAL SHELL/FOAM WIDTH	9"	11"	13"	15"	17"

		Seat Width				
		10	12	14	16	18
Frame Depth	12" (8" - 12" Range)	☐ JFITS1012	☐ JFITS1212			
	14" (9" - 14" Range)	☐ JFITS1014	☐ JFITS1214	☐ JFITS1414		
	16" (11" - 16" Range)		☐ JFITS1216	☐ JFITS1416	☐ JFITS1616	☐ JFITS1816
	18" (13" - 18" Range)			☐ JFITS1418	☐ JFITS1618	☐ JFITS1818
	20" (15" - 20" Range)				☐ JFITS1620	☐ JFITS1820

The JAY ZIP Backrest has been dynamically tested and meets all criteria for use in motor vehicles according to Section 5.1 of ISO 16840-4 and Section 20, ANSI/RESNA WC/Volume 4.
Please visit www.sunparts.us for spare parts information.

JAY Fit Seat - Seat Mounting Hardware (Optional)

Note: Std Quick Mount Hardware is include with the seat but if Snaptite or Ht Adj Quick Mount Hardware is selected Std Quick Mnt hardware will automatically be omitted.

☐	Included in Seat	Standard Quick Mount	NC	☐	JYW-JFTSHDW-OPT02	Snaptite Hardware	\$77.00
☐	JYW-JFTSHDW-OPT01	Height Adjustable Quick Mount					

JAY Fit Seat - Seat Cover - Required

Oversized covers are required when Lateral Pelvic and/or Lateral Thigh Positioning Supports are selected below.

☐	JYW-JFTSCV-OPT01	3DX with Black Stretch	NC	☐	JYW-JFTSCV-OPT03	Oversized - 3DX with Black Stretch	NC
☐	JYW-JFTSCV-OPT02	Reverse Dartex	\$74.00	☐	JYW-JFTSCV-OPT04	Oversized - Reverse Dartex	\$74.00

JAY Fit Seat - Well Inserts - Required

Top Layer			Middle Layer		
☐	JYW-JFTSTF-OPT01	1/2" Soft Poly Foam	☐	JYW-JFTSMF-OPT01	1/2" Soft Poly
☐	JYW-JFTSTF-OPT02	1/2" Med/Soft Poly Foam	☐	JYW-JFTSMF-OPT02	1/2" Med/Soft Poly
☐	JYW-JFTSTF-OPT03	1/2" Medium Poly Foam	☐	JYW-JFTSMF-OPT03	1/2" Medium Poly
☐	JYW-JFTSTF-OPT04	1/2" Firm Poly Foam	☐	JYW-JFTSMF-OPT04	1/2" Firm Poly
☐	JYW-JFTSTF-OPT05	1/2" Soft Conform Foam	☐	JYW-JFTSMF-OPT05	1/2" Soft Conform
☐	JYW-JFTSTF-OPT06	1/2" Med/Soft Conform Foam	☐	JYW-JFTSMF-OPT06	1/2" Med/Soft Conform
☐	JYW-JFTSTF-OPT07	1/2" Medium Conform Foam	☐	JYW-JFTSMF-OPT07	1/2" Medium Conform
☐	JYW-JFTSTF-OPT08	1/2" Firm Conform Foam	☐	JYW-JFTSMF-OPT08	1/2" Firm Conform
☐	JYW-JFTSTF-OPT09	JAY Fluid Insert - Well Only	☐	JYW-JFTSMF-OPT09	Omit Middle Layer Foam
☐	JYW-JFTSTF-OPT10	JAY Fluid Insert - Entire Seat			NC

Omit Middle Layer Foam must and can only be selected when JAY Fluid insert is selected.

Bottom Layer		
☐	JYW-JFTSBF-OPT01	1/2" Soft Poly
☐	JYW-JFTSBF-OPT02	1/2" Med/Soft Poly
☐	JYW-JFTSBF-OPT03	1/2" Medium Poly
☐	JYW-JFTSBF-OPT04	1/2" Firm Poly
☐	JYW-JFTSBF-OPT05	1/2" Soft Conform
☐	JYW-JFTSBF-OPT06	1/2" Med/Soft Conform
☐	JYW-JFTSBF-OPT07	1/2" Medium Conform
☐	JYW-JFTSBF-OPT08	1/2" Firm Conform
☐	JYW-JFTSBF-OPT09	Omit Bottom Layer Foam

Omit Bottom Layer Foam must and can only be selected when JAY Fluid insert is selected.

JAY Fit Seat - Lateral Well Reducers (Optional)

Comes Standard with 1 Right and 1 Left Well Reducer

☐	JYW-JFTSWR-OPT01	1 Left Well Reducer	☐	JYW-JFTSWR-OPT07	1 Pair Well Reducer
☐	JYW-JFTSWR-OPT02	2 Left Well Reducer	☐	JYW-JFTSWR-OPT08	2 Pair Well Reducer
☐	JYW-JFTSWR-OPT03	3 Left Well Reducer	☐	JYW-JFTSWR-OPT09	3 Pair Well Reducer
☐	JYW-JFTSWR-OPT04	1 Right Well Reducer			
☐	JYW-JFTSWR-OPT05	2 Right Well Reducer			
☐	JYW-JFTSWR-OPT06	3 Right Well Reducer			

JAY Fit Seat - Anterior Well Reducers (Optional)

Comes Standard with 2 Well Reducers on 12" Seat depths, 3 on 14 & 16, and 4 on 18 and 20.

☐	JYW-JFTSAR-OPT01	1 Anterior Well Reducer
☐	JYW-JFTSAR-OPT02	2 Anterior Well Reducer
☐	JYW-JFTSAR-OPT03	3 Anterior Well Reducer

JAY Fit Seat - Adjustable Lateral Pelvic Supports (Optional)

Oversized covers are required when Lateral Pelvic Positioning Supports are selected.

☐	JYW-JFTSPS-OPT01	0"-1/2" inset Brkts w/2" Pad-Reg*(pr)	☐	JYW-JFTSPS-OPT05	0"-1/2" inset Brkts w/2" Pad-Long**(pr)
☐	JYW-JFTSPS-OPT02	0"-1/2" inset Brkts w/3" Pad-Reg*(pr)	☐	JYW-JFTSPS-OPT06	0"-1/2" inset Brkts w/3" Pad-Long**(pr)
☐	JYW-JFTSPS-OPT03	1/2"-1" inset Brkts w/2" Pad-Reg*(pr)	☐	JYW-JFTSPS-OPT07	1/2"-1" inset Brkts w/2" Pad-Long**(pr)
☐	JYW-JFTSPS-OPT04	1/2"-1" inset Brkts w/3" Pad-Reg*(pr)	☐	JYW-JFTSPS-OPT08	1/2"-1" inset Brkts w/3" Pad-Long**(pr)

*Note: Regular length pads measure 5-3/4" in Length

**Note: Long length pads measure 7-1/2" in Length

JAY Fit Seat - Adjustable Lateral Thigh Supports

Oversized covers are required when Lateral Thigh Positioning Supports are selected.

☐	JYW-JFTSLTS-OPT01	0"-1" inset Brkts w/2" Pad-Reg*(pr)	☐	JYW-JFTSLTS-OPT05	0"-1" inset Brkts w/2" Pad-Long**(pr)
☐	JYW-JFTSLTS-OPT02	0"-1" inset Brkts w/3" Pad-Reg*(pr)	☐	JYW-JFTSLTS-OPT06	0"-1" inset Brkts w/3" Pad-Long**(pr)
☐	JYW-JFTSLTS-OPT03	1/2" out-1/2" in Brkts w/2" Pad-Reg*(pr)	☐	JYW-JFTSLTS-OPT07	1/2" out-1/2" in Brkts w/2" Pad-Lg**(pr)
☐	JYW-JFTSLTS-OPT04	1/2" out-1/2" in Brkts w/3" Pad-Reg*(pr)	☐	JYW-JFTSLTS-OPT08	1/2" out-1/2" in Brkts w/3" Pad-Lg**(pr)

*Note: Regular length pads measure 9-1/4" in Length

**Note: Long length pads measure 11-1/4" in Length

JAY Fit Seat - Flip Down Medial Thigh Support

Medial Thigh Support includes flip down bracket, pads, and cover.

☐	JYW-JFTSMTS-OPT01	3" H x 4" L Pad Size - Pair	☐	JYW-JFTSMTS-OPT03	4" H x 6" L Pad Size - Pair
☐	JYW-JFTSMTS-OPT02	3" H x 6" L Pad Size - Pair	☐	JYW-JFTSMTS-OPT04	4" H x 8" L Pad Size - Pair

JAY FIT - ACCESSORIES

Anterior Pelvic Belts

Standard 2 Point Belt - Push Button 1"

☐	JFT5110	No Pad	\$64.00
☐	JFT5115	4.5" Pad	\$124.00
☐	JFT5116	5.5" Pad	\$124.00
☐	JFT5117	7.5" Pad	\$124.00

Standard 2 Point Belt - Push Button 1 1/2"

☐	JFT7210	No Pad	\$64.00
☐	JFT7212	7.5" Pad	\$119.00
☐	JFT7213	9.5" Pad	\$119.00

Medial 4 Point Belt - Push Button 1 1/2"

☐	JFT7242	4.5" Pad	\$148.00
☐	JFT7243	5.5" Pad	\$148.00

V-Contour Pelvic 4 Point Belt - Push Button Varies

☐	JFT3716	Small Pad - 1" Push Button	\$232.00
☐	JFT3727	Medium Pad - 1 1/2" Push Button	\$232.00
☐	JFT3738	Large Pad - 2" Push Button	\$232.00

Shoe Holders

Padded Shoe Holders (EA) - Must Select Left and Right

☐	TST6001	Padded Shoe Holder Small - Right	\$122.00
☐	TST6001	Padded Shoe Holder Small - Left	\$122.00
☐	TST6101	Padded Shoe Holder Medium - Right	\$122.00
☐	TST6101	Padded Shoe Holder Medium - Left	\$122.00
☐	TST6201	Padded Shoe Holder Large - Right	\$122.00
☐	TST6201	Padded Shoe Holder Large - Left	\$122.00
☐	TST6301	Padded Shoe Holder X-Large - Right	\$122.00
☐	TST6301	Padded Shoe Holder X-Large - Left	\$122.00

Tendon Relief Shoe Holder (EA) - Must Select Left and Right

☐	TST6061	Tendon Relief Shoe Holder Sml-Right	\$114.00
☐	TST6061	Tendon Relief Shoe Holder Sml-Left	\$114.00
☐	TST6161	Tendon Relief Shoe Holder Med-Right	\$114.00
☐	TST6161	Tendon Relief Shoe Holder Med-Left	\$114.00
☐	TST6261	Tendon Relief Shoe Holder Lrg-Right	\$114.00
☐	TST6261	Tendon Relief Shoe Holder Lrg-Left	\$114.00
☐	TST6361	Tendon Relief Shoe Holder XL-Right	\$114.00
☐	TST6361	Tendon Relief Shoe Holder XL-Left	\$114.00

Ankle Positioners

Ankle Positioners (Foot & Ankle) (EA) - Must Select Left and Right

☐	TST6510R	Small - Right	\$155.00
☐	TST6510L	Small - Left	\$155.00
☐	TST6520R	Medium - Right	\$155.00
☐	TST6520L	Medium - Left	\$155.00
☐	TST6540R	Large - Right	\$155.00
☐	TST6540L	Large - Left	\$155.00
☐	TST6560R	X-Large - Right	\$155.00
☐	TST6560L	X-Large - Left	\$155.00

Flexi Feet (PR)

☐	32600XS-RT	X-Small	\$257.00
☐	32600S-RT	Small	\$257.00
☐	32600M-RT	Medium	\$257.00
☐	32600L-RT	Large	\$257.00



Sunrise Medical Canada Inc.

1000 Creditstone Road Unit 2 · Concord, ON L4K 4P8