



ORDER FORM

# PAL Classroom Seat

# EILK03

The perfect balance of postural support and classroom practicality.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

## PAL Classroom Seat

Size 1 Weight Capacity: 66 lbs | Size 2 Weight Capacity: 77 lbs | Size 3 Weight Capacity: 88 lbs | Size 4 Weight Capacity: 110 lbs

✓ OPTION	DESCRIPTION	RETAIL	
<b>Seat</b>		<b>\$</b>	
Includes Lap Strap			
<input type="checkbox"/>	151ST92	Size 1 PAL Classroom Seat - Blue	\$1,075
<input type="checkbox"/>	151ST93	Size 1 PAL Classroom Seat - Green	\$1,075
<input type="checkbox"/>	151ST94	Size 2 PAL Classroom Seat - Blue	\$1,211
<input type="checkbox"/>	151ST95	Size 2 PAL Classroom Seat - Green	\$1,211
<input type="checkbox"/>	151ST96	Size 3 PAL Classroom Seat - Blue	\$1,348
<input type="checkbox"/>	151ST97	Size 3 PAL Classroom Seat - Green	\$1,348
<input type="checkbox"/>	151ST98	Size 4 PAL Classroom Seat - Blue	\$1,614
<input type="checkbox"/>	151ST99	Size 4 PAL Classroom Seat - Green	\$1,614
<b>Side Pad</b>		<b>\$</b>	
Must match seat size.			
<input type="checkbox"/>	151ZV34	Size 1 PAL Side Pads 1" Thick - Blue (Pair)	\$139
<input type="checkbox"/>	151ZV35	Size 1 PAL Side Pads 1" Thick - Green (Pair)	\$139
<input type="checkbox"/>	151ZV36	Size 2 PAL Side Pads 1" Thick - Blue (Pair)	\$139
<input type="checkbox"/>	151ZV37	Size 2 PAL Side Pads 1" Thick - Green (Pair)	\$139
<input type="checkbox"/>	151ZV38	Size 3 and 4 PAL Side Pads 1" Thick - Blue (Pair)	\$139
<input type="checkbox"/>	151ZV39	Size 3 and 4 PAL Side Pads 1" Thick - Green (Pair)	\$139
<input type="checkbox"/>	151ZV40	Size 1 PAL Side Pads 2" Thick - Blue (Pair)	\$139
<input type="checkbox"/>	151ZV41	Size 1 PAL Side Pads 2" Thick - Green (Pair)	\$139
<input type="checkbox"/>	151ZV42	Size 2 PAL Side Pads 2" Thick - Blue (Pair)	\$139
<input type="checkbox"/>	151ZV43	Size 2 PAL Side Pads 2" Thick - Green (Pair)	\$139
<input type="checkbox"/>	151ZV44	Size 3 and 4 PAL Side Pads 2" Thick - Blue (Pair)	\$139
<input type="checkbox"/>	151ZV45	Size 3 and 4 PAL Side Pads 2" Thick - Green (Pair)	\$139
<b>Ramped Seat Cushion</b>		<b>\$</b>	
<input type="checkbox"/>	151CU115	Size 1 PAL Ramped Seat Cushion - Blue	\$93
<input type="checkbox"/>	151CU116	Size 1 PAL Ramped Seat Cushion - Green	\$93
<input type="checkbox"/>	151CU117	Size 2 PAL Ramped Seat Cushion - Blue	\$93
<input type="checkbox"/>	151CU118	Size 2 PAL Ramped Seat Cushion - Green	\$93
<input type="checkbox"/>	151CU119	Size 3 PAL Ramped Seat Cushion - Blue	\$93
<input type="checkbox"/>	151CU120	Size 3 PAL Ramped Seat Cushion - Green	\$93
<input type="checkbox"/>	151CU121	Size 4 PAL Ramped Seat Cushion - Blue	\$93
<input type="checkbox"/>	151CU122	Size 4 PAL Ramped Seat Cushion - Green	\$93
<b>Tray and Tray Pad</b>		<b>\$</b>	
Must match seat size.			
<input type="checkbox"/>	151A222	Size 1 PAL Tray	\$253
<input type="checkbox"/>	151A223	Size 2 PAL Tray	\$253
<input type="checkbox"/>	151A224	Size 3 PAL Tray	\$253
<input type="checkbox"/>	151A225	Size 4 PAL Tray	\$253
<input type="checkbox"/>	151A226	Size 1 PAL Tray Pad	\$116
<input type="checkbox"/>	151A227	Size 2 PAL Tray Pad	\$116
<input type="checkbox"/>	151A228	Size 3 PAL Tray Pad	\$116
<input type="checkbox"/>	151A229	Size 4 PAL Tray Pad	\$116

## PAL Classroom Seat

✓ OPTION	DESCRIPTION	RETAIL
<b>Pommel</b>		\$
Size 1 Pommel suitable for Size 1 seat, Size 2 Pommel for Size 2, 3 and 4 seats.		
<input type="checkbox"/> 1518A114	Size 1 PAL Pommel - Blue	\$185
<input type="checkbox"/> 1518A115	Size 1 PAL Pommel - Green	\$185
<input type="checkbox"/> 1518A116	Size 2 PAL Pommel - Blue	\$185
<input type="checkbox"/> 1518A117	Size 2 PAL Pommel - Green	\$185
<b>Footrest</b>		\$
Must match seat size		
<input type="checkbox"/> 151F378	Size 1 PAL Footrest	\$141
<input type="checkbox"/> 151F379	Size 2 PAL Footrest	\$141
<input type="checkbox"/> 151F380	Size 3 PAL Footrest	\$141
<input type="checkbox"/> 151F381	Size 4 PAL Footrest	\$141
<b>Rocker / Stabilizer</b>		\$
Not available with PAL Ski with Casters		
<input type="checkbox"/> 151RX19	Size 1 PAL Rocker/Stabilizer (Pair)	\$159
<input type="checkbox"/> 151RX20	Size 2 PAL Rocker/Stabilizer (Pair)	\$159
<input type="checkbox"/> 151RX21	Size 3 PAL Rocker/Stabilizer (Pair)	\$159
<input type="checkbox"/> 151RX22	Size 4 PAL Rocker/Stabilizer (Pair)	\$159
<b>Caster</b>		\$
Not available with Rocker/Stabilizer		
<input type="checkbox"/> 151SY5	Size 1 PAL Ski with Casters (Pair)	\$277
<input type="checkbox"/> 151SY6	Size 2 PAL Ski with Casters (Pair)	\$277
<input type="checkbox"/> 151SY7	Size 3 PAL Ski with Casters (Pair)	\$277
<input type="checkbox"/> 151SY8	Size 4 PAL Ski with Casters (Pair)	\$277



### Technical Specifications - Imperial

Frame Sizing Specification	Size 1	Size 2	Size 3	Size 4
Age	1 - 3 years	2 - 5 years	4 - 8 years	7 - 12 years
Child Height	27.6" - 37.4"	35.4" - 41.3"	39.4" - 51.2"	47" - 63"
Max User Weight	66 lbs	77 lbs	88 lbs	110 lbs
Seat to Footrest	5.5" - 9.4"	6.5" - 11.2"	9.1" - 13.8"	11.2" - 16"
Seat Depth	7.1" - 9.4"	8.7" - 11.4"	9.8" - 13"	12" - 16"
Backrest Height	10.4"	12.4"	15"	17.7"
Backrest Angle	95° - 105°	95° - 105°	95° - 105°	95° - 105°
Seat Width	5.9" - 9.8"	6.7" - 10.6"	9.1" - 13"	11.1" - 15.2"
Seat Weight	9.6 lbs	13.2 lbs	17.9 lbs	22 lbs
Tray Size	20"W x 18"D	21.3"W x 19.7"D	23.6"W x 21.7"D	28.7"W x 23.6"D

### Technical Specifications - Metric

Frame Sizing Specification	Size 1	Size 2	Size 3	Size 4
Age	1 - 3 years	2 - 5 years	4 - 8 years	7 - 12 years
Child Height	700 - 950 mm	900 - 1050 mm	1000 - 1300 mm	1200 - 1600 mm
Max User Weight	30 kg	35 kgs	40 kgs	50 kgs
Seat to Footrest	40 - 240 mm	165 - 285 mm	230 - 350 mm	280 - 400 mm
Seat Depth	80 - 240 mm	220 - 290 mm	250 - 330 mm	300 - 400 mm
Backrest Height	265 mm	315mm	380 mm	450 mm
Backrest Angle	95° - 105°	95° - 105°	95° - 105°	95° - 105°
Seat Width	50 - 250 mm	170 - 270 mm	230 - 330 mm	280 - 380 mm
Seat Weight	4.35kg	6 kgs	8.15 kgs	10 kgs
Tray Size	508mm W x 457mm D	541mm W x 500mm D	599mm W x 551mm D	729mmW x 599mmD

Sunrise Medical Canada Inc.  
1000 Creditstone Road, Unit 2  
Concord, ON L4K 4P8



MK-129852-CAN Rev. 9 © 04/2025

