



Zippie X'Cape
May 2025
EIFZ1

| | | | |
|--|--|--|--|
| Date: _____ | | SUBMITTING FOR: <input type="checkbox"/> Quote <input type="checkbox"/> Order | |
| Dealer Acct #: _____ | | PO#: _____ | |
| Dealer: _____ | | Ship To: _____ | |
| Dealer Contact: _____ | | Attention: _____ | |
| Dealer Address: _____ | | Address: _____ | |
| Dealer City: _____ Prov: _____ PC: _____ | | Address: _____ | |
| Dealer Phone: () _____ Fax: () _____ | | Ship To City: _____ Prov: _____ PC: _____ | |
| Confirm Via: <input type="checkbox"/> Fax <input type="checkbox"/> Email | | Ship To Phone: () _____ Fax: () _____ | |
| Mark For: _____ | | | |
| | | | |

| ✓ | Option # | Definition | Price | ✓ | Option # | Definition | Price |
|---|----------|--|---------|---|----------|---|---------|
| 1. FRAME | | | | | | | |
| A. Zippie X'CAPE - TRANSIT INCLUDED 165 lbs. Weight Capacity | | | | B. Zippie X'CAPE - TRANSIT OPTIONAL 165 lbs. Weight Capacity | | | |
| A. FRAME STYLE (Includes True Fit) (N/A w/ Xtender) | | | | B. FRAME STYLE (Includes True Fit) | | | |
| <input type="checkbox"/> | 121WM87 | X'CAPE - Transit Incl. without seating (JAY seating opts available) | \$3,494 | <input type="checkbox"/> | 121WM83 | X'CAPE - Transit Opt. without seating (JAY seating opts available) | \$3,387 |
| <input type="checkbox"/> | 121WM84 | X'CAPE - Transit Incl. w/ Seating (Sling and Back Uph.) | \$3,494 | <input type="checkbox"/> | 121WM85 | X'CAPE - Transit Opt. w/ Seating (Sling and Back Uph.) | \$3,387 |

| TRANSIT SYSTEM WC-19 | | | | | | | |
|--------------------------|---------|-------------------------------|----|--------------------------|---------|--------------|----|
| <input type="checkbox"/> | 121TO48 | Transit Opt. (N/A w/ Xtender) | NC | <input type="checkbox"/> | 121TO40 | Omit Transit | NC |

| FRAME TYPE | | | | | | | |
|--------------------------|---------|-------------------------|----|--------------------------|---------|------------------------------|----|
| <input type="checkbox"/> | 121FT64 | Zippie-Swing Away Frame | NC | <input type="checkbox"/> | 121FT65 | Zippie-80° Fixed Front Frame | NC |

| 2. SIZE/SEATING | | | | | | | |
|--|--------|----------------------------------|----|--------------------------|--------|----------------------------------|----|
| 2A. FRAME WIDTH | | | | | | | |
| 8", 9" and 10" widths N/A with Xtender | | | | | | | |
| <input type="checkbox"/> | 121W48 | 8" (Frame Width Growth 9"-10") | NC | <input type="checkbox"/> | 121W53 | 13" (Frame Width Growth 14"-15") | NC |
| <input type="checkbox"/> | 121W49 | 9" (Frame Width Growth 10"-11") | NC | <input type="checkbox"/> | 121W54 | 14" (Frame Width Growth 15"-16") | NC |
| <input type="checkbox"/> | 121W50 | 10" (Frame Width Growth 11"-12") | NC | <input type="checkbox"/> | 121W55 | 15" (Frame Width Growth 16"-17") | NC |
| <input type="checkbox"/> | 121W51 | 11" (Frame Width Growth 12"-13") | NC | <input type="checkbox"/> | 121W56 | 16" (Frame Width Growth 17"-18") | NC |
| <input type="checkbox"/> | 121W52 | 12" (Frame Width Growth 13"-14") | NC | | | | |

| 2B/C. SLING DEPTH, FRAME DEPTH AND BACK POST POSITION. See matrix on page 2. | | | | | | | |
|---|--|--|--|--|--|--|--|
| 8"-11" Frame Depth with 20", 22" or 24" wheels and Personalized Seating must have hublock or Fixed Front frame. When selecting Frame Depth with JAY seating, add 1" to accommodate seating. Example: need 12" depth, order 13" depth. | | | | | | | |

| 2B. Sling Depth when choosing with Seating Sling must be up to 3" less than or equal to the Frame Depth (C.) | | | | | | | |
|--|--------|-----|----|--------------------------|---------|-----|----|
| <input type="checkbox"/> | 121GJ1 | 8" | NC | <input type="checkbox"/> | 121GJ7 | 14" | NC |
| <input type="checkbox"/> | 121GJ2 | 9" | NC | <input type="checkbox"/> | 121GJ8 | 15" | NC |
| <input type="checkbox"/> | 121GJ3 | 10" | NC | <input type="checkbox"/> | 121GJ9 | 16" | NC |
| <input type="checkbox"/> | 121GJ4 | 11" | NC | <input type="checkbox"/> | 121GJ10 | 17" | NC |
| <input type="checkbox"/> | 121GJ5 | 12" | NC | <input type="checkbox"/> | 121GJ11 | 18" | NC |
| <input type="checkbox"/> | 121GJ6 | 13" | NC | | | | |

| 2C. Frame Depth and Back Post Position with Growth Note: Must be equal or greater than Sling Depth if chosen with seat Sling | | | | | | | |
|--|---------|---|----|--------------------------|---------|----------------------------------|----|
| <input type="checkbox"/> | 121FD40 | 8" (Frame depth growth 9"-11", N/A w/ 24" wheels) | NC | <input type="checkbox"/> | 121FD46 | 14" (Frame depth growth 15"-17") | NC |
| <input type="checkbox"/> | 121FD41 | 9" (Frame depth growth 10"-12", N/A w/ 24" wheels) | NC | <input type="checkbox"/> | 121FD47 | 15" (Frame depth growth 16"-18") | NC |
| <input type="checkbox"/> | 121FD42 | 10" (Frame depth growth 11"-13") | NC | <input type="checkbox"/> | 121FD48 | 16" (Frame depth growth 17"-19") | NC |
| <input type="checkbox"/> | 121FD43 | 11" (Frame depth growth 12"-14") | NC | <input type="checkbox"/> | 121FD49 | 17" (Frame depth growth 18"-20") | NC |
| <input type="checkbox"/> | 121FD44 | 12" (Frame depth growth 13"-15") | NC | <input type="checkbox"/> | 121FD50 | 18" (Frame depth growth 19"-20") | NC |
| <input type="checkbox"/> | 121FD45 | 13" (Frame depth growth 14"-16") | NC | | | | |

SEATING

Seat Sling

| | | | |
|--------------------------|---------|---|-----|
| <input type="checkbox"/> | 121OP01 | Standard Seat Sling | STD |
| <input type="checkbox"/> | 121OP04 | Omit Seat Sling (available when base is ordered with seating) | NC |

JAY Cushions (Optional)

To order any other JAY cushions on-chair, please see the JAY price list.
Cushion name or part number must be filled out w/ desired width and depth

| | | | |
|--------------------------|---------|-----------------------|--------------|
| <input type="checkbox"/> | 121JC20 | JAY Bundled Cushion | Price List |
| | | Cushion Name or Part# | |
| | | Width (in.): | Depth (in.): |

Seat Sling to Frame Depth Matrix

| | | Sling Length | | | | | | | | | | |
|--------------|-----|--|----|----|----|----|----|----|----|----|----|----|
| | | Difference from front of frame to front of sling | | | | | | | | | | |
| Frame Length | 8" | 0" | | | | | | | | | | |
| | 9" | 0" | 0" | | | | | | | | | |
| | 10" | 0" | 0" | 0" | | | | | | | | |
| | 11" | 1" | 1" | 1" | 0" | | | | | | | |
| | 12" | | 2" | 2" | 1" | 0" | | | | | | |
| | 13" | | | 3" | 2" | 1" | 0" | | | | | |
| | 14" | | | | 3" | 2" | 1" | 0" | | | | |
| | 15" | | | | | 3" | 2" | 1" | 0" | | | |
| | 16" | | | | | | 2" | 1" | 0" | 0" | | |
| | 17" | | | | | | | 1" | 0" | 0" | 0" | |
| | 18" | | | | | | | | 0" | 0" | 0" | 0" |

3. CASTERS/SEAT HEIGHT

Caster Wheel

| | | | |
|--------------------------|---------|--------------------------------|-------|
| <input type="checkbox"/> | 121CW30 | 3" Micro Light Up | NC |
| <input type="checkbox"/> | 121CW3 | 4" Low-Profile Polyurethane | NC |
| <input type="checkbox"/> | 121CW31 | 4" Micro Light Up | \$110 |
| <input type="checkbox"/> | 121CW64 | 4" x 1-1/2" Aluminum Soft Roll | \$110 |
| <input type="checkbox"/> | 121CW5 | 5" Low-Profile Polyurethane | NC |
| <input type="checkbox"/> | 121CW51 | 5" Micro Light Up | \$110 |
| <input type="checkbox"/> | 121CW6 | 5" Soft Roll (Grey) | \$126 |
| <input type="checkbox"/> | 121CW7 | 5"x1-1/2" Semi-Pneumatic | \$110 |
| <input type="checkbox"/> | 121CW66 | 5"x1-1/2" Aluminum Soft Roll | \$110 |
| <input type="checkbox"/> | 121CW34 | 6" Poly | NC |
| <input type="checkbox"/> | 121CW35 | 6"x1-1/2" Semi-Pneumatic | \$110 |

Caster Forks

| | | | |
|--------------------------|---------|-------------------|----|
| <input type="checkbox"/> | 121CF65 | Multiple Position | NC |
|--------------------------|---------|-------------------|----|

Front Seat Heights Matrix

| Rear Wheel Size | Front Caster | | | |
|--|--------------------------------|--------------------------------|--------------------------------|--------------------------------|
| | <input type="checkbox"/> 3" | <input type="checkbox"/> 4" | <input type="checkbox"/> 5" | <input type="checkbox"/> 6" |
| <input type="checkbox"/> 12" <input type="checkbox"/> 16" <input type="checkbox"/> 20" <input type="checkbox"/> 22" | <input type="checkbox"/> 13.0" | <input type="checkbox"/> 13.0" | <input type="checkbox"/> 14.0" | <input type="checkbox"/> 15.0" |
| | <input type="checkbox"/> 13.5" | <input type="checkbox"/> 13.5" | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15.5" |
| | <input type="checkbox"/> 14.0" | <input type="checkbox"/> 14.0" | <input type="checkbox"/> 15.0" | <input type="checkbox"/> 16.0" |
| | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16.5" |
| | <input type="checkbox"/> 15.0" | <input type="checkbox"/> 15.0" | <input type="checkbox"/> 16.0" | <input type="checkbox"/> 17.0" |
| | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17.5" |
| | <input type="checkbox"/> 16.0" | <input type="checkbox"/> 16.0" | <input type="checkbox"/> 17.0" | <input type="checkbox"/> 18.0" |
| | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17.5" | |
| | <input type="checkbox"/> 17.0" | <input type="checkbox"/> 17.0" | | |
| | N/A | <input type="checkbox"/> 17.5" | | |
| <input type="checkbox"/> 24" | N/A | N/A | <input type="checkbox"/> 14.0" | <input type="checkbox"/> 15.0" |
| | | | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15.5" |
| | <input type="checkbox"/> 14.0" | <input type="checkbox"/> 14.0" | <input type="checkbox"/> 15.0" | <input type="checkbox"/> 16.0" |
| | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16.5" |
| | <input type="checkbox"/> 15.0" | <input type="checkbox"/> 15.0" | <input type="checkbox"/> 16.0" | <input type="checkbox"/> 17.0" |
| | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17.5" |
| | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17.0" | <input type="checkbox"/> 18.0" |
| | <input type="checkbox"/> 17.0" | <input type="checkbox"/> 17.0" | <input type="checkbox"/> 17.5" | |
| | N/A | <input type="checkbox"/> 17.5" | N/A | |
| | | | | N/A |

Standard Frame

Hemi Fame (N/A w/ 90° hanger, 80° hanger footplate in lowest position may cause interference w/caster)

Depends on Caster setting

FRONT SEAT HEIGHTS ARE SHOWN IN THE BOXES IN THE FRONT SEAT HEIGHTS MATRIX. DEPENDING ON CG SETTING, SEAT DEPTH, AND FRAME SIZE SOME OF THE SIZES MAY NOT BE AVAILABLE

4.REAR WHEELS/SEAT HEIGHT

XTENDER POWER ASSIST OPTION

22" & 24" wheels only. RSH 14" to 16" for 22", RSH 15.0" to 17" for 24". 0° camber, pneumatic tire. Anti-Tips no charge. Dual drive mode available through B4Me.

Xtender Power Assist

| | | | |
|--------------------------|--------|--|----------|
| <input type="checkbox"/> | 121XT6 | 22" Xtender Power Assist JWX-2 NICKEL | \$9,598 |
| <input type="checkbox"/> | 121XT7 | 22" Xtender Power Assist JWX-2 LITHIUM | \$11,143 |
| <input type="checkbox"/> | 121XT4 | 24" Xtender Power Assist JWX-2 NICKEL | \$9,598 |
| <input type="checkbox"/> | 121XT5 | 24" Xtender Power Assist JWX-2 LITHIUM | \$11,143 |

Please select one of the following handrim options

| | | | |
|--------------------------|--------|-----------------|----|
| <input type="checkbox"/> | 121XH5 | Plastic Coated | NC |
| <input type="checkbox"/> | 121XH6 | Stainless Steel | NC |

Xtender - Other Options

| | | | |
|--------------------------|---------|-----------------------------------|---------|
| <input type="checkbox"/> | 121XA10 | Spare NiMH Power Pack - NICKEL | \$1,848 |
| <input type="checkbox"/> | 121XA11 | Spare Li-ion Power Pack - LITHIUM | \$2,379 |

Xtender Wheel Options

| | | | |
|--------------------------|--------|-----------------------------------|-------|
| <input type="checkbox"/> | 121XA8 | SPARE 22" Lite Spoke Pneu Alum | NC |
| <input type="checkbox"/> | 121XA9 | SPARE 22" Lite Spoke Pneu Plastic | \$141 |
| <input type="checkbox"/> | 121XA2 | SPARE 24" Lite Spoke Pneu Alum | NC |
| <input type="checkbox"/> | 121XA4 | SPARE 24" Lite Spoke Pneu Plastic | NC |

Rear Wheels

N/A w/ Xtender & 5 Spoke Mag Wheels. N/A w/ Hub Lock. Std Mag Wheels only available with One Arm Drive. 24" wheels. N/A with 8" and 9" frame depths. 8"-12" Frame Depth with 20-24" Rear Wheels and ELR's is ONLY available with Hub Lock Wheels.

| | | | |
|--------------------------|----------|--|---------|
| <input type="checkbox"/> | 121RW48 | 12" Mag (Not available with handrims or one arm drive) | \$47 |
| <input type="checkbox"/> | 121RW14 | 16" Mag (Not available with handrims or one arm drive) | NC |
| <input type="checkbox"/> | 121RW3 | 20" Mag Wheel | NC |
| <input type="checkbox"/> | 121RW8 | 22" Mag Wheel | NC |
| <input type="checkbox"/> | 121RW1 | 24" Mag Wheel | NC |
| <input type="checkbox"/> | 121RW103 | 20" Mag Wheel 5 Spoke | NC |
| <input type="checkbox"/> | 121RW104 | 22" Mag Wheel 5 Spoke | NC |
| <input type="checkbox"/> | 121RW102 | 24" Mag Wheel 5 Spoke | NC |
| <input type="checkbox"/> | 121RW74 | 20" Lite Spoke Wheel | NC |
| <input type="checkbox"/> | 121RW75 | 22" Lite Spoke Wheel | NC |
| <input type="checkbox"/> | 121RW76 | 24" Lite Spoke Wheel | NC |
| <input type="checkbox"/> | 121RW06 | 22" Ultra light wt. Spoke Wheel 2.0 | \$357 |
| <input type="checkbox"/> | 121RW10 | 24" Ultra light wt. Spoke Wheel 2.0 | \$357 |
| <input type="checkbox"/> | 121RW28 | 22" Spinergy SPOX Black Spokes | \$864 |
| <input type="checkbox"/> | 121RW95 | 24" Spinergy SPOX Black Spokes | \$864 |
| <input type="checkbox"/> | 121RW32 | 22" Spinergy LX Black Spokes | \$1,111 |
| <input type="checkbox"/> | 121RW98 | 24" Spinergy LX Black Spokes | \$1,111 |

Rear Wheel Tire

| | | | |
|--------------------------|---------|----------------------------------|-------|
| <input type="checkbox"/> | 121RT3 | Full Poly Tire | NC |
| <input type="checkbox"/> | 121RT6 | High Pressure Clincher | \$105 |
| <input type="checkbox"/> | 121RT27 | Knobby V-TRAK by Kenda | \$121 |
| <input type="checkbox"/> | 121RT38 | SHOX Performance Poly | \$146 |
| <input type="checkbox"/> | 121RT1 | Pneumatic | NC |
| <input type="checkbox"/> | 121RT2 | Pneumatic Airless Insert | \$247 |
| <input type="checkbox"/> | 121RT21 | Schwalbe Marathon Plus Evolution | \$157 |
| <input type="checkbox"/> | 121RT28 | Schwalbe Speed | \$162 |
| <input type="checkbox"/> | 121RT4 | Low Poly | NC |

| Tire | Wheel | 5 Spoke Mag | Std. Mag | Lite Spoke | Ultra Light-weight | Spinergy |
|---------------------------|-------|-------------|---------------------|-------------|--------------------|----------|
| Full Poly Tire | | 20",22",24" | 20",22",24" | 20",22",24" | | |
| High Pressure Clincher | | | | 22",24" | 24" | 22",24" |
| SHOX Performance Poly | | 20",22",24" | | 20",22",24" | | 22",24" |
| Knobby V-TRAK | | | | 24" | | 24" |
| Pneumatic | | 20",22",24" | 12",16",20",22",24" | 20",22",24" | | 22",24" |
| Pneumatic Airless Insert | | 20",22",24" | 12",16",20",22",24" | 20",22",24" | | 22",24" |
| Schwalbe Marathon Plus | | | | 24" | 22",24" | 22",24" |
| Schwalbe Speed Run | | | | 24" | 24" | 24" |
| Low Poly | | | 12" | 24" | | |
| 22" = 489 mm 24" = 540 mm | | | | | | |

Handrim

N/A w/ Xtender. Natural Fit N/A w/ 20" and 22" 5 Spoke Mag Wheel

| | | | |
|--------------------------|---------|--|-------|
| <input type="checkbox"/> | 121HR1 | Aluminum Anodized | NC |
| <input type="checkbox"/> | 121HR4 | Plastic Coated (N/A on Ultra and 22" Spinergy LX,SPOX) | \$146 |
| <input type="checkbox"/> | 121HR22 | Natural Fit -Standard Grip (N/A on Ultra Lite or 20" Lite Spoke) | \$374 |
| <input type="checkbox"/> | 121HR23 | Natural Fit - Super Grip (N/A on Ultra Lite or 20" Lite Spoke) | \$434 |
| <input type="checkbox"/> | 121HR30 | Natural Fit LT - Standard Grip (N/A on Ultra Lite or 20" Lite Spoke) | \$521 |
| <input type="checkbox"/> | 121HR31 | Natural Fit LT - Super Grip (N/A on Ultra Lite or 20" Lite Spoke) | \$592 |
| <input type="checkbox"/> | 121HR24 | Spinergy Flex Rim (Only on 24" Spinergy LX) | \$701 |
| <input type="checkbox"/> | 121HR7 | Omit Handrim (Option may not be available on all rear wheels) | NC |

Axle Plate

| | | | |
|--------------------------|---------|----------|----|
| <input type="checkbox"/> | 121AL33 | Standard | NC |
|--------------------------|---------|----------|----|

4. REAR WHEELS/SEAT HEIGHT Cont.

One Arm Drive

NA w/ Xtender. N/A w/ 12"-16" wheels. Only avail. On 12"-16" widths and 20", 22" and 24" mag wheels. N/A with Natural Fit handrims. N/A with SPHA arms

| | | | |
|--------------------------|---------|--|---------|
| <input type="checkbox"/> | 121AL16 | One Arm Drive | \$1,382 |
| <input type="checkbox"/> | 121OA2 | Standard Handrim | NC |
| <input type="checkbox"/> | 121HM1 | <input type="checkbox"/> Right Hand Mount | NC |
| <input type="checkbox"/> | 121HM2 | <input type="checkbox"/> Left Hand Mount | NC |
| <input type="checkbox"/> | 121OA1 | Small Handrim (N/A w/ 20" rear wheels) | NC |
| <input type="checkbox"/> | 121HM1 | <input type="checkbox"/> Right Hand Mount | NC |
| <input type="checkbox"/> | 121HM2 | <input type="checkbox"/> Left Hand Mount | NC |

Axle

| | | | |
|--------------------------|---------|---|-------|
| <input type="checkbox"/> | 121AX1 | Quick Release (N/A w/ One Arm Drive) | NC |
| <input type="checkbox"/> | 121AX11 | Quad Quick Release (Only Available with 12" and 16" Mag Wheels) | \$182 |

Center of Gravity Preset (dependent on wheels, casters & depth, From standard frame depth +2" forward and - depth growth up to -3")

| | | | |
|--------------------------|---------|--------------------------------|----|
| <input type="checkbox"/> | 121CC28 | -0.5" | NC |
| <input type="checkbox"/> | 121CC26 | 0 (Flush w/ backrest) | NC |
| <input type="checkbox"/> | 121CC16 | 0.5" (Only avail 9"-18" deep) | NC |
| <input type="checkbox"/> | 121CC20 | 1" (Only avail 10"-18" deep) | NC |
| <input type="checkbox"/> | 121CC15 | 1.5" (Only avail 10"-18" deep) | NC |
| <input type="checkbox"/> | 121CC17 | 2" (Only avail 16"-18" deep) | NC |
| <input type="checkbox"/> | 121CC18 | 2.5" (Only avail 16"-18" deep) | NC |

Camber

| | | | |
|--------------------------|--------|--------------------------------------|----|
| <input type="checkbox"/> | 121CI1 | 0° | NC |
| <input type="checkbox"/> | 121CI2 | 3° (N/A w/ Single Post Ht. Adj. Arm) | NC |

Rear Seat Height (see Allowable Squeeze chart)

| | | | |
|--------------------------|---------|-------|----|
| <input type="checkbox"/> | 121RH32 | 11" | NC |
| <input type="checkbox"/> | 121RH33 | 11.5" | NC |
| <input type="checkbox"/> | 121RH9 | 12" | NC |
| <input type="checkbox"/> | 121RH10 | 12.5" | NC |
| <input type="checkbox"/> | 121RH11 | 13" | NC |
| <input type="checkbox"/> | 121RH12 | 13.5" | NC |
| <input type="checkbox"/> | 121RH13 | 14" | NC |
| <input type="checkbox"/> | 121RH14 | 14.5" | NC |
| <input type="checkbox"/> | 121RH15 | 15" | NC |
| <input type="checkbox"/> | 121RH16 | 15.5" | NC |
| <input type="checkbox"/> | 121RH17 | 16" | NC |
| <input type="checkbox"/> | 121RH18 | 16.5" | NC |
| <input type="checkbox"/> | 121RH19 | 17" | NC |
| <input type="checkbox"/> | 121RH26 | 17.5" | NC |
| <input type="checkbox"/> | 121RH20 | 18" | NC |

| Rear Wheel | 11" | 11.5" | 12" | 12.5" | 13" | 13.5" | 14" | 14.5" | 15" | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" |
|--------------|-----|-------|-----|-------|-----|-------|-----|-------|-----|-------|-----|-------|-----|-------|-----|
| 12" low poly | | | | | | | | | | | | | | | |
| 12" pneu | | | | | | | | | | | | | | | |
| 16" | | | | | | | | | | | | | | | |
| 20" | | | | | | | | | | | | | | | |
| 22" | | | | | | | | | | | | | | | |
| 24" | | | | | | | | | | | | | | | |
| 24" pneu | | | | | | | | | | | | | | | |

When tires are not listed default to the wheel size. All Heights are subject to vary slightly

Allowable Squeeze (may vary +/- 1/2")

| Frame Depth | Max Squeeze |
|-------------|-------------|
| 8"-9" | 1.0" |
| 10"-12" | 1.5" |
| 13"-15" | 2.0" |
| 16"-18" | 2.5" |

5. WHEEL LOCK

Hub Lock

Only available with 12" or 16" mag wheels or Lite Spoke wheels.
Not available with One Arm Drive.

| | | | |
|--------------------------|---------|-----------------------------|-------|
| <input type="checkbox"/> | 121LK01 | 12" Hub Lock | \$162 |
| <input type="checkbox"/> | 121LK02 | 16" Hub Lock | \$162 |
| <input type="checkbox"/> | 121LK04 | 20"-24" Hub Lock Lite Spoke | \$162 |

Hub Lock Mount Position

Foot lock may limit folding width.

| | | | |
|--------------------------|---------|--------------------------------------|----|
| <input type="checkbox"/> | 121HH01 | Foot lock Mount | NC |
| <input type="checkbox"/> | 121HH02 | Seat Mount | NC |
| <input type="checkbox"/> | 121HH03 | Back Mount (back ht. may be limited) | NC |

Wheel Lock

High Mount will be mounted in low position on 12 and 16 rear wheels.
Some may not be available depending on wheel size and position.

| | | | |
|--------------------------|---------|-----------------|-------|
| <input type="checkbox"/> | 121WL1 | High Mount Push | NC |
| <input type="checkbox"/> | 121WL2 | High Mount Pull | NC |
| <input type="checkbox"/> | 121WL30 | Compact | \$162 |
| <input type="checkbox"/> | 121WL10 | Omit | NC |

Extension Handle (Only avail w/ High Mount Wheel Locks)

| | | | |
|--------------------------|--------|---------|------|
| <input type="checkbox"/> | 121EH1 | 6" Pair | \$73 |
|--------------------------|--------|---------|------|

Swing Away 8"-12" Frame Depth

Wheel Lock go/no-go matrix

| Frame Depth | C.G. | 20" wheel | 22" wheel | 24" wheel | Frame Depth | C.G. | 20" wheel | 22" wheel | 24" wheel |
|-------------|-------|-----------|-----------|-----------|-------------|-------------------------------|-----------|-----------|-----------|
| 8" | -0.5" | | | | 11 | -0.5" | | | |
| | 0" | | | | | 0" | | | |
| | 0.5" | | | | | 0.5" | | | |
| | 1" | | | | | 1" | | | |
| 9" | -0.5" | | | | 12 | 1.5" | | | |
| | 0" | | | | | -0.5" | | | |
| | 0.5" | | | | | 0" | | | |
| | 1" | | | | | 0.5" | | | |
| 10" | -0.5" | | | | | 1" | | | |
| | 0" | | | | | 1.5" | | | |
| | 0.5" | | | | | No-go. User must use hublock. | | | |
| | 1" | | | | | Go | | | |
| | 1.5" | | | | | | | | |

6. HANGER AND FOOTREST

Z-finity Hangers/Footrest System

N/A on Fixed Frame.

Some interference may exist in certain configurations.

| | | | |
|--------------------------|---------|---|-------|
| <input type="checkbox"/> | 121H162 | Z-finity Hanger | STD |
| <input type="checkbox"/> | 121H160 | Elevating Hanger - Ext Mount <i>Includes ELR pad. N/A w/ Z-finity. Some limitations may exist.</i> | \$147 |
| <input type="checkbox"/> | 121H144 | Omit Hanger | NC |

Z-finity Footrest Mount (Must Choose see Illustration)

| | | | |
|--------------------------|---------|-------------------------------------|----|
| <input type="checkbox"/> | 121DI22 | Direct Hanger Mt/Fixed Frame - High | NC |
| <input type="checkbox"/> | 121DI23 | Direct Hanger Mt/Fixed Frame - Low | NC |
| <input type="checkbox"/> | 121DI24 | Extension Tube Mount | NC |

Z-finity Footrest Pivot Arm (Must Choose, use matrix for leg length)

| | | | |
|--------------------------|---------|-------------|----|
| <input type="checkbox"/> | 121DI20 | Short 2" | NC |
| <input type="checkbox"/> | 121DI21 | Long 3-1/4" | NC |

Footplate

Some interference may exist in certain configurations.

| | | | |
|--------------------------|---------|---|-------|
| <input type="checkbox"/> | 121F320 | Multi-Angle Adjustable - Pair <i>(N/A w/ 8"-9" widths)</i> | \$247 |
| <input type="checkbox"/> | 121F321 | Angle Adjustable Platform Footboard | \$231 |
| <input type="checkbox"/> | 121F6 | Omit Footplate | NC |

Extension Tube Length (When Choosing Ext. Tube Mount)

| | | | |
|--------------------------|--------|--------|----|
| <input type="checkbox"/> | 121E62 | Hemi | NC |
| <input type="checkbox"/> | 121E1 | Short | NC |
| <input type="checkbox"/> | 121E2 | Medium | NC |
| <input type="checkbox"/> | 121E5 | Long | NC |

Hanger Mounted Footrest (Footplates Mount to Hanger Tubing)

N/A with 8" and 9" width. Footplate mounted in lowest position may cause interference w/ caster.

| | | | |
|--------------------------|---------|------------------------------|----|
| <input type="checkbox"/> | 121H157 | 80° Swing In/Out Front Mount | NC |
| <input type="checkbox"/> | 121H158 | 90° Swing In/Out Front Mount | NC |

Extensions (may interfere with some seat heights & casters)

| | | | |
|--------------------------|--------|------------------|------|
| <input type="checkbox"/> | 121E23 | 2" (not mounted) | \$53 |
| <input type="checkbox"/> | 121E24 | 4" (not mounted) | \$53 |


Footplate

| | | | |
|--------------------------|-------|--|-------|
| <input type="checkbox"/> | 121F9 | Kids Angle Adjustable (only available w/ 10"-16" frame widths) | NC |
| <input type="checkbox"/> | 121F8 | Adult Angle Adj. (only available with 10"-16" widths) | \$176 |

Swing In/Out Hanger Leg Lengths

(Seat sling to footrest measured)


| 80° Hanger | Min | Max |
|--------------|-------|-------|
| No Ext. Tube | 5.7" | 9.7" |
| 2" Ext. Tube | 11.2" | 12.2" |
| 4" Ext. Tube | 11.2" | 14.1" |
| 90° Hanger | Min | Max |
| No Ext. Tube | 6.0" | 10.0" |
| 2" Ext. Tube | 11.5" | 12.5" |
| 4" Ext. Tube | 11.5" | 14.5" |



Z-finity Matrix

Direct Hanger Mount

Extension Tube Mount



Direct Hanger Mounted

| Pivot Arm | Hanger Mounting Position | Min | Max |
|-----------|--------------------------|-------|------|
| Short | High Mount | 0.5" | 4.5" |
| | Low Mount | 3.5" | 8.0" |
| Long | High Mount | -1.0" | 5.5" |
| | Low Mount | 2.5" | 9.0" |

Extension Tube Mounted

| Pivot Arm | Extension Tube Length | Min | Max |
|-----------|-----------------------|-------|-------|
| Short | Hemi Short | 5.8" | 14.1" |
| | Short | 6.5" | 15.1" |
| | Medium | 8.4" | 17.0" |
| | Long | 10.3" | 19.0" |
| Long | Hemi Short | 4.6" | 15.3" |
| | Short | 5.2" | 16.3" |
| | Medium | 7.2" | 18.3" |
| | Long | 9.2" | 20.2" |

ELR Ext Mt.(MED Ext Tube)

| Pivot Arm | Min | Max |
|-----------|------|-------|
| Short | 7.0" | 18.0" |
| Long | 6.0" | 19.0" |

Footrest Options

| | | | |
|--------------------------|---------|--|------|
| <input type="checkbox"/> | 121HL8 | Leg Strap | \$61 |
| <input type="checkbox"/> | 121HL11 | Heel Loops (N/A w/ Kids Angle Adj. Z-finity. N/A on widths below 14") | NC |
| <input type="checkbox"/> | 121IG2 | Neoprene Impact Guards (N/A with footrest in highest mounted position) | \$61 |

7. BACKRESTS

BACKREST

| | | | | | |
|----------------------------------|--|----|-----------------------------------|--------------------------------------|---------|
| <input type="checkbox"/> 121BT10 | Standard Backrest (Angle adj./folding) | NC | <input type="checkbox"/> 121BT101 | Dynamic Backrest (N/A w/ Fixed Back) | \$1,558 |
|----------------------------------|--|----|-----------------------------------|--------------------------------------|---------|

BACK TYPE (Choose Height Adjustable or Fixed Height)

| | | | | | |
|-----------------------------------|--------------------------------------|----|----------------------------------|--------------|----|
| <input type="checkbox"/> 121BT120 | Height Adjustable (see matrix below) | NC | <input type="checkbox"/> 121BT53 | Fixed Height | NC |
|-----------------------------------|--------------------------------------|----|----------------------------------|--------------|----|

Height Adjustable (See Backrest Height Range matrix)

Must choose Lower Base Tube, Backrest Style, and Upper Back Tube.

Lower Base Tube (See Backrest Height Matrix Below)

| | | |
|----------------------------------|----------------------------------|----|
| <input type="checkbox"/> 121BP58 | Short | NC |
| <input type="checkbox"/> 121BP59 | Medium | NC |
| <input type="checkbox"/> 121BP60 | Tall (N/A with Dynamic Backrest) | NC |

Backrest Styles (See Images Below)

| | | |
|----------------------------------|---|-------|
| <input type="checkbox"/> 121PH16 | Straight Non Push Handle | NC |
| <input type="checkbox"/> 121PH6 | Integral Push handle | NC |
| <input type="checkbox"/> 121PH17 | Folding Push Handle with Ht. Adj. Back (8"-10" widths, stroller handle will face outward) | \$310 |
| <input type="checkbox"/> 121PH15 | Omit Back Tube (See omit matrix) | NC |

Upper Back Tube (See Backrest Height Matrix Below)

| | | |
|----------------------------------|--------|----|
| <input type="checkbox"/> 121B245 | Short | NC |
| <input type="checkbox"/> 121B5 | Medium | NC |

Height Adjustable Backrest


Preset Back height (Choose for Set Up)
No Upholstery w/ 19"-24"

| | | |
|------------------------------|------------------------------|------------------------------|
| <input type="checkbox"/> 9" | <input type="checkbox"/> 15" | <input type="checkbox"/> 20" |
| <input type="checkbox"/> 11" | <input type="checkbox"/> 16" | <input type="checkbox"/> 21" |
| <input type="checkbox"/> 12" | <input type="checkbox"/> 17" | <input type="checkbox"/> 22" |
| <input type="checkbox"/> 13" | <input type="checkbox"/> 18" | <input type="checkbox"/> 23" |
| <input type="checkbox"/> 14" | <input type="checkbox"/> 19" | <input type="checkbox"/> 24" |

Fixed Height

Fixed Height back N/A with upholstery


| | | |
|----------------------------------|-------------------------------------|-------|
| <input type="checkbox"/> 121BP9 | Standard (16" Height) | NC |
| <input type="checkbox"/> 121BP10 | Medium (20" Height) | NC |
| <input type="checkbox"/> 121BP11 | Tall (23" Height) | NC |
| <input type="checkbox"/> 121BP61 | Standard w/ Folding PH (16" Height) | \$310 |
| <input type="checkbox"/> 121BP62 | Medium w/ Folding PH (20" Height) | \$310 |
| <input type="checkbox"/> 121BP63 | Tall w/ Folding PH (23" Height) | \$310 |



Backrest Height Ranges

Upper Back Tube (Choose One)

| | | | Straight non-push handle Short | Straight non-push handle Medium | Integral push handle Short | Integral push handle Medium | Folding Push Handle Short | Folding Push Handle Medium | Omit |
|------------------------------|----------|--------|--------------------------------|---------------------------------|----------------------------|-----------------------------|---------------------------|----------------------------|------|
| Lower Base Tube (Choose one) | Standard | Short | 11"-15" | 15"-19" | 11"-15" | 15"-19" | 12"-16" | 16"-20" | 9" |
| | | Medium | 15"-19" | 18"-23" | 15"-19" | 18"-23" | 16"-20" | 19"-24" | 13" |
| | | Tall | 18"-22" | 18"-24" | 18"-22" | 18"-24" | 19"-23" | 19"-24" | 16" |
| | Dynamic | Short | 14"-18" | 18"-22" | 14"-18" | 18"-22" | 15"-19" | 19"-23" | 12" |
| | | Medium | 18"-22" | 21"-24" | 18"-22" | 21"-24" | 19"-23" | 22"-24" | 16" |



7. BACKRESTS Cont.

Back Upholstery

Uph N/A with Fixed Height back. Std Uph N/A w/Flip-Up. Adj. Locking Flip-Up and Length Adj Locking Flip-up . 3DX N/A on 8"-11" widths. Ballistic N/A on 8"-9" widths.

| | | | |
|--------------------------|---------|----------------------------------|-------|
| <input type="checkbox"/> | 121BU1 | Standard Upholstery | STD |
| <input type="checkbox"/> | 121BU71 | Backrest Tension Adj. 3DX Vented | \$471 |
| <input type="checkbox"/> | 121BU70 | Backrest Tension Adj. Ballistic | \$190 |
| <input type="checkbox"/> | 121BU66 | Omit Back Upholstery | NC |

JAY Back (Optional)

To order any other JAY back on-chair, please See the JAY Price Lst. Requires Part #

| | | | |
|--------------------------|---------|-------------------------|---------------|
| <input type="checkbox"/> | 121JB02 | JAY Bundled Backrest | Price List |
| | | Backrest Name or Part # | |
| | | Width (in.): | Height (in.): |

Stabilizer Bar

N/A w/ Dynamic back and Fixed Back. Only available w/ Integral push handle.

| | | | |
|--------------------------|--------|----------------|-------|
| <input type="checkbox"/> | 121SZ1 | Stabilizer Bar | \$168 |
|--------------------------|--------|----------------|-------|

8. ARMRESTS/SIDEGUARDS

Armrest and Armpad - Height Adjustable

Flip-Up. Adj. Locking Flip-Up and Length Adj. Locking Flip-up N/A with Std Back Upholstery. Low RSH in combination with large wheels may result in arm interference. SPHA arms N/A with One Arm Drive or Xtender.

| | | | |
|--------------------------|---------|---|-------|
| <input type="checkbox"/> | 121AA58 | T-Post Xtra Tall(with sideguard) | \$278 |
| <input type="checkbox"/> | 121AA59 | T-Post 2X Tall(with sideguard) | \$278 |
| <input type="checkbox"/> | 121AA3 | Flip-Up | NC |
| <input type="checkbox"/> | 121AA1 | Adjustable Locking Flip-Up | \$339 |
| <input type="checkbox"/> | 121AA7 | Length Adjust. Locking Flip-Up | \$339 |
| <input type="checkbox"/> | 121AA11 | Single-Post Ht Adj. Std. (N/A w/ 3" camber or 8"-9" frame depths) | \$339 |
| <input type="checkbox"/> | 121AA29 | Single-Post Height Adj. Low (N/A w/ 3" camber or 8"-9" frame depths) | \$339 |
| <input type="checkbox"/> | 121AA14 | Omit Armrest | NC |

Side Guard

N/A with height adjustable armrests and T-Posts

| | | | |
|--------------------------|---------|-------------|-------|
| <input type="checkbox"/> | 121SG33 | Plastic Std | \$188 |
| <input type="checkbox"/> | 121SG29 | Plastic Low | \$188 |


Armpad

N/A with T-Post and Adj. Locking Flip-Up and Foam only available with Length Adj. Locking Flip-Up.


| | | | |
|--------------------------|---------|-------------------------------------|----|
| <input type="checkbox"/> | 121AM23 | Desk Classic | NC |
| <input type="checkbox"/> | 121AM24 | Full Classic (N/A w/ 10"-13" Depth) | NC |
| <input type="checkbox"/> | 121AM9 | Foam Short (N/A w/ T Post) | NC |
| <input type="checkbox"/> | 121AM10 | Foam Long (N/A w/ T Post) | NC |

Single Post Height Adjustable Armrest

| Pad | Armrest Size | Transfer Tube Length (D) | Range (H) | Pad Length (L) | |
|---------|--------------|--------------------------|-------------------|----------------|------|
| | | | | Desk | Full |
| Classic | STD | 8" | 9-1/2" to 14-1/2" | 10-1/4" | 14" |
| | LOW | 5-1/4" | 8" to 12" | | |



| Armrest Size | T-POST | |
|--------------|------------------|----------------|
| | Range (H) | Pad Length (L) |
| XTALL | 6" to 7" | 7" |
| 2X TALL | 7-1/2" to 8-1/2" | |



| 9. ACCESSORIES | | | |
|---|----------|--|-------|
| Back Pack <i>(one Backpack Black/Charcoal already included)</i> | | Seat Pouch | |
| <input type="checkbox"/> | 121A1 | Black Utility Backpack | \$112 |
| <input type="checkbox"/> | 121A2 | Backpack Black/Charcoal Kids | \$105 |
| Miscellaneous | | | |
| <input type="checkbox"/> | 1218A03 | Caddy <i>(Only avail w/ Ext. Mt. hangers)</i> | \$89 |
| <input type="checkbox"/> | 1213A01 | Tool Kit | \$43 |
| <input type="checkbox"/> | 121A2100 | Crutch Holder | \$115 |
| Positioning Belt | | | |
| <input type="checkbox"/> | 121A823 | Belt 1" Dual Adj. Sqz Blk Pad | \$105 |
| <input type="checkbox"/> | 121A824 | Belt 2" Aircraft Buckle | \$139 |
| <input type="checkbox"/> | 121A825 | Belt 2" Aircraft Buckle Padded | \$204 |
| <input type="checkbox"/> | 121A826 | Belt 1.5" Auto Buckle | \$89 |
| <input type="checkbox"/> | 121A814 | 1" Comfort Fit APB5116 | \$146 |
| Upper Extremity Support Tray | | | |
| <i>N/A with 8"-9" widths</i> | | | |
| <input type="checkbox"/> | 121A29 | X-Small <i>(Suggested 10"-12" chair width)</i> | \$381 |
| <input type="checkbox"/> | 121A30 | Small <i>(Suggested 13"-14" chair width)</i> | \$381 |
| <input type="checkbox"/> | 121A31 | Medium <i>(Suggested 15"-17" chair width)</i> | \$381 |
| <input type="checkbox"/> | 121A32 | Large <i>(Suggested 18"-20" chair width)</i> | \$381 |
| Canopy (SunShade) | | | |
| <i>Need 2" back post clearance to mount canopy hardware</i> | | | |
| <input type="checkbox"/> | 1218A06 | Canopy (Sunshade) 10"-14" Width | \$428 |
| <input type="checkbox"/> | 1218A07 | Canopy (Sunshade) 15"-20" Width | \$428 |

| 10. COLOR <i>(Standard with frame and cross brace color choices)</i> | | | |
|--|---------|---|-------|
| Frame Color (Must choose) | | Cross Brace Color | |
| Standard Colors | | Standard Colors | |
| <input type="checkbox"/> | 121C211 | Celestial Sapphire | NC |
| <input type="checkbox"/> | 121C212 | Desert Sage | NC |
| <input type="checkbox"/> | 121C213 | Desert Tan | NC |
| <input type="checkbox"/> | 121C214 | Gloss Black | NC |
| <input type="checkbox"/> | 121C215 | Matte Black | NC |
| <input type="checkbox"/> | 121C216 | Metallic Shadow | NC |
| <input type="checkbox"/> | 121C217 | Midnight Blue | NC |
| <input type="checkbox"/> | 121C218 | Mystic Purple | NC |
| <input type="checkbox"/> | 121C219 | Petal Pink | NC |
| <input type="checkbox"/> | 121C220 | Pink Power | NC |
| <input type="checkbox"/> | 121C221 | Red Pepper | NC |
| <input type="checkbox"/> | 121C222 | Royal Blue | NC |
| <input type="checkbox"/> | 121C223 | Slaute Green | NC |
| <input type="checkbox"/> | 121C224 | Scarlet Noir | NC |
| <input type="checkbox"/> | 121C225 | Silver Ore | NC |
| <input type="checkbox"/> | 121C226 | Solar Orange | NC |
| <input type="checkbox"/> | 121C227 | Starlight White | NC |
| <input type="checkbox"/> | 121C228 | Storm Gray | NC |
| <input type="checkbox"/> | 121C229 | Turbo Green | NC |
| <input type="checkbox"/> | 121C230 | Velocity Yellow | NC |
| Personalized Colors (Must choose cross brace std. color) | | | |
| <input type="checkbox"/> | 121C43 | Mossy Oak Camouflage | \$368 |
| <input type="checkbox"/> | 121C42 | Stars and Stripes | \$368 |
| <input type="checkbox"/> | 121C108 | Zebra <i>Choose base colour. N/A w/ Gloss Black or Matte Black.</i> | \$368 |
| <input type="checkbox"/> | 121C627 | Celestial Sapphire - Cross Brace | NC |
| <input type="checkbox"/> | 121C628 | Desert Sage - Cross Brace | NC |
| <input type="checkbox"/> | 121C629 | Desert Tan - Cross Brace | NC |
| <input type="checkbox"/> | 121C630 | Gloss Black - Cross Brace | NC |
| <input type="checkbox"/> | 121C631 | Matte Black - Cross Brace | NC |
| <input type="checkbox"/> | 121C632 | Metallic Shadow - Cross Brace | NC |
| <input type="checkbox"/> | 121C633 | Midnight Blue - Cross Brace | NC |
| <input type="checkbox"/> | 121C634 | Mystic Purple - Cross Brace | NC |
| <input type="checkbox"/> | 121C635 | Petal Pink - Cross Brace | NC |
| <input type="checkbox"/> | 121C636 | Pink Power - Cross Brace | NC |
| <input type="checkbox"/> | 121C637 | Red Pepper - Cross Brace | NC |
| <input type="checkbox"/> | 121C638 | Royal Blue - Cross Brace | NC |
| <input type="checkbox"/> | 121C639 | Salute Green - Cross Brace | NC |
| <input type="checkbox"/> | 121C640 | Scarlet Noir - Cross Brace | NC |
| <input type="checkbox"/> | 121C641 | Silver Ore - Cross Brace | NC |
| <input type="checkbox"/> | 121C642 | Solar Orange - Cross Brace | NC |
| <input type="checkbox"/> | 121C643 | Starlight White - Cross Brace | NC |
| <input type="checkbox"/> | 121C644 | Storm Gray - Cross Brace | NC |
| <input type="checkbox"/> | 121C645 | Turbo Green - Cross Brace | NC |
| <input type="checkbox"/> | 121C246 | Velocity Yellow - Cross Brace | NC |

Additional Options Through Built-4-Me

To inquire about any unique personalization for Zippie products, contact: Phone: (800) 263-3390, Email: cscanada@sunmed.com

Commonly Requested Options Include:

- ✓ Solid Seat
- ✓ Grade Aid
- ✓ Drum Brake



Sunrise Medical Inc. Canada
1000 Creditstone Road, Unit 2 · Concord, ON L4K 4P8 · Canada