



ORDER FORM

Squiggles Stander

EILK02

The prone, upright, or adjustable 3-in-1 stander designed for growing children.

Customer			
Account			
Account #		Date	
☐ Quote	☐ Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

Squiggles Stander

Weight Capacity: 48.4 lbs

✓ OPTION	DESCRIPTION	RETAIL
Multi-Position Packages		\$
Package includes Squiggles Standard Model- Multi-Position, Pivot Base, Positioning Pads, Flat Headrest Hardware, Cushion and Laterals, Tray, Small Sandals with Straps. Substitutions are not available for Packages.		
☐ 141QE19	Squiggles Stander Package - Blue	\$6,061
☐ 141QE20	Squiggles Stander Package - Blue Vinyl	\$6,061
☐ 141QE21	Squiggles Stander Package - Green	\$6,061
☐ 141QE22	Squiggles Stander Package - Pink	\$6,061
Standard Model		\$

Selection Required.

Includes hip and chest support with sternum prompt; chest and hip lateral supports; size 2 knee cups, footplate and tray attachment.

☐ 141BS34	Standard Support Shell - Multi-Position (Prone, Supine, and Upright)	\$2,156
☐ 141BS23	Squiggles Pivot Base	\$1,742
Positioning Pads		\$


Selection Required. Includes chest and hip pad covers and knee support covers.

☐ 141ZI15	Squiggles Stander Positioning Pads - Blue	\$559
☐ 141ZI16	Squiggles Stander Positioning Pads - Blue Vinyl	\$559
☐ 141ZI17	Squiggles Stander Positioning Pads - Green	\$559
☐ 141ZI18	Squiggles Stander Positioning Pads - Pink	\$559

Accessories		\$
☐ 141XX31	Size 1 Knee Cup	\$88
☐ 141XX32	Size 3 Knee Cup	\$88
☐ 141A208	Tray - Gray	\$402
☐ 141A216	Padded Tray Insert - Black	\$331
☐ 1418A106	Grab Rail / Toy Bar for Tray	\$326
☐ 1418A107	Small Sandals w/ Straps (Pair)	\$273
☐ 1418A108	Medium Sandals w/ Straps (Pair)	\$291
☐ 1419A24	Footplate Raiser 1" (Each)	\$97



Sandals



Grab Rail

Squiggles Stander

✓ OPTION	DESCRIPTION	RETAIL
Headrest		\$

Cushion selection required with hardware.

Lateral Bracket Assembly Requires Flat Headrest Support Hardware and Cushion.

Headrest selection NOT available with Prone Configuration.

☐	141AW24	Contoured Headrest Support Hardware	\$396
☐	141Y040	Contoured Headrest Cushion - Blue	\$194
☐	141Y041	Contoured Headrest Cushion - Blue Vinyl	\$194
☐	141Y042	Contoured Headrest Cushion - Green	\$194
☐	141Y043	Contoured Headrest Cushion - Pink	\$194
☐	141Y027	Contoured Headrest Cushion - Black	\$188
☐	141AW22	Flat Headrest Support Hardware	\$247
☐	141Y044	Flat Headrest Cushion - Blue	\$124
☐	141Y045	Flat Headrest Cushion - Blue Vinyl	\$124
☐	141Y046	Flat Headrest Cushion - Green	\$124
☐	141Y047	Flat Headrest Cushion - Pink	\$124
☐	141Y028	Flat Headrest Cushion - Black	\$124
☐	141AW23	Lateral Bracket Assembly Hardware (Pair) Requires Flat Headrest Support Hardware and Cushion	\$194
☐	141Y048	Flat Headrest Lateral Cushions - Blue (Pair)	\$194
☐	141Y049	Flat Headrest Lateral Cushions - Blue Vinyl (Pair)	\$194
☐	141Y050	Flat Headrest Lateral Cushions - Green (Pair)	\$194
☐	141Y051	Flat Headrest Lateral Cushions - Pink (Pair)	\$194
☐	141Y052	Flat Headrest Lateral Cushions - Black (Pair)	\$302



Headrest Support



Flat Headrest



Lateral Cushions



Contoured Headrest

Technical Specifications

Frame Sizing Specification	Imperial	Metric
Age (approximate)	1 - 5 years	1 - 5 years
Maximum User Weight	48.4 lbs	22 kgs
User Standing Height	28.5" - 43.3"	730 mm - 1100 mm
Chest Width Between Laterals	6.3" - 9"	160 mm - 230 mm
Chest Support Height	21.7" - 32"	550 mm - 820 mm
Chest Support Angle Adjustment	+/- 30°	+/- 30°
Chest Support Depth Adjustment	2"	50 mm
Hip Width Between Hip Guides	6.3" - 9"	160 mm - 230 mm
Hip Support Height from Footplate to Mid Buttocks	11.8" - 20.8"	300 mm - 530 mm
Footplate Angle	10° Plantarflexion / 10° Dorsiflexion	10° Plantarflexion / 10° Dorsiflexion
Small Sandal Size	6.3" L x 3" W	175 mm L x 75 mm W
Medium Sandal Size	8.3" L x 3.2" W	210 mm L x 80 mm W
Tray Height: Prone	21.9" - 32.5"	555 mm - 825 mm
Tray Height: Supine	18.9" - 32.5"	480 mm - 825 mm
Tray Angle Adjustment	Prone 30° / Supine 40°	Prone 30° / Supine 40°
Stander Support Frame Weight	15.4 lbs	7 kgs
Stander Support Frame Dimensions	26.8" L x 15.4" W x 15.4" H	680 x 390 x 390 mm
Pivot Base Weight	12.1 lbs	5.5 kgs
Pivot Base Dimensions (Unfolded)	32.3" L x 21.7" W x 20.3" H	820 x 550 x 515 mm
Pivot Base Dimensions (Folded)	27.6" L x 21.7" W x 11.8" H	700 x 550 x 300 mm
Knee Cup Opening Width	Size 1: 2.5" Size 3: 3.5"	Size 1: 63 mm Size 3: 89 mm

Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2
Concord, ON L4K 4P8



MK-129857-CAN Rev. 11 © 04/2025

