
IN MINI

Q300//Mini

Imagine the Possibilities

The sky's the limit with the ZIPPIE Mini's ultra-compact base and small turning radius.

ZIPPIE



zippie.

Q30C/JMini

Imagine the places they will go



Imagine Independence

Your child can move shoulder to shoulder with family or friends while driving 3 MPH (walking speed). A smooth riding drive suspension is always active, no matter what seat elevation you drive from. Your child will be able to steal cookies from cabinets, write on the chalkboard, and give you smooches at your level.

Tight places, such as classrooms or busy living spaces and narrow doorways won't slow down the exploration and fun. The TRUE mid-wheel drive technology allows the ZIPPIE Mini to turn on the spot for intuitive and easy driving, providing the best indoor maneuverability.







Imagine your style

Personalize the ZIPPIE Mini with accent colors or choose from a shroud color for a bit more color.

Available shroud colors



Available accent colors





















ZIPPIE Q300 M Mini Specifications

HCPCS codes	Group 3 — Standard Duty No Power Seating: K0848 Single Power Seating: K0856
Max. user weight	165 lbs.
Seat width	10" to 16"
Seat depth	10" to 19"
Seat-to-floor height	16.5" to 19"
Back height	12" to 20"
Total width	20.5″
Total length (without footplates)	33″
Min. turning radius (without footplates)	17"
Max. speed	6 mph.
Battery size	Group 22 –50 Ah
Overall weight	No Power Seating: starting at 250 lbs. Single Power Seating: starting at 274 lbs.
Drive wheel position	Mid
Max. ground clearance	3″
Wheel options	12" or 14" drive wheels, 7" casters
Max. range	Up to 17 miles
Transit securement	ISO 7176-19 2008 Transit Approved

See order form as specifications may vary based upon configuration.



Imagine Freedom

The ZIPPIE Mini seating system supports your child so they are prepared for all the big and small adventures in everyday life. The chair can also adapt and grow based on your child's needs.



Order Form



Web Page





800.333.4000 www.SunriseMedical.com