

JAY[®]

X2™ CUSHION

- ≡ Lightweight, reduced profile skin protection for active users



JAY® X2™ CUSHION

≡ The JAY X2 is a lightweight, maintenance-free, reduced profile fluid cushion designed to provide stability and pressure relief to active clients. Using elements from the JAY Xtreme and Active cushions, the X2 combines JAY's time-tested designs with its latest technologies.

Lightweight, Stable Foam Base with Reduced Profile

- 2.25" cushion height provides support without sacrificing overall seat-to-floor height
- Flat base allows client to sit at the front of the cushion comfortably during daily tasks and activities
- Notch cut-out for additional protection over the coccyx bone
- Optional 1" lateral thigh supports provide additional support for the active client



Maintenance-Free JAY Flow™ Fluid Insert

- JAY Flow™ fluid insert conforms to each individual and ensures proper fluid placement to protect skin from breakdown
- Zero-maintenance design delivers superior skin protection throughout the life of the cushion



JAY's Innovative Dual-Cover Technologies

- Stretch outer cover with zipper enclosure for secure placement
- Water-resistant inner cover with AquaGuard® zipper and anti-wicking seam to protect foam base from moisture build-up



SPECIFICATIONS

JAY X2 Cushion	
HCPCS Code	E2622
Weight	3.8 lbs. (16 x 16)
Product Width	14" to 20"
Product Depth	14" to 20"
Product Height	2.25"
Weight Capacity	300 lbs.
Outer Cover	Stretch